



Simplified Edge 3 TB Flex Order Form

BASE MODEL

300 lb. weight capacity with standard 6mph programmable motors and 8 amp off board charger. Dimensions of base without legrests: 24"Wx35.5"L.

Quantum Edge 3.....\$8,450
Standard with Switch Operated Fender Lights

PATIENT WEIGHT

0 - 200 lbs

251 - 300 lbs

200 - 250 lbs

BASE COLOR

Red Corvette

Raspberry Beret

Lemon Crush

Back in Black

Evolution Orange

Green Machine

Night Train

Midnight Blue

Military Green

Grey Street

Electric Blue

Pretty in Pink

Deep Purple

Blue Jean

DRIVE MOTORS

6 mph Drive Motors Standard

Accu-Trac Motor Technology \$995

BASE MOUNTED FOOT PLATFORM

Base Mounted Foot Platform Standard

Omit base mounted foot platform No Charge

DRIVE WHEELS

3"x14" Flat Free Black Tires Standard

3"x14" Pneumatic Black Tires No Charge

Drive Wheel Accent Color

MatchingStandard

ELECTRONICS

QLNESStandard

Q-Logic 3 E\$620

Q-Logic 3\$900

Standard with Colour Match Shroud

JOYSTICK MOUNTING BRACKETS

Fixed, In-line, Bracket Standard

Adjustable Swing Away \$265

Offset Swing Away \$265

Basic Swing Away \$265

Left Side Mount

Right Side Mount

JOYSTICK HANDLES

Mushroom Handle (A) \$120

Soft Ball Handle (E) \$120

Chin Cup (B) \$90

Extended I-Handle (F) \$140

T-Bar (C) \$160

3" Goal Post Handle \$180

Medium Ball Handle (D) \$70

4" Goal Post Handle \$180

SEATING SYSTEM

TB Flex Seat \$650

With TB Flex 24"H Quick Angle Adjustable, Flip Forward back canes

Seat Dimension - Width

15" W (#1005) 17" W (#1007) 19" W (#1009) 21" W (#1011)

16" W (#1006) 18" W (#1008) 20" W (#1010) 22" W (#1012)

Seat Dimension - Depth

14" D 15" D 16" D 17" D 18" D 19" D 20" D

Seat to Floor

Lowest (16.5") Middle (17.5") Highest (18.25")

Account Number: _____ Date: _____ Marked for: _____

SEATING SYSTEM (CONTINUED)

1b. Single Post, Height Adjustable, Drop In Armrests

☐ L ☐ R Full Length, Drop In Armrests No Charge
Only available with Standard Armpad (no armpad selection required)

☐ L ☐ R Desk Length, Drop In Armrests No Charge
Only available with Standard Armpad (no armpad selection required)

1b. Two-Post, Flip Back, Height Adjustable, Removable Armrests

☐ L ☐ R Full Length, Two-Post Armrests.....\$145 ea.

☐ L ☐ R Desk Length, Two-Post Armrests\$145 ea.

Straight Armpads..... Standard

Waterfall Armpads..... No Charge

Standard Lap Belts

☐ 38" Lap Belt

☐ 50" Lap Belt

☐ 60" Lap Belt

☐ 70" Lap Belt

☐ 80" Lap Belt

☐ 90" Lap Belt

LEGRESTS

Center Mount Foot Platform

Foot Platform Assembly

Manually Adjustable Center Mount Foot Platform... \$265

Lower Extension and Footplates

10"-16" Extension, 5"x6" Footplates

10"-16" Extension, 5"x11" Footplates

13"-19" Extension, 5"x8" Footplates

10"-16" Extension, 5"x8" Footplates

13"-19" Extension, 5"x6" Footplates

13"-19" Extension, 5"x11" Footplates

Center Mount Calf Support

Small Calf Pads (pair) 6"Wx6"H \$210*

Large Calf Pads (pair) 8" taper to 6"Wx8"H \$210*

Modular Footplate Extension Components

Width Extension (pair) \$55

Available with 8"D and 11"D footplates

Outside Edge (pair) \$55

Available with 8"D and 11"D footplates

60°-90° Swing-Away Legrests

☐ P ☐ L ☐ R 60°-90° Swing-Away Legrest..... \$285

☐ P ☐ L ☐ R 60°-90° Swing-Away Legrest, Short \$285.

2. Lower Extension

☐ P ☐ L ☐ R 11"-13" Length

☐ P ☐ L ☐ R 13"-16.5" Length

☐ P ☐ L ☐ R 15.5"-19.5" Length

3b. Angle Adjustable Footplate(s)

¹The 6"x11" footplates may cause interference with seat widths less than 18". ²The 8"x11" footplates may cause interference with seat widths less than 20".

☐ P ☐ L ☐ R 5"Wx8"D \$87.50 ea.

☐ P ☐ L ☐ R 6"Wx11"D¹ \$87.50 ea.

☐ P ☐ L ☐ R 6"Wx8"D¹ \$87.50 ea.

☐ P ☐ L ☐ R 8"Wx11"D² \$87.50 ea.

BATTERIES

NF-22, 55 Ah Gel Cell, Interceptor Battery \$300ea.
Up to 3 extra miles additional range!

NF-22, 50 Ah Gel Cell, MK Battery \$400 ea.