



Simplified Edge 3 Static Seating Order Form

BASE MODEL

300 lb. weight capacity with standard 6mph programmable motors and 8 amp off board charger. Dimensions of base without legrests: 24"Wx35.5"L. XLR USB Charger is included with Edge 3.

Quantum Edge 3.....\$8,450
Standard with Switch Operated Fender Lights

PATIENT WEIGHT

0 - 200 lbs

251 - 300 lbs

200 - 250 lbs

BASE COLOR

Red Corvette

Raspberry Beret

Lemon Crush

Back in Black

Evolution Orange

Green Machine

Night Train

Midnight Blue

Military Green

Grey Street

Electric Blue

Pretty in Pink

Deep Purple

Blue Jean

DRIVE MOTORS

6 mph Drive Motors Standard

Accu-Trac Motor Technology \$995

BASE MOUNTED FOOT PLATFORM

Base Mounted Foot Platform Standard

Omit base mounted foot platform No Charge

DRIVE WHEELS

3"x14" Flat Free Black Tires Standard

3"x14" Pneumatic Black Tires No Charge

Drive Wheel Accent Color

Matching Standard

ELECTRONICS

QLNES Standard

Q-Logic 3 E \$620

Q-Logic 3 \$900

Standard with Colour Match Shroud

JOYSTICK MOUNTING BRACKETS

Fixed, In-line, Bracket Standard

Adjustable Swing Away \$265

Offset Swing Away \$265

Basic Swing Away \$265

Left Side Mount

Right Side Mount

JOYSTICK HANDLES

Mushroom Handle (A) \$120

Soft Ball Handle (E) \$120

Chin Cup (B) \$90

Extended I-Handle (F) \$140

T-Bar (C) \$160

3" Goal Post Handle \$180

Medium Ball Handle (D) \$70

4" Goal Post Handle \$180

SEATING SYSTEM

TB3 Static Seat \$1,045

Standard Towers

Tru Balance 3 12"
Power Adjustable
Seat Lift with
iLevel.....\$2,995

Seat Dimension - Width

14" W
(#1004)

16" W
(#1006)

18" W
(#1008)

20" W
(#1010)

22" W
(#1012)

15" W
(#1005)

17" W
(#1007)

19" W
(#1009)

21" W
(#1011)

Seat Dimension - Width Package

12"-16" Width Package

16"-20" Width Package

20"-24" Width Package

(Maximum width is 22")

Account Number: _____ Date: _____ Marked for: _____

SEATING SYSTEM (CONTINUED)

Seat Dimension - Depth

14" D 15" D 16" D 17" D 18" D 19" D 20" D

Seat Dimension - Depth Package

12"-16" Depth Package 16"-20" Depth Package 20"-24" Depth Package (Maximum depth is 22")

Seat to Floor

Lowest Middle Highest

Back Cane Height

16" High Canes 20" High Canes 24" High Canes
17" High Canes 21" High Canes 25" High Canes

1b. Two-Post, Flip Back, Height Adjustable, Removable Armrests

☐ L ☐ R Full Length, Two-Post Armrests.....\$145 ea. ☐ L ☐ R Desk Length, Two-Post Armrests\$145 ea.

1c. Single-Post, Heavy Duty, 8"-12" Height Adjustable, Removable Armrests

☐ L ☐ R Full Length, Heavy Duty Armrests, 8"-12" H\$130 ea. ☐ L ☐ R Desk Length, Heavy Duty Armrests, 8"-12"H.....\$130 ea.

1d. Single-Post, Heavy Duty, 10"-14" Height Adjustable, Removable Armrests

☐ L ☐ R Full Length, Heavy Duty Armrests, 10"-14"H.....\$130 ea. ☐ L ☐ R Desk Length, Heavy Duty Armrests, 10"-14"H\$130 ea.
Straight Armpads..... Standard Waterfall Armpads..... No Charge

Standard Lap Belts

☐ 38" Lap Belt ☐ 50" Lap Belt ☐ 60" Lap Belt ☐ 70" Lap Belt ☐ 80" Lap Belt ☐ 90" Lap Belt

LEGRESTS

Center Mount Foot Platform

Foot Platform Assembly

Manually Adjustable Center Mount Foot Platform... \$265

Lower Extension and Footplates

10"-16" Extension, 5"x6" Footplates 10"-16" Extension, 5"x11" Footplates 13"-19" Extension, 5"x8" Footplates
10"-16" Extension, 5"x8" Footplates 13"-19" Extension, 5"x6" Footplates 13"-19" Extension, 5"x11" Footplates

Center Mount Calf Support

Small Calf Pads (pair) 6"Wx6"H \$210* Large Calf Pads (pair) 8" taper to 6"Wx8"H \$210*

Modular Footplate Extension Components

Width Extension (pair) \$55
Available with 8"D and 11"D footplates
Outside Edge (pair) \$55
Available with 8"D and 11"D footplates

60°-90° Swing-Away Legrests

☐ P ☐ L ☐ R 60°-90° Swing-Away Legrest..... \$285 ☐ P ☐ L ☐ R 60°-90° Swing-Away Legrest, Short \$285.

2. Lower Extension

☐ P ☐ L ☐ R 11"-13" Length ☐ P ☐ L ☐ R 13"-16.5" Length ☐ P ☐ L ☐ R 15.5"-19.5" Length

3b. Angle Adjustable Footplate(s)

¹The 6"x11" footplates may cause interference with seat widths less than 18". ²The 8"x11" footplates may cause interference with seat widths less than 20".

☐ P ☐ L ☐ R 5"Wx8"D \$87.50 ea. ☐ P ☐ L ☐ R 6"Wx11"D¹ \$87.50 ea.
☐ P ☐ L ☐ R 6"Wx8"D¹ \$87.50 ea. ☐ P ☐ L ☐ R 8"Wx11"D² \$87.50 ea.

BATTERIES

NF-22, 55 Ah Gel Cell, Interceptor Battery \$300ea. NF-22, 50 Ah Gel Cell, MK Battery \$400 ea.
Up to 3 extra miles additional range!