

# QUANTUM® 4FRONT<sup>2</sup>

## 4FRONT<sup>2</sup>

shown with TRU-Balance<sup>®</sup> 3 seating system, Q-Logic 3 EX joystick on a TRU-Balance<sup>®</sup> 3 swing-away mount, TRU-Comfort 2 seat cushion, TRU-Comfort 4 back rest and Stealth Comfort Plus UniLink<sup>™</sup> headrest

  
**PROUDLY  
BUILT IN  
CANADA**  
WITH IMPORTED PARTS



AVAILABLE WITH 3.5 MPH  
10"  
SAFE SEAT ELEVATION



The 4Front<sup>®</sup> 2 is the next generation of the 4Front<sup>®</sup> Power Chair! New Smart Traction Control (STC), along with Smooth Ride Suspension (SRS), provide stability for smooth obstacle transitions and unmatched driving performance. Additional new features include 4-pole motors and redesigned caster arms for improved durability and aesthetics and reduced weight.

  
**QUANTUM<sup>®</sup>**  
#1 FOR REHAB POWER

**SMART TRACTION CONTROL**

 | [quantumrehab.ca](http://quantumrehab.ca)

## Specifications

<b>DRIVE WHEELS</b>	Standard: 14" flat free Optional: 14" pneumatic 14" pneumatic knobby 14" flat free knobby
<b>CASTER WHEELS</b>	Rear: 9" flat free
<b>MAXIMUM SPEED</b> <sup>1</sup>	Up to 6.25 mph (10 km/h) Up to 3.5 mph (5.6 km/h) elevated
<b>GROUND CLEARANCE</b>	3.15" at battery tray
<b>TURNING RADIUS</b> <sup>2</sup>	25.9"
<b>OVERALL LENGTH</b> <sup>2</sup>	39.9" without front riggings
<b>BASE WIDTH</b>	24"
<b>SEATING SIZES</b>	
<b>TRU-Balance® 3 Static Seat</b>	<b>Width:</b> 12-22" <b>Depth:</b> 12-22"
<b>TRU-Balance® 3 Power Tilt</b>	<b>Width:</b> 12-22" <b>Depth:</b> 12-22"
<b>TRU-Balance® 3 Power Lift &amp; Tilt</b>	<b>Width:</b> 12-22" <b>Depth:</b> 12-22"
<b>SEAT-TO-FLOOR HEIGHTS</b> <sup>9</sup>	
<b>TRU-Balance® 3 Static Seat</b> <sup>5</sup>	15.2-17.25"
<b>TRU-Balance® 3 Power Tilt</b> <sup>5</sup>	17.5-19.5"
<b>TRU-Balance® 3 Tilt &amp; iLevel</b> <sup>5</sup>	17.75-19.75"
<b>BATTERY SIZE</b> <sup>3</sup>	Group 24 or Group 34
<b>BATTERY WEIGHT</b> <sup>6</sup>	Group 24: 51.5 lb Group 34: 44.8 lb
<b>AVAILABLE ELECTRONICS</b> <sup>7</sup>	120A Q-Logic 3e (standard) 120A Q-Logic 3 EX (optional)
<b>BATTERY CHARGER</b>	8A Off-board
<b>MOTOR PACKAGES</b>	4-pole, front-wheel drive
<b>WEIGHT CAPACITY</b> <sup>4</sup>	300 lb
<b>BASE WEIGHT</b>	130 lb

1) Speed varies with user weight, terrain type, battery charge, battery condition and tire pressure.

2) Without front riggings

3) Battery capacity may change when adding power positioning systems.

4) Weight capacity is dictated by the selection of seating.

5) All measurements +/- 1/8" with 18" x 18" seat, fixed towers, 300 lbs. weight capacity and solid tires (add 3/8" with pneumatic tires). Widths may increase seat-to-floor height. Contact Customer Service for details. Listed measurements are in the lowest position; does not include elevating system unless specified.

6) Battery weight may vary +/- 3 lbs. based on manufacturer.

7) All Quantum Rehab® units can be ordered with specialty controls.

8) Call for availability

9) Seat-to-floor heights are measured with the unit unloaded. Due to the 4Front 2's suspension system, seat-to-floor heights are approximately 0.25" lower when the unit is loaded.

## Standard Features

- 4-pole motors
- Smart Traction Control (STC)
- Superior, automotive-grade SRS (Smooth Ride Suspension)
- Unmatched driving performance and ride quality
- Advanced stability for smooth obstacle transitions
- Redesigned caster arms for improved durability and aesthetics and reduced weight
- Q-Logic 3e controller
- Speeds up to 6.25 mph (10 km/h)
- LED front and rear fender lights
- Quiet ride technology

## Options\*


\* See order form for complete list of options.

- TB4 Seating System with Anterior Tilt up to 30 degrees
- 10" of power adjustable safe seat elevation at speeds of up to 3.5 mph (5.6 km/h)
- Quantum backup camera
- Trilogy vent tray <sup>8</sup>
- Compact vent tray <sup>8</sup>
- TRU-Balance® 3 static and power positioning systems
- TB-Flex seating system
- Q-Captain's seat
- Q-Logic 3 EX electronics
- Specialty drive controls
- Stealth Products® iDrive®
- Stealth Products® Seating and Positioning
- Q-Back 2.0

## Colors



## Hydrographics



**QUANTUM**®

#1 FOR REHAB POWER

quantumrehab.ca

