



**Pride Mobility Products Co.**

5096 South Service Rd, Beamsville, ON, L0R 1B3

Phone: 888-570-1113 | Fax: 866-514-1303 | Email: [productorders@pridemobility.com](mailto:productorders@pridemobility.com)

# Quantum 4Front<sup>®</sup> 2 Order Form with TRU-Balance 3 Positioning

Account Number: \_\_\_\_\_ Date: \_\_\_\_\_

Provider Name: \_\_\_\_\_

Contact: \_\_\_\_\_

Phone: \_\_\_\_\_ Fax: \_\_\_\_\_

Email: \_\_\_\_\_

PO Number: \_\_\_\_\_

Marked for: \_\_\_\_\_

Ship to Address: \_\_\_\_\_

City: \_\_\_\_\_ Prov: \_\_\_\_\_ PC: \_\_\_\_\_

## QUANTUM 4FRONT<sup>2</sup>



### INTRODUCTION

This form is interactive when viewed with [Adobe Acrobat Reader](#). Complete the form by placing checks in the desired boxes and provide information in the interactive fields. Buttons shown at the bottom of the form may be utilized to print or submit the order form through a desktop email application. To email via a web-based application, please 'Save As' and attach the PDF to your email.

This order form contains a large variety of options to fulfill various patient needs. Descriptions and section notes, such as optional or required, are included to help you complete your order. Other part numbers can be found on the quote acknowledgment after this order form is submitted. Please contact Quantum<sup>®</sup> Sales at 888-570-1113 if assistance is needed. Send the completed order form by fax (866-514-1303) or email ([productorders@pridemobility.com](mailto:productorders@pridemobility.com)). Incomplete forms may delay the quote or order. Customer service will contact you if the order is incomplete or if there are compatibility issues. If special order requests are needed, be sure that the Patient Information section is completed

All prices are MSRP. Prices, specifications, part numbers, and availability are subject to change without notice. Prices and part numbers as shown are only available when configured as a complete power chair. Please contact technical service for accurate parts ordering.

**Thank you for choosing Quantum!**


### PATIENT INFORMATION

Weight: \_\_\_\_\_ Height: \_\_\_\_\_

A. Shoulder Width: \_\_\_\_\_ D. Max Sitting Height: \_\_\_\_\_ G. Thigh Depth: \_\_\_\_\_

B. Chest Width: \_\_\_\_\_ E. Shoulder Height: \_\_\_\_\_ H. Lower Leg Length: \_\_\_\_\_

C. Hip Width: \_\_\_\_\_ F. Axilla Height: \_\_\_\_\_ I. Elbow Height: \_\_\_\_\_



### 1. BASE MODEL

**REQUIRED.** Select a group 3 power base. 300 lb. weight capacity with standard 6 mph programmable motors and 8 amp off board charger. Dimensions of base without legrests: 24.2"Wx39.9"L

Quantum 4Front 2 .....\$8,450  
Part: 4FRONT 2

### 1A. AUTO CORRECTION SYSTEM


















**REQUIRED.** Smart Traction Control is auto-selected with every 4Front 2.

Smart Traction Control ..... \$400  
Part: Various

**2. BASE COLOR**

**REQUIRED.** Select one shroud color. Colors shown here may differ from actual product depending on monitor/printer calibration used.

**2.1a. Standard Shrouds**

- |   |  |  |  |
|---|--|--|--|
|  <b>Red Corvette</b><br>Part: DGN2003028-PRD1  |  <b>White Satin</b><br>Part: DGN2003028-PWT2      |  <b>Midnight Blue</b><br>Part: DGN2003028-PBU10 |  <b>Military Green</b><br>Part: DGN2003028-PGN8 |
|  <b>Cherry Bomb</b><br>Part: DGN2003028-PRD4   |  <b>Deep Purple</b><br>Part: DGN2003028-PPR5      |  <b>Electric Blue</b><br>Part: DGN2003028-PBU9  |  <b>Sugar Plum</b><br>Part: DGN2003028-PPR7     |
|  <b>Back in Black</b><br>Part: DGN2003028-PBK1 |  <b>Raspberry Beret</b><br>Part: DGN2003028-PPK4  |  <b>Blue Jean</b><br>Part: DGN2003028-PBU11     |  |
|  <b>Night Train</b><br>Part: DGN2003028-PBK3   |  <b>Pretty in Pink</b><br>Part: DGN2003028-PPK3   |  <b>Lemon Crush</b><br>Part: DGN2003028-PYL4    |  |
|  <b>Grey Street</b><br>Part: DGN2003028-PSL3   |  <b>Evolution Orange</b><br>Part: DGN2003028-POR7 |  <b>Green Machine</b><br>Part: DGN2003028-PGN11 |  |

**2.1b. Hydrographic Shrouds**

-  **Military Digital Camo** ..... \$325  
Part: DGN2003028-WCM7

**3. BASE MOUNTED LIGHTS**

**Optional.** Select one option. Seat Mounted, LED, Full Light Package option can be found later on this form.

- Switch Operated Fender Lights**..... No Charge  
Part: DGN2004686  
Fender lights operate by pressing included single switch button. Fender lighting is not available with Single Actuator Control through Toggle Box.
- Base Mounted, LED, Full Light Package** ..... \$770  
Part: ELE1700035  
Requires Q-Logic 3 EX Joystick for Full Lights or Q-Logic 3e Joystick.
- Omit Fender Lights** ..... No Charge  
Part: DGN1700032

**4. QUANTUM BACKUP CAMERA**

**Optional.** The backup camera may interfere with armrest mounted accessories such as media device holder.

- Quantum Backup Camera** ..... \$475  
Part: MEC1910365

**1. Optional Mounting Bracket**

- Backup Camera Headrest Bracket**..... \$30  
Part: MEC1910373  
Works with .5" x .5" vertical supports only.



**5. DRIVE WHEELS**

**REQUIRED.** Select a type of drive wheels and accent color.

**5.1. Drive Wheels**

- |  |  |
|--|--|
| <b>3"x14" Flat Free Black Tires with Accent</b> ..... Standard<br>Part: PTOINDV3525    | <b>4"x14" Flat Free Knobby Black Tires</b> ..... \$60 ea.<br>Part: WHL2006647 x2 |
| <b>3"x14" Pneumatic Black Tires with Accent</b> ..... No Charge<br>Part: WHL1704059 x2 | <b>4"x14" Pneumatic Knobby Black Tires</b> ..... \$60 ea.<br>Part: WHL2006822 x2 |

**5.2. Drive Wheel Accent Color**

Not available with knobby tires.

- |  |  |   |
|--|--|---|
| <b>Back in Black</b> ..... Standard<br>Part: DGN164245-PBK1 x2 | <b>Raspberry Beret</b> ..... No Charge<br>Part: DGN164245-PPK4 x2  | <b>Green Machine</b> ..... No Charge<br>Part: DGN164245-PGN11 x2          |
| <b>Red Corvette</b> ..... No Charge<br>Part: DGN164245-PRD1 x2 | <b>Pretty in Pink</b> ..... No Charge<br>Part: DGN164245-PPK3 x2   | <b>Military Green</b> ..... No Charge<br>Part: DGN164245-PGN8 x2          |
| <b>Cherry Bomb</b> ..... No Charge<br>Part: DGN164245-PRD4 x2  | <b>Evolution Orange</b> ..... No Charge<br>Part: DGN164245-POR7 x2 | <b>Sugar Plum</b> ..... No Charge<br>Part: DGN164245-PPR7 x2              |
| <b>Night Train</b> ..... No Charge<br>Part: DGN164245-PBK3 x2  | <b>Midnight Blue</b> ..... No Charge<br>Part: DGN164245-PBU10 x2   | <b>Military Digital Camo</b> ..... \$42.50 ea.<br>Part: DGN164245-WCM7 x2 |
| <b>Grey Street</b> ..... No Charge<br>Part: DGN164245-PSL3 x2  | <b>Electric Blue</b> ..... No Charge<br>Part: DGN164245-PBU9 x2    |   |
| <b>White Satin</b> ..... No Charge<br>Part: DGN164245-PWH2 x2  | <b>Blue Jean</b> ..... No Charge<br>Part: DGN164245-PBU11 x2       |   |
| <b>Deep Purple</b> ..... No Charge<br>Part: DGN164245-PPR5 x2  | <b>Lemon Crush</b> ..... No Charge<br>Part: DGN164245-PYL4 x2      |   |

**6. BATTERIES**

**Optional.** The powerchair requires 2 flat top batteries and will be installed if ordered.

- Group 24, AGM Cell Battery (MK)**.....\$510 ea.  
Part: BAT1712207.
- Group 34, Gel Cell Battery (MK)** ..... \$420 ea.  
Part: BATGEL1006.  
Not available with TRU-Balance 3 Power Lift.

**7. PATIENT'S WEIGHT**

**REQUIRED.** The patient's weight is required to quote and order the power chair.

**Up to 200 lb User Weight**  
Part: PTOINDV3837

**201 lb-300 lb User Weight**  
Part: PTOINDV3838

**8. ELECTRONICS**

**REQUIRED.** Select a type of drive control from the following sections. The No Joystick option is used to omit the joystick and utilize another control, such as a specialty control at the end of this form or found on the [Alternative Controls and Electronics Order Form](#).

**8.1a. Omit Joystick**

**No Joystick, Setup for Q-Logic Electronics**  
Part: XRF000010

**8.1c. Q-Logic 3 Expandable Capacity Controls**

Q-Logic 3 expandable capacity controls allow use of specialty controls and multiple actuators.

**1. Q-Logic 3 Expandable Capacity Joysticks**

Q-Logic 3 Expandable Capacity Joystick .....\$900  
Part: CTL167724

Q-Logic 3 Expandable Capacity Joystick .....\$900  
for Full Light Package  
Part: CTL167725  
Required with full light packages.

Q-Logic 3e Joystick .....\$620  
Part: CTL1910600  
Can operate full light package. Not available with alternative or specialty controls.



**8.2. Q-Logic 3 Joystick Armor**

**REQUIRED with Q-Logic 3 Joystick selections above. Not available with Q-Logic 3e Joystick.** Reference colors shown on page 2.

**Black, Soft Touch** ..... No Charge  
Part: DGN124307

**Match Joystick Armor to Base Color** ..... Standard  
Part: PTOINDV3421  
Part below will be selected with system configurator.

**Non-matching Joystick Armor** ..... No Charge  
Part: PTOINDV3422  
Select color below

**Red Corvette**  
Part: DGN124336-PRD1

**Grey Street**  
Part: DGN124336-PSL3

**Pretty in Pink**  
Part: DGN124336-PPK3

**Blue Jean**  
Part: DGN124336-PBU11

**Sugar Plum**  
Part: DGN124336-PPR7

**Cherry Bomb**  
Part: DGN124336-PRD4

**White Satin**  
Part: DGN124336-PWH2

**Evolution Orange**  
Part: DGN124336-POR7

**Lemon Crush**  
Part: DGN124336-PYL4

**Military Digital Camo**  
Part: DGN124336-WCM7

**Back in Black**  
Part: DGN124336-PBK1

**Deep Purple**  
Part: DGN124336-PPR5

**Midnight Blue**  
Part: DGN124336-PBU10

**Green Machine**  
Part: DGN124336-PGN11

**Night Train\***  
\*See Black, Soft Touch

**Raspberry Beret**  
Part: DGN124336-PPK4

**Electric Blue**  
Part: DGN124336-PBU9

**Military Green**  
Part: DGN124336-PGN8

**9. JOYSTICK HANDLES**

**Optional.** The standard Q-Logic 3 joystick knob can be seen in an image in the previous section. An optional handle selected below will be mounted to the joystick and the standard handle will be shipped with the unit.

**Mushroom Handle (B)** ..... \$120  
Part: ACCASMB1860.

**Extended I-Handle (G)** ..... \$140  
Part: ACCASMB1857.

**Chin Cup (C)** ..... \$85  
Part: ACCASMB186190

**Bodypoint 3" Goal Post Handle** ..... \$180  
Part: PLSKNCP1057.

**T-Bar (D)** ..... \$160  
Part: ACCASMB1858.

**Bodypoint 4" Goal Post Handle** ..... \$180  
Part: PLSKNCP1058.

**Medium Ball Handle (E)** ..... \$70  
Part: ACCASMB1859.

**Soft Ball Handle (F)** ..... \$120  
Part: ACCASMB1856.



**10. JOYSTICK MOUNTING BRACKETS**

**Optional.** Select a mount to be used with a joystick selected in previous section.

**10.1a. Fixed, In-line Joystick Bracket**


Fixed, In-line, Bracket ..... Standard  
 Part: FRMASMB10617



**10.1b. Swing-Away (Retractable) Joystick Bracket**

Allows the joystick to retract to the side of the armrest while maintaining the same joystick orientation.

Swing-away, Q-Logic 3 ..... \$265  
 Part: ACC1908434.



**10.1c. Specialty Joystick Mounts**

Select a bracket and side to mount it.

Retract4 Mount ..... \$450  
 Part: XRF000003.

Full 180° of movement that translates into multi-position use, straight, right offset, left offset, or fold away. Joystick orientation varies as bracket is adjusted.

Multi-Axis Mount ..... \$100  
 Part: XRF000004  
 Allows height adjustment of joystick.

Stealth Swing-Away/Flip Down Bracket ..... \$300  
 Part: XRF000005.

Allows the joystick to pivot to side of armrest, or flip-down if bracket is adjusted.



**10.2. Joystick Mount Location**

Left Side Mount  
 Part: PTOINDV1225

Right Side Mount  
 Part: PTOINDV1226

**11. TRU-BALANCE 3® SEATING - QUANTUM INTELLIGENT SYSTEMS™**

**11.1. TRU-Balance 3 Configuration**

**REQUIRED.** Select one seating system. Be sure to complete the seat dimensions, seat to floor, and back selections.

TRU-Balance 3 Static Seat ..... \$1,045  
 Part: PTOINDV3442  
 20° of limited recline adjustment (-5/ +15) via gas strut within back cover.

TRU-Balance 3 Power Tilt Only ..... \$6,245  
 Part: PTOINDV3443.  
 50° Center of Gravity (CG) power tilt seating system.

TRU-Balance 3 Power Recline Only ..... \$7,875  
 Part: PTOINDV3444.  
 175° of power recline.

TRU-Balance 3 Power Tilt & Recline ..... \$13,125  
 Part: PTOINDV3446.  
 50° CG power tilt with 175° of power recline. Full drive lockout at 20°.

**11.2. 10" Power Lift**

**Optional.** Select one option below.

TRU-Balance 3 Power 10" Lift ..... \$2,495  
 Part: PTOINDV3775.  
 10" of scissor lift. 22"W x 22"D maximum seat size.

**11.3. Occupied Transit**

**Optional.** Available in seat dimensions of 16" to 20" wide and 16"-20" deep. Requires TRU-Balance 3 armrests and 60" Sure-lok lap belt. Only available with swing-away legrests and center mount foot platforms.

Occupied Transit ..... \$400  
 Part: KIT152047

To determine if your chair is WC-19 compliant (occupied), please refer to this website: <http://wc-transportation-safety.umtri.umich.edu/crash-tested-product-lists/wheelchairs>

**11.4. Seat Width**

**REQUIRED with seat selection.** 15" width is only available with non-matching back width. Minimum seat width is 12". Maximum seat width is 22".

12" W (#1002)	13" W (#1003)	14" W (#1004)	15" W (#1005)	16" W (#1006)	17" W (#1007)	18" W (#1008)	19" W (#1009)	20" W (#1010)	21" W (#1011)	22" W (#1012)
------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------

**11.5. Seat Width Package**

**REQUIRED with seat selection.**

12"-16" Width Package  
 Part: XRF000006  
 Not available with Power Lift.

14"-16" Width Package  
 Part: XRF000007  
 Only available with Power Lift.

16"-20" Width Package  
 Part: XRF000008

20"-24" Width Package  
 Part: XRF000158  
 Maximum width is 22".

**11.6. Seat Depth**

**REQUIRED with seat selection.** Minimum seat depth is 12". Maximum seat depth is 22".

12" D (#1002)	13" D (#1003)	14" D (#1004)	15" D (#1005)	16" D (#1006)	17" D (#1007)	18" D (#1008)	19" D (#1009)	20" D (#1010)	21" D (#1011)	22" D (#1012)
------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------

**11.7. Seat Depth Package**

**REQUIRED with seat selection.**

12"-16" Depth Package  
 Part: XRF000011  
 Not available with Power Lift.

14"-16" Depth Package  
 Part: XRF000012  
 Only available with Power Lift.

16"-20" Depth Package  
 Part: XRF000013

20"-24" Depth Package  
 Part: XRF000159  
 Maximum depth is 22".

TRU-Balance 3 seating section continued on following page.

**11. TRU-BALANCE 3 SEATING CONTINUED**

**11.8. Seat to Floor Height**

**REQUIRED.**

Lowest Seat to Floor  
Part: PTOINDV3558

Middle Seat to Floor  
Part: PTOINDV3559

Highest Seat to Floor  
Part: PTOINDV3560

Low	Med	High
Static Seat/ Recline Only		
15"	16"	17"
Tilt/Tilt & Recline/Lift		
17.75"	18.75"	19.75"
TRU-Balance 3 Seat to Floor Heights		


**11.9. Matching Seat Cushion**

**Optional. Specific cushion/component sizes available, part numbers, weight capacities, MSRP, and special order cushions can be found on the respective complete Stealth Cushion Order Forms.** This section allows a quick selection of a single cushion that matches the seat pan width and depth. See specific order form for non-matching sizes, alternate covers and liners, or additional covers and liners.

**Stealth Simplicity G Cushion**

General use cushion with Stealth Tek 2, 2-way stretch, breathable cover.  
Part series: ST-SIMGXXYY


Simplicity G, Matching Seat Size  
Part: XRF000111



**Stealth TRU-Comfort 2 SPP Cushion**

Skin protection and positioning with Coolcore, 4-way stretch, breathable cover.  
Part series: ST-TC2SPPXXYY


TRU-Comfort 2 SPP, Matching Seat Size  
Part: XRF000116



**Stealth Spectrum Foam SPP Cushion**

Skin protection and positioning with Coolcore, 4-way stretch, breathable cover. Includes urethane liner.  
Part series: ST-SPCFSPXXYY


Spectrum Foam SPP, Matching Seat Size  
Part: XRF000112



**Stealth Spectrum Gel SPP Cushion**

Skin protection and positioning with Coolcore, 4-way stretch, breathable cover. Includes urethane liner.  
Part series: ST-SPCGSPPXXYY


Spectrum Gel SPP, Matching Seat Size  
Part: XRF000117



**Stealth Solution SPP Cushion**

Skin protection and positioning with Coolcore, 4-way stretch, breathable cover. Includes urethane liner.  
Part series: ST-SOLSPPXXYY

Solution SPP, Matching Seat Size  
Part: XRF000113



**Stealth Essence SPP Cushion**

Skin protection and positioning with Coolcore, 4-way stretch, breathable cover.  
Part series: ST-ESSSPPXXYY


Essence SPP, No Wedge, Matching Seat Size  
Part: XRF000118

Essence SPP, Small Wedge, Matching Seat Size  
Part: XRF000119

Essence SPP, Medium Wedge, Matching Seat Size  
Part: XRF000120

Essence SPP, Large Wedge, Matching Seat Size  
Part: XRF000121


Essence Plus SPP, Large Wedge, Matching Seat Size  
Part: XRF000121



**Stealth Glacial SP Cushion**

Skin protection with Coolcore, 4-way stretch, breathable cover.  
Part series: ST-GLSPXXYY


Glacial SP, Matching Seat Size  
Part: XRF000114



**Stealth Zen SP Cushion**

Skin protection with Coolcore, 4-way stretch, breathable cover.  
Part series: ST-ZENSPXXYY

Zen SP, Matching Seat Size  
Part: XRF000115



**Stealth Premiere P Cushion**

Positioning with Coolcore, 4-way stretch, breathable cover.  
Part series: ST-PREMPXXYY

Premiere P, Matching Seat Size  
Part: XRF000122



**11.10. Back Type Selection**

**REQUIRED.** Match a choice in this section with the options in the following back cushion configuration section.

Redesigned, Static Type Back..... \$735  
Part: XRF000014

12" and 24" back width will use the TRU-Balance 3 Back assembly. Select a back cushion configuration in the subsequent section.

Back Cane Option..... No Charge  
Part: XRF000015

Available only with TRU-Balance 3 Static Seat and Power Tilt (1 actuator) selections. Select cane type, back cane height, and back cushion configuration below.

Redesigned, Flip Forward Back..... \$735  
Part: XRF000016

12" and 24" back width will use the TRU-Balance 3 Back assembly. Not available with manual or power recline systems. Only available with 2-post, flip back or pediatric armrests. Not available with 12"-14" width or with attendant manual recline. Select a back cushion type in the subsequent section.

**11.11. Back Cane Selection**

**REQUIRED when Back Cane Option is selected above.** Select type and height of back canes.

**1. Back Cane Type**

Angle Adjustable Back Cane..... No Charge  
Part: MEC140266

Flip Forward Back Cane ..... \$25  
Part: XRF000273

2" Endomorph (1" per side) Back Cane ..... \$70  
Part: XRF000269

4" Endomorph (2" per side) Back Cane..... \$70  
Part: XRF000271

**2. Back Cane Height**

16"H Cane  
Part: FRMASMB15834

17"H Cane  
Part: FRMASMB15835

20"H Cane  
Part: FRMASMB15838

21"H Cane  
Part: FRMASMB15839

24"H Cane  
Part: FRMASMB15842

25"H Cane  
Part: FRMASMB15843

TRU-Balance 3 seating section continued on following page.

**11. TRU-BALANCE 3 SEATING CONTINUED**

**11.12. Back Cushion Configuration**

**REQUIRED.** Sport Back and TRU-Comfort Back selections up to 22" back width will be the Sport 4 or TRU-Comfort 4 back cushions. Sport back and TRU-Comfort back selections of 24" back width will be the Sport or TRU-Comfort 2 back cushions. Select a configuration type and applicable dimensions. See the Sport, Q-Back, and TRU-Comfort back availability chart for stock order back size ranges. Other size backs may be available with special order and extended lead time.

**1a. Sport Back/Omit/After-Market Interface**

Sport Back ..... No Charge  
 Part: XRF000017  
 Complete following width package and size.

Omit Back Pan/Cushion ..... No Charge  
 Part: PTOINDV1287  
 Select a following width package.

After-market Back Interface ..... No Charge  
 Part: ST-SSAPT4-Q. Select a set width in following section, range of 12"-21". Overall seat depth may need to be adjusted to accommodate certain backs. May interfere with gel/waterfall armpads and arm through. Does not include brackets to mount after-market backs. A Special Order is required on back widths less than 16".

**1b. TRU-Comfort Back**

TRU-Comfort Back (Static/Flip Forward), Up to 15" Width..... \$835  
 Part: SET123622. HCPCS: E2620  
 Complete following width package and size.

TRU-Comfort Back (Static/Flip Forward), 16" to 20" Width... \$1,035  
 Part: TC4V10W21W. HCPCS: E2620  
 Complete following width package and size.

TRU-Comfort Back (Static/Flip Forward), 22" Width ..... \$1,235  
 Part: TC4V22W32W. HCPCS: E2621  
 Complete following width package and size.

TRU-Comfort Back, Back Cane Option ..... \$835  
 Part: TC4V22W32WBC. HCPCS: E2621  
 Complete following width package and size.

**1c. ADI Back, AL Series, 4-Point, Active Contour**

**Available with Back Cane Option only.** Select a back type and the specific overall back height in subsequent section. Previously selected back cane height must be equal to the greater of overall height range.

Active Contour, 14" Wide x 13" or 14" Overall Height ..... \$865  
 Part: ALBR14LA4.

Active Contour, 14" Wide x 16" or 17" Overall Height ..... \$885  
 Part: ALBR14MA4.

Active Contour, 16" Wide x 19" or 20" Overall Height ..... \$905  
 Part: ALBR16SA4.

Active Contour, 16" Wide x 21" or 22" Overall Height ..... \$925  
 Part: ALBR16TA4.

Active Contour, 16" Wide x 23" or 24" Overall Height ..... \$945  
 Part: ALBR16XA4.

Active Contour, 18" Wide x 19" or 20" Overall Height ..... \$905  
 Part: ALBR18SA4.

Active Contour, 18" Wide x 21" or 22" Overall Height ..... \$925  
 Part: ALBR18TA4.

Active Contour, 18" Wide x 23" or 24" Overall Height ..... \$945  
 Part: ALBR18XA4.

Active Contour, 20" Wide x 19" or 20" Overall Height ..... \$905  
 Part: ALBR20SA4.

Active Contour, 20" Wide x 21" or 22" Overall Height ..... \$925  
 Part: ALBR20TA4.

Active Contour, 20" Wide x 23" or 24" Overall Height ..... \$945  
 Part: ALBR20XA4.

**1d. ADI Back, AL Series, 4-Point, Deep Contour**

**Available with Back Cane Option only.** Select a back type and the specific overall back height in subsequent section. Previously selected back cane height must be equal to the greater of overall height range.

Deep Contour, 14" Wide x 13" or 14" Overall Height..... \$965  
 Part: ALBR14LD4.

Deep Contour, 14" Wide x 16" or 17" Overall Height..... \$985  
 Part: ALBR14MD4.

Deep Contour, 16" Wide x 19" or 20" Overall Height..... \$1,005  
 Part: ALBR16SD4.

Deep Contour, 16" Wide x 21" or 22" Overall Height..... \$1,025  
 Part: ALBR16TD4.

Deep Contour, 16" Wide x 23" or 24" Overall Height..... \$1,045  
 Part: ALBR16XD4.

Deep Contour, 18" Wide x 19" or 20" Overall Height..... \$1,005  
 Part: ALBR18SD4.

Deep Contour, 18" Wide x 21" or 22" Overall Height..... \$1,025  
 Part: ALBR18TD4.

Deep Contour, 20" Wide x 19" or 20" Overall Height ..... \$1,005  
 Part: ALBR20SD4.

Deep Contour, 20" Wide x 21" or 22" Overall Height..... \$1,025  
 Part: ALBR20TD4.

Deep Contour, 20" Wide x 23" or 24" Overall Height ..... \$1,045  
 Part: ALBR20XD4.

**1e. ADI Back, AL Series, 4-Point, Extra Deep Contour**

**Available with Back Cane Option only.** Select a back type and the specific overall back height in subsequent section. Previously selected back cane height must be equal to the greater of overall height range.

Extra Deep Contour, 16" Wide x 21" or 22" Overall Height .... \$1,125  
 Part: ALBR16T64.

Extra Deep Contour, 16" Wide x 23" or 24" Overall Height .... \$1,145  
 Part: ALBR16X64.

Extra Deep Contour, 18" Wide x 19" or 20" Overall Height .... \$1,105  
 Part: ALBR18S64.

Extra Deep Contour, 18" Wide x 21" or 22" Overall Height .... \$1,125  
 Part: ALBR18T64.

Extra Deep Contour, 18" Wide x 23" or 24" Overall Height.... \$1,145  
 Part: ALBR18X64.

Extra Deep Contour, 20" Wide x 19" or 20" Overall Height.... \$1,105  
 Part: ALBR20S64.

Extra Deep Contour, 20" Wide x 21" or 22" Overall Height.... \$1,125  
 Part: ALBR20T64.

Extra Deep Contour, 20" Wide x 23" or 24" Overall Height.... \$1,145  
 Part: ALBR20X64.

TRU-Balance 3 seating section continued on following page.

**11. TRU-BALANCE 3 SEATING CONTINUED**

**1f. Q-Back 2.0, AL Series, Active Contour**

**Not available with Back Cane Option.** Select a back type and the specific overall back height in subsequent section.

Active Contour, 14"Wx19-20"H Tilt/22"-23" H Recline..... \$985 Part: QB2AL14AS. E2613	Active Contour, 18"Wx19-20"H Tilt/22"-23" H Recline..... \$1,145 Part: QB2AL18AS. E2613
Active Contour, 14"Wx21-22"H Tilt/24"-25" H Recline ..... \$985 Part: QB2AL14AT. E2613	Active Contour, 18"Wx21-22"H Tilt/24"-25" H Recline ..... \$1,145 Part: QB2AL18AT. E2613
Active Contour, 16"Wx19-20"H Tilt/22"-23" H Recline ..... \$1,145 Part: QB2AL16AS. E2613	Active Contour, 20"Wx19-20"H Tilt/22"-23" H Recline..... \$1,145 Part: QB2AL20AS. E2613
Active Contour, 16"Wx21-22"H Tilt/24"-25" H Recline ..... \$1,145 Part: QB2AL16AT. E2613	Active Contour, 20"Wx21-22"H Tilt/24"-25" H Recline..... \$1,145 Part: QB2AL20AT. E2613

**1g. Q-Back 2.0, AL Series, Deep Contour**

**Not available with Back Cane Option.** Select a back type and the specific overall back height in subsequent section.

Deep Contour, 14"Wx19-20"H Tilt/22"-23" H Recline ..... \$1,145 Part: QB2AL14DS.	Deep Contour, 18"Wx19-20"H Tilt/22"-23" H Recline ..... \$1,325 Part: QB2AL18DS.
Deep Contour, 14"Wx21-22"H Tilt/24"-25" H Recline ..... \$1,145 Part: QB2AL14DT.	Deep Contour, 18"Wx21-22"H Tilt/24"-25" H Recline ..... \$1,325 Part: QB2AL18DT.
Deep Contour, 16"Wx19-20"H Tilt/22"-23" H Recline ..... \$1,325 Part: QB2AL16DS.	Deep Contour, 20"Wx19-20"H Tilt/22"-23" H Recline ..... \$1,325 Part: QB2AL20DS.
Deep Contour, 16"Wx21-22"H Tilt/24"-25" H Recline ..... \$1,325 Part: QB2AL16DT.	Deep Contour, 20"Wx21-22"H Tilt/24"-25" H Recline ..... \$1,325 Part: QB2AL20DT.

**11.14. TRU-Balance 3 Back Shroud**

**REQUIRED.** The back shroud is a feature of the TRU-Balance 3 seating system and is available in matching shroud colors. The back shroud is not available when Omit Back Pan/Cushion is selected. Reference colors shown on page 2.

<b>Red Corvette</b> Part: XRF000019	<b>Grey Street</b> Part: XRF000023	<b>Pretty in Pink</b> Part: XRF000027	<b>Blue Jean</b> Part: XRF000031	<b>Sugar Plum</b> Part: XRF000035
<b>Cherry Bomb</b> Part: XRF000020	<b>White Satin</b> Part: XRF000024	<b>Evolution Orange</b> Part: XRF000028	<b>Lemon Crush</b> Part: XRF000032	
<b>Back in Black</b> Part: XRF000021	<b>Deep Purple</b> Part: XRF000025	<b>Midnight Blue</b> Part: XRF000029	<b>Green Machine</b> Part: XRF000033	
<b>Night Train</b> Part: XRF000022	<b>Raspberry Beret</b> Part: XRF000026	<b>Electric Blue</b> Part: XRF000030	<b>Military Green</b> Part: XRF000034	

**11.15. TRU-Balance 3 Power Tilt Side Shrouds**

**REQUIRED with TRU-Balance 3 Power Tilt selections.** Reference colors shown on page 2.

<b>Red Corvette</b> Part: DGN1709093-PRD1	<b>Grey Street</b> Part: DGN1709093-PSL3	<b>Pretty in Pink</b> Part: DGN1709093-PPK3	<b>Blue Jean</b> Part: DGN1709093-PBU11	<b>Sugar Plum</b> Part: DGN1709093-PPR7
<b>Cherry Bomb</b> Part: DGN1709093-PRD4	<b>White Satin</b> Part: DGN1709093-PWT2	<b>Evolution Orange</b> Part: DGN1709093-POR7	<b>Lemon Crush</b> Part: DGN1709093-PYL4	
<b>Back in Black</b> Part: DGN1709093-PBK1	<b>Deep Purple</b> Part: DGN1709093-PPR5	<b>Midnight Blue</b> Part: DGN1709093-PBU10	<b>Green Machine</b> Part: DGN1709093-PGN11	
<b>Night Train</b> Part: DGN1709093-PBK3	<b>Raspberry Beret</b> Part: DGN1709093-PPK4	<b>Electric Blue</b> Part: DGN1709093-PBU9	<b>Military Green</b> Part: DGN1709093-PGN8	

**11.16. TRU-Balance 3 Power Lift Side Shrouds**

**REQUIRED with TRU-Balance 3, 10" Power Lift selection.** The power lift side shrouds are located at the base of the TRU-Balance 3 system and is available in matching shroud colors. Reference colors shown on page 2.

<b>Red Corvette</b> Part: DGN1800070	<b>Grey Street</b> Part: DGN1800071	<b>Pretty in Pink</b> Part: DGN1800068	<b>Blue Jean</b> Part: DGN1905060	<b>Sugar Plum</b> Part: DGN1905063
<b>Cherry Bomb</b> Part: DGN1905045	<b>White Satin</b> Part: DGN1800073	<b>Evolution Orange</b> Part: DGN1905051	<b>Lemon Crush</b> Part: DGN1905054	
<b>Back in Black</b> Part: DGN1800061	<b>Deep Purple</b> Part: DGN1800069	<b>Midnight Blue</b> Part: DGN1800063	<b>Green Machine</b> Part: DGN1905057	
<b>Night Train</b> Part: DGN1800062	<b>Raspberry Beret</b> Part: DGN1905048	<b>Electric Blue</b> Part: DGN1800064	<b>Military Green</b> Part: DGN1800066	

**11.17. Back Width Package**

**REQUIRED** except when Omit Back Pan/Cushion is selected. Available with +/-4" of seat width.

<b>Match Seat Width Selections</b> Part: PTOINDV3468	<b>14"W-16"W Package</b> Part: XRF000037	<b>18"W-22"W Package</b> Part: XRF000039
<b>10"W-14"W Package</b> Part: XRF000036	<b>16"W-20"W Package</b> Part: XRF000038	<b>20"W-24"W Package</b> Part: XRF000040

TRU-Balance 3 seating section continued on following page.

**11. TRU-BALANCE 3 SEATING CONTINUED**

**11.18. Back Width**

**REQUIRED** except when Omit Back Pan/Cushion is selected. Part number is PTOINDV plus number following # below. Available with +/-4" of seat width.

12" W (#3562)	14" W (#3563)	16" W (#3564)	17" W (#3565)	18" W (#3566)	19" W (#3567)	20" W (#3568)	21" W (#3569)	22" W (#3570)	24" W (#3571)
------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------

**11.19. Overall Back Height**

**Static Seat and Power Tilt available height range: 14"-25"H. Odd heights include 4" gap and even heights include 3" gap between seat pan and bottom of back cushion. Power Recline available height range: 19"-28"H. Odd heights include 6" gap and even heights include 7" gap. See the Sport and TRU-Comfort back availability chart for stock order back size ranges.**

<b>14" Height</b> Part: SETHGHT1027	<b>17" Height</b> Part: SETHGHT1030	<b>20" Height</b> Part: SETHGHT1033	<b>23" Height</b> Part: SETHGHT1036	<b>26" Height</b> Part: SETHGHT1039
<b>15" Height</b> Part: SETHGHT1028	<b>18" Height</b> Part: SETHGHT1031	<b>21" Height</b> Part: SETHGHT1034	<b>24" Height</b> Part: SETHGHT1037	<b>27" Height</b> Part: SETHGHT1040
<b>16" Height</b> Part: SETHGHT1029	<b>19" Height</b> Part: SETHGHT1032	<b>22" Height</b> Part: SETHGHT1035	<b>25" Height</b> Part: SETHGHT1038	<b>28" Height</b> Part: SETHGHT1041

Overall Height	Sport Back	TRU-Comfort Back	Sport or TRU-Comfort Back	Sport, Q-Back, or TRU-Comfort Back	Tilt or Static Width									Recline Width									
					12"	14"	15"	16"	17"	18"	19"	20"	21"	22"	24"	12"	14"	15"	16"	17"	18"	19"	20"
14"	T	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
15"	T	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
16"	S	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T
17"	S	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T
18"	-	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T
19"	-	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T
20"	-	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T
21"	-	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T
22"	-	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T
23"	-	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T
24"	-	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T
25"	-	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T
26"	-	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T
27"	-	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T
28"	-	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T
29"	-	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T
30"	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
31"	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
32"	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-

**Sport, Q-Back, and TRU-Comfort Back Availability Chart**

**11.20. Power Tilt Privacy Flap**

**Optional.** A privacy flap is standard with power recline systems. This selection is to add a flap for power tilt only systems.

Privacy Flap, 9.75" ..... \$35 Part: DWR1044E456B	Privacy Flap, 12.75" ..... \$35 Part: DWR1044E456A
--	---

**12. TRU-BALANCE 3 ARMRESTS**

**REQUIRED.** Select the desired type of armrest below. Please read notes to ensure seating compatibility. Be sure to also select an armrest size and armpad type.

**12.1a. Flip-Back, 8"-16" H Armrests**

Choose a left and right armrest and set height in the following section. Not available with waterfall armpads. Not available with Back Cane Option.

<b>L R</b> Flip-Back, 8"-16"H Armrests ..... \$145 ea. Left part: ARM1806662. Right part: ARM1806663
---

**Armrest Support**

**Optional.** Available up to 2" difference between seat and back width via endomorph or mesomorph configuration. Not available less than 16" width seat or back.

<b>L R</b> Flip-Back Armrest Downpost ..... \$115 ea. Left part: XRF000265. Right part: XRF000266
--

**12.1b. Two-Post, Height Adjustable, Removable Armrests**

Select the armrest and set the left and right heights in the following section. These armrests may interfere with offset backs/ back canes.

<b>L R</b> Full Length, Two-Post Armrest ..... \$145 ea. Left part: XRF000047. Right part: XRF000048	<b>L R</b> Desk Length, Two-Post Armrest ..... \$145 ea. Left part: XRF000049. Right part: XRF000050
---	---

*TRU-Balance 3 armrests section continued on following page.*



**12. TRU-BALANCE 3 ARMRESTS CONTINUED**

**12.1c. Heavy Duty, Removable Armrests**

Select a type of HD removable armrests and the set height in the following section. Only one height package is selectable as a normal configuration. Payer coding, coverage, and bundling guidelines may apply.

8"-12"H, Single Post, Height Adjustable, Removable, .....\$130 ea.  
Heavy Duty (Sold in pairs)  
Part: XRF000055.

10"-14"H, Single Post, Height Adjustable, Removable,.....\$130 ea.  
Heavy Duty (Sold in pairs)  
Part: XRF000056.

**12.2. Armrest Height**

**REQUIRED.** Set left and right armrest height. Note the available height range in the option description.

**1. Left Armrest Height**

8" H PTOINDV3535	10" H PTOINDV3537	12" H PTOINDV3539	14" H PTOINDV3541	16" H PTOINDV3543
9" H PTOINDV3536	11" H PTOINDV3538	13" H PTOINDV3540	15" H PTOINDV3542	

**2. Right Armrest Height**

8" H PTOINDV3450	10" H PTOINDV3452	12" H PTOINDV3454	14" H PTOINDV3456	16" H PTOINDV3458
9" H PTOINDV3451	11" H PTOINDV3453	13" H PTOINDV3455	15" H PTOINDV3457	

**12.3. Armrest Size**

**REQUIRED.** Select a left (L) and right (R) armrest size. This also affects armpad selection.

<input type="checkbox"/> L <input type="checkbox"/> R 10" Length Left part: XRF000057. Right part: XRF000123	<input type="checkbox"/> L <input type="checkbox"/> R 14" Length Left part: XRF000058. Right part: XRF000124	<input type="checkbox"/> L <input type="checkbox"/> R 17" Length Left part: XRF000059. Right part: XRF000125 Only available with Flip-Back, 8"-16" H Armrests and straight arm pads.
--	--	--

**13. TRU-BALANCE 3 ARMPADS**

**REQUIRED.** Select to omit an armpad selection or select a left and right armpad type. Please note that waterfall and gel armpads may cause interference with the seat back.

**13.1a. Omit Armpad**

Select one below if aftermarket armpads or handpads are required.

No Armpad Selection..... No Charge Part: XRF000060 Not available with Flip-back, 8"-16" H armrests.	Interface for Aftermarket Pad ..... No Charge Part: XRF000308 Only available with Flip-back, 8"-16" armrests.
---	---

**13.1b. Foam Armpads**

<input type="checkbox"/> L <input type="checkbox"/> R Straight, Armpad ..... Standard Left part: XRF000061. Right part: XRF000126	<input type="checkbox"/> L <input type="checkbox"/> R Waterfall, Armpad ..... No Charge Left part: XRF000062. Right part: XRF000127 Not available with Flip-back, 8"-16" armrests.
--	--

**13.1c. Gel Armpads**

Sold only in pairs. Armpads are 3" wide. Manufacturer part number for full length is 314G and desk length is 310G.

**a. Falcon Gel Armpads for Flip-Back Armrests**

Falcon Gel, Full Length Left and Right ..... \$300 Part: ARM1806807	Falcon Gel, Desk Length Left, Full Length Right..... \$300 Part: ARM1806810
Falcon Gel, Full Length Left, Desk Length Right..... \$250 Part: ARM1806813	Falcon Gel, Desk Length Left and Right ..... \$250 Part: ARM1806806

**b. Falcon Gel Armpads for Tube Style Armrests**

Falcon Gel, Full Length Left and Right ..... \$300 Part: ARM136438	Falcon Gel, Desk Length Left, Full Length Right..... \$300 Part: ARM136439
Falcon Gel, Full Length Left, Desk Length Right..... \$250 Part: ARM136442	Falcon Gel, Desk Length Left and Right ..... \$200 Part: ARM136437

**14. TRU-BALANCE 3 JOYSTICK BRACKETS**

**REQUIRED.**

Standard Fixed Receiver .....Standard Part: XRF000063	Stealth Height Adjustable Receiver ..... \$165 Part: XRF000064
--	---

**15. TRU-BALANCE 3 HEADRESTS**

**Optional.** Select a headrest pad and bracket.

**15.1a. Headrests, Comfort Plus with Cool Core**

**1. Headrest Pads**

Stealth 8" Comfort Plus, Cool Core (CP08CCFPM) ..... \$190  
 Part: ST-CP08CCFPM-Q. E0955

Stealth 10" Comfort Plus, Cool Core (CP10CCFPM)..... \$190  
 Part: ST-CP10CCFPM-Q. E0955

Stealth 14" Comfort Plus, Cool Core (CP14CCFPM)..... \$190  
 Part: ST-CP14CCFPM-Q. E0955

Stealth 6", Comfort Plus, Cool Core (CP06CCFPM)..... \$190  
 Part: ST-CP06CCFPM-Q. E0955

**2. Headrest Brackets**

Comfort Plus Headrest Mounting Bracket, UniLink.....\$255  
 Part: ST-UHB2D33200-Q. E1028

**15.1b. Headrests, Comfort Plus**

**1. Headrest Pads**

Stealth 8" Comfort Plus (CP450)..... \$190  
 Part: POS148484. HCPCS: E0955

Stealth 10" Comfort Plus (CP250)..... \$190  
 Part: SETHEAD1057. HCPCS: E0955

Stealth 14" Comfort Plus (CP150) ..... \$190  
 Part: INDPART2327. HCPCS: E0955

Pediatric Comfort Plus (6"Wx4.5"H) (CPN350)..... \$190  
 Part: SETHEAD1069. HCPCS: E0955

**2. Headrest Brackets**

Multi-Axis Removable Headrest Bracket (HMO475P) ..... \$200  
 Part: HDWASMB1084. HCPCS: E1028

**16. LAP BELTS**

**REQUIRED.** Select a lap belt below.

**16.1a. Standard Lap Belts**

38" Lap Belt  
 Part: ACCBELT1006

50" Lap Belt  
 Part: ACCBELT1000

60" Lap Belt  
 Part: ACCBELT1002

70" Lap Belt  
 Part: ACCBELT1003

80" Lap Belt  
 Part: ACCBELT1004

90" Lap Belt  
 Part: ACCBELT1001

**16.1b. Optional Lap Belts**

Stealth Retractable Lap Belt (IFRB5350).....\$175  
 Part: ST-IFRB5350

Stealth Positioning, 60" Single Pull, Push Button, Medium..... \$65  
 (7.5") Padded, 1.5" Belt (IFPB5213)  
 Part: ST-IFPB5213-Q

Stealth Positioning, 60" Single Pull, Push Button, Large ..... \$65  
 Padded (9.5"), 2" Belt (IFPB5314)  
 Part: ST-IFPB5314-Q

Stealth Positioning, 60" Single Pull, Push Button, Extra Large... \$65  
 Padded (11.5"), 2" Belt (IFPB5315)  
 Part: ST-IFPB5315-Q

60" Sure-Lok Lap Belt (Non-retractable)..... \$200  
 Required with ISO 7176-19 Compliant Transit Option.  
 Part: ACCBELT1016

**17. TRU-BALANCE 3 LEGREST CONTROLS**

**Optional.** Specific power legrest options are selected later in this form. Select the control type here, if power legrests are desired, to assist the selection of the power positioning electronics.

Independent Legrests (Requires 2 actuator controls)  
 Part: XRF000108

Articulating Foot Platform (Requires 1 Actuator Control)  
 Part: XRF000110

Combined Legrests (Requires 2 actuator control)  
 Part: XRF000109

**18. LEGRESTS**

**Optional.** All efforts will be made to accommodate legrest length selection, but minor changes may need to be made to meet shipping guidelines. Some adjustments may be needed upon receipt of the unit. Reference the section's compatibility matrix for legrest interference.

**18.1a. Center Mount Foot Platform**

If power articulating foot platform is selected, be sure power legrest control has been selected in earlier section. Select legrest length in following section.

**1. Foot Platform Assembly**

Power Articulating Foot Platform .....\$3,145  
 Part: PTOINDV3754. E1012+K0040

Manually Adjustable Center Mount Foot Platform... \$265  
 Part: XRF000155. K0040



**2. Lower Extension and Angle Adjustable Footplates**

8"-12" Extension, 5"x6" Footplates  
 Part: RIG152577  
 Not available with power articulating platform.

8"-12" Extension, 5"x8" Footplates  
 Part: RIG152578  
 Not available with power articulating platform.

10"-16" Extension, 5"x6" Footplates  
 Part: RIG150206

10"-16" Extension, 5"x8" Footplates with  
 Width Extension  
 Part: MEC1707469

10"-16" Extension, 5"x12" Footplates with  
 Width Extension  
 Part: MEC1707465

13"-19" Extension, 5"x6" Footplates  
 Part: RIG151505

13"-19" Extension, 5"x8" Footplates with  
 Width Extension  
 Part: MEC1707474

13"-19" Extension, 5"x12" Footplates with  
 Width Extension  
 Part: MEC1707468

Legrests section continued on following pages.

**18. LEGRESTS CONTINUED**

Center Mount Foot Platform Compatibility Matrix for TRU-Balance® 3	Power Recline or Static Seat			Power Tilt			Power Lift		
	Low	Med.	High	Low	Med.	High	Low	Med.	High
10"-16" Length Extension	!	✓	✓	✓	✓	✓	✓	✓	✓
13"-19" Length Extension	!	!	!	!	!	!	!	!	✓

- No interference issues at this seat-to-floor height.
- May have less than 1.75" clearance when fully extended and may interfere with the ground or reduce obstacle climbing ability.

**3a. Center Mount Calf Support**

\*Calf pads are standard with a power articulating foot platform. Select calf pads covers for additional charge. Calf pad and panel covers include 1/2" of soft visco foam.

Small Calf Pads, AFP Only, 6"Wx6"H (pair) ..... Standard  
Part: RIGASMB7120018

Large Calf Pads, AFP Only, 8" taper to 6"Wx8"H (pair) ... No Charge  
Part: RIGASMB7120019

Small Calf Pads, Manually Adjustable Center Mount Only, ..... \$210  
6"Wx6"H (pair)  
Part: XRF000156

Large Calf Pads, Manually Adjustable Center Mount Only, ..... \$210  
8" taper to 6"Wx8"H (pair)  
Part: XRF000157

**3b. Center Mount Calf Panel**

Calf pads are standard with a power articulating foot platform.

Small Calf Panel 18"Wx8"H..... \$210  
Part: FRMASMB7110012

Large Calf Panel 24"Wx8"H ..... \$265  
Part: RIG123824  
Only available on widths of 20" or greater.

**4. Modular Footplate Extension Components**

Optional. Assembly required. May require seating configuration adjustments to avoid front caster interference.

Width Extension for 8"D Footplates ..... \$55  
Part: RIG152517

4" Outside Edge (pair) ..... \$55  
Part: RIG152510

Width Extension for 11"D Footplates ..... \$55  
Part: RIG152518

Only available with 6"D footplates and 8"D footplate extension.

2" Bell-Shaped Length Extension ..... \$55  
Part: RIG152519

6" Outside Edge (pair) ..... \$55  
Part: RIG152516  
Available with 8"D and 11"D footplates and 11"D footplate extension.

Not compatible with width extensions. Can be used with 6" outside edge with 11" footplates.

**5. Additional Accessories**

Small Reversed Dartex Calf Pad Cover (pair) ..... \$150  
Part: RIG133717

Large Reversed Dartex Panel Cover (pair) ..... \$150  
Part: RIG134744

Large Reversed Dartex Calf Pad Cover (pair) ..... \$150  
Part: RIG133719

Footplate Connector ..... \$55  
Part: RIGASMB7120023

Small Reversed Dartex Panel Cover (pair)..... \$150  
Part: RIG134673

**6. Stealth Gel Footplate Covers**

5"x8" Tapered Gel Footplate Covers..... \$275  
Part: ST-FPG5008T-Q, K0108

6.5"x8" Tapered Gel Footplate Covers for Width Extension..... \$275  
Part: ST-FPG6508T-Q, K0108

5"x11" Tapered Gel Footplate Covers..... \$275  
Part: ST-FPG5011T-Q, K0108

7"x11" Tapered Gel Footplate Cover for Width Extension ..... \$275  
Part: ST-FPG7011T-Q, K0108

**18.1b. 70° Power Articulating, Swing Away Legrests**

4" of articulating range. Be sure power legrest control has been selected in previous section. Select an upper extension type and footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes. Select legrest length in following section.

**1. Articulating Legrest Upper Assembly**

70° Articulating, Swing-Away..... \$1,432.50 ea.  
Legrest with 14.5"-18" Extension  
Pair Part: FRMASMB9073. Left Part: FRMASMB9074.  
Right Part: FRMASMB9075. Pair E1010+K0040+E0951

70° Articulating, Swing-Away ..... \$1,432.50 ea.  
Legrest with 11"-13" Extension  
Pair Part: FRMASMB11444. Left Part: FRMASMB11366.  
Right Part: FRMASMB11367. Pair E1010+K0040+E0951

70° Power Articulating, Swing-Away Legrest Compatibility Matrix	Power Recline or Static Seat			Power Tilt			Power Lift		
	Low	Med.	High	Low	Med.	High	Low	Med.	High
Angle Adjustable, 11"-13"	✓	✓	✓	✓	✓	✓	✓	✓	✓
Angle Adjustable, 14.5"-18"	!	!	!	!	!	!	!	!	!

- This rigging configuration causes no interference issues at this seat-to-floor height.
- This rigging configuration may have less than 1.75" clearance when fully extended and may interfere with the ground or reduce obstacle climbing ability.

**2. Angle Adjustable Footplate with Heel Loop and Wheel Bumper**

\*The 6"x11" footplates may cause interference with seat widths less than 18". \*The 8"x11" footplates may cause interference with seat widths less than 20".

4"Wx6"D  
Pair part: FRMASMB7448. Left part: RIG133543. Right part: RIG133566.

5"Wx8"D  
Pair part: FRMASMB7450. Left part: RIG133545. Right part: RIG133568.

6"Wx11"D<sup>1</sup>  
Pair part: FRMASMB7452. Left part: RIG133548. Right part: RIG133570.

5"Wx6"D  
Pair part: FRMASMB7449. Left part: RIG133544. Right part: RIG133567.

6"Wx8"D<sup>1</sup>  
Pair part: FRMASMB7451. Left part: RIG133547. Right part: RIG133569.

8"Wx11"D<sup>2</sup>  
Pair part: FRMASMB7453. Left part: RIG133549. Right part: RIG133571.

**3. Power Articulating, Swing-Away Legrest Accessories**

Calf Pad (pair) 6"Wx6"H..... No Charge  
Part: MEC153735

Legrests section continued on following pages.

**18. LEGRESTS CONTINUED**

**18.1c. 60°-90° Swing-Away Legrests**

Adjustable in 5° increments. Select a legrest upper, a legrest extension, and type of footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes. Select legrest length in following section.

**1a. 60°-90° Swing-Away Legrest Upper**

**P L R** 60°-90° Swing-Away Legrest..... \$285  
 Pair part: RIG120062. Left: FRMASMB17558.  
 Right: FRMASMB17559  
 Not available with 11"-13" extensions.

**P L R** 60°-90° Swing-Away Legrest, Short..... \$285  
 Pair part: RIG136481. Left: RIG132607.  
 Right: RIG132609.  
 For use only with 11"-13" lower extensions.

**2. Swing-Away Legrests Lower Extension**

**P L R** 11"-13" Length  
 Pair part: XRF000088.  
 Left: XRF000089.  
 Right: XRF000090.  
 N/A with Multi-axis.

**P L R** 13"-16.5" Length  
 Pair part: XRF000094.  
 Left: XRF000095.  
 Right: XRF000096.  
 N/A with Multi-Axis.

**P L R** 11"-16.5" Length  
 Pair part: XRF000091.  
 Left: XRF000092.  
 Right: XRF000093.  
 N/A with Angle Adjustable.

**P L R** 15.5"-19.5" Length  
 Pair part: XRF000097.  
 Left: XRF000098.  
 Right: XRF000099.

60°-90° Swing-Away Legrest Compatibility Matrix	Power Recline or Static Seat			Power Tilt			Power Lift		
	Low	Med.	High	Low	Med.	High	Low	Med.	High
Extruded, 11"-13"	!	✓	✓	✓	✓	✓	✓	✓	✓
Extruded, 13"-16.5"	!	✓	✓	✓	✓	✓	✓	✓	✓
Extruded, 15.5"-19.5"	!	!	!	!	!	!	!	!	✓
Angle Adjustable, 13"-16.5"	!	✓	✓	✓	✓	✓	✓	✓	✓
Angle Adjustable, 15.5"-19.5"	!	!	!	!	!	!	!	!	✓
Multi-Axis, 11"-16.5"	!	✓	✓	✓	✓	✓	✓	✓	✓
Multi-Axis, 15.5"-19.5"	!	!	!	!	!	!	!	!	✓

11"-13" extension causes no interference at all seat-to-floor heights.

This rigging configuration causes no interference issues at this seat-to-floor height.

This rigging configuration may have less than 1.75" clearance when fully extended and may interfere with the ground or reduce obstacle climbing ability.

**3a. Extruded Footplate(s)**

**P L R** 7.25"Wx6"D Extruded Footplate..... Standard  
 Part pair: FRMASMB7847. Left: FRMASMB7861  
 Right: FRMASMB7862

**P L R** 8.5"Wx6"D Extruded Footplate ..... Standard  
 Part pair: FRMASMB7848. Left: FRMASMB7863  
 Right: FRMASMB7864  
 May cause interference on seat widths less than 18."

**3b. Angle Adjustable Footplate(s) and Heel Loops**

<sup>1</sup>The 6"x11" footplates may cause interference with seat widths less than 18". <sup>2</sup>The 8"x11" footplates may cause interference with seat widths less than 20".

**P L R** 4"Wx6"D ..... \$87.50 ea.  
 Pair part: FRMASMB10564. Left: RIG133537. Right: RIG133560

**P L R** 6"Wx8"D<sup>1</sup> ..... \$87.50 ea.  
 Pair part: FRMASMB10567. Left: RIG133540. Right: RIG133563

**P L R** 5"Wx6"D..... \$87.50 ea.  
 Pair part: FRMASMB10565. Left: RIG133538. Right: RIG133561

**P L R** 6"Wx11"D<sup>1</sup> ..... \$87.50 ea.  
 Pair part: FRMASMB10568. Left: RIG133541. Right: RIG133564

**P L R** 5"Wx8"D ..... \$87.50 ea.  
 Pair part: FRMASMB10566. Left: RIG133539. Right: RIG133562

**P L R** 8"Wx11"D<sup>2</sup> ..... \$87.50 ea.  
 Pair part: FRMASMB10569. Left: RIG133542. Right: RIG133565

**3c. Multi-Axis Footplate(s) and Heel Loops**

<sup>1</sup>The 6"x11" footplates may cause interference with seat widths less than 18". <sup>2</sup>The 8"x11" footplates may cause interference with seat widths less than 20".

**P L R** 4"Wx6"D ..... \$260 ea.  
 Pair part: FRMASMB9456. Left: RIG133531. Right: RIG133552

**P L R** 6"Wx8"D<sup>1</sup> ..... \$260 ea.  
 Pair part: FRMASMB9459. Left: RIG133534. Right: RIG133555

**P L R** 5"Wx6"D..... \$260 ea.  
 Pair part: FRMASMB9457. Left: RIG133532. Right: RIG133553

**P L R** 6"Wx11"D<sup>1</sup> ..... \$260 ea.  
 Pair part: FRMASMB9460. Left: RIG133535. Right: RIG133558

**P L R** 5"Wx8"D ..... \$260 ea.  
 Pair part: FRMASMB9458. Left: RIG133533. Right: RIG133554

**P L R** 8"Wx11"D<sup>2</sup> ..... \$260 ea.  
 Pair part: FRMASMB9461. Left: RIG133536. Right: RIG133559

**4. 60°-90° Swing-Away Legrest Accessories**

Gel Neoprene Footrest Sleeve (pair) ..... \$80  
 Part: INDPART2243. HCPCS: K0108

Legrests section continued on following pages.

**18. LEGRESTS CONTINUED**

**18.1d. Elevating Legrests**

Select a legrest upper, a legrest extension, and type of footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes. Select legrest length in following section.

**1. Elevating Legrests Upper**

P  L  R Elevating Legrest Assembly ..... \$275 ea.  
 Pair part: FRMASMB4755. Left: RIG138012  
 Right: RIG138013. E0990  
 Select lower extension and footplates

**2. Lower Extension**

P  L  R 11"-14" Length  
 Pair: FRMASMB4756. Left: RIG138026. Right: RIG138030  
 For reinforced only

P  L  R 8.5"-12" HD Drop-Ins/12.5"-16.5" ELR, Ext Length  
 Pair part: XRF000101. Left: XRF000102. Right: XRF000103  
 For angle adjustable only

P  L  R 14"-17" Length  
 Pair: FRMASMB4756. Left: RIG138033. Right: RIG138035  
 For reinforced only

P  L  R 13"-18" HD Drop-Ins/15.5"-19.5" ELR, Ext Length  
 Pair part: XRF000104. Left: XRF000105. Right: XRF000106  
 For angle adjustable only

**3a. Reinforced Cast Footplate(s)**

P  L  R 7.25"Wx6"D Reinforced Footplate..... Standard  
 Pair part: FRMASMB4748. Left: RIG138055. Right: RIG138056

P  L  R 8.5"Wx6"D Reinforced Footplate ..... Standard  
 Pair part: FRMASMB4747. Left: RIG138060. Right: RIG138062  
 May cause interference on seat widths less than 18."

**3b. Angle Adjustable Footplate(s)**

<sup>1</sup>The 6"x11" footplates may cause interference with seat widths less than 18". <sup>2</sup>The 8"x11" footplates may cause interference with seat widths less than 20". 4"x6" and 5"x6" footplates are only available in pairs.

P  L  R 4"Wx6"D ..... \$87.50 ea.  
 Pair parts: FRMASMB8812 x2.

P  L  R 6"Wx8"D<sup>1</sup> Footplates with Heel Loops ..... \$87.50 ea.  
 Pair part: FRMASMB4738. Left: RIG137793. Right: RIG137795

P  L  R 5"Wx6"D..... \$87.50 ea.  
 Pair parts: FRMASMB8811 x2.

P  L  R 6"Wx11"D<sup>1</sup> Footplates with Heel Loops ..... \$87.50 ea.  
 Pair part: FRMASMB4741. Left: RIG137796. Right: RIG137797

P  L  R 5"Wx8"D Footplates with Heel Loops ..... \$87.50 ea.  
 Pair part: FRMASMB4735. Left: RIG137791. Right: RIG137792

P  L  R 8"Wx11"D<sup>2</sup> Footplates with Heel Loops ..... \$87.50 ea.  
 Pair part: FRMASMB4744. Left: RIG137798. Right: RIG137799

Elevating Legrest Compatibility Matrix	Power Recline or Static Seat			Power Tilt			Power Lift		
	Low	Med.	High	Low	Med.	High	Low	Med.	High
Reinforced, 11"-14"	✓	✓	✓	✓	✓	✓	✓	✓	✓
Reinforced, 14"-17"	!	!	✓	!	✓	✓	✓	✓	✓
Angle Adjustable, 12"-16.5"	✓	✓	✓	✓	✓	✓	✓	✓	✓
Angle Adjustable, 15.5"-19.5"	!	!	!	!	!	✓	!	✓	✓

This rigging configuration causes no interference issues at this seat-to-floor height.  
 This rigging configuration may have less than 1.75" clearance when fully extended and may interfere with the ground or reduce obstacle climbing ability.

**18.2. Legrest Accessories**

Heel Loops (pair)..... \$35  
 Part: FRMASMB7873. HCPCS: E0951  
 Standard on angle adjustable and multi axis footplates when noted.

Wheel Bumper..... \$15 ea.  
 Part: XRF000100

**18.3. Set Legrest Extension Length**

**REQUIRED with matching left and right legrest configuration. Part number is PTOINDV plus number following # below. Non-matching or individual legrest configurations lengths are not set by the manufacturer.**

7" L (#3470)	8.5" L (#3473)	10" L (#3476)	11.5" L (#3479)	13" L (#3482)	14.5" L (#3485)	16" L (#3488)	17.5" L (#3491)	19" L (#3494)	20.5" L (#3497)	22" L (#3669)
7.5" L (#3471)	9" L (#3474)	10.5" L (#3477)	12" L (#3480)	13.5" L (#3483)	15" L (#3486)	16.5" L (#3489)	18" L (#3492)	19.5" L (#3495)	21" L (#3498)	22.5" L (#3670)
8" L (#3472)	9.5" L (#3475)	11" L (#3478)	12.5" L (#3481)	14" L (#3484)	15.5" L (#3487)	17" L (#3490)	18.5" L (#3493)	20" L (#3496)	21.5" L (#3600)	

**19. TRU-BALANCE 3 POWER POSITIONING ELECTRONICS**

**REQUIRED** with power seating or legrest options.

**19.1a. Q-Logic 3 Controls through Joystick**

Single Actuator Through Q-Logic 3 ..... No Charge  
 Controller and Interface (AM1)  
 Part: XRF000067.

Actuator Through Q-Logic 3 ..... \$1,675  
 Controller and Interface (AAM)  
 Part: XRF000068.



**1. Optional Switch Control**

Not available with AM1 electronics. 4-Way switches will be mounted to same side of joystick.

4-Function Push Button.....\$225  
 Part: CTL123218  
 See next section for switch mounting.

4-Way Toggle.....\$575  
 Part: ELEASMB5925

4-Way Button Toggle.....\$575  
 Part: ELEASMB5926

iAccess Flat Button.....\$575  
 Part: CTL144009  
 See next section for switch mounting.

iAccess Toggle Button.....\$575  
 Part: CTL144010  
 See next section for switch mounting.

**19.2. 4-Function Push Button/iAccess Mount**

**REQUIRED** with 4-Function Push Button or iAccess options.

Joystick Mount ..... No Charge  
 Part: MEC1907500

Attendant Mount ..... \$35  
 Part: MEC1906564

Under Arm Mount, Flip-Back Armrest..... No Charge  
 Part: MEC1906531

Under Arm Mount, Flip-Back Armrest, plus iAccess Guard ..... \$35  
 Part: MEC2001400

Under Arm Mount, Tube Style Armrest ..... No Charge  
 Part: MEC1906530

Under Arm Mount, Tube Style Armrest, plus iAccess Guard ..... \$35  
 Part: MEC2001401

**1. Under Arm Mount Side Selection**

Left Side Under Arm Mount  
 Part: PTOINDV3678

Right Side Under Arm Mount  
 Part: PTOINDV3679

**20. TRU-BALANCE 3 ATTENDANT CONTROL**

**Optional.** \*Attendant control will be at no charge when ordered with Tilt and Recline, alternative drive controls, and a Q-Logic Display or SCIM module. Mounting bracket is included.

Attendant Control ..... \$950  
 Part: XRF000078. E2331

**21. SEAT MOUNTED LED LIGHTS**

**Optional.** Requires Q-Logic 3 Joystick for Full Light Package or Q-Logic 3e Joystick.

Seat Mounted, LED, Full Light Package ..... \$770  
 Part: XRF000130  
 Not available with Base Mounted, LED, Full Light Package.




**22. PROVIDER INSTALLED ACCESSORIES**

**Optional.** Multiple accessories can be mounted along the various tracks but some interference may occur.

- Privacy Bag ..... \$80  
Part: ACC1907570
- TRU-Balance 3 Cane/Crutch Holder ..... \$100  
Part: ACC152515. HCPCS: E2207
- Quantum Backpack ..... \$200  
Part: ACC1905853
- Transfer Handle (D)..... \$305  
Part: MEC1809930. HCPCS: E0705  
Not available with Heavy Duty armrests.
- TRU-Balance 3 Cup Holder..... \$20  
Part: ACC125003  
For use with TRU-Balance 3 (track style) armrests.
- Cup Holder (Tube-style Armrest)..... \$20  
Part: ACCASMB2634  
For use with tube-style armrests.
- TRU-Balance 3 Personal Item Hook..... \$20  
Part: ACC125008  
Not available with back cane option.
- TRU-Balance 3 Reflector Kit..... \$35  
Part: ACC124670
- TRU-Balance 3 Backpack Holder (Pair) (A) ..... \$42  
Part: ACC147774  
Not available with back cane option.
- TRU-Balance 3 Front Mounted ..... \$95  
Cell Phone/ Media Holder, X-Grip (E)  
Part: ACC155835  
Not available with Heavy Duty armrests.
- TRU-Balance 3 Rear View Mirror ..... \$35  
Part: ACC155833  
Not available with Heavy Duty armrests.
- TRU-Balance 3 Oxygen Holder ..... \$110  
Part: ACC125006. HCPCS: E2208  
Includes accessory bar, a credit will be applied PTOIN-DV3769.
- TRU-Balance 3 Clothing Guard (Pair) ..... \$275  
Part: ACC125007

- Not available with Heavy Duty armrests.
- TRU-Balance 3 Glove Box (C)..... \$85  
Part: ACC125011
- Stealth Large Tablet Holder ..... \$565  
Part: ACC141818  
Not available with Heavy Duty armrests.
- Trilogy Vent Tray ..... \$420  
Part: ACC136791. HCPCS: E1029  
Includes accessory bar, a credit will be applied PTOIN-DV3769. Contact custom for models other than 100 and 200.
- Compact Vent Tray ..... \$1,395  
Part: ACC133889. HCPCS: E1029
- VOCSN Vent Tray ..... \$1,000  
Part: ACC2000162
- Push Handles (pair)..... \$140  
Part: ACC2005027  
Not available with oxygen holder, back cane option, or vent tray.
- Accessory Bar ..... \$90  
Part: ACC1909634  
Horizontal rail attached to rear of seat. Standard with the Trilogy vent tray and oxygen holder. Not available with back cane option.
- Large, Black, Weather Cover ..... \$240  
Part: ACCCOVR1021
- Q-Logic EX/NE Programmer ..... \$625  
Part: ELEASMB6976



**24. POSITIONING COMPONENTS**

Optional. More options are available on the [Positioning Components Order Form](#).

**24.1. Stealth Thoracic Laterals**

Swing-away laterals may interfere with TRU-Balance 3 armrests and 2-post, flip up armrests. These laterals are not available with the back cane option.

**1. Bracket**

Swing-away, UniLink, Redesign Back, pair (ULS4042000-Q).... \$590  
 Part: ST-ULS4042000-Q. E1028  
 Only available with Redesigned back assemblies.

Swing-away, Rail Mounts, pair (TWBTLTC) ..... \$590  
 Part: INDPART2677. E1028  
 Not available with Redesign back assemblies.

**2. Pads**

Curved Pads, 3"x5", pair (TWBL3X5)..... \$140  
 Part: INDPART2663. E0956

Flat Pads, 3"x5", pair (TWBLP3X5) ..... \$140  
 Part: INDPART2676. E0956

Curved Pads, 4"x6", pair (TWBL4X6)..... \$140  
 Part: INDPART2664. E0956

Flat Pads, 4"x6", pair (TWBLP4X6) ..... \$140  
 Part: INDPART2675. E0956

**3. Installation**

Install Laterals  
 Part: PTOINDV3499

Do Not Install  
 Part: PTOINDV3500

**24.2. Stealth Pelvic/Thigh Guides**

Pelvic/thigh guides may interfere with the joystick receiver block and down posts. Please provide quantity next to check box.

**1. Pelvic/Thigh Guide Bracket**

Flip Down Bracket (TWBTLADD-HW) ..... \$300 ea.  
 Part: INDPART2804. : E1028

Flip Down with LEFT Full Surface Contact ..... \$400 ea.  
 Feature (TWBADD-L & LPHW407)  
 Part: POS140428.

Removable Bracket (LP430)..... \$330 ea.  
 Part: INDPART2805. : E1028

Flip Down with RIGHT Full Surface Contact..... \$400 ea.  
 Feature (TWBADD-R & LPHW407)  
 Part: POS140429.

Removable Bracket with 1" Offset (LP430 & LHW-121)... \$430 ea.  
 Part: INDPART2854. : E1028

**2. Pelvic/Thigh Guide Adapter**

2" Height Adapter Plate (LHW-122).....\$40  
 Part: INDPART2821.  
 Height adapter is needed with TRU-Comfort seat cushions. Not available with the flip down bracket.

2 1/4" Adapter Plate (LHW-129) .....\$40  
 Part: POS140430.  
 Only available with the flip down full surface bracket.

**3. Pelvic/Thigh Guide Pad**

Short Pad, 4"x4" (TWBADD-4X4)..... \$70 ea.  
 Part: INDPART2801.

Long Pad, 4"x12" (TWBADD-4X12) ..... \$70 ea.  
 Part: INDPART2803.

Medium Pad, 4"x8" (TWBADD-4X8)..... \$70 ea.  
 Part: INDPART2802.

**4. Pelvic/Thigh Guide Installation**

Install on Left  
 Part: PTOINDV3585

Install on Right  
 Part: PTOINDV3586

Do Not Install  
 Part: PTOINDV3719

**24.3. Stealth UniLink Pelvic/Thigh Guide**

**1. UniLink Pelvic/Thigh Hardware**

UniLink, 5" Height Adjustable, ..... \$320 ea.  
 Removable Mounting Hardware  
 Part: ST-UTTP5-Q

UniLink, 7" Height Adjustable, ..... \$320 ea.  
 Removable Mounting Hardware  
 Part: ST-UTTP7-Q

**2. UniLink Pelvic/Thigh Pads**

^Please verify the appropriate billing code to be used for this item with all third-party payors.

UniLink, 4" x 4" Pad ..... \$125 ea.  
 Part: ST-UTP544056100-Q

UniLink, 4" x 12" Pad ..... \$165 ea.  
 Part: ST-UTP54D056100-Q

UniLink, 4" x 8" Pad..... \$135 ea.  
 Part: ST-UTP548056100-Q

**3. UniLink Pelvic/Thigh Additional Links**

UniLink, 1.5" Link ..... \$40 ea.  
 Part: ST-UIL001150A000-Q

UniLink, 2" Link..... \$40 ea.  
 Part: ST-UIL001200A000-Q

**4. UniLink Pelvic/Thigh Guide Installation**

Install on Left  
 Part: PTOINDV3585

Install on Right  
 Part: PTOINDV3586

Do Not Install  
 Part: PTOINDV3719



Positioning components section continued on following page.



**24. POSITIONING COMPONENTS CONTINUED**

**24.4. Stealth Stealth Positioning Chest Harness**

Upper Classic, Medium (IFPS1103) ..... \$190 ea.  
 Part: ST-IFPS1103-Q.  
 Upper Classic, Large (IFPS1104) ..... \$190 ea.  
 Part: ST-IFPS1104-Q.

Upper Contour, Medium (IFPS2103)..... \$190 ea.  
 Part: ST-IFPS2103-Q.  
 Upper Contour, Large (IFPS2104) ..... \$190 ea.  
 Part: ST-IFPS2104-Q.



Contoured and Classic Chest Harness

**24.5. Stealth Stealth Positioning Chest Straps**

Chest Strap, Medium w/ Belt Anchor ..... \$120 ea.  
 (IFCS4113)  
 Part: ST-IFCS4113-Q.

Chest Strap, Large w/ Belt Anchor ..... \$120 ea.  
 (IFCS4114)  
 Part: ST-IFCS4114-Q.



Chest Strap

**24.6. Stealth Swing-Away Headrest Mounts**

These swing-away headrest mounts are only available with the Stealth Comfort Series headrests.

**1. Facial Lateral Hardware**

Facial Lateral for Comfort Series.....\$430 ea.  
 Right Side (CPS9-R-HW)  
 Part: INDPART2647.

Facial Laterals for Comfort Series .....\$430 ea.  
 Left and Right Side (CPS9-B-HW)  
 Part: INDPART2649.

Facial Lateral for Comfort Series.....\$430 ea.  
 Left Side (CPS9-L-HW)  
 Part: INDPART2648.



Left and Right Facial Lateral with Pads

**2. Facial Lateral Pads and Accessories**

1" Ball Kit for Facial Laterals (TWB-1BK) (D) ..... \$100  
 Part: POS150532  
 Required with facial laterals to mount to comfort headrest brackets.


Facial Pad, 4"x2" (900-Q)..... \$70 ea.  
 Part: ST-900-Q. : K0108

Egg Switch Mount (ESM) (A) ..... \$70 ea.  
 Part: INDPART2815

Facial Pad, 6"x2" (920-Q)..... \$70 ea.  
 Part: ST-920-Q. : K0108

Figure 8 Bracket with 13" Rod..... \$315 ea.  
 and Egg Switch Mount (B)  
 Part: ST-F8675-13-Q

Swing-Away Joystick Mount, 18" Rod (C) .... \$295 ea.  
 Part: ST-SUS18-Q  
 Allows mounting of facial pads and attaches to the headrest



Headrest Accessories

**24.7. Stealth Elbow Support**


**1. Elbow Support Bracket**

Left, Swing-Away, Elbow Support ..... \$230 ea.  
 Bracket Tube Style Arms (ARMS-ES-L)  
 Part: POS124064.

Right, UniLink, Fixed, Elbow Support Bracket ..\$360 ea.  
 for Flip-Back, 8"-16"H Armrest (UETB4021FR-Q)  
 Part: ST-UETB4021FR-Q

Left, UniLink, Fixed, Elbow Support Bracket ..\$360 ea.  
 for Flip-Back, 8"-16"H Armrests (UETB4021FL-Q)  
 Part: ST-UETB4021FL-Q

Right, Swing-Away, Elbow Support..... \$230 ea.  
 Bracket Tube Style Arms (ARMS-ES-R)  
 Part: POS124066.



Elbow Supports

**2. Elbow Support Pad**

3"Wx3"H Elbow Pad (ARMS-ES-3X3) ..... \$70 ea.  
 Part: POS124065. : K0108

4"Wx4"H Elbow Pad (ARMS-ES-4X4)..... \$70 ea.  
 Part: POS124067. : K0108

**24.8. Stealth Foot Boxes**

Stealth Small Foot Boxes, 5"Wx8"Dx7"H ..... \$625  
 (SPLE107) (pair)  
 Part: INDPART2786.

Stealth Large Foot Boxes, 7"Wx12"Dx11"H ..... \$625  
 (SPLE109) (pair)  
 Part: INDPART2788.

Stealth Medium Foot Boxes, 6"Wx10"Dx9"H ..... \$625  
 (SPLE108) (pair)  
 Part: INDPART2787.

**24.9. Stealth Padded Shoe Holder**

Stealth Small Padded Shoe Holder 3.75"Wx6.25"D ..... \$220  
 (SPLE111) (pair)  
 Part: INDPART2783.

Stealth X-Large Padded Shoe Holder, 5"Wx10.75"D ..... \$220  
 (SPLE114) (pair)  
 Part: INDPART2785.

Stealth Medium Padded Shoe Holder, 4"Wx8.25"D ..... \$220  
 (SPLE112) (pair)  
 Part: INDPART2784.

**25. ALTERNATIVE CONTROLS AND ELECTRONICS**

**Optional.** Select a type of specialty controls below. More options are available on the [Alternative Controls and Electronics order form](#).

**25.1. Enhanced Display**

**Q-Logic 3 Enhanced Display Kit**..... \$2,045  
 Part: CTL167726

Select a color option below. Includes enhanced display, flexible gooseneck and mounting bracket. Features: 3.5" LCD display; up to 5 profiles to control drive, seat, and auxiliary functions; Bluetooth connectivity for mobile devices and mouse emulation; integrated single switch scanner; IR controls.



Enhanced Display

**Q-Logic 3 Enhanced Display Shroud**

Required with Q-Logic 3 Enhanced Display selection above. Reference colors shown on page 2.

**Black, Soft Touch** .....Standard      **Non-matching Shroud**..... \$50      **Match Shroud to Base Color** .... \$50  
 Part: DGN137728      Select color below      Part: See color below

<b>Red Corvette</b> Part: DGN137729-PRD1	<b>Grey Street</b> Part: DGN137729-PSL3	<b>Pretty in Pink</b> Part: DGN137729-PPK3	<b>Electric Blue</b> Part: DGN137729-PBU9	<b>Green Machine</b> Part: DGN137729-PGN11
<b>Back in Black</b> Part: DGN137729-PBK1	<b>White Satin</b> Part: DGN137729-PWH2	<b>Evolution Orange</b> Part: DGN137729-POR7	<b>Lemon Crush</b> Part: DGN137729-PYL4	<b>Military Green</b> Part: DGN137729-PGN8
<b>Night Train</b> Part: DGN137729-PBK3	<b>Deep Purple</b> Part: DGN137729-PPR5	<b>Midnight Blue</b> Part: DGN137729-PBU10		

**25.2. Q-Logic SCIM Module**

**Optional.**

**SCIM Module for Q-Logic 3** ..... \$850  
 Part: CTL168055

Allows use of joystick display for specialty devices. Available with Q-Logic 3 EX Joystick selection above. Up to 5 profiles to control drive, seat, and auxiliary functions, excluding sip-n-puff and single switch functions. A switch is required to operate the input device attached to the SCIM module.



SCIM Module

**25.3. Q-Logic 3 SAJ**

**Q-Logic 3 EX SAJ** ..... \$850  
 Part: CTL143980

Must be used in conjunction with Enhanced Display or Q-Logic 3 Joystick.



EX SAJ

**25.4. Head Array**

**Adult, i-Drive Tri-Array Series, 3-Switch** ..... \$5,855  
**Head Array (IDHBT200-1)**

Part: ST-IDHBT200-1-Q.

Features bluetooth, head support, 3 non-proportional proximity sensors for drive control, an egg switch for reset mode change, and i-Drive smart box.

**Adult, i-Drive Pro Series Head Array** ..... \$5,855  
**(IDHBT300-1)**

Part: ST-IDHBT300-1-Q.

Features bluetooth, head support with occipital pad and 2 swing-away spot pads for left and right side, 3 non-proportional proximity sensors, an egg switch for reset mode change, and i-Drive smart box.



Adult Pro-Series

**25.5. Sip and Puff Drive**

**1. Module**

**Sip and Puff Module for Q-Logic 3**.....\$1,695  
 Part: CTL168056 ADP: WEX5

Includes module and mounting hardware.



Sip & Puff Module

**2. Straw and Mount**

**The straw and mount do not include electrical components to operate the Sip&Puff with a head array. This component requires the selection of a Q-Logic Sip&Puff module or additional module from Stealth Products for operation.**

**Stealth Sip and Puff for Comfort Series, Left Mount (SP625-CL)** ..... \$410

Part: ST-SP625-CL-Q.

Includes straw and mounting hardware. Use with Stealth Comfort Plus headrests.

**Stealth Sip and Puff for Stealth i-Drive, Head Array (SP625-U)** ..... \$410

Part: ST-SP625-U-Q.

Includes straw and mounting hardware. Use with Stealth i-Drive Head Arrays.

**Stealth Sip and Puff for Comfort Series, Right Mount (SP622-CR)** ..... \$410

Part: ST-SP625-CR-Q.

Includes straw and mounting hardware. Use with Stealth Comfort Plus headrests.



Stealth Sip & Puff

**25.6. Sip and Puff Drive**

**Red Egg Switch (ESRED)**..... \$100  
 Part: SWTMCR01077

**Blue Egg Switch (ESBLU)** ..... \$100  
 Part: SWTMCR01079

**Black Egg Switch (ESBLK)**..... \$100  
 Part: SWTMCR01081

**Green Egg Switch (ESGRN)** ..... \$100  
 Part: SWTMCR01078

**Yellow Egg Switch (ESYEL)** ..... \$100  
 Part: SWTMCR01080



Egg Switch

Alternative controls and electronics section continued on following page.

**25. ALTERNATIVE CONTROLS AND ELECTRONICS CONTINUED**

**25.7. Switch Adapter**

2 to 1 Switch Adapter ..... \$135  
 Part: DWR1361H037  
 Connects 2 single switches to a stereo jack to act as a dual switch.

**25.8. Environmental Controls**

Q-Logic 3 ECU Module Kit ..... \$1,780  
 Part: CTL168057

**26. TRU-BALANCE 3 MULTIPLIER HARNESS**

**REQUIRED.** Multiple controllers or power accessories may require an additional multiplier harness or special order. Reference the chart below to determine the appropriate selection. The rows refer to the total amount of arm mounted devices and the columns refer to the total amount of back mounted devices. The online order configurator will automatically make this selection under specialty controls.

Secondary Multiplier (SM) ..... \$250  
 Part: XRF000107

Arm Mounts	Joystick (Hand Control)		1 Arm	2 Arm	3 Arm	4+ Arm
	Enhanced Display	0 Back	NA	SM	SM	SOF
	iAccess	1 Back	NA	SM	SM	SOF
	USB Charger	2 Back	SM	SM	SM	SOF
	SAJ	3+ Back	SOF	SOF	SOF	SOF

Power Positioning Chart

Back Mounts	Attendant Control		1 Arm	2 Arm	3 Arm	4+ Arm
	SCIM	0 Back	NA	SM	SOF	SOF
	Sip n Puff	1 Back	NA	SM	SOF	SOF
	ECU	2 Back	SM	SM	SOF	SOF
		3+ Back	SOF	SOF	SOF	SOF

Static Seat Chart

**Q-Logic 3 Multiplier Charts**