

Cruiser X4 Wheelchair

HCPCS: K0004,K0004



FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E



FIGURE F



FIGURE G



FIGURE H

- Weighs 34 pounds or less when properly equipped in the basic configuration (excluding front riggings)
- Easy to propel and transport
- Frame style eliminates seat rail guides and allows for custom backs and accessories
- Flip-back, removable arm is height adjustable from top of seat rail 10" to 14" and features a front mounted easy-to-use lever release mechanism
- Padded armrests provide added user comfort
- Height adjustable back
- Aluminum and carbon steel blended frame with combination powder coated and silver vein finish is attractive, chip resistant and easy-to-maintain
- Durable, flame retardant nylon upholstery resists mildew and bacteria
- Patented adjustable depth seat form 16" and 18"
- Three position quick release rear axles mounted on composite mag-style wheels with composite hand rims are lightweight and maintenance free
- Polyurethane rear tires and composite casters offer superior performance, smoother ride and are lightweight
- Dual over center cross frame assists provide added rigidity
- Includes swing-away footrests or elevating leg rests with tool-free adjustable length riggings
- Composite footplates with heel loops are durable and lightweight
- Precision sealed wheel bearings ensure durability and reliability
- 8" front casters with three height adjustments and are angle adjustable
- Triple rear axle positions provides easy transition of seat height from standard to hemi position
- Comes with push-to-lock wheel locks

PRODUCTS

Description

CX416ADDA-ELR	CruiserX4, 16" Adj DeskArm, ElvLg, 1/ea
CX416ADDA-SF	CruiserX4, 16" Adj DeskArm, SwgFt, 1/ea
CX416ADFA-ELR	CruiserX4, 16" Adj FullArm, ElvLg, 1/ea
CX416ADFA-SF	CruiserX4, 16" Adj FullArm, SwgFt, 1/ea
CX418ADDA-ELR	CruiserX4, 18" Adj DeskArm, ElvLg, 1/ea
CX418ADDA-SF	CruiserX4, 18" Adj DeskArm, SwgFt, 1/ea
CX418ADFA-ELR	CruiserX4, 18" Adj FullArm, ElvLg, 1/ea
CX418ADFA-SF	CruiserX4, 18" Adj FullArm, SwgFt, 1/ea
CX420ADDA-ELR	CruiserX4, 20" Adj DeskArm, ElvLg, 1/ea
CX420ADDA-SF	CruiserX4, 20" Adj DeskArm, SwgFt, 1/ea
CX420ADFA-ELR	CruiserX4, 20" Adj FullArm, ElvLg, 1/ea
CX420ADFA-SF	CruiserX4, 20" Adj FullArm, SwgFt, 1/ea

SPECIFICATIONS	Description
Warranty	Limited Lifetime
Seat Dimensions	16" (W) x 16" - 18" (D)
Weight Capacity	300 lb
Seat To Floor Height	17" 18" 19"
Carton Shipping Weight	46 lb 48 lb 50 lb
Folded Dimensions	43" (L) w/ riggings x 12" (W) x 36"
Back Height	16" - 17.5"
ACCESSORIES	Description
STDSSHK-K3	Super Hemi Kit, Cruiser 3 & X4, 1/ea, EA
STDS851	Seat Belt, Velcro Close,1/ea, EA
STDS855	Seat Belt,60"" Auto-type,1/ea, EA
STDS820	Telescoping IV Pole,W/C, 1/bx, BX
STDS856	Seat Belt,60"Velcro-type,1/ea, EA
STDS806	Anti-Fold/Theft Device for w/c,1/ea, BX
STDS850	Seat Belt, Clasp Type,1/ea, EA
STDS834	OverheadAnti-Theft Bar,Single,1/ea, EA
STDS1034	Wheelchair Cane/Crutch Holder,1/ea, EA
WASL	Limb Support,Left,1/cs, CS
WASR	Limb Support, Right,1/cs, CS
STDS601	Brake Extension,6",1/ea, EA
STDS1002	Caster,8in Sentra,Viper,Astaire,1/ea, EA
STDS860	Forwrd Stablzrs/Univ.Fit ,1pr/bx, BX
PARTS	Description
STDS4SC4X-03	LeftFrntCasterHousingW/Bearings&HrdwrCX4, EA
STDS4SC4X-04	RightFrntCasterHousingW/Bearings&HrdwCX4, EA
STDS4SC4X-17	Fork, Cruiser 4X, 1/ea, EA
STDS3J24AFL	Adj. Full Arm, Cruiser 3 / X4, Left, 1e, EA
STDS3J24AFR	Adj. Full Arm, Cruiser 3 / X4,Right,1ea, EA
STDS3J24AFR-UFR	UpprSectOnlyW/Pad,AdjFullArm,Rt,Crsr3/4, EA
STDS3J24AL-UDL	UpprSectOnlyW/Pad,AdjDskArm,Lft,Crsr3/4, EA
STDS4SC4X-08	HeightAdjustmentLever,Left,CruiserX4,1e, EA
STDS4SC4X-09	HeightAdjustmentLever,Rght,CruiserX4,1e, EA
STDS4SC4X-15	PushHandleW/Grip&AttachHrdwr,CruiserX4, EA
STDS4SC4X-10	Back Upholstery, 16", CruiserX4, 1/ea, EA
STDS4SC4X-11	Back Upholstery, 18", CruiserX4, 1/ea, EA
STDS4SC4X-12	Back Upholstery, 20", CruiserX4, 1/ea, EA
STDS4SC4X-01	RearWheelComplW/Bearings&Hndrm,Cruiser4, EA
STDS4SC4X-02	Quick Release Rear Axle, CruiserX4, 1/e, EA
STDS4SC4X-16	RearWheelAxleReceptacle, CruiserX4, 1/e, EA
STDS4SC4X-05	Side Frame, Left, CruiserX4, 1/ea, EA
STDS4SC4X-06	Side Frame, Right, CruiserX4, 1/ea, EA
STDS4S1069	HangerPins,CS/SilverSport/Cruiser, 4/pkg, PKG
STDS3J24X6-AD	Cross Frame,Cruiser3,Adj Depth5D,16",1ea, EA