

Rebel Wheelchair



HCPCS:




FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E



FIGURE F

- Quick-release rear wheels allow for easy storage or transport (Figure A)
- Fold-down back and folding frame allows transport in almost any trunk (Figure B)
- Easily removable swing-away footrests are also convenient for storage
- Tool-free removable front casters
- Desk-style armrests make transfers easy and accessibility better
- Lightweight and extremely maneuverable, easy to navigate tight areas
- Large accommodating seating surface
- Carton size makes the Rebel perfect for retail

PRODUCTS

Description

RTLREB18DDA-SF	Rebel W/C,18",FlameRd,DetDsk,SwgFt,1/ea
----------------	---

SPECIFICATIONS

Description

Back Height	17"
Carton Shipping Weight	42 lbs
Overall Height	35"
Overall Length w/ Riggings	43"
Overall Width	25.5" (open) 10.5" (closed)
Seat Dimensions	18" (W) x 16" (D)
Seat To Floor Height	19.5"
Seat Width	18"
Warranty	Limited Lifetime
Weight	37 lbs
Weight Capacity	300 lbs

ACCESSORIES

Description

STDS1040S	Cup Holder, Swivel,1 BX,6/CS, BX
STDS850	Seat Belt, Clasp Type,1/ea, EA
STDS855	Seat Belt,60" Auto-type,1/ea, EA
STDS851	Seat Belt, Velcro Close,1/ea, EA
STDS856	Seat Belt,60"Velcro-type,1/ea, EA
STDS806	Anti-Fold/Theft Device for w/c,1/ea, BX
STDS801GT	WheelLockExtension,ViperGT,SSP2,1/ea, EA
STDS831	Heel Loops, 1pr, BX
STDS1034	Wheelchair Cane/Crutch Holder,1/ea, EA
STDS803	Oxygn Tank Holdr for W/C,1/cs, CS

PARTS	Description
RTLREB-01	BackUpholsteryWithScrewsRTLREB,1 ea, EA
RTLREB-02	BRCaneW/PuHandle&AttachingBoltRTLREB,1e, EA
RTLREB-03	FoldingCaneHingeRTLREB,1ea, EA
RTLREB-05	QuickReleaseAxleRTLREB,1ea, EA
RTLREB-04	SeatUpholsteryW/ScrewsRTLREB,1ea, EA
RTLREB-06	RearWheelRTLREB,1ea, EA
RTLREB-07	FrontAxleBoltRTLREB,1ea, EA
RTLREB-09	FrontFork/CasterAssyRTLREB,1ea, EA
RTLREB-10	FrontForkFrameKnobRTLREB,1pr, PR
RTLREB-11	ForkLockingClipRTLREB,1pr, PR
RTLREB-16	SlideTubeRTLREB,1ea, EA
RTLREB-17	FramePlugSetRTLREB,1pr, PR
RTLREB-18	CalfStrapRTLREB,1ea, EA
RTLREB-19	LeftDeskArmAssyRTLREB,1ea, EA
RTLREB-20	RightDeskArmAssyRTLREB,1ea, EA
STDS8M4603	BrakeTipSleeve,Cruiser,SSP, 1/ea, EA
STDS4Y4729	Rear Wheel Bearing, SSPRB, CX4, 2ea/pk, PKG
STDS4020	Tipping Sleeve 1/ea, EA
STDSDSKVB	UniversalWheelChairVinylDeskArmPad,1/ea, EA
STDS8M241L	Brake Fixed Arm,Left,SSP,1ea, EA