

Company:

Ship To:

Account Number:

Date:

Address:

Quote: PO Number:

City:

Province:

Clinician / Sales Rep:

Postal Code:

Phone:

Marked For:

Email/Fax:

Associated Quote/Order#:

Submit Form

Effective: 5/22/2023

CR45 - Transit Included _____ **\$3150**

300 lb weight capacity, Tilt range: 0° anterior to 45° posterior

ADP Code: WM0001935

Frame

Frame Width

Select box corresponding to desired Frame Width

Frame Width								
14" N/C	15" N/C	16" N/C	17" N/C	18" N/C	19" WAMA N/C	20" WAMA N/C	21" WAMA N/C	22" WAMA N/C

Seating Setup

Make one selection from either Option 1 at right or Option 2 below based on desired seating setup

Setup By Ki Mobility (Option 1) _____ N/C

Select Seat Depth corresponding to desired Frame Depth Range row and Backrest Offset column. Note: All Seat Depths in a column and within the same frame depth range are achievable with future adjustment.

Setup by Frame Depth (Option 2) _____ N/C

Select Frame Depth within desired Frame Depth Range

Option 2: Setup By Frame Depth

Frame Depth Range	Frame Depth				
Short	14" N/C	15" N/C	16" N/C	17" N/C	18" N/C
Medium	16" N/C	17" N/C	18" N/C	19" WAMB N/C	20" WAMB N/C
Long	18" N/C	19" WAMB N/C	20" WAMB N/C	21" WAMB N/C	22" WAMB N/C

Seat Pan

Growing Adjustable Aluminum Seat Pan _____ \$145


Omit Seat Pan _____ N/C

Transit

Omit Transit _____ -\$50

Option 1: Setup By Ki Mobility

Frame Depth Range	Backrest Offset A						
	-2"	-1"	Flush	+1"	+2"	+3"	+4"
Short	Seat Depth						
	16"	15"	14"				
	17"	16"	15"	14"			
	18"	17"	16"	15"	14"		
Medium	19"	18"	17"	16"	15"	14"	
	20" ¹	19"	18"	17"	16"	15"	14"
	21" ¹	20" ¹	19"	18"	17"	16"	15"
	22" ²	21" ¹	20" ¹	19"	18"	17"	16"
Long	20" ¹	19"	18"	17"	16"	15"	14"
	21" ¹	20" ¹	19"	18"	17"	16"	15"
	22" ²	21" ¹	20" ¹	19"	18"	17"	16"
		22" ²	21" ¹	20" ¹	19"	18"	17"
		22" ²	21" ¹	20" ¹	19"	18"	
			22" ²	21" ¹	20" ¹	19"	



Tilt Mechanism

Hand Tilt - Dual	N/C
Foot Tilt - Right	N/C
Foot Tilt - Left	N/C

Backrest

Backrest

Stroller Handle Back Post (Short 20")	N/C
Stroller Handle Back Post (Tall 24")	N/C
Fixed Height Back with Adjustable Handle (Short 22")	N/C
Fixed Height Back with Adjustable Handle (Tall 26")	N/C
Adjustable Height Back with Adjustable Handle (Short 18" - 25")	N/C
<i>Not available with O2 Holder</i>	
Adjustable Height Back with Adjustable Handle (Tall 20" - 27")	\$36
Reclining Back with Adjustable Handle (Short 22") (Non-Folding)	\$1232
<i>Not available with O2 Holder. If no Gas Spring selection is made, Standard Spring will be used.</i>	
Reclining Back with Adjustable Handle (Tall 26") (Non-Folding)	\$1232
<i>Not available with O2 Holder. If no Gas Spring selection is made, Standard Spring will be used.</i>	

Backrest Height

18" 19" 20" 21" 22" 23" 24" 25" 26" 27" N/C

Stroller Handle Back Post Option

Available only with Stroller Handle Back Post

Removable Stroller Handle Extension	\$173
-------------------------------------	-------

Other Backrest Options

Non Folding Back	N/C
<i>Standard with Reclining Back</i>	
Adjustable Height Rigidizer Bar	\$75

Reclining Back Options

Select option below if Reclining Back is selected on chart above. If no selection is made, Standard Spring will be used.

Standard Spring	N/C
Heavy Duty Spring	N/C

Footrest

Hanger / Front Frame Bend

Hanger	60°	70°	80°	90°	Retail
Extension Tube Mount			1		N/C
Front Mount			1		N/C
Pro Elevating Leg Rests WAMK					\$376
Omit Hangers					-\$25

1 7" and 8" Casters may interfere with less than 5° tilt

PRO ELR Cover

Available only on PRO ELR

Gel Pad Covers (Pair)	\$557
-----------------------	-------

3" Extension

Available only with Front Mount Hangers

3" Extension (Provided as Parts)	\$44
----------------------------------	------

Seat to Footrest Length

Select Extension and Seat to Footrest Length in row corresponding to desired Footrest. Footrest Length of at least 2.5" shorter than Front Seat Height recommended. Extension Tube Hanger measurements based on 70° Hanger Bend.

Footrest	Extension	Seat to Footrest Length											
		10"	10.5"	11"	11.5"	12"	12.5"	13"	13.5"	15"	18"		
Extension Tube Hangers (60°, 70°, 80°)	Short (10" - 13.5")	10"	10.5"	11"	11.5"	12"	12.5"	13"	13.5"	15"	18"		
	Medium (11.5" - 15.5")	11.5"	12"	12.5"	13"	13.5"	14"	14.5"	15"	15.5"			
	Long (14" - 18")	14"	14.5"	15"	15.5"	16"	16.5"	17"	17.5"	18"			
Front Mount ¹ Hangers	Short Front Mount Hanger	Standard Selections									Requires 3" Extension for \$44		
		4" ²	4.5" ²	5"	5.5"	6"					8.5"	9"	9.5"
	Medium Front Mount Hanger	6.5"	7"	7.5"	8"						10"	10.5"	11"
		4" ²	4.5" ²	5"	5.5"	6"	6.5"	7"			10.5"	11"	11.5"
	Long Front Mount Hanger	7.5"	8"	8.5"	9"	9.5"	10"				12"	12.5"	13"
		5"	5.5"	6"	6.5"	7"	7.5"	8"	8.5"		12.5"	13"	13.5"
		9"	9.5"	10"	10.5"	11"	11.5"	12"		14"	14.5"	15"	
PRO ELR Low Pivot Setting	PRO ELR Low Pivot Setting Short (10.5" - 15.5")	10.5" - 15.5"											
	PRO ELR Low Pivot Setting Long (13" - 19")	13" - 19"											
PRO ELR High Pivot Setting	PRO ELR High Pivot Setting Short (8.5" - 17")	8.5" - 17"											
	PRO ELR High Pivot Setting Long (11.5" - 17.5")	11.5" - 17.5"											

¹ See Minimum Front Seat Height chart below

² Available only with 90° Front Mount Hangers

Minimum Front Seat Height with Front Mount Hangers						
Hanger Length	Short		Medium		Long	
	70°	80/90°	70°	80/90°	70°	80°
Hanger Bend	70°	80/90°	70°	80/90°	70°	80°
Minimum Front Seat Height with Standard Selection	13"	13"	13"	13.5"	14.5"	15.5"
Minimum Front Seat Height with 3" Extension	13.5"	14"	15.5"	16"	17"	18"

Hanger Release

If no selection is made 4-Way - In line Hanger Release will be used

Classic	N/C
4-Way - In-Line	N/C
Not available with Pro ELR	
4-Way - Out	N/C
Not available with Pro ELR	

Heel Loops

Not available with One Piece Flip Up Angle Adjustable Footplate

Heel Loops	\$65
Adjustable Heel Loops	\$45
Omit Heel Loops	N/C

Calf Straps

Velcro Adjustable	\$52
Padded Velcro Adjustable	\$85
Bodypoint Padded AeroMesh	\$134

Footplates

Composite	N/C
Not available with Front Mount Hanger	
Composite Angle Adjustable WAMJ	\$195
Not available with 14" Frame Width	
Aluminum Angle Adjustable WAMJ	\$172
One Piece Flip Up Angle Adjustable WAMJ	\$238
Not available with Front Mount Hanger or PRO ELR	

Residual Limb Support

	Limb Support	Price	Omit	Cover		Pad Size					
				Smooth Cover Left	Air Mesh Cover Left	5X8 Pad Left Side	7X10 Pad Left Side	9x10 Pad Left Side	9x14 Pad Left Side	11x10 Pad Left Side	11x14 Pad Left Side
Left	Left Limb Support	\$462	Omit Left Hanger	Smooth Cover Left	Air Mesh Cover Left	5X8 Pad Left Side	7X10 Pad Left Side	9x10 Pad Left Side	9x14 Pad Left Side	11x10 Pad Left Side	11x14 Pad Left Side
Right	Right Limb Support	\$462	Omit Right Hanger	Smooth Cover Right	Air Mesh Cover Right	5X8 Pad Right Side	7X10 Pad Right Side	9x10 Pad Right Side	9x14 Pad Right Side	11x10 Pad Right Side	11x14 Pad Right Side

Seat Height, Casters and Rear Wheels

Select Seat Height corresponding to desired Rear Wheel and Caster Size

Rear Wheel Size	Caster Size					
	4"	5"	6"	6" x 2"	7"	8"
12" with Pneumatic or Pneumatic with Airless Insert Tire	Seat Height					
		15.5"				
	16"	16"	16"	16"		
	16.5"	16.5"	16.5"	16.5"		
		17"	17"	17"	17"	
12" with Full Poly Tire		17.5"	17.5"	17.5"	17.5"	
	15"	15"				
	15.5"	15.5"				
	16"	16"	16"	16"		
	16.5"	16.5"	16.5"	16.5"		
16" with Pneumatic, Pneumatic with Airless Insert or Full Poly Tire	14"					
	14.5"					
	15"	15"				
	15.5"	15.5" ¹				
		17"		17" ²	17" ²	
		17.5"	17.5"	17.5"	17.5"	
		18"	18"	18"	18"	18"
		18.5"	18.5"	18.5"	18.5"	18.5"
16" with Low Poly Tire	19"	19"	19"	19"	19"	
	14"					
	14.5"					
		16"				
	16.5"	16.5"	16.5"	16.5"		
		17"	17"	17"	17"	
		17.5"	17.5"	17.5"	17.5"	
	18"	18"	18"	18"	18"	

Rear Wheel Size	Caster Size						
	4"	5"	6"	6" x 2"	7"	8"	
20"	Seat Height						
	16"	16"	16"	16"			
	16.5"	16.5"	16.5"	15.5"			
		17"	17"	17"	17"		
		17.5"	17.5"	17.5"	17.5"		
			19" ¹				
		19.5"	19.5"	19.5"	19.5"	19.5"	
		20"	20"	20"	20"	20"	
		20.5"	20.5"	20.5"	20.5"	20.5"	
	22"		16.5"		16.5"		
		17"	17"	17"	17"		
		17.5"	17.5"	17.5"	17.5"		
		18"	18"	18"	18"	18"	
		18.5"	18.5"	18.5"	18.5"	18.5"	
		20.5"	20.5"	20.5"	20.5"	20.5"	
24"				17" ¹	17" ¹		
		17.5"	17.5"	17.5"	17.5"		
		18"	18"	18"	18"	18"	
		18.5"	18.5"	18.5"	18.5"	18.5"	
	19"	19"	19"	19"	19"		

¹ Available only with Full Poly tire
² Not available with Full Poly tire

Casters and Rear Wheels

Casters

		Caster Size					
		Price	4"	5"	6"	7"	8"
Caster Type	1" Poly	N/C					
	1.5" Poly	\$90					
	2" Poly	\$109					

Axle

Fixed	N/C
Quick Release	N/C

Wheels

	Rear Wheel Size					
	Price	12"	16"	20"	22"	24"
Mag	N/C					

Wheel Locks

Wheel Lock

Push to Lock	N/C
Pull to Lock	N/C
<i>Not available with 12" or 16" Full Poly tires</i>	
Attendant Foot Lock	\$102
Omit Wheel Locks	N/C

Rear Tire Type

Full Poly	N/C
Pneumatic	N/C
Pneumatic with Airless Insert	\$207
Low Poly	N/C
<i>Available only with 16" wheel size</i>	

Handrims

Handrims not available on 12" or 16" wheel sizes

Standard Handrims

Aluminum Anodized	N/C
Plastic Coated	\$151
Projection (8 Vertical Projections) WAMQ	\$325
<i>Available only with full profile tires in combination with 22" and 24" Mag wheels</i>	
Omit Handrim	N/C

Armrests

Armrest Type	Length	Price
Height Adjustable T-Arm Tall (11.5"-15")	Desk	\$314
	Full	
Height Adjustable T-Arm (8.5"-13.5")	Desk	\$314
	Full	
Height Adjustable T-Arm Low (6.5"-10.5") ¹	Desk	\$314
	Full	

¹ If requested configuration is unavailable, Standard Height T-Arm will be used
² Minimum armrest height with Reclining Backrest is 12"

Wheel Lock Options

Extension Handles	\$91
<i>Available only with Push to Lock or Pull to Lock Wheel Locks</i>	

Armrest Type	Length		Price
	Desk	Full	
Angle Adjustable Locking Extendable Flip Up ²			\$312
Omit Armrests			N/C

Arm Pads

Standard Arm Pad	N/C
Waterfall Arm Pad	\$35
<i>Not available with Angle Adjustable Locking Extendable Flip Up Armrests</i>	
Foam Grip	N/C
<i>Available only with Angle Adjustable Locking Extendable Flip Up Armrest</i>	

Accessories

Removable Side Guard

Composite Side Guards	\$175
<i>Not available with T-Arm Armrests</i>	

Anti-tips

Rear Anti-Tips	\$145
<i>Required with Reclining Backrest</i>	
Omit	N/C

Pelvic Positioning Belt

1 1/2" Auto Buckle Pelvic Positioning Belt	\$80
1 1/2" Padded Auto Buckle Pelvic Positioning Belt	\$185
2" Airline Buckle Pelvic Positioning Belt	\$150
Small Bodypoint Evoflex (16-24")	\$225
Medium Bodypoint Evoflex (19-29")	\$225
Large Bodypoint Evoflex (25-37")	\$225
Belt Mounting Band Clamp	\$25


O2 Holder

O2 Holder WAMX	\$365
----------------	-------

Other Accessories

Tool Kit	\$60
Removable Underseat Pouch	\$75
Ki Mobility Backpack	\$90
Black Neoprene Impact Guards	\$102
Tilt Stop (Pair)	\$40
IV Pole	\$188
Canopy with Hardware	\$335

Not available on frame widths greater than 18". See image below.



Frame Color

Standard Colors

Black	N/C
Military Green	N/C
Candy Blue	N/C
Charcoal Gray	N/C
Candy Purple	N/C
Iridescent Burgundy	N/C
Candy Red	N/C
Iridescent Navy	N/C
Copper	N/C
Shadow Black	N/C

Special Order Colors

Extended lead time

Blue Skies	\$90
Emerald Green	\$90
Bubblegum Pink	\$90
Retro Red	\$90
Royal Blue	\$90
Granny Smith	\$90
Smooth White	\$90
Metallic Orange	\$90
Teal	\$90
Yellow	\$90

Axiom Backs

The following is an abbreviated list of Seating Options. Refer to full Seating and Accessories Order Forms within the Ki Mobility website for more selections.

Axiom Posterior Lateral Back (PL)

Select back that is the same or up to 2" narrower than the width of the frame. Hardware is included with back. Weight capacity: 300 lbs.

		Retail	Width		
			16"	18"	20"
Height	Mid	14"	AXPL1614	AXPL1814	AXPL2014
		16"	AXPL1616	AXPL1816	AXPL2016
		18"	AXPL1618	AXPL1818	AXPL2018
	Tall	20"	AXPL1620	AXPL1820	AXPL2020
		22"	AXPL1622	AXPL1822	AXPL2022
			\$1006	\$1006	\$1006

Back Hardware

DEX Hardware BTX Hardware FXT Hardware _____ N/C

Axiom Deep Lateral Back (DL)

Select back that is the same or up to 2" narrower than the width of the frame. Hardware is included with back. Weight capacity: 300 lbs.

		Retail	Width		
			16"	18"	20"
Height	Mid	14"	AXDL1614	AXDL1814	AXDL2014
		16"	AXDL1616	AXDL1816	AXDL2016
		18"	AXDL1618	AXDL1818	AXDL2018
	Tall	20"	AXDL1620	AXDL1820	AXDL2020
		22"	AXDL1622	AXDL1822	AXDL2022
			\$1216	\$1216	\$1216

Back Hardware

DEX Hardware BTX Hardware FXT Hardware _____ N/C

Axiom Back Accessories

Headrest

Choose headrest and hardware combination in chart below. Not available with Low Backs.

Headrest Pad	Description	Item Code	Hardware	
			Kio Code	Retail
			Multi Axis Removable Hardware with Horizontal Bar	
			HEDHDW201 / HEDHDW217	
	Comfort Plus Headrest 8"	CS08MA	CS08MA	\$514.29
	Comfort Plus Headrest 10"	CS10MA	CS10MA	\$514.29
	Comfort Plus Headrest 14"	CS14MA	CS14MA	\$514.29

Axiom SP Cushions

The following is an abbreviated list of Seating Options. Refer to full Seating and Accessories Order Forms within the Ki Mobility website for more selections.

Axiom SP - Visco

ADP Code for Cushions Sized 14" - 20" Wide is SE0001397 Cushions sized 21" - 22" Wide is SE0001398

		Width																														
		14"			15"			16"			17"			18"			19"			20"			21"			22"						
		Cushion	Cover		Cushion	Cover		Cushion	Cover		Cushion	Cover		Cushion	Cover		Cushion	Cover		Cushion	Cover		Cushion	Cover								
	Inner	Outer	Inner	Outer	Inner	Outer	Inner	Outer	Inner	Outer	Inner	Outer	Inner	Outer	Inner	Outer	Inner	Outer	Inner	Outer	Inner	Outer	Inner	Outer								
Retail	\$745	\$157			\$745	\$157			\$745	\$157			\$745	\$157			\$745	\$157			\$745	\$157			\$906	\$157			\$906	\$157		
Height	14"																															
	15"																															
	16"																															
	17"																															
	18"																															
	19"																															
	20"																															
	21"																															
22"																																

Item Code for Axiom SP-Visco Cushions are XSP(Cushion Width and Depth)V example XSP1415V. For the Extra Outer Cover XSPC(Cushion Width and Depth)V example XSPC1415V. For the Extra Inner Cover XSPI(Cushion Width and Depth)V example XSPI1414V.

Star Cushions

Star Stabil-Air Cushion (4" Tall Cells) and Extra Standard / Incontinent Cover

ADP Code for Stabil-Air Cushions, All Sizes: SE0001562

		Width																	
		13.5"			15.5"			17.25"			19.25"			21"			22.75"		
		Cushion w/ Inco Cover	Extra Cover		Cushion w/ Inco Cover	Extra Cover		Cushion w/ Inco Cover	Extra Cover		Cushion w/ Inco Cover	Extra Cover		Cushion w/ Inco Cover	Extra Cover		Cushion w/ Inco Cover	Extra Cover	
	Standard	Incontinent	Standard	Incontinent	Standard	Incontinent	Standard	Incontinent	Standard	Incontinent	Standard	Incontinent	Standard	Incontinent	Standard	Incontinent	Standard	Incontinent	
Retail	\$1195	\$89	\$103	\$1195	\$89	\$103	\$1195	\$89	\$103	\$1195	\$89	\$103	\$1195	\$89	\$103	\$1360	\$117	\$132	
Height	13.5"																		
	15.5"																		
	17.25"																		
	19.25"																		
	21"																		

Item Code for the Star Stabil-Air Guard Cushion (4" Tall Cells) Is ESSA4 (Cushion width and Depth)-1 example ESSA41517-1. Item code for the extra standard covers are ESZC4(width and depth) example ESZC42117. Or for the Incontinent cover is ESIZC4(width and depth) ESIZC42117