



# Simplified Edge 3 and Q6 Edge 2.0 TB Flex Order Form

## BASE MODEL

300 lb. weight capacity with standard 6mph programmable motors and 8 amp off board charger. Dimensions of base without legrests: 24"Wx35.5"L. XLR USB Charger is included with Edge 3.

Quantum Edge 3.....\$8,450      Quantum Q6 Edge 2.0..... \$8,170  
*Standard with Switch Operated Fender Lights*

## PATIENT WEIGHT

0 - 200 lbs      201 - 300 lbs

## BASE COLOR

Red Corvette	White Satin	Electric Blue	Sugar Plum (E3 only)
Cherry Bomb (E3 only)	Deep Purple	Blue Jean (E3 only)	Pretty in Pink
Back in Black	Raspberry Beret (E3 only)	Lemon Crush	
Night Train	Evolution Orange	Green Machine	
Grey Street	Midnight Blue	Military Green	

## DRIVE MOTORS

6 mph Drive Motors ..... Standard      Accu-Trac Motor Technology ..... \$995

## BASE MOUNTED FOOT PLATFORM

Base Mounted Foot Platform ..... Standard      Omit base mounted foot platform ..... No Charge

## DRIVE WHEELS

3"x14" Flat Free Black Tires ..... Standard      3"x14" Pneumatic Black Tires ..... No Charge

## Drive Wheel Accent Color

Matching ..... Std/\$85  
 Standard on Edge3, \$85 on Edge 2.0

## ELECTRONICS

NE 4-Key Joystick.....Standard      NE+ 6-Key Joystick ..... No Charge

Q-Logic 3 .....\$900      Q-Logic 3 E .....\$620  
*Standard with Colour Match Shroud*

## JOYSTICK MOUNTING BRACKETS

Fixed, In-line, Bracket ..... Standard

Adjustable Swing Away ..... \$265

Offset Swing Away ..... \$265

Basic Swing Away ..... \$265

Left Side Mount      Right Side Mount

## JOYSTICK HANDLES

Mushroom Handle (A) ..... \$120	Soft Ball Handle (E)..... \$120
Chin Cup (B)..... \$90	Extended I-Handle (F) ..... \$140
T-Bar (C) ..... \$160	3" Goal Post Handle ..... \$180
Medium Ball Handle (D)..... \$70	4" Goal Post Handle ..... \$180

## SEATING SYSTEM

TB Flex Seat..... \$650  
 With TB Flex 24"H Quick Angle Adjustable, Flip Forward back canes

### Seat Dimension - Width

15" W (#1005)	17" W (#1007)	19" W (#1009)	21" W (#1011)
16" W (#1006)	18" W (#1008)	20" W (#1010)	22" W (#1012)

### Seat Dimension - Depth

14" D      15" D      16" D      17" D      18" D      19" D      20" D

### Seat to Floor

Lowest (16.5")      Middle (17.5")      Highest (18.25")

**SEATING SYSTEM (CONTINUED)**

**1b. Single Post, Height Adjustable, Drop In Armrests**

<input type="checkbox"/> L <input type="checkbox"/> R	Full Length, Drop In Armrests ..... No Charge <i>Only available with Standard Armpad (no armpad selection required)</i>	<input type="checkbox"/> L <input type="checkbox"/> R	Desk Length, Drop In Armrests ..... No Charge <i>Only available with Standard Armpad (no armpad selection required)</i>
---	--	---	--

**1b. Two-Post, Flip Back, Height Adjustable, Removable Armrests**

<input type="checkbox"/> L <input type="checkbox"/> R	Full Length, Two-Post Armrests.....\$145 ea.	<input type="checkbox"/> L <input type="checkbox"/> R	Desk Length, Two-Post Armrests .....\$145 ea.
	Straight Armpads..... Standard		Waterfall Armpads..... No Charge

**Standard Lap Belts**

38" Lap Belt     50" Lap Belt     60" Lap Belt     70" Lap Belt     80" Lap Belt     90" Lap Belt

**LEGRESTS**

**Center Mount Foot Platform**

**Foot Platform Assembly**

Manually Adjustable Center Mount Foot Platform... \$265

**Lower Extension and Footplates**

10"-16" Extension, 5"x6" Footplates	10"-16" Extension, 5"x11" Footplates	13"-19" Extension, 5"x8" Footplates
10"-16" Extension, 5"x8" Footplates	13"-19" Extension, 5"x6" Footplates	13"-19" Extension, 5"x11" Footplates

**Center Mount Calf Support**

Small Calf Pads (pair) 6"Wx6"H..... \$210*	Large Calf Pads (pair) 8" taper to 6"Wx8"H..... \$210*
--	--

**Modular Footplate Extension Components**

Width Extension for 8"D Footplates ..... \$55	4" Outside Edge (pair) ..... \$55 <i>Only available with 6"D footplates and 8"D footplate extension.</i>
Width Extension for 11"D Footplates ..... \$55	6" Outside Edge (pair) ..... \$55 <i>Available with 8"D and 11"D footplates and 11"D footplate extension.</i>
2" Bell-Shaped Length Extension ..... \$55 <i>Not compatible with width extensions. Can be used with 6" outside edge with 11" footplates.</i>	

**60°-90° Swing-Away Legrests**

<input type="checkbox"/> P <input type="checkbox"/> L <input type="checkbox"/> R 60°-90° Swing-Away Legrest..... \$285	<input type="checkbox"/> P <input type="checkbox"/> L <input type="checkbox"/> R 60°-90° Swing-Away Legrest, Short..... \$285.
--	--

**2. Lower Extension**

<input type="checkbox"/> P <input type="checkbox"/> L <input type="checkbox"/> R 11"-13" Length	<input type="checkbox"/> P <input type="checkbox"/> L <input type="checkbox"/> R 13"-16.5" Length	<input type="checkbox"/> P <input type="checkbox"/> L <input type="checkbox"/> R 15.5"-19.5" Length
---	---	---

**3b. Angle Adjustable Footplate(s)**

<sup>1</sup>The 6"x11" footplates may cause interference with seat widths less than 18". <sup>2</sup>The 8"x11" footplates may cause interference with seat widths less than 20".

<input type="checkbox"/> P <input type="checkbox"/> L <input type="checkbox"/> R 4"Wx6"D ..... \$87.50 ea.	<input type="checkbox"/> P <input type="checkbox"/> L <input type="checkbox"/> R 6"Wx8"D <sup>1</sup> ..... \$87.50 ea.
<input type="checkbox"/> P <input type="checkbox"/> L <input type="checkbox"/> R 5"Wx6"D..... \$87.50 ea..	<input type="checkbox"/> P <input type="checkbox"/> L <input type="checkbox"/> R 6"Wx11"D <sup>1</sup> ..... \$87.50 ea.
<input type="checkbox"/> P <input type="checkbox"/> L <input type="checkbox"/> R 5"Wx8"D ..... \$87.50 ea.	<input type="checkbox"/> P <input type="checkbox"/> L <input type="checkbox"/> R 8"Wx11"D <sup>2</sup> ..... \$87.50 ea.

**BATTERIES**

NF-22, 55 Ah Gel Cell, Interceptor Battery ..... \$300ea. <i>Up to 3 extra miles additional range!</i>	NF-22, 50 Ah Gel Cell, MK Battery ..... \$400 ea.
---	---