



Pride Mobility Products Co.

5096 South Service Rd, Beamsville, ON, L0R 1B3

Phone: 888-570-1113 | Fax: 866-514-1303 | Email: productorders@pridemobility.com

Q6 Edge® 3 Order Form with TRU-Balance® 4 Positioning

Account Number: _____ Date: _____

Provider Name: _____

Contact: _____

Phone: _____ Fax: _____

Email: _____

PO Number: _____

Marked for: _____

Ship to Address: _____

City: _____ State: _____ Zip: _____

EDGE 3



INTRODUCTION

This form is interactive when viewed with [Adobe Acrobat Reader](#). Complete the form by placing checks in the desired boxes and provide information in the interactive fields. Buttons shown at the bottom of the form may be utilized to print or submit the order form through a desktop email application. To email via a web-based application, please 'Save As' and attach the PDF to your email.

This order form contains a large variety of options to fulfill various patient needs. Descriptions and section notes, such as optional or required, are included to help you complete your order. Other part numbers can be found on the quote acknowledgment after this order form is submitted. Please contact Quantum® Sales at 888-570-1113 if assistance is needed. Send the completed order form by fax (866-514-1303) or email (productorders@pridemobility.com). Incomplete forms may delay the quote or order. Customer service will contact you if the order is incomplete or if there are compatibility issues. If special order requests are needed, be sure that the Patient Information section is completed

All prices are MSRP. Prices, specifications, part numbers, and availability are subject to change without notice. Prices and part numbers as shown are only available when configured as a complete power chair. Please contact technical service for accurate parts ordering.

Thank you for choosing Quantum!


PATIENT INFORMATION

Weight: _____ Height: _____

A. Shoulder Width: _____ D. Max Sitting Height: _____ G. Thigh Depth: _____

B. Chest Width: _____ E. Shoulder Height: _____ H. Lower Leg Length: _____

C. Hip Width: _____ F. Axilla Height: _____ I. Elbow Height: _____



1. BASE MODEL

REQUIRED. 300 lb. weight capacity with standard 6mph programmable motors and 8 amp off board charger. XLR USB charger is included. Dimensions of base without legrests: 24"Wx35.5"L

Quantum Q6 Edge 3 3MP-SS \$8,450
Part: Q6 EDGE 3. ADP: WP0001575 (ADULT), WP0001576 (PEDIATRIC)

2. DRIVE MOTORS

REQUIRED.

6 mph 2 Pole Drive Motors..... Standard
Part: MOT147433

6 mph, 2 Pole, Accu-Trac Motor Technology..... \$995
Part: MOT147881.
Requires Q-Logic 3 electronics. Delivers enhanced traction and tracking, especially with switched drive controls.

7 mph, 4 Pole Drive Motors..... \$700
Part: MOT1902455

Requires Q-Logic 3 electronics.

7 mph, 4 Pole, Accu-Trac Motor Technology..... \$1,695
Part: MOT2008501
Requires Q-Logic 3 electronics. Delivers enhanced traction and tracking, especially with switched drive controls.

3. BASE COLOR

REQUIRED. Select one shroud color. Colors shown here may differ from actual product depending on monitor/printer calibration used.

3.1a. Standard Shrouds

 **Red Corvette**
Part: DGN1900199-PRD1

 **Back in Black**
Part: DGN1900199-PBK1

 **Night Train**
Part: DGN1900199-PBK3

 **Grey Street**
Part: DGN1900199-PSL3

 **Deep Purple**
Part: DGN1900199-PPR5

 **Raspberry Beret**
Part: DGN1900199-PPK4

 **Pretty in Pink**
Part: DGN1900199-PPK3

 **Evolution Orange**
Part: DGN1900199-POR7

 **Midnight Blue**
Part: DGN1900199-PBU10

 **Electric Blue**
Part: DGN1900199-PBU9

 **Blue Jean**
Part: DGN1900199-PBU11

 **Lemon Crush**
Part: DGN1900199-PYL4

 **Green Machine**
Part: DGN1900199-PGN11

 **Military Green**
Part: DGN1900199-PGN8

3.1b. Hydrographic Shrouds

 **Military Digital Camo**..... \$325
Part: DGN1801101-WCM7

 **Faded Glory**..... \$325
Part: DGN1801101-WFL2

 **Carbon Fiber**..... \$325
Part: DGN1801101-WCB1

4. BASE MOUNTED LIGHTS

Optional. Select one option. Seat Mounted, LED, Full Light Package option can be found later on this form.

Switch Fender Lights..... No Charge
Part: DGN154117
Fender lights operate by pressing included single switch button. Fender lighting is not available with Single Actuator Control through Toggle Box.

Base Mounted, LED, Full Light Package \$770
Part: XRF000129
Requires Q-Logic 3 EX Joystick for Full Lights or Q-Logic 3e Joystick.

5. QUANTUM BACKUP CAMERA

Optional. The backup camera may interfere with armrest mounted accessories such as media device holder.

Quantum Backup Camera \$475
Part: MEC1910365

1. Optional Mounting Bracket

Backup Camera Headrest Bracket..... \$30
Part: MEC1910373
Works with .5" x .5" vertical supports only.

Backup Camera Track Mount Bracket \$35
Part: MEC2104227

2. Optional Display Mounting Bracket

Backup Camera Q-Logic Joystick Mount Bracket \$145
Part: MEC2106123



Backup Camera

6. DRIVE WHEELS

REQUIRED. Select a type of drive wheels and accent color. Knobby tires increase overall base width by 1.5".

6.1. Drive Wheels

3"x14" Flat Free Black Tires Standard
Part: PTOINDV3525 (pair)

4"x14" Knobby, Flat Free Black Tires \$60 ea.
Part: WHL2006647 x2

3"x14" Pneumatic Black Tires No Charge
Part: WHL1713925 x2

4"x14" Knobby, Pneumatic Black Tires \$60 ea.
Part: WHL2006822 x2

6.2. Drive Wheel Accent Color

Optional. Not available with knobby tires.

Back in Black.....Standard
Part: DGN164245-PBK1 x2

Red Corvette..... No Charge
Part: DGN164245-PRD1 x2

Night Train..... No Charge
Part: DGN164245-PBK3 x2

Grey Street No Charge
Part: DGN164245-PSL3 x2

Deep Purple No Charge
Part: DGN164245-PPR5 x2

Raspberry Beret..... No Charge
Part: DGN164245-PPK4 x2

Pretty in Pink No Charge
Part: DGN164245-PPK3 x2

Evolution Orange..... No Charge
Part: DGN164245-POR7 x2

Midnight Blue No Charge
Part: DGN164245-PBU10 x2

Electric Blue No Charge
Part: DGN164245-PBU9 x2

Blue Jean No Charge
Part: DGN164245-PBU11 x2

Lemon Crush No Charge
Part: DGN164245-PYL4 x2

Green Machine No Charge
Part: DGN164245-PGN11 x2

Military Green No Charge
Part: DGN164245-PGN8 x2

Military Digital Camo \$50 ea.
Part: DGN164245-WCM7 x2

Drive wheels section continued on following page.

6. DRIVE WHEELS CONTINUED

- Faded Glory \$50 ea.
Part: DGN164245-WFL2 x2
- Carbon Fiber \$50 ea.
Part: DGN164245-WCB1 x2

7. BATTERIES

Optional. The powerchair requires 2 batteries and will be installed if ordered.

- NF-22, Gel Cell, Interceptor 2 Battery \$300 ea.
Part: BAT2004376 x2. HCPCS: E2361
- NF-22, Gel Cell, MK Battery \$400 ea.
Part: BATGEL1008 x2. HCPCS: E2361

8. PATIENT'S WEIGHT

REQUIRED. The patient's weight is required to quote and order the power chair.

- Up to 200 lb User Weight
Part: PTOINDV3837
- 201 lb-300 lb User Weight
Part: PTOINDV3838

9. ELECTRONICS

REQUIRED. Select a type of drive control from the following sections. The No Joystick option is used to omit the joystick and utilize another control, such as a specialty control at the end of this form or found on the [Alternative Controls and Electronics Order Form](#).

9.1. Expandable Controls - Quantum Intelligent Systems™

Select expandable controls or omit selection if another controller type is desired.

2. Q-Logic 3 Joystick

- Q-Logic 3 EX Joystick \$900
Part: CTL167724
- Q-Logic 3e Joystick \$620
Part: CTL1910600
Can operate full light package. Not available with alternative or specialty controls.
- Q-Logic 3 EX Joystick for Full Lights \$900
Part: CTL167725
- No Joystick, Setup for Q-Logic
Part: XRF000010



9.2. Q-Logic 3 EX Joystick Armor

REQUIRED with Q-Logic 3 EX Joystick selections above. Not available with Q-Logic 3e Joystick. Reference colors shown on page 2.

- | | | | | |
|---------------------------------------|--|--|---|--------------------------------------|
| Red Corvette
Part: DGN124336-PRD1 | Deep Purple
Part: DGN124336-PPR5 | Midnight Blue
Part: DGN124336-PBU10 | Green Machine
Part: DGN124336-PGN11 | Carbon Fiber
Part: DGN124336-WCB1 |
| Back in Black
Part: DGN124336-PBK1 | Raspberry Beret
Part: DGN124336-PPK4 | Electric Blue
Part: DGN124336-PBU9 | Military Green
Part: DGN124336-PGN8 | |
| Night Train*
See Black, Soft Touch | Pretty in Pink
Part: DGN124336-PPK3 | Blue Jean
Part: DGN124336-PBU11 | Military Digital Camo
Part: DGN124336-WCM7 | |
| Grey Street
Part: DGN124336-PSL3 | Evolution Orange
Part: DGN124336-POR7 | Lemon Crush
Part: DGN124336-PYL4 | Faded Glory
Part: DGN124336-WFL2 | |

10. JOYSTICK HANDLES

Optional. The optional handle will be mounted to the joystick and the standard handle will be shipped with the unit.

- Mushroom Handle (B) \$120
Part: ACCASMB1860.
- Chin Cup (C) \$85
Part: ACCASMB186190
- T-Bar (D) \$160
Part: ACCASMB1858.
- Medium Ball Handle (E) \$70
Part: ACCASMB1859.
- Soft Ball Handle (F) \$120
Part: ACCASMB1856.
- Extended I-Handle (G) \$140
Part: ACCASMB1857.
- Bodypoint 3" Goal Post Handle \$180
Part: PLSKNCP1057.
- Bodypoint 4" Goal Post Handle \$180
Part: PLSKNCP1058.



11. JOYSTICK MOUNTING BRACKETS

Optional. Select a mount to be used with a joystick selected in previous section.

11.1a. Fixed, In-line Joystick Bracket

Fixed, In-line, BracketStandard
Part: Various




11.1b. Swing-Away (Retractable) Joystick Bracket

Allows the joystick to retract to the side of the armrest while maintaining the same joystick orientation.

Adjustable Swing-away \$265
Part: Various.
As pictured to the right "Swing Away Joystick"

Offset Swing-away \$265
Part: Various.
Magnetic latching mechanism.
As pictured to the right "Swing Away Joystick Offset Mount (Right)"




Basic Swing-away \$265
Part: Various.
Magnetic latching mechanism.
Stealth Height Adjustable Receiver recommended.

11.1c. Specialty Joystick Mounts

Select a bracket and side to mount it.

Retract4 Mount \$450
Part: XRF000003.
Full 180° of movement that translates into multi-position use, straight, right offset, left offset, or fold away. Joystick orientation varies as bracket is adjusted.

Stealth Swing-Away/Flip Down Bracket \$300
Part: XRF000005.
Allows the joystick to pivot to side of armrest, or flip-down if bracket is adjusted.



Multi-Axis Mount \$100
Part: XRF000004
Allows height adjustment of joystick.



11.2. Joystick Mount Location

Left Side Mount
Part: PTOINDV1225

Right Side Mount
Part: PTOINDV1226

12. TRU-BALANCE® 4 SEATING

REQUIRED. Maximum seat size is 22"W x 22"D.

12.1. TRU-Balance 4 Configuration

REQUIRED.

✓ TRU-Balance 4 Power Tilt and Recline \$14,995
Part: PTOINDV3954, ADP: WP0001925
50° of posterior tilt and 165° of recline. Contains positional feedback.

12.2. 12" Power Lift

REQUIRED.

✓ Power Lift with iLevel \$2,995
Part: PTOINDV3956
12" scissor lift with positional feedback.

12.3. Anterior Tilt

Chest support is required with anterior tilt 20 and 30 degree functions. See chest support section. Omitting a chest support selection requires the provider to facilitate aftermarket chest supports that best fits consumer needs.

1. Anterior Tilt Angle

REQUIRED.

10° Anterior Tilt Angle Standard
Part: MEC2009356

30° Anterior Tilt Angle \$3,195
Part: MEC2009400

20° Anterior Tilt Angle \$2,695
Part: MEC2009399
Requires chest support and knee block.

Requires chest support and knee block.

2. Disable Anterior Tilt

Optional.

Disable Anterior Tilt..... No Charge
Part: PTOINDV3980
Only available with 10° Anterior Tilt Angle. Anterior tilt is restricted with programming.

12.4. Knee Block Assembly

REQUIRED with 20° or 30° anterior tilt. **Optional** with 10° anterior tilt.

Anterior Knee Block Assembly (ST-KBATB4ATC-1-Q) \$700
Part: ST-KBATB4ATC-1-Q

12.5. Occupied Transit Kit

Optional. The power base is standard with Unoccupied Transit Loops. 60" Sure-Lok Lap Belt is required with Occupied Transit Option.

WC-19 Occupied Transit Option..... \$400
Part: KIT152047

To determine if your chair is WC-19 compliant (occupied), please refer to this website: <http://wc-transportation-safety.umtri.umich.edu/crash-tested-product-lists/wheelchairs>

TRU-Balance 4 Seating section continued on following pages.

12. TRU-BALANCE 4 SEATING CONTINUED

12.6. Seat Width

REQUIRED with seat selection. Maximum seat width is 22". Part number is SETWDTH plus number following # below.

16" W (#1006)	17" W (#1007)	18" W (#1008)	19" W (#1009)	20" W (#1010)	21" W (#1011)	22" W (#1012)
------------------	------------------	------------------	------------------	------------------	------------------	------------------

12.7. Seat Width Package

REQUIRED with seat selection.

16"-20" Width Package
Part: XRF000008

20"-24" Width Package
Part: XRF000158

12.8. Seat Depth

REQUIRED with seat selection. Maximum seat depth is 22". Part number is SETDPHT plus number following # below.

14" D (#1004)	15" D (#1005)	16" D (#1006)	17" D (#1007)	18" D (#1008)	19" D (#1009)	20" D (#1010)	21" D (#1011)	22" D (#1012)
------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------

12.9. Seat Depth Package

REQUIRED with seat selection.

16"-22" Depth Package
Part: XRF000309

Lowest	Low	Med	High
Q6 Edge 3			
17.50"	18.25"	18.75"	19.25"
TRU-Balance 4 Seat to Floor Heights			

12.10. Seat to Floor Height

REQUIRED.

Lowest Seat to Floor
Part: PTOINDV3558

Low Seat to Floor
Part: PTOINDV3976

Middle Seat to Floor
Part: PTOINDV3559

Highest Seat to Floor
Part: PTOINDV3560


12.11. Matching Seat Cushion

Optional. Specific cushion/component sizes available, part numbers, weight capacities, MSRP, and special order cushions can be found on the respective order forms. This section allows a quick selection of a single cushion that matches the seat pan width and depth. See specific order form for non-matching sizes, alternate covers and liners, or additional covers and liners.

Stealth Simplicity G Cushion

General use cushion with Stealth Tek 2, 2-way stretch, breathable cover.
Part series: ST-SIMGXXYY

Simplicity G, Matching Seat Size
Part: XRF000111



Stealth Spectrum Foam SPP Cushion

Skin protection and positioning with Coolcore, 4-way stretch, breathable cover. Includes urethane liner.
Part series: ST-SPCFSPXXYY


Spectrum Foam SPP, Matching Seat Size
Part: XRF000112



Stealth Solution SPP Cushion

Skin protection and positioning with Coolcore, 4-way stretch, breathable cover. Includes urethane liner.
Part series: ST-SOLSPPXXYY


Solution SPP, Matching Seat Size
Part: XRF000113



Stealth Glacial SP Cushion

Skin protection with Coolcore, 4-way stretch, breathable cover.
Part series: ST-GLSPXXYY

Glacial SP, Matching Seat Size
Part: XRF000114



Stealth Zen SP Cushion

Skin protection with Coolcore, 4-way stretch, breathable cover.
Part series: ST-ZENSPXXYY


Zen SP, Matching Seat Size
Part: XRF000115



Stealth TRU-Comfort 2 SPP Cushion

Skin protection and positioning with Coolcore, 4-way stretch, breathable cover.
Part series: ST-TC2SPPXXYY

TRU-Comfort 2 SPP, Matching Seat Size
Part: XRF000116



Stealth Spectrum Gel SPP Cushion

Skin protection and positioning with Coolcore, 4-way stretch, breathable cover. Includes urethane liner.
Part series: ST-SPCGSPXXYY

Spectrum Gel SPP, Matching Seat Size
Part: XRF000117



Stealth Essence SPP Cushion

Skin protection and positioning with Coolcore, 4-way stretch, breathable cover.
Part series: ST-ESSSPXXYY

Essence SPP, No Wedge, Matching Seat Size
Part: XRF000118

Essence SPP, Small Wedge, Matching Seat Size
Part: XRF000119

Essence SPP, Medium Wedge, Matching Seat Size
Part: XRF000120

Essence SPP, Large Wedge, Matching Seat Size
Part: XRF000121

Essence Plus SPP, Large Wedge, Matching Seat Size
Part: XRF000121



Stealth Premiere P Cushion

Positioning with Coolcore, 4-way stretch, breathable cover.
Part series: ST-PREMPXXYY

Premiere P, Matching Seat Size
Part: XRF000122



12.12. Back Cushion Configuration

REQUIRED. Select a configuration type and applicable dimensions. See the Sport, Q-Back, and TRU-Comfort back availability chart for stock order back size ranges. Other size backs may be available with special order and extended lead time.

1a. Sport Back/Omit/After-Market Interface

Sport Back..... No Charge
Part: XRF000017
Complete following width package and size.

Omit Back Pan/Cushion..... No Charge
Part: PTOINDV1287

Select a following width package.

After-market Back Interface No Charge
Part: ST-SSAPT4-Q. Select a set width in following section, range of 12"-21". Overall seat depth may need to be adjusted to accommodate certain backs. May interfere with gel/waterfall armpads and arm

TRU-Balance 4 Seating section continued on following pages.

12. TRU-BALANCE 4 SEATING CONTINUED

1b. TRU-Comfort Back

TRU-Comfort Back (Static/Flip Forward), Up to 15" Width..... \$835
 Part: SET123622. HCPCS: E2620
 Complete following width package and size.

TRU-Comfort Back (Static/Flip Forward), 16" to 20" Width... \$1,035
 Part: TC4V10W21W. HCPCS: E2620
 Complete following width package and size.

TRU-Comfort Back (Static/Flip Forward), 22" Width \$1,235
 Part: TC4V22W32W. HCPCS: E2621
 Complete following width package and size.

1c. Q-Back 2.0, AL Series, Active Contour

Select a back type and the specific overall back height in subsequent section.

Active Contour, 14"Wx19-20"H Tilt/22"-23" H Recline \$985
 Part: QB2AL14AS. HCPCS: E2613

Active Contour, 14"Wx21-22"H Tilt/24"-25" H Recline \$985
 Part: QB2AL14AT. HCPCS: E2613

Active Contour, 16"Wx19-20"H Tilt/22"-23" H Recline \$1,145
 Part: QB2AL16AS. HCPCS: E2613

Active Contour, 16"Wx21-22"H Tilt/24"-25" H Recline \$1,145
 Part: QB2AL16AT. HCPCS: E2613

Active Contour, 16"Wx23-24"H Tilt/26"-27" H Recline \$1,145
 Part: QB2AL16AX. HCPCS: E2613

Active Contour, 16"Wx25-26"H Tilt/28"-29" H Recline \$1,145
 Part: QB2AL16AXX. HCPCS: E2613

Active Contour, 18"Wx19-20"H Tilt/22"-23" H Recline \$1,145
 Part: QB2AL18AS. HCPCS: E2613

Active Contour, 18"Wx21-22"H Tilt/24"-25" H Recline \$1,145
 Part: QB2AL18AT. HCPCS: E2613

Active Contour, 18"Wx23-24"H Tilt/26"-27" H Recline \$1,145
 Part: QB2AL18AX. HCPCS: E2613

Active Contour, 18"Wx25-26"H Tilt/28"-29" H Recline \$1,145
 Part: QB2AL18AXX. HCPCS: E2613

Active Contour, 20"Wx19-20"H Tilt/22"-23" H Recline \$1,145
 Part: QB2AL20AS. HCPCS: E2613

Active Contour, 20"Wx21-22"H Tilt/24"-25" H Recline \$1,145
 Part: QB2AL20AT. HCPCS: E2613

Active Contour, 20"Wx23-24"H Tilt/26"-27" H Recline \$1,145
 Part: QB2AL20AX. HCPCS: E2613

Active Contour, 20"Wx25-26"H Tilt/28"-29" H Recline \$1,145
 Part: QB2AL20AXX. HCPCS: E2613

1d. Q-Back 2.0, AL Series, Deep Contour

Select a back type and the specific overall back height in subsequent section.

Deep Contour, 14"Wx19-20"H Tilt/22"-23" H Recline \$1,145
 Part: QB2AL14DS. HCPCS: E2620

Deep Contour, 14"Wx21-22"H Tilt/24"-25" H Recline \$1,145
 Part: QB2AL14DT. HCPCS: E2620

Deep Contour, 16"Wx19-20"H Tilt/22"-23" H Recline \$1,325
 Part: QB2AL16DS. HCPCS: E2620

Deep Contour, 16"Wx21-22"H Tilt/24"-25" H Recline \$1,325
 Part: QB2AL16DT. HCPCS: E2620

Deep Contour, 16"Wx23-24"H Tilt/26"-27" H Recline \$1,325
 Part: QB2AL16DX. HCPCS: E2620

Deep Contour, 16"Wx25-26"H Tilt/28"-29" H Recline \$1,325
 Part: QB2AL16DXX. HCPCS: E2620

Deep Contour, 18"Wx19-20"H Tilt/22"-23" H Recline \$1,325
 Part: QB2AL18DS. HCPCS: E2620

Deep Contour, 18"Wx21-22"H Tilt/24"-25" H Recline \$1,325
 Part: QB2AL18DT. HCPCS: E2620

Deep Contour, 18"Wx23-24"H Tilt/26"-27" H Recline \$1,325
 Part: QB2AL18DX. HCPCS: E2620

Deep Contour, 18"Wx25-26"H Tilt/28"-29" H Recline \$1,325
 Part: QB2AL18DXX. HCPCS: E2620

Deep Contour, 20"Wx19-20"H Tilt/22"-23" H Recline \$1,325
 Part: QB2AL20DS. HCPCS: E2620

Deep Contour, 20"Wx21-22"H Tilt/24"-25" H Recline \$1,325
 Part: QB2AL20DT. HCPCS: E2620

Deep Contour, 20"Wx23-24"H Tilt/26"-27" H Recline \$1,325
 Part: QB2AL20DX. HCPCS: E2620

Deep Contour, 20"Wx25-26"H Tilt/28"-29" H Recline \$1,325
 Part: QB2AL20DXX. HCPCS: E2620

1e. Q-Back 2.0, AL Series, Extra Deep Contour

Select a back type and the specific overall back height in subsequent section.

Extra Deep Contour, 14"Wx19-20"H Tilt/22"-23" H Recline \$1,225
 Part: QB2AL14ES. HCPCS: E2620

Extra Deep Contour, 14"Wx21-22"H Tilt/24"-25" H Recline \$1,225
 Part: QB2AL14ET. HCPCS: E2620

Extra Deep Contour, 16"Wx19-20"H Tilt/22"-23" H Recline \$1,410
 Part: QB2AL16ES. HCPCS: E2620

Extra Deep Contour, 16"Wx21-22"H Tilt/24"-25" H Recline \$1,410
 Part: QB2AL16ET. HCPCS: E2620

Extra Deep Contour, 16"Wx23-24"H Tilt/26"-27" H Recline \$1,410
 Part: QB2AL16EX. HCPCS: E2620

Extra Deep Contour, 16"Wx25-26"H Tilt/28"-29" H Recline \$1,410
 Part: QB2AL16EXX. HCPCS: E2620

Extra Deep Contour, 18"Wx19-20"H Tilt/22"-23" H Recline \$1,410
 Part: QB2AL18ES. HCPCS: E2620

Extra Deep Contour, 18"Wx21-22"H Tilt/24"-25" H Recline \$1,410
 Part: QB2AL18ET. HCPCS: E2620

Extra Deep Contour, 18"Wx23-24"H Tilt/26"-27" H Recline \$1,410
 Part: QB2AL18EX. HCPCS: E2620

Extra Deep Contour, 18"Wx25-26"H Tilt/28"-29" H Recline \$1,410
 Part: QB2AL18EXX. HCPCS: E2620

Extra Deep Contour, 20"Wx19-20"H Tilt/22"-23" H Recline \$1,410
 Part: QB2AL20ES. HCPCS: E2620

Extra Deep Contour, 20"Wx21-22"H Tilt/24"-25" H Recline \$1,410
 Part: QB2AL20ET. HCPCS: E2620

Extra Deep Contour, 20"Wx23-24"H Tilt/26"-27" H Recline \$1,410
 Part: QB2AL20EX. HCPCS: E2620

Extra Deep Contour, 20"Wx25-26"H Tilt/28"-29" H Recline \$1,410
 Part: QB2AL20EXX. HCPCS: E2620

TRU-Balance 4 Seating section continued on following page.

12. TRU-BALANCE 4 SEATING CONTINUED

12.13. TRU-Balance 4 Back Shroud

REQUIRED. The back shroud is a feature of the TRU-Balance 4 seating system and is available in matching shroud colors. The back shroud is not available when Omit Back Pan/Cushion is selected. Reference colors shown on page 2.

- | | | | |
|---|--|---|--|
| Red Corvette
Part: XRF000019 | Deep Purple
Part: XRF000025 | Midnight Blue
Part: XRF000029 | Green Machine
Part: XRF000033 |
| Back in Black
Part: XRF000021 | Raspberry Beret
Part: XRF000026 | Electric Blue
Part: XRF000030 | Military Green
Part: XRF000034 |
| Night Train
Part: XRF000022 | Pretty in Pink
Part: XRF000027 | Blue Jean
Part: XRF000031 | Faded Glory
Part: DGN1714085-WFL2 |
| Grey Street
Part: XRF000023 | Evolution Orange
Part: XRF000028 | Lemon Crush
Part: XRF000032 | Digital Camo
Part: DGN1714085-WCM7 |

12.14. Lateral Positioning Wedge

Optional. Available with Sport and TRU-Comfort backs. Not available on power recline back 19"-20"H with limited adjustability on 21"H.

Lateral Positioning Wedge \$30
Part: SET123465 x2

12.15. Back Width Package

REQUIRED except when Omit Back Pan/Cushion is selected. Available with +/-4" of seat width.

- Match Seat Width Selections**
Part: PTOINDV3468
- 14"W-16"W Package**
Part: XRF000037
- 16"W-20"W Package**
Part: XRF000038
- 18"W-22"W Package**
Part: XRF000039

- S Sport Back
- T TRU-Comfort Back
- ST Sport or TRU-Comfort Back
- SQT Sport, Q-Back, or TRU-Comfort Back

		Back Width												
		12"	14"	15"	16"	17"	18"	19"	20"	21"	22"			
Overall Height	17"	T	-	-	-	-	-	-	-	-	-	-	-	-
	18"	T	-	-	-	-	-	-	-	-	-	-	-	-
	19"	S	ST	ST	ST	ST	ST	ST	ST	ST	ST	ST	ST	ST
	20"	S	ST	ST	ST	ST	ST	ST	ST	ST	ST	ST	ST	ST
	21"	-	ST	ST	SQT	ST	ST	ST	ST	SQT	ST	ST	ST	ST
	22"	-	ST	ST	SQT	ST	SQT	ST	SQT	ST	SQT	ST	ST	ST
	23"	-	ST	ST	SQT	ST	SQT	ST	SQT	ST	SQT	ST	ST	ST
	24"	-	ST	ST	SQT	ST	SQT	ST	SQT	ST	SQT	ST	ST	ST
	25"	-	ST	ST	SQT	ST	SQT	ST	SQT	ST	SQT	ST	ST	ST
	26"	-	ST	ST	ST	ST	ST	ST	ST	ST	ST	ST	ST	ST
	27"	-	ST	ST	ST	ST	ST	ST	ST	ST	ST	ST	ST	ST
	28"	-	ST	ST	ST	ST	ST	ST	ST	ST	ST	ST	ST	ST
	29"	-	ST	ST	ST	ST	ST	ST	ST	ST	ST	ST	ST	ST
30"	-	ST	ST	ST	ST	ST	ST	ST	ST	ST	ST	ST	ST	
31"	-	ST	ST	ST	ST	ST	ST	ST	ST	ST	ST	ST	ST	
32"	-	ST	ST	ST	ST	ST	ST	ST	ST	ST	ST	ST	ST	

Sport, Q-Back, and TRU-Comfort Back Availability Chart

12.16. Back Width

REQUIRED except when Omit Back Pan/Cushion is selected. Available with +/-4" of seat width.

- | | | | | |
|--|--|--|--|--|
| 12" Back Width
Part: PTOINDV3562 | 16" Back Width
Part: PTOINDV3564 | 18" Back Width
Part: PTOINDV3566 | 20" Back Width
Part: PTOINDV3568 | 22" Back Width
Part: PTOINDV3570 |
| 14" Back Width
Part: PTOINDV3563 | 17" Back Width
Part: PTOINDV3565 | 19" Back Width
Part: PTOINDV3567 | 21" Back Width
Part: PTOINDV3569 | |

12.17. Overall Back Height

Power Recline available height range: 19"-32"H. Odd heights include 6" gap and even heights include 7" gap between seat pan and bottom of back cushion. See the availability chart for stock order back size ranges.

- | | | | | |
|--|--|--|--|--|
| 19" Height
Part: SETHGHT1032 | 22" Height
Part: SETHGHT1035 | 25" Height
Part: SETHGHT1038 | 28" Height
Part: SETHGHT1041 | 31" Height
Part: SETHGHT1102 |
| 20" Height
Part: SETHGHT1033 | 23" Height
Part: SETHGHT1036 | 26" Height
Part: SETHGHT1039 | 29" Height
Part: SETHGHT1042 | 32" Height
Part: SETHGHT1103 |
| 21" Height
Part: SETHGHT1034 | 24" Height
Part: SETHGHT1037 | 27" Height
Part: SETHGHT1040 | 30" Height
Part: SETHGHT1043 | |

12.18. Optional Back Height Extension

Optional. Adds 2" of seat back pan height. 25"H selection includes a back height extension.

Optional Back Height Extension \$48
Part: XRF000041

13. ARMRESTS

REQUIRED. Select the desired type of armrest below. Please read notes to ensure seating compatibility. Be sure to also select an armrest size and armpad type.

13.1. Flip-Back, 8"-16" H Armrests

Choose a left and right armrest and set height in the following section. Not available with waterfall armpads. Not available with Back Cane Option.

L R Flip-Back, 8"-16"H Armrests \$145 ea.
Left part: ARM1806662. Right part: ARM1806663

Armrest Support

Optional. Available up to 2" difference between seat and back width via endomorph or mesomorph configuration. Not available less than 16" width seat or back.

L R Flip-Back Armrest Downpost \$115 ea.
Left part: XRF000265. Right part: XRF000266

Armrests section continued on following page.

13. ARMRESTS CONTINUED

13.2. Armrest Height

REQUIRED. Set left and right armrest height. Note the available height range in the option description.

1. Left Armrest Height

8" H PTOINDV3535	10" H PTOINDV3537	12" H PTOINDV3539	14" H PTOINDV3541	16" H PTOINDV3543
9" H PTOINDV3536	11" H PTOINDV3538	13" H PTOINDV3540	15" H PTOINDV3542	

2. Right Armrest Height

8" H PTOINDV3450	10" H PTOINDV3452	12" H PTOINDV3454	14" H PTOINDV3456	16" H PTOINDV3458
9" H PTOINDV3451	11" H PTOINDV3453	13" H PTOINDV3455	15" H PTOINDV3457	

13.3. Armrest Size

REQUIRED. Select a left (L) and right (R) armrest size. This also affects armpad selection.

<input type="checkbox"/> L <input type="checkbox"/> R	10" Length L part: XRF000057. R part: XRF000123	<input type="checkbox"/> L <input type="checkbox"/> R	14" Length L part: XRF000058. R part: XRF000124	<input type="checkbox"/> L <input type="checkbox"/> R	17" Length L part: XRF000059. R part: XRF000125
---	---	---	---	---	---

14. ARMPADS

REQUIRED.

14.1a. Aftermarket Armpad

Not available with Chest Bar. Select one below if aftermarket armpads or handpads are required.

Interface for Aftermarket Pad No Charge
Part: XRF000308

14.1b. Foam Armpads

L R Straight, Armpad Standard
Left part: XRF000061. Right part: XRF000126

14.1c. Gel Armpads

Sold only in pairs. Armpads are 3" wide. Manufacturer part number for full length is 314G and desk length is 310G.

Falcon Gel, Full Length Left and Right \$300 Part: ARM1806807	Falcon Gel, Desk Length Left, Full Length Right..... \$250 Part: ARM1806810
Falcon Gel, Full Length Left, Desk Length Right..... \$250 Part: ARM1806813	Falcon Gel, Desk Length Left and Right \$200 Part: ARM1806806

15. JOYSTICK BRACKETS

REQUIRED.

Standard Fixed Receiver Standard Part: XRF000063	Stealth Height Adjustable Receiver.....\$165 Part: XRF000064
---	---

16. HEADRESTS

Optional. Select a headrest pad and bracket.

16.1a. Headrests, Comfort Plus with Cool Core

1. Headrest Pads

Stealth 8" Comfort Plus, Cool Core (CP08CCFPM) \$190 Part: ST-CP08CCFPM-Q. E0955	Stealth 14" Comfort Plus, Cool Core (CP14CCFPM) \$190 Part: ST-CP14CCFPM-Q. E0955
Stealth 10" Comfort Plus, Cool Core (CP10CCFPM)..... \$190 Part: ST-CP10CCFPM-Q. E0955	Stealth 6", Comfort Plus, Cool Core (CP06CCFPM) \$190 Part: ST-CP06CCFPM-Q. E0955

2. Headrest Brackets

Comfort Plus Headrest Mounting Bracket, UniLink.....\$255
Part: ST-UHB2D33200-Q. E1028

16.1b. Headrests, Comfort Plus

1. Headrest Pads

Stealth 8" Comfort Plus (CP450)..... \$190 Part: POS148484. HCPCS: E0955	Stealth 14" Comfort Plus (CP150) \$190 Part: INDPART2327. HCPCS: E0955
Stealth 10" Comfort Plus (CP250)..... \$190 Part: SETHEAD1057. HCPCS: E0955	Pediatric Comfort Plus (6"Wx4.5"H) (CPN350) \$190 Part: SETHEAD1069. HCPCS: E0955

2. Headrest Brackets

Multi-Axis Removable Headrest Bracket (HMO475P) \$200
Part: HDWASMB1084. HCPCS: E1028

17. LAP BELTS

REQUIRED. Select a lap belt below.

17.1a. Standard Lap Belts

- 38" Lap Belt Part: ACCBELT1006 50" Lap Belt Part: ACCBELT1000 60" Lap Belt Part: ACCBELT1002 70" Lap Belt Part: ACCBELT1003 80" Lap Belt Part: ACCBELT1004 90" Lap Belt Part: ACCBELT1001

17.1b. Optional Lap Belts

- Stealth Positioning, 60" Single Pull, Push Button, Medium..... \$65
(7.5") Padded, 1.5" Belt (IFPB5213)
Part: ST-IFPB5213-Q
- Stealth Positioning, 60" Single Pull, Push Button, Extra Large, Padded (11.5"), 2" Belt (IFPB5315)
Part: ST-IFPB5315-Q
- Stealth Positioning, 60" Single Pull, Push Button, Large \$65
Padded (9.5"), 2" Belt (IFPB5314)
Part: ST-IFPB5314-Q
- 60" Sure-Lok Lap Belt (Non-retractable)..... \$200
Part: ACCBELT1016
Required with ISO 7176-19 Compliant Transit Option.

18. LEGRESTS

All efforts will be made to accommodate legrest length selection, but minor changes may need to be made to meet shipping guidelines. Some adjustments may be needed upon receipt of the unit.

18.1. Center Mount Foot Platform

REQUIRED.

1. Foot Platform Assembly

- ✓ Dual Actuator, Power Articulating Foot Platform..... \$4,295
Part: RIG2009246. HCPCS: E1012+K0040, ADP: WPO001926
Positional feedback with independent control of articulation.

2. Lower Extension and Angle Adjustable Footplates

Medium length extension range is 11"-17" and long extension range is 14"-20". See compatibility matrix for possible interference.

- Medium Extension, 5"x6" Footplates Part: RIG150206 Long Extension, 5"x6" Footplates Part: RIG151505
- Medium Extension, 5"x8" Footplates Part: RIG150203 Long Extension, 5"x8" Footplates Part: RIG151506
- Medium Extension, 5"x11" Footplates Part: RIG150212 Long Extension, 5"x11" Footplates Part: RIG151508

Legrest Compatibility Matrix Articulating Foot Platform TRU-Balance 4	Seat to Floor Option			
	Lowest	Low	Middle	High
Medium Length Extension	11"-15"	11"-16"	11"-16.5"	11"-17"
Long Length Extension	14"-15"	14"-16"	14"-16.5"	14"-17"

3a. Center Mount Calf Pads

Calf pads are standard with a power articulating foot platform.

- Small Calf Pads, AFP Only, 6"Wx6"H (pair) Standard Part: RIGASMB7120018
- Large Calf Pads, AFP Only, 8" taper to 6"Wx8"H (pair) Standard Part: RIGASMB7120019

4. Center Mount Modular Footplate Extension Components

Optional. Assembly required. May require seating configuration adjustments to avoid front caster interference.

- Width Extension for 8"D Footplates \$55
Part: RIG152517
- Width Extension for 11"D Footplates \$55
Part: RIG152518
- 2" Bell-Shaped Length Extension \$55
Part: RIG152519
Not compatible with width extensions. Can be used with 6" outside edge with 11" footplates.
- 4" Outside Edge (pair) \$55
Part: RIG152510
Only available with 6"D footplates and 8"D footplate extension.
- 6" Outside Edge (pair) \$55
Part: RIG152516
Available with 8"D and 11"D footplates and 11"D footplate extension.

5. Center Mount Additional Accessories

- Small Reversed Dartex Calf Pad Cover (pair) \$150
Part: RIG133717
- Large Reversed Dartex Calf Pad Cover (pair) \$150
Part: RIG133719
- Footplate Connector \$55
Part: RIGASMB7120023

18.2. Legrest Accessories

- Wheel Bumper \$15 ea.
Part: XRF000100

18.3. Legrest Length

REQUIRED with matching left and right legrest configuration. Part number is PTOINDV plus number following # below. All efforts will be made to accommodate legrest length selection, but minor changes may need to be made to meet shipping guidelines. Some adjustments may be needed upon receipt of the unit. Reference the section's compatibility matrix for legrest interference. Non-matching or individual legrest configurations lengths are not set by the manufacturer.

- | | | | | | | |
|-----------------|-----------------|-----------------|-----------------|-----------------|-----------------|---------------|
| 11" L (#3478) | 12" L (#3480) | 13" L (#3482) | 14" L (#3484) | 15" L (#3486) | 16" L (#3488) | 17" L (#3490) |
| 11.5" L (#3479) | 12.5" L (#3481) | 13.5" L (#3483) | 14.5" L (#3485) | 15.5" L (#3487) | 16.5" L (#3489) | |

19. POWER POSITIONING ELECTRONICS

REQUIRED.

19.1a. Q-Logic 3 Controls through Joystick

✓ Multiple Actuator Through Controller and Interface..... \$2,175
 Part: XRF000336. HCPCS: E2311, ADP: SECX15005

19.1b. Q-Logic 3 Controls with iAccess Switch Control

iAccess Flat Button.....No Charge Part: CTL144009 See next section for switch mounting.	iAccess Toggle Button.....No Charge Part: CTL144010 See next section for switch mounting.
---	---

19.2. iAccess Mount

REQUIRED with iAccess options.

1. iAccess Mount Selection

Joystick Mount No Charge Part: MEC1907500	Attendant Mount \$35 Part: MEC1906564
Under Arm Mount, Flip-Back Armrest..... No Charge Part: MEC1906531	Under Arm Mount, Flip-Back Armrest, plus iAccess Guard \$35 Part: MEC2001392

2. Under Arm Mount Side

Left Side Under Arm Mount Part: PTOINDV3678	Right Side Under Arm Mount Part: PTOINDV3679
--	---

20. ATTENDANT CONTROL

Optional. Mounting bracket is included.

Attendant Control No Charge
Part: XRF000078. HCPCS: E2331

21. SEAT MOUNTED LED LIGHTS


Optional. Requires Q-Logic 3 Joystick for Full Light Package or Q-Logic 3e Joystick.

Seat Mounted, LED, Full Light Package\$770
 Part: XRF000130
 Not available with Base Mounted, LED, Full Light Package.

22. ACCESSORIES

Optional. Multiple accessories can be mounted along the various tracks but some interference may occur.

Accessory Bar \$95 Part: ACC1909634 Horizontal rail attached to rear of seat. Standard with the Trilogy vent tray and oxygen holder.	Personal Item Hook..... \$20 Part: ACC125008
Backpack Holder (Pair) (A)..... \$42 Part: ACC147774	Privacy Bag..... \$80 Part: ACC1907570
Cane/Crutch Holder \$100 Part: ACC152515. HCPCS: E2207	Push Handles (pair)..... \$140 Part: ACC2005027 Not available with oxygen holder or vent tray.
Clothing Guard (Pair)..... \$275 Part: ACC125007	Q-Logic EX/NE Programmer \$625 Part: ELEASMB6976
Compact Vent Tray\$1,395 Part: ACC133889. HCPCS: E1029	Quantum Backpack \$200 Part: ACC1905853
Cup Holder \$20 Part: ACC125003	Rear View Mirror \$35 Part: ACC155833
Flag Kit \$125 Part: ACC1700443	Reflector Kit \$35 Part: ACC124670
Flip-Back Armrest, Side Track Accessory Bracket .. \$25 Part: MEC1809178	Stealth Large Tablet Holder \$565 Part: ACC141818
Front Mounted, Cell Phone/ Media Holder, \$95 X-Grip (E) Part: ACC155835	Transfer Handle (D)..... \$152.50 ea. Part: MEC1809930 x2. HCPCS: E0705
Glove Box (C)..... \$85 Part: ACC125011	Trilogy Vent Tray \$420 Part: ACC136791. HCPCS: E1029 Includes accessory bar, a credit will be applied PTOIN-DV3769. Contact custom for models other than 100 and 200.
Large, Black, Weather Cover \$240 Part: ACCCOVR1021	VOCSN Vent Tray for Accessory Bar Mount..... \$1,000 Part: ACC2000166. HCPCS: E1029 Includes accessory bar, a credit will be applied PTOIN-DV3769.
Oxygen Holder..... \$110 Part: ACC125006. HCPCS: E2208 Includes accessory bar, a credit will be applied PTOIN-DV3769.	



24. POSITIONING COMPONENTS

Optional. More options are available on the [Positioning Components Order Form](#).

24.1. Stealth Unilink Pelvic/Thigh Guide

1. UniLink Pelvic/Thigh Hardware

UniLink, 5" Height Adjustable, \$320 ea.
 Removable Mounting Hardware
 Part: ST-UTTP5-Q
 QTY

UniLink, 7" Height Adjustable, \$320 ea.
 Removable Mounting Hardware
 Part: ST-UTTP7-Q
 QTY



2. UniLink Pelvic/Thigh Pads

UniLink, 4" x 4" Pad \$125 ea.
 Part: ST-UTP544056100-Q
 QTY

UniLink, 4" x 12" Pad \$165 ea.
 Part: ST-UTP54D056100-Q
 QTY

UniLink, 4" x 8" Pad..... \$135 ea.
 Part: ST-UTP548056100-Q
 QTY

3. UniLink Pelvic/Thigh Guide Installation

Install on Left
 Part: PTOINDV3585

Install on Right
 Part: PTOINDV3586

Do Not Install
 Part: PTOINDV3719

24.2. Stealth Thoracic Laterals

Swing-away laterals may interfere with flip back armrests. These laterals are not available with the back cane option.

1. Bracket

Swing-away, UniLink, Redesign Back, pair (ULS4042000-Q).... \$590
 Part: ST-ULS4042000-Q. E1028
 Only available with Redesigned back assemblies.

Swing-away, Rail Mounts, pair (TWBTLTC) \$590
 Part: INDPART2677. E1028
 Not available with Redesign back assemblies.

2. Pads

Curved Pads, 3"x5", pair (TWBL3X5)..... \$140
 Part: INDPART2663. E0956

Flat Pads, 3"x5", pair (TWBLP3X5) \$140
 Part: INDPART2676. E0956

Curved Pads, 4"x6", pair (TWBL4X6)..... \$140
 Part: INDPART2664. E0956

Flat Pads, 4"x6", pair (TWBLP4X6) \$140
 Part: INDPART2675. E0956

3. Installation

Install Laterals
 Part: PTOINDV3499

Do Not Install
 Part: PTOINDV3500

24.3. Stealth Pelvic/Thigh Guides

Pelvic/thigh guides may interfere with the joystick receiver block and down posts. Please provide quantity next to check box.

1. Pelvic/Thigh Guide Bracket

Flip Down Bracket (TWBTLADD-HW) \$300 ea.
 Part: INDPART2804. : E1028
 QTY

Flip Down with LEFT Full Surface Contact \$400 ea.
 Feature (TWBADD-L & LPHW407)
 Part: POS140428.

Removable Bracket (LP430)..... \$330 ea.
 Part: INDPART2805. : E1028
 QTY

Flip Down with RIGHT Full Surface Contact..... \$400 ea.
 Feature (TWBADD-R & LPHW407)
 Part: POS140429.

Removable Bracket with 1" Offset (LP430 & LHW-121)... \$430 ea.
 Part: INDPART2854. : E1028
 QTY

2. Pelvic/Thigh Guide Adapter

2" Height Adapter Plate (LHW-122).....\$40
 Part: INDPART2821.
 Height adapter is needed with TRU-Comfort seat cushions. Not available with the flip down bracket.

2 1/4" Adapter Plate (LHW-129) \$40
 Part: POS140430.
 Only available with the flip down full surface bracket.

3. Pelvic/Thigh Guide Pad

Short Pad, 4"x4" (TWBADD-4X4)..... \$70 ea.
 Part: INDPART2801.
 QTY

Long Pad, 4"x12" (TWBADD-4X12) \$70 ea.
 Part: INDPART2803.
 QTY

Medium Pad, 4"x8" (TWBADD-4X8)..... \$70 ea.
 Part: INDPART2802.
 QTY

4. Pelvic/Thigh Guide Installation

Install on Left
 Part: PTOINDV3585

Install on Right
 Part: PTOINDV3586

Do Not Install
 Part: PTOINDV3719

Positioning components section continued on following page.

24. POSITIONING COMPONENTS CONTINUED

24.4. Chest Bar Assembly

Small Chest Bar, 14" Foam Pad (ST-CBATB4ATT-S-Q)..... \$300
 Part: ST-IFPS1103-Q

Medium Chest Bar, 16" Foam Pad (ST-CBATB4ATT-M-Q)..... \$300
 Part: ST-CBATB4ATT-M-Q

24.5. Stealth Positioning Chest Harness

Upper Classic, Medium (IFPS1103) \$190 ea.
 Part: ST-IFPS1103-Q

Upper Contour, Medium (IFPS2103)..... \$190 ea.
 Part: ST-IFPS2103-Q

Upper Classic, Large (IFPS1104) \$190 ea.
 Part: ST-IFPS1104-Q

Upper Contour, Large (IFPS2104) \$190 ea.
 Part: ST-IFPS2104-Q



24.6. Stealth Positioning Chest Straps

Chest Strap, Medium w/ Belt Anchor \$120 ea.
 (IFCS4113)
 Part: ST-IFCS4113-Q

Chest Strap, Large w/ Belt Anchor \$120 ea.
 (IFCS4114)
 Part: ST-IFCS4114-Q



24.7. Stealth Swing-Away Headrest Mounts

These swing-away headrest mounts are only available with the Stealth Comfort Series headrests.

1. Facial Lateral Hardware

Facial Lateral for Comfort Series..... \$430 ea.
 Right Side (CPS9-R-HW)
 Part: INDPART2647.

Facial Laterals for Comfort Series \$430 ea.
 Left and Right Side (CPS9-B-HW)
 Part: INDPART2649.

Facial Lateral for Comfort Series..... \$430 ea.
 Left Side (CPS9-L-HW)
 Part: INDPART2648.



2. Facial Lateral Pads and Accessories

1" Ball Kit for Facial Laterals (TWB-1BK) (D) \$100
 Part: POS150532
 Required with facial laterals to mount to comfort headrest brackets.


Facial Pad, 4"x2" (900-Q)..... \$70 ea.
 Part: ST-900-Q : K0108

Egg Switch Mount (ESM) (A) \$70 ea.
 Part: INDPART2815

Facial Pad, 6"x2" (920-Q)..... \$70 ea.
 Part: ST-920-Q : K0108

Figure 8 Bracket with 13" Rod..... \$315 ea.
 and Egg Switch Mount (B)
 Part: ST-F8675-13-Q

Swing-Away Joystick Mount, 18" Rod (C) \$295 ea.
 Part: ST-SUS18-Q
 Allows mounting of facial pads and attaches to the headrest



24.8. Stealth Elbow Support


1. Elbow Support Bracket

Left, Swing-Away, Elbow Support \$230 ea.
 Bracket Tube Style Arms (ARMS-ES-L)
 Part: POS124064.

Right, Swing-Away, Elbow Support..... \$230 ea.
 Bracket Tube Style Arms (ARMS-ES-R)
 Part: POS124066.

Left, UniLink, Fixed, Elbow Support Bracket ..\$360 ea.
 for Flip-Back, 8"-16"H Armrests (UETB4021FL-Q)
 Part: ST-UETB4021FL-Q

Right, UniLink, Fixed, Elbow Support Bracket . \$360 ea.
 for Flip-Back, 8"-16"H Armrest (UETB4021FR-Q)
 Part: ST-UETB4021FR-Q



2. Elbow Support Pad

3"Wx3"H Elbow Pad (ARMS-ES-3X3) \$70 ea.
 Part: POS124065. HCPCS: K0108

4"Wx4"H Elbow Pad (ARMS-ES-4X4)..... \$70 ea.
 Part: POS124067. HCPCS: K0108

24.9. Stealth Foot Boxes

Stealth Small Foot Boxes, 5"Wx8"Dx7"H..... \$625
 (SPLE107) (pair)
 Part: INDPART2786.

Stealth Large Foot Boxes, 7"Wx12"Dx11"H..... \$625
 (SPLE109) (pair)
 Part: INDPART2788.

Stealth Medium Foot Boxes, 6"Wx10"Dx9"H \$625
 (SPLE108) (pair)
 Part: INDPART2787.

24.10. Stealth Padded Shoe Holder

Stealth Small Padded Shoe Holder 3.75"Wx6.25"D \$220
 (SPLE111) (pair)
 Part: INDPART2783.

Stealth X-Large Padded Shoe Holder, 5"Wx10.75"D \$220
 (SPLE114) (pair)
 Part: INDPART2785.

Stealth Medium Padded Shoe Holder, 4"Wx8.25"D \$220
 (SPLE112) (pair)
 Part: INDPART2784.

25. ALTERNATIVE CONTROLS AND ELECTRONICS

Optional. Select a type of specialty controls below. More options are available on the [Alternative Controls and Electronics order form](#).

25.1. Enhanced Display

Q-Logic 3 Enhanced Display Kit..... \$2,045
 Part: CTL167726
 Select a color option below. Includes enhanced display, flexible gooseneck and mounting bracket. Features: 3.5" LCD display; up to 5 profiles to control drive, seat, and auxiliary functions; Bluetooth connectivity for mobile devices and mouse emulation; integrated single switch scanner; IR controls.



Enhanced Display

Q-Logic 3 Enhanced Display Shroud

Required with Q-Logic 3 Enhanced Display selection above. Reference colors shown on page 2.

- | | | |
|--|---|--|
| Black, Soft Touch
Part: DGN137728 | Raspberry Beret
Part: DGN137729-PPK4 | Lemon Crush
Part: DGN137729-PYL4 |
| Red Corvette
Part: DGN137729-PRD1 | Pretty in Pink
Part: DGN137729-PPK3 | Green Machine
Part: DGN137729-PGN11 |
| Back in Black
Part: DGN137729-PBK1 | Evolution Orange
Part: DGN137729-POR7 | Military Green
Part: DGN137729-PGN8 |
| Night Train
Part: DGN137729-PBK3 | Midnight Blue
Part: DGN137729-PBU10 | Military Digital Camo
Part: DGN137729-WCM7 |
| Grey Street
Part: DGN137729-PSL3 | Electric Blue
Part: DGN137729-PBU9 | Faded Glory
Part: DGN137729-WFL2 |
| Deep Purple
Part: DGN137729-PPR5 | Blue Jean
Part: DGN137729-PBU11 | Carbon Fiber
Part: DGN137729-WCB1 |

25.2. SCIM Module

Optional.
SCIM Module for Q-Logic 3 \$850
 Part: CTL2001963
 Allows use of joystick display for specialty devices. Available with Q-Logic 3 EX Joystick selection above. Up to 5 profiles to control drive, seat, and auxiliary functions, excluding sip-n-puff and single switch functions. A switch is required to operate the input device attached to the SCIM module.



SCIM Module

25.3. Q-Logic 3 SAJ

Q-Logic 3 EX SAJ \$850
 Part: CTL143980
 Must be used in conjunction with Enhanced Display or Q-Logic 3 Joystick.



EX SAJ

25.4. Head Array

- | | |
|---|---|
| Adult, i-Drive Tri-Array Series, 3-Switch \$5,855
Head Array (IDHBT200-1)
Part: ST-IDHBT200-1-Q. HCPCS: E2330
Features bluetooth, head support, 3 non-proportional proximity sensors for drive control, an egg switch for reset mode change, and i-Drive smart box. | Adult, i-Drive Pro Series Head Array \$5,855
(IDHBT300-1)
Part: ST-IDHBT300-1-Q. HCPCS: E2330
Features bluetooth, head support with occipital pad and 2 swing-away spot pads for left and right side, 3 non-proportional proximity sensors, an egg switch for reset mode change, and i-Drive smart box. |
|---|---|



Adult Pro-Series

25.5. Sip and Puff Drive

1. Sip and Puff Module

Sip and Puff Module for Q-Logic 3 \$1,695
 Part: CTL168056. HCPCS: E2325
 Includes module and mounting hardware.



Sip & Puff Module

2. Sip and Puff Straw and Mount

The straw and mount do not include electrical components to operate the Sip&Puff with a head array. This component requires the selection of a Q-Logic Sip&Puff module or additional module from Stealth Products for operation.

- | | |
|---|--|
| Stealth Sip and Puff for Comfort Series, Left Mount (SP625-CL) \$410
Part: ST-SP625-CL-Q. HCPCS: E2326
Includes straw and mounting hardware. Use with Stealth Comfort Plus headrests. | Stealth Sip and Puff for Stealth i-Drive, Head Array (SP625-U) \$410
Part: ST-SP625-U-Q. HCPCS: E2326
Includes straw and mounting hardware. Use with Stealth i-Drive Head Arrays. |
| Stealth Sip and Puff for Comfort Series, Right Mount (SP622-CR) \$410
Part: ST-SP625-CR-Q. HCPCS: E2326
Includes straw and mounting hardware. Use with Stealth Comfort Plus headrests. | |



Stealth Sip & Puff

Alternative controls and electronics section continued on following page.

25. ALTERNATIVE CONTROLS AND ELECTRONICS CONTINUED

25.6. Egg Switches

Red Egg Switch (ESRED)..... \$100 <i>Part: SWTMCRO1077</i>	Green Egg Switch (ESGRN) \$100 <i>Part: SWTMCRO1078</i>
Blue Egg Switch (ESBLU) \$100 <i>Part: SWTMCRO1079</i>	Yellow Egg Switch (ESYEL) \$100 <i>Part: SWTMCRO1080</i>
Black Egg Switch (ESBLK)..... \$100 <i>Part: SWTMCRO1081</i>	



25.7. Switch Adapter

2 to 1 Switch Adapter..... \$135
Part: DWR1361H037
Connects 2 single switches to a stereo jack to act as a dual switch.

25.8. Environmental Controls

Q-Logic 3 ECU Module Kit..... \$1,780
Part: CTL168057

25.9. USB Mobile Device Charger

USB Mobile Device Charger \$185
Part: ELE168069
Can charge 1 device. A multiplier harness is required and must be added when a Static Seat or Power Tilt Only seating system is chosen. The selection of more than 3 modules requires an additional multiplier.



26. TRU-BALANCE 3 MULTIPLIER HARNESS

REQUIRED. Multiple controllers or power accessories may require an additional multiplier harness or special order. Reference the chart below to determine the appropriate selection. The rows refer to the total amount of arm mounted devices and the columns refer to the total amount of back mounted devices. The online order configurator will automatically make this selection under specialty controls.

Secondary Multiplier (SM)..... \$250
Part: XRF000107

Arm Mounts	Joystick (Hand Control)	1 Arm	2 Arm	3 Arm	4+ Arm
	Enhanced Display	NA	SM	SM	SOF
	iAccess	NA	SM	SM	SOF
	USB Charger	SM	SM	SM	SOF
	SAJ	SOF	SOF	SOF	SOF

Power Positioning Chart

Back Mounts	Attendant Control	1 Arm	2 Arm	3 Arm	4+ Arm
	SCIM	NA	SM	SOF	SOF
	Sip n Puff	NA	SM	SOF	SOF
	ECU	SM	SM	SOF	SOF
		SOF	SOF	SOF	SOF

Static Seat Chart

NA No additional Multiplier
SM Secondary Multiplier
SOF Special Order Required

Q-Logic 3 Multiplier Charts