

ETHOS

Ki* Mobility





THE RIDE DOESN'T HAVE TO BE.

Navigate city streets, cruise the backroads or tackle cobblestones. No matter where you go, let Ethos take you there. Ethos delivers **unrivaled ride quality and performance** packaged in a strikingly sculpted design that is unlike anything you have ever experienced.

Ki Mobility's patent pending **Isolation Technology** separates you from the road, dramatically reducing the amount of vibration you feel. Vibration that published research has linked to pain and fatigue.

ETHOS A NEW ERA IN WHEELCHAIR DESIGN

ISO TECH

Iso Tech (Isolation Technology) is a patent pending method of isolating the seat frame from the base frame using vibration damping polymers. The result is a ride so smooth cracks and bumps disappear in the wake of every push.

SEAT FRAME

For the first time ever the seat frame is completely isolated from the base frame. The casters are no longer attached to the frame you sit on. As a result, vibrations from the casters are no longer transferred to the footplate and your body.

VIBRATION DAMPING POLYMERS

Interposed between the seat frame and base frame are four polymers. These act in unison to significantly reduce vibration that has been linked to pain and fatigue in studies. Ethos offers four different polymers allowing you to create your ultimate ride experience.

BASE FRAME

The base frame has been engineered with incredible rigidity for optimum performance. Its front caster design provides just the right amount of flex to smoothly navigate bumps and cracks. The perfect balance between performance and ride quality.



SEAT FRAME

BASE FRAME

The graph above shows actual vibration damping from using an instrumented Ethos on real world surfaces.



UNRIVALED PERFORMANCE

Easily optimize Ethos to create a perfect balance of agility, efficiency and responsiveness. Enhance the stability and quality of your ride experience with the first ever independently adjustable wheelbase. The ultra-rigid base frame has just the right amount of flexion in the front to keep the casters on track when navigating all environments.



MAKE IT YOURS

Whether your style preference is minimalistic or edgy you can personalize your Ethos with a combination of 21 frame colors, 8 highlights and 8 anodized options. Make it as bright and bold as your personality or make a power statement with a more subtle and sleek design. No other chair allows you to express yourself like Ethos.



DON'T JUST TAKE OUR WORD FOR IT



Karah S.

Working professional and wheelchair user since 2009

"In my old wheelchair my back and hips were so sore, and my feet were always bouncing off my footplate. Ethos has reduced my spasticity and improved how I feel throughout the day."

Micah Fowler

Actor on ABC's "Speechless" and wheelchair user since 2003

"Ethos is a game changer. It's so much easier for me to roll this wheelchair. When I use it outside it is so smooth and nice! It is the best wheelchair I've ever owned."

Rick H.

Disability advocate and wheelchair user since 1976

"It rides smoother than any other chair I have been in before. Ethos significantly reduces the vibration I feel from the road. I come back from rides feeling less fatigued."

SPECIFICATIONS

Transport Weight*	Seat Width	Seat Depth	Front Seat Height	Rear Seat Height
11.6 lbs. 5.26 kg	12-20" 30-51 cm	14-20" 36-51 cm	13.5-20.5" 34-52 cm	13.5-20" 34-51 cm
Camber	Tubing Size	Weight Capacity	Ansi Resna WC-19 Transit Approved	
0°, 2°, 4°, 6°, 8°	1 1/8" 2.86 cm	275 lbs. 124.7 kg	Yes	

* 16"x16" Standard Length Frame, Fixed Height Back, 4" Poly Casters, Tubular Open Footplate and Standard Upholstery. Rear Wheels Removed.

