



# Golden Scooters & Power Wheelchairs

SUPERIOR DESIGNS | FIRST CLASS LUXURY



**BUZZ**<sup>™</sup>  
AROUND  
*Carry On*

*Carry On with Life*

 **GOLDEN**



## **Relax & Ride in Style!**

**Our scooters and power wheelchairs offer excellent maneuverability and performance to get you where you want to go in comfort and style. Explore our line of mobility products and choose the one that's perfect for you! Plus, all models are backed by the best warranties in the industry!**

Highly Dependable • Maximum Maneuverability  
Unparalleled Stability • Greatest Functionality  
Superior Designs • First Class Luxury

*Golden scooters and power wheelchairs  
are Class II medical devices.*



## **Luxury Meets Rugged in TWO NEW Powerful Scooters!**

**INTRODUCING America's ONLY Full-time Luxury Travel Scooter! The GB119 3-wheel and GB149 4-wheel Buzzaround LX models offer comfort and convenience in a stylish new design.**

Enjoy upgraded features in a whole new class of scooters  
in our Buzzaround Series!



Full-size comfort comes with the convenience of true portability on the Buzzaround LX, along with standard features that include a stylish captain's seat, an infinitely adjustable tiller at your fingertips, and a USB charging port right under the dashboard. With a 375 lb. weight capacity, comfort suspension springs, and up to 18 mile range on two U1 batteries, the Buzzaround LX offers ample foot room to take you where you want to go in comfort and style!

**BEST IN CLASS  
RANGE!**

**BEST IN CLASS  
WEIGHT CAPACITY!**



GB119

**WARRANTIES:**

**FRAME:** Lifetime Warranty

**DRIVE TRAIN:** 2 Year Warranty

**ELECTRONICS:** 2 Year Warranty

**BATTERIES:** 1 Year Warranty



**BUZZAROUND LX  
FEATURES:**

- Stylish Captain's Seat
- 375 lb. Weight Capacity
- Infinitely Adjustable Tiller with EZ Reach Adjustment Lever
- Up to 18 Mile Range with two U1 Batteries
- Ample Foot Room
- USB Charging Port
- Comfort Spring Suspension for a Smoother Ride
- Red & Blue Color Shrouds Standard



GB149

# BUZZ<sup>TM</sup> AROUND EX

EXTREME

## FULL-SIZE PORTABLE

The all new Golden BuzzAround Extreme Scooter features the portability and easy disassembly of a lightweight scooter with the amazing benefits of a full size scooter. Both the three and four-wheel BuzzAround Extreme Scooters come standard with 2-U1 batteries - industry leading and the largest battery in this category - for maximum operating range up to 18 miles. The Extreme also features a brand new ultrabright LED headlight and taillight, and comfort suspension springs for the smoothest ride.

Go the extra mile on a BuzzAround Extreme Scooter!

### BUZZAROUND EX FEATURES:

- 350 lb. weight capacity
- 4" ground clearance
- Full size & portable, yet easy to disassemble
- Larger 90 amp controller
- Larger 420 watt motor
- Larger 9" tires
- 2-U1 batteries for maximum range
- Ultrabright LED headlight and taillight
- Comfort suspension springs
- Operating range up to 18 miles per charge

Optional color  
shrouds available.  
\*Upcharge applies.



SCAN FOR VIDEO



**BEST IN CLASS  
RANGE!**





GB118



GB148



### OPTIONAL BUMPER CONVERSION KIT

Can be ordered as a KIT or pre-installed on models GB118, GB148, GB119 & GB149.

### WARRANTIES:

**FRAME:** Lifetime Warranty

**DRIVE TRAIN:** 2 Year Warranty

**ELECTRONICS:** 2 Year Warranty

**BATTERIES:** 1 Year Warranty



### OPTIONAL BATTERY DOCKING STATION

The Battery Docking Station for models GB118, GB148, GB119 and GB149 is portable. It connects and charges the U1 batteries. Perfect for charging batteries off the scooter, especially when visiting family or friends, or at a vacation destination. Enjoy being on the go!



# BUZZAROUND XL™



GB117D

Optional color shrouds available.  
\*Upcharge applies.



Optional color shrouds available.  
\*Upcharge applies.



SCAN FOR VIDEO



GB147D

## BUZZAROUND XL FEATURES:

- 300 lb. weight capacity
- Ultra Bright LED headlight standard
- Charging port below the control panel
- 3 inch ground clearance
- Red and Blue Shroud panels included free
- Disassembles quickly and easily for transport

## WARRANTIES:

- FRAME:** Lifetime Warranty
- DRIVE TRAIN:** 2 Year Warranty
- ELECTRONICS:** 2 Year Warranty
- BATTERIES:** 1 Year Warranty



# BUZZAROUND XL<sup>TM</sup> AROUND HD

Optional color  
shrouds available.  
\*Upcharge applies.



GB117H

Optional color  
shrouds available.  
\*Upcharge applies.



GB147H

## BUZZAROUND XLHD FEATURES:

- Increased weight capacities up to 325 lbs.
- Ultra Bright LED headlight standard
- Charging port below the control panel
- 3 inch ground clearance
- Red & Blue Shroud panels included free
- Disassembles quickly and easily for transport
- 22Amp batteries standard
- Operating range up to 14 miles per charge

## WARRANTIES:

**FRAME:** Lifetime Warranty  
**DRIVE TRAIN:** 2 Year Warranty  
**ELECTRONICS:** 2 Year Warranty  
**BATTERIES:** 1 Year Warranty



# BUZZAROUND XLS<sup>TM</sup> HD



Optional color shrouds available.  
\*Upcharge applies.



GB117Z



Optional color shrouds available.  
\*Upcharge applies.



GB147Z

## BUZZAROUND XLSHD FEATURES:

- Increased weight capacity up to 325 lbs.
- Ultra Bright LED headlight standard
- Charging port below the control panel
- 3 inch ground clearance
- Red and Blue Shroud panels included free
- Disassembles quickly and easily for transport
- 22 Amp batteries standard
- Operating range up to 14 miles per charge
- NEW Comfort Spring Suspension for a smoother ride

## **WARRANTIES:**

- FRAME:** Lifetime Warranty
- DRIVE TRAIN:** 2 Year Warranty
- ELECTRONICS:** 2 Year Warranty
- BATTERIES:** 1 Year Warranty





# BUZZ AROUND CarryOn™

The patented\* Buzzaround CarryOn foldable scooter is small in stature and big on fun! This super portable, sporty new scooter will take you where you want to go in comfort and style! It features an infinitely adjustable tiller, full front and rear light package, lithium-ion battery (travel battery sold separately\*\*), and 300 lb. weight capacity. The heaviest piece is 47.5 lbs. (without seat and batteries) and the height when folded is only 12". It disassembles easily into two pieces and is designed to store in the smallest compartments. Pack your dreams and go anywhere with the CarryOn!

**ALL-NEW Sporty Foldable Scooter! Airline Friendly!\*\*  
USB Port Standard!**



## BUZZAROUND CARRYON FEATURES:

- Designed for minimal storage space
- 4-point carry handles for easy transport
- 12" height when folded (without seat)
- 4" ground clearance
- 300 lb. weight capacity
- Flip up armrests
- Infinitely adjustable tiller
- USB port & water bottle holder on tiller
- Disassembles easily
- Long range coverage up to 18 miles
- Easy access 15 AH lithium-ion battery standard
- Optional airline friendly 6.5 AH lithium-ion battery removable for carry on (Available Late Summer 2020)



GB120

SCAN FOR VIDEO



## WARRANTIES:

**FRAME:** Lifetime Warranty  
**DRIVE TRAIN:** 2 Year Warranty  
**ELECTRONICS:** 2 Year Warranty  
**BATTERIES:** 1 Year Warranty

\* Patent Pending.

\*\* Optional 6.5 AH lithium-ion battery comes with a stylish bag and can be removed and taken as a carry on while traveling by air, sea or train.

Optional color shrouds available.  
\*Upcharge applies.



GB107

# BUZZ AROUND<sup>LT</sup>

## BUZZAROUND LT FEATURES:

- 300 lb. weight capacity
- Increased foot room
- Color through chip resistant plastic (red only)
- Charger port in the tiller
- Quickly disassembles into five, lightweight pieces for easy transport
- Comes standard in matte red
- Pink, yellow, silver and camouflage color shrouds available
- Accepts new interchangeable color shrouds

## WARRANTIES:

**FRAME:** Lifetime Warranty

**DRIVE TRAIN:** 2 Year Warranty

**ELECTRONICS:** 2 Year Warranty

**BATTERIES:** 1 Year Warranty



(Buzz XL  
model shown)

# LiteRider<sup>®</sup> GOLDEN

## LARGE PORTABLE

The unique tubular design of the LiteRider scooters gets you wherever you need to go in a style all their own. The only scooters to come with two cargo baskets standard, the LiteRiders are easy to disassemble, and maneuver smoothly indoors and out.



### **WARRANTIES:**

**FRAME:** Lifetime Warranty

**DRIVE TRAIN:** 2 Year Warranty

**ELECTRONICS:** 2 Year Warranty

**BATTERIES:** 1 Year Warranty



### **GL111 & GL141 FEATURES:**

- 300 lb. weight capacity
- Wraparound Delta Tiller standard
- Easy to read control panel with backlit battery gauge
- Two storage baskets standard
- Comfortable 17"x16" stadium style seat standard
- Adjustable armrests
- Optional 20" x 18" stadium style seat
- 2-22 AH Batteries; increased operating range up to 12 miles
- Completely wireless disassembly
- 9" front and rear tires
- New Ultrabright LED adjustable headlight and rugged plastic basket standard!

# GOLDEN COMPANION

FULL-SIZE LUXURY

**ALL NEW Rugged Steel Frame with Weight Capacity up to 400 lbs.!**

The Golden Companion scooters are your ticket to ride first class every day! Allow the luxurious, full-size Companions to enhance your life and take you to the places you've always had trouble getting to. Shop in the mall, attend sporting events and enjoy the great outdoors with the unmatched comfort of your favorite new Companion! All Companion models feature:



*Companion GC340 shown with optional power lift seat.*

## **WARRANTIES:**

**FRAME:** Lifetime Warranty  
**DRIVE TRAIN:** 2 Year Warranty  
**ELECTRONICS:** 2 Year Warranty  
**BATTERIES:** 1 Year Warranty

SCAN FOR VIDEO



## **COMPANIONS FEATURE:**

- All new design for easy disassembly
- Wraparound Delta Tiller standard
- Adjustable LED headlight
- EZ Reach tiller adjustment lever
- Auto lock tiller
- Easy to read control panel
- Easy access charging port
- Black, non-marking, low profile tires on sturdy steel rims
- LED battery gauge
- Comfortable, fully adjustable padded armrests
- Under floor mat troubleshooting guide
- Removable shopping basket with handle
- Four mobility transportation tie downs



# golden **Avenger™**



GA541  
(Freight Supplement applies.)

When you want rugged, heavy duty scooters to tackle the challenging terrain of the great outdoors, look no further than the Golden Avenger™ and Patriot™. Both feature oversized tires with sporty mag wheels, a full-size captain's seat for all day riding comfort, and an impressive operating range of up to 18 miles.

## **AVENGER FEATURES:**

- 500 lb. weight capacity
- Adjustable tiller for riding comfort
- Deluxe, reclining high back captain's seat
- Large, 13" front and rear tires
- Adjustable front headlight and four-way flashers
- Large front cargo basket standard
- Operating range up to 18 miles per charge
- Batteries not included: 2-NF22 batteries optional

# Golden **Patriot™**

## **PATRIOT FEATURES:**

- Full front and rear lighting package with tail lights, headlights and flashers
- Four-wheel independent suspension for smoothly traversing outdoor terrain
- Rugged rubber mat and anti-skid material on the floor deck
- Convenient water bottle holders for two cold beverage bottles
- Full control panel featuring safety horn, flashers, lights, turn signals, speed knob and backlit battery gauge
- Black, pneumatic tires with brushed alloy wheels
- 13-Inch Pneumatic Tire and Rims
- Batteries not included: 2-NF22 or 2-24 batteries optional

## **WARRANTIES:**

- FRAME:** 5 Year Warranty
- DRIVE TRAIN:** 2 Year Warranty
- ELECTRONICS:** 2 Year Warranty
- BATTERIES:** 1 Year Warranty

## **WARRANTIES:**

- FRAME:** 5 Year Warranty
- DRIVE TRAIN:** 2 Year Warranty
- ELECTRONICS:** 2 Year Warranty
- BATTERIES:** 1 Year Warranty



GR575  
(Freight Supplement applies.)

# Live the Golden Life

Test drive the industry's most comprehensive line of consumer power wheelchairs today!

## THE NEW PORTABLE LITERIDER ENVY'S DESIGN MAKES OUR POWER WHEELCHAIRS THE EASIEST TO TRANSPORT!

### LITERIDER ENVY FEATURES:

- Lightweight & Convenient
- Disassembles into 4 pieces
- 22 Amp Hour battery for all-day range
- Off board battery charging for on the go

### WARRANTIES:

- FRAME:** 5 Year Warranty
- DRIVE TRAIN:** 13 Month Warranty
- ELECTRONICS:** 13 Month Warranty
- BATTERIES:** 12 Month Warranty

SCAN FOR VIDEO



**ONLY  
36 lbs.!**





# LiteRider<sup>®</sup> GOLDEN ENVY

*Available in 6 Vibrant Colors  
with a Larger Foot Plate and  
New Shroud Design!*

**Golden's RETAIL CASH POWER  
WHEELCHAIR features:**

- 300 lb. weight capacity
- Operating range up to 15.5 miles
- Turning radius: 28.5"
- NEW premium 17" Seat Option

## SOLAR FLARE ORANGE

Shown with NEW high back  
premium seat option

GP162



## ENVY GREEN



## CHERRY RED



## PARADISE BLUE



## SUNBURST YELLOW



## SATIN SILVER



# Experience Constant Contact



Inclines and thresholds are no match for the Golden Compass Sport's true center wheel drive design with our patented springless articulating chassis!

This unique suspension system keeps all six wheels on the ground simultaneously, eliminating the rocking sensation when traversing inclines, giving both Compass™ models unparalleled curb and threshold climbing ability. No other power base can match the constant wheel to ground contact of a Compass!



golden  
**COMPASS™  
SPORT**

GP605

## **COMPASS SPORT FEATURES:**

- NEW All Black Seat
- NEW aluminum wheels provide increased durability
- Medicare Code: K0823
- 300 lb. Weight Capacity
- Turning Radius: 19.5"
- Batteries 2-U1
- Single Post Seat Mount
- 40 Amp Controller Standard

## **WARRANTIES:**

**FRAME:** Lifetime Warranty

**DRIVE TRAIN:** 13 Month Warranty

**ELECTRONICS:** 13 Month Warranty

**BATTERIES:** 12 Month Warranty

## **THE COMPASS HD HAS THE MOST STABLE BASE IN ITS CLASS!**

Compass™ models provide superior forward stability making them the safest power wheelchairs to mount or dismount. The Compass line has fewer moving parts than a power base with springs for the best forward stability and no tipping!



golden  
**COMPASS™  
HD**

GP620

## **COMPASS HD FEATURES:**

- NEW All Black Seat
- Medicare Code: K0825
- 450 lb. Weight Capacity
- Turning Radius: 20"
- Four Post Seat Mount with Slide Seat Bracket
- 75 Amp Controller Standard
- Batteries not included: 2-NF22 batteries optional



*You can't stand like this on any other power wheelchair base.*

## **WARRANTIES:**

**FRAME:** 5 Year Warranty

**DRIVE TRAIN:** 13 Month Warranty

**ELECTRONICS:** 13 Month Warranty

**BATTERIES:** 12 Month Warranty



# Alanté SPORT

- Low-maintenance Color Through Shroud Design!
- Larger, Lightweight Composite Footplate!
- Generous High Back Captain's Seat



## Alanté SPORT

GP208



### ALANTÉ SPORT FEATURES:

- Low center of gravity
- Excellent maneuverability
- 28" turning radius
- 300 lb weight capacity
- Upgraded user friendly motor connections
- Auto shut-off charger system
- Electromagnetic braking system
- Padded flip up arms standard
- 9" driving wheels
- Larger 6" rear caster wheels

### WARRANTIES:

- FRAME:** 5 Year Warranty
- DRIVE TRAIN:** 13 Month Warranty
- ELECTRONICS:** 13 Month Warranty
- BATTERIES:** 12 Month Warranty

# Ride in Style

Secure Your Personal Items with Scooter and Power Wheelchair Accessories by Golden!



Travel Case

Optional BuzzAround Shroud Cover Colors: Pink, Silver, Yellow & Camouflage  
(Red & Blue come standard on most models. GB107 standard in red only.)



Cane Holder



Scooter Cover



Quad Cane Holder



Rear View Mirror



Crutch Holder



Rear Basket



Oxygen Tank Holder



Power Elevating Seat\*



Walker Holder



Safety Flag

*(Included with most models.)*



Cup Holder



Pack'n'Go Armbag



Cell Phone Holder

**PACK BUNDLE**

\*Only available on Companion models.

# SPECIFICATIONS\*

Item	Description	Length			Width			Front Tire	Rear Tire	Center Tire	Casters	Weight Capacity (lbs.) <sup>1</sup>	Maximum Speed (mph)	Ground Clearance	Operating Range <sup>2</sup> (miles)	Turning Radius	Batteries (2)**	Weight w/ Batteries	Front Half Weight	Rear Half Weight	Seat Weight w/ Arms	Heaviest Piece	Base Weight
		40.5"	21.5"	8"	N/A	N/A	300																
GB107	Buzz LT (3-wheel)	40.5"	21.5"	8"	N/A	N/A	300	4	3"	8	37"	12 AH	97	33	25	20	33	N/A				N/A	
GB117D	Buzz XL (3-wheel)	42.5"	21.5"	8"	N/A	N/A	300	4	3"	8	37"	12 AH	102	35	25	23	35	N/A				N/A	
GB117H	Buzz XLHD (3-wheel)	42.5"	21.5"	8"	N/A	N/A	325	4	3"	14	37"	22 AH	112	35	25	23	35	N/A				N/A	
GB117Z	Buzz XLSHD (3-wheel)	42.5"	21.5"	8"	N/A	N/A	325	4	3"	14	37"	22 AH	115	35	28	23	35	N/A				N/A	
GB147D	Buzz XL (4-wheel)	41.5"	21.5"	8"	N/A	N/A	300	4	3"	8	47"	12 AH	110	43	25	23	43	N/A				N/A	
GB147H	Buzz XLHD (4-wheel)	41.5"	21.5"	8"	N/A	N/A	325	4	3"	14	47"	22 AH	120	43	25	23	43	N/A				N/A	
GB147Z	Buzz XLSHD (4-wheel)	41.5"	21.5"	8"	N/A	N/A	325	4	3"	14	47"	22 AH	123	43	28	23	43	N/A				N/A	
GB118	Buzz EX (3-wheel)	45"	22"	9"	N/A	N/A	350	5	4"	18	40.5"	U1	154	46	36	23	46	N/A				N/A	
GB148	Buzz EX (4-wheel)	44"	22"	9"	N/A	N/A	350	5	4"	18	56"	U1	161	53	36	23	53	N/A				N/A	
GB119	Buzz LX (3-wheel)	45"	22"	9"	N/A	N/A	375	5	4"	18	40.5"	U1	168	50	36	32	50	N/A				N/A	
GB149	Buzz LX (4-wheel)	44"	22"	9"	N/A	N/A	375	5	4"	18	54"	U1	174	57	36	32	57	N/A				N/A	
GB120	BuzzCarryOn (4-wheel)	38"	19.5"	6"	N/A	N/A	300	3.6	4"	18"	48.75"	15 AH	65.5	47.5	N/A	12	47.5	N/A				N/A	
GL111	LiteRider (3-wheel)	41.5"	22"	9"	N/A	N/A	300	5	3.75"	12	39"	22 AH	125	44	28	21	44	N/A				N/A	
GL141	LiteRider (4-wheel)	43.25"	22"	9"	N/A	N/A	300	5	3.75"	12.5	48"	22 AH	133	52	28	21	52	N/A				N/A	
GC240	Companion (3-wheel)	42.25"	23.5"	9"	10.75"	N/A	350	5.3	3.75"	19.9	37"	U1	163	45	48	26	48	N/A				N/A	
GC340	Companion (3-wheel)	45.5"	23.5"	9"	10.75"	N/A	400	5.3	4.5"	20.9	41"	U1	176	56	48	26	56	N/A				N/A	
GC440	Companion (4-wheel)	47.5"	23.5"	10.75"	10.75"	N/A	400	5.3	4.5"	18.9	56.5"	U1	194	74	48	26	74	N/A				N/A	
GA541	Avenger	57.5"	24.5"	13"	13"	N/A	500	7	4"	18 <sup>2</sup>	67"	NF22**	265	68	55	68	N/A					N/A	
GR575	Patriot	57"	28"	13"	13"	N/A	400	10	6.75"	16	76.5"	75 AH**	302	N/A	N/A	44	152	152				152	
GP162	LiteRider Envy	36"	22.5"	N/A	9"	N/A	300	3.5	2.5"	15.5	28.5"	22 AH	123	37	36	23	37	N/A				N/A	
GP208	Alanté Sport	41.5"	25"	9"	N/A	N/A	300	4.3	2.25"	15 <sup>2</sup>	28"	U1	146.5	N/A	N/A	37	66	66				66	
GP605	Compass Sport	39"	24"	N/A	N/A	10"	300	4	3.25"	24	19.5"	U1	185	N/A	N/A	46	92	92				92	
GP620	Compass HD	40"	25"	N/A	N/A	10"	450	4	3.5"	12 <sup>2</sup>	20"	NF22**	226	N/A	N/A	51	101	101				101	

\* All weight shown in lbs.  
 \*\* Batteries not included.

All products feature Freewheel Mode, Electronic Speed Control and Off-Board Charger.

**NOTE:** Literature is current at time of printing. Golden Technologies reserves the right to make changes to the product or literature at any time.

<sup>1</sup> Battery range at 200 lbs., and will vary due to rider weight, drive surface, terrain and battery type.  
<sup>2</sup> Battery range at maximum weight capacity, and will vary due to rider weight, drive surface, terrain and battery type.

Brakes: All Golden Technologies scooters come standard with Dynamic/Electromechanical braking systems.  
 Rear wheel drive: All Golden Technologies scooters come standard with a sealed transaxle, 24 volt motor.

# Choose from a Complete Line of Lifestyle Solution Products by Golden



## SCOOTERS

- BuzzAround Series
- Lite Rider Series
- Companion Series
- Heavy Duty Series



## POWER WHEELCHAIRS

- Compass Sport
- Compass HD
- LiteRider Envy
- Alanté Sport



## POWER LIFT RECLINERS

- Value Series
- Signature Series
- Comforter Series
- Comforter Wide Series
- MaxiComforter Series
- PATENTED Twilight Technology Positioning Features!

Visit Your Local Authorized Golden Dealer.