

Elevate & Circulate



PR761 EZ Sleeper
shown in Brisa® Shiitake

 **GOLDEN**
POSITIONING YOU FOR LIFE



A MESSAGE FROM RICH GOLDEN, CEO

Golden is known for comfort, quality and innovation. Since 1985, we have produced power lift recliners that push the boundaries of expectation, providing seating options with unique body positions.

In 2018, we launched the MaxiComfort® Cloud with Twilight Technology, a luxury power recliner capable of providing never before seen or felt levels of comfort, thanks to our patented three Comfort Zone positioning system. In 2019, we launched the next generation of patented Twilight Technology offering five Comfort Zones. It offers you a new range of rejuvenating positions to comfort the body and ease the spirit.

Our chairs are handcrafted with the highest quality components available, engineered to give you maximum service with minimum maintenance, and backed not only with a lifetime warranty on our lift frame and scissor mechanism, but our wooden frame as well. We hope your Golden chair provides years of comfort and enjoyment while sitting down and relaxing as well as standing up when lift assist is needed most.

We are proud to manufacture America's favorite power lift recliners - for your family, and ours.

OUR VISION...

To build the most innovative and comfortable products that put you and your family's lives in motion.



Features that make Golden Power Lift Recliners America's #1 brand!

Golden Technologies is the world's leading manufacturer of power lift recliners. Founded in 1985, Golden has become America's most trusted and reliable brand with the most innovative products and features. All chairs are handcrafted by American workers with domestic and imported materials (with the exception of the DeLuna Series). Each chair undergoes a 42+ point quality control inspection before the product is approved for shipping.

Golden lift chairs are Class II Medical Devices.



TWILIGHT

Experience the ultimate in personalized comfort with the EZ Sleeper. This luxury power lift recliner offers a modern look and features our exclusive 5-Zone positioning system with patented technologies for hours of unmatched comfort.



PROUDLY MADE

 BY AMERICANS
 with Domestic & Imported Materials

PR761 shown in
 Brisa® Shiitake*

*Upcharge applies


We Patented Comfort™

Golden holds multiple patents on our technology, allowing us to offer an EXCLUSIVE comfort experience.




PATENTED Twilight Technology

The patented Twilight Technology includes an extra adjustable zone that provides a tilting motion, cradling your body and offering more positioning than any other recliner on the market.



PATENTED MaxiComfort® Positioning

Our patented MaxiComfort® Positioning provides ergonomic support thanks to the four-way motion of the seat combined with the independent movement of the back. This unique movement provides extra support to your lower back and can provide relief to muscle and joint pain.



PATENTED SmartTek® Diagnostic System

SmartTek® is a patented technology standard on all Golden lift recliners. SmartTek uses three LED lights to help diagnose the electrical health of the chair.



com · fort zone

noun | a place or situation where one feels safe or at ease without stress



Comfort Zone™ by GOLDEN

The Comfort Zone by Golden offers a stress-free reclining experience by providing patented, adjustable and supportive positioning zones that can increase your body's wellness and optimize the feeling of comfort.

Zone 1



LOWER BODY

The four-way motion of the seat helps provide relief to muscle and joint pain by placing your lower body in a supported ergonomic position.

Zone 2



UPPER BODY

The independent movement of the backrest, combined with the four-way motion seat box in Zone 1, creates our patented positioning that allows stress and strain to be released from your upper body by offering adjustable support for your back and shoulders.

Zone 3



HEADREST

The upward motion of the adjustable headrest allows for tension to be released from your neck and shoulders by supporting your head placement, while encouraging the natural shape of your spine.

Zone 4



LUMBAR

The custom and firm support of the adjustable lumbar can take pressure off your spine while giving your back muscles the chance to rest and rejuvenate.

Zone 5



TWILIGHT TILT

The tilting motion of our patented Twilight Technology provides a feeling of weightlessness and offers even more positioning options not available with other recliners. Encourage circulation by relaxing with your legs elevated or promote the release of muscle tension and joint pain in a comforting neutral position.

Comfort Zone™

by GOLDEN



Scan with your smartphone camera to learn more!

WE PATENTED COMFORT®

FIVE-ZONE SYSTEM

Comfort Zones 1 - 5 This positioning system includes all five adjustable comfort zones to create the ultimate wellness comfort experience with our patented Twilight Technology. Each Comfort Zone can be adjusted independently by using the exclusive AutoDrive™ hand control.



FOUR-ZONE SYSTEM

Comfort Zones 1 - 4 provide an ergonomic comfort experience with an emphasis on head and back support. Each Comfort Zone can be adjusted independently by using the exclusive AutoDrive™ hand control.



TWO-ZONE SYSTEM

Comfort Zones 1 & 2 work together to create individualized wellness positioning and comfort with our patented four-way motion seat box and independent backrest. Each Comfort Zone can be adjusted independently by using the exclusive AutoDrive™ hand control.



ONE-ZONE SYSTEM

Comfort Zone 1 provides ergonomic placement of your lower body allowing your muscles to relax for long-lasting comfort. Comfort Zone 1 can be adjusted using the two buttons on the hand control.





Shown in
Brisa® Bridle*

PROUDLY MADE

BY AMERICANS
with Domestic & Imported Materials




TWILIGHT
★★★★★

Kim K. | Google Review

"The power lift chair from Golden Technologies with Twilight Technology was the best decision. The 5 different zones with endless amounts of positioning made it possible for me to recover from an abdominal surgery. Comfort is one of the best features about this chair."



TWILIGHT



The Cloud™
with Twilight Positioning
PR515-SME shown in Sterling



The Cloud™
with Twilight Positioning
PR515-MLA shown in Brisa® Coffee Bean*



EZ Sleeper
with Twilight Positioning
PR761-MLA shown in Bourbon



These chairs offer our unique, ergonomic seating system for long-term comfort and support. They feature our custom bucket seat and plush backrest. Positioning options are endless with our 5-zone positioning system, which includes our patented Twilight Technology, patented MaxiComfort® powered positioning, adjustable headrest and lumbar support.

The EZ Sleeper features a modern look with a tailored back and stylish track arm design, while the plush chaise pad and pocketed coil spring seat provide hours of unmatched comfort. Positioning options are endless with our 5-zone positioning system, which includes our patented Twilight Technology, patented MaxiComfort® powered positioning, adjustable headrest and lumbar support.

Standard Fabrics for PR515



Calypso

Copper

Hazelnut

Shiraz

Sterling

Standard Fabrics for PR761



Hazelnut

Bourbon

Smoke


*Upcharge applies

Deborahannette B | Google Review

"Amazing family owned and run business that takes pride in what they do. The chairs built are not only beautiful, but provide comfort and so many health features for all ages."

MaxiComfort Day Dreamer PR632-MED shown in Brisa® Coffee Bean*

This chair is a medium-sized power lift recliner offering our exclusive PowerPillow™ articulating headrest and power lumbar support that comes standard with The Cloud™ seating system and features our custom bucket seat design. Adjust your head and neck position at the touch of a button.



Standard Fabrics for PR632



Calypso



Copper



Hazelnut



Shiraz



Sterling

*Upcharge applies





GOLDEN AutoDrive™ 3.0

CUSTOMIZED COMFORT

Designed with the consumer in mind, every feature of the AutoDrive™ hand control has been purposefully added to ensure personalized comfort and optimal wellness positioning options, available through easy-to-use features, making operating your recliner a breeze.

Easy-Grip Design – A tapered design with a finger grip indentation on the back and “soft touch” finish helps keep the controller firmly in place in anyone’s grasp.

Button Programmability – Four programmable buttons help the user find their perfect position again and again without hassle.

SmartTek® Diagnostic System – The electrical health of the chair is communicated using LED lit buttons, allowing the user to know immediately if power to the chair has been lost.

Quick Disconnect – This feature allows for the hand control to be easily replaced, if needed.

In addition to the physical features, all AutoDrive™ hand controls offer the following functional features:

- **Full Lockout** – Disables all the buttons on the hand control
- **Partial Lockout** – Disables all buttons except “Rise” and “Recline”
- **Lights Out** – Turns off the light that illuminates the buttons
- **Factory Default Reset** – Resets the Position Control buttons back to the original factory settings
- **Calibration** – Fine-tunes your positioning and keeps your recliner movements accurate

(See Owner’s Manual for more details.)

AVAILABLE IN THE THREE VERSIONS, THIS UNIQUE HAND CONTROL OPERATES OUR 2, 4 & 5-ZONE POSITIONING SYSTEMS, ALLOWING FOR BASIC, CUSTOM OR POSITION CONTROL OPTIONS.



AUTODRIVE



Scan with your smartphone camera to learn more!

BASIC CONTROL – The top two buttons take you to the basic lift and recline functions of the chair while still allowing you to experience our PATENTED MaxiComfort® positioning.

POSITION CONTROL – The labeled buttons in the middle are preset, but can also be programmed to your personalized positions, offering a customized reclining experience every time.

ZG – Zero Gravity offers a supportive feeling of weightlessness.

Sit – Provides an easy and stress-free way to return your chair to the proper seated position.

TV – Takes you to the perfect position to relax while watching TV or reading a book.

Rejuv – Allows you to elevate your feet and legs comfortably for a long period of time. *(Position only available on 5-Zone recliners).*

Sleep – Takes your body to a traditional lay flat position, which can be comfortable for napping during the day or sleeping at night. *(Button only available on 2 and 4-Zone recliners).*

CUSTOM CONTROL – The arrow buttons toward the bottom of the hand control allow you to adjust each zone independently which is the best way to truly find your perfect position to program.



PROUDLY MADE

BY AMERICANS
with Domestic & Imported Materials

PR510 shown in Brisa®
Coffee Bean*



 GOLDEN



Cheryl S. | Google Review

"It has been a life changer for my husband. Living with Parkinson's is challenging enough, but the MaxiComfort lift chair has added to his independence and comfort immensely. It is his home base."

MaxiComfort SERIES

The Cloud™ offers our unique, ergonomic Cloud seating system for long-term comfort and support. They feature our custom bucket seat, plush backrest, and solid molded high-density arm pillows for the ultimate in relaxation and rejuvenation.



The Cloud™
PR510-SME *shown in Calypso*



The Cloud™
PR510-MLA *shown in Brisa® Coffee Bean**



Upgrade with
Deluxe Therapeutic
Heat & Massage
on these models!



The Relaxer is the original MaxiComfort® power lift recliner. These chairs feature a luxurious three-pillow waterfall back with individual zippers for customizable comfort to relieve pressure and improve circulation. A full chaise pad provides lower body and leg support.



Relaxer
PR766-SMA
shown in Calypso




Relaxer
PR766-MED
*shown in Brisa® Coffee Bean**



Relaxer
PR766-LAR
shown in Hazelnut

Standard Fabrics for PR510 & PR766



Calypso



Copper



Hazelnut



Shiraz



Sterling

*Upcharge applies

SMALL TO TALL, WE BUILD THEM ALL!

Our long-standing, small-to-tall sized collection of power lift recliners now features a contemporary design and plush fabric offerings! You'll notice our popular MaxiComforter Series in a whole new light.



**America's Favorite Lift Chair
for Over 20 Years!**

**An Old Favorite
Meets New
Upgrades in the
MaxiComforter
Power Lift Recliner!**



*PR535 shown
in Anchor*

MaxiComforter SERIES

small^{to}TALL
WE BUILD THEM ALL



Upgrade with
Deluxe Therapeutic
Heat & Massage
on these models!



MaxiComforter Petite/Small
PR535-PSA shown in Sandstorm



MaxiComforter Medium
PR535-MED shown in Port



MaxiComforter Large
PR535-LAR shown in Anchor



MaxiComforter Tall
PR535-TAL shown in Bittersweet




MaxiComforter Medium Wide
PR535-M26 shown in Oxford
AutoDrive™ not available on this model.



Standard Fabrics for PR535



Anchor



Bittersweet



Oxford



Port



Sandstorm

The Comforter Series is all about finding the perfect size for the perfect fit. These newly upgraded power lift recliners were designed to fit individuals small to tall with a weight capacity up to 500 lbs. They feature our signature seat, which provides the utmost in comfort and support, and a new arm style with additional plush padding for more comfort and value.



The Comforter Series has been reinvigorated with a design refresh and plush fabric offering!



PR531 shown in Sandstorm



Bubba D. | Google Review

"THE BEST Chair I have ever had and in 25 years I have had 5. My only complaint is that it is TOO comfortable! So much so that I rarely use my bed. Worth every single penny."

Comforter SERIES

small to TALL
WE BUILD THEM ALL



Upgrade with
Deluxe Therapeutic
Heat & Massage
on these models!



Comforter Petite/Small
PR531-PSA *shown in Oxford*



Comforter Medium
PR531-MED *shown in Sandstorm*



Comforter Large
PR531-LAR *shown in Bittersweet*



Comforter Tall
PR531-TAL *shown in Anchor*

Comforter Wides

The Comforter Wide is perfect for those who need more space between the armrests. Constructed using the highest quality components for reliable, steady lifting and reclining, these chairs offer a weight capacity of up to 500 lbs.



Comforter Wide Small
PR531-S23 *shown in Port*



Comforter Wide Medium
PR531-M26 *shown in Anchor*




Comforter Wide Tall
PR531-T28 *shown in Bittersweet*




Standard Fabrics for PR531



Anchor



Bittersweet



Oxford



Port



Sandstorm

Signature SERIES

The popular Cambridge power lift recliner features our exclusive hybrid sizing for a better fit. The Small/Medium size fits people 5'1" to 5'6" and the Medium/Large size fits people 5'7" to 6'2". The Cambridge features a supportive biscuit backrest design and a soft comfortable seat for hours of relaxation.



Upgrade with
Deluxe Therapeutic
Heat & Massage
on these models!




Cambridge
PR401-SME *shown in Calypso*



Cambridge
PR401-MLA *shown in Copper*



Standard Fabrics for PR401




Calypso




Copper



Hazelnut



Shiraz



Sterling



GOLDEN OFFERS THE HIGHEST QUALITY RECLINER FABRICS

When it comes to performance, comfort, and style, Golden Technologies' array of power lift recliners is unmatched in the industry. It's no wonder that the fabric we use gives you the same high quality and range of choices to please every taste and need – including our exclusive luxury Brisa® and high-performance fabric with Alta™ technology, plus a full palette of plush fabrics and colors. At Golden, we understand how important your lift chair is to you, and we want you to relish in its appearance and durability, as well as its functionality.

THE BRISA® COLLECTION




Cream Buckskin Saddle Bridle Coffee Bean Shiitake Iron Black Onyx Beet Root Night Navy

Supple Brisa® fabric offers the refined look of leather, but in a luxuriously soft, breathable, and temperate design. Brisa fabric comes in an extensive selection of lovely colors to choose from and will never stick to you or feel too hot or cold. This beautiful fabric is scratch resistant and easy to clean and disinfect, so it always stays looking new! These fabrics are double rubbed 300,000 times to ensure the highest level of quality and durability. We're so sure of the quality of this fabric, that Golden offers a Limited Lifetime Warranty* on all Brisa fabrics.

** Limited Lifetime Warranty, as long as normal use is maintained. Restrictions apply.*

NEW HIGH-PERFORMANCE FABRIC with Alta™ Technology

Our high-performance fabric with Alta™ technology is perfect for heavy-duty or high-traffic seating. Alta ramps up performance with a total liquid barrier, antimicrobial protection, and high repellency to liquids and oil-based stains. Exceptional cleanability reduces the need for harsh cleaning agents. Durablock™ liquid barrier provides breathable protection. Golden offers a 10-year warranty on all high-performance fabric with Alta technology.



NOTE: Please view all Golden Fabrics in person for accuracy.



Scan with your smartphone camera to learn more about Brisa®



Scan with your smartphone camera to learn more about Alta®





DeLuna™

Restful nights, energizing days.

Inspired by the moon, the DeLuna series of luxury power lift recliners is designed to offer you the most calming comfort and restorative sleep.

Just like the many moons, each DeLuna recliner has unique visual characteristics. Enjoy the powerful energy and take time to better connect with yourself in your DeLuna recliner by Golden!

*PR448 shown
in Maple*



GOLDEN



Keith F. | Google Review

"I received my chair as a gift and I have loved it ever since. So durable and very comfortable. I look forward each day to using it."

** These models were designed and engineered in Old Forge, PA, and assembled in China using Golden components.*




DeLuna™ SERIES

The **Dione** and **Titan** luxury power lift recliners offer stylish, overstuffed seamed back designs, plush chaise pads, and numerous standard features providing hours of unmatched comfort. Enjoy endless positioning options with our 4-Zone Comfort Zone positioning system, which includes our patented MaxiComfort® positioning, adjustable headrest and adjustable lumbar all operated by our programmable AutoDrive™ hand control.

Standard Fabrics for PR446 & PR448

No crossover fabrics available.



Graphite

Maple



Dione

PR446-MED* shown in Maple



Dione

PR446-LAR* shown in Graphite



Titan

PR448-MED* shown in Maple



The **Dione** and **Titan** models feature a convenient USB charging port for your phone or tablet located on the side of the chair. The **Titan** offers even more extras, such as a built-in cup holder, wireless smartphone charger and an accessory bracket with a removable table conveniently located in the open-arm design.

The **Elara** is medium-sized with a back style that provides maximum comfort in a tasteful seam design. It features a manually adjustable three-position headrest, supportive chaise pad, and a USB port on the hand control. Available in luxurious new fabrics to match any décor.

Standard Fabrics for PR118 & PR202

No crossover fabrics available.



Antler

Elk

Lakefront

Chestnut



Elara

PR118-MED* shown in Antler



The **Eirene** is medium-sized with a back style that provides maximum comfort in a tasteful button design. It features a supportive seating pad, an easy grip armrest design, and a USB port on the hand control for your phone or tablet.



Eirene

PR202-MED* shown in Lakefront



Specifications

	Name	Model	Comfort Zones	Weight Capacity	Overall Width	Overall Height	Seat to Top of Back	Width Between Arms	Seat Depth	Floor to Top of Seat	Extension - While Reclined from Back of Seat to End of Footrest	Heat & Massage Option Available	Footrest Extension Option Available*	Straight Lift Option Available
MaxiComfort	EZ Sleeper w/ Twilight	PR761-MLA	5	375 lbs.	34"	46"	30"	22"	21"	21.5"	43"	No	No	No
	Cloud w/ Twilight	PR515-SME	5	375 lbs.	35.5"	44"	28"	20"	18.5"	20"	43"	No	No	No
	Cloud w/ Twilight	PR515-MLA	5	375 lbs.	39.5"	46.5"	31"	20"	20"	21.5"	43"	No	No	No
	Cloud SME	PR510-SME	2	375 lbs.	33.5"	41"	27"	18.5"	19.5"	19"	43"	Yes	No	Yes
	Cloud MLA	PR510-MLA	2	375 lbs.	38"	45"	29"	20"	21.5"	20"	44"	Yes	No	Yes
	Day Dreamer	PR632-MED	4	375 lbs.	38"	45"	31"	22.5"	20.5"	20"	44"	No	Yes	Yes
	Relaxer SMA	PR766-SMA	2	375 lbs.	34"	41.5"	26"	18.5"	19.5"	19.5"	43"	Yes	No	Yes
	Relaxer MED	PR766-MED	2	375 lbs.	36"	44.5"	28.5"	20.5"	20.5"	20"	44"	Yes	Yes	Yes
	Relaxer LAR	PR766-LAR	2	375 lbs.	38.5"	48"	29.5"	22.5"	22.5"	22"	47"	Yes	Yes	Yes
	Comforter PSA	PR535-PSA	2	375 lbs.	32"	40"	26"	21"	19"	20"	40"	Yes	Yes	Yes
	Comforter Md.	PR535-MED	2	375 lbs.	32"	42"	27"	21"	19"	20"	41"	Yes	Yes	Yes
	Comforter Lg.	PR535-LAR	2	375 lbs.	34"	43"	29"	23"	21"	21"	42"	Yes	Yes	Yes
	Comforter Tall	PR535-TAL	2	375 lbs.	34"	48"	31"	23"	21"	23"	45"	Yes	Yes	Yes
	Comforter Med-Wide	PR535-M26	2	500 lbs.	37"	45"	30"	27"	20"	20"	43"	Yes	Yes	Yes
Comforter	Petite/Small	PR531-PSA	1	375 lbs.	31.5"	40"	27"	21"	19"	19"	41"	Yes	Yes	Yes
	Medium	PR531-MED	1	375 lbs.	32"	41"	27"	21"	19"	20"	41"	Yes	Yes	Yes
	Large	PR531-LAR	1	375 lbs.	34"	44"	28.5"	23"	20.5"	21"	42"	Yes	Yes	Yes
	Tall	PR531-TAL	1	375 lbs.	34.5"	47"	30.5"	23"	22"	22.5"	47"	Yes	Yes	Yes
	Small-Wide	PR531-S23	1	375 lbs.	34"	40.5"	27"	23"	20"	19.5"	40"	Yes	Yes	Yes
	Medium-Wide	PR531-M26	1	500 lbs.	38"	43.5"	29"	27"	20.5"	20.5"	42"	Yes	Yes	Yes
	Tall-Wide	PR531-T28	1	500 lbs.	40"	47.5"	32"	29"	22"	22"	46"	Yes	Yes	Yes
Signature	Cambridge	PR401-SME	1	375 lbs.	33"	41.5"	27"	20"	19"	19"	41.5"	Yes	Yes	Yes
	Cambridge	PR401-MLA	1	375 lbs.	35"	44.5"	31"	22"	22"	22"	44"	Yes	Yes	Yes
DeLuna	Eirene	PR202-MED	1	375 lbs.	32"	43"	28"	21"	19.5"	20.5"	43"	No	No	Yes
	Elara	PR118-MED	1	375 lbs.	32"	45"	29"	21"	20"	21"	42"	No	No	Yes
	Dione	PR446-MED	4	375 lbs.	34"	43.5"	28"	22"	21"	20.5"	42.5"	No	No	Yes
	Dione	PR446-LAR	4	375 lbs.	35"	46"	30"	22.5"	21.5"	22"	44"	No	No	Yes
	Titan	PR448-MED	4	375 lbs.	38.5"	46"	31"	20.5"	21"	20.5"	43"	No	No	Yes

*May operate differently depending on model.

Golden is the world's leading manufacturer of power lift recliners. Founded in 1985 and family owned, Golden has become America's most trusted and reliable brand. Each and every power lift recliner must pass a 42+ point inspection at our manufacturing facility located in Old Forge, PA, to ensure the highest quality before packaging.




<https://m.youtube.com/user/GoldenTechChannel>



www.facebook.com/goldentechnologies



www.goldentech.com



Visit your Authorized Golden Dealer:

Note: Literature is current at time of printing. Golden Technologies reserves the right to make changes to the product or literature at any time.

BR-LIFT042022