



Helio C2

ORDER FORM - CANADA

1- Purchasing information

P.O. number:* _____ Date: _____

Company Name:* _____

Purchasing contact:* _____

Phone:* _____

Email:* _____

Contact preference:* Phone Email

Marked for: _____

3- Shipping information

Ship to:* _____

Name: _____

Address: _____

City: _____ Province: _____

Postal Code: _____

Phone: _____

2- Dealer/Representative information

Dealer/Representative:* _____

Phone:* _____

Email:* _____

Contact preference:* Phone Email

* Required		MA080 <input type="checkbox"/>	MA081 <input type="checkbox"/>	Taxes not included					
		English Manual (STD)	Manuel Français						
MODEL									
WC05	<input type="checkbox"/>	Ultralight Folding Carbon Fiber Helio C2 Wheelchair (265 lb weight capacity)		\$3,995					
WC09	<input type="checkbox"/>	HD Heavy Duty package (350 lb weight capacity)		\$180					
WC10	<input type="checkbox"/>	E-Motion adaptor Kit (2 adaptors) - Only necessary with 6° camber		\$35					
Frame									
FREE Growth Kit within the first 5 years of the Helio C2 wheelchair ownership The growth kit includes carbon crossbraces, crossbrace links, seat sling, back upholstery and footplates.									
SEAT WIDTH ³									
SW	<input type="checkbox"/>	14" <input type="checkbox"/>	15" <input type="checkbox"/>	16" <input type="checkbox"/>	17" <input type="checkbox"/>	18" ¹ <input type="checkbox"/>			NC
	<input type="checkbox"/>	19" ¹ <input type="checkbox"/>	20" ¹ <input type="checkbox"/>	21" ² <input type="checkbox"/>	22" ² <input type="checkbox"/>			WAMA	\$360
SEAT DEPTH									
SD	<input type="checkbox"/>	14" <input type="checkbox"/>	15" <input type="checkbox"/>	16" <input type="checkbox"/>	17" <input type="checkbox"/>	18" <input type="checkbox"/>			NC
	<input type="checkbox"/>	19" <input type="checkbox"/>	20" <input type="checkbox"/>					WAMB	\$275
TRANSIT TIE DOWN									
AC04	<input type="checkbox"/>	Transit Tie Down (Unoccupied)							NC
AC05-S	<input type="checkbox"/>	WC-19 Transit Tie Down w/Q'straint belt ⁴		Short belt 20"-41"					\$365
AC05-L	<input type="checkbox"/>	WC-19 Transit Tie Down w/Q'straint belt ⁴		Long belt 20"-46"					\$365
AC35	<input type="checkbox"/>	WC-19 Transit Tie Down without belt ⁴							\$335
1- Available with HD option; 2- Only available with HD option 3. Width calculator 3- May interfere with swing away armrest and Clamp on height adjustable push handle; Not available with HD option									
NEWTON ULTRALIGHT PACKAGE									
PK01	<input type="checkbox"/>	Includes: Carbon fiber angle adjustable footplates + Newton Gravity Ultralight rear wheels			Save \$245				\$530

C2 FRONT Seat-to-floor Height - Please select your choice

Forks / Position		3" caster			4" caster			5" caster			6" caster			7" caster			8" caster		
		Stem bolt Length			Stem bolt Length			Stem bolt Length			Stem bolt Length			Stem bolt Length			Stem bolt Length		
		STD	+1"	+2"	STD	+1"	+2"	STD	+1"	+2"	STD	+1"	+2"	STD	+1"	+2"	STD	+1"	+2"
4"	P1	13"	14"	15"	13 1/2"	14 1/2"	15 1/2"												
	P2	13 1/2"	14 1/2"	15 1/2"	14"	15"	16"	14 1/2"	15 1/2"	16 1/2"									
5"	P1																		
	P2				14"	15"	16"	14 1/2"	15 1/2"	16 1/2"									
	P3	14"	15"	16"	14 1/2"	15 1/2"	16 1/2"	15"	16"	17"	15 1/2"	16 1/2"	17 1/2"						
	P4	14 1/2"	15 1/2"	16 1/2"	15"	16"	17"	15 1/2"	16 1/2"	17 1/2"	16"	17"	18"	16 1/2"	17 1/2"	18 1/2"			
7"	P1										16 1/2"	17 1/2"	18 1/2"	17"	18"	19"			
	P2							16 1/2"	17 1/2"	18 1/2"	17"	18"	19"	17 1/2"	18 1/2"	19 1/2"			
	P3							17"	18"	19"	17 1/2"	18 1/2"	19 1/2"	18"	19"	20"	18 1/2"	19 1/2"	20 1/2"
	P4	16 1/2"	17 1/2"	18 1/2"	17"	18"	19"	17 1/2"	18 1/2"	19 1/2"	18"	19"	20"	18 1/2"	19 1/2"	20 1/2"	19"	20"	21"
	P5	17"	18"	19"	17 1/2"	18 1/2"	19 1/2"	18"	19"	20"	18 1/2"	19 1/2"	20 1/2"	19"	20"	21"	19 1/2"	20 1/2"	21 1/2"

Height available +/- 1/4" (see stem bolt and fork combination at the end of the order form)

CASTERS						
CA01-03	<input type="checkbox"/>	Composite wheel w/ PU tire	3" x 1"			NC
CA01-04	<input type="checkbox"/>		4" x 1"			NC
CA01-05	<input type="checkbox"/>		5" x 1"			NC
CA01-06	<input type="checkbox"/>		6" x 1"			NC
CA01-07	<input type="checkbox"/>		7" x 1"			NC
CA01-08	<input type="checkbox"/>		8" x 1"			NC
CA10	<input type="checkbox"/>	Pneumatic casters	6" x 1 1/4"			\$170
CA11	<input type="checkbox"/>		8" x 1 1/4"			\$170
CA14	<input type="checkbox"/>		4" x 1 1/2"			\$100
CA15	<input type="checkbox"/>	Newton UltraCaster Composite	5" x 1 1/2"			\$100
CA16	<input type="checkbox"/>		6" x 1 1/2"			\$100
CA17	<input type="checkbox"/>	Newton UltraCaster Aluminum	4" x 1 1/2"			\$185
CA18	<input type="checkbox"/>		5" x 1 1/2"			\$195
CA19	<input type="checkbox"/>		6" x 1 1/2"			\$205
CASTER OPTIONS						
CA12	<input type="checkbox"/>	Caster Pin Locks (Available with +1" or +2" Stem bolt ONLY) ¹			WAMS	\$275
CA13	<input type="checkbox"/>	Carbon fiber Frog Legs suspension forks (Not available with caster pin locks)				\$510

1 - Not available with WC-19

C2 - REAR Seat-to-floor Height

WHEEL SIZES	Height available	Specify desired height: _____
20" (451) wheel <input type="checkbox"/>	12 1/2 to 17 1/4	
22" (501) wheel <input type="checkbox"/>	13 1/2 to 18 1/4	
24" (540) wheel <input type="checkbox"/>	14 1/2 to 19 1/4	
25" (559) wheel <input type="checkbox"/>	15 to 19 3/4	
26" (590) wheel <input type="checkbox"/>	15 1/2 to 20 1/4	

Note: Tire size, center of gravity setting and seat slope may result in a variation of +/- 1/2" of selected seat to floor height"

Footrest

SWING-IN/SWING-OUT FOOTREST						
FR03-01	<input type="checkbox"/>	60° Swing-in/swing-out footrest - 13.5" - 17"		Set at: *		NC
FR03-02	<input type="checkbox"/>	60° Swing-in/swing-out footrest - 14.5" - 19.5"		Set at: *		NC
FR03-03	<input type="checkbox"/>	60° Swing-in/swing-out footrest - 15.5" - 20.5"		Set at: *		NC
FR04-01	<input type="checkbox"/>	70° Swing-in/swing-out footrest - 13.5" - 17"		Set at: *		NC
FR04-02	<input type="checkbox"/>	70° Swing-in/swing-out footrest - 14.5" - 19.5"		Set at: *		NC
FR04-03	<input type="checkbox"/>	70° Swing-in/swing-out footrest - 15.5" - 20.5"		Set at: *		NC
FR15-01	<input type="checkbox"/>	80° Swing-in/swing-out footrest - 13.5" - 16"		Set at: *		NC
FR15-02	<input type="checkbox"/>	80° Swing-in/swing-out footrest - 14.5" - 19.5"		Set at: *		NC
FR15-03	<input type="checkbox"/>	80° Swing-in/swing-out footrest - 15.5" - 20.5"		Set at: *		NC
FR06	<input type="checkbox"/>	90° swing-in/swing-out with front mount - Standard (no extension) ¹	7.5" - 11"	Set at: *		NC
FR07	<input type="checkbox"/>	90° swing-in/swing-out with front mount - Short extension 2" ^{1,2}	9.5" - 13"	Set at: *		NC
FR08	<input type="checkbox"/>	90° swing-in/swing-out with front mount - Medium extension 4" ^{1,2}	9.5" - 14.5"	Set at: *		NC
FOOTREST OPTIONS ³						
FR09	<input type="checkbox"/>	High mount (footplate attached to the hanger with clamps) ³	5" - 13.5"	Set at: *		\$100
FR10	<input type="checkbox"/>	Elevating legrest with calf pad, Length adjustable ⁴	Left <input type="checkbox"/> Right <input type="checkbox"/> 14" - 22"	Set at:	WAMK	\$295
FR99	<input type="checkbox"/>	Omit Footrest				NC

1- Available only with two piece adjustable footplates; 2- May interfere with front casters; 3- Not available with front mount; * +/- 1/2"

4- The elevating legrest replace the footrest side selected; Not available with one piece footplate

FOOTPLATES						
FP01	<input type="checkbox"/>	Standard - Composite ¹				STD
FP02	<input type="checkbox"/>	Newton adjustable angle - Composite			WAMJ	\$155
FP03	<input type="checkbox"/>	Newton adjustable angle - Aluminum			WAMJ	\$175
FP04	<input type="checkbox"/>	Newton adjustable angle - Carbon fiber (with rubber inserts)			WAMJ	\$395
FP06	<input type="checkbox"/>	One piece adjustable angle flip-up - Aluminum		WAMJ	SEFND0015	\$215
FP07	<input type="checkbox"/>	One piece adjustable angle flip-up - Carbon fiber (with rubber inserts)		WAMJ	SEFND0015	\$480
FOOTPLATE ACCESSORIES						
FP08	<input type="checkbox"/>	Length adjustable heel loop (not available with one piece footplate)				NC
FP09	<input type="checkbox"/>	Calf strap			SEFND0010	\$57
FP30	<input type="checkbox"/>	Bodypoint Aeromesh Padded calf strap			SEFND0010	\$140

1- Not available with front mount

SWING AWAY RESIDUAL LIMB SUPPORT						
FR27	<input type="checkbox"/>	Residual limb support - flat 6"x8"	Left <input type="checkbox"/> Right <input type="checkbox"/>		SEFND0005	\$505
FR28	<input type="checkbox"/>	Residual limb support - flat 8"x10"	Left <input type="checkbox"/> Right <input type="checkbox"/>		SEFND0005	\$505
FR29	<input type="checkbox"/>	Residual limb support - flat 10"x10"	Left <input type="checkbox"/> Right <input type="checkbox"/>		SEFND0005	\$505
FR30	<input type="checkbox"/>	Residual limb support - flat 10"x14"	Left <input type="checkbox"/> Right <input type="checkbox"/>		SEFND0005	\$505
FR32	<input type="checkbox"/>	Residual limb support - contour 6"x8"	Left <input type="checkbox"/> Right <input type="checkbox"/>		SEFND0005	\$505
FR33	<input type="checkbox"/>	Residual limb support - contour 8"x10"	Left <input type="checkbox"/> Right <input type="checkbox"/>		SEFND0005	\$505

Backrest

FIXED ANGLE, HEIGHT ADJUSTABLE BACK CANES WITH PUSH HANDLE OPTIONS (MUST select ONE push handle options)							
BC06-912		Straight Carbon fiber back cane	9" - 12"	increment of 1/2"	Set at ² :		NC
BC06-1215		Straight Carbon fiber back cane	12" - 15"	increment of 1/2"	Set at ² :		NC
BC06-1518		Straight Carbon fiber back cane	15" - 18"	increment of 1/2"	Set at ² :		NC
BC06-1821		Straight Carbon fiber back cane	18" - 21"	increment of 1/2"	Set at ² :		NC
BC08-1315		8° Bend Carbon fiber back cane	13.5" - 15"	increment of 1/2"	Set at ² :		\$165
BC08-1518		8° Bend Carbon fiber back cane	15" - 18"	increment of 1/2"	Set at ² :		\$165
BC08-1821		8° Bend Carbon fiber back cane	18" - 21"	increment of 1/2"	Set at ² :		\$165
ANGLE AND HEIGHT ADJUSTABLE BACK CANE WITH PUSH HANDLE OPTIONS (MUST select ONE push handle options)							
BC09-1215	<input type="checkbox"/>	Straight Carbon fiber back cane	12" - 15"	increment of 1/2"	Set at ² :		\$330
BC09-1518	<input type="checkbox"/>	Straight Carbon fiber back cane	15" - 18"	increment of 1/2"	Set at ² :		\$330
BC09-1821	<input type="checkbox"/>	Straight Carbon fiber back cane	18" - 21"	increment of 1/2"	Set at ² :		\$330
BC10-1315		8° bend Carbon fiber back cane	13.5" - 15"	increment of 1/2"	Set at ² :		\$380
BC10-1518		8° bend Carbon fiber back cane	15" - 18"	increment of 1/2"	Set at ² :		\$380
BC10-1821		8° bend Carbon fiber back cane	18" - 21"	increment of 1/2"	Set at ² :		\$380
		Select back cane angle:	85° ¹ <input type="checkbox"/>	90° <input type="checkbox"/>	95° STD <input type="checkbox"/>	100° <input type="checkbox"/>	105° <input type="checkbox"/> 110° <input type="checkbox"/>
PUSH HANDLES OPTIONS (ONE option at a time)							
PH01	<input type="checkbox"/>	Push handles					STD
PH05	<input type="checkbox"/>	Fold down push handles ^{3,5}					\$160
PH06-920	<input type="checkbox"/>	Clamp on height adjustable push handles 9"-20" ^{3,4,5}					\$345
PH06-923	<input type="checkbox"/>	Clamp on height adjustable push handles 9"-23" ^{3,4,5}					\$345
PH02	<input type="checkbox"/>	Omit push handles					NC
STROLLER HANDLE AND STABILIZER BAR (only available with push handle)							
PH04	<input type="checkbox"/>	Stroller handle (not available with the HD option)				WAMW	\$220
PH03	<input type="checkbox"/>	Newton folding stabilizer bar (not available with stroller handle)					\$275
1- Not available with 8 degree bend cane; 2- If not indicated, default height setting is at the lowest position; measured from seat sling to top of back upholstery; 3- 265 lb weight limit; 4- Not compatible with swing away armrests, 8° bend back canes and amputee axle plates; 5- Not compatible with stroller handle, stabilizer bar, headrest and attendant lock							
BACK UPHOLSTERY							
BH07	<input type="checkbox"/>	Slip-on Active-Mesh back upholstery					NC
BH08	<input type="checkbox"/>	Adjustable tension Active-Mesh back upholstery				WAMD/WAME	\$260
BH99	<input type="checkbox"/>	Omit Back upholstery					STD
	<input type="checkbox"/>	NXT back Support			See pricelist		
SEAT CUSHION							
SC01	<input type="checkbox"/>	2" Cushion					\$185
SC02	<input type="checkbox"/>	3" Cushion					\$205
SC99	<input type="checkbox"/>	Omit seat cushion					STD
	<input type="checkbox"/>	NXT cushion			See price list		

Rear wheels

ONE ARM DRIVE ¹		Only available with Newton One rear wheels with Aluminum, Plastic coated Handrims, permanent axles, 0° camber				Available		
HR22	<input type="checkbox"/>	One arm drive - aluminum inner handrim ^{2,3}	24"	Left <input type="checkbox"/>	Right <input type="checkbox"/>		WAML	\$1,480
HR60	<input type="checkbox"/>	One arm drive - plastic coated inner handrim ^{2,3}	24"	Left <input type="checkbox"/>	Right <input type="checkbox"/>		WAML	\$1,480
HR98	<input type="checkbox"/>	Omit handrim on opposite side of the One Arm Drive (OAD)						NC
REAR WHEELS		Available						
RW01	<input type="checkbox"/>	Mag ⁴						\$140
RW02	<input type="checkbox"/>	Newton One - spoke wheel ⁵						NC
RW03	<input type="checkbox"/>	Newton Gravity - ultralight wheel				20", 22", 24", 25", 26"		\$380
RW04	<input type="checkbox"/>	Spinergy Spox	Blk STD <input type="checkbox"/>	Wht ⁶ <input type="checkbox"/>	Blue ⁶ <input type="checkbox"/>	Red ⁶ <input type="checkbox"/>	Yel ⁶ <input type="checkbox"/>	\$850
RW05	<input type="checkbox"/>	Spinergy LX	Blk STD <input type="checkbox"/>	Pur ⁶ <input type="checkbox"/>	Wht ⁶ <input type="checkbox"/>	Blue ⁶ <input type="checkbox"/>	Red ⁶ <input type="checkbox"/>	Yel ⁶ <input type="checkbox"/>
RW07	<input type="checkbox"/>	Spinergy CLX ⁷	Blk STD <input type="checkbox"/>	Pur ⁶ <input type="checkbox"/>	Wht ⁶ <input type="checkbox"/>	Blue ⁶ <input type="checkbox"/>	Red ⁶ <input type="checkbox"/>	Yel ⁶ <input type="checkbox"/>
								\$1,055
								\$1,635
TIRES		Available						
TI01	<input type="checkbox"/>	Soft Urethane - 1 3/8"				20", 22" 24"		STD
TI02	<input type="checkbox"/>	Soft Urethane - knobby 1 3/8 "				22", 24", 25", 26"		\$120
TI04	<input type="checkbox"/>	Soft Urethane - Shox 1"				22", 24", 25", 26"		\$140
TI05	<input type="checkbox"/>	Pneumatic - 1 3/8"				20", 22", 24", 25" 26"		NC
TI08	<input type="checkbox"/>	Pneumatic with Airless insert - 1 3/8" ⁸				22", 24", 26"		\$225
TI06	<input type="checkbox"/>	Pneumatic HP - 1"				20", 22", 24", 25", 26"		\$120
TI07	<input type="checkbox"/>	Pneumatic HP - Puncture resistant, black 1"				24", 25", 26"		\$200

1- Offered with handrims on each side of the wheelchair; select the outer handrim in the handrims section; 2- Only standard handrim position available; 3- Not available with HD option
 4- Not compatible with HP tires ; Only available with aluminum or plastic coated handrim, standard position; 5- Not compatible with Shox tires; 6- Colour spokes require an additional 10 days lead-time
 7- Only available with 1" tires; 8- Only available with Mag, Newton One and Newton Gravity rear wheels

Handrims

HANDRIM		Available						
HR01	<input type="checkbox"/>	Aluminum anodized				20", 22", 24", 25", 26"		STD
HR61	<input type="checkbox"/>	Aluminum black hard anodized ²				20", 22", 24", 25"		\$80
HR02	<input type="checkbox"/>	Plastic coated ¹				20", 22", 24", 25", 26"	WAMO	\$170
HR03	<input type="checkbox"/>	High friction coated ^{1,2}				22", 24", 25", 26"	WAMO	\$460
HR07 to HR18	<input type="checkbox"/>	Newton Air Grip ²	<i>See on page 6 for color choice and price</i>			20", 22", 24", 25", 26"	WAMO	
HR28	<input type="checkbox"/>	Surge ^{1,2} (1 1/2" Oval handrim w/ Gripton high friction strip)				22", 24", 25", 26"		\$485
HR05	<input type="checkbox"/>	Surge LT ^{1,2} (1 1/8" Oval handrim w/ Gripton high friction strip)				22", 24", 25", 26"		\$515
HR04	<input type="checkbox"/>	Natural Fit ^{1,2}				22", 24", 25", 26"		\$320
HR63	<input type="checkbox"/>	Nova H ^{1,2} (1 1/4" Oval handrim with high friction strip)				22", 24", 25", 26"		\$515
HR64	<input type="checkbox"/>	Optimum H ^{1,2} (1 1/8" heart-shaped handrim)				22", 24", 25", 26"		\$320
HR99	<input type="checkbox"/>	Omit handrim						NC

1- Only standard and narrow handrim position. By default, the tab won't be cut with the narrow position; 2- Not available with Mag wheels

Handrims option

ASSEMBLY POSITION		Available						
HR29	<input type="checkbox"/>	Standard Handrim position	Full tab					STD
HR030	<input type="checkbox"/>	Narrow Handrim position no cut	Full tab					NC
HR30	<input type="checkbox"/>	Narrow Handrim position with cut ³	Short tab 2 positions					\$27
HR31	<input type="checkbox"/>	Super Narrow Handrim ³	Short tab 2 positions STD <input type="checkbox"/>			Short tab 1 position <input type="checkbox"/>		\$27
THUMB PIECE SELECTION FOR NATURAL FIT, SURGE, SURGE LT HANDRIMS		Available						
HR23	<input type="checkbox"/>	Standard grip						\$200
HR27	<input type="checkbox"/>	Super grip						\$230
HR099	<input type="checkbox"/>	Omit thumb piece						NC


3- Only available with Aluminum and Air Grip handrims;

AXLE ¹		Available						
AX01	<input type="checkbox"/>	Quick release axle					WAMU	STD
AX02	<input type="checkbox"/>	Permanent axle						NC
AX03	<input type="checkbox"/>	Quad release axle					WAMU	\$180
AXLE PLATE		Available						
AX04	<input type="checkbox"/>	Amputee adjustable axle plate ²					WAMT	\$270

1- Not required with one arm drive 2 - May interfere with swing away armrest and limit seat to floor height available

REAR WHEEL CAMBER			
CM00	<input type="checkbox"/>	0° camber	STD
CM03	<input type="checkbox"/>	3° camber	NC
CM06	<input type="checkbox"/>	6° camber (not available with amputee axle plate)	\$85

CENTER OF GRAVITY



Choose position of rear wheel axle plate:

Measure distance X from front of back post to center of axle receiver.

Most stable													
X =	1"	1/4"	1/2"	3/4"	2"	1/4"	1/2"	3/4"	3"	1/4"	1/2"	3/4"	4"
	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

If not indicated, default setting is at 1.5"

Possible conflict between a seat depth under 16" and a center of gravity over 2.5".

Possible conflict between a 26" rear wheel and a center of gravity over 2.5".

Motion Composites will select the closest possible center of gravity (COG) if an issue arises with the COG

WHEEL LOCK			
WL08	<input type="checkbox"/>	Newton Wheel lock - push to lock	STD
WL16	<input type="checkbox"/>	Aluminum push to lock with extension	\$73
WL02	<input type="checkbox"/>	Pull to lock	NC
WL04	<input type="checkbox"/>	Newton grade aid push to lock (suggested with pneumatic or high pressure tire)	WAMR \$185
WL13 WL14	<input type="checkbox"/>	6" removable extension handle	\$73
WL05	<input type="checkbox"/>	Scissors Lock	\$185
WL06	<input type="checkbox"/>	Attendant Lock (additional set of wheel lock mounted to the back)	\$240
WL15	<input type="checkbox"/>	Unilateral wheel lock - pull to lock ¹ Left <input type="checkbox"/> Right <input type="checkbox"/>	WAMN \$335
WL99	<input type="checkbox"/>	Omit wheel locks	NC

1- For more information about incompatibilities and restriction, refer to the table at page 8 and the "How-to documents" page of our website:

https://www.motioncomposites.com/en_ca/support-and-education/how-to-documents

Armrests

HEIGHT ADJUSTABLE FLIP BACK ARMREST (Convertible to "T" armrest) ^{1,4}			
AR01-812	<input type="checkbox"/>	10" Desk length armpad 8" - 12" Locked <input type="checkbox"/> Unlocked STD <input type="checkbox"/>	\$395
AR01-1014	<input type="checkbox"/>	10" Desk length armpad 10" - 14" Locked <input type="checkbox"/> Unlocked STD <input type="checkbox"/>	\$395
AR02-812	<input type="checkbox"/>	14" Full length armpad 8" - 12" Locked <input type="checkbox"/> Unlocked STD <input type="checkbox"/>	\$395
AR02-1014	<input type="checkbox"/>	14" Full length armpad 10" - 14" Locked <input type="checkbox"/> Unlocked STD <input type="checkbox"/>	\$395
HEIGHT ADJUSTABLE "T" ARMREST ¹			
AR13-610	<input type="checkbox"/>	10" Desk length armpad ⁵ 6" - 10"	\$380
AR03-812	<input type="checkbox"/>	10" Desk length armpad 8" - 12"	\$380
AR03-1014	<input type="checkbox"/>	10" Desk length armpad 10" - 14"	\$380
AR14-610	<input type="checkbox"/>	14" Full length armpad ⁵ 6" - 10"	\$380
AR04-812	<input type="checkbox"/>	14" Full length armpad 8" - 12"	\$380
AR04-1014	<input type="checkbox"/>	14" Full length armpad 10" - 14"	\$380
HEIGHT ADJUSTABLE SWING AWAY ARMREST ^{2,4}			
AR05	<input type="checkbox"/>	Tubular armrest - short pad 8" - 11"	\$315
AR23	<input type="checkbox"/>	Tubular armrest - long pad 8" - 11"	\$315
AR06	<input type="checkbox"/>	10" Desk length flat armpad 8" - 11"	\$315
AR17	<input type="checkbox"/>	14" Full length flat armpad 8" - 11"	\$315
ARMRESTS OPTIONS			
AR12	<input type="checkbox"/>	14" Full Length armpad - Gel Ovation (pair) ³	SEAND0030 \$160
AR99	<input type="checkbox"/>	Omit armrest	NC

1- May have 1/2" adjust. restriction on both limits of the range; 2- Could interfere with some rigid backrests; 3- Not available with swing away & flip back armrest; 4- Not available with attendant lock

5- May have some restriction with seat-to-floor height and rear wheel

SIDE GUARDS (not available with Flip Back and "T" Armrest)			
SG01	<input type="checkbox"/>	Plastic side guards (removable)	WAMV \$235
SG02	<input type="checkbox"/>	Carbon fiber side guards (fixed)	WAMV \$275
SG03	<input type="checkbox"/>	Carbon fiber side guards (removable)	WAMV \$315
SG06	<input type="checkbox"/>	Plastic fender side guards - straight bracket (removable)	WAMV \$340
SG08	<input type="checkbox"/>	Plastic fender side guards - offset bracket (removable) ¹	WAMV \$340
SG04	<input type="checkbox"/>	Carbon fiber fender side guards (removable)	WAMV \$510
SG99	<input type="checkbox"/>	Omit side guard	NC
SIDE GUARDS OPTION			
SG07	<input type="checkbox"/>	Side guards clip (pair) ²	NC

1- only compatible with 0° camber; Increase distance between fender side guards by 1 1/8"; 2- Not available with fixed carbon fiber side guard

POSITIONING BELT			
SB01-S	<input type="checkbox"/>	Velcro adjustable - Standard 35"	\$50
SB01-L	<input type="checkbox"/>	Velcro adjustable - Long 60"	\$50
SB02-S	<input type="checkbox"/>	Auto buckle - Standard 35"	\$60
SB02-L	<input type="checkbox"/>	Auto buckle - Long 60"	\$60
SB04-M	<input type="checkbox"/>	BodyPoint Hip Belt - 2 point not padded - Medium 56"	\$90
SB04-L	<input type="checkbox"/>	BodyPoint Hip Belt - 2 point not padded - Long 61"	\$90
SB05-S	<input type="checkbox"/>	BodyPoint Hip Belt - 2 point padded - Small 51"	SERND0030 \$180
SB05-M	<input type="checkbox"/>	BodyPoint Hip Belt - 2 point padded - Medium 56"	SERND0030 \$180
SB05-L	<input type="checkbox"/>	BodyPoint Hip Belt - 2 point padded - Long 61"	SERND0030 \$180
SB06-S	<input type="checkbox"/>	BodyPoint Hip Belt - 4 point padded - Small 51"	SERND0020 \$205
SB06-M	<input type="checkbox"/>	BodyPoint Hip Belt - 4 point padded - Medium 56"	SERND0020 \$205
SB06-L	<input type="checkbox"/>	BodyPoint Hip Belt - 4 point padded - Long 61"	SERND0020 \$205
SB03	<input type="checkbox"/>	Hardware attachment for 2 extra points (1 pair)	\$15
SB08-M	<input type="checkbox"/>	BodyPoint EVOFLEX® Positioning Belt Rehab latch and Band Clamps - Medium 18"	\$345
SB08-L	<input type="checkbox"/>	BodyPoint EVOFLEX® Positioning Belt Rehab latch and Band Clamps - Long 24"	\$345
SB99	<input type="checkbox"/>	Omit positioning belt	STD

HEADREST ¹			
AC01-S	<input type="checkbox"/>	Profiled Headrest - Small (3.5"X6")	\$215
AC01-M	<input type="checkbox"/>	Profiled Headrest - Medium (4.5"X7")	\$215
AC01-L	<input type="checkbox"/>	Profiled Headrest - Large (5.5"X8")	\$215
AC28	<input type="checkbox"/>	Headrest mounting hardware with removable non folding stabilizer bar ¹	\$360

1- Only available with push handles - Not available with a stroller handle or a folding stabilizer bar

ACCESSORIES			
AC03	<input type="checkbox"/>	Anti-tippers	\$120
AC29	<input type="checkbox"/>	Swing-away anti-tipper ¹ Left <input type="checkbox"/> Right <input type="checkbox"/>	\$180
AC27	<input type="checkbox"/>	Tip assist ²	\$65
AC06	<input type="checkbox"/>	Cane and crutch holder Left <input type="checkbox"/> Right <input type="checkbox"/>	\$115
AC07	<input type="checkbox"/>	Oxygen tank holder ⁴ - D and E Type - installed on the left	\$695
AC08	<input type="checkbox"/>	Spoke guards ³ Available: 20", 22", 24", 26"	\$155
AC09	<input type="checkbox"/>	Footrest neoprene impact guard (pair) - Black	\$85
AC10	<input type="checkbox"/>	Seat pouch	\$100
AC30	<input type="checkbox"/>	Newton multi-purpose pouch	\$75
AC26	<input type="checkbox"/>	BodyPoint Mobility Bag	\$145
AC23	<input type="checkbox"/>	Newton Backpack	\$130
AC12	<input type="checkbox"/>	Velcro Closing Strap	\$55
AC13	<input type="checkbox"/>	Tool kit	\$55
AC16	<input type="checkbox"/>	Motion Composites Polo Shirt W/S <input type="checkbox"/> W/M <input type="checkbox"/> W/L <input type="checkbox"/> M/S <input type="checkbox"/> M/M <input type="checkbox"/> M/L <input type="checkbox"/> M/XL <input type="checkbox"/> M/XXL <input type="checkbox"/> M/XXXL <input type="checkbox"/>	\$36

1- Must select the pair (left and right) with HD option² Only available with rear seat to floor height between 13 1/2" and 19 3/4"

2 - Cannot be combined with Anti Tippers

3 - Not available with Mag wheels and Spinergy; 4- Not compatible with swing away armrest or amputee axle plate and WC-19 transit tie-down


Colors																					
Options ¹	Black	High Gloss Carbon Fiber	Sunkissed Orange	Burgundy	Sapphire Blue	Monster Green	Charcoal	Fuchsia	Steel Blue	Ferrari Red	Acid Green	Electric Purple	Caribbean Blue	White	Pink Pearl	Plum Purple	Sky Blue	Black Cherry	Khaki Green (matte)	Desert Sand (matte)	
AirGrip	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	\$350
Frame ²	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	NC
FC99	<input type="checkbox"/>	Custom Color Specify PPG Paintit color: additional 10 days lead-time																	\$422		

1- Color representation may vary from one monitor to another. Therefore, the actual paint color on your wheelchair might differ slightly from what you see on your screen.;

Finish between Airgrip and frame might be different;

2- With black carbon fiber crossbrace;

Fork and Stem bolt combinations




Unilateral wheel lock restrictions


Seat depth: shall be 16 in and greater		
Rear wheel: shall be 24 in and less		
"T" or flip-back armrest only compatible with 0° Camber		
One arm drive: Not compatible with "T" or flip-back armrest		
"T" or flip-back armrest may require to move the rear wheels outwards		
Not compatible with wheel lock lever extension		
Rear wheel diameter	Seat depth (in)	COG Compatibility Range (in)
20 in	16+	1.5 to 3.75
	16	1.0 to 2.75
	17	1.0 to 3.25
22 in	18+	1.0 to 3.75
	16	1.0 to 2.0
24 in	17	1.0 to 2.75
	18+	1.0 to 3.75
	16	1.0 to 2.0

Notes :

FOLDING ON CHAIR - CDN ORDER FORM

BioFit™ CUSHION	WIDTH											
	12"	13"	14"	15"	16"	17"	18"	19"	20"	21"	22"	
 NXCABF _____ (W" x D")	<input type="checkbox"/> \$625	<input type="checkbox"/> \$625	<input type="checkbox"/> \$625									
	<input type="checkbox"/> \$625	<input type="checkbox"/> \$625	<input type="checkbox"/> \$625									
		<input type="checkbox"/> \$625	<input type="checkbox"/> \$625	<input type="checkbox"/> \$625	<input type="checkbox"/> \$625	<input type="checkbox"/> \$625						
		<input type="checkbox"/> \$625	<input type="checkbox"/> \$625	<input type="checkbox"/> \$625	<input type="checkbox"/> \$625	<input type="checkbox"/> \$625	<input type="checkbox"/> \$625	<input type="checkbox"/> \$625	<input type="checkbox"/> \$655	<input type="checkbox"/> \$655	<input type="checkbox"/> \$655	
				<input type="checkbox"/> \$625	<input type="checkbox"/> \$625	<input type="checkbox"/> \$625	<input type="checkbox"/> \$625	<input type="checkbox"/> \$655	<input type="checkbox"/> \$655			
				<input type="checkbox"/> \$625	<input type="checkbox"/> \$625	<input type="checkbox"/> \$625	<input type="checkbox"/> \$625	<input type="checkbox"/> \$655	<input type="checkbox"/> \$655	<input type="checkbox"/> \$655	<input type="checkbox"/> \$655	
				<input type="checkbox"/> \$625	<input type="checkbox"/> \$625	<input type="checkbox"/> \$625	<input type="checkbox"/> \$625	<input type="checkbox"/> \$655	<input type="checkbox"/> \$655			
				<input type="checkbox"/> \$625	<input type="checkbox"/> \$625	<input type="checkbox"/> \$625	<input type="checkbox"/> \$625	<input type="checkbox"/> \$655	<input type="checkbox"/> \$655	<input type="checkbox"/> \$655	<input type="checkbox"/> \$655	
Option: Change cover for incontinent / infection control outer cover											<input type="checkbox"/> N/C	


SEATING AND CUSHION ACCESSORIES

	Model #	Description	Fits Back Width	\$	QTY
	RSSI 1600	Solid Seat Insert	Fits cushions 16" & 17" wide (40.6 & 43.2 cm) and 14"-18" deep (35.6-45.7 cm)	80	
	RSSI 1800	Solid Seat Insert	Fits cushions 18" & 19" wide (45.7 & 48.3 cm) and 14"-18" deep (35.6-45.7 cm)	80	
	RSSI 2000	Solid Seat Insert	Fits cushions 20"-22" wide (50.8-55.9 cm) and 18"-20" deep (45.7-50.8 cm)	80	


Optima™ THORACIC BACK SUPPORT

	Optima™	Fits Width	Height	✓	\$
	NOACB1613	16"-18"	13"	<input type="checkbox"/>	\$725
	NOACB1616	16"-18"	16"	<input type="checkbox"/>	\$725
	NOACB1619	16"-18"	19"	<input type="checkbox"/>	\$725
	NOACB1813	18"-20"	13"	<input type="checkbox"/>	\$725
	NOACB1816	18"-20"	16"	<input type="checkbox"/>	\$725
	NOACB1819	18"-20"	19"	<input type="checkbox"/>	\$725
	NOACB2016	20"-22"	16"	<input type="checkbox"/>	\$725
	NOACB2019	20"-22"	19"	<input type="checkbox"/>	\$725
	NOACB2213	22"-24"	13"	<input type="checkbox"/>	\$755
	NOACB2216	22"-24"	16"	<input type="checkbox"/>	\$755
	NOACB2219	22"-24"	19"	<input type="checkbox"/>	\$755


Optima™ DEEP THORACIC BACK SUPPORT

	Optima™ DEEP	Fits Width	Height	✓	\$
	NOPCB1616	16"-18"	16"	<input type="checkbox"/>	\$825
	NOPCB1619	16"-18"	19"	<input type="checkbox"/>	\$825
	NOPCB1816	18"-20"	16"	<input type="checkbox"/>	\$825
	NOPCB1819	18"-20"	19"	<input type="checkbox"/>	\$825
	NOPCB2016	20"-22"	16"	<input type="checkbox"/>	\$825
	NOPCB2019	20"-22"	19"	<input type="checkbox"/>	\$825
	NOPCB2213	22"-24"	13"	<input type="checkbox"/>	\$860
	NOPCB2216	22"-24"	16"	<input type="checkbox"/>	\$860
	NOPCB2219	22"-24"	19"	<input type="checkbox"/>	\$860

Optima™ CARBON THORACIC BACK SUPPORT

	Optima™ Carbon	Fits Width	Height	✓	\$
	NOCFACB1413	14"-16"	13"	<input type="checkbox"/>	\$1150
	NOCFACB1416	14"-16"	16"	<input type="checkbox"/>	\$1150
	NOCFACB1513	15"-17"	13"	<input type="checkbox"/>	\$1150
	NOCFACB1516	15"-17"	16"	<input type="checkbox"/>	\$1150
	NOCFACB1613	16"-18"	13"	<input type="checkbox"/>	\$1150
	NOCFACB1616	16"-18"	16"	<input type="checkbox"/>	\$1150
	NOCFACB1713	17"-19"	13"	<input type="checkbox"/>	\$1150
	NOCFACB1716	17"-19"	16"	<input type="checkbox"/>	\$1150
	NOCFACB1813	18"-20"	13"	<input type="checkbox"/>	\$1150
	NOCFACB1816	18"-20"	16"	<input type="checkbox"/>	\$1150
	NOCFACB2013	20"-22"	13"	<input type="checkbox"/>	\$1150
	NOCFACB2016	20"-22"	16"	<input type="checkbox"/>	\$1150

Optima™ CARBON DEEP THORACIC BACK SUPPORT

	Optima™ Carbon DEEP	Fits Width	Height	✓	\$
	NOCFPCB1413	14"-16"	13"	<input type="checkbox"/>	\$1250
	NOCFPCB1416	14"-16"	16"	<input type="checkbox"/>	\$1250
	NOCFPCB1513	15"-17"	13"	<input type="checkbox"/>	\$1250
	NOCFPCB1516	15"-17"	16"	<input type="checkbox"/>	\$1250
	NOCFPCB1613	16"-18"	13"	<input type="checkbox"/>	\$1250
	NOCFPCB1616	16"-18"	16"	<input type="checkbox"/>	\$1250
	NOCFPCB1713	17"-19"	13"	<input type="checkbox"/>	\$1250
	NOCFPCB1716	17"-19"	16"	<input type="checkbox"/>	\$1250
	NOCFPCB1813	18"-20"	13"	<input type="checkbox"/>	\$1250
	NOCFPCB1816	18"-20"	16"	<input type="checkbox"/>	\$1250
	NOCFPCB2013	20"-22"	13"	<input type="checkbox"/>	\$1250
	NOCFPCB2016	20"-22"	16"	<input type="checkbox"/>	\$1250

THORACIC SHAPE

FOLDING ON CHAIR - CDN ORDER FORM

Xtend™ LOW HEIGHT ADJUSTABLE THORACIC BACK SUPPORT

Xtend™ DEEP HEIGHT ADJUSTABLE THORACIC BACK SUPPORT

THORACIC SHAPE		Xtend™	Fits Width	Height	✓	\$
		**NHAACB (L)1213/17	12"-14"	13"-17"	<input type="checkbox"/>	\$725
		NHAACB (L)1413/17	14"-16"	13"-17"	<input type="checkbox"/>	\$725
		NHAACB (L)1613/17	16"-18"	13"-17"	<input type="checkbox"/>	\$725
		NHAACB (L)1813/17	18"-20"	13"-17"	<input type="checkbox"/>	\$755
		NHAACB (L)2013/17	20"-22"	13"-17"	<input type="checkbox"/>	\$755
		NHAACB (L)2213/17	22"-24"	13"-17"	<input type="checkbox"/>	\$755

	Xtend™ DEEP	Fits Width	Height	✓	\$
	NHAPCB1616/20	16"-18"	16"-20"	<input type="checkbox"/>	\$725
	NHAPCB1816/20	18"-20"	16"-20"	<input type="checkbox"/>	\$725
	NHAPCB2016/20	20"-22"	16"-20"	<input type="checkbox"/>	\$725
	NHAPCB2216/20	22"-24"	16"-20"	<input type="checkbox"/>	\$755

Classic™ PELVIC BACK SUPPORT

Classic™ DEEP PELVIC BACK SUPPORT

PELVIC SHAPE		Classic™	Fits Width	Height	✓	\$
		NACB1616	16"-18"	16"	<input type="checkbox"/>	\$720
		NACB1618	16"-18"	18"	<input type="checkbox"/>	\$720
		NACB1816	18"-20"	16"	<input type="checkbox"/>	\$720
		NACB1818	18"-20"	18"	<input type="checkbox"/>	\$720
		NACB1820	18"-20"	20"	<input type="checkbox"/>	\$720
		NACB2020	20"-22"	20"	<input type="checkbox"/>	\$720
		NACB2216	22"-24"	16"	<input type="checkbox"/>	\$755
		NACB2218	22"-24"	18"	<input type="checkbox"/>	\$755
		NACB2220	22"-24"	20"	<input type="checkbox"/>	\$755

	Classic DEEP™	Fits Width	Height	✓	\$
	NPCB1616	16"-18"	16"	<input type="checkbox"/>	\$825
	NPCB1618	16"-18"	18"	<input type="checkbox"/>	\$825
	NPCB1818	18"-20"	18"	<input type="checkbox"/>	\$825
	NPCB1820	18"-20"	20"	<input type="checkbox"/>	\$825
	NPCB2020	20"-22"	20"	<input type="checkbox"/>	\$825
	NPCB2216	22"-24"	16"	<input type="checkbox"/>	\$860
	NPCB2218	22"-24"	18"	<input type="checkbox"/>	\$860
	NPCB2220	22"-24"	20"	<input type="checkbox"/>	\$860

BACK SUPPORT ACCESSORIES - PRIVACY FLAP

	Model #	Description	Fits Back Width	\$	QTY
	NIB-1938	Extra Small Privacy Flap for 11" - 13" wide back supports	11" - 13" (27.9 - 33.0)	70	
	NIB-1078	Small Privacy Flap for 14"-16" wide back supports	14" - 16" (35.6-40.6 cm)	70	
	NIB-1079	Medium Flap for 17"-19" wide back supports	17" - 19" (43.2-48.3 cm)	70	
	NIB-1080	Large Privacy Flap for 20"-22" wide back supports	20" - 22" (50.8-55.9 cm)	70	
	NIB-1081	Extra Large Privacy Flap for 23" - 26" wide back supports	23" - 26" (58.4 - 66.0 cm)	70	

MultiFit™ HEAD SUPPORT - MULTI-ADJUSTABLE PADS

MultiFit™ HEAD SUPPORT - FIXTURES

	Model # NIB-0784	✓	\$
10" x 5" Large Rectangular	<input type="checkbox"/>	\$205	

	Model # NIB-0818	✓	\$
14" x 5" Extra-Large Rectangular	<input type="checkbox"/>	\$205	

	Model # NIB-0747	✓	\$
Multi-Axis Hardware	<input type="checkbox"/>	\$150	

	Model # NIB-0823	✓	\$
Head Support Receiver Bracket	<input type="checkbox"/>	\$60	

ALL PRICES QUOTED IN CDN. PRICES ARE SUBJECT TO CHANGE WITHOUT NOTICE.

For additional NXT models, sizes and hardware not shown, please refer to the complete NXT Seating Order Form at www.motioncomposites.com