



MOVE
ORDER FORM - CANADA

1- Purchasing information

P.O. number:* _____ Date: _____
 Company Name:* _____
 Purchasing contact:* _____
 Phone:* _____
 Email:* _____
 Contact preference:* Phone Email
 Marked for: _____

3- Shipping information

Ship to:* _____
 Name: _____
 Address: _____
 City: _____ Province: _____
 Postal Code: _____
 Phone: _____

2- Dealer/Representative information

Dealer/Representative:* _____
 Phone:* _____
 Email:* _____
 Contact preference:* Phone Email

MA080 English Manual (STD) MA081 Manuel Français Taxes not included

MODEL							
WC01	<input type="checkbox"/>	Ultralight MOVE wheelchair aluminum. (265 lb weight capacity)			WM0001196	\$1 425	
WC09	<input type="checkbox"/>	Heavy duty kit (350 lb weight capacity)			WAMF	\$800	
SEAT WIDTH ³							
SW	<input type="checkbox"/>	14" <input type="checkbox"/>	15" <input type="checkbox"/>	16" <input type="checkbox"/>	17" <input type="checkbox"/>	18" ¹ <input type="checkbox"/>	NC
	<input type="checkbox"/>	19" ¹ <input type="checkbox"/>	20" ¹ <input type="checkbox"/>	22" ² <input type="checkbox"/>			WAMA \$305
SEAT DEPTH							
SD	<input type="checkbox"/>	16" <input type="checkbox"/>			18" <input type="checkbox"/>	NC	
	<input type="checkbox"/>	20" <input type="checkbox"/>					WAMB \$230

1 - Available with HD option 2 - Only available with HD option [3. Width calculator](#)

Move FRONT Seat-to-floor Height - Please select your choice

Fork/Position		4" caster			5" caster			6" caster			7" caster			8" caster		
		Stem bolt length			Stem bolt length			Stem bolt length			Stem bolt length			Stem bolt length		
		STD	+1"	+2"	STD	+1"	+2"	STD	+1"	+2"	STD	+1"	+2"	STD	+1"	+2"
4"	P1	12 1/2"	13 1/2"	14 1/2"												
	P2	13"	14"	15"	13 1/2"	14 1/2"	15 1/2"									
5"	P1															
	P2	13"	14"	15"	13 1/2"	14 1/2"	15 1/2"									
	P3	13 1/2"	14 1/2"	15 1/2"	14"	15"	16"	14 1/2"	15 1/2"	16 1/2"						
	P4	14"	15"	16"	14 1/2"	15 1/2"	16 1/2"	15"	16"	17"	15 1/2"	16 1/2"	17 1/2"			
7"	P1							15 1/2"	16 1/2"	17 1/2"	16"	17"	18"			
	P2				15 1/2"	16 1/2"	17 1/2"	16"	17"	18"	16 1/2"	17 1/2"	18 1/2"			
	P3				16"	17"	18"	16 1/2"	17 1/2"	18 1/2"	17"	18"	19"	17 1/2"	18 1/2"	19 1/2"
	P4	16"	17"	18"	16 1/2"	17 1/2"	18 1/2"	17"	18"	19"	17 1/2"	18 1/2"	19 1/2"	18"	19"	20"
	P5	16 1/2"	17 1/2"	18 1/2"	17"	18"	19"	17 1/2"	18 1/2"	19 1/2"	18"	19"	20"	18 1/2"	19 1/2"	20 1/2"

Height available +/- 1/4" (see stem bolt and fork combination at the end of the order form)

CASTERS						
CA01-04	<input type="checkbox"/>	Composite wheel w/ PU tire	4" x 1"			NC
CA01-05	<input type="checkbox"/>		5" x 1"			NC
CA01-06	<input type="checkbox"/>		6" x 1"			NC
CA01-07	<input type="checkbox"/>		7" x 1"			NC
CA01-08	<input type="checkbox"/>		8" x 1"			NC
CA10	<input type="checkbox"/>	Pneumatic casters	6" x 1 1/4"			\$165
CA11	<input type="checkbox"/>		8" x 1 1/4"			\$165
CA14	<input type="checkbox"/>	Newton UltraCaster Composite	4" x 1 1/2"			\$98
CA15	<input type="checkbox"/>		5" x 1 1/2"			\$98
CA16	<input type="checkbox"/>		6" x 1 1/2"			\$98
CA17	<input type="checkbox"/>		4" x 1 1/2"			\$180
CA18	<input type="checkbox"/>	Newton UltraCaster Aluminum	5" x 1 1/2"			\$190
CA19	<input type="checkbox"/>		6" x 1 1/2"			\$200

Move REAR Seat-to-floor Height - Please circle your choice

20" (451)	22" (501)	24" (540)	25" (559)	26" (590)
12"	13"	14"	14 1/2"	15"
13"	14"	15"	15 1/2"	16"
14"	15"	16"	16 1/2"	17"
15"	16"	17"	17 1/2"	18"
16"	17"	18"	18 1/2"	19"
17"	18"	19"	19 1/2"	20"

Note: Tire size, center of gravity setting and seat slope may result in a variation of +/- 1/2" of selected seat to floor height"

Footrest						
SWING-AWAY FOOTREST						
FR01-01	<input type="checkbox"/>	60° swing away footrest	13.5 - 17.5"		Set at:*	\$260
FR01-02	<input type="checkbox"/>	60° swing away footrest	15.5 - 19.5"		Set at:*	\$260
FR01-03	<input type="checkbox"/>	60° swing away footrest	17 - 21"		Set at:*	\$260
FR02-01	<input type="checkbox"/>	70° swing away footrest	13.5 - 17.5"		Set at:*	\$260
FR02-02	<input type="checkbox"/>	70° swing away footrest	15.5 - 19.5"		Set at:*	\$260
FR02-03	<input type="checkbox"/>	70° swing away footrest	17 - 21"		Set at:*	\$260
FOOTREST OPTIONS						
FR09	<input type="checkbox"/>	High mount (footplate attached to the hanger with clamps) ³	5" - 13.5"		Set at:*	\$95
FR10	<input type="checkbox"/>	Elevating legrest with calf pad, Length adjustable ¹	Left <input type="checkbox"/> Right <input type="checkbox"/>	14" - 22"	Set at:	WAMK \$390
FR99	<input type="checkbox"/>	Omit Footrest				NC

¹ - The elevating legrest replace the footrest side selected

* +/- 1/2"

FOOTPLATES						
FP01	<input type="checkbox"/>	Standard - Composite				STD
FP02	<input type="checkbox"/>	Newton adjustable angle - Composite			WAMJ	\$150
FP03	<input type="checkbox"/>	Newton adjustable angle - Aluminum			WAMJ	\$170
FOOTPLATE ACCESSORIES						
FP08	<input type="checkbox"/>	Length adjustable heel loop				\$25
FP09	<input type="checkbox"/>	Calf strap			SEFND0010	\$55
FP30	<input type="checkbox"/>	Bodypoint Aeromesh Padded calf strap			SEFND0010	\$135

SWING AWAY RESIDUAL LIMB SUPPORT						
FR27	<input type="checkbox"/>	Residual limb support - flat 6"x8"	Left <input type="checkbox"/>	Right <input type="checkbox"/>		SEFND0005 \$490
FR28	<input type="checkbox"/>	Residual limb support - flat 8"x10"	Left <input type="checkbox"/>	Right <input type="checkbox"/>		SEFND0005 \$490
FR29	<input type="checkbox"/>	Residual limb support - flat 10"x10"	Left <input type="checkbox"/>	Right <input type="checkbox"/>		SEFND0005 \$490
FR30	<input type="checkbox"/>	Residual limb support - flat 10"x14"	Left <input type="checkbox"/>	Right <input type="checkbox"/>		SEFND0005 \$490
FR32	<input type="checkbox"/>	Residual limb support - contour 6"x8"	Left <input type="checkbox"/>	Right <input type="checkbox"/>		SEFND0005 \$490
FR33	<input type="checkbox"/>	Residual limb support - contour 8"x10"	Left <input type="checkbox"/>	Right <input type="checkbox"/>		SEFND0005 \$490

Backrest

FIXED ANGLE BACK CANES WITH INTEGRATED PUSH HANDLES						
BC01-16	<input type="checkbox"/>	8° bend aluminum back cane	16"			NC
BC01-18	<input type="checkbox"/>	8° bend aluminum back cane	18"			NC
BC01-20	<input type="checkbox"/>	8° bend aluminum back cane	20"			NC
BC59	<input type="checkbox"/>	8° bend aluminum back cane	21" - 24"	increment of 1"	Set at ² :	New \$95
ADJUSTABLE ANGLE BACK WITH INTEGRATED PUSH HANDLES						
BC02-16	<input type="checkbox"/>	8° bend aluminum back cane	16"			\$320
BC02-20	<input type="checkbox"/>	8° bend aluminum back cane	20"			\$320
Select back cane angle:		90° <input type="checkbox"/>	95° STD <input type="checkbox"/>	100° <input type="checkbox"/>	105° <input type="checkbox"/>	110° <input type="checkbox"/>

STROLLER HANDLE AND STABILIZER BAR ¹						
PH04	<input type="checkbox"/>	Stroller handle (not available with the HD option)			WAMW	\$215
PH03	<input type="checkbox"/>	Newton folding stabilizer bar (not available with stroller handle)				\$265

1 - Stroller handle and stabilizer bar can't be combined 2- If not indicated, will be at the lowest position

BACK UPHOLSTERY						
BH05	<input type="checkbox"/>	Slip-on Nylon back upholstery				NC
BH06	<input type="checkbox"/>	Adjustable tension Nylon back upholstery			WAMD/WAME	250 \$
BH99	<input type="checkbox"/>	Omit back upholstery				STD
	<input type="checkbox"/>	NXT back Support		See pricelist		

SEAT CUSHIONS						
SC01	<input type="checkbox"/>	2" Cushion				180 \$
SC02	<input type="checkbox"/>	3" Cushion				200 \$
SC99	<input type="checkbox"/>	Omit seat cushion				NC
	<input type="checkbox"/>	NXT cushion		See pricelist		

Rear wheels

ONE ARM DRIVE ¹ Only available with Newton One rear wheels with Aluminum, Plastic coated Handrims						
HR22	<input type="checkbox"/>	One arm drive ^{2,3} 24" Left <input type="checkbox"/>	Right <input type="checkbox"/>			\$1435
HR98	<input type="checkbox"/>	Omit handrim on opposite side of the One Arm Drive (OAD)				NC

REAR WHEELS				Available		
RW01	<input type="checkbox"/>	Mag ⁴		20", 22" 24"		NC
RW02	<input type="checkbox"/>	Newton One - spoke wheel ⁵		20", 22" 24"		NC
RW03	<input type="checkbox"/>	Newton Gravity - ultralight wheel		20", 22", 24", 25", 26"		\$370

TIRES				Available		
TI01	<input type="checkbox"/>	Soft Urethane - 1 3/8"		20", 22" 24"		STD
TI02	<input type="checkbox"/>	Soft Urethane - knobby 1 3/8 "		22", 24", 25", 26"		\$135
TI04	<input type="checkbox"/>	Soft Urethane - Shox ¹ "		22", 24", 25", 26"		\$160
TI05	<input type="checkbox"/>	Pneumatic - 1 3/8"		20", 22", 24", 26"		\$135
TI08	<input type="checkbox"/>	Pneumatic with Airless insert - 1 3/8" ⁷		22", 24", 26"		\$220
TI06	<input type="checkbox"/>	Pneumatic HP - 1"		20", 22", 24", 25", 26"		\$155
TI07	<input type="checkbox"/>	Pneumatic HP - Puncture resistant, black 1"		24", 25", 26"		\$235

1- Offered with handrims on each side of the wheelchair 2- Only standard handrim position available 3- Not available with HD option
4 - Not compatible with HP tires ; Only available with aluminum or plastic coated handrim in standard position 5- Not compatible with Shox tires

Handrims

HANDRIM				Available		
HR01	<input type="checkbox"/>	Aluminum anodized		20", 22", 24", 25", 26"		STD
HR61	<input type="checkbox"/>	Aluminum black hard anodized ²		20", 22", 24", 25"		\$77
HR02	<input type="checkbox"/>	Plastic coated ¹		20", 22", 24", 25", 26"	WAMO	\$165
HR03	<input type="checkbox"/>	Q-Grip High friction Neoprene coated ^{1,2}		22", 24", 25", 26"	WAMO	\$445
HR07	<input type="checkbox"/>	Newton Air Grip ²	Black <input type="checkbox"/>	20", 22", 24", 25", 26"	WAMO	\$155
HR	<input type="checkbox"/>	Newton Air Grip ²	(08) Sunkissed Orange <input type="checkbox"/> (09) Burgundy <input type="checkbox"/> (10) Sapphire Blue <input type="checkbox"/> (11) Monster Green <input type="checkbox"/> (12) Charcoal <input type="checkbox"/> (13) Fuchsia <input type="checkbox"/> (14) Steel Blue <input type="checkbox"/> (15) Ferrari Red <input type="checkbox"/> (16) Acid Green <input type="checkbox"/> (17) White <input type="checkbox"/> (18) Sky Blue <input type="checkbox"/>	20", 22", 24", 25", 26"	WAMO	\$340
HR28	<input type="checkbox"/>	Surge ^{1,2} (1-1/2" Oval handrim w/ Gripton high friction strip)		**Limited availability	22", 24", 25", 26"	\$485
HR05	<input type="checkbox"/>	Surge LT ^{1,2} (1-1/8" Oval handrim w/ Gripton high friction strip)		**Limited availability	22", 24", 25", 26"	\$515
HR04	<input type="checkbox"/>	Natural Fit ^{1,2}		**Limited availability	22", 24", 25", 26"	\$320
HR99	<input type="checkbox"/>	Omit handrim				NC

1- Only standard handrim position available 2- Not available with Mag wheels or OAD

Handrims option

Assembly position for aluminum and Newton AirGrip handrims						
HR29	<input type="checkbox"/>	Standard Handrim position	Full tab			STD
HR30	<input type="checkbox"/>	Narrow Handrim position	Short tab 2 positions			\$26
HR31	<input type="checkbox"/>	Super Narrow Handrim	Short tab 2 positions STD <input type="checkbox"/>	Short tab 1 position <input type="checkbox"/>		\$26
Thumb piece ³ selection Natural Fit handrims						
HR23	<input type="checkbox"/>	Standard grip		**Limited availability		\$195
HR27	<input type="checkbox"/>	Super grip		**Limited availability		\$225
HR099	<input type="checkbox"/>	Omit thumb piece				NC


3- Compatible with Surge, Surge LT handrims

AXLES (0° camber) ¹			
AX01	<input type="checkbox"/>	Quick release axle	WAMU STD
AX02	<input type="checkbox"/>	Permanent axle	NC
AX03	<input type="checkbox"/>	Quad release axle	WAMU \$175
AXLE PLATE			
AX04	<input type="checkbox"/>	Amputee adjustable axle plate ²	WAMT \$260

1- Not required with one arm drive

2- May limit seat to floor height available

CENTER OF GRAVITY



Choose position of rear wheel axle plate:

Measure distance from front of back post to center of axle receiver.

Default setting is at 3/4"

Most stable		Most responsive
AMPUTEE ¹	3/4"	2"
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

1- Only available with amputee axle plate

WHEEL LOCKS			
WL08	<input type="checkbox"/>	Newton Wheel lock - push to lock	STD
WL02	<input type="checkbox"/>	Pull to lock	NC
WL04	<input type="checkbox"/>	Newton grade aid push to lock (suggested with pneumatic or high pressure tire)	WAMR \$180
WL13 WL14	<input type="checkbox"/>	6" removable extension handle	\$70
WL06	<input type="checkbox"/>	Attendant Lock (additional set of wheel lock mounted to the back)	\$235
WL15	<input type="checkbox"/>	Unilateral wheel lock - pull to lock ¹ Left <input type="checkbox"/> Right <input type="checkbox"/>	New \$325
WL99	<input type="checkbox"/>	Omit wheel locks	NC

1- Only available with: seat depths of ≥ 16"; rear wheels of 20", 22" and 24"; Not available with the following combinations: T or U armrests with One Arm Drive; May interfere with some COG and other configurations. For more details, refer to the table at page 6 and the "How-to documents" page of our website:

https://www.motioncomposites.com/en_ca/support-and-education/how-to-documents

Armrests

HEIGHT ADJUSTABLE FLIP BACK ARMREST (Convertible to "T" armrest) (improved) ^{1,2}			
AR01-812	<input type="checkbox"/>	10" Desk length armpad 8" - 12" Locked <input type="checkbox"/> Unlocked STD <input type="checkbox"/>	\$390
AR01-1014	<input type="checkbox"/>	10" Desk length armpad 10" - 14" Locked <input type="checkbox"/> Unlocked STD <input type="checkbox"/>	\$390
AR02-812	<input type="checkbox"/>	14" Full length armpad 8" - 12" Locked <input type="checkbox"/> Unlocked STD <input type="checkbox"/>	\$390
AR02-1014	<input type="checkbox"/>	14" Full length armpad 10" - 14" Locked <input type="checkbox"/> Unlocked STD <input type="checkbox"/>	\$390
HEIGHT ADJUSTABLE "T" ARMREST (improved) ¹			
AR13-610	<input type="checkbox"/>	10" Desk length armpad ⁵ 6" - 10" New	\$370
AR03-812	<input type="checkbox"/>	10" Desk length armpad 8" - 12" New	\$370
AR03-1014	<input type="checkbox"/>	10" Desk length armpad 10" - 14" New	\$370
AR14-610	<input type="checkbox"/>	14" Full length armpad ⁵ 6" - 10" New	\$370
AR04-812	<input type="checkbox"/>	14" Full length armpad 8" - 12" New	\$370
AR04-1014	<input type="checkbox"/>	14" Full length armpad 10" - 14" New	\$370
ARMREST OPTIONS			
AR12	<input type="checkbox"/>	14" Full Length armpad - Gel Ovation (pair) ³	SEAND0030 \$155
AR99	<input type="checkbox"/>	Omit armrest	NC

1- May have 1/2" adjust. restriction on both limits of the range; 2- Not available with attendant lock; 3- Not available with flip back armrest; 4- May have restriction with seat-to-floor height and rear wheel

POSITIONING BELTS						
SB01-S	<input type="checkbox"/>	Velcro adjustable - Standard 35"				\$47
SB01-L	<input type="checkbox"/>	Velcro adjustable - Long 60"				\$47
SB02-S	<input type="checkbox"/>	Auto buckle - Standard 35"			Improved	\$57
SB02-L	<input type="checkbox"/>	Auto buckle - Long 60"			Improved	\$57
SB07	<input type="checkbox"/>	Hardware attachment for 2 point seatbelt (1 pair)				NC
SB03	<input type="checkbox"/>	Additional hardware for 2 extra attachment points (needed for the 4 point seat belt) (1 pair)				\$14
SB04-M	<input type="checkbox"/>	BodyPoint Hip Belt - 2 point not padded - Medium 56"				\$85
SB04-L	<input type="checkbox"/>	BodyPoint Hip Belt - 2 point not padded - Long 61"				\$85
SB05-S	<input type="checkbox"/>	BodyPoint Hip Belt - 2 point padded - Small 51"			SERND0030	\$175
SB05-M	<input type="checkbox"/>	BodyPoint Hip Belt - 2 point padded - Medium 56"			SERND0030	\$175
SB05-L	<input type="checkbox"/>	BodyPoint Hip Belt - 2 point padded - Long 61"			SERND0030	\$175
SB06-S	<input type="checkbox"/>	BodyPoint Hip Belt - 4 point padded - Small 51"			SERND0020	\$200
SB06-M	<input type="checkbox"/>	BodyPoint Hip Belt - 4 point padded - Medium 56"			SERND0020	\$200
SB06-L	<input type="checkbox"/>	BodyPoint Hip Belt - 4 point padded - Long 61"			SERND0020	\$200
SB08-M	<input type="checkbox"/>	BodyPoint EVOFLEX® Positioning Belt Rehab latch and Band Clamps - Medium 18"				\$335
SB08-L	<input type="checkbox"/>	BodyPoint EVOFLEX® Positioning Belt Rehab latch and Band Clamps - Long 24"				\$335
SB99	<input type="checkbox"/>	Omit positioning belt				STD

HEADREST ¹						
AC01-S	<input type="checkbox"/>	Profiled Headrest - Small (3.5"X6")				\$210
AC01-M	<input type="checkbox"/>	Profiled Headrest - Medium (4.5"X7")				\$210
AC01-L	<input type="checkbox"/>	Profiled Headrest - Large (5.5"X8")				\$210
AC28	<input type="checkbox"/>	Headrest mounting hardware with removable non folding stabilizer bar ¹				\$350

1- Not available with a stroller handle or a folding stabilizer bar

ACCESSORIES						
AC03	<input type="checkbox"/>	Anti-tippers				\$120
AC29	<input type="checkbox"/>	Swing-away anti-tipper ¹	Left <input type="checkbox"/>	Right <input type="checkbox"/>		\$170
AC27	<input type="checkbox"/>	Tip assist ²				\$62
AC04	<input type="checkbox"/>	Transit Tie Down (Unoccupied)				NC
AC05-S	<input type="checkbox"/>	WC-19 Transit Tie Down w/Q'straint belt ³	Short belt 20"-41"			\$355
AC05-L	<input type="checkbox"/>	WC-19 Transit Tie Down w/Q'straint belt ³	Long belt 20"-46"			\$355
AC35	<input type="checkbox"/>	WC-19 Transit Tie Down without belt ³			New	\$325
AC06	<input type="checkbox"/>	Cane and crutch holder	Left <input type="checkbox"/>	Right <input type="checkbox"/>		\$110
AC07	<input type="checkbox"/>	Oxygen tank holder - D and E Type Not compatible with swing away armrest or amputee axle plate	Left			\$385
AC08	<input type="checkbox"/>	Spoke guards ⁴	Available: 20", 22", 24", 26"			\$150
AC09	<input type="checkbox"/>	Footrest neoprene impact guard (pair) - Black				\$82
AC10	<input type="checkbox"/>	Seat pouch				\$98
AC30	<input type="checkbox"/>	Newton multi-purpose pouch				\$72
AC26	<input type="checkbox"/>	BodyPoint Mobility Bag				\$140
AC23	<input type="checkbox"/>	Newton Backpack				\$125
AC12	<input type="checkbox"/>	Velcro Closing Strap				\$52
AC13	<input type="checkbox"/>	Tool kit				\$52
AC16	<input type="checkbox"/>	Motion Composites Polo Shirt	W/S <input type="checkbox"/>	W/M <input type="checkbox"/>	W/L <input type="checkbox"/>	\$35
			M/S <input type="checkbox"/>	M/M <input type="checkbox"/>	M/L <input type="checkbox"/>	
				M/XL <input type="checkbox"/>	M/XXL <input type="checkbox"/>	
					M/XXXL <input type="checkbox"/>	
					M/XXXXL <input type="checkbox"/>	

1- Must select the pair (left and right) with HD option; Only available with rear seat to floor height between 13 1/2" and 19 3/4"


2 - Cannot be combined with Anti Tippers

3 - May interfere with swing away armrest and Clamp on height adjustable push handle; Not available with HD option

4 - Not available with Mag wheels and Spinergy

FRAME COLOR (with black aluminum crossbrace)						
FC01	<input type="checkbox"/>	Black				NC
FC	<input type="checkbox"/>	(03) Sunkissed Orange <input type="checkbox"/>	(04) Burgundy <input type="checkbox"/>	(05) Sapphire Blue <input type="checkbox"/>	(06) Monster Green <input type="checkbox"/>	\$270
		(07) Charcoal <input type="checkbox"/>	(08) Fuchsia <input type="checkbox"/>	(09) Steel Blue <input type="checkbox"/>	(17) Ferrari Red <input type="checkbox"/>	(18) Acid Green <input type="checkbox"/>
FC99	<input type="checkbox"/>	Custom Color Specify PPG Paintit color:			additional 10 days lead-time	\$410

Fork and Stem bolt combinations




Unilateral wheel lock restrictions

Seat depth must be 16in and greater		
Rear wheel diameter 24 in and less		
"T" or flip-back armrest only available with 0° Camber		
One arm drive: Not compatible with "T" or flip-back armrest		
"T" or flip-back armrest may require to move the rear wheels outwards		
Wheel lock lever extension not available		
Rear wheel diameter	Seat depth (in)	COG Compatibility Range (in)
20 in	16+	2
	16	0.75 to 2
22 in	17+	0.75 to 2
	16	0.75
24 in	17	0.75 to 2
	18+	0.75 to 2

Notes :

BioFit™ CUSHION	WIDTH											
	12"	13"	14"	15"	16"	17"	18"	19"	20"	21"	22"	
	12"	<input type="checkbox"/> \$595	<input type="checkbox"/> \$595	<input type="checkbox"/> \$595								
14"	<input type="checkbox"/> \$595	<input type="checkbox"/> \$595	<input type="checkbox"/> \$595									
15"		<input type="checkbox"/> \$595	<input type="checkbox"/> \$595	<input type="checkbox"/> \$595	<input type="checkbox"/> \$595							
16"		<input type="checkbox"/> \$595	<input type="checkbox"/> \$595	<input type="checkbox"/> \$595	<input type="checkbox"/> \$595	<input type="checkbox"/> \$595	<input type="checkbox"/> \$595	<input type="checkbox"/> \$625	<input type="checkbox"/> \$625	<input type="checkbox"/> \$625	<input type="checkbox"/> \$625	
17"				<input type="checkbox"/> \$595	<input type="checkbox"/> \$595	<input type="checkbox"/> \$595	<input type="checkbox"/> \$595	<input type="checkbox"/> \$625	<input type="checkbox"/> \$625			
18"				<input type="checkbox"/> \$595	<input type="checkbox"/> \$595	<input type="checkbox"/> \$595	<input type="checkbox"/> \$595	<input type="checkbox"/> \$625	<input type="checkbox"/> \$625	<input type="checkbox"/> \$625	<input type="checkbox"/> \$625	
19"				<input type="checkbox"/> \$595	<input type="checkbox"/> \$595	<input type="checkbox"/> \$595	<input type="checkbox"/> \$595	<input type="checkbox"/> \$625	<input type="checkbox"/> \$625			
20"				<input type="checkbox"/> \$595	<input type="checkbox"/> \$595	<input type="checkbox"/> \$595	<input type="checkbox"/> \$625	<input type="checkbox"/> \$625	<input type="checkbox"/> \$625	<input type="checkbox"/> \$625	<input type="checkbox"/> \$625	

Optional Incontinent / Infection Control Outer Cover NXCAF _ _ _ _ V N/C



NXCABF _ _ _ _
(W" x D")

HCPCS CODE: E2607-E2608

Optima™ THORACIC BACK SUPPORT

THORACIC SHAPE	Optima™	Fits Width	Height	✓	\$
	NOACB1613	16"-18"	13"	<input type="checkbox"/>	\$690
	NOACB1616	16"-18"	16"	<input type="checkbox"/>	\$690
	NOACB1619	16"-18"	19"	<input type="checkbox"/>	\$690
	NOACB1813	18"-20"	13"	<input type="checkbox"/>	\$690
	NOACB1816	18"-20"	16"	<input type="checkbox"/>	\$690
	NOACB1819	18"-20"	19"	<input type="checkbox"/>	\$690
	NOACB2016	20"-22"	16"	<input type="checkbox"/>	\$690
	NOACB2019	20"-22"	19"	<input type="checkbox"/>	\$690
	NOACB2213	22"-24"	13"	<input type="checkbox"/>	\$720
	NOACB2216	22"-24"	16"	<input type="checkbox"/>	\$720
	NOACB2219	22"-24"	19"	<input type="checkbox"/>	\$720

HCPCS CODE: E2620-E2621

Optima™ DEEP THORACIC BACK SUPPORT

Optima™ DEEP	Fits Width	Height	✓	\$
NOPCB1616	16"-18"	16"	<input type="checkbox"/>	\$790
NOPCB1619	16"-18"	19"	<input type="checkbox"/>	\$790
NOPCB1816	18"-20"	16"	<input type="checkbox"/>	\$790
NOPCB1819	18"-20"	19"	<input type="checkbox"/>	\$790
NOPCB2016	20"-22"	16"	<input type="checkbox"/>	\$790
NOPCB2019	20"-22"	19"	<input type="checkbox"/>	\$790
NOPCB2213	22"-24"	13"	<input type="checkbox"/>	\$820
NOPCB2216	22"-24"	16"	<input type="checkbox"/>	\$820
NOPCB2219	22"-24"	19"	<input type="checkbox"/>	\$820

HCPCS CODE: E2620-E2621



Xtend™ HEIGHT ADJUSTABLE THORACIC BACK SUPPORT

Xtend™	Fits Width	Height	✓	\$
NHAACB1616/20	16"-18"	16"-20"	<input type="checkbox"/>	\$690
NHAACB1816/20	18"-20"	16"-20"	<input type="checkbox"/>	\$690
NHAACB2016/20	20"-22"	16"-20"	<input type="checkbox"/>	\$690
NHAACB2216/20	22"-24"	16"-20"	<input type="checkbox"/>	\$720

HCPCS CODE: E2620-E2621

Xtend™ DEEP HEIGHT ADJUSTABLE THORACIC BACK SUPPORT

Xtend™ DEEP	Fits Width	Height	✓	\$
NHAPCB1616/20	16"-18"	16"-20"	<input type="checkbox"/>	\$690
NHAPCB1816/20	18"-20"	16"-20"	<input type="checkbox"/>	\$690
NHAPCB2016/20	20"-22"	16"-20"	<input type="checkbox"/>	\$690
NHAPCB2216/20	22"-24"	16"-20"	<input type="checkbox"/>	\$720

HCPCS CODE: E2620-E2621



Classic™ PELVIC BACK SUPPORT

Classic™	Fits Width	Height	✓	\$
NACB1616	16"-18"	16"	<input type="checkbox"/>	\$690
NACB1618	16"-18"	18"	<input type="checkbox"/>	\$690
NACB1816	18"-20"	16"	<input type="checkbox"/>	\$690
NACB1818	18"-20"	18"	<input type="checkbox"/>	\$690
NACB1820	18"-20"	20"	<input type="checkbox"/>	\$690
NACB2020	20"-22"	20"	<input type="checkbox"/>	\$690
NACB2216	22"-24"	16"	<input type="checkbox"/>	\$720
NACB2218	22"-24"	18"	<input type="checkbox"/>	\$720
NACB2220	22"-24"	20"	<input type="checkbox"/>	\$720

HCPCS CODE: E2620-E2621

Classic™ DEEP PELVIC BACK SUPPORT

Classic DEEP™	Fits Width	Height	✓	\$
NPCB1616	16"-18"	16"	<input type="checkbox"/>	\$790
NPCB1618	16"-18"	18"	<input type="checkbox"/>	\$790
NPCB1818	18"-20"	18"	<input type="checkbox"/>	\$790
NPCB1820	18"-20"	20"	<input type="checkbox"/>	\$790
NPCB2020	20"-22"	20"	<input type="checkbox"/>	\$790
NPCB2216	22"-24"	16"	<input type="checkbox"/>	\$820
NPCB2218	22"-24"	18"	<input type="checkbox"/>	\$820
NPCB2220	22"-24"	20"	<input type="checkbox"/>	\$820

HCPCS CODE: E2620-E2621



MultiFit™ HEAD SUPPORT - MULTI-ADJUSTABLE PADS

	HCPCS CODE: E0955
Model #NIB-0784	✓ \$
10" x 5" Large Rectangular	<input type="checkbox"/> \$195

	HCPCS CODE: E0955
Model # NIB-0818	✓ \$
14" x 5" Extra-Large Rectangular	<input type="checkbox"/> \$195

MultiFit™ HEAD SUPPORT - FIXTURES

	Model #NIB-0747	✓ \$
Multi-Axis Hardware	<input type="checkbox"/> \$145	

	Model # NIB-0823	✓ \$
Head Support Receiver Bracket	<input type="checkbox"/> \$55	

ALL PRICES QUOTED IN CDN. PRICES ARE SUBJECT TO CHANGE WITHOUT NOTICE.