



Maple Leaf WHEELCHAIR

SUPER TILT



The **SUPERTILT** is Maple Leaf Wheelchair's flagship product. It is a steel manual tilt wheelchair designed perfectly for those long term users in need of pressure relief. It is engineered with the most consistent and smooth tilting function in the market.




Wide range of colour selection.
Also custom colour matching available.

**available on all chairs*

◆ Angle adjustable front casters



◆ Width Adjustable Frame



◆ Two-way locking gas struts for the smoothest most consistent tilting feature in the market



45° of Tilt for pressure relief

◆ 12" Wheel Assembly Available
**custom inside frame mount to narrow overall width also available with 12" wheels*



Maple Leaf WHEELCHAIR

LOWRIDER

The **LOWRIDER** is our dynamic tilt wheelchair. Its design allows up to 35 degrees of pressure relieving tilt while maintaining the seat to floor height, perfect for foot propelling.



- ◆ Self-Tilting Option (Power Tilt)
*Also available on the SuperTilt



- ◆ Front pivot point with recessed frame for easy foot propulsion

Feature Item

- ◆ New Wheelchair Fender Guards available



Standard Features SuperTilt & LowRider

- ◆ Angle Adjustable Back Posts
- ◆ Pin & Plug-In Footrests available
- ◆ Angle Adjustable Stroller Handle



Specifications SuperTilt & LowRider

- ◆ 14" to 36" chair widths
- ◆ 16" to 24" chair depths
- ◆ 13.5" to 18.25" seat heights
- ◆ 250 lbs. Weight Capacity
- ◆ Available heavy duty packages
- ◆ Customize to any specifications
- ◆ Will accommodate any seating system
- ◆ 5 year warranty on frame
- ◆ 2 year warranty on additional components

Visit our website for a full list of standard and additional options

www.mapleleafwheelchair.ca