

COMPANY _____ SHIP TO _____

DATE _____ NAME _____ PHONE _____

CONTACT NAME _____ ADDRESS _____

EMAIL _____ CITY _____ PROVINCE _____ POSTAL CODE _____



Optima™ THORACIC BACK SUPPORT

CARBON FIBRE SHELL = WEIGHT SAVINGS UP TO 2.4 LBS WITH FIXED MICRO-MOUNT HARDWARE OPTION

Has a 3" (7.6 cm) mild thoracic contour. Weight capacity for backs up to 20" (50.8 cm) is 300lbs (136 kg), for 22" (55.9 cm) and larger is 400lbs (181.5 kg). When using 4-Point Fixed Mount Hardware, 23" (58.4 cm) and larger have a weight capacity of 550lbs.



Optima™ DEEP THORACIC BACK SUPPORT

CARBON FIBRE SHELL = WEIGHT SAVINGS UP TO 2.4 LBS WITH FIXED MICRO-MOUNT HARDWARE OPTION

Has a depth adjustable deep thoracic contour of 5.5" - 7" (14-17.8 cm). Weight capacity for backs up to 20" (50.8 cm) is 300lbs (136 kg), for 22" (55.9 cm) and larger is 400lbs (181.5 kg). When using four point fixed mount hardware, 23" (58.4 cm) and larger have a weight capacity of 550lbs.

THORACIC SHAPE

	Model No. Aluminum	Width	Height	\$ Carbon	QTY	\$ Aluminum	QTY
	NOACB1307	13"	7"			725	
	NOACB1309	13"	9"			725	
	NOACB1313	13"	13"			725	
	NOACB1316	13"	16"			725	
	NOACB1407	14"	7"			725	
NOCFACB1409	NOACB1409	14"	9"	1150		725	
NOCFACB1413	NOACB1413	14"	13"	1150		725	
NOCFACB1416	NOACB1416	14"	16"	1150		725	
	NOACB1507	15"	7"			725	
NOCFACB1509	NOACB1509	15"	9"	1150		725	
NOCFACB1513	NOACB1513	15"	13"	1150		725	
NOCFACB1516	NOACB1516	15"	16"	1150		725	
	NOACB1607	16"	7"			725	
NOCFACB1609	NOACB1609	16"	9"	1150		725	
NOCFACB1613	NOACB1613	16"	13"	1150		725	
NOCFACB1616	NOACB1616	16"	16"	1150		725	
	NOACB1619	16"	19"			725	
	NOACB1707	17"	7"			725	
NOCFACB1709	NOACB1709	17"	9"	1150		725	
NOCFACB1713	NOACB1713	17"	13"	1150		725	
NOCFACB1716	NOACB1716	17"	16"	1150		725	
	NOACB1719	17"	19"			725	
	NOACB1807	18"	7"			725	
NOCFACB1809	NOACB1809	18"	9"	1150		725	
NOCFACB1813	NOACB1813	18"	13"	1150		725	
NOCFACB1816	NOACB1816	18"	16"	1150		725	
	NOACB1819	18"	19"			725	
	NOACB1916	19"	16"			725	
	NOACB1919	19"	19"			725	
	NOACB2007	20"	7"			725	
NOCFACB2009	NOACB2009	20"	9"	1150		725	
NOCFACB2013	NOACB2013	20"	13"	1150		725	
NOCFACB2016	NOACB2016	20"	16"	1150		725	
	NOACB2019	20"	19"			725	

	Model No. Carbon Fibre	Model No. Aluminum	Width	Height	\$ Carbon	QTY	\$ Aluminum	QTY
		NOPCB1307	13"	7"			825	
		NOPCB1309	13"	9"			825	
		NOPCB1313	13"	13"			825	
		NOPCB1316	13"	16"			825	
		NOPCB1407	14"	7"			825	
NOCFPCB1409	NOPCB1409		14"	9"	1250		825	
NOCFPCB1413	NOPCB1413		14"	13"	1250		825	
NOCFPCB1416	NOPCB1416		14"	16"	1250		825	
		NOPCB1507	15"	7"			825	
NOCFPCB1509	NOPCB1509		15"	9"	1250		825	
NOCFPCB1513	NOPCB1513		15"	13"	1250		825	
NOCFPCB1516	NOPCB1516		15"	16"	1250		825	
		NOPCB1607	16"	7"			825	
NOCFPCB1609	NOPCB1609		16"	9"	1250		825	
NOCFPCB1613	NOPCB1613		16"	13"	1250		825	
NOCFPCB1616	NOPCB1616		16"	16"	1250		825	
		NOPCB1619	16"	19"			825	
		NOPCB1707	17"	7"			825	
NOCFPCB1709	NOPCB1709		17"	9"	1250		825	
NOCFPCB1713	NOPCB1713		17"	13"	1250		825	
NOCFPCB1716	NOPCB1716		17"	16"	1250		825	
		NOPCB1719	17"	19"			825	
		NOPCB1807	18"	7"			825	
NOCFPCB1809	NOPCB1809		18"	9"	1250		825	
NOCFPCB1813	NOPCB1813		18"	13"	1250		825	
NOCFPCB1816	NOPCB1816		18"	16"	1250		825	
		NOPCB1819	18"	19"			825	
		NOPCB1916	19"	16"			825	
		NOPCB1919	19"	19"			825	
		NOPCB2007	20"	7"			825	
NOCFPCB2009	NOPCB2009		20"	9"	1250		825	
NOCFPCB2013	NOPCB2013		20"	13"	1250		825	
NOCFPCB2016	NOPCB2016		20"	16"	1250		825	
		NOPCB2019	20"	19"			825	

ITEM	WIDTH X DEPTH	\$	QTY
Replacement Cover & Foam	___ X ___	130	

ITEM	WIDTH X DEPTH	\$	QTY
Replacement Cover & Foam	___ X ___	130	

PLEASE NOTE: Each back will fit widths up to 2" wider.

Optima™ THORACIC BACK SUPPORT

Optima™ DEEP THORACIC BACK SUPPORT

Model No. Carbon Fibre	Model No. Aluminum	Width	Height	\$ Carbon	QTY	\$ Aluminum	QTY
N/A	NOACB2116	21"	16"	N/A		755	
N/A	NOACB2119	21"	19"	N/A		755	
N/A	NOACB2213	22"	13"	N/A		755	
N/A	NOACB2216	22"	16"	N/A		755	
N/A	NOACB2219	22"	19"	N/A		755	
N/A	NOACB2316	23"	16"	N/A		755	
N/A	NOACB2319	23"	19"	N/A		755	
N/A	NOACB2416	24"	16"	N/A		755	
N/A	NOACB2419	24"	19"	N/A		755	
N/A	NOACB2516	25"	16"	N/A		755	
N/A	NOACB2519	25"	19"	N/A		755	
N/A	NOACB2616	26"	16"	N/A		755	
N/A	NOACB2619	26"	19"	N/A		755	
N/A	NOACB2716	27"	16"	N/A		755	
N/A	NOACB2719	27"	19"	N/A		755	
N/A	NOACB2816	28"	16"	N/A		755	
N/A	NOACB2819	28"	19"	N/A		755	
N/A	NOACB2916	29"	16"	N/A		755	
N/A	NOACB2919	29"	19"	N/A		755	

Model No. Carbon Fibre	Model No. Aluminum	Width	Height	\$ Carbon	QTY	\$ Aluminum	QTY
N/A	NOPCB2116	21"	16"	N/A		860	
N/A	NOPCB2119	21"	19"	N/A		860	
N/A	NOPCB2213	22"	13"	N/A		860	
N/A	NOPCB2216	22"	16"	N/A		860	
N/A	NOPCB2219	22"	19"	N/A		860	
N/A	NOPCB2316	23"	16"	N/A		860	
N/A	NOPCB2319	23"	19"	N/A		860	
N/A	NOPCB2416	24"	16"	N/A		860	
N/A	NOPCB2419	24"	19"	N/A		860	
N/A	NOPCB2516	25"	16"	N/A		860	
N/A	NOPCB2519	25"	19"	N/A		860	
N/A	NOPCB2616	26"	16"	N/A		860	
N/A	NOPCB2619	26"	19"	N/A		860	
N/A	NOPCB2716	27"	16"	N/A		860	
N/A	NOPCB2719	27"	19"	N/A		860	
N/A	NOPCB2816	28"	16"	N/A		860	
N/A	NOPCB2819	28"	19"	N/A		860	
N/A	NOPCB2916	29"	16"	N/A		860	
N/A	NOPCB2919	29"	19"	N/A		860	

ITEM	WIDTH X DEPTH	\$	QTY
Replacement Cover & Foam	__ x __	130	

ITEM	WIDTH X DEPTH	\$	QTY
Replacement Cover & Foam	__ x __	130	

ADP NOTES: ADP Non-Standard Sizes - for non-standard sizes please use ADP code SEMCF202L and apply an up-charge of \$86.00 * One Hand Release does not fit on backs smaller than 14" wide

Xtend™ HEIGHT ADJUSTABLE THORACIC BACK SUPPORT

.Step-less height adjustable range from 16"-20" (40.6-50.8 cm) with a mild, 3.5" (8.9 cm) thoracic lateral contour.
 .Weight capacity for backs up to 20" (40.6 cm) is 300lbs (136 kg), for 22" (55.9 cm) and larger is 400lbs (181.5 kg)

MODEL NO.	WIDTH	HEIGHT	\$	QTY
NHAACB1616/20	16"	16"-20"	725	
NHAACB1816/20	18"	16"-20"	725	
NHAACB2016/20	20"	16"-20"	725	
NHAACB2216/20	22"	16"-20"	755	
NHAACB2416/20	24"	16"-20"	755	

Optional: Change cover for incontinent/infection control cover

ITEM	WIDTH X DEPTH	\$	QTY
Replacement Cover & Foam	__ x __	130	

Xtend™ DEEP HEIGHT ADJUSTABLE THORACIC BACK SUPPORT

.Step-less height adjustable range from 16"-20" (40.6-50.8 cm) with a 5.5"-7" (14-17.8 cm) depth adjustable thoracic lateral contour.
 .Weight capacity for backs up to 20" (40.6 cm) is 300lbs (136 kg), for 22" (55.9 cm) and larger is 400lbs (181.5 kg).

MODEL NO.	WIDTH	HEIGHT	\$	QTY
NHAPCB1616/20	16"	16"-20"	725	
NHAPCB1816/20	18"	16"-20"	725	
NHAPCB2016/20	20"	16"-20"	725	
NHAPCB2216/20	22"	16"-20"	755	
NHAPCB2416/20	24"	16"-20"	755	

Optional: Change cover for incontinent/infection control cover

ITEM	WIDTH X DEPTH	\$	QTY
Replacement Cover & Foam	__ x __	130	

ADP NOTES: ADP code SEMCF202L can be used for 22' wide or wider back widths. ADP PRICE \$86.00

Xtend™ LOW HEIGHT ADJUSTABLE THORACIC BACK SUPPORT

.Step-less height adjustable range from 13"-17" (33-43.2 cm) with a mild, 3.5" (8.9 cm) thoracic lateral contour.
 .Weight capacity for backs up to 20" (50.8 cm) is 300lbs (136 kg), for 22" (55.9 cm) and larger is 400lbs (181.5 kg).

MODEL NO.	WIDTH	HEIGHT	\$	QTY
**NHAACB (L) 1213/17	12"	13"-17"	725	
NHAACB (L) 1413/17	14"	13"-17"	725	
NHAACB (L) 1613/17	16"	13"-17"	725	
NHAACB (L) 1813/17	18"	13"-17"	725	
NHAACB (L) 2013/17	20"	13"-17"	755	
NHAACB (L) 2213/17	22"	13"-17"	755	

Optional: Change cover for incontinent/infection control cover

ITEM	WIDTH X DEPTH	\$	QTY
Replacement Cover & Foam	__ x __	130	

Xtend™ LOW DEEP HEIGHT ADJUSTABLE THORACIC BACK SUPPORT

.Step-less height adjustable range from 13"-17" (33-43.2 cm) with a 5.5"-7" (14-17.8 cm) depth adjustable thoracic lateral contour.
 .Weight capacity for backs up to 20" (40.6 cm) is 300lbs (136 kg), for 22" (55.9 cm) and larger is 400lbs (181.5 kg).

MODEL NO.	WIDTH	HEIGHT	\$	QTY
**NHAPCB (L) 1213/17	12"	13"-17"	725	
NHAPCB (L) 1413/17	14"	13"-17"	725	
NHAPCB (L) 1613/17	16"	13"-17"	725	
NHAPCB (L) 1813/17	18"	13"-17"	725	
NHAPCB (L) 2013/17	20"	13"-17"	755	
NHAPCB (L) 2213/17	22"	13"-17"	755	

Optional: Change cover for incontinent/infection control cover

ITEM	WIDTH X DEPTH	\$	QTY
Replacement Cover & Foam	__ x __	130	

** 12" Wide Xtend LOW does not come with one hand release

ADP NOTES: ADP code SEMCF202L can be used for 22' wide or wider back widths. ADP PRICE \$86.00

PLEASE NOTE: Each back will fit widths up to 2" wider.


All prices are quoted in Canadian Dollars. Prices are subject to change without notice. JW-07-02-24

THORACIC SHAPE

Classic™ PELVIC BACK SUPPORT

.Has a 3" (7.6 cm) mild pelvic contour.
 .Weight capacity for backs up to 20" (50.8 cm) is 300lbs (136 kg), for 22" (55.9 cm) and larger is 400lbs (181.5 kg).


PELVIC SHAPE

	MODEL NO.	WIDTH	HEIGHT	\$	QTY
	NACB1616	16"	16"	725	
	NACB1618	16"	18"	725	
	NACB1620	16"	20"	725	
	NACB1816	18"	16"	725	
	NACB1818	18"	18"	725	
	NACB1820	18"	20"	725	
	NACB2016	20"	16"	725	
	NACB2018	20"	18"	725	
	NACB2020	20"	20"	725	
	NACB2216	22"	16"	755	
	NACB2218	22"	18"	755	
	NACB2220	22"	20"	755	

ITEM	WIDTH X DEPTH	\$	QTY
Replacement Cover & Foam	__ __ x __ __	130	

Classic™ DEEP PELVIC BACK SUPPORT

.Has a depth adjustable deep pelvic contour of 5.5"- 7" (14-17.8 cm).
 .Weight capacity for backs up to 20" (50.8 cm) is 300lbs (136 kg), for 22" (55.9 cm) and larger is 400lbs (181.5 kg).

	MODEL NO.	WIDTH	HEIGHT	\$	QTY
	NPCB1616	16"	16"	830	
	NPCB1618	16"	18"	830	
	NPCB1620	16"	20"	830	
	NPCB1816	18"	16"	830	
	NPCB1818	18"	18"	830	
	NPCB1820	18"	20"	830	
	NPCB2016	20"	16"	830	
	NPCB2018	20"	18"	830	
	NPCB2020	20"	20"	830	
	NPCB2216	22"	16"	860	
	NPCB2218	22"	18"	860	
	NPCB2220	22"	20"	860	


ITEM	WIDTH X DEPTH	\$	QTY
Replacement Cover & Foam	__ __ x __ __	130	

ADP NOTES: ADP code SEMCF202L can be used for 22' wide or wider back widths. ADP PRICE \$86.00
 PLEASE NOTE: Each back will fit widths up to 2" wider.

üFit™ CUSTOMIZABLE PELVIC BACK SUPPORT

.Has a 6.25" (15.9 cm) contour (deduct 3.25" (8.3 cm) for foam cushioning).
 .Designed for users with positioning needs due to spinal deformities including Kyphosis, Scoliosis and Lordosis
 .Weight capacity for backs up to 20" (50.8 cm) is 300lbs (136 kg), for 22" (55.9 cm) and larger is 400lbs (181.5 kg).


CUSTOMIZABLE PELVIC SHAPE

	MODEL NO.	WIDTH	HEIGHT	\$	QTY
	NOUFACB1416	14"	16"	940	
	NOUFACB1418	14"	18"	940	
	NOUFACB1420	14"	20"	940	
	NOUFACB1616	16"	16"	940	
	NOUFACB1618	16"	18"	940	
	NOUFACB1620	16"	20"	940	
	NOUFACB1816	18"	16"	940	
	NOUFACB1818	18"	18"	940	
	NOUFACB1820	18"	20"	940	
	NOUFACB2016	20"	16"	940	
	NOUFACB2018	20"	18"	940	
	NOUFACB2020	20"	20"	940	
	NOUFACB2216	22"	16"	1045	
	NOUFACB2218	22"	18"	1045	
	NOUFACB2220	22"	20"	1045	

ITEM	WIDTH X DEPTH	\$	QTY
Replacement Cover & Foam	__ __ x __ __	130	

üFit™ DEEP CUSTOMIZABLE PELVIC BACK SUPPORT

.Has a 8.75" (22.2 cm) pelvic contour (deduct 3.25" (8.3 cm) for foam cushioning)
 .Designed for users with positioning needs due to spinal deformities including Kyphosis, Scoliosis and Lordosis
 .Weight capacity for backs up to 20" (50.8 cm) is 300lbs (136 kg), for 22" (55.9 cm) and larger is 400lbs (181.5 kg).


	MODEL NO.	WIDTH	HEIGHT	\$	QTY
	NOUFPCB1416	14"	16"	1040	
	NOUFPCB1418	14"	18"	1040	
	NOUFPCB1420	14"	20"	1040	
	NOUFPCB1616	16"	16"	1040	
	NOUFPCB1618	16"	18"	1040	
	NOUFPCB1620	16"	20"	1040	
	NOUFPCB1816	18"	16"	1040	
	NOUFPCB1818	18"	18"	1040	
	NOUFPCB1820	18"	20"	1040	
	NOUFPCB2016	20"	16"	1040	
	NOUFPCB2018	20"	18"	1040	
	NOUFPCB2020	20"	20"	1040	
	NOUFPCB2216	22"	16"	1255	
	NOUFPCB2218	22"	18"	1255	
	NOUFPCB2220	22"	20"	1255	

ITEM	WIDTH X DEPTH	\$	QTY
Replacement Cover & Foam	__ __ x __ __	130	

ADP NOTES: ADP code SEMCF202L can be used for 22' wide or wider back widths. ADP PRICE \$86.00


üFit™ II CUSTOMIZABLE THORACIC BACK SUPPORT

. Has a 7.75" (19.6 cm) thoracic contour (deduct 3.25" (8.3 cm) for foam cushioning)
 . Designed for users with positioning needs due to spinal deformities including Kyphosis, Scoliosis and Lordosis
 . Weight capacity for backs up to 20" (50.8 cm) is 300lbs (136 kg), for 22" (55.9 cm) and larger is 400lbs (181.5 kg).

	MODEL NO.	WIDTH	HEIGHT	\$	QTY
	NOUFXACB1416	14"	16"	955	
	NOUFXACB1418	14"	18"	955	
	NOUFXACB1420	14"	20"	955	
	NOUFXACB1616	16"	16"	955	
	NOUFXACB1618	16"	18"	955	
	NOUFXACB1620	16"	20"	955	
	NOUFXACB1816	18"	16"	955	
	NOUFXACB1818	18"	18"	955	
	NOUFXACB1820	18"	20"	955	
	NOUFXACB2016	20"	16"	955	
	NOUFXACB2018	20"	18"	955	
	NOUFXACB2020	20"	20"	955	
	NOUFXACB2216	22"	16"	1045	
NOUFXACB2218	22"	18"	1045		
NOUFXACB2220	22"	20"	1045		
ITEM	WIDTH X DEPTH	\$	QTY		
Replacement Cover & Foam	__ x __	130			

üFit™ II DEEP CUSTOMIZABLE THORACIC BACK SUPPORT

. Has a 9.55" (23.5 cm) thoracic contour (deduct 3.25" (8.3 cm) for foam cushioning)
 . Designed for users with positioning needs due to spinal deformities including Kyphosis, Scoliosis and Lordosis
 . Weight capacity for backs up to 20" (50.8 cm) is 300lbs (136 kg), for 22" (55.9 cm) and larger is 400lbs (181.5 kg).


	MODEL NO.	WIDTH	HEIGHT	\$	QTY
	NOUFXPCB1416	14"	16"	1095	
	NOUFXPCB1418	14"	18"	1095	
	NOUFXPCB1420	14"	20"	1095	
	NOUFXPCB1616	16"	16"	1095	
	NOUFXPCB1618	16"	18"	1095	
	NOUFXPCB1620	16"	20"	1095	
	NOUFXPCB1816	18"	16"	1095	
	NOUFXPCB1818	18"	18"	1095	
	NOUFXPCB1820	18"	20"	1095	
	NOUFXPCB2016	20"	16"	1095	
	NOUFXPCB2018	20"	18"	1095	
	NOUFXPCB2020	20"	20"	1095	
	NOUFXPCB2216	22"	16"	1255	
NOUFXPCB2218	22"	18"	1255		
NOUFXPCB2220	22"	20"	1255		
ITEM	WIDTH X DEPTH	\$	QTY		
Replacement Cover & Foam	__ x __	130			

üFIT™ BACK SUPPORT OPTIONS

MODEL NO.	FITS WHEELCHAIR WIDTH	BACK HEIGHT	\$	QTY
UF Custom-Visco	Optional Visco gel mid section foam overlay	Specify Back Model # when ordering _____	160	
NIB-1935	Pack of Polyurethane replacement 1" üCubes™ (Includes 10, 1"x2" dual üCube sets)	All ÜFit Back Supports	70	
NIB-1531	Pack of Polyurethane replacement 2" üCubes™ (Includes 10, 2"x2" dual üCube sets)	All ÜFit Back Supports	70	
NIB-1936	Pack of Polyurethane replacement 3" üCubes™ (includes 10, 3"x2" dual üCube sets)	All ÜFit Back Supports	70	

Armadillo™ BACK SUPPORT

. Three-section modular shell designed to adjust to fit the unique contour of the user's spine
 . Designed for users with kyphosis or other spinal deformities
 . Weight capacity for backs up to 20" (50.8 cm) is 300lbs (136 kg), for 22" (55.9 cm) and larger is 400lbs (181.5 kg)
 . Armadillo does not come with the One Hand Quick Release

	ARMADILLO™	WIDTH	HEIGHT	\$	QTY
	NMAPSB14	14"	16"	1045	
	NMAPSB16	16"	16"	1045	
	NMAPSB18	18"	18"	1045	
	NMAPSB18-TALL	18"	20"	1045	
	NMAPSB20	20"	18"	1045	
	NMAPSB20-TALL	20"	20"	1045	
	NMAPSB22	22"	18"	1195	
	NMAPSB22-TALL	22"	20"	1195	
	NMAPSB24	24"	18"	1195	
	NMAPSB24-TALL	24"	20"	1195	

Armadillo™ ADDITIONAL LATERAL CONTOUR ADD-ON

MODEL NO.	\$	QTY
NIB-1937 - 2" Lateral Contour Add-On (Pair)	94	




PLEASE NOTE: Each back will fit widths up to 2" wider.

All prices are quoted in Canadian Dollars. Prices are subject to change without notice. JW-07-02-24

BACK SUPPORT HARDWARE OPTIONS

.All Backs come standard with NIB-1004 Right and NIB-1178 Left, Quick Release Mounting Hardware Below are fixed hardware options available at no charge.

.Both the 2-point and 4-point fixed hardware have 1" depth adjustment. The 2-point fixed mount has 25 degrees of posterior & anterior angle adjustment, the 4-point hardware has an average of 15 degrees of posterior & anterior angle adjustment (varies, depending on spacing between the brackets).

ITEM	MODEL #	DESCRIPTION	OPTIMA™	OPTIMA™ CARBON	XTEND™	XTEND™ LOW	CLASSIC™	ÜFIT™	ARMADILLO	QTY
2-Point Fixed Hardware 	NIB-1207	2 Point fixed hardware	Optional 1	Optional	N/A	N/A	Optional 1	Optional 1	N/A	
2-Point Fixed Micro Mount Hardware  *Order as a pair	NIB-1642*	Micro Mount 2 Point fixed hardware	Optional 2 & 3	Optional 2 & 3	N/A	N/A	N/A	N/A	N/A	
4-Point Fixed Hardware 	NIB-1969	4 Point fixed hardware	Optional 3	Optional 3	Optional 3 & 4	Optional 3 & 4	Optional 3	Optional 3	N/A	

NOTE: * Must be ordered 2x to get the pair

1 - Use only on back supports with 16" (40.6cm) back height or lower. 2 - Use only on back supports with 13" (33cm) back height or lower. 3 - Not width adjustable - fits one back/chair width. 4 - Limits height adjustment to 3"

BACK SUPPORT OPTIONS AND ACCESSORIES

MODEL	DESCRIPTION	FITS BACK WIDTH	\$	QTY
Custom - Visco	Custom Visco Gel Foam Option (Dual density 0.5" (1.3 cm) Polyurethane foam backing and 1" (2.5 cm) top layer of Visco gel foam)	Specify back model # when ordering	160	
Custom - Cover	Optional Infection Control Wipe Down Back Support Cover	Specify back model # when ordering	205	
NIB-0190	Pelvic-Lumbar Support Pads (PAIR)	All	32	
NIB-XXXXXX	Back Support Inner Incontinent Cover. PLEASE NOTE: Specify model and size in part number when ordering: i.e. For a Classic DEEP Back 1820 the part number is - NIB-CD1820		70	
NIB-0191	Small Lateral Support Pads (PAIR)	14"-16" (35.6-40.6 cm)	32	
NIB-0192	Large Lateral Support Pads (PAIR)	17"-22" (43.2-55.9 cm)	32	
NIB-1004	Mounting Hardware Clamp Kit, right	N/A	100	
NIB-1178	Mounting Hardware Clamp Kit, left	N/A	100	
NIB-1207	2-Point Fixed Mounting Hardware (PAIR). Use only on back supports with 16" (40.6 cm) back height or lower	14"-20" (35.5-50.8 cm)	100	
NIB-1642	2-Point Micro Mount Fixed Hardware (Kit of 4). Use only on back supports with 13" (33 cm) back height or lower. 2 each required when ordering for a complete back support.	14"-20" (35.5-50.8 cm)	50	
NIB-1969	4-Point Mount Fixed Hardware	N/A	100	
NIB-1077	Mounting Hardware Clamp Kit, right for Armadillo	N/A	100	
NIB-1189	Mounting Hardware Clamp Kit, left for Armadillo	N/A	100	
NIB-1938	Extra Small Privacy Flap for 11" - 13" wide back supports	11" - 13" (27.9 - 33.0)	70	
NIB-1078	Small Privacy Flap for 14"-16" wide back supports	14"-16" (35.6-40.6 cm)	70	
NIB-1079	Medium Flap for 17"-19" wide back supports	17"-19" (43.2-48.3 cm)	70	
NIB-1080	Large Privacy Flap for 20"-22" wide back supports	20"-22" (50.8-55.9 cm)	70	
NIB-1081	Extra Large Privacy Flap for 23" - 26" wide back supports	23" - 26" (58.4 - 66.0 cm)	70	
NIB-1174	Replacement One Hand Release Kit	14"-16" (35.6-40.6 cm)	45	
NIB-1051	Replacement One Hand Release Kit	16"-22" (40.6-55.9 cm)	47	
NIB-1053	Replacement One Hand Release Kit	22"-24" (55.9-61cm)	47	

MultiFit™ HEAD SUPPORT

HEAD SUPPORT




MODEL NO.	MULTI-ADJUSTABLE PADS	\$	QTY
NIB-0785	9" x 5" CONTOURED	205	
NIB-0784	10" x 5" LARGE RECTANGULAR	205	
NIB-0818	14" x 5" EXTRA-LARGE RECTANGULAR	205	
MODEL NO.	FEATURES	\$	QTY
NIB-0747	MULTI-AXIS HARDWARE	150	
NIB-0823	HEAD SUPPORT RECEIVER BRACKET	60	
NIB-1306	FLIP BACK MULTI-AXIS HARDWARE	265	
NIB-1210	FLIP BACK KIT	110	
NIB-1200	UNIVERSAL HEAD SUPPORT ADAPTER	70	

PRIVACY FLAP



MODEL NO.	DESCRIPTION	\$	QTY
NIB-1938	Extra Small Privacy Flap for 11" - 13" wide back supports	70	
NIB-1078	Small Privacy Flap for 14"-16" wide back supports	70	
NIB-1079	Medium Flap for 17"-19" wide back supports	70	
NIB-1080	Large Privacy Flap for 20"-22" wide back supports	70	
NIB-1081	Extra Large Privacy Flap for 23" - 26" wide back supports	70	



BioFit™

- .Multi-layer Cushions that utilize Dynamic's smartGEL™ technology
- .Multi-layer design provides deep immersion and support for medium to high skin protection
- .Includes smarTX^{3D} outer cover and fluid proof inner cover.
- .Weight capacity for 10"-20" (25.4-50.8 cm) cushions is 300lbs (136 kg), 22"-26" (55.9-61 cm) cushions is 400lbs (181.5 kg).
- .Optional Incontinent / Infection Control Cover is a no charge alternative option when ordered with a cushion. Specify model number when ordering.

NXT CUSHIONS



BIOFIT™	WIDTH	DEPTH	\$	QTY	BIOFIT	WIDTH	DEPTH	\$	QTY
NXCABF1010	10"	10"	625		NXCABF1818	18"	18"	655	
NXCABF1012	10"	12"	625		*NXCABF1819	18"	19"	655	
NXCABF1014	10"	14"	625		NXCABF1820	18"	20"	655	
NXCABF1210	12"	10"	625		*NXCABF1916	18"	22"	655	
NXCABF1212	12"	12"	625		*NXCABF1917	19"	16"	655	
NXCABF1214	12"	14"	625		*NXCABF1918	19"	17"	655	
*NXCABF1312	13"	12"	625		*NXCABF1919	19"	18"	655	
*NXCABF1314	13"	14"	625		*NXCABF1920	19"	19"	655	
*NXCABF1315	13"	15"	625		NXCABF2016	19"	20"	655	
*NXCABF1316	13"	16"	625		*NXCABF2017	20"	16"	655	
NXCABF1412	14"	12"	625		NXCABF2018	20"	17"	655	
NXCABF1414	14"	14"	625		*NXCABF2019	20"	18"	655	
*NXCABF1415	14"	15"	625		NXCABF2020	20"	20"	655	
NXCABF1416	14"	16"	625		*NXCABF2022	20"	22"	655	
*NXCABF1515	15"	15"	625		*NXCABF2116	21"	16"	655	
*NXCABF1516	15"	16"	625		*NXCABF2118	21"	18"	655	
*NXCABF1517	15"	17"	625		*NXCABF2120	21"	20"	655	
*NXCABF1518	15"	18"	625		NXCABF2216	22"	16"	655	
*NXCABF1519	15"	19"	625		NXCABF2218	22"	18"	655	
*NXCABF1520	15"	20"	625		NXCABF2220	22"	20"	655	
*NXCABF1615	16"	15"	625		*NXCABF2222	22"	22"	655	
NXCABF1616	16"	16"	625		*NXCABF2316	23"	16"	655	
*NXCABF1617	16"	17"	625		*NXCABF2318	23"	18"	655	
NXCABF1618	16"	18"	625		*NXCABF2320	23"	20"	655	
*NXCABF1619	16"	19"	625		NXCABF2416	24"	16"	655	
NXCABF1620	16"	20"	625		NXCABF2418	24"	18"	655	
*NXCABF1716	17"	16"	625		NXCABF2420	24"	20"	655	
*NXCABF1717	17"	17"	625		*NXCABF2422	24"	22"	730	
*NXCABF1718	17"	18"	625		*NXCABF2616	26"	16"	730	
*NXCABF1719	17"	19"	625		*NXCABF2618	26"	18"	730	
*NXCABF1720	17"	20"	625		*NXCABF2620	26"	20"	730	
NXCABF1816	18"	16"	625		*NXCABF2622	26"	22"	730	
*NXCABF1817	18"	17"	625		Optional: Change cover for incontinent/infection control cover				

* These are our non-standard/custom cushion sizes. ADP code SESND1040 can be used for non-standard / custom sizing. Width or depth or both. ADP price \$69. / customization
 ** ADP code SESND1050 can be used for extra replacement covers-standard sizes. ADP price \$65.
 *** ADP code SESND1055 can be used for extra covers - non-standard sizes. ADP price \$99.
 **** ADP code SESND1060 can be used for extra covers for Bariatric sizes. ADP price \$134.

ITEM	WIDTH X DEPTH	\$	QTY
Replacement Incontinent Cover	___ x ___	130	For a cover up to 18" x 19"
		150	For a cover 18" x 20" and above
Replacement Smartx ^{3D} Cover	___ X ___	130	For a cover up to 18" x 19"
		150	For a cover 18" x 20" and above

NüFit™

.Multi-layer Contoured Cushions designed for users who need a comfortable cushion with good positioning and are at low to medium risk of skin breakdown
 .Includes smarTX^{3D} outer cover and fluid proof inner cover
 .Weight capacity 300lbs (136 kg).
 .Optional Incontinent / Infection Control Cover is a no charge alternative option when ordered with a cushion. Specify model number when ordering.




NXT CUSHIONS

NüFit™	Width	Depth	\$	QTY
NXCANF1616	16"	16"	475	
NXCANF1618	16"	18"	475	
NXCANF1620	16"	20"	475	
NXCANF1816	18"	16"	475	
NXCANF1818	18"	18"	475	
*NXCANF1820	18"	20"	475	
*NXCANF2016	20"	16"	475	
*NXCANF2018	20"	18"	475	
*NXCANF2020	20"	20"	475	
*NXCANF2216	22"	16"	475	
*NXCANF2218	22"	18"	475	
*NXCANF2220	22"	20"	475	

Optional: Change cover for incontinent/infection control cover

ITEM	WIDTH X DEPTH	\$		QTY
Replacement Incontinent Cover	___ x ___	130	For a cover up to 18" x 18"	
		150	For a cover 18" x 20" and above	
Replacement Smartx ^{3D} Cover	___ X ___	130	For a cover up to 18" x 18"	
		150	For a cover 18" x 20" and above	



Solid Seat Insert



Pelvic Obliquity Wedges

SEATING AND CUSHION ACCESSORIES

MODEL NO.	DESCRIPTION	FITS BACK WIDTH	\$	QTY
NIB-0711	Pair of 0.5" (1.2 cm) pelvic build up wedges	All Sizes	80	
RSSI 1600	Solid Seat Insert	Fits cushions 16" & 17" wide (40.6 & 43.2 cm) and 14"-18" deep (35.6-45.7 cm)	80	
RSSI 1800	Solid Seat Insert	Fits cushions 18" & 19" wide (45.7 & 48.3 cm) and 14"-18" deep (35.6-45.7 cm)	80	
RSSI 2000	Solid Seat Insert	Fits cushions 20"-22" wide (50.8-55.9 cm) and 18"-20" deep (45.7-50.8 cm)	80	

BIOFIT AND NÜFIT™ SEAT CUSHIONS CUSTOMIZATION

- > When ordering custom cushions: Please specify which cushion model you wish to purchase and then the modification.
- > Note: The modification will be made to the specific cushion size selected. i.e. "NXCABF1818 with NSDR2."

MODEL NO.	DESCRIPTION	DIMENSIONS OF CHANGE	\$	QTY
NSDR1	Seat depth reduction	Cushion depth less 1"	\$160	
NSDR2	Seat depth reduction	Cushion depth less 2"	\$160	
NSDR3	Seat depth reduction	Cushion depth less 3"	\$160	
NLLDR1	Right leg length discrepancy	Cushion depth users right less 1"	\$160	
NLLDR2	Right leg length discrepancy	Cushion depth users right less 2"	\$160	
NLLDL1	Left leg length discrepancy	Cushion depth users left less 1"	\$160	
NLLDL2	Left leg length discrepancy	Cushion depth users left less 2"	\$160	
NCWR	Cushion width reduction	Special request - requires approval	\$160	