



**Pride Mobility Products Co.**

5096 South Service Rd, Beamsville, ON, L0R 1B3

Phone: 888-570-1113 | Fax: 866-514-1303 | Email: productorders@pridemobility.com

## Q6 Edge® 3 Order Form with Contoured Seating

Account Number: \_\_\_\_\_ Date: \_\_\_\_\_

Provider Name: \_\_\_\_\_

Contact: \_\_\_\_\_

Phone: \_\_\_\_\_ Fax: \_\_\_\_\_

Email: \_\_\_\_\_

PO Number: \_\_\_\_\_

Marked for: \_\_\_\_\_

Ship to Address: \_\_\_\_\_

City: \_\_\_\_\_ Prov: \_\_\_\_\_ PC: \_\_\_\_\_



### INTRODUCTION

This form is interactive when viewed with [Adobe Acrobat Reader](#). Complete the form by placing checks in the desired boxes and provide information in the interactive fields. Buttons shown at the bottom of the form may be utilized to print or submit the order form through a desktop email application. To email via a web-based application, please 'Save As' and attach the PDF to your email.

This order form contains a large variety of options to fulfill various patient needs. Descriptions and section notes, such as optional or required, are included to help you complete your order. Other part numbers can be found on the quote acknowledgment after this order form is submitted. Please contact Quantum® Sales at 888-570-1113 if assistance is needed. Send the completed order form by fax (866-514-1303) or email ([productorders@pridemobility.com](mailto:productorders@pridemobility.com)). Incomplete forms may delay the quote or order. Customer service will contact you if the order is incomplete or if there are compatibility issues. If special order requests are needed, be sure that the Patient Information section is completed

All prices are MSRP. Prices, specifications, part numbers, and availability are subject to change without notice. Prices and part numbers as shown are only available when configured as a complete power chair. Please contact technical service for accurate parts ordering.

Thank you for choosing Quantum!


### PATIENT INFORMATION

Weight: \_\_\_\_\_ Height: \_\_\_\_\_

A. Shoulder Width: \_\_\_\_\_ D. Max Sitting Height: \_\_\_\_\_ G. Thigh Depth: \_\_\_\_\_

B. Chest Width: \_\_\_\_\_ E. Shoulder Height: \_\_\_\_\_ H. Lower Leg Length: \_\_\_\_\_

C. Hip Width: \_\_\_\_\_ F. Axilla Height: \_\_\_\_\_ I. Elbow Height: \_\_\_\_\_



### 1. BASE MODEL

**REQUIRED.** Select a group 3 power base. 300 lb. weight capacity with standard 6 mph programmable motors and 8 amp off board charger. XLR USB charger is included. Dimensions of base without legrests: 24"Wx35.5"L

Quantum Q6 Edge 3 3S-SS ..... \$8,450  
Part: Q6 EDGE 3. ADP: WP0001575 (ADULT), WP0001576 (PEDIATRIC)

**2. DRIVE MOTORS**


**REQUIRED.**

- |   |   |
|---|---|
| 6 mph 2 Pole Drive Motors..... Standard<br>Part: MOT147433  | 7 mph, 4 Pole Drive Motors..... \$700<br>Part: MOT1902455<br>Requires Q-Logic 3 electronics.  |
| 6 mph, 2 Pole, Accu-Trac Motor Technology..... \$995<br>Part: MOT147881.<br>Requires Q-Logic 3 electronics. Delivers enhanced traction and tracking, especially with switched drive controls. | 7 mph, 4 Pole, Accu-Trac Motor Technology..... \$1,695<br>Part: MOT2008501<br>Requires Q-Logic 3 electronics. Delivers enhanced traction and tracking, especially with switched drive controls. |

**3. BASE COLOR**

**REQUIRED.** Select one shroud color. Colors shown here may differ from actual product depending on monitor/printer calibration used.

**3.1a. Standard Shrouds**

- |   |  |  |
|---|--|--|
|  <b>Red Corvette</b><br>Part: DGN1900199-PRD1  |  <b>Raspberry Beret</b><br>Part: DGN1900199-PPK4  |  <b>Blue Jean</b><br>Part: DGN1900199-PBU11     |
|  <b>Back in Black</b><br>Part: DGN1900199-PBK1 |  <b>Pretty in Pink</b><br>Part: DGN1900199-PPK3   |  <b>Lemon Crush</b><br>Part: DGN1900199-PYL4    |
|  <b>Night Train</b><br>Part: DGN1900199-PBK3   |  <b>Evolution Orange</b><br>Part: DGN1900199-POR7 |  <b>Green Machine</b><br>Part: DGN1900199-PGN11 |
|  <b>Grey Street</b><br>Part: DGN1900199-PSL3   |  <b>Midnight Blue</b><br>Part: DGN1900199-PBU10   |  <b>Military Green</b><br>Part: DGN1900199-PGN8 |
|  <b>Deep Purple</b><br>Part: DGN1900199-PPR5   |  <b>Electric Blue</b><br>Part: DGN1900199-PBU9    |  |

**3.1b. Hydrographic Shrouds**

- |  |  |   |
|--|--|---|
| <input type="checkbox"/>  <b>Military Digital Camo</b> ..... \$325<br>Part: DGN1801101-WCM7 | <input type="checkbox"/>  <b>Faded Glory</b> ..... \$325<br>Part: DGN1801101-WFL2 | <input type="checkbox"/>  <b>Carbon Fiber</b> ..... \$325<br>Part: DGN1801101-WCB1 |
|--|--|---|

**4. POWER ADJUSTABLE SEAT HEIGHT**

**Optional.** The seat to floor height will be increased. Lift/elevate will be operated through the joystick. Power elevating seat includes a 1/4 speed inhibit.

- 12" Power Lift for Captain Seat with iLevel..... \$2,995  
Part: PTOINDV3857. HCPCS: E2300  
For use with K0857 base code. Must select [Power Lift Captain Seat](#) in contoured seat section. Includes armrests and lap belt. Available with swing away legrests or manually adjustable center mount foot platform.

**5. BASE MOUNTED LIGHTS**

**Optional.** Select one option.

- |  |  |
|--|--|
| Switch Fender Lights..... No Charge<br>Part: DGN154117<br>Fender lights operate by pressing included single switch button. Fender lighting is not available with Single Actuator Control through Toggle Box. | Base Mounted, LED, Full Light Package ..... \$770<br>Part: XRF000129<br>Requires Q-Logic 3 Joystick for Full Light Package or Q-Logic 3e Joystick. |
|--|--|

**6. DRIVE WHEELS**

**REQUIRED.** Select a type of drive wheels and accent color. Knobby tires increase overall base width by 1.5".

**6.1. Drive Wheels**

- |   |  |
|---|--|
| 3"x14" Flat Free Black Tires ..... Standard<br>Part: PTOINDV3525 (pair) | 3"x14" Pneumatic Black Tires ..... No Charge<br>Part: WHL1713925 x2    |
| 4"x14" Knobby, Flat Free Black Tires ..... \$120<br>Part: WHL2006647 x2 | 4"x14" Knobby, Pneumatic Black Tires..... \$120<br>Part: WHL2006822 x2 |

**6.2. Drive Wheel Accent Color**

Not available with knobby tires.

- |  |   |  |
|--|---|--|
| <b>Back in Black</b> .....Standard<br>Part: DGN164245-PBK1 x2    | <b>Evolution Orange</b> .....No Charge<br>Part: DGN164245-POR7 x2 | <b>Military Digital Camo</b> ..... \$50 ea.<br>Part: DGN164245-WCM7 x2 |
| <b>Red Corvette</b> .....No Charge<br>Part: DGN164245-PRD1 x2    | <b>Midnight Blue</b> .....No Charge<br>Part: DGN164245-PBU10 x2   | <b>Faded Glory</b> ..... \$50 ea.<br>Part: DGN164245-WFL2 x2           |
| <b>Night Train</b> .....No Charge<br>Part: DGN164245-PBK3 x2     | <b>Electric Blue</b> .....No Charge<br>Part: DGN164245-PBU9 x2    | <b>Carbon Fiber</b> ..... \$50 ea.<br>Part: DGN164245-WCB1 x2          |
| <b>Grey Street</b> .....No Charge<br>Part: DGN164245-PSL3 x2     | <b>Blue Jean</b> .....No Charge<br>Part: DGN164245-PBU11 x2       |  |
| <b>Deep Purple</b> .....No Charge<br>Part: DGN164245-PPR5 x2     | <b>Lemon Crush</b> .....No Charge<br>Part: DGN164245-PYL4 x2      |  |
| <b>Raspberry Beret</b> .....No Charge<br>Part: DGN164245-PPK4 x2 | <b>Green Machine</b> .....No Charge<br>Part: DGN164245-PGN11 x2   |  |
| <b>Pretty in Pink</b> .....No Charge<br>Part: DGN164245-PPK3 x2  | <b>Military Green</b> .....No Charge<br>Part: DGN164245-PGN8 x2   |  |

**7. BATTERIES**

**Optional.** The powerchair requires 2 batteries and will be installed if ordered.

NF-22, Gel Cell, Interceptor Battery..... \$300 ea.  
Part: BAT1704613.

NF-22, Gel Cell, MK Battery ..... \$400 ea.  
Part: BATGEL1008.

**8. PATIENT'S WEIGHT**

**REQUIRED.**

Up to 200 lb User Weight  
Part: PTOINDV3837

201 lb-300 lb User Weight  
Part: PTOINDV3838

**10. ELECTRONICS**

**REQUIRED.** Select a type of drive control from the following sections. The No Joystick option is used to omit the joystick and utilize another control, such as a specialty control at the end of this form or found on the [Alternative Controls and Electronics Order Form](#).

**10.1a. Non-Expandable Controls**

Select a joystick below or omit selection if another controller type is desired. Non-expandable does not allow supplementary drive controllers.

NE 4-Key Joystick.....Standard  
Part: CTLDC1554

NE+ 8-Key Joystick (2 Actuators) ..... No Charge  
Part: CTLDC1561

NE+ 6-Key Joystick..... No Charge  
Part: CTLDC1560

Control up to 2 actuators through joystick and only Seat Mounted, LED, full light package.

No Joystick, Setup for NE+ Joystick ..... No Charge  
Part: XRF000009



NE+ 6 Key Joystick

**10.1b. Expandable Controls - Quantum Intelligent Systems™**

Select expandable controls or omit selection if another controller type is desired. Expandable controls allow multiple or specialized controllers.

**Q-Logic 3 Joystick**

Q-Logic 3 EX Joystick .....\$900  
Part: CTL167724

Q-Logic 3 EX Joystick for Full Lights .....\$900  
Part: CTL167725

Q-Logic 3e Joystick .....\$620  
Part: CTL1910600

No Joystick, Setup for Q-Logic  
Part: XRF000010

Can operate full light package. Not available with alternative or specialty controls.



Q-Logic 3 EX Joystick

**10.2. Q-Logic 3 EX Joystick Armor**

**REQUIRED with Q-Logic 3 EX Joystick selections above. Not available with Q-Logic 3e Joystick.** Reference colors shown on page 2. \*Night train joystick armor will be the black, soft touch.

Black, Soft Touch ..... No Charge  
Part: DGN124307  
Also matches Night Train shrouds.

Non-matching Armor ..... No Charge  
Part: PTOINDV3422  
Select color below.

Match Armor to Base Color ..... Standard  
Part: PTOINDV3421

Red Corvette  
Part: DGN124336-PRD1

Deep Purple  
Part: DGN124336-PPR5

Midnight Blue  
Part: DGN124336-PBU10

Green Machine  
Part: DGN124336-PGN11

Carbon Fiber  
Part: DGN124336-WCB1

Back in Black  
Part: DGN124336-PBK1

Raspberry Beret  
Part: DGN124336-PPK4

Electric Blue  
Part: DGN124336-PBU9

Military Green  
Part: DGN124336-PGN8

Night Train\*  
See Black, Soft Touch

Pretty in Pink  
Part: DGN124336-PPK3

Blue Jean  
Part: DGN124336-PBU11

Military Digital Camo  
Part: DGN124336-WCM7

Grey Street  
Part: DGN124336-PSL3

Evolution Orange  
Part: DGN124336-POR7

Lemon Crush  
Part: DGN124336-PYL4

Faded Glory  
Part: DGN124336-WFL2

**11. JOYSTICK HANDLES**

**Optional.** The NE/NE Plus Joysticks are standard with a carrot handle (See A in image). The optional handle will be mounted to the joystick and the standard handle will be shipped with the unit.

Mushroom Handle (B) ..... \$120  
Part: ACCASMB1860.

Extended I-Handle (G) ..... \$140  
Part: ACCASMB1857.

Chin Cup (C)..... \$85  
Part: ACCASMB186190

Bodypoint 3" Goal Post Handle ..... \$180  
Part: PLSKNCP1057.

T-Bar (D) ..... \$160  
Part: ACCASMB1858.

Bodypoint 4" Goal Post Handle ..... \$180  
Part: PLSKNCP1058.

Medium Ball Handle (E)..... \$70  
Part: ACCASMB1859.

Soft Ball Handle (F)..... \$120  
Part: ACCASMB1856.



Joystick Handles

## 12. JOYSTICK MOUNTING BRACKETS

**Optional.** Select a mount to be used with a joystick selected in previous section.

### 12.1a. Fixed, In-line Joystick Bracket

Fixed, In-line, Bracket .....Standard  
 Part: Various




### 12.1b. Swing-Away (Retractable) Joystick Bracket

Allows the joystick to retract to the side of the armrest while maintaining the same joystick orientation.

Adjustable Swing-away ..... \$265  
 Part: Various.  
 As pictured to the right "Swing Away Joystick"

Offset Swing-away ..... \$265  
 Part: Various.  
 Magnetic latching mechanism.  
 As pictured to the right "Swing Away Joystick Offset Mount (Right)"




Basic Swing-away ..... \$265  
 Part: Various.  
 Magnetic latching mechanism.  
 Stealth Height Adjustable Receiver recommended.

### 12.1c. Specialty Joystick Mounts

Select a bracket and side to mount it.

Retract4 Mount ..... \$450  
 Part: XRF000003.  
 Full 180° of movement that translates into multi-position use, straight, right offset, left offset, or fold away. Joystick orientation varies as bracket is adjusted.

Stealth Swing-Away/Flip Down Bracket ..... \$300  
 Part: XRF000005.  
 Allows the joystick to pivot to side of armrest, or flip-down if bracket is adjusted.



Multi-Axis Mount ..... \$100  
 Part: XRF000004  
 Allows height adjustment of joystick.



### 12.2. Joystick Mount Location

Left Side Mount  
 Part: PTOINDV1225

Right Side Mount  
 Part: PTOINDV1226

## 13. CONTOURED SEATING

**Optional.** Select only one seating system. Special order seat sizes can be requested by checking the seat type and writing in the seat size in the special order seat size.

### 13.2a. Quantum Captain Seat

16"W x 16"D ..... No Charge  
 Part: SET2005145  
 16"W x 18"D ..... No Charge  
 Part: SET2005146  
 16"W x 20"D ..... No Charge  
 Part: SET2005147  
 18"W x 16"D ..... No Charge  
 Part: SET2005149  
 18"W x 18"D ..... No Charge  
 Part: SET2005150

18"W x 20"D ..... No Charge  
 Part: SET2005151  
 18"W x 22"D ..... No Charge  
 Part: SET2005908  
 20"W x 16"D ..... No Charge  
 Part: SET2005152  
 20"W x 18"D ..... No Charge  
 Part: SET2005153  
 20"W x 20"D ..... No Charge  
 Part: SET2005154

20"W x 22"D ..... No Charge  
 Part: SET2005909  
 22"W x 18"D ..... No Charge  
 Part: SET2005155  
 22"W x 20"D ..... No Charge  
 Part: SET2005156  
 22"W x 22"D ..... No Charge  
 Part: SET2005154

### 13.2b. Quantum Captain Seat with Solid Seat Pan

16"W x 16"D ..... No Charge  
 Part: SET2005394  
 16"W x 18"D ..... No Charge  
 Part: SET2005395  
 16"W x 20"D ..... No Charge  
 Part: SET2005396

18"W x 16"D ..... No Charge  
 Part: SET2005397  
 18"W x 18"D ..... No Charge  
 Part: SET2005398  
 18"W x 20"D ..... No Charge  
 Part: SET2005399

18"W x 22"D ..... No Charge  
 Part: SET2005910  
 20"W x 16"D ..... No Charge  
 Part: SET2005400  
 20"W x 18"D ..... No Charge  
 Part: SET2005401

Contoured seating section continued on following page.

**13. CONTOURED SEATING CONTINUED**

20"W x 20"D..... No Charge Part: SET2005402	22"W x 18"D..... No Charge Part: SET2005403	22"W x 22"D..... No Charge Part: SET2005405
20"W x 22"D..... No Charge Part: SET2005911	22"W x 20"D..... No Charge Part: SET2005404	

**13.3. Matching Seat Cushion**

Optional with Quantum Captain Seat with Solid Seat Pan. Specific cushion/component sizes available, part numbers, weight capacities, MSRP, HCPCS code, and special order cushions can be found on the respective complete [Stealth Cushion Order Forms](#). This section allows a quick selection of a single cushion that matches the seat pan width and depth. See specific order form for non-matching sizes, alternate covers and liners, or additional covers and liners.

**Stealth Simplicity G Cushion**

General use cushion with Stealth Tek 2, 2-way stretch, breathable cover.  
Part series: ST-SIMGXXYY. HCPCS: E2601/E2602


Simplicity G, Matching Seat Size  
Reference part: XRF000111. MSRP: \$110/\$165



**Stealth Glacial SP Cushion**

Skin protection with Coolcore, 4-way stretch, breathable cover.  
Part series: ST-GLSPXXYY. HCPCS: E2603/E2604

Glacial SP, Matching Seat Size  
Reference part: XRF000114. MSRP: \$330/\$440



**Stealth Zen SP Cushion**

Skin protection with Coolcore, 4-way stretch, breathable cover.  
Part series: ST-ZENSPXXYY. HCPCS: E2603/E2604

Zen SP, Matching Seat Size  
Reference part: XRF000115. MSRP: \$315/\$425



**Stealth Premiere P Cushion**

Positioning with Coolcore, 4-way stretch, breathable cover.  
Part series: ST-PREMPXXYY. HCPCS: E2605/E2606

Premiere P, Matching Seat Size  
Reference part: XRF000122. MSRP: \$325/\$475



**Stealth Spectrum Foam SPP Cushion**

Skin protection and positioning with Coolcore, 4-way stretch, breathable cover. Includes urethane liner.

Part series: ST-SPCFSPXXYY. HCPCS: E2607/E2608

Spectrum Foam SPP, Matching Seat Size  
Reference part: XRF000112. MSRP: \$460/\$530



**Stealth Essence SPP Cushion**

Skin protection and positioning with Coolcore, 4-way stretch, breathable cover.  
Part series: ST-ESSSPXXYY. HCPCS: E2607/E2608

Essence SPP, No Wedge, Matching Seat Size  
Reference part: XRF000118. MSRP: \$400/\$500

Essence SPP, Small Wedge, Matching Seat Size  
Reference part: XRF000119. MSRP: \$435/\$535

Essence SPP, Medium Wedge, Matching Seat Size  
Reference part: XRF000120. MSRP: \$465/\$565

Essence SPP, Large Wedge, Matching Seat Size  
Reference part: XRF000121. MSRP: \$495/\$595




**Stealth Solution SPP Cushion**

Skin protection and positioning with Coolcore, 4-way stretch, breathable cover. Includes urethane liner.

Part series: ST-SOLSPPXXYY. HCPCS: E2607/E2608

Solution SPP, Matching Seat Size  
Reference part: XRF000113. MSRP: \$375/\$430



**13.4. Armrests**

Optional. Configured in pairs. Left (L) and right (R) boxes may be checked to configure different size arms within one section arm type.

**1a. Drop-In Armrests**

Includes flat armpads.

Full Length, Drop-In Armrests..... No Charge  
L R Left part: ARM1803253. Right part: ARM156742

Desk Length, Drop-In Armrests ..... No Charge  
L R Left part: ARM1803254. Right part: ARM156743

**1b. Flip-Up, Height Adjustable Armrests**

Includes full length, flat armpads.

Full Length, Left Armrest ..... No Charge  
Part: MEC1906767

Full Length, Right Armrest..... No Charge  
Part: MEC1906768

**1c. Single Post, Height Adjustable, Removable Armrests**

Height adjustment from 11"-17" in half-inch increments.

Full Length, Single Post Armrests .....\$145 ea.  
L R Left part: ARM2000906. Right part: ARM2000907.

Desk Length, Single Post Armrests .....\$145 ea.  
L R Left part: ARM2000908. Right part: ARM2000909.

**1d. Bariatric, Height Adjustable, Removable Armrests**

Select armrests and armpads. Armrests are sold in same-size pairs.

Full Length, Bariatric Armrests, 8"-12" H .....\$130 ea.  
L R Left part: FRMASMB6728. Right part: FRMASMB6729

Full Length, Bariatric Armrests, 10"-14"H .....\$130 ea.  
L R Left part: FRMASMB6730. Right part: FRMASMB6731

**1e. Heavy Duty, Height Adjustable, Removable Armrests**

Select armrests and armpads. Armrests are sold in same-size pairs.

Full Length, Heavy Duty Armrests, 8"-12" H .....\$130 ea.  
L R Left part: XRF000253. Right part: XRF000254.

Full Length, Heavy Duty Armrests, 10"-14"H.....\$130 ea.  
L R Left Part: XRF000257. Right part: XRF000258.

**1f. Pediatric, Removable, Height Adjustable Armrests**

Height adjustment from 7.75" to 8.25". Includes desk length, straight armpads.

Pediatric, Removable Armrests ..... \$145 ea.  
L R Left Part: ARM142876. Right part: ARM142877. HCPCS: E0973

Contoured seating section continued on following pages.

**13. CONTOURED SEATING CONTINUED**

**13.5. Armrest Height**

**REQUIRED with Bariatric, Heavy Duty, and Pediatric Height Adjustable armrests.**

**2. Left Armrest Height**

7.75" H PTOINDV3722	8.25" H PTOINDV3724	10" H PTOINDV3537	12" H PTOINDV3539	14" H PTOINDV3541	16" H PTOINDV3543
8" H PTOINDV3535	9" H PTOINDV3536	11" H PTOINDV3538	13" H PTOINDV3540	15" H PTOINDV3542	

**3. Right Armrest Height**

7.75" H PTOINDV3721	8.25" H PTOINDV3723	10" H PTOINDV3452	12" H PTOINDV3454	14" H PTOINDV3456	16" H PTOINDV3458
8" H PTOINDV3450	9" H PTOINDV3451	11" H PTOINDV3453	13" H PTOINDV3455	15" H PTOINDV3457	

**13.6. Armpads**

**REQUIRED except with Flip-Up and Pediatric armrests.**

**1a. Straight Armpads for Drop-In, or Single Post Armrests**

Straight, Full Length Left and Right ..... No Charge Part: ARM134275	Straight, Desk Length Left, Full Length Right ..... No Charge Part: ARM134280
Straight, Full Length Left, Desk Length Right ..... No Charge Part: ARM134281	Straight, Desk Length Left and Right ..... No Charge Part: ARM134254

**1b. Waterfall Armpads for Drop-In, or Single Post Armrests**

Waterfall, Full Length Left and Right ..... No Charge Part: ARM134285	Waterfall, Desk Length Left, Full Length Right ..... No Charge Part: ARM134288
Waterfall, Full Length Left, Desk Length Right ..... No Charge Part: ARM134289	Waterfall, Desk Length Left and Right ..... No Charge Part: ARM134282

**1c. Falcon Gel Armpads for Drop-In, or Single Post Armrests**

Falcon Gel, Full Length Left and Right ..... \$300 Part: ARM136438	Falcon Gel, Desk Length Left, Full Length Right ..... \$300 Part: ARM136439
Falcon Gel, Full Length Left, Desk Length Right ..... \$250 Part: ARM136442	Falcon Gel, Desk Length Left and Right ..... \$200 Part: ARM136437

**1d. Armpads for Bariatric or Heavy Duty Armrests**

Bariatric Armpads ..... \$40 ea. Part: UPSASMB1027 Only for use with Bariatric Armrests.	Waterfall Armpads ..... No Charge Part: UPSASMB1024 Only for use with Heavy Duty Armrests
Straight Armpads ..... No Charge Part: UPSASMB1021 Only for use with Heavy Duty Armrests	

**14. LAP BELTS**

**REQUIRED.** Select one lap belt below.

**14.1a. Standard Lap Belts**

<input type="checkbox"/> 38" Lap Belt Part: ACCBELT1006	<input type="checkbox"/> 50" Lap Belt Part: ACCBELT1000	<input type="checkbox"/> 60" Lap Belt Part: ACCBELT1002	<input type="checkbox"/> 70" Lap Belt Part: ACCBELT1003	<input type="checkbox"/> 80" Lap Belt Part: ACCBELT1004	<input type="checkbox"/> 90" Lap Belt Part: ACCBELT1001
--	--	--	--	--	--

**14.1b. Optional Lap Belts**

Stealth Retractable Lap Belt (IFRB5350).....\$175 Part: ST-IFRB5350	Stealth Positioning, 60" Single Pull, Push Button, Extra Large... \$65 Padded (11.5"), 2" Belt (IFPB5315) Part: ST-IFPB5315-Q
Stealth Positioning, 60" Single Pull, Push Button, Medium..... \$65 (7.5") Padded, 1.5" Belt (IFPB5213) Part: ST-IFPB5213-Q	60" Sure-Lok Lap Belt (Non-retractable)..... \$200 Required with ISO 7176-19 Compliant Transit Option. Part: ACCBELT1016
Stealth Positioning, 60" Single Pull, Push Button, Large ..... \$65 Padded (9.5"), 2" Belt (IFPB5314) Part: ST-IFPB5314-Q	

**15. LEGRESTS**

**Optional.** All efforts will be made to accommodate legrest length selection, but minor changes may need to be made to meet shipping guidelines. Some adjustments may be needed upon receipt of the unit. Reference the section's compatibility matrix for legrest interference.

**15.1a. Base Mounted Foot Platform**

May interfere with seat mounted legrests. Not available with center mount foot platforms.

- Base Mounted Foot Platform ..... No Charge  
Part: FRMASMB16077
- Omit base mounted foot platform ..... No Charge  
Part: PTOINDV3674
- Bell Shaped, Base Mounted Foot Platform ..... \$235  
Part: RIGASMB7120041

**15.1b. Basic Center Mount Foot Platform**

**1. Basic Center Mount Foot Platform Assembly**

- Basic Center Mount..... \$200  
Part: XRF000323

Short	Medium	Long
8"-10"	10"-14"	14"-20"
Basic Center Mount Extension Range		

**2. Basic Center Mount Extension and Angle Adjustable Footplates**

- Short Extension, 5"x9" Footplates .....  
Part: XRF000312. HCPCS: K0040
- Long Extension, 5"x9" Footplates .....  
Part: XRF000314. HCPCS: K0040
- Medium Extension, 5"x9" Footplates .....  
Part: XRF000313. HCPCS: K0040

**3. Basic Center Mount Calf Pads**

- Small Calf Pads, Basic Center Mount..... \$210  
Part: RIG2104369
- Large Calf Pads, Basic Center Mount ..... \$210  
Part: RIG2202997  
Only available with Long Extension.

**4. Basic Center Mount Additional Accessory**

- Basic Center Mount Footplate Connector ..... \$55  
Part: RIG2104240

**15.1c. Center Mount Foot Platform**

If power articulating foot platform is selected, be sure power legrest control has been selected in earlier section. Select legrest length in following section.

**1. Center Mount Foot Platform Type**

- Power Articulating Foot Platform ..... \$3,145  
Part: PTOINDV3754. HCPCS: E1012+K0040
- Quick Adjustable Center Mount Foot Platform..... \$265  
Part: XRF000311. HCPCS: K0040

**2. Extension and Angle Adjustable Footplates**

Pediatric extensions are not compatible with Power Articulating Foot Platform.

- Pediatric Extension, 5"x6" Footplates .....  
Part: RIG152577
- Medium Extension, 5"x11" Footplates .....  
Part: RIG150212
- Pediatric Extension, 5"x8" Footplates .....  
Part: RIG152578
- Long Extension, 5"x6" Footplates .....  
Part: RIG151505
- Medium Extension, 5"x6" Footplates .....  
Part: RIG150206
- Long Extension, 5"x8" Footplates .....  
Part: RIG151506
- Medium Extension, 5"x8" Footplates .....  
Part: RIG150203
- Long Extension, 5"x11" Footplates .....  
Part: RIG151508

Medium	Long
11"-17"	14"-20"
Power AFP Extension Range	

Pediatric	Medium	Long
8.5"-12.5"	11"-16.5"	14"-20"
Quick Adj Center Mount Extension Range		

**3a. Center Mount Calf Pads**

Calf pads are standard with a power articulating foot platform.

- Small Calf Pads, AFP Only, 6"Wx6"H (pair) ..... Standard  
Part: RIGASMB7120018
- Large Calf Pads, AFP Only, 8" taper to 6"Wx8"H (pair) ... No Charge  
Part: RIGASMB7120019
- Small Calf Pads, Manually Adjustable Center Mount Only, ..... \$210  
6"Wx6"H (pair)  
Part: XRF000156
- Large Calf Pads, Manually Adjustable Center Mount Only, ..... \$210  
8" taper to 6"Wx8"H (pair)  
Part: XRF000157

**3b. Center Mount Calf Panel**

Calf pads are standard with a power articulating foot platform.

- Small Calf Panel 18"Wx8"H..... \$210  
Part: FRMASMB7110012
- Large Calf Panel 24"Wx8"H ..... \$265  
Part: RIG123824  
Only available on widths of 20" or greater.

**4. Center Mount Modular Footplate Extension Components**

Optional. Assembly required. May require seating configuration adjustments to avoid front caster interference.

- Width Extension for 8"D Footplates ..... \$55  
Part: RIG152517
- 4" Outside Edge (pair) ..... \$55  
Part: RIG152510  
Only available with 6"D footplates and 8"D footplate extension.
- Width Extension for 11"D Footplates ..... \$55  
Part: RIG152518
- 6" Outside Edge (pair) ..... \$55  
Part: RIG152516  
Available with 8"D and 11"D footplates and 11"D footplate extension.
- 2" Bell-Shaped Length Extension ..... \$55  
Part: RIG152519  
Not compatible with width extensions. Can be used with 6" outside edge with 11" footplates.

Legrests section continued on following pages.

**15. LEGRESTS CONTINUED**

**5. Center Mount Additional Accessories**

- Small Reversed Dartex Calf Pad Cover (pair) ..... \$150  
Part: RIG133717
- Large Reversed Dartex Calf Pad Cover (pair) ..... \$150  
Part: RIG133719
- Small Reversed Dartex Panel Cover (pair)..... \$150  
Part: RIG134673
- Large Reversed Dartex Panel Cover (pair) ..... \$150  
Part: RIG134744
- Small Foam Calf Pad Cover (pair)..... \$150  
Part: ST-CSCS0775-Q
- Large Foam Calf Pad Cover (pair) ..... \$150  
Part: ST-CSCS0985-Q
- Small Foam Panel Cover (pair) ..... \$150  
Part: ST-CSCP0920-Q
- Large Foam Panel Cover (pair)..... \$150  
Part: ST-CSCP0927-Q
- Footplate Connector ..... \$55  
Part: RIGASMB7120023
- Footplate Riser ..... \$150  
Part: XRF000348  
Not available with Long Extension. Select legrest length between 2" to 9" at 1" increments.

**6. Stealth Gel Footplate Covers**

- 5"x8" Tapered Gel Footplate Covers ..... \$275  
Part: ST-FPG5008T-Q. HCPCS: K0108
- 5"x11" Tapered Gel Footplate Covers..... \$275  
Part: ST-FPG5011T-Q. HCPCS: K0108
- 5"x8" Tapered Gel Footplate Covers for Width Extension ..... \$275  
Part: ST-FPG6508T-Q. HCPCS: K0108
- 5"x11" Tapered Gel Footplate Cover for Width Extension ..... \$275  
Part: ST-FPG7011T-Q. HCPCS: K0108

**15.1c. Swing-Away Legrests**

Select a legrest upper, a lower extension, and type of footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes.

**1a. 60°-90° Swing-Away Legrest Uppers**

Adjustable in 5° increments.

- 60°-90° Swing-Away Legrest** ..... \$285  
Pair part: RIG120062. Left: FRMASMB17558.  
Right: FRMASMB17559  
Not available with 11"-13" extensions.
- 60°-90° Swing-Away Legrest, Short** ..... \$285  
Pair part: RIG136481. Left: RIG132607.  
Right: RIG132609.  
For use only with 11"-13" lower extensions.

**2. Lower Extension**

- 11"-13" Length**  
Pair part: XRF000088.  
Left: XRF000089.  
Right: XRF000090.  
N/A with Multi-axis.
- 11"-16.5" Length**  
Pair part: XRF000091.  
Left: XRF000092.  
Right: XRF000093.  
N/A with Angle Adjustable.
- 13"-16.5" Length**  
Pair part: XRF000094.  
Left: XRF000095.  
Right: XRF000096.  
N/A with Multi-Axis.
- 15.5"-19.5" Length**  
Pair part: XRF000097.  
Left: XRF000098.  
Right: XRF000099.

**3a. Extruded Footplate(s)**

- 7.25"Wx6"D Extruded Footplate**..... Standard  
Part pair: FRMASMB7847. Left: FRMASMB7861  
Right: FRMASMB7862
- 8.5"Wx6"D Extruded Footplate** ..... Standard  
Part pair: FRMASMB7848. Left: FRMASMB7863  
Right: FRMASMB7864  
May cause interference on seat widths less than 18."

60°-90° Swing-Away Legrest Compatibility Matrix	Contoured/Synergy Seat			
	Low	Medium	High	
Extruded, 13" - 16.5"	!	✓	✓	<input checked="" type="checkbox"/> This rigging configuration causes no interference issues at this seat-to-floor height. <input type="checkbox"/> This rigging configuration may have less than 1.75" clearance when fully extended and may interfere with the ground or reduce obstacle climbing ability.
Extruded, 15.5" - 19.5"	!	!	!	
Angle Adjustable, 13" - 16.5"	!	✓	✓	<input checked="" type="checkbox"/> This rigging configuration does not work with this seat-to-floor height.
Angle Adjustable, 15.5" - 19.5"	!	!	!	
Multi-Axis, 15.5" - 19.5"	!	!	!	

**3b. Angle Adjustable Footplate(s) with Heel Loops**

<sup>1</sup>The 6"x11" footplates may cause interference with seat widths less than 18". <sup>2</sup>The 8"x11" footplates may cause interference with seat widths less than 20".

- 4"Wx6"D** ..... \$87.50 ea.  
Pair part: FRMASMB10564. Left: RIG133537. Right: RIG133560
- 5"Wx6"D**..... \$87.50 ea.  
Pair part: FRMASMB10565. Left: RIG133538. Right: RIG133561
- 5"Wx8"D** ..... \$87.50 ea.  
Pair part: FRMASMB10566. Left: RIG133539. Right: RIG133562
- 6"Wx8"D<sup>1</sup>** ..... \$87.50 ea.  
Pair part: FRMASMB10567. Left: RIG133540. Right: RIG133563
- 6"Wx11"D<sup>1</sup>** ..... \$87.50 ea.  
Pair part: FRMASMB10568. Left: RIG133541. Right: RIG133564
- 8"Wx11"D<sup>2</sup>** ..... \$87.50 ea.  
Pair part: FRMASMB10569. Left: RIG133542. Right: RIG133565

**3c. Multi-Axis Footplate(s) with Heel Loops**

<sup>1</sup>The 6"x11" footplates may cause interference with seat widths less than 18". <sup>2</sup>The 8"x11" footplates may cause interference with seat widths less than 20".

- 4"Wx6"D** ..... \$260 ea.  
Pair part: FRMASMB9456. Left: RIG133531. Right: RIG133552
- 5"Wx6"D**..... \$260 ea.  
Pair part: FRMASMB9457. Left: RIG133532. Right: RIG133553
- 5"Wx8"D** ..... \$260 ea.  
Pair part: FRMASMB9458. Left: RIG133533. Right: RIG133554
- 6"Wx8"D<sup>1</sup>** ..... \$260 ea.  
Pair part: FRMASMB9459. Left: RIG133534. Right: RIG133555
- 6"Wx11"D<sup>1</sup>** ..... \$260 ea.  
Pair part: FRMASMB9460. Left: RIG133535. Right: RIG133558
- 8"Wx11"D<sup>2</sup>** ..... \$260 ea.  
Pair part: FRMASMB9461. Left: RIG133536. Right: RIG133559

Legrests section continued on following page.



**15. LEGRESTS CONTINUED**

**15.1d. Pediatric, Swing-Away Legrests**

Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes. Select legrest length in following section.

**1. Pediatric Legrests Upper**

**90°, Pediatric, Swing-Away Legrest (7"-10" Length Extension) ..... \$285**  
 Pair part: FRMASMB8808. Left: FRMASMB8809  
 Right: FRMASMB8810

**2. Angle Adjustable Footplate(s)**

**4"Wx6"D Angle Adjustable ..... \$87.50 ea.**         **5"Wx6"D Angle Adjustable ..... \$87.50 ea.**  
 Pair part: FRMASMB8812. Left: RIG133537. Right: RIG133560      Pair part: FRMASMB8811. Left: RIG133538. Right: RIG133561  
 HCPCS: K0040      HCPCS: K0040

**15.1e. Elevating Legrests**

Select a legrest upper, a legrest extension, and type of footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes. Select legrest length in following section.

**1. Elevating Legrests Upper**

**Elevating Legrest Assembly ..... \$275 ea.**  
 Pair: FRMASMB6483. Left: RIG138012  
 Right: RIG138013. HCPCS: E0990  
 Select lower extension and footplates

**11"-14" Length**  
 Pair: FRMASMB4755  
 Left: RIG138026  
 Right: RIG138030  
 For reinforced only

**14"-17" Length**  
 Pair: FRMASMB4756  
 Left: RIG138033  
 Right: RIG138035  
 For reinforced only

**12"-14.5" Length**  
 Pair part: XRF000101.  
 Left: XRF000102.  
 Right: XRF000103.  
 For angle adjustable only

**14.5"-18" Length**  
 Pair part: XRF000104.  
 Left: XRF000105.  
 Right: XRF000106  
 For angle adjustable only

Elevating Legrest Compatibility Matrix	Contoured/Synergy Seat		
	Low	Medium	High
Standard, 11"-14"	✓	✓	✓
Standard, 14"-17"	!	!	!
Angle Adjustable, 12"-14.5"	✓	✓	✓
Angle Adjustable, 14.5"-18"	!	!	!

This rigging configuration causes no interference issues at this seat-to-floor height.  
 This rigging configuration may have less than 1.75" clearance when fully extended and may interfere with the ground or reduce obstacle climbing ability.

**3a. Reinforced Cast Footplate(s)**

**7.25"Wx6"D Reinforced Footplate ..... Standard**  
 Pair part: FRMASMB4748. Left: RIG138055. Right: RIG138056

**8.5"Wx6"D Reinforced Footplate ..... Standard**  
 Pair part: FRMASMB4747. Left: RIG138060. Right: RIG138062  
 May cause interference on seat widths less than 18."

**3b. Angle Adjustable Footplate(s)**

\*The 6"x11" footplates may cause interference with seat widths less than 18". \*The 8"x11" footplates may cause interference with seat widths less than 20". 4"x6" and 5"x6" footplates are only available in pairs.

**4"Wx6"D ..... \$87.50 ea.**  
 Pair parts: FRMASMB8812 x2.

**6"Wx8"D<sup>1</sup> Footplates with Heel Loops ..... \$87.50 ea.**  
 Pair part: FRMASMB4738. Left: RIG137793. Right: RIG137795

**5"Wx6"D ..... \$87.50 ea.**  
 Pair parts: FRMASMB8811 x2.

**6"Wx11"D<sup>1</sup> Footplates with Heel Loops ..... \$87.50 ea.**  
 Pair part: FRMASMB4741. Left: RIG137796. Right: RIG137797

**5"Wx8"D Footplates with Heel Loops ..... \$87.50 ea.**  
 Pair part: FRMASMB4735. Left: RIG137791. Right: RIG137792

**8"Wx11"D<sup>2</sup> Footplates with Heel Loops ..... \$87.50 ea.**  
 Pair part: FRMASMB4744. Left: RIG137798. Right: RIG137799

**15.2. Legrest Accessories**

**Heel Loops (pair) ..... \$35**  
 Part: FRMASMB7873. HCPCS: E0951  
 Standard on angle adjustable and multi axis footplates when noted.

**Wheel Bumper ..... \$15 ea.**  
 Part: XRF000100

**15.3. Set Legrest Extension Length**

**REQUIRED with matching left and right legrest configuration. Part number is PTOINDV plus number following # below.** All efforts will be made to accommodate legrest length selection, but minor changes may need to be made to meet shipping guidelines. Some adjustments may be needed upon receipt of the unit. Reference the section's compatibility matrix for legrest interference. Non-matching or individual legrest configurations lengths are not set by the manufacturer.

7" L (#3470)	8.5" L (#3473)	10" L (#3476)	11.5" L (#3479)	13" L (#3482)	14.5" L (#3485)	16" L (#3488)	17.5" L (#3491)	19" L (#3494)	20.5" L (#3497)	22" L (#3669)
7.5" L (#3471)	9" L (#3474)	10.5" L (#3477)	12" L (#3480)	13.5" L (#3483)	15" L (#3486)	16.5" L (#3489)	18" L (#3492)	19.5" L (#3495)	21" L (#3498)	22.5" L (#3670)
8" L (#3472)	9.5" L (#3475)	11" L (#3478)	12.5" L (#3481)	14" L (#3484)	15.5" L (#3487)	17" L (#3490)	18.5" L (#3493)	20" L (#3496)	21.5" L (#3600)	

**16. POWER LIFT ELECTRONICS**

**REQUIRED with 12" Power Lift for Captain Seat with iLevel.**

Single Actuator Control through NE+ Joystick ..... No Charge  
 Part: ELE2006331

Single Actuator Control through Q-Logic Joystick ..... No Charge  
 Part: ELE2006330

**17. ACCESSORIES**

**Optional.** Multiple accessories can be mounted along the various tracks but some interference may occur. Accessories may need to be installed by the provider.

<b>Cup Holder, C Style</b> ..... \$20 <i>Part: ACCASMB2288</i> For use with Flip-Up armrests.	<b>Quantum Backpack</b> ..... \$200 <i>Part: ACC1905853</i>
<b>Cup Holder, Bariatric Arms or Gel Armpads</b> ..... \$20 <i>Part: ACCASMB2632</i>	<b>Reflector Kit</b> ..... \$35 <i>Part: ACC124670</i>
<b>Cup Holder, Tube Style Armrests</b> ..... \$20 <i>Part: XRF000221</i> For HD Removable, Single Post, or Pediatric Arms.	<b>Clothing Guard (Pair)</b> ..... \$275 <i>Part: ACC125007</i> Available with Flip-Up armrests only. Not available with Elevating Legrests.
<b>Oxygen Holder</b> ..... \$110 <i>Part: ACC2006353. HCPCS: E2208</i>	<b>Transfer Handles</b> ..... \$305 <i>Part: MEC1809930. HCPCS: E0705</i> Available with Flip-Up armrests only. Not available with Elevating Legrests.
<b>Cane/Crutch Holder</b> ..... \$195 <i>Part: ACC2006352. HCPCS: E2207</i>	<b>Privacy Bag</b> ..... \$80 <i>Part: ACC1907570</i>
<b>Rear Basket</b> ..... \$200 <i>Part: ACC2006351</i>	
<b>Saddle Bag</b> ..... \$15 <i>Part: DWR1020D010</i>	
<b>Large, Black, Weather Cover</b> ..... \$240 <i>Part: ACCCOVR1021</i>	

**18. PROGRAMMER**

**Optional.**

<input type="checkbox"/> <b>Q-Logic EX/NE Programmer</b> ..... \$625 <i>Part: ELEASMB6976</i>
--