



**Pride Mobility Products Co.**

5096 South Service Rd, Beamsville, ON, L0R 1B3

Phone: 888-570-1113 | Fax: 866-514-1303 | Email: productorders@pridemobility.com

# Quantum 4Front<sup>®</sup> 2 Order Form with TRU-Balance Flex

Account Number: \_\_\_\_\_ Date: \_\_\_\_\_

Provider Name: \_\_\_\_\_

Contact: \_\_\_\_\_

Phone: \_\_\_\_\_ Fax: \_\_\_\_\_

Email: \_\_\_\_\_

PO Number: \_\_\_\_\_

Marked for: \_\_\_\_\_

Ship to Address: \_\_\_\_\_

City: \_\_\_\_\_ Prov: \_\_\_\_\_ PC: \_\_\_\_\_



## INTRODUCTION

This form is interactive when viewed with [Adobe Acrobat Reader](#). Complete the form by placing checks in the desired boxes and provide information in the interactive fields. Buttons shown at the bottom of the form may be utilized to print or submit the order form through a desktop email application. To email via a web-based application, please 'Save As' and attach the PDF to your email.

This order form contains a large variety of options to fulfill various patient needs. Descriptions and section notes, such as optional or required, are included to help you complete your order. Other part numbers can be found on the quote acknowledgment after this order form is submitted. Please contact Quantum<sup>®</sup> Sales at 888-570-1113 if assistance is needed. Send the completed order form by fax (866-514-1303) or email ([productorders@pridemobility.com](mailto:productorders@pridemobility.com)). Incomplete forms may delay the quote or order. Customer service will contact you if the order is incomplete or if there are compatibility issues. If special order requests are needed, be sure that the Patient Information section is completed

All prices are MSRP. Prices, specifications, part numbers, and availability are subject to change without notice. Prices and part numbers as shown are only available when configured as a complete power chair. Please contact technical service for accurate parts ordering.


**Thank you for choosing Quantum!**

## PATIENT INFORMATION

Weight: \_\_\_\_\_ Height: \_\_\_\_\_ A. Shoulder Width: \_\_\_\_\_ D. Max Sitting Height: \_\_\_\_\_ G. Thigh Depth: \_\_\_\_\_

B. Chest Width: \_\_\_\_\_ E. Shoulder Height: \_\_\_\_\_ H. Lower Leg Length: \_\_\_\_\_

C. Hip Width: \_\_\_\_\_ F. Axilla Height: \_\_\_\_\_ I. Elbow Height: \_\_\_\_\_



### 1. BASE MODEL

**REQUIRED.** Select a group 3 power base. 300 lb. weight capacity with standard 6 mph programmable motors and 8 amp off board charger. Dimensions of base without legrests: 24.2"Wx39.9"L

Quantum 4Front 2 .....\$8,450  
Part: 4FRONT 2

### 1A. AUTO CORRECTION SYSTEM







**REQUIRED.** Smart Traction Control is auto-selected with every 4Front 2.

Smart Traction Control ..... \$400  
Part: Various

**2. BASE COLOR**

**REQUIRED.** Select one shroud color. Colors shown here may differ from actual product depending on monitor/printer calibration used.

**2.1a. Standard Shrouds**

- |   |  |  |
|---|--|--|
|  <b>Red Corvette</b><br>Part: DGN1900199-PRD1  |  <b>Raspberry Beret</b><br>Part: DGN1900199-PPK4  |  <b>Blue Jean</b><br>Part: DGN1900199-PBU11     |
|  <b>Back in Black</b><br>Part: DGN1900199-PBK1 |  <b>Pretty in Pink</b><br>Part: DGN1900199-PPK3   |  <b>Lemon Crush</b><br>Part: DGN1900199-PYL4    |
|  <b>Night Train</b><br>Part: DGN1900199-PBK3   |  <b>Evolution Orange</b><br>Part: DGN1900199-POR7 |  <b>Green Machine</b><br>Part: DGN1900199-PGN11 |
|  <b>Grey Street</b><br>Part: DGN1900199-PSL3   |  <b>Midnight Blue</b><br>Part: DGN1900199-PBU10   |  <b>Military Green</b><br>Part: DGN1900199-PGN8 |
|  <b>Deep Purple</b><br>Part: DGN1900199-PPR5   |  <b>Electric Blue</b><br>Part: DGN1900199-PBU9    |  |

**2.1b. Hydrographic Shrouds**



**Military Digital Camo** ..... \$325  
Part: DGN2003028-WCM7

**3. BASE MOUNTED LIGHTS**

**Optional.** Select one option. Seat Mounted, LED, Full Light Package option can be found later on this form.

**Switch Operated Fender Lights**..... No Charge  
Part: DGN2004686  
Fender lights operate by pressing included single switch button. Fender lighting is not available with Single Actuator Control through Toggle Box.

**Base Mounted, LED, Full Light Package** ..... \$770  
Part: ELE1700035  
Requires Q-Logic 3 EX Joystick for Full Lights or Q-Logic 3e Joystick.

**Omit Fender Lights** ..... No Charge  
Part: DGN1700032

**4. QUANTUM BACKUP CAMERA**

**Optional.** The backup camera may interfere with armrest mounted accessories such as media device holder.

**Quantum Backup Camera** ..... \$475  
Part: MEC1910365

**1. Optional Mounting Bracket**

**Backup Camera Headrest Bracket**..... \$30  
Part: MEC1910373  
Works with .5" x .5" vertical supports only.

**Backup Camera Track Mount Bracket** ..... \$35  
Part: MEC2104227

**2. Optional Display Mounting Bracket**

**Backup Camera Q-Logic Joystick Mount Bracket** ....\$145  
Part: MEC2106123



Backup Camera

**5. DRIVE WHEELS**

**REQUIRED.** Select a type of drive wheels and accent color.

**5.1. Drive Wheels**

**3"x14" Flat Free Black Tires with Accent** ..... Standard  
Part: PTOINDV3525

**3"x14" Pneumatic Black Tires with Accent** ..... No Charge  
Part: WHL1704059 x2

**4"x14" Flat Free Knobby Black Tires**..... \$60 ea.  
Part: WHL2006647 x2

**4"x14" Pneumatic Knobby Black Tires**..... \$60 ea.  
Part: WHL2006822 x2

**5.2. Drive Wheel Accent Color**

Not available with knobby tires.

**Back in Black** .....Standard  
Part: DGN164245-PBK1 x2

**Red Corvette** ..... No Charge  
Part: DGN164245-PRD1 x2

**Night Train** ..... No Charge  
Part: DGN164245-PBK3 x2

**Grey Street** ..... No Charge  
Part: DGN164245-PSL3 x2

**Deep Purple** ..... No Charge  
Part: DGN164245-PPR5 x2

**Raspberry Beret** ..... No Charge  
Part: DGN164245-PPK4 x2

**Pretty in Pink** ..... No Charge  
Part: DGN164245-PPK3 x2

**Evolution Orange** ..... No Charge  
Part: DGN164245-POR7 x2

**Midnight Blue** ..... No Charge  
Part: DGN164245-PBU10 x2

**Electric Blue** ..... No Charge  
Part: DGN164245-PBU9 x2

**Blue Jean** ..... No Charge  
Part: DGN164245-PBU11 x2

**Lemon Crush** ..... No Charge  
Part: DGN164245-PYL4 x2

**Green Machine** ..... No Charge  
Part: DGN164245-PGN11 x2

**Military Green** ..... No Charge  
Part: DGN164245-PGN8 x2

**Military Digital Camo** ..... \$50 ea.  
Part: DGN164245-WCM7 x2

**6. BATTERIES**

**Optional.** The powerchair requires 2 flat top batteries and will be installed if ordered.

**Group 24, AGM Cell Battery (MK)**.....\$510 ea.  
Part: BAT1712207.

**Group 34, Gel Cell Battery (MK)** ..... \$400 ea.  
Part: BATGEL1006.  
Not available with TRU-Balance 3 Power Lift.

**7. PATIENT'S WEIGHT**

**REQUIRED.** The patient's weight is required to quote and order the power chair.

Up to 200 lb User Weight  
Part: PTOINDV3837

201 lb-300 lb User Weight  
Part: PTOINDV3838

**8. ELECTRONICS**

**REQUIRED.** Select a type of drive control from the following sections. The No Joystick option is used to omit the joystick and utilize another control, such as a specialty control at the end of this form or found on the [Alternative Controls and Electronics Order Form](#).

**8.1a. Omit Joystick**

No Joystick, Setup for Q-Logic Electronics  
Part: XRF000010

**8.1c. Q-Logic 3 Expandable Capacity Controls**

Q-Logic 3 expandable capacity controls allow use of specialty controls and multiple actuators.

**1. Q-Logic 3 Expandable Capacity Joysticks**

Q-Logic 3 Expandable Capacity Joystick ..... \$900  
Part: CTL167724

Q-Logic 3 Expandable Capacity Joystick ..... \$900  
for Full Light Package  
Part: CTL167725  
Required with full light packages.

Q-Logic 3e Joystick ..... \$620  
Part: CTL1910600

Can operate full light package. Not available with alternative or specialty controls.



**8.2. Q-Logic 3 Joystick Armor**

**REQUIRED with Q-Logic 3 Joystick selections above. Not available with Q-Logic 3e Joystick.** Reference colors shown on page 2.

Black, Soft Touch ..... No Charge  
Part: DGN124307

Match Joystick Armor to Base Color ..... Standard  
Part: PTOINDV3421

Non-matching Joystick Armor ..... No Charge  
Part: PTOINDV3422  
Select color below

Part below will be selected with system configurator.

Red Corvette  
Part: DGN124336-PRD1

Deep Purple  
Part: DGN124336-PPR5

Midnight Blue  
Part: DGN124336-PBU10

Green Machine  
Part: DGN124336-PGN11

Back in Black  
Part: DGN124336-PBK1

Raspberry Beret  
Part: DGN124336-PPK4

Electric Blue  
Part: DGN124336-PBU9

Military Green  
Part: DGN124336-PGN8

Night Train\*  
See Black, Soft Touch

Pretty in Pink  
Part: DGN124336-PPK3

Blue Jean  
Part: DGN124336-PBU11

Military Digital Camo  
Part: DGN124336-WCM7

Grey Street  
Part: DGN124336-PSL3

Evolution Orange  
Part: DGN124336-POR7

Lemon Crush  
Part: DGN124336-PYL4

Carbon Fiber  
Part: DGN124336-WCB1

**9. JOYSTICK HANDLES**

**Optional.** The standard Q-Logic 3 joystick knob can be seen in an image in the previous section. An optional handle selected below will be mounted to the joystick and the standard handle will be shipped with the unit.

Mushroom Handle (B) ..... \$120  
Part: ACCASMB1860.

Extended I-Handle (G) ..... \$140  
Part: ACCASMB1857.

Chin Cup (C) ..... \$85  
Part: ACCASMB186190

Bodypoint 3" Goal Post Handle ..... \$180  
Part: PLSKNCP1057.

T-Bar (D) ..... \$160  
Part: ACCASMB1858.

Bodypoint 4" Goal Post Handle ..... \$180  
Part: PLSKNCP1058.

Medium Ball Handle (E) ..... \$70  
Part: ACCASMB1859.

Soft Ball Handle (F) ..... \$120  
Part: ACCASMB1856.



**10. JOYSTICK MOUNTING BRACKETS**

**Optional.** Select a mount to be used with a joystick selected in previous section.

**10.1a. Fixed, In-line Joystick Bracket**

Fixed, In-line, Bracket ..... Standard  
Part: FRMASMB10617



**10.1b. Swing-Away (Retractable) Joystick Bracket**

Allows the joystick to retract to the side of the armrest while maintaining the same joystick orientation.

Adjustable Swing-away ..... \$265  
Part: Various.  
As pictured to the right "Swing Away Joystick"

Offset Swing-away ..... \$265  
Part: Various.

Basic Swing-away ..... \$265  
Part: Various.  
Magnetic latching mechanism.  
Stealth Height Adjustable Receiver recommended.

As pictured to the right "Swing Away Joystick Offset Mount (Right)"

Left Side Mount  
Part: PTOINDV1225

Right Side Mount  
Part: PTOINDV1226



**10. JOYSTICK MOUNTING BRACKETS**

**Optional.** Select a mount to be used with a joystick selected in previous section.

**10.1c. Specialty Joystick Mounts**

Select a bracket and side to mount it.

**Retract4 Mount**.....\$450

Part: XRF000003.

Full 180° of movement that translates into multi-position use, straight, right offset, left offset, or fold away. Joystick orientation varies as bracket is adjusted.

**Multi-Axis Mount**.....\$100


Part: XRF000004

Allows height adjustment of joystick.

**Stealth Swing-Away/Flip Down Bracket**.....\$300

Part: XRF000005.

Allows the joystick to pivot to side of armrest, or flip-down if bracket is adjusted.



Retract4



Multi-Axis Joystick Mount

**10.2. Joystick Mount Location**

**Left Side Mount**  
Part: PTOINDV1225

**Right Side Mount**  
Part: PTOINDV1226

**11. TRU-BALANCE FLEX SEATING SYSTEM**

**Select one option from each section.** Seat width and depth can be adjusted in 1" increments from 15"W-22"W and 16"D-20"D. Maximum seat size 22"W x 20"D.

**11.1. TRU-Balance Flex Frame Selection**

**TRU-Balance Flex Seating System**.....\$650  
Part: PTOINDV3853

**11.2. Occupied Transit Kit**

**Optional.** The power base is standard with Unoccupied Transit Loops. 60" Sure-Lok Lap Belt is required with Occupied Transit Option. To determine if your chair is WC-19 compliant (occupied), please refer to this website: <http://wc-transportation-safety.umtri.umich.edu/crash-tested-product-lists/wheelchairs>

**WC-19 Occupied Transit Option**.....\$400  
Part: KIT152047

**11.3. Width**

Maximum seat width is 22". Part number is SETWDTH plus number following # below.

15" W (#1005)	16" W (#1006)	17" W (#1007)	18" W (#1008)	19" W (#1009)	20" W (#1010)	21" W (#1011)	22" W (#1012)
------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------

**11.4. Depth**

Maximum seat depth is 20". Part number is SETDPTH plus number following # below.

15" D (#1005)	16" D (#1006)	17" D (#1007)	18" D (#1008)	19" D (#1009)	20" D (#1010)
------------------	------------------	------------------	------------------	------------------	------------------

High	18.5" at front with 4" posterior angle
Med	17.5" at front with 2" posterior angle
Low	16.5" at front and rear
TRU-Balance Flex Seat to Floor Heights	

**11.5. Seat to Floor Height**

**Lowest Seat to Floor**  
Part: PTOINDV3558

**Middle Seat to Floor**  
Part: PTOINDV3559

**Highest Seat to Floor**  
Part: PTOINDV3560

**11.6. Matching Seat Cushion**

**Optional. Specific cushion/component sizes available, part numbers, weight capacities, MSRP, HCPCS code, and special order cushions can be found on the respective complete [Stealth Cushion Order Forms](#).** This section allows a quick selection of a single cushion that matches the seat pan width and depth. See specific order form for non-matching sizes, alternate covers and liners, or additional covers and liners.

**Stealth Simplicity G Cushion**

General use cushion with Stealth Tek 2, 2-way stretch, breathable cover.  
Part series: ST-SIMGXXYY. HCPCS: E2601/E2602


**Simplicity G, Matching Seat Size**  
Reference part: XRF000111. MSRP: \$110/\$165



**Stealth Glacial SP Cushion**

Skin protection with Coolcore, 4-way stretch, breathable cover.  
Part series: ST-GLSPXXYY. HCPCS: E2603/E2604

**Glacial SP, Matching Seat Size**  
Reference part: XRF000114. MSRP: \$330/\$440



**Stealth Zen SP Cushion**

Skin protection with Coolcore, 4-way stretch, breathable cover.  
Part series: ST-ZENSPXXYY. HCPCS: E2603/E2604

**Zen SP, Matching Seat Size**  
Reference part: XRF000115. MSRP: \$315/\$425



**Stealth Premiere P Cushion**

Positioning with Coolcore, 4-way stretch, breathable cover.  
Part series: ST-PREMPXXYY. HCPCS: E2605/E2606

**Premiere P, Matching Seat Size**  
Reference part: XRF000122. MSRP: \$325/\$475




TRU-Balance Flex Seating System section continued on following pages.

**11. TRU-BALANCE FLEX SEATING SYSTEM CONTINUED**

**Stealth Spectrum Foam SPP Cushion**

Skin protection and positioning with Coolcore, 4-way stretch, breathable cover. Includes urethane liner.  
Part series: ST-SPCFSPXXYY. HCPCS: E2607/E2608


Spectrum Foam SPP, Matching Seat Size  
Reference part: XRF000112. MSRP: \$460/\$530



**Stealth Solution SPP Cushion**

Skin protection and positioning with Coolcore, 4-way stretch, breathable cover. Includes urethane liner.  
Part series: ST-SOLSPPXXYY. HCPCS: E2607/E2608

Solution SPP, Matching Seat Size  
Reference part: XRF000113. MSRP: \$375/\$430



**Stealth Essence SPP Cushion**

Skin protection and positioning with Coolcore, 4-way stretch, breathable cover. Part series: ST-ESSSPXXYY. HCPCS: E2607/E2608

Essence SPP, No Wedge, Matching Seat Size  
Reference part: XRF000118. MSRP: \$400/\$500

Essence SPP, Small Wedge, Matching Seat Size  
Reference part: XRF000119. MSRP: \$435/\$535

Essence SPP, Medium Wedge, Matching Seat Size  
Reference part: XRF000120. MSRP: \$465/\$565

Essence SPP, Large Wedge, Matching Seat Size  
Reference part: XRF000121. MSRP: \$495/\$595

**11.7. TRU-Balance Flex Back Cane Selection**

Angle Adjustable Back Cane .....\$195  
Part: MEC140266

Flip Forward Back Cane .....\$195  
Part: MEC140272

Not available with oxygen holder, or cane mounted and cantilever armrests.

2" Mesomorph (1" per side) Back Cane ..... \$205  
Part: MEC1804736

4" Mesomorph (2" per side) Back Cane..... \$205  
Part: MEC1804739

2" Endomorph (1" per side) Back Cane ..... \$205  
Part: MEC141455

4" Endomorph (2" per side) Back Cane..... \$205  
Part: MEC141456

24"H Flip Forward Flex Cane ..... Standard

**11.8. TRU-Balance Flex Back Cane Height Selection**

16" and 17" back canes are not available with a Synergy Back and cane mounted, flip up armrests.

14" High Canes  
Part: FRMASMB15832

17" High Canes  
Part: FRMASMB15835

20" High Canes  
Part: FRMASMB15838

23" High Canes  
Part: FRMASMB15841

15" High Canes  
Part: FRMASMB15833

18" High Canes  
Part: FRMASMB15836

21" High Canes  
Part: FRMASMB15839

24" High Canes  
Part: FRMASMB15842

16" High Canes  
Part: FRMASMB15834

19" High Canes  
Part: FRMASMB15837

22" High Canes  
Part: FRMASMB15840

25" High Canes  
Part: FRMASMB15843

**11.9. TRU-Balance Flex Back Cushion Options**

Select one type of back. See the seat back availability chart for stock order back size ranges. Other size backs may be available with special order and extended lead time. Select a back type and overall back height for configured part number.

		Back Width										
		12"	13"	14"	15"	16"	17"	18"	19"	20"	21"	22"
Overall Height	13"	-	-	S	-	-	-	-	-	-	-	-
	14"	P T	-	S P	-	-	-	-	-	-	-	-
	15"	S P T	-	S P	-	-	-	-	-	-	-	-
	16"	S P T	-	S P T	-	P T	P T	P T	P T	P T	P T	P T
	17"	P T	S	S P T	S	S P T	S P T	S P T	S P T	S P T	P T	S P T
	18"	-	S	S	S	S P T	S P T	S P T	S P T	S P T	P T	S P T
	19"	-	-	-	S	S P T	P T	S P T	S P T	S P T	S P T	S P T
	20"	-	-	-	S	S P T	P T	S P T	S P T	S P T	S P T	S P T
	21"	-	-	-	S	S P T	S P T	S P T	S P T	S P T	S P T	S P T
	22"	-	-	-	S	S	S P T	S P T	S P T	S P T	S P T	S P T
	23"	-	-	S	S	S P T	S P T	S P T	S P T	S P T	S P T	S P T
24"	-	-	S	S	S P	S P	S P T	S P	S P T	S P	S P	
25"	-	-	S	S	S P	S P	S P T	S P	S P T	S P	S P	

S Synergy Back   
 P Sport Back   
 T TRU-Comfort Back

**TRU-Balance Flex Seat Back Availability Chart**

TRU-Balance Flex Seating System section continued on following pages.

**11. TRU-BALANCE FLEX SEATING SYSTEM CONTINUED**

**1a. Synergy Back Selection**

Synergy Back Cushion ..... \$940  
 Select size below. Height is overall. Part number is SETASMB plus number following #.

12"Wx15"H Part: #6448	14"Wx21"H Part: #6491	15"Wx22"H Part: #6668	16"Wx23"H Part: #6512	17"Wx24"H Part: #6689	18"Wx25"H Part: #6533	20"Wx17"H Part: #6467	21"Wx20"H Part: #6655	22"Wx21"H Part: #6499
12"Wx16"H Part: #6624	14"Wx22"H Part: #6667	15"Wx23"H Part: #6511	16"Wx24"H Part: #6688	17"Wx25"H Part: #6532	19"Wx17"H Part: #6466	20"Wx18"H Part: #6643	21"Wx21"H Part: #6498	22"Wx22"H Part: #6675
13"Wx17"H Part: #6460	14"Wx23"H Part: #6510	15"Wx24"H Part: #6687	16"Wx25"H Part: #6531	18"Wx17"H Part: #6465	19"Wx18"H Part: #6642	20"Wx19"H Part: #6478	21"Wx22"H Part: #6674	22"Wx23"H Part: #6518
13"Wx18"H Part: #6636	14"Wx24"H Part: #6686	15"Wx25"H Part: #6530	17"Wx17"H Part: #6464	18"Wx18"H Part: #6641	19"Wx19"H Part: #6477	20"Wx20"H Part: #6654	21"Wx23"H Part: #6517	22"Wx24"H Part: #6694
14"Wx13"H Part: #6447	14"Wx25"H Part: #6529	16"Wx17"H Part: #6463	17"Wx18"H Part: #6640	18"Wx19"H Part: #6476	19"Wx20"H Part: #6653	20"Wx21"H Part: #6497	21"Wx24"H Part: #6693	22"Wx25"H Part: #6537
14"Wx14"H Part: #6623	15"Wx17"H Part: #6462	16"Wx18"H Part: #6639	17"Wx19"H Part: #6475	18"Wx20"H Part: #6652	19"Wx21"H Part: #6496	20"Wx22"H Part: #6673	21"Wx25"H Part: #6536	
14"Wx15"H Part: #6449	15"Wx18"H Part: #6638	16"Wx19"H Part: #6474	17"Wx20"H Part: #6651	18"Wx21"H Part: #6495	19"Wx22"H Part: #6672	20"Wx23"H Part: #6516	22"Wx17"H Part: #6469	
14"Wx16"H Part: #6625	15"Wx19"H Part: #6473	16"Wx20"H Part: #6650	17"Wx21"H Part: #6494	18"Wx22"H Part: #6671	19"Wx23"H Part: #6515	20"Wx24"H Part: #6692	22"Wx18"H Part: #6645	
14"Wx17"H Part: #6461	15"Wx20"H Part: #6649	16"Wx21"H Part: #6493	17"Wx22"H Part: #6670	18"Wx23"H Part: #6514	19"Wx24"H Part: #6691	20"Wx25"H Part: #6535	22"Wx19"H Part: #6480	
14"Wx18"H Part: #6637	15"Wx21"H Part: #6492	16"Wx22"H Part: #6669	17"Wx23"H Part: #6513	18"Wx24"H Part: #6690	19"Wx25"H Part: #6534	21"Wx19"H Part: #6479	22"Wx20"H Part: #6656	

**1b. Sport Back Selection**

Sport Back Cushion, 20" Back Width or Less ..... No Charge  
 Select size below. Height is overall. Part number is SET plus number following #.

12" Wide, 10"-14" Width Package: Part: #130206	12"Wx14"-15"H Part: #130206	12"Wx16"-17"H Part: #130246		
14" Wide, 10"-14" Width Package: Part: #130247	14"Wx14"-15"H Part: #130247	14"Wx16"-17"H Part: #130248		
16" Wide, 16"-20" Width Package: Part: #130249	16"Wx16"-17"H Part: #130249	16"Wx18"-19"H Part: #130250	16"Wx20"-21"H Part: #130251	16"Wx22"-23"H Part: #130252
17" Wide, 16"-20" Width Package: Part: #130255	17"Wx18"-19"H Part: #130255	17"Wx22"-23"H Part: #130257	17"Wx24"-25"H Part: #130258	
18" Wide, 16"-20" Width Package: Part: #130259	18"Wx16"-17"H Part: #130259	18"Wx18"-19"H Part: #130260	18"Wx20"-21"H Part: #130261	18"Wx22"-23"H Part: #130262
19" Wide, 16"-20" Width Package: Part: #130267	19"Wx22"-23"H Part: #130267	19"Wx24"-25"H Part: #130268		
20" Wide, 16"-20" Width Package: Part: #130269	20"Wx16"-17"H Part: #130269	20"Wx18"-19"H Part: #130271	20"Wx22"-23"H Part: #130275	20"Wx24"-25"H Part: #130277
20" Wide, 20"-24" Width Package: Part: #130270	20"Wx16"-17"H Part: #130270	20"Wx18"-19"H Part: #130272	20"Wx22"-23"H Part: #130276	20"Wx24"-25"H Part: #130278

Sport Back Cushion, 22" Back Width or Greater ..... No Charge  
 Select size below. Height is overall. Available when selected with Mesomorph back canes. Part number is SET plus number following #.

21" Wide, 20"-24" Width Package: Part: #130292	21"Wx22"-23"H Part: #130292	21"Wx24"-25"H Part: #130293		
22" Wide, 20"-24" Width Package: Part: #130294	22"Wx16"-17"H Part: #130294	22"Wx22"-23"H Part: #130297	22"Wx24"-25"H Part: #130298	
24" Wide, 20"-24" Width Package: Part: #130299	24"Wx16"-17"H Part: #130299	24"Wx20"-21"H Part: #130303	24"Wx22"-23"H Part: #130305	24"Wx24"-25"H Part: #130308
24" Wide, 24"-28" Width Package: Part: #130300	24"Wx16"-17"H Part: #130300	24"Wx20"-21"H Part: #130304	24"Wx22"-23"H Part: #130306	24"Wx24"-25"H Part: #130309

**1c. TRU-Comfort Back Selection**

TRU-Comfort Back, Back Cane Option, Up to 20" Width ..... \$835  
 Part: TC4V10W21WBC. HCPCS: E2620  
 Back size based on seat width, back cane type, and overall height option selections.

TRU-Comfort Back, Back Cane Option, 22" Width ..... \$835  
 Part: TC4V22W32WBC. HCPCS: E2621  
 Back size based on seat width, back cane type, and overall height option selections.

TRU-Balance Flex Seating System section continued on following page.

**11. TRU-BALANCE FLEX SEATING SYSTEM CONTINUED**

**1d. ADI Back, AL Series, 4-Point, Active Contour**

Select a back type and the specific overall back height in subsequent section. Previously selected back cane height must be equal to the greater of overall height range.

- |  |  |
|--|--|
| Active Contour, 14" Wide x 13" or 14" Overall Height ..... \$865<br>Part: ALBR14LA4. | Active Contour, 18" Wide x 21" or 22" Overall Height ..... \$925<br>Part: ALBR18TA4. |
| Active Contour, 14" Wide x 16" or 17" Overall Height ..... \$885<br>Part: ALBR14MA4. | Active Contour, 18" Wide x 23" or 24" Overall Height ..... \$945<br>Part: ALBR18XA4. |
| Active Contour, 16" Wide x 19" or 20" Overall Height ..... \$905<br>Part: ALBR16SA4. | Active Contour, 20" Wide x 19" or 20" Overall Height ..... \$905<br>Part: ALBR20SA4. |
| Active Contour, 16" Wide x 21" or 22" Overall Height ..... \$925<br>Part: ALBR16TA4. | Active Contour, 20" Wide x 21" or 22" Overall Height ..... \$925<br>Part: ALBR20TA4. |
| Active Contour, 16" Wide x 23" or 24" Overall Height ..... \$945<br>Part: ALBR16XA4. | Active Contour, 20" Wide x 23" or 24" Overall Height ..... \$945<br>Part: ALBR20XA4. |
| Active Contour, 18" Wide x 19" or 20" Overall Height ..... \$905<br>Part: ALBR18SA4. |  |

**1e. ADI Back, AL Series, 4-Point, Deep Contour**

Select a back type and the specific overall back height in subsequent section. Previously selected back cane height must be equal to the greater of overall height range.

- |  |  |
|--|--|
| Deep Contour, 14" Wide x 13" or 14" Overall Height ..... \$965<br>Part: ALBR14LD4.   | Deep Contour, 18" Wide x 19" or 20" Overall Height ..... \$1,005<br>Part: ALBR18SD4. |
| Deep Contour, 14" Wide x 16" or 17" Overall Height ..... \$985<br>Part: ALBR14MD4.   | Deep Contour, 18" Wide x 21" or 22" Overall Height ..... \$1,025<br>Part: ALBR18TD4. |
| Deep Contour, 16" Wide x 19" or 20" Overall Height ..... \$1,005<br>Part: ALBR16SD4. | Deep Contour, 20" Wide x 19" or 20" Overall Height ..... \$1,005<br>Part: ALBR20SD4. |
| Deep Contour, 16" Wide x 21" or 22" Overall Height ..... \$1,025<br>Part: ALBR16TD4. | Deep Contour, 20" Wide x 21" or 22" Overall Height ..... \$1,025<br>Part: ALBR20TD4. |
| Deep Contour, 16" Wide x 23" or 24" Overall Height ..... \$1,045<br>Part: ALBR16XD4. | Deep Contour, 20" Wide x 23" or 24" Overall Height ..... \$1,045<br>Part: ALBR20XD4. |

**1f. ADI Back, AL Series, 4-Point, Extra Deep Contour**

Select a back type and the specific overall back height in subsequent section. Previously selected back cane height must be equal to the greater of overall height range.

- |   |   |
|---|---|
| Extra Deep Contour, 16" Wide x 21" or 22" Overall Height .... \$1,125<br>Part: ALBR16T64. | Extra Deep Contour, 18" Wide x 23" or 24" Overall Height .... \$1,145<br>Part: ALBR18X64. |
| Extra Deep Contour, 16" Wide x 23" or 24" Overall Height .... \$1,145<br>Part: ALBR16X64. | Extra Deep Contour, 20" Wide x 19" or 20" Overall Height .... \$1,105<br>Part: ALBR20S64. |
| Extra Deep Contour, 18" Wide x 19" or 20" Overall Height .... \$1,105<br>Part: ALBR18S64. | Extra Deep Contour, 20" Wide x 21" or 22" Overall Height .... \$1,125<br>Part: ALBR20T64. |
| Extra Deep Contour, 18" Wide x 21" or 22" Overall Height .... \$1,125<br>Part: ALBR18T64. | Extra Deep Contour, 20" Wide x 23" or 24" Overall Height .... \$1,145<br>Part: ALBR20X64. |

**1h. Omit Back Selection**

- Omit Back ..... No Charge  
Part: PTOINDV1287

TRU-Balance Flex Seating System section continued on following page.

**11. TRU-BALANCE FLEX SEATING SYSTEM CONTINUED**

**11.10. Overall Back Height**

<b>16" Overall Back Height</b> Part: SETHGHT1029	<b>19" Overall Back Height</b> Part: SETHGHT1032	<b>22" Overall Back Height</b> Part: SETHGHT1035	<b>25" Overall Back Height</b> Part: SETHGHT1038
<b>17" Overall Back Height</b> Part: SETHGHT1030	<b>20" Overall Back Height</b> Part: SETHGHT1033	<b>23" Overall Back Height</b> Part: SETHGHT1036	
<b>18" Overall Back Height</b> Part: SETHGHT1031	<b>21" Overall Back Height</b> Part: SETHGHT1034	<b>24" Overall Back Height</b> Part: SETHGHT1037	

**11.11. TRU-Balance Flex Armrests**

Left (L) and right (R) boxes may be checked to configure different size arms within one section arm type.

**1a. Drop-In, Height Adjustable, Removable, T-Arm Armrests**

Includes flat armpads.

<input type="checkbox"/> <input type="checkbox"/> Full Length, Drop-In Armrests..... No Charge L R Left part: ARM1803253. Right part: ARM156742	<input type="checkbox"/> <input type="checkbox"/> Desk Length, Drop-In Armrests ..... No Charge L R Left part: ARM1803254. Right part: ARM156743
--	---

**1b. Flip Back, Cane Mounted Armrests**

May interfere with swing-away laterals.

<input type="checkbox"/> <input type="checkbox"/> Full Length, Flip Back Armrests..... No Charge L R Left part: FRMASMB7038. Right part: FRMASMB7039	<input type="checkbox"/> <input type="checkbox"/> Desk Length, Flip Back Armrests ..... No Charge L R Left part: FRMASMB7040. Right part: FRMASMB7041
---	--

**1c. Two-Post, Flip Back, Height Adjustable, Removable Armrests**

Not available with a depth less than 16". Required with flip forward back and laptrays.

<input type="checkbox"/> <input type="checkbox"/> Full Length, Two-Post Armrests.....\$145 ea. L R Left part: XRF000051. Right part: XRF000052. HCPCS: E0973	<input type="checkbox"/> <input type="checkbox"/> Desk Length, Two-Post Armrests .....\$145 ea. L R Left part: XRF000053. Right part: XRF000054. HCPCS: E0973
---	--

**1d. Single Post, Height Adjustable, Removable Armrests**

Height adjustment from 11"-17" in half-inch increments.

<input type="checkbox"/> <input type="checkbox"/> Full Length, Single Post Armrests .....\$145 ea. L R Left part: ARM2000906. Right part: ARM2000907. HCPCS: E0973	<input type="checkbox"/> <input type="checkbox"/> Desk Length, Single Post Armrests .....\$145 ea. L R Left part: ARM2000908. Right part: ARM2000909. HCPCS: E0973
---	---

**11.12. TRU-Balance Flex Armpads**

Waterfall armpads may cause interference with the seat back. Armpads are configured as pairs.

**1a. Straight Armpads**

Straight, Full Length Left and Right ..... No Charge Part: ARM134275	Straight, Desk Length Left, Full Length Right ..... No Charge Part: ARM134280
Straight, Full Length Left, Desk Length Right ..... No Charge Part: ARM134281	Straight, Desk Length Left and Right..... No Charge Part: ARM134254

**1b. Waterfall Armpads**

Waterfall, Full Length Left and Right..... No Charge Part: ARM134285	Waterfall, Desk Length Left, Full Length Right ..... No Charge Part: ARM134288
Waterfall, Full Length Left, Desk Length Right ..... No Charge Part: ARM134289	Waterfall, Desk Length Left and Right ..... No Charge Part: ARM134282

**1c. Falcon Gel Armpads**

Falcon Gel, Full Length Left and Right ..... \$300 Part: ARM136438	Falcon Gel, Desk Length Left, Full Length Right..... \$250 Part: ARM136439
Falcon Gel, Full Length Left, Desk Length Right..... \$250 Part: ARM136442	Falcon Gel, Desk Length Left and Right ..... \$200 Part: ARM136437

**12. HEADRESTS**

**Optional.** Select a headrest pad and bracket.

**12.1a. Headrests, Comfort Plus with Cool Core**

**1. Headrest Pads**

Stealth 8" Comfort Plus, Cool Core (CP08CCFPM) ..... \$279 Part: ST-CP08CCFPM-Q. HCPCS: E0955	Stealth 14" Comfort Plus, Cool Core (CP14CCFPM) ..... \$279 Part: ST-CP14CCFPM-Q. HCPCS: E0955
Stealth 10" Comfort Plus, Cool Core (CP10CCFPM)..... \$279 Part: ST-CP10CCFPM-Q. HCPCS: E0955	Stealth 6", Comfort Plus, Cool Core (CP06CCFPM) ..... \$279 Part: ST-CP06CCFPM-Q. HCPCS: E0955

**2. Headrest Brackets**

Comfort Plus Headrest Mounting Bracket, UniLink..... \$485 Part: ST-UHB2D33200-Q. HCPCS: E1028	Comfort Plus Headrest Mounting Bracket, HMO475.....\$271 Part: ST-HMO475-17M-Q.
---	--

Headrests section continued on following pages.



**13. LAP BELTS**

**REQUIRED.** Select a lap belt below.

**13.1a. Standard Lap Belts**

<b>38" Lap Belt</b> Part: ACCBELT1006	<b>50" Lap Belt</b> Part: ACCBELT1000	<b>60" Lap Belt</b> Part: ACCBELT1002	<b>70" Lap Belt</b> Part: ACCBELT1003	<b>80" Lap Belt</b> Part: ACCBELT1004	<b>90" Lap Belt</b> Part: ACCBELT1001
--	--	--	--	--	--

**14. LEGRESTS**

**Optional.** All efforts will be made to accommodate legrest length selection, but minor changes may need to be made to meet shipping guidelines. Some adjustments may be needed upon receipt of the unit. Reference the section's compatibility matrix for legrest interference.

**14.1a. Basic Center Mount Foot Platform**

Select legrest length in following section.

Short	Medium	Long
8"-10"	10"-14"	14"-20"
Basic Center Mount Extension Range		

**1. Center Mount Foot Platform Type**

Basic Center Mount..... \$200  
Part: XRF000323

**2. Basic Center Mount Foot Platform Extension and Angle Adjustable Footplates**

See legrest length matrix for compatibility.

<b>Short Extension, 5"x9" Footplates</b> Part: MEC1906675. HCPCS: K0040	<b>Medium Lower Extension, 6.25"W Tapered 5"W x 12"D</b> Part: MEC2104853. HCPCS: K0040
<b>Medium Extension, 5"x9" Footplates</b> Part: MEC1906680. HCPCS: K0040	<b>Long Lower Extension, 6.25"W Tapered 5"W x 12"D</b> Part: MEC2104854. HCPCS: K0040
<b>Long Extension, 5"x9" Footplates</b> Part: MEC1906677. HCPCS: K0040	

**3. Basic Center Mount Calf Pads**

Calf pads are standard with a power articulating foot platform.

<b>Small Calf Pads, Basic Center Mount..... \$210</b> Part: RIG2104369	<b>Large Calf Pads, Basic Center Mount ..... \$210</b> Part: RIG2202997 Only available with Long Length Extension.
---	--

**4. Basic Center Mount Additional Accessories**

Basic Center Mount Footplate Connector ..... \$55  
Part: RIG2104240

*Legrests section continued on following pages.*

**14. LEGRESTS CONTINUED**

**14.1b. Center Mount Foot Platform**

If power articulating foot platform is selected, be sure power legrest control has been selected in previous section. Select legrest length in following section.

**1. Center Mount Foot Platform Type**

Quick Adjustable Center Mount Foot Platform..... \$265  
Part: XRF000311. HCPCS: K0040

Pediatric	Medium	Long
8.5"-12.5"	11"-16.5"	14"-20"
Quick Adj Center Mount Extension Range		

**2. Extension and Angle Adjustable Footplates**

See legrest length matrix for compatibility.

Pediatric Extension, 5"x6" Footplates  
Part: RIG152577

Not available with power articulating platform.

Pediatric Extension, 5"x8" Footplates  
Part: RIG152578

Not available with power articulating platform.

Medium Extension,  
5"x6" Footplates  
Part: RIG150206

Medium Extension,  
5"x8" Footplates with Width Extension  
Part: MEC1707469

Medium Extension,  
5"x12" Footplates with Width Extension  
Part: MEC1707465

Long Extension,  
5"x6" Footplates  
Part: RIG151505

Long Extension,  
5"x8" Footplates with Width Extension  
Part: MEC1707474

Long Extension,  
5"x12" Footplates with Width Extension  
Part: MEC1707468

**3a. Center Mount Calf Support**

\*Calf pads are standard with a power articulating foot platform. Select calf pads covers for additional charge. Calf pad and panel covers include 1/2" of soft visco foam.

Small Calf Pads, Quick Adjustable Center Mount Only, ..... \$210  
6"Wx6"H (pair)  
Part: XRF000156

Large Calf Pads, Quick Adjustable Center Mount Only, ..... \$210  
8" taper to 6"Wx8"H (pair)  
Part: XRF000157

**3b. Center Mount Calf Panel**

Calf pads are standard with a power articulating foot platform.

Small Calf Panel, 18"Wx8"H..... \$210  
Part: FRMASMB7110007

Large Calf Panel, 24"Wx8"H ..... \$265  
Part: RIG123824  
Only available on widths of 20" or greater.

**4. Modular Footplate Extension Components**

Optional. Assembly required. May require seating configuration adjustments to avoid front caster interference.

4" Outside Edge (pair) ..... \$55  
Part: RIG152510  
Only available with 6"D footplates and 8"D footplate extension.

2" Bell-Shaped Length Extension ..... \$55  
Part: RIG152519  
Only available with the 5x6 footplates.

**5. Additional Accessories**

Small Reversed Dartex Calf Pad Cover (pair) ..... \$150  
Part: RIG133717

Large Reversed Dartex Calf Pad Cover (pair) ..... \$150  
Part: RIG133719

Small Reversed Dartex Panel Cover (pair) ..... \$150  
Part: RIG134673

Large Reversed Dartex Panel Cover (pair) ..... \$150  
Part: RIG134744

Small Foam Calf Pad Cover (pair) ..... \$150  
Part: ST-CSCS0775-Q

Large Foam Calf Pad Cover (pair) ..... \$150  
Part: ST-CSCS0985-Q

Small Foam Panel Cover (pair) ..... \$150  
Part: ST-CSCP0920-Q

Large Foam Panel Cover (pair) ..... \$150  
Part: ST-CSCP0927-Q

Footplate Connector ..... \$55  
Part: RIGASMB7120023

Footplate Riser ..... \$150  
Part: XRF000348  
Not available with Long Extension. Select legrest length between 2" to 9" at 1" increments.

**6. Stealth Gel Footplate Covers**

5"x8" Tapered Gel Footplate Covers ..... \$420  
Part: ST-FPG5008T-Q. HCPCS: K0108

5"x11" Tapered Gel Footplate Covers ..... \$420  
Part: ST-FPG5011T-Q. HCPCS: K0108

6.5"x8" Tapered Gel Footplate Covers for Width Extension..... \$420  
Part: ST-FPG6508T-Q. HCPCS: K0108

5"x11" Tapered Gel Footplate Cover for Width Extension ..... \$420  
Part: ST-FPG7011T-Q. HCPCS: K0108

Legrests section continued on following pages.

**14. LEGRESTS CONTINUED**

**14.1b. 60°-90° Swing-Away Legrests**

Adjustable in 5° increments. Select a legrest upper, a legrest extension, and type of footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes. Select legrest length in following section.

**1a. 60°-90° Swing-Away Legrest Upper**

**60°-90° Swing-Away Legrest** ..... \$285  
 Pair part: RIG120062. Left: FRMASMB17558.  
 Right: FRMASMB17559  
 Not available with 11"-13" extensions.

**60°-90° Swing-Away Legrest, Short** ..... \$285  
 Pair part: RIG136481. Left: RIG132607.  
 Right: RIG132609.  
 For use only with 11"-13" lower extensions.

**2. Swing-Away Legrests Lower Extension**

**11"-13" Length**  
 Pair part: XRF000088.  
 Left: XRF000089.  
 Right: XRF000090.  
 N/A with Multi-axis.

**13"-16.5" Length**  
 Pair part: XRF000094.  
 Left: XRF000095.  
 Right: XRF000096.  
 N/A with Multi-Axis.

**11"-16.5" Length**  
 Pair part: XRF000091.  
 Left: XRF000092.  
 Right: XRF000093.  
 N/A with Angle Adjustable.

**15.5"-19.5" Length**  
 Pair part: XRF000097.  
 Left: XRF000098.  
 Right: XRF000099.

Legrest Compatibility Matrix 60°-90° Swing-Away Legrest TRU-Balance Flex	Seat to Floor Option		
	Low	Middle	High
Extruded/Angle A., 11"-13"	Full	Full	Full
Extruded/Angle A., 13"-16.5"	13"	13"-14"	13"-15"
Extruded/Angle A., 15.5"-19.5"	N/A	N/A	N/A
Multi-Axis, 11"-16.5"	Full	Full	Full
Multi-Axis, 15.5"-19.5"	N/A	N/A	N/A

**3a. Extruded Footplate(s)**

**7.25"Wx6"D Extruded Footplate** ..... Standard  
 Part pair: FRMASMB7847. Left: FRMASMB7861  
 Right: FRMASMB7862

**8.5"Wx6"D Extruded Footplate** ..... Standard  
 Part pair: FRMASMB7848. Left: FRMASMB7863  
 Right: FRMASMB7864  
 May cause interference on seat widths less than 18."

**3b. Angle Adjustable Footplate(s) and Heel Loops**

<sup>1</sup>The 6"x11" footplates may cause interference with seat widths less than 18". <sup>2</sup>The 8"x11" footplates may cause interference with seat widths less than 20".

**4"Wx6"D** ..... \$87.50 ea.  
 Pair part: FRMASMB10564. Left: RIG133537. Right: RIG133560

**6"Wx8"D<sup>1</sup>** ..... \$87.50 ea.  
 Pair part: FRMASMB10567. Left: RIG133540. Right: RIG133563

**5"Wx6"D** ..... \$87.50 ea.  
 Pair part: FRMASMB10565. Left: RIG133538. Right: RIG133561

**6"Wx11"D<sup>1</sup>** ..... \$87.50 ea.  
 Pair part: FRMASMB10568. Left: RIG133541. Right: RIG133564

**5"Wx8"D** ..... \$87.50 ea.  
 Pair part: FRMASMB10566. Left: RIG133539. Right: RIG133562

**8"Wx11"D<sup>2</sup>** ..... \$87.50 ea.  
 Pair part: FRMASMB10569. Left: RIG133542. Right: RIG133565

**3c. Multi-Axis Footplate(s) and Heel Loops**

<sup>1</sup>The 6"x11" footplates may cause interference with seat widths less than 18". <sup>2</sup>The 8"x11" footplates may cause interference with seat widths less than 20".

**4"Wx6"D** ..... \$260 ea.  
 Pair part: FRMASMB9456. Left: RIG133531. Right: RIG133552

**6"Wx8"D<sup>1</sup>** ..... \$260 ea.  
 Pair part: FRMASMB9459. Left: RIG133534. Right: RIG133555

**5"Wx6"D** ..... \$260 ea.  
 Pair part: FRMASMB9457. Left: RIG133532. Right: RIG133553

**6"Wx11"D<sup>1</sup>** ..... \$260 ea.  
 Pair part: FRMASMB9460. Left: RIG133535. Right: RIG133558

**5"Wx8"D** ..... \$260 ea.  
 Pair part: FRMASMB9458. Left: RIG133533. Right: RIG133554

**8"Wx11"D<sup>2</sup>** ..... \$260 ea.  
 Pair part: FRMASMB9461. Left: RIG133536. Right: RIG133559

**4. 60°-90° Swing-Away Legrest Accessories**

Gel Neoprene Footrest Sleeve (pair) ..... \$130  
 Part: INDPART2243. HCPCS: K0108

Legrests section continued on following pages.

**14. LEGRESTS CONTINUED**

**14.1c. Elevating Legrests**

Select a legrest upper, a legrest extension, and type of footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes. Select legrest length in following section.

**1. Elevating Legrests Upper**

P  L  R Elevating Legrest Assembly ..... \$275 ea.  
 Pair part: FRMASMB6483. Left: RIG138012  
 Right: RIG138013. HCPCS: E0990  
 Select lower extension and footplates

Legrest Compatibility Matrix Elevating Legrest TRU-Balance Flex	Seat to Floor Option		
	Low	Middle	High
Reinforced, 11"-14"	11"-13"	Full	Full
Reinforced, 14"-17"	N/A	14"	14"-15"
Angle Adjustable, 12"-16.5"	12"-13"	12.5"-14"	12.5"-15"
Angle Adjustable, 15.5"-19.5"	N/A	N/A	N/A

**2. Lower Extension**

P  L  R 11"-14" Length  
 Pair: FRMASMB4755. Left: RIG138026. Right: RIG138030  
 For reinforced only

P  L  R 14"-17" Length  
 Pair: FRMASMB4756. Left: RIG138033. Right: RIG138035  
 For reinforced only

P  L  R 8.5"-12" HD Drop-Ins/12.5"-16.5" ELR, Ext Length  
 Pair part: XRF000101. Left: XRF000102. Right: XRF000103  
 For angle adjustable only

P  L  R 13"-18" HD Drop-Ins/15.5"-19.5" ELR, Ext Length  
 Pair part: XRF000104. Left: XRF000105. Right: XRF000106  
 For angle adjustable only

**3a. Reinforced Cast Footplate(s)**

P  L  R 7.25"Wx6"D Reinforced Footplate..... Standard  
 Pair part: FRMASMB4748. Left: RIG138055. Right: RIG138056

P  L  R 8.5"Wx6"D Reinforced Footplate ..... Standard  
 Pair part: FRMASMB4747. Left: RIG138060. Right: RIG138062  
 May cause interference on seat widths less than 18."

**3b. Angle Adjustable Footplate(s)**

\*The 6"x11" footplates may cause interference with seat widths less than 18". \*The 8"x11" footplates may cause interference with seat widths less than 20". 4"x6" and 5"x6" footplates are only available in pairs.

P  L  R 4"Wx6"D ..... \$87.50 ea.  
 Pair parts: FRMASMB8812 x2.

P  L  R 6"Wx8"D<sup>1</sup> Footplates with Heel Loops ..... \$87.50 ea.  
 Pair part: FRMASMB4738. Left: RIG137793. Right: RIG137795

P  L  R 5"Wx6"D..... \$87.50 ea.  
 Pair parts: FRMASMB8811 x2.

P  L  R 6"Wx11"D<sup>1</sup> Footplates with Heel Loops ..... \$87.50 ea.  
 Pair part: FRMASMB4741. Left: RIG137796. Right: RIG137797

P  L  R 5"Wx8"D Footplates with Heel Loops ..... \$87.50 ea.  
 Pair part: FRMASMB4735. Left: RIG137791. Right: RIG137792

P  L  R 8"Wx11"D<sup>2</sup> Footplates with Heel Loops ..... \$87.50 ea.  
 Pair part: FRMASMB4744. Left: RIG137798. Right: RIG137799

**14.2. Legrest Accessories**

Heel Loops (pair)..... \$35  
 Part: FRMASMB7873. HCPCS: E0951  
 Standard on angle adjustable and multi axis footplates when noted.

Wheel Bumper..... \$15 ea.  
 Part: XRF000100

**14.3. Set Legrest Extension Length**


REQUIRED with matching left and right legrest configuration. Part number is PTOINDV plus number following # below. Non-matching or individual legrest configurations lengths are not set by the manufacturer.

7" L (#3470)	8.5" L (#3473)	10" L (#3476)	11.5" L (#3479)	13" L (#3482)	14.5" L (#3485)	16" L (#3488)	17.5" L (#3491)	19" L (#3494)	20.5" L (#3497)	22" L (#3669)
7.5" L (#3471)	9" L (#3474)	10.5" L (#3477)	12" L (#3480)	13.5" L (#3483)	15" L (#3486)	16.5" L (#3489)	18" L (#3492)	19.5" L (#3495)	21" L (#3498)	22.5" L (#3670)
8" L (#3472)	9.5" L (#3475)	11" L (#3478)	12.5" L (#3481)	14" L (#3484)	15.5" L (#3487)	17" L (#3490)	18.5" L (#3493)	20" L (#3496)	21.5" L (#3600)	

**15. SEAT MOUNTED LED LIGHTS**

Optional. Requires Q-Logic 3 Joystick for Full Light Package or Q-Logic 3e Joystick.

Seat Mounted, LED, Full Light Package ..... \$770  
 Part: XRF000130  
 Not available with Base Mounted, LED, Full Light Package.



**16. ACCESSORIES**

**Optional.** Multiple accessories can be mounted along the various tracks but some interference may occur. Accessories may need to be installed by the provider.

**Flex Seat Cup Holder** ..... No Charge  
 Part: ACCASMB2634  
 Not available on TRU-Balance Flex seating without two-post armrests.

**Oxygen Holder for Flex Seat** ..... \$110  
 Part: ACC125006. HCPCS: E2208  
 Not available with Manual Recline.

**Cane/Crutch Holder** ..... \$198  
 Flex part: ACC152515. HCPCS: E2207

**Flag Holder** ..... \$70  
 Flex part: ACCASMB2415.

**Saddle Bag** ..... \$15  
 Part: DWR1020D010

**Quantum Backpack** ..... \$200  
 Part: ACC1905853

**Accessory Bar** ..... \$95  
 Part: XRF000128

**Large, Black, Weather Cover** ..... \$240  
 Part: ACCCOVR1021

**VOCSN Vent Tray, Accessory Bar Mounted** ..... \$1,000  
 Part: ACC2000166  
 Accessory bar is required. Credit for bar may be applied.

**Privacy Bag** ..... \$80  
 Part: ACC1907570

**Transfer Handle** ..... \$305  
 Part: MEC1809930. HCPCS: E0705