



Pride Mobility Products Co.

5096 South Service Rd, Beamsville, ON, L3J 1V4

Phone: 888-570-1113 | Fax: 866-514-1303 | Email: productorders@pridemobility.com

# Quantum 4Front<sup>®</sup> 2 Order Form with TRU-Balance 3 Positioning

Account Number: \_\_\_\_\_ Date: \_\_\_\_\_

Provider Name: \_\_\_\_\_

Contact: \_\_\_\_\_

Phone: \_\_\_\_\_ Fax: \_\_\_\_\_

Email: \_\_\_\_\_

PO Number: \_\_\_\_\_

Marked for: \_\_\_\_\_

Ship to Address: \_\_\_\_\_

City: \_\_\_\_\_ Prov: \_\_\_\_\_ PC: \_\_\_\_\_

## QUANTUM 4FRONT<sup>2</sup>



### INTRODUCTION

This form is interactive when viewed with [Adobe Acrobat Reader](#). Complete the form by placing checks in the desired boxes and provide information in the interactive fields. Buttons shown at the bottom of the form may be utilized to print or submit the order form through a desktop email application. To email via a web-based application, please 'Save As' and attach the PDF to your email.

This order form contains a large variety of options to fulfill various patient needs. Descriptions and section notes, such as optional or required, are included to help you complete your order. Other part numbers can be found on the quote acknowledgment after this order form is submitted. Please contact Quantum<sup>®</sup> Sales at 888-570-1113 if assistance is needed. Send the completed order form by fax (866-514-1303) or email ([productorders@pridemobility.com](mailto:productorders@pridemobility.com)). Incomplete forms may delay the quote or order. Customer service will contact you if the order is incomplete or if there are compatibility issues. If special order requests are needed, be sure that the Patient Information section is completed

All prices are MSRP. Prices, specifications, part numbers, and availability are subject to change without notice. Prices and part numbers as shown are only available when configured as a complete power chair. Please contact technical service for accurate parts ordering.


**Thank you for choosing Quantum!**

### PATIENT INFORMATION

Weight: \_\_\_\_\_ Height: \_\_\_\_\_ A. Shoulder Width: \_\_\_\_\_ D. Max Sitting Height: \_\_\_\_\_ G. Thigh Depth: \_\_\_\_\_

B. Chest Width: \_\_\_\_\_ E. Shoulder Height: \_\_\_\_\_ H. Lower Leg Length: \_\_\_\_\_

C. Hip Width: \_\_\_\_\_ F. Axilla Height: \_\_\_\_\_ I. Elbow Height: \_\_\_\_\_



\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

### 1. BASE MODEL

**REQUIRED.** Select a group 3 power base. 300 lb. weight capacity with standard 6 mph programmable motors and 8 amp off board charger. Dimensions of base without legrests: 24.2"Wx39.9"L

Quantum 4Front 2 ..... \$8,450  
Part: 4FRONT 2















**REQUIRED.** Smart Traction Control is auto-selected with every 4Front 2.

Smart Traction Control ..... \$400

**2. BASE COLOR**

**REQUIRED.** Select one shroud color. Colors shown here may differ from actual product depending on monitor/printer calibration used.

**2.1a. Standard Shrouds**

- |   |  |  |  |
|---|--|--|--|
|  <b>Red Corvette</b><br>Part: DGN2003028-PRD1  |  <b>Deep Purple</b><br>Part: DGN2003028-PPR5      |  <b>Midnight Blue</b><br>Part: DGN2003028-PBU10 |  <b>Green Machine</b><br>Part: DGN2003028-PGN11 |
|  <b>Back in Black</b><br>Part: DGN2003028-PBK1 |  <b>Raspberry Beret</b><br>Part: DGN2003028-PPK4  |  <b>Electric Blue</b><br>Part: DGN2003028-PBU9  |  <b>Military Green</b><br>Part: DGN2003028-PGN8 |
|  <b>Night Train</b><br>Part: DGN2003028-PBK3   |  <b>Pretty in Pink</b><br>Part: DGN2003028-PPK3   |  <b>Blue Jean</b><br>Part: DGN2003028-PBU11     |  |
|  <b>Grey Street</b><br>Part: DGN2003028-PSL3   |  <b>Evolution Orange</b><br>Part: DGN2003028-POR7 |  <b>Lemon Crush</b><br>Part: DGN2003028-PYL4    |  |

**2.1b. Hydrographic Shrouds**

- |   |   |
|---|---|
|  <b>Military Digital Camo</b> ..... \$325<br>Part: DGN2003028-WCM7 |  <b>Faded Glory</b> ..... \$325<br>Part: DGN2003028-WFL2 |
|---|---|

**3. BASE MOUNTED LIGHTS**

**Optional.** Select one option. Seat Mounted, LED, Full Light Package option can be found later on this form.

- |  |   |
|--|---|
| <b>Switch Operated Fender Lights</b> ..... No Charge<br>Part: DGN2004686 | <b>Omit Fender Lights</b> ..... No Charge<br>Part: DGN1700032 |
|--|---|

Fender lights operate by pressing included single switch button. Fender lighting is not available with Single Actuator Control through Toggle Box.

**4. QUANTUM BACKUP CAMERA**

**Optional.** The backup camera may interfere with armrest mounted accessories such as media device holder.

- Quantum Backup Camera** ..... \$475  
Part: MEC1910365

**1. Optional Camera Mounting Bracket**

- |  |   |
|--|---|
| <b>Backup Camera Headrest Bracket</b> ..... \$30<br>Part: MEC1910373 | <b>Backup Camera Track Mount Bracket</b> ..... \$35<br>Part: MEC2104227 |
|--|---|

Works with .5" x .5" vertical supports only.

**2. Optional Display Mounting Bracket**

- Backup Camera Q-Logic Joystick Mount Bracket** .... \$145  
Part: MEC2106123



**5. DRIVE WHEELS**

**REQUIRED.** Select a type of drive wheels and accent color.

**5.1. Drive Wheels**

- |  |   |
|--|---|
| <b>3"x14" Flat Free Black Tires with Accent</b> ..... Standard<br>Part: PTOINDV3525    | <b>4"x14" Flat Free Knobby Black Tires</b> ..... \$120<br>Part: WHL2006647 x2 |
| <b>3"x14" Pneumatic Black Tires with Accent</b> ..... No Charge<br>Part: WHL1704059 x2 | <b>4"x14" Pneumatic Knobby Black Tires</b> ..... \$120<br>Part: WHL2006822 x2 |

**5.2. Drive Wheel Accent Color**

Not available with knobby tires.

- |   |  |  |
|---|--|--|
| <b>Back in Black</b> ..... Standard<br>Part: DGN164245-PBK1 x2    | <b>Pretty in Pink</b> ..... No Charge<br>Part: DGN164245-PPK3 x2   | <b>Green Machine</b> ..... No Charge<br>Part: DGN164245-PGN11 x2       |
| <b>Red Corvette</b> ..... No Charge<br>Part: DGN164245-PRD1 x2    | <b>Evolution Orange</b> ..... No Charge<br>Part: DGN164245-POR7 x2 | <b>Military Green</b> ..... No Charge<br>Part: DGN164245-PGN8 x2       |
| <b>Night Train</b> ..... No Charge<br>Part: DGN164245-PBK3 x2     | <b>Midnight Blue</b> ..... No Charge<br>Part: DGN164245-PBU10 x2   | <b>Military Digital Camo</b> ..... \$50 ea.<br>Part: DGN164245-WCM7 x2 |
| <b>Grey Street</b> ..... No Charge<br>Part: DGN164245-PSL3 x2     | <b>Electric Blue</b> ..... No Charge<br>Part: DGN164245-PBU9 x2    | <b>Faded Glory</b> ..... \$50 ea.<br>Part: DGN164245-WFL2 x2           |
| <b>Deep Purple</b> ..... No Charge<br>Part: DGN164245-PPR5 x2     | <b>Blue Jean</b> ..... No Charge<br>Part: DGN164245-PBU11 x2       |  |
| <b>Raspberry Beret</b> ..... No Charge<br>Part: DGN164245-PPK4 x2 | <b>Lemon Crush</b> ..... No Charge<br>Part: DGN164245-PYL4 x2      |  |

**6. BATTERIES**

**Optional.** The powerchair requires 2 flat top batteries and will be installed if ordered.

Group 24, AGM Cell Battery (MK).....\$510 ea.  
Part: BAT1712207. HCPCS: E2363

Group 34, Gel Cell Battery (MK) ..... \$400 ea.  
Part: BATGEL1006. HCPCS: E2359  
Not available with TRU-Balance 3 Power Lift.

**7. PATIENT'S WEIGHT**

**REQUIRED.** The patient's weight is required to quote and order the power chair.

Up to 200 lb User Weight  
Part: PTOINDV3837

201 lb-300 lb User Weight  
Part: PTOINDV3838

**8. ELECTRONICS**

**REQUIRED.** Select a type of drive control from the following sections. The No Joystick option is used to omit the joystick and utilize another control, such as a specialty control at the end of this form or found on the [Alternative Controls and Electronics Order Form](#).

**8.1a. Omit Joystick**

No Joystick, Setup for Q-Logic Electronics  
Part: XRF000010  
Requires selection of Enhanced Display

**8.1c. Q-Logic 3 Expandable Capacity Controls**

Q-Logic 3 expandable capacity controls allow use of specialty controls and multiple actuators.

**1. Q-Logic 3 Expandable Capacity Joysticks**

Q-Logic 3 Expandable Capacity Joystick .....\$900  
Part: CTL167724

Q-Logic 3 Expandable Capacity Joystick .....\$900  
for Full Light Package  
Part: CTL167725  
Required with full light packages.

Q-Logic 3e Joystick .....\$620  
Part: CTL1910600

Can operate full light package. Not available with alternative or specialty controls.



**8.2. Q-Logic 3 Joystick Armor**

**REQUIRED with Q-Logic 3 Joystick selections above. Not available with Q-Logic 3e Joystick. Reference colors shown on page 2.**

Red Corvette  
Part: DGN124336-PRD1

Deep Purple  
Part: DGN124336-PPR5

Midnight Blue  
Part: DGN124336-PBU10

Green Machine  
Part: DGN124336-PGN11

Carbon Fiber  
Part: DGN124336-WCB1

Back in Black  
Part: DGN124336-PBK1

Raspberry Beret  
Part: DGN124336-PPK4

Electric Blue  
Part: DGN124336-PBU9

Military Green  
Part: DGN124336-PGN8

Night Train\*  
See Black, Soft Touch

Pretty in Pink  
Part: DGN124336-PPK3

Blue Jean  
Part: DGN124336-PBU11

Military Digital Camo  
Part: DGN124336-WCM7

Grey Street  
Part: DGN124336-PSL3

Evolution Orange  
Part: DGN124336-POR7

Lemon Crush  
Part: DGN124336-PYL4

Faded Glory  
Part: DGN124336-WFL2

**9. JOYSTICK HANDLES**

**Optional.** The standard Q-Logic 3 joystick knob can be seen in an image in the previous section. An optional handle selected below will be mounted to the joystick and the standard handle will be shipped with the unit.

Mushroom Handle (B) ..... \$120  
Part: ACCASMB1860.

Chin Cup (C)..... \$85  
Part: ACCASMB186190

T-Bar (D) ..... \$160  
Part: ACCASMB1858.

Medium Ball Handle (E)..... \$70  
Part: ACCASMB1859.

Soft Ball Handle (F)..... \$120  
Part: ACCASMB1856.

Extended I-Handle (G) ..... \$140  
Part: ACCASMB1857.

Bodypoint 3" Goal Post Handle ..... \$180  
Part: PLSKNCP1057.

Bodypoint 4" Goal Post Handle ..... \$180  
Part: PLSKNCP1058.



**10. JOYSTICK MOUNTING BRACKETS**

**Optional.** Select a mount to be used with a joystick selected in previous section.

**10.1a. Fixed, In-line Joystick Bracket**

Fixed, In-line, Bracket ..... Standard  
Part: FRMASMB10617



**10.1b. Swing-Away (Retractable) Joystick Bracket**

Allows the joystick to retract to the side of the armrest while maintaining the same joystick orientation.

Swing-away, Q-Logic 3 ..... \$265  
Part: ACC1908434. HCPCS: E1028

**10.1c. Basic Swing Away Joystick Mount**


Basic Swing Away, Q-Logic, Left Side ..... \$265  
Part: ACC2201041. HCPCS: E1028

Basic Swing Away, Q-Logic, Right Side ..... \$265  
Part: ACC2201042. HCPCS: E1028

**10.1d. Offset Swing Away Joystick Mount**

Offset Swing Away, Q-Logic, Left Side ..... \$265  
Part: ACC142742. HCPCS: E1028

Offset Swing Away, Q-Logic, Right Side ..... \$265  
Part: ACC142743. HCPCS: E1028



**10.1d. Specialty Joystick Mounts**

Select a bracket and side to mount it.

Retract4 Mount ..... \$450  
Part: XRF000003. HCPCS: E1028  
Full 180° of movement that translates into multi-position use, straight, right offset, left offset, or fold away. Joystick orientation varies as bracket is adjusted.

Stealth Swing-Away/Flip Down Bracket ..... \$300  
Part: XRF000005. HCPCS: E1028  
Allows the joystick to pivot to side of armrest, or flip-down if bracket is adjusted.

Multi-Axis Mount ..... \$100  
Part: XRF000004  
Allows height adjustment of joystick.



**10.2. Joystick Mount Location**

Left Side Mount  
Part: PTOINDV1225

Right Side Mount  
Part: PTOINDV1226

**11. TRU-BALANCE 3® SEATING - QUANTUM INTELLIGENT SYSTEMS™**

**11.1. TRU-Balance 3 Configuration**

**REQUIRED.** Select one seating system. Be sure to complete the seat dimensions, seat to floor, and back selections.

TRU-Balance 3 Static Seat ..... \$1,045  
Part: PTOINDV3442  
20° of limited recline adjustment (-5/ +15) via gas strut within back cover.

TRU-Balance 3 Power Recline Only ..... \$7,875  
Part: PTOINDV3444. HCPCS: E1004  
175° of power recline.

TRU-Balance 3 Power Tilt Only ..... \$6,245  
Part: PTOINDV3443. HCPCS: E1002  
50° Center of Gravity (CG) power tilt seating system.

TRU-Balance 3 Power Tilt & Recline ..... \$13,125  
Part: PTOINDV3446. HCPCS: E1007  
50° CG power tilt with 175° of power recline. Full drive lockout at 20°.

**11.2. 10" Power Lift**

**Optional.** Select one option below. ^Please verify the appropriate billing code to be used for this item with all third-party payors.

TRU-Balance 3 Power 10" Lift ..... \$2,995  
Part: PTOINDV3775. HCPCS: E2298/E2300^  
10" of scissor lift. 22"W x 22"D maximum seat size.

**11.3. Occupied Transit**

**Optional.** Available in seat dimensions of 16" to 20" wide and 16"-20" deep. Requires TRU-Balance 3 armrests and 60" Sure-lok lap belt. Only available with swing-away legrests and center mount foot platforms.

Occupied Transit ..... \$400  
Part: KIT152047

To determine if your chair is WC-19 compliant (occupied), please refer to this website: <https://wc-transportation-safety.umtri.umich.edu/>

**11.4. Seat Width**

**REQUIRED with seat selection.** 15" width is only available with non-matching back width. Minimum seat width is 12". Maximum seat width is 22".

12" W (#1002)	13" W (#1003)	14" W (#1004)	15" W (#1005)	16" W (#1006)	17" W (#1007)	18" W (#1008)	19" W (#1009)	20" W (#1010)	21" W (#1011)	22" W (#1012)
------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------

**11.5. Seat Width Package**

**REQUIRED with seat selection.**

12"-16" Width Package  
Part: XRF000006  
Not available with Power Lift.

14"-16" Width Package  
Part: XRF000007  
Only available with Power Lift.

16"-20" Width Package  
Part: XRF000008

20"-24" Width Package  
Part: XRF000158  
Maximum width is 22".

**11.6. Seat Depth**

**REQUIRED with seat selection.** Minimum seat depth is 12". Maximum seat depth is 22".

12" D (#1002)	13" D (#1003)	14" D (#1004)	15" D (#1005)	16" D (#1006)	17" D (#1007)	18" D (#1008)	19" D (#1009)	20" D (#1010)	21" D (#1011)	22" D (#1012)
------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------

TRU-Balance 3 seating section continued on following pages.

**11. TRU-BALANCE 3 SEATING CONTINUED**

**11.7. Seat Depth Package**

**REQUIRED with seat selection.**

12"-16" Depth Package  
Part: XRF000011

Not available with Power Lift.

14"-16" Depth Package  
Part: XRF000012

Only available with Power Lift.

16"-20" Depth Package  
Part: XRF000013

20"-24" Depth Package  
Part: XRF000159

Maximum depth is 22".

**11.8. Seat to Floor Height**

**REQUIRED.**

Lowest Seat to Floor  
Part: PTOINDV3558

Middle Seat to Floor  
Part: PTOINDV3559

Highest Seat to Floor  
Part: PTOINDV3560

Low	Med	High
Static Seat/ Recline Only 15"	16"	17"
Tilt/Tilt & Recline/Lift 17.75"	18.75"	19.75"
TRU-Balance 3 Seat to Floor Heights		


**11.9. Matching Seat Cushion**

**Optional.**

**Stealth TRU-Comfort 2 SPP Cushion**

Skin protection and positioning with Coolcore, 4-way stretch, breathable cover. Includes urethane liner.  
Part series: ST-TC2SPPXXYY. HCPCS: E2607/E2608


TRU-Comfort 2 SPP, Matching Seat Size  
Reference part: XRF000116.



**Stealth Solution SPP Cushion**

Skin protection and positioning with Coolcore, 4-way stretch, breathable cover. Includes urethane liner.  
Part series: ST-SOLSPPXXYY. HCPCS: E2607/E2608


Solution SPP, Matching Seat Size  
Reference part: XRF000113.



**Stealth Spectrum Foam SPP Cushion**

Skin protection and positioning with Coolcore, 4-way stretch, breathable cover. Includes urethane liner.  
Part series: ST-SPCFSPXXYY. HCPCS: E2607/E2608

Spectrum Foam SPP, Matching Seat Size  
Reference part: XRF000112.



**Stealth Premiere P Cushion**

Positioning with Coolcore, 4-way stretch, breathable cover.  
Part series: ST-PREMPXXYY. HCPCS: E2605/E2606


Premiere P, Matching Seat Size  
Reference part: XRF000122.



**Stealth Zen SP Cushion**

Skin protection with Coolcore, 4-way stretch, breathable cover.  
Part series: ST-ZENSPXXYY. HCPCS: E2603/E2604


Zen SP, Matching Seat Size  
Reference part: XRF000115.



**Stealth Spectrum Gel SPP Cushion**

Skin protection and positioning with Coolcore, 4-way stretch, breathable cover. Includes urethane liner.  
Part series: ST-SPCGSPPXXYY. HCPCS: E2607/E2608


Spectrum Gel SPP, Matching Seat Size  
Reference part: XRF000117. MSRP: \$507/\$585



**Stealth Essence SPP Cushion**

Skin protection and positioning with Coolcore, 4-way stretch, breathable cover.  
Part series: ST-ESSSPPXXYY. HCPCS: E2607/E2608

Essence SPP, No Wedge, Matching Seat Size  
Reference part: XRF000118.



Essence SPP, Small Wedge, Matching Seat Size  
Reference part: XRF000119.

Essence SPP, Medium Wedge, Matching Seat Size  
Reference part: XRF000120.


Essence SPP, Large Wedge, Matching Seat Size  
Reference part: XRF000121.

Essence Plus SPP, Large Wedge, Matching Seat Size  
Reference part: XRF000121.

**Stealth Glacial SP Cushion**

Skin protection with Coolcore, 4-way stretch, breathable cover.  
Part series: ST-GLSPXXYY. HCPCS: E2603/E2604


Glacial SP, Matching Seat Size  
Reference part: XRF000114.



**Stealth Simplicity G Cushion**

General use cushion with Stealth Tek 2, 2-way stretch, breathable cover.  
Part series: ST-SIMGXXYY. HCPCS: E2601/E2602

Simplicity G, Matching Seat Size  
Reference part: XRF000111.



**11.10. Back Type Selection**

**REQUIRED.** Match a choice in this section with the options in the following back cushion configuration section.

Redesigned, Static Type Back..... \$735  
Part: XRF000014

12" and 24" back width will use the TRU-Balance 3 Back assembly. Select a back cushion configuration in the subsequent section.

Back Cane Option..... Standard  
Part: XRF000015

Available only with TRU-Balance 3 Static Seat and Power Tilt (1 actuator) selections. Select cane type, back cane height, and back cushion configuration below.

Redesigned, Flip Forward Back..... \$735  
Part: XRF000016

12" and 24" back width will use the TRU-Balance 3 Back assembly. Not available with power recline systems. Only available with 2-post, flip back or pediatric armrests. Not available with 12"-14" width. Select a back cushion type in the subsequent section.

**11.11. Back Cane Selection**

**REQUIRED when Back Cane Option is selected above.** Select type and height of back canes.

**1. Back Cane Type**

Angle Adjustable Back Cane..... Standard  
Part: MEC140266

Flip Forward Back Cane ..... No Charge  
Part: XRF000273

2" Endomorph (1" per side) Back Cane..... No Charge  
Part: XRF000269

4" Endomorph (2" per side) Back Cane..... No Charge  
Part: XRF000271

TRU-Balance 3 seating section continued on following pages.

**11. TRU-BALANCE 3 SEATING CONTINUED**

**2. Back Cane Height**

16"H Cane Part: FRMASMB15834	20"H Cane Part: FRMASMB15838	24"H Cane Part: FRMASMB15842
17"H Cane Part: FRMASMB15835	21"H Cane Part: FRMASMB15839	25"H Cane Part: FRMASMB15843

**11.12. Back Cushion Configuration**

**REQUIRED.** Sport Back and TRU-Comfort Back selections up to 22" back width will be the Sport 4 or TRU-Comfort 4 back cushions. Sport back and TRU-Comfort back selections of 24" back width will be the Sport or TRU-Comfort 2 back cushions. Select a configuration type and applicable dimensions. See the Sport, Q-Back, and TRU-Comfort back availability chart for stock order back size ranges. Other size backs may be available with special order and extended lead time.

**1a. Sport Back/Omit/After-Market Interface**

Sport Back ..... No Charge Part: XRF000017 Complete following width package and size.	After-market Back Interface ..... No Charge Part: ST-SSAPT4-Q. Select a set width in following section, range of 12"-21". Overall seat depth may need to be adjusted to accommodate certain backs. May interfere with gel/waterfall armpads and arm through. Does not include brackets to mount after-market backs.
Omit Back Pan/Cushion ..... No Charge Part: PTOINDV1287 Select a following width package.	

**1b. TRU-Comfort Back**

TRU-Comfort Back (Static/Flip Forward), Up to 15" Width ..... \$835 Part: SET123622. HCPCS: E2620 Complete following width package and size.	TRU-Comfort Back, Back Cane Option ..... \$835 Part: TC4V10W21WBC. HCPCS: E2620 Complete following width package and size.
TRU-Comfort Back (Static/Flip Forward), 16" to 20" Width... \$1,035 Part: TC4V10W21W. HCPCS: E2620 Complete following width package and size.	
TRU-Comfort Back (Static/Flip Forward), 22" Width ..... \$1,235 Part: TC4V22W32W. HCPCS: E2621 Complete following width package and size.	

**1c. ADI Back, AL Series, 4-Point, Active Contour**

**Available with Back Cane Option only.** Select a back type and the specific overall back height in subsequent section. Previously selected back cane height must be equal to the greater of overall height range.

Active Contour, 14" Wide x 13" or 14" Overall Height ..... \$865 Part: ALBR14LA4. HCPCS: E2613	Active Contour, 18" Wide x 21" or 22" Overall Height ..... \$925 Part: ALBR18TA4. HCPCS: E2613
Active Contour, 14" Wide x 16" or 17" Overall Height ..... \$885 Part: ALBR14MA4. HCPCS: E2613	Active Contour, 18" Wide x 23" or 24" Overall Height ..... \$945 Part: ALBR18XA4. HCPCS: E2613
Active Contour, 16" Wide x 19" or 20" Overall Height ..... \$905 Part: ALBR16SA4. HCPCS: E2613	Active Contour, 20" Wide x 19" or 20" Overall Height ..... \$905 Part: ALBR20SA4. HCPCS: E2613
Active Contour, 16" Wide x 21" or 22" Overall Height ..... \$925 Part: ALBR16TA4. HCPCS: E2613	Active Contour, 20" Wide x 21" or 22" Overall Height ..... \$925 Part: ALBR20TA4. HCPCS: E2613
Active Contour, 16" Wide x 23" or 24" Overall Height ..... \$945 Part: ALBR16XA4. HCPCS: E2613	Active Contour, 20" Wide x 23" or 24" Overall Height ..... \$945 Part: ALBR20XA4. HCPCS: E2613
Active Contour, 18" Wide x 19" or 20" Overall Height ..... \$905 Part: ALBR18SA4. HCPCS: E2613	

**1d. ADI Back, AL Series, 4-Point, Deep Contour**

**Available with Back Cane Option only.** Select a back type and the specific overall back height in subsequent section. Previously selected back cane height must be equal to the greater of overall height range.

Deep Contour, 14" Wide x 13" or 14" Overall Height ..... \$965 Part: ALBR14LD4. HCPCS: E2615	Deep Contour, 18" Wide x 19" or 20" Overall Height ..... \$1,005 Part: ALBR18SD4. HCPCS: E2615
Deep Contour, 14" Wide x 16" or 17" Overall Height ..... \$985 Part: ALBR14MD4. HCPCS: E2615	Deep Contour, 18" Wide x 21" or 22" Overall Height ..... \$1,025 Part: ALBR18TD4. HCPCS: E2615
Deep Contour, 16" Wide x 19" or 20" Overall Height ..... \$1,005 Part: ALBR16SD4. HCPCS: E2615	Deep Contour, 20" Wide x 19" or 20" Overall Height ..... \$1,005 Part: ALBR20SD4. HCPCS: E2615
Deep Contour, 16" Wide x 21" or 22" Overall Height ..... \$1,025 Part: ALBR16TD4. HCPCS: E2615	Deep Contour, 20" Wide x 21" or 22" Overall Height ..... \$1,025 Part: ALBR20TD4. HCPCS: E2615
Deep Contour, 16" Wide x 23" or 24" Overall Height ..... \$1,045 Part: ALBR16XD4. HCPCS: E2615	Deep Contour, 20" Wide x 23" or 24" Overall Height ..... \$1,045 Part: ALBR20XD4. HCPCS: E2615

**1e. ADI Back, AL Series, 4-Point, Extra Deep Contour**

**Available with Back Cane Option only.** Select a back type and the specific overall back height in subsequent section. Previously selected back cane height must be equal to the greater of overall height range.

Extra Deep Contour, 16" Wide x 21" or 22" Overall Height .... \$1,125 Part: ALBR16T64. HCPCS: E2615	Extra Deep Contour, 18" Wide x 23" or 24" Overall Height .... \$1,145 Part: ALBR18X64. HCPCS: E2615
Extra Deep Contour, 16" Wide x 23" or 24" Overall Height .... \$1,145 Part: ALBR16X64. HCPCS: E2615	Extra Deep Contour, 20" Wide x 19" or 20" Overall Height .... \$1,105 Part: ALBR20S64. HCPCS: E2615
Extra Deep Contour, 18" Wide x 19" or 20" Overall Height .... \$1,105 Part: ALBR18S64. HCPCS: E2615	Extra Deep Contour, 20" Wide x 21" or 22" Overall Height .... \$1,125 Part: ALBR20T64. HCPCS: E2615
Extra Deep Contour, 18" Wide x 21" or 22" Overall Height .... \$1,125 Part: ALBR18T64. HCPCS: E2615	Extra Deep Contour, 20" Wide x 23" or 24" Overall Height .... \$1,145 Part: ALBR20X64. HCPCS: E2615

TRU-Balance 3 seating section continued on following pages.

**11. TRU-BALANCE 3 SEATING CONTINUED**

**1f. Q-Back 2.0, AL Series, Active Contour**

**Not available with Back Cane Option.** Select a back type and the specific overall back height in subsequent section.

Active Contour, 14"Wx19-20"H Tilt/22"-23" H Recline ..... \$9,85 Part: QB2AL14AS. HCPCS: E2613	Active Contour, 18"Wx21-22"H Tilt/24"-25" H Recline ..... \$1,145 Part: QB2AL18AT. HCPCS: E2613
Active Contour, 14"Wx21-22"H Tilt/24"-25" H Recline ..... \$985 Part: QB2AL14AT. HCPCS: E2613	Active Contour, 18"Wx23-24"H Tilt/26"-27" H Recline..... \$1,145 Part: QB2AL18AX. HCPCS: E2613
Active Contour, 16"Wx19-20"H Tilt/22"-23" H Recline ..... \$1,145 Part: QB2AL16AS. HCPCS: E2613	Active Contour, 18"Wx25-26"H Tilt/28"-29" H Recline..... \$1,145 Part: QB2AL18AXX. HCPCS: E2613
Active Contour, 16"Wx21-22"H Tilt/24"-25" H Recline ..... \$1,145 Part: QB2AL16AT. HCPCS: E2613	Active Contour, 20"Wx19-20"H Tilt/22"-23" H Recline..... \$1,145 Part: QB2AL20AS. HCPCS: E2613
Active Contour, 16"Wx23-24"H Tilt/26"-27" H Recline ..... \$1,145 Part: QB2AL16AX. HCPCS: E2613	Active Contour, 20"Wx21-22"H Tilt/24"-25" H Recline..... \$1,145 Part: QB2AL20AT. HCPCS: E2613
Active Contour, 16"Wx25-26"H Tilt/28"-29" H Recline ..... \$1,145 Part: QB2AL16AXX. HCPCS: E2613	Active Contour, 20"Wx23-24"H Tilt/26"-27" H Recline..... \$1,145 Part: QB2AL20AX. HCPCS: E2613
Active Contour, 18"Wx19-20"H Tilt/22"-23" H Recline ..... \$1,145 Part: QB2AL18AS. HCPCS: E2613	Active Contour, 20"Wx25-26"H Tilt/28"-29" H Recline..... \$1,145 Part: QB2AL20AXX. HCPCS: E2613

**1g. Q-Back 2.0, AL Series, Deep Contour**

**Not available with Back Cane Option.** Select a back type and the specific overall back height in subsequent section.

Deep Contour, 14"Wx19-20"H Tilt/22"-23" H Recline ..... \$1,145 Part: QB2AL14DS. HCPCS: E2620	Deep Contour, 18"Wx21-22"H Tilt/24"-25" H Recline ..... \$1,325 Part: QB2AL18DT. HCPCS: E2620
Deep Contour, 14"Wx21-22"H Tilt/24"-25" H Recline ..... \$1,145 Part: QB2AL14DT. HCPCS: E2620	Deep Contour, 18"Wx23-24"H Tilt/26"-27" H Recline ..... \$1,325 Part: QB2AL18DX. HCPCS: E2620
Deep Contour, 16"Wx19-20"H Tilt/22"-23" H Recline ..... \$1,325 Part: QB2AL16DS. HCPCS: E2620	Deep Contour, 18"Wx25-26"H Tilt/28"-29" H Recline ..... \$1,325 Part: QB2AL18DXX. HCPCS: E2620
Deep Contour, 16"Wx21-22"H Tilt/24"-25" H Recline ..... \$1,325 Part: QB2AL16DT. HCPCS: E2620	Deep Contour, 20"Wx19-20"H Tilt/22"-23" H Recline ..... \$1,325 Part: QB2AL20DS. HCPCS: E2620
Deep Contour, 16"Wx23-24"H Tilt/26"-27" H Recline ..... \$1,325 Part: QB2AL16DX. HCPCS: E2620	Deep Contour, 20"Wx21-22"H Tilt/24"-25" H Recline ..... \$1,325 Part: QB2AL20DT. HCPCS: E2620
Deep Contour, 16"Wx25-26"H Tilt/28"-29" H Recline ..... \$1,325 Part: QB2AL16DXX. HCPCS: E2620	Deep Contour, 20"Wx23-24"H Tilt/26"-27" H Recline ..... \$1,325 Part: QB2AL20DX. HCPCS: E2620
Deep Contour, 18"Wx19-20"H Tilt/22"-23" H Recline ..... \$1,325 Part: QB2AL18DS. HCPCS: E2620	Deep Contour, 20"Wx25-26"H Tilt/28"-29" H Recline ..... \$1,325 Part: QB2AL20DXX. HCPCS: E2620

**1h. Q-Back 2.0, AL Series, Extra Deep Contour**

**Not available with Back Cane Option.** Select a back type and the specific overall back height in subsequent section.

Extra Deep Contour, 14"Wx19-20"H Tilt/22"-23" H Recline ..... \$1,225 Part: QB2AL14ES. HCPCS: E2620	Extra Deep Contour, 18"Wx21-22"H Tilt/24"-25" H Recline ..... \$1,410 Part: QB2AL18ET. HCPCS: E2620
Extra Deep Contour, 14"Wx21-22"H Tilt/24"-25" H Recline ..... \$1,225 Part: QB2AL14ET. HCPCS: E2620	Extra Deep Contour, 18"Wx23-24"H Tilt/26"-27" H Recline ..... \$1,410 Part: QB2AL18EX. HCPCS: E2620
Extra Deep Contour, 16"Wx19-20"H Tilt/22"-23" H Recline..... \$1,410 Part: QB2AL16ES. HCPCS: E2620	Extra Deep Contour, 18"Wx25-26"H Tilt/28"-29" H Recline ..... \$1,410 Part: QB2AL18EXX. HCPCS: E2620
Extra Deep Contour, 16"Wx21-22"H Tilt/24"-25" H Recline..... \$1,410 Part: QB2AL16ET. HCPCS: E2620	Extra Deep Contour, 20"Wx19-20"H Tilt/22"-23" H Recline ..... \$1,410 Part: QB2AL20ES. HCPCS: E2620
Extra Deep Contour, 16"Wx23-24"H Tilt/26"-27" H Recline ..... \$1,410 Part: QB2AL16EX. HCPCS: E2620	Extra Deep Contour, 20"Wx21-22"H Tilt/24"-25" H Recline ..... \$1,410 Part: QB2AL20ET. HCPCS: E2620
Extra Deep Contour, 16"Wx25-26"H Tilt/28"-29" H Recline ..... \$1,410 Part: QB2AL16EXX. HCPCS: E2620	Extra Deep Contour, 20"Wx23-24"H Tilt/26"-27" H Recline ..... \$1,410 Part: QB2AL20EX. HCPCS: E2620
Extra Deep Contour, 18"Wx19-20"H Tilt/22"-23" H Recline ..... \$1,410 Part: QB2AL18ES. HCPCS: E2620	Extra Deep Contour, 20"Wx25-26"H Tilt/28"-29" H Recline ..... \$1,410 Part: QB2AL20EXX. HCPCS: E2620

**11.13. Lateral Positioning Wedge**

**Optional.**

Lateral Positioning Wedge, pair ..... \$32 Part: SET123465 Available with Sport and TRU-Comfort backs. Not available on power tilt backs 16"-17"H with limited adjustability on 18" H. Not available on power recline back 19"-20"H with limited adjustability on 21"H.
---

**11.14. TRU-Balance 3 Back Shroud**

**REQUIRED.** The back shroud is a feature of the TRU-Balance 3 seating system and is available in matching shroud colors. The back shroud is not available when Omit Back Pan/Cushion is selected. Reference colors shown on page 2.

Red Corvette Part: XRF000019	Grey Street Part: XRF000023	Pretty in Pink Part: XRF000027	Electric Blue Part: XRF000030	Green Machine Part: XRF000033
Back in Black Part: XRF000021	Deep Purple Part: XRF000025	Evolution Orange Part: XRF000028	Blue Jean Part: XRF000031	Military Green Part: XRF000034
Night Train Part: XRF000022	Raspberry Beret Part: XRF000026	Midnight Blue Part: XRF000029	Lemon Crush Part: XRF000032	

TRU-Balance 3 seating section continued on following page.

**11. TRU-BALANCE 3 SEATING CONTINUED**

**11.15. Back Width Package**

**REQUIRED** except when Omit Back Pan/Cushion is selected. Available with +/-4" of seat width.

Match Seat Width Selections  
Part: PTOINDV3468

14"W-16"W Package  
Part: XRF000037

18"W-22"W Package  
Part: XRF000039

10"W-14"W Package  
Part: XRF000036

16"W-20"W Package  
Part: XRF000038

20"W-24"W Package  
Part: XRF000040

**11.16. Back Width**

**REQUIRED** except when Omit Back Pan/Cushion is selected. Part number is PTOINDV plus number following # below. Available with +/-4" of seat width.

12" W (#3562)	14" W (#3563)	15" W (#3856)	16" W (#3564)	17" W (#3565)	18" W (#3566)	19" W (#3567)	20" W (#3568)	21" W (#3569)	22" W (#3570)	24" W (#3571)
------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------

**11.17. Overall Back Height**

**Static Seat and Power Tilt available height range: 14"-25"H. Odd heights include 4" gap and even heights include 3" gap between seat pan and bottom of back cushion. Power Recline available height range: 19"-28"H. Odd heights include 6" gap and even heights include 7" gap. See the Sport and TRU-Comfort back availability chart for stock order back size ranges.**

14" Height Part: SETHGHT1027	17" Height Part: SETHGHT1030	20" Height Part: SETHGHT1033	23" Height Part: SETHGHT1036	26" Height Part: SETHGHT1039
15" Height Part: SETHGHT1028	18" Height Part: SETHGHT1031	21" Height Part: SETHGHT1034	24" Height Part: SETHGHT1037	27" Height Part: SETHGHT1040
16" Height Part: SETHGHT1029	19" Height Part: SETHGHT1032	22" Height Part: SETHGHT1035	25" Height Part: SETHGHT1038	28" Height Part: SETHGHT1041

Overall Height	Tilt or Static Width												Recline Width											
	12"	14"	15"	16"	17"	18"	19"	20"	21"	22"	24"	12"	14"	15"	16"	17"	18"	19"	20"	21"	22"	24"		
14"	T	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	14"	
15"	T	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	15"
16"	S	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	16"
17"	S	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	17"
18"	-	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	18"
19"	-	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	19"
20"	-	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	20"
21"	-	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	21"
22"	-	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	22"
23"	-	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	23"
24"	-	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	24"
25"	-	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	25"
26"	-	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	26"
27"	-	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	27"
28"	-	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	28"
29"	-	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	29"
30"	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	30"
31"	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	31"
32"	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	32"

Sport Back     Q-Back     TRU-Comfort Back

**Sport, Q-Back, and TRU-Comfort Back Availability Chart**

**11.18. Power Tilt Privacy Flap**

**Optional.** A privacy flap is standard with power recline systems. This selection is to add a flap for power tilt only systems.

Privacy Flap ..... \$40  
Part: XRF000351

**12. TRU-BALANCE 3 ARMRESTS**

**REQUIRED.** Select the desired type of armrest below. Please read notes to ensure seating compatibility. Be sure to also select an armrest size and armpad type.

**12.1a. Flip-Back, 8"-16" H Armrests**

Choose a left and right armrest and armrest height in the following section. Not available with waterfall armpads. Not available with Back Cane Option.

Flip-Back, 8"-16"H Armrests .....\$145 ea.  
L R Left part: ARM1806662. Right part: ARM1806663

**Armrest Support**

**Optional.** Available up to 2" difference between seat and back width via endomorph or mesomorph configuration. Not available less than 16" width seat or back.

Flip-Back Armrest Downpost ..... \$115 ea.  
L R Left part: XRF000265. Right part: XRF000266

**12.1b. Drop-In, Height Adjustable, Removable, T-Arm Armrests**

**Only available with TRU-Balance 3 Static Seat.** Includes flat armpads.

Full Length, Drop-In Armrests ..... No Charge  
L R Left part: ARM1803253. Right part: ARM156742

Desk Length, Drop-In Armrests ..... No Charge  
L R Left part: ARM1803254. Right part: ARM156743



**12. TRU-BALANCE 3 ARMRESTS CONTINUED**

**12.1c. Two-Post, Height Adjustable, Removable Armrests**

Select the armrest and the left and right armrest heights in the following section. These armrests may interfere with offset backs/ back canes.

**1a. Two-Post Armrest**

L  R Full Length, Two-Post Armrest.....\$145 ea.  
Left part: XRF000051. Right part: XRF000052. HCPCS: E0973

L  R Desk Length, Two-Post Armrest .....\$145 ea.  
Left part: XRF000053. Right part: XRF000054. HCPCS: E0973

**12.1d. Heavy Duty, Removable Armrests**

HD removable armrests require armrest height and armrest size selections. Only one height package is selectable as a normal configuration. Payer coding, coverage, and bundling guidelines may apply.

8"-12"H, Single Post, Height Adjustable, Removable, .....\$130 ea.  
Heavy Duty (Sold as pair)  
Part: XRF000055. HCPCS w/ Static Seat: E0973

10"-14"H, Single Post, Height Adjustable, Removable,.....\$130 ea.  
Heavy Duty (Sold as pair)  
Part: XRF000056. HCPCS w/ Static Seat: E0973

**12.2. Armrest Height**

Set left and right armrest height. Note the available height range in the option description.

**1. Left Armrest Height**

8" H  
PTOINDV3535

10" H  
PTOINDV3537

12" H  
PTOINDV3539

14" H  
PTOINDV3541

16" H  
PTOINDV3543

9" H  
PTOINDV3536

11" H  
PTOINDV3538

13" H  
PTOINDV3540

15" H  
PTOINDV3542

**2. Right Armrest Height**

8" H  
PTOINDV3450

10" H  
PTOINDV3452

12" H  
PTOINDV3454

14" H  
PTOINDV3456

16" H  
PTOINDV3458

9" H  
PTOINDV3451

11" H  
PTOINDV3453

13" H  
PTOINDV3455

15" H  
PTOINDV3457

**12.3. Armrest Size**

Select a left (L) and right (R) armrest size. This also affects armpad selection.

L  R 10" Length  
Left part: XRF000057. Right part:  
XRF000123

L  R 14" Length  
Left part: XRF000058. Right part:  
XRF000124

L  R 17" Length  
Left part: XRF000059. Right part:  
XRF000125  
Only available with Flip-Back, 8"-16" H  
Armrests and straight arm pads.

**13. TRU-BALANCE 3 ARMPADS**

**REQUIRED.** Select to omit an armpad selection or select a left and right armpad type. Please note that waterfall and gel armpads may cause interference with the seat back.

**13.1a. Omit Armpad**

Select one below if aftermarket armpads or handpads are required.

No Armpad Selection..... No Charge  
Part: XRF000060  
Not available with Flip-back, 8"-16" H armrests.

Interface for Aftermarket Pad ..... No Charge  
Part: XRF000308  
Only available with Flip-back, 8"-16" armrests.

**13.1b. Foam Armpads**

L  R Straight, Armpad ..... Standard  
Left part: XRF000061. Right part: XRF000126

L  R Waterfall, Armpad ..... No Charge  
Left part: XRF000062. Right part: XRF000127  
Not available with Flip-back, 8"-16" armrests.

**13.1c. Gel Armpads**

Sold only in pairs. Armpads are 3" wide. Manufacturer part number for full length is 314G and desk length is 310G.

**a. Falcon Gel Armpads for Flip-Back Armrests**

Falcon Gel, Full Length Left and Right ..... \$300  
Part: ARM1806807

Falcon Gel, Desk Length Left, Full Length Right..... \$250  
Part: ARM1806810

Falcon Gel, Full Length Left, Desk Length Right..... \$250  
Part: ARM1806813

Falcon Gel, Desk Length Left and Right ..... \$250  
Part: ARM1806806

**b. Falcon Gel Armpads for Tube Style Armrests**

Falcon Gel, Full Length Left and Right ..... \$300  
Part: ARM136438

Falcon Gel, Desk Length Left, Full Length Right..... \$250  
Part: ARM136439

Falcon Gel, Full Length Left, Desk Length Right..... \$250  
Part: ARM136442

Falcon Gel, Desk Length Left and Right ..... \$200  
Part: ARM136437

**14. TRU-BALANCE 3 JOYSTICK BRACKETS**

**REQUIRED.**

Standard Fixed Receiver ..... Standard  
Part: XRF000063

Stealth Height Adjustable Receiver ..... \$194  
Part: XRF000064. HCPCS: K0108

**15. TRU-BALANCE 3 HEADRESTS**

**Optional.** Select a headrest pad and bracket.

**15.1a. Headrests, Comfort Plus with Cool Core**

**1. Headrest Pads**

- Stealth 8" Comfort Plus, Cool Core (CP08CCFPM) ..... \$279  
Part: ST-CP08CCFPM-Q, HCPCS: E0955
- Stealth 10" Comfort Plus, Cool Core (CP10CCFPM)..... \$279  
Part: ST-CP10CCFPM-Q, HCPCS: E0955

- Stealth 14" Comfort Plus, Cool Core (CP14CCFPM)..... \$279  
Part: ST-CP14CCFPM-Q, HCPCS: E0955
- Stealth 6", Comfort Plus, Cool Core (CP06CCFPM)..... \$279  
Part: ST-CP06CCFPM-Q, HCPCS: E0955

**2. Headrest Brackets**

- Comfort Plus Headrest Mounting Bracket, UniLink..... \$485  
Part: ST-UHB2D33200-Q, HCPCS: E1028
- Comfort Plus Headrest Mounting Bracket, HMO475.....\$271  
Part: ST-HMO475-17M-Q.

**16. LAP BELTS**

**REQUIRED.** Select a lap belt below.

**16.1a. Standard Lap Belts**

- |                                   |                                   |                                   |                                   |                                   |                                   |
|-----------------------------------|-----------------------------------|-----------------------------------|-----------------------------------|-----------------------------------|-----------------------------------|
| 38" Lap Belt<br>Part: ACCBELT1006 | 50" Lap Belt<br>Part: ACCBELT1000 | 60" Lap Belt<br>Part: ACCBELT1002 | 70" Lap Belt<br>Part: ACCBELT1003 | 80" Lap Belt<br>Part: ACCBELT1004 | 90" Lap Belt<br>Part: ACCBELT1001 |
|-----------------------------------|-----------------------------------|-----------------------------------|-----------------------------------|-----------------------------------|-----------------------------------|

**16.1b. Optional Lap Belts**

- Stealth Retractable Lap Belt (IFRB5350).....\$175  
Part: ST-IFRB5350
- Stealth Positioning, 60" Single Pull, Push Button, Medium..... \$134 (7.5") Padded, 1.5" Belt (IFPB5213)  
Part: ST-IFPB5213-Q
- Stealth Positioning, 60" Single Pull, Push Button, Large ..... \$134 Padded (9.5"), 2" Belt (IFPB5314)  
Part: ST-IFPB5314-Q
- Stealth Positioning, 60" Single Pull, Push Button, Extra Large. \$134 Padded (11.5"), 2" Belt (IFPB5315)  
Part: ST-IFPB5315-Q
- 60" Sure-Lok Lap Belt (Non-retractable)..... \$200  
Required with ISO 7176-19 Compliant Transit Option.  
Part: ACCBELT1016

**17. POWER LEGREST CONTROLS**

**Optional.** Specific power legrest options are selected later in this form. Select the control type here, if power legrests are desired, to assist the selection of the power positioning electronics.

- |  |   |
|--|---|
| Independent Legrests (Requires 2 actuator controls)<br>Part: XRF000108 | Articulating Foot Platform (Requires 1 Actuator Control)<br>Part: XRF000110 |
| Combined Legrests (Requires 2 actuator control)<br>Part: XRF000109     |   |

**18. LEGRESTS**

**Optional.** All efforts will be made to accommodate legrest length selection, but minor changes may need to be made to meet shipping guidelines. Some adjustments may be needed upon receipt of the unit.

**18.1a. Basic Center Mount Foot Platform**

Select legrest length in following section.

Short	Medium	Long
8"-10"	10"-14"	14"-19"
Basic Center Mount Extension Range		

**1. Center Mount Foot Platform Type**

- Basic Center Mount..... \$200  
Part: XRF000323

**2. Basic Center Mount Foot Platform Extension and Angle Adjustable Footplates**

See legrest length matrix for compatibility.

- |   |  |
|---|--|
| Short Extension, 5"x9" Footplates .....No Charge<br>Part: MEC1906675  | Medium Lower Extension, 6.25"W Tapered 5"W x 12"D... No Charge<br>Part: MEC2104853 |
| Medium Extension, 5"x9" Footplates .....No Charge<br>Part: MEC1906680 | Long Lower Extension, 6.25"W Tapered 5"W x 12"D..... No Charge<br>Part: MEC2104854 |
| Long Extension, 5"x9" Footplates .....No Charge<br>Part: MEC1906677   |  |

**3. Basic Center Mount Calf Pads**

Calf pads are standard with a power articulating foot platform.

- Small Calf Pads, Basic Center Mount..... \$210  
Part: RIG2104369
- Large Calf Pads, Basic Center Mount ..... \$210  
Part: RIG2202997  
Only available with Long Length Extension.

*Legrests section continued on following page.*

**18. LEGRESTS CONTINUED**

**4. Basic Center Mount Additional Accessories**

Basic Center Mount Footplate Connector ..... \$55  
 Part: RIG2104240

**18.1b. Center Mount Foot Platform**

If power articulating foot platform is selected, be sure power legrest control has been selected in previous section. Select legrest length in following section.

**1. Center Mount Foot Platform Type**

Power Articulating Foot Platform ..... \$3,145  
 Part: PTOINDV3754. HCPCS: E1012+K0040

Quick Adjustable Center Mount Foot Platform ..... \$265  
 Part: XRF000311. HCPCS: K0040

Medium	Long
11"-17"	14"-20"
Power AFP Extension Range	

Pediatric	Medium	Long
8.5"-12.5"	11"-16.5"	14"-20"
Quick Adj Center Mount Extension Range		

**2. Extension and Angle Adjustable Footplates**

See legrest length matrix for compatibility.

Pediatric Extension, 5"x6" Footplates Part: RIG152577 Not available with power articulating platform.	Medium Extension, 5"x8" Footplates with Width Extension Part: MEC1707469	Long Extension, 5"x8" Footplates with Width Extension Part: MEC1707474
Pediatric Extension, 5"x8" Footplates Part: RIG152578 Not available with power articulating platform.	Medium Extension, 5"x12" Footplates with Width Extension Part: MEC1707465	Long Extension, 5"x12" Footplates with Width Extension Part: MEC1707468
Medium Extension, 5"x6" Footplates Part: RIG150206	Long Extension, 5"x6" Footplates Part: RIG151505	

**3a. Center Mount Calf Support**

\*Calf pads are standard with a power articulating foot platform. Select calf pads covers for additional charge. Calf pad and panel covers include 1/2" of soft visco foam.

Small Calf Pads, 6"Wx6"H (pair) ..... \$210  
 Part: XRF000156. No Charge with Power AFP.

Large Calf Pads, 8" taper to 6"Wx8"H (pair) ..... \$210  
 Part: XRF000157. No Charge with Power AFP.

**3b. Center Mount Calf Panel**

Calf pads are standard with a power articulating foot platform.

Small Calf Panel, 18"Wx8"H ..... \$210  
 Part: FRMASMB7110007

Large Calf Panel, 24"Wx8"H ..... \$265  
 Part: RIG123824  
 Only available on widths of 20" or greater.

**4. Modular Footplate Extension Components**

Optional. Assembly required. May require seating configuration adjustments to avoid front caster interference.

4" Outside Edge (pair) ..... \$69  
 Part: RIG152510  
 Only available with 6"D footplates and 8"D footplate extension.

**5. Additional Accessories**

Small Reversed Dartex Calf Pad Cover (pair) ..... \$150 Part: RIG133717	Small Foam Panel Cover (pair) ..... \$150 Part: ST-CSCP0920-Q
Large Reversed Dartex Calf Pad Cover (pair) ..... \$150 Part: RIG133719	Large Foam Panel Cover (pair) ..... \$150 Part: ST-CSCP0927-Q
Small Reversed Dartex Panel Cover (pair) ..... \$150 Part: RIG134673	Footplate Connector ..... \$55 Part: RIGASMB7120023
Large Reversed Dartex Panel Cover (pair) ..... \$150 Part: RIG134744	Footplate Riser ..... \$150 Part: XRF000348 Not available with Long Extension. Select legrest length between 2" to 9" at 1" increments.
Small Foam Calf Pad Cover (pair) ..... \$150 Part: ST-CSCS0775-Q	
Large Foam Calf Pad Cover (pair) ..... \$150 Part: ST-CSCS0985-Q	

**6. Stealth Gel Footplate Covers**

5"x8" Tapered Gel Footplate Covers ..... \$420 Part: ST-FPG5008T-Q. HCPCS: K0108	6.5"x8" Tapered Gel Footplate Covers for Width Extension ..... \$420 Part: ST-FPG6508T-Q. HCPCS: K0108
5"x11" Tapered Gel Footplate Covers ..... \$420 Part: ST-FPG5011T-Q. HCPCS: K0108	5"x11" Tapered Gel Footplate Cover for Width Extension ..... \$420 Part: ST-FPG7011T-Q. HCPCS: K0108

Legrests section continued on following pages.

**18. LEGRESTS CONTINUED**

**18.1c. Power Elevating Legrests**

Be sure power legrest control has been selected in previous section. Select an upper extension type and footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes. Select legrest length in following section.

**1. Power Elevating Legrest Upper Assembly**

Includes calf pads.

Power Elevating Legrests, Pair..... \$2,865  
 Part: XRF000341. HCPCS: E1010+K0040+E0951

Legrest Length Matrix	Seat to Floor Option			Legrest Length Matrix	Seat to Floor Option		
	Low	Middle	High		Low	Middle	High
Power ELR: Static Seat, Recline Only				Power ELR: Tilt Only, Tilt and Recline, Any Power Lift			
8.5"-14.5" Lower Length	8.5"-12.5"	8.5"-13.5"	8.5"-14.5"	8.5"-14.5" Lower Length	8.5"-14.5"	8.5"-14.5"	8.5"-14.5"
12.5"-21" Lower Length	N/A	12.5"-13.5"	12.5"-14.5"	12.5"-21" Lower Length	12.5"-15.5"	12.5"-15.5"	12.5"-17"

**2. Lower Extension**

8.5" to 14.5" Lower Length                      12.5" to 21" Lower Length  
 Part: XRF000344                                      Part: XRF000347

**3. Angle Adjustable Footplate with Heel Loop and Wheel Bumper**

<sup>1</sup>The 6"W footplates may cause interference with seat widths less than 16". <sup>2</sup>The 8"W footplates may cause interference with seat widths less than 19".

4"Wx6"D, Pair Part: FRMASMB7448	5"Wx8"D, Pair Part: FRMASMB7450	6"Wx11"D <sup>1</sup> , Pair Part: FRMASMB7452
5"Wx6"D, Pair Part: FRMASMB7449	6"Wx8"D <sup>1</sup> , Pair Part: FRMASMB7451	8"Wx11"D <sup>2</sup> , Pair Part: FRMASMB7453

**18.1d. 60°-90° Swing-Away Legrests**

Adjustable in 5° increments. Select a legrest upper, a legrest extension, and type of footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes. Select legrest length in following section.

**1a. 60°-90° Swing-Away Legrest Upper**

P  L  R 60°-90° Swing-Away Legrest..... \$285  
 Pair part: RIG120062. Left: FRMASMB17558.  
 Right: FRMASMB17559  
 Available with Power Positioning only. Not available with 11"-11" or 11"-16.5" extensions.

P  L  R 60°-90° Swing-Away Legrest, Short..... \$285  
 Pair part: RIG136481. Left: RIG132607.  
 Right: RIG132609.  
 Available with Power Positioning only. For use only with 11"-13" or 11"-16.5" lower extensions.

**2. Swing-Away Legrests Lower Extension**

<input type="checkbox"/> P <input type="checkbox"/> L <input type="checkbox"/> R 11"-13" Length Pair part: XRF000088. Left: XRF000089. Right: XRF000090. N/A with Multi-axis.	<input type="checkbox"/> P <input type="checkbox"/> L <input type="checkbox"/> R 13"-16.5" Length Pair part: XRF000094. Left: XRF000095. Right: XRF000096. N/A with Multi-Axis.
---	---

<input type="checkbox"/> P <input type="checkbox"/> L <input type="checkbox"/> R 11"-16.5" Length Pair part: XRF000091. Left: XRF000092. Right: XRF000093. N/A with Angle Adjustable.	<input type="checkbox"/> P <input type="checkbox"/> L <input type="checkbox"/> R 15.5"-19.5" Length Pair part: XRF000097. Left: XRF000098. Right: XRF000099.
---	---

**3a. Extruded Footplate(s)**

P  L  R 7.25"Wx6"D Extruded Footplate..... Standard  
 Part pair: FRMASMB7847. Left: FRMASMB7861  
 Right: FRMASMB7862

P  L  R 8.5"Wx6"D Extruded Footplate ..... Standard  
 Part pair: FRMASMB7848. Left: FRMASMB7863  
 Right: FRMASMB7864  
 May cause interference on seat widths less than 18."

60°-90° Swing-Away Legrest Compatibility Matrix	Power Recline or Static Seat			Power Tilt			Power Lift		
	Low	Med.	High	Low	Med.	High	Low	Med.	High
Extruded, 11"-13"	!	✓	✓	✓	✓	✓	✓	✓	✓
Extruded, 13"-16.5"	!	✓	✓	✓	✓	✓	✓	✓	✓
Extruded, 15.5"-19.5"	!	!	!	!	!	!	!	!	✓
Angle Adjustable, 13"-16.5"	!	✓	✓	✓	✓	✓	✓	✓	✓
Angle Adjustable, 15.5"-19.5"	!	!	!	!	!	!	!	!	✓
Multi-Axis, 11"-16.5"	!	✓	✓	✓	✓	✓	✓	✓	✓
Multi-Axis, 15.5"-19.5"	!	!	!	!	!	!	!	!	✓

- This rigging configuration causes no interference issues at this seat-to-floor height.
- This rigging configuration may have less than 1.75" clearance when fully extended and may interfere with the ground or reduce obstacle climbing ability.

Legrests section continued on following pages.

**18. LEGRESTS CONTINUED**

**3b. Angle Adjustable Footplate(s) and Heel Loops**

<sup>1</sup>The 6"x11" footplates may cause interference with seat widths less than 18". <sup>2</sup>The 8"x11" footplates may cause interference with seat widths less than 20".

- 4"Wx6"D ..... \$87.50 ea.  
Pair part: FRMASMB10564. Left: RIG133537. Right: RIG133560
- 5"Wx6"D..... \$87.50 ea.  
Pair part: FRMASMB10565. Left: RIG133538. Right: RIG133561
- 5"Wx8"D ..... \$87.50 ea.  
Pair part: FRMASMB10566. Left: RIG133539. Right: RIG133562

- 6"Wx8"D<sup>1</sup> ..... \$87.50 ea.  
Pair part: FRMASMB10567. Left: RIG133540. Right: RIG133563
- 6"Wx11"D<sup>1</sup> ..... \$87.50 ea.  
Pair part: FRMASMB10568. Left: RIG133541. Right: RIG133564
- 8"Wx11"D<sup>2</sup> ..... \$87.50 ea.  
Pair part: FRMASMB10569. Left: RIG133542. Right: RIG133565

**3c. Multi-Axis Footplate(s) and Heel Loops**

<sup>1</sup>The 6"x11" footplates may cause interference with seat widths less than 18". <sup>2</sup>The 8"x11" footplates may cause interference with seat widths less than 20".

- 4"Wx6"D ..... \$260 ea.  
Pair part: FRMASMB9456. Left: RIG133531. Right: RIG133552
- 5"Wx6"D..... \$260 ea.  
Pair part: FRMASMB9457. Left: RIG133532. Right: RIG133553
- 5"Wx8"D ..... \$260 ea.  
Pair part: FRMASMB9458. Left: RIG133533. Right: RIG133554

- 6"Wx8"D<sup>1</sup> ..... \$260 ea.  
Pair part: FRMASMB9459. Left: RIG133534. Right: RIG133555
- 6"Wx11"D<sup>1</sup> ..... \$260 ea.  
Pair part: FRMASMB9460. Left: RIG133535. Right: RIG133558
- 8"Wx11"D<sup>2</sup> ..... \$260 ea.  
Pair part: FRMASMB9461. Left: RIG133536. Right: RIG133559

**4. 60°-90° Swing-Away Legrest Accessories**

- Gel Neoprene Footrest Sleeve (pair) ..... \$130  
Part: INDPART2243. HCPCS: K0108

**18.1e. Elevating Legrests**

Select a legrest upper, a legrest extension, and type of footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes. Select legrest length in following section.

**1. Elevating Legrests Upper**

- Elevating Legrest Assembly ..... \$275 ea.  
Pair part: FRMASMB6483. Left: RIG138012  
Right: RIG138013. HCPCS: E0990  
Select lower extension and footplates

**2a. Lower Extension**

- 11"-14" Length  
Pair: FRMASMB4755. Left: RIG138026. Right: RIG138030  
For reinforced only
- 8.5"-12" HD Drop-Ins/12.5"-16.5" ELR, Ext Length  
Pair part: XRF000101. Left: XRF000102. Right: XRF000103  
For angle adjustable only
- 14"-17" Length  
Pair: FRMASMB4756. Left: RIG138033. Right: RIG138035  
For reinforced only
- 13"-18" HD Drop-Ins/15.5"-19.5" ELR, Ext Length  
Pair part: XRF000104. Left: XRF000105. Right: XRF000106  
For angle adjustable only

**3a. Reinforced Cast Footplate(s)**

- 7.25"Wx6"D Reinforced Footplate..... Standard  
Pair part: FRMASMB4748. Left: RIG138055. Right: RIG138056

- 8.5"Wx6"D Reinforced Footplate ..... Standard  
Pair part: FRMASMB4747. Left: RIG138060. Right: RIG138062  
May cause interference on seat widths less than 18."

**3b. Angle Adjustable Footplate(s)**

<sup>1</sup>The 6"x11" footplates may cause interference with seat widths less than 18". <sup>2</sup>The 8"x11" footplates may cause interference with seat widths less than 20". 4"x6" and 5"x6" footplates are only available in pairs.

- 4"Wx6"D ..... \$87.50 ea.  
Pair parts: FRMASMB8812 x2.
- 5"Wx6"D..... \$87.50 ea.  
Pair parts: FRMASMB8811 x2.
- 5"Wx8"D Footplates with Heel Loops ..... \$87.50 ea.  
Pair part: FRMASMB4735. Left: RIG137791. Right: RIG137792

- 6"Wx8"D<sup>1</sup> Footplates with Heel Loops ..... \$87.50 ea.  
Pair part: FRMASMB4738. Left: RIG137793. Right: RIG137795
- 6"Wx11"D<sup>1</sup> Footplates with Heel Loops ..... \$87.50 ea.  
Pair part: FRMASMB4741. Left: RIG137796. Right: RIG137797
- 8"Wx11"D<sup>2</sup> Footplates with Heel Loops ..... \$87.50 ea.  
Pair part: FRMASMB4744. Left: RIG137798. Right: RIG137799

Elevating Legrest Compatibility Matrix	Power Recline or Static Seat			Power Tilt			Power Lift		
	Low	Med.	High	Low	Med.	High	Low	Med.	High
Reinforced, 11"-14"	✓	✓	✓	✓	✓	✓	✓	✓	✓
Reinforced, 14"-17"	!	!	✓	!	✓	✓	✓	✓	✓
Angle Adjustable, 12"-16.5"	✓	✓	✓	✓	✓	✓	✓	✓	✓
Angle Adjustable, 15.5"-19.5"	!	!	!	!	!	✓	!	✓	✓

This rigging configuration causes no interference issues at this seat-to-floor height.  
 This rigging configuration may have less than 1.75" clearance when fully extended and may interfere with the ground or reduce obstacle climbing ability.

Legrests section continued on following page.

**18. LEGRESTS CONTINUED**

**18.2. Legrest Accessories**

Heel Loops (pair)..... \$35      Wheel Bumper..... \$15 ea.  
 Part: FRMASMB7873. HCPCS: E0951      Part: XRF000100  
 Standard on angle adjustable and multi axis footplates when noted.

**18.3. Set Legrest Extension Length**

**REQUIRED with matching left and right legrest configuration. Part number is PTOINDV plus number following # below. Non-matching or individual legrest configurations lengths are not set by the manufacturer.**

**1a. Centermount Extension Length**

2" L (#3989)	5" L (#3992)	8" L (#3680)	9.5" L (#3683)	11" L (#3686)	12.5" L (#3689)	14" L (#3692)	15.5" L (#3695)	17" L (#3698)	18.5" L (#3701)
3" L (#3990)	6" L (#3993)	8.5" L (#3681)	10" L (#3684)	11.5" L (#3687)	13" L (#3690)	14.5" L (#3693)	16" L (#3696)	17.5" L (#3699)	19" L (#3702)
4" L (#3991)	7" L (#3994)	9" L (#3682)	10.5" L (#3685)	12" L (#3688)	13.5" L (#3691)	15" L (#3694)	16.5" L (#3697)	18" L (#3700)	

**1b. Front Rigging Extension Length**

Non-matching or individual legrest configurations lengths are not set by the manufacturer. See previous section for center mount extension lengths.

7" L (#3470)	8.5" L (#3473)	10" L (#3476)	11.5" L (#3479)	13" L (#3482)	14.5" L (#3485)	16" L (#3488)	17.5" L (#3491)	19" L (#3494)	20.5" L (#3497)	22" L (#3669)
7.5" L (#3471)	9" L (#3474)	10.5" L (#3477)	12" L (#3480)	13.5" L (#3483)	15" L (#3486)	16.5" L (#3489)	18" L (#3492)	19.5" L (#3495)	21" L (#3498)	22.5" L (#3670)
8" L (#3472)	9.5" L (#3475)	11" L (#3478)	12.5" L (#3481)	14" L (#3484)	15.5" L (#3487)	17" L (#3490)	18.5" L (#3493)	20" L (#3496)	21.5" L (#3600)	

**19. TRU-BALANCE 3 POWER POSITIONING ELECTRONICS**

**REQUIRED** with power seating or legrest options.

**19.1a. Tilt Only through Toggle**

Single Actuator Control, Power Tilt Only, Through Toggle (NE 4-Key Joystick Only)..... No Charge  
 Part: ELE1905907  
 Select toggle side. May require extended lead time when selected with Fender Lights.

**1. Toggle Mount Side Selection**

Left Side Toggle      Right Side Toggle  
 Part: PTOINDV1839      Part: PTOINDV1838

**19.1b. NE+ 6 Key Controls through Joystick**

Single or Dual Actuator Control Through NE+ Joystick..... No Charge  
 Part: XRF000065

**19.1c. Q-Logic 3 Controls through Joystick**

Single Actuator Through Q-Logic 3 Controller and Interface (AM3)..... No Charge

2 Actuators Through Q-Logic 3 Controller and Interface (AM3).....\$1,250  
 Not funded by ADP

3 Actuators Through Q-Logic 3 Controller and Interface (AM3)..... \$1,500  
 Not funded by ADP

Up to 5 Actuators Through Q-Logic 3 Controller and Interface (AAM).....\$1,675  
 ADP: SECX15005

**1. Optional Switch Control**

Available with Q-Logic 3 only. 4-Way switches will be mounted to same side of joystick.

4-Function Push Button.....\$225 Part: CTL123218 See next section for switch mounting. Only available with AAM.	4-Way Toggle.....\$575 Part: ELEASMB5925 Only available with AAM.	4-Way Button Toggle.....\$575 Part: ELEASMB5926 Only available with AAM.
iAccess Flat Button.....\$575 Part: CTL144009 See next section for switch mounting.	iAccess Toggle Button.....\$575 Part: CTL144010 See next section for switch mounting.	



TRU-Balance 3 power positioning electronics section continued on following page.

**19. TRU-BALANCE 3 POWER POSITIONING ELECTRONICS CONTINUED**

**19.2. 4-Function Push Button/iAccess Mount**

**REQUIRED with 4-Function Push Button or iAccess options.** Under arm mount is not available with TRU-Balance 3 HD armrests.

Joystick Mount ..... No Charge Part: MEC1907500	Attendant Mount ..... \$35 Part: MEC1906564
Under Arm Mount, Flip-Back Armrest..... No Charge Part: MEC1906531	Under Arm Mount, Flip-Back Armrest, plus iAccess Guard ..... \$35 Part: MEC2001400
Under Arm Mount, Tube Style Armrest ..... No Charge Part: MEC1906530	Under Arm Mount, Tube Style Armrest, plus iAccess Guard ..... \$35 Part: MEC2001401

**1. Under Arm Mount Side Selection**

Left Side Under Arm Mount Part: PTOINDV3678	Right Side Under Arm Mount Part: PTOINDV3679
--	---

**20. TRU-BALANCE 3 ATTENDANT CONTROL**

**Optional.** \*Attendant control will be at no charge when ordered with Tilt and Recline. Mounting bracket is included.

Attendant Control ..... \$950 Part: XRF000078. HCPCS: E2331
--

**21. SEAT MOUNTED LED LIGHTS**


**Optional.** Requires Q-Logic 3 Joystick for Full Light Package or Q-Logic 3e Joystick.

Seat Mounted, LED, Full Light Package ..... \$770 Part: XRF000130 Not available with Base Mounted, LED, Full Light Package.
---

**22. PROVIDER INSTALLED ACCESSORIES**

**Optional.** Multiple accessories can be mounted along the various tracks but some interference may occur.

Accessory Bar ..... \$90 Part: XRF000128 Horizontal rail attached to rear of seat. Standard with the Trilogy vent tray and oxygen holder.	Personal Item Hook..... \$20 Part: ACC125008
Backpack Holder (Pair) (A)..... \$42 Part: ACC147774	Personal Item Pouch ..... \$175 Part: ACC2101017
Cane/Crutch Holder ..... \$100 Part: ACC152515. HCPCS: E2207	Privacy Bag..... \$80 Part: ACC2005723
Cell Phone/ Media Holder, X-Grip (E)..... \$95 Part: ACC155835 Not available with Heavy Duty armrests.	Push Handles (pair)..... \$140 Part: ACC2005027 Not available with oxygen holder or vent tray.
Clothing Guard (Pair)..... \$275 Part: ACC125007	Q-Logic EX/NE Programmer.....\$625 Part: ELEASMB6976
Compact Vent Tray .....\$1,395 Part: ACC133889. HCPCS: E1029	Quantum Backpack.....\$200 Part: ACC1905853
Cup Holder ..... \$20 Part: ACC125003 For use with track-style armrests.	Rear View Mirror.....\$35 Part: ACC155833
Cup Holder (Tube-style Armrest)..... \$20 Part: ACCASMB2634 For use with tube-style armrests.	Reflector Kit.....\$35 Part: ACC124670
Flag Kit ..... \$125 Part: ACC1700443	Stealth Large Tablet Holder.....\$565 Part: ACC141818
Flip-Back Armrest, Side Track Accessory Bracket .. \$25 Part: MEC1809178	Transfer Handle (D).....\$152.50 ea. Part: MEC2205680 x2. HCPCS: E0705
Glove Box (C) \$85 Part: ACC125011	Trilogy Vent Tray.....\$420 Part: ACC136791. HCPCS: E1029 Includes accessory bar, a credit will be applied PTOINDV3769. Contact custom for models other than 100 and 200.
Large, Black, Weather Cover .....\$240 Part: ACCCOVR1021	VOCSN Vent Tray for Accessory Bar Mount.....\$1,000 Part: ACC2000166. HCPCS: E1029 Includes accessory bar, a credit will be applied PTOINDV3769.
Oxygen Holder..... \$110 Part: ACC125006. HCPCS: E2208 Includes accessory bar, a credit will be applied PTOINDV3769.	Depth Adjustable Vent Tray.....\$1,000 Part: ACC2302500. HCPCS: E1029 Includes accessory bar, a credit will be applied PTOINDV3769.



**24. POSITIONING COMPONENTS**

**Optional.** More options are available on the [Positioning Components Order Form](#).

**24.1. Stealth Thoracic Laterals**

Swing-away laterals may interfere with TRU-Balance 3 armrests and 2-post, flip up armrests. These laterals are not available with the back cane option.

**1. Bracket**

Swing-away, UniLink, Redesign Back, pair (ULS4042000-Q).. \$1,032  
 Part: ST-ULS4042000-Q. HCPCS: E1028  
 Only available with Redesigned back assemblies.

**2. Pads**

Curved Pads, 3"x5", pair (TWBL3X5)..... \$164  
 Part: INDPART2663. HCPCS: E0956

Flat Pads, 3"x5", pair (TWBLP3X5) ..... \$164  
 Part: INDPART2676. HCPCS: E0956

Curved Pads, 4"x6", pair (TWBL4X6)..... \$164  
 Part: INDPART2664. HCPCS: E0956

Flat Pads, 4"x6", pair (TWBLP4X6) ..... \$164  
 Part: INDPART2675. HCPCS: E0956

**3. Installation**

Install Laterals  
 Part: PTOINDV3499

Do Not Install  
 Part: PTOINDV3500

**24.2. Stealth Pelvic/Thigh Guides**

Pelvic/thigh guides may interfere with the joystick receiver block and down posts. Please provide quantity next to check box.

**1. Pelvic/Thigh Guide Bracket**

Flip Down Bracket (TWBTLADD-HW) ..... \$343 ea.  
 Part: INDPART2804. HCPCS: E1028

Flip Down with LEFT Full Surface Contact ..... \$459  
 Feature (TWBADD-L & LPHW407)  
 Part: POS140428. HCPCS: E1028

Removable Bracket (LP430)..... \$383 ea.  
 Part: INDPART2805. HCPCS: E1028

Flip Down with RIGHT Full Surface Contact..... \$459  
 Feature (TWBADD-R & LPHW407)  
 Part: POS140429. HCPCS: E1028

Removable Bracket with 1" Offset ..... \$383 ea.  
 (LP430 & LHW-121)  
 Part: INDPART2854. HCPCS: E1028

**2. Pelvic/Thigh Guide Adapter**

2" Height Adapter Plate (LHW-122)..... \$35  
 Part: INDPART2821. HCPCS: K0108  
 Height adapter is needed with TRU-Comfort seat cushions. Not available with the flip down bracket.

2 1/4" Adapter Plate (LHW-129) ..... \$82  
 Part: POS140430. HCPCS: K0108  
 Only available with the flip down full surface bracket.

**3. Pelvic/Thigh Guide Pad**

Short Pad, 4"x4" (TWBADD-4X4)..... \$85 ea.  
 Part: INDPART2801. HCPCS: E0953/E0956

Long Pad, 4"x12" (TWBADD-4X12) ..... \$85 ea.  
 Part: INDPART2803. HCPCS: E0953/E0956

Medium Pad, 4"x8" (TWBADD-4X8)..... \$85 ea.  
 Part: INDPART2802. HCPCS: E0953/E0956

**4. Pelvic/Thigh Guide Installation**

Install on Left  
 Part: PTOINDV3585

Install on Right  
 Part: PTOINDV3586

Do Not Install  
 Part: PTOINDV3719

**24.3. Stealth UniLink Pelvic/Thigh Guide**

**1. UniLink Pelvic/Thigh Hardware**

UniLink, 5" Height Adjustable, ..... \$405 ea.  
 Removable Mounting Hardware  
 Part: ST-UTTP5-Q. HCPCS: E1028

UniLink, 7" Height Adjustable, ..... \$423 ea.  
 Removable Mounting Hardware  
 Part: ST-UTTP7-Q. HCPCS: E1028

**2. UniLink Pelvic/Thigh Pads**

UniLink, 4" x 4" Pad ..... \$164 ea.  
 Part: ST-UTP544056100-Q. HCPCS: E0953/E0956

UniLink, 4" x 12" Pad ..... \$164 ea.  
 Part: ST-UTP54D056100-Q. HCPCS: E0953/E0956

UniLink, 4" x 8" Pad..... \$164 ea.  
 Part: ST-UTP548056100-Q. HCPCS: E0953/E0956

**3. UniLink Pelvic/Thigh Additional Links**

UniLink, 1.5" Link ..... \$121 ea.  
 Part: ST-UIL001150A000-Q

UniLink, 2" Link..... \$121 ea.  
 Part: ST-UIL001200A000-Q

**4. UniLink Pelvic/Thigh Guide Installation**

Install on Left  
 Part: PTOINDV3585

Install on Right  
 Part: PTOINDV3586

Do Not Install  
 Part: PTOINDV3719



Positioning components section continued on following page.



**24. POSITIONING COMPONENTS CONTINUED**

**24.4. Stealth Positioning Chest Harness**

Upper Classic, Medium (IFPS1103) ..... \$221  
 Part: ST-IFPS1103-Q. HCPCS: E0960  
 Upper Classic, Large (IFPS1104) ..... \$221  
 Part: ST-IFPS1104-Q. HCPCS: E0960

Upper Contour, Medium (IFPS2103)..... \$221  
 Part: ST-IFPS2103-Q. HCPCS: E0960  
 Upper Contour, Large (IFPS2104) ..... \$221  
 Part: ST-IFPS2104-Q. HCPCS: E0960



Contoured and Classic Chest Harness

**24.5. Stealth Positioning Chest Straps**

Chest Strap, Medium w/ Belt Anchor, ..... \$140 (IFCS4113)  
 Part: ST-IFCS4113-Q. HCPCS: E0960

Chest Strap, Large w/ Belt Anchor, ..... \$140 (IFCS4114)  
 Part: ST-IFCS4114-Q. HCPCS: E0960



Chest Strap

**24.6. Stealth Swing-Away Headrest Mounts**

These swing-away headrest mounts are only available with the Stealth Comfort Series headrests.

**1. Facial Lateral Hardware**

UniLink Facial Lateral, Swing Away, Right..... \$496  
 Part: ST-SAR2A-R-Q. HCPCS: E1028  
 UniLink Facial Lateral, Swing Away, Left ..... \$496  
 Part: ST-SAR2A-L-Q. HCPCS: E1028  
 UniLink Facial Lateral, Dual Swing Away ..... \$765  
 2 Upper & 1 Lower Right  
 Part: ST-SAR4A-DBU-RL-Q. HCPCS: E1028

UniLink Facial Lateral, Dual Swing Away ..... \$765  
 2 Upper & 1 Lower Left  
 Part: ST-SAR4A-DBU-LL-Q. HCPCS: E1028  
 UniLink Facial Lateral, Dual Swing Away ..... \$674  
 2 Upper  
 Part: ST-SAR4A-DBU-Q. HCPCS: E1028




Left and Right Facial Lateral with Pads

**2. Facial Lateral Pads and Accessories**

Egg Switch Mount (ESM) (A) ..... \$70 ea.  
 Part: INDPART2815  
 Figure 8 Bracket with 13" Rod..... \$366 ea.  
 and Egg Switch Mount (B)  
 Part: ST-F8675-13-Q  
 Facial Pad, 4".5x2", Cool Core (900M-Q) .... \$96 ea.  
 Part: ST-900M-Q. HCPCS: K0108

Facial Pad, 6"x2", Cool Core (920M-Q)..... \$96 ea.  
 Part: ST-920M-Q. HCPCS: K0108  
 Swing-Away Joystick Mount, 18" Rod (C)..... \$100 ea.  
 Part: ST-SUS18-Q  
 Allows mounting of facial pads and attaches to the headrest




Headrest Accessories

**24.7. Stealth Elbow Support**

**1. Elbow Support Bracket**

Left, Swing-Away, Elbow Support ..... \$279  
 Bracket Tube Style Arms (ARMS-ES-L)  
 Part: POS124064. HCPCS: E1028  
 Left, UniLink, Fixed, Elbow Support Bracket..... \$309  
 for Flip-Back, 8"-16"H Armrests (UETB4021FL-Q)  
 Part: ST-UETB4021FL-Q

Right, Swing-Away, Elbow Support..... \$279  
 Bracket Tube Style Arms (ARMS-ES-R)  
 Part: POS124066. HCPCS: E1028  
 Right, UniLink, Fixed, Elbow Support Bracket..... \$309  
 for Flip-Back, 8"-16"H Armrests (UETB4021FR-Q)  
 Part: ST-UETB4021FR-Q



Elbow Supports

**2. Elbow Support Pad**

3"Wx3"H Elbow Pad (ARMS-ES-3X3) ..... \$57 ea.  
 Part: POS124065. HCPCS: K0108

4"Wx4"H Elbow Pad (ARMS-ES-4X4)..... \$57 ea.  
 Part: POS124067. HCPCS: K0108

**24.8. Stealth Foot Boxes**

Stealth Small Foot Boxes, 5"Wx8"Dx7"H..... \$725 (SPLE107) (pair)  
 Part: INDPART2786. HCPCS: E0954  
 Stealth Medium Foot Boxes, 6"Wx10"Dx9"H ..... \$725 (SPLE108) (pair)  
 Part: INDPART2787. HCPCS: E0954

Stealth Large Foot Boxes, 7"Wx12"Dx11"H..... \$725 (SPLE109) (pair)  
 Part: INDPART2788. HCPCS: E0954

**24.9. Stealth Positioning Shoe Holder**

Small Shoe Holder, 3/4" Webbing (IFSH4211)..... \$93 ea.  
 Part: ST-IFSH4211-Q. HCPCS: E0951+E0952  
 Medium Shoe Holder, 3/4" Webbing (IFSH4212) ..... \$93 ea.  
 Part: ST-IFSH4212-Q. HCPCS: E0951+E0952

Large Shoe Holder, 3/4" Webbing (IFSH4213) ..... \$93 ea.  
 Part: ST-IFSH4213-Q. HCPCS: E0951+E0952

**25. ALTERNATIVE CONTROLS AND ELECTRONICS**

**Optional.** Select a type of specialty controls below. More options are available on the [Alternative Controls and Electronics order form](#).

**25.1. Enhanced Display**

**Q-Logic 3 Enhanced Display Kit**..... \$2,045

Part: CTL167726

Select a color option below. Includes enhanced display, flexible gooseneck and mounting bracket. Features: 3.5" LCD display; up to 5 profiles to control drive, seat, and auxiliary functions; Bluetooth connectivity for mobile devices and mouse emulation; integrated single switch scanner; IR controls.



Enhanced Display

**Q-Logic 3 Enhanced Display Shroud**

Required with Q-Logic 3 Enhanced Display selection above. Reference colors shown on page 2.

**Black, Soft Touch** ..... Standard

Part: DGN137728

**Deep Purple**

Part: DGN137729-PPR5

**Electric Blue**

Part: DGN137729-PBU9

**Red Corvette**

Part: DGN137729-PRD1

**Raspberry Beret**

Part: DGN137729-PPK4

**Blue Jean**

Part: DGN137729-PBU11

**Back in Black**

Part: DGN137729-PBK1

**Pretty in Pink**

Part: DGN137729-PPK3

**Lemon Crush**

Part: DGN137729-PYL4

**Night Train**

Part: DGN137729-PBK3

**Evolution Orange**

Part: DGN137729-POR7

**Green Machine**

Part: DGN137729-PGN11

**Grey Street**

Part: DGN137729-PSL3

**Midnight Blue**

Part: DGN137729-PBU10

**Military Green**

Part: DGN137729-PGN8

**Military Digital Camo**..... \$50

Part: DGN137729-WCM7

**Faded Glory** ..... \$50

Part: DGN137729-WFL2

**Carbon Fiber** ..... \$50

Part: DGN137729-WCB1

**25.2. Q-Logic SCIM Module**

**Optional.**

**SCIM Module for Q-Logic 3** ..... \$850

Part: CTL2001963

Allows use of joystick display for specialty devices. Available with Q-Logic 3 EX Joystick selection above. Up to 5 profiles to control drive, seat, and auxiliary functions, excluding sip-n-puff and single switch functions. A switch is required to operate the input device attached to the SCIM module.



SCIM Module

**25.3. Q-Logic 3 SAJ**

**Q-Logic 3 EX SAJ**..... \$850

Part: CTL143980

Must be used in conjunction with Enhanced Display or Q-Logic 3 Joystick.



EX SAJ

**25.4. Head Array**

**i-Drive 4.0 Bluetooth Tri-array Series**,..... \$7,406

**3 Switch with Egg Switch**

Part: ST-IDHBT200-U-Q. HCPCS: E2330

**Adult, Curved, i-Drive Tri-Array with** ..... \$7,406

**Mounting Hardware**

Part: ST-IDHBT400-U-Q. HCPCS: E2330

**i-Drive 4.0 Bluetooth Pro Series**,..... \$7,486

**3 Switch with Egg Switch**

Part: ST-IDHBT300-U-Q. HCPCS: E2330

**Junior, Curved, i-Drive Tri-Array with** ..... \$7,406

**Mounting Hardware**

Part: ST-IDHBT450-U-Q. HCPCS: E2330



Adult Pro-Series

**25.5. Sip and Puff Drive**

**1. Module**

**Sip and Puff Module for Q-Logic 3**..... \$1,695

Part: CTL168056. HCPCS: E2325

Includes module and mounting hardware.



Sip & Puff Module

**2. Straw and Mount**

The straw and mount do not include electrical components to operate the Sip&Puff with a head array. This component requires the selection of a Q-Logic Sip&Puff module or additional module from Stealth Products for operation.

**Stealth Sip and Puff for Stealth i-Drive**, ..... \$496

**Head Array (SP625-U)**

Part: ST-SP625-U-Q. HCPCS: E2326

Includes straw and mounting hardware. Use with Stealth i-Drive Head Arrays.



Stealth Sip & Puff

**3. Sip & Puff Accessories**

**Sip & Puff Straws, 12 pack (SP625-S12)**..... \$35

Part: ST-SP625-S12-Q

Alternative controls and electronics section continued on following page.

**25. ALTERNATIVE CONTROLS AND ELECTRONICS CONTINUED**

**25.6. Egg Switches**

Red Egg Switch (ESRED)..... \$116  
*Part: SWTMCRO1077*  
 \_\_\_\_\_  
 QTY

Blue Egg Switch (ESBLU)..... \$116  
*Part: SWTMCRO1079*  
 \_\_\_\_\_  
 QTY

Black Egg Switch (ESBLK)..... \$116  
*Part: SWTMCRO1081*  
 \_\_\_\_\_  
 QTY

Green Egg Switch (ESGRN) ..... \$116  
*Part: SWTMCRO1078*  
 \_\_\_\_\_  
 QTY

Yellow Egg Switch (ESYEL) ..... \$116  
*Part: SWTMCRO1080*  
 \_\_\_\_\_  
 QTY



**25.7. Switch Adapter**

2 to 1 Switch Adapter..... \$135  
*Part: DWR1361H037*  
 Connects 2 single switches to a stereo jack to act as a dual switch.

**25.8. Environmental Controls**

Q-Logic 3 ECU Module Kit..... \$1,780  
*Part: CTL168057*