

Prime Plus Care Bed Model P750

HCPCS:



- 750 lbs weight capacity
- 700 lbs Net weight capacity
- 11" low height 30" high working height
- Push button integrated width expansion from 39" to 48" width
- Integrated 80", 84" and 88" sleep surface length helps comfortably accommodate the tallest residents
- One touch floor lock
- Quiet, smooth operation/Easy-to-use Hand Pendant
- Cardiac Comfort Chair
- 4 heavy duty locking casters roll easy at any height and lock secure
- Fixed adjustable wall bumper
- Warranty: 15 years on frame lifetime on welds 5 years on electronics
- Options:
- Deluxe Assist Rails, Assist Bars
- Bed ends available in 2 styles and 3 colors
- Trapeze
- Battery Backup
- Bed Transporter

PRODUCTS	Description
P750	LTC Bariatric Bed 750 lbs

SPECIFICATIONS	Description
Built-in Width Expansion	39", 42", 48"
Height range	11"-30"
One Step Locking	Yes
Sleep Surface	80" - 84" - 88"
Warranty	15 Year on Frame, 5 Year on Electrical and Mechanical. Lifetime on All Welds
Weight Capacity	750 lbs

ACCESSORIES	Description
PBASR	Fold Rails f/1750/52;2000/02 BariBed,PR, PR

PARTS	Description
SP01-TA31-1060-007	P750 leg actuator, 1/EA, EA
SP01-P750-10016E	P750 brake pedal assy, 1/EA, EA
SP01-P907-10E	P750 wall-bumper assy top level, 1/EA, EA
SP01-P750-10001E	P750 legs assy 1-1, 1/ea, EA

PARTS	Description
SP01-P907-05E	P750 head board assy, 1/PR, PR
SP01-P903-10006E	Head Board Bracket Knobs , P903, 1/pr, PR
SP01-P907-06E	P750 foot board assy, 1/PR, PR
SP01-P750-10002E	P750 legs assy 2-1, 1/EA, EA
SP01-P750-10003E	P750 legs assy 2-3, 1/EA, EA
SP01-P750-10004E	P750 head deck, 1/EA, EA
SP01-P750-10007E	P750 head/thigh deck bolt/washers, 1/pkg, EA
SP01-P750-10005E	P750 thigh deck, 1/EA, EA
SP01-P750-10006E	P750 Foot Deck, 1/EA, EA
SP01-P750-10008E	P750 foot deck bolt/washers, 1/PKG, EA
SP01-P902-10018E	P750 bearing spare part,1/PR, PR
SP01-P750-10014E	P750 bed ends plastic assy, 1/PKG, EA
SP01-P907-13E	P750 Brake Assy, 1/EA, EA
SP01-PO10040209E-G	Brake Pedal, Green, P903, 1/ea, EA
SP01-PO10040209E-R	Brake Pedal, Red, P903, 1/ea, EA
SP01-P750-10017E	P750foot end truck assembly RT, 1/ea, EA