




Prevention Through  
Innovation.

### Addressing *Shear* from the Inside Out



Low friction material separates two layers of alternative foam - absorbs and dissipates shear forces.




## Align Mattress

with *Shear Reduction Technology*

The **Align** reactive foam replacement surface with patented Shear Reduction Technology (SRT) is an innovative solution designed to address shearing forces present during client movement in bed or when changing the position of the bed frame, at the same time ensuring proper pressure redistribution. The Align is a cost effective choice for clients at risk of experiencing skin breakdown. The side bolsters complete with accordion cuts provide an excellent transfer edge.

The surface is designed to absorb forces that damage skin, thereby reducing the risk for skin breakdown. A third party clinical trial was completed by Durham College Department of Nursing and Department of Engineering to validate SRT. The trial provided non-biased support for the efficacy and use in populations where shear forces can damage vulnerable skin.

*5 Minds Mobility received the Colleges of Ontario Network for Industry Innovations Award for SRT.*



**5 MINDS MOBILITY GROUP**

ggarland@5mindsmobility.com

orders@5mindsmobility.com

555 Washington Ave., Suite 310

St. Louis, MO 63101

Ph: 866-709-1332 Fax: (877)-426-5697




Prevention Through  
Innovation.

## Product Features:

- 360 degrees of movement offering reduction of shear in all directions, essential when reducing the risk of skin breakdown related to shear forces
- Head and foot sections of foam are fixed in place to ensure that the component parts of the mattress 're-set' themselves
- Designed to reduce need for re-positioning
- Bolstered sides for ease of transfer
- Multi-zoned to address pressure redistribution requirements at head, torso and feet
- Higher cut middle section to improve immersion without compromising transfers
- 5 way stretch incontinent cover
- Two-year warranty

## Addressing *Shear* from the Inside Out



**The Align is available in a Raised Bolster version.** 3" are added to the head and foot section portions of the bolster allowing proper egress through the middle section. This minimizes the falls risk and meets/exceeds all restraint protocols.

Weight Capacity: 350lbs  
Standard Size: 36"x7"x80"  
Warranty: 2 years  
Custom Sizes Available

**5 MINDS MOBILITY GROUP**  
ggarland@5mindsmobility.com  
orders@5mindsmobility.com  
555 Washington Ave., Suite 310  
St. Louis, MO 63101

Ph: 866-709-1332 Fax: (877)-426-5697