

Hälsa Bed & Mattress Options
Retail Price List



*Permobil's ideal home care bed now available
with Optima and 5 Minds Mobility mattresses!*



Hälsa Bed features:

- exceptionally low positioning option
- the Trendelenburg capability comes standard
- upgraded large caster option available
- easy disassembly for transportation and storage

Hälsa Bed is available in regular and plus sizes, and with Maple (Nature) or Mahogany (Wine) coloured head and foot boards.



The Hälsa Bed from Permobil is an electric, adjustable bed that pairs well with the durable foam Optima mattress or a variety of 5 Minds Mobility mattresses (Adapt, Allay, or Align Bariatric).

Hälsa Bed has unique qualities that separate it from other electrical beds on the market. Most home health care beds have Trendelenburg capabilities as an option, but with the Hälsa Bed, the Trendelenburg feature comes standard.

Additionally, most beds don't offer the capability of moving both the head and foot sections simultaneously - the Hälsa Bed does.



Table of Contents

Hälsa Bed Comparison 3

Hälsa Bed Packages 3

Hälsa Bed and Components..... 4

Optima Mattress..... 4

5MM Mattresses 5

Hälsa Replacement Parts 6

Warranty 6

Discounts and Freight..... 6

Ordering Info 6

Hälsa Bed Comparison

Hälsa Bed

- » Dual posts on each side increases stability at all heights
- » Removable composite deck sections for easy cleaning
- » **Weight Capacity: 375 lb**
- » **Dimensions: 36" x 80"**
- » Low Position: 9.875" (12.125" with large caster)
- » High Position: 25.625" (27.875" with large caster)
- » Trendelenburg/Reverse Trendelenburg
- » Rolls at any height
- » Emergency battery back-up
- » Virtually noise free with low electrical consumption DC motor
- » Cover for headboard and footboard
- » Hooks for incontinence bag
- » 3 locking casters & 1 directional caster
- » Side and foot mattress retainers
- » Receptacle for trapeze bar
- » 2 receptacles for IV holder stems
- » Disassembles into 4 sections for transport/storage
- » Ideal for home care patients

Hälsa Plus Bariatric Bed

- » Dual posts on each side increases stability at all heights
- » **Weight Capacity: 500 lb**
- » **Width expansion: 4 options available**
- » **Dimensions: 36", 39", 42", 48" x 80"**
- » Low Position: 9.875" (12.125" with large caster)
- » High Position: 25.625" (27.875" with large caster)
- » Trendelenburg/Reverse Trendelenburg
- » Rolls at any height
- » Emergency battery back-up
- » Virtually noise free with low electrical consumption DC motor
- » Cover for headboard and footboard
- » Hooks for incontinence bag
- » 3 locking casters & 1 directional caster
- » Side and foot mattress retainers
- » Receptacle for trapeze bar
- » 2 receptacles for IV holder stems
- » Disassembles into 4 sections for transport/storage
- » Ideal for home care patients

Hälsa Bed Packages

The Hälsa Bed HSR package includes: head and foot board with cover, casters, motor, optional battery and a set of 1/2 rails.

The Hälsa Bed FSR package includes: head and foot board with cover, casters, motor, optional battery and a set of full size rails.

The Hälsa Plus Bariatric Bed HSR package includes: head and foot board with cover, casters, motor, optional battery and a set of 1/2 rails.



Part Number	Description	Retail Price
HASLA-STD-HSR-PKG	HALSA STANDARD BED w/ HALF SIDE RAILS	\$2,995.00
HASLA-STD-FSR-PKG	HALSA STANDARD BED w/ FULL SIDE RAILS	\$2,995.00
HASLA-PLUS-HSR-PKG	HALSA PLUS BARIATRIC BED w/ HALF SIDE RAILS	\$3,995.00
With Battery Backup		
HASLA-STD-HSR-PKG-BAT	HALSA STANDARD BED w/ HALF SIDE RAILS + BATTERY BACKUP	\$3,145.00
HASLA-STD-FSR-PKG-BAT	HALSA STANDARD BED w/ FULL SIDE RAILS + BATTERY BACKUP	\$3,145.00
HASLA-PLUS-HSR-PKG-BAT	HALSA PLUS BARIATRIC BED w/ HALF SIDE RAILS + BATTERY BACKUP	\$4,145.00

★ NEW fixed freight pricing for Hälsa beds & packages

# of beds	ON	QC	Vancouver, BC	Other cities, BC	AB	SK	MB	NS/NB/PEI	NL
3 or more beds (each)	\$59.00	\$59.00	\$99.00	\$129.00	\$89.00	\$89.00	\$79.00	\$79.00	\$179.00
1 or 2 beds (total cost)	\$118.00	\$118.00	\$198.00	\$258.00	\$178.00	\$178.00	\$158.00	\$158.00	\$358.00


Hälsa Bed and Components

Part Number	Description	Retail Price
HALSA-STD	HALSA STANDARD BED WITH HEAD/FOOT BOARDS	\$2,795.00
HALSA-PLUS	HALSA PLUS BARIATRIC BED WITH HEAD/FOOT BOARDS	\$3,795.00
HALSA-BAT	HALSA BED BATTERY BACKUP	\$150.00
HALSA-HSR	HALSA BED HALF SIDE RAILS (SET OF 2)	\$235.00
HALSA-FSR	HALSA BED FULL SIDE RAILS (SET OF 2)	\$280.00
HALSA-TB	HALSA BED TRAPEZE BAR STANDARD (150 LB WEIGHT CAPACITY)	\$200.00
HALSA-TBB	HALSA BED TRAPEZE BAR BARIATRIC (300 LB WEIGHT CAPACITY)	\$225.00
HALSA-CASTUPG-5"	HALSA CASTERS ASSEMBLY 5" (SET OF 4)	\$255.00

Optima Mattress


Optima mattresses are specifically engineered to provide lasting performance using better pressure reducing materials.

Comfort, safety, durability and budget are addressed by the anti-microbial, foam core constructed of 100% medical grade green friendly non-allergenic foam. Mattress is also flippable and reversible for longer product life. *Custom sizing is available.*



Part Number	Description	Retail Price
310	OPTIMA MATTRESS - 70D LT NYLON COVER, 36" x 80" x 5.5", 350 LB WEIGHT CAPACITY	\$460.00
314	OPTIMA MATTRESS - 70D LT NYLON COVER, 39" x 80" x 5.5", 350 LB WEIGHT CAPACITY	\$485.00
316	OPTIMA XT MATTRESS - 70D LT NYLON COVER, 42" x 80" x 5.5", 650 LB WEIGHT CAPACITY	\$625.00
318	OPTIMA XT MATTRESS - 70D LT NYLON COVER, 48" x 80" x 5.5", 650 LB WEIGHT CAPACITY	\$665.00

The **Optima Symmetry** premium mattress provides a multi-zoned foam support with: 2 layer, 1 zone reversible design, firm internal foam side rails, anti-shear liner and 70D Nylon fabric cover. *Custom sizing is available.*



Part Number	Description	Retail Price
311	OPTIMA SYMMETRY MATTRESS w/ FIRM SR - 70D LT NYLON COVER, 36" x 80" x 5.5", 350 LB WEIGHT CAPACITY	\$535.00
323	OPTIMA SYMMETRY MATTRESS w/ FIRM SR - 70D LT NYLON COVER, 39" x 80" x 5.5", 350 LB WEIGHT CAPACITY	\$560.00
325	OPTIMA SYMMETRY XT MATTRESS w/ FIRM SR - 70D LT NYLON COVER, 42" x 80" x 5.5", 650 LB WEIGHT CAPACITY	\$660.00
327	OPTIMA SYMMETRY XT MATTRESS w/ FIRM SR - 70D LT NYLON COVER, 48" x 80" x 5.5", 650 LB WEIGHT CAPACITY	\$750.00

* Additional 5% discount on Optima mattresses only applicable when purchased with a Hälsa bed.

Optima Key Terms

70D Nylon Cover - A durable, snug fitting low denier 70D Nylon is soft, moisture vapor permeable, and resistant to staining. It also includes a low slip vinyl bottom cover. Will not allow fluids to pass through and can be easily cleaned with standard cleaning products.

Non-Allergenic - The entire foam core is constructed of 100% medical grade green friendly non-allergenic foam.

Perimeter Foam (on the SYMMETRY MATTRESS) - Constructed of 2.0 lb density and IFD 86-94. Provides firm edge support and helps to keep the patient in a central position and from falling off the mattress.

Anti-Shear Liner (on the SYMMETRY MATTRESS) - Dissipates heat and reduces friction.


5MM Adapt, Allay and Align Bariatric Mattresses

The **Adapt** modifiable foam surface with SRT addresses shearing forces present during client movement in bed or when changing the position of the bed frame, at the same time ensuring proper pressure redistribution. This mattress comes complete with three removable GelFlex™ inserts, but allows for the addition of ROHO® air mattress sections. *Custom sizing is available.*



Part Number	Description	Retail Price
AD3680	ADAPT w/ SRT - 3 MOVABLE GELFLEX SECTIONS, 36" x 80" x 7", 350 LB WEIGHT CAPACITY	\$1,295.00
AD3980	ADAPT w/ SRT - 3 MOVABLE GELFLEX SECTIONS, 39" x 80" x 7", 350 LB WEIGHT CAPACITY	\$1,395.00
AD4280	ADAPT w/ SRT - 3 MOVABLE GELFLEX SECTIONS, 42" x 80" x 7", 350 LB WEIGHT CAPACITY	\$1,495.00
AD4880	ADAPT w/ SRT - 3 MOVABLE GELFLEX SECTIONS, 48" x 80" x 7", 350 LB WEIGHT CAPACITY	\$1,695.00

The **Allay** reactive GelFlex™ foam replacement surface with SRT is our premium option for the client wanting optimal pressure redistribution and comfort. The cooling benefit of the gel-infused visco foam top combined with a supportive core make this the best option for users not requiring the option of ROHO® mattress sections. *Custom sizing is available.*



Part Number	Description	Retail Price
AL3680	ALLAY w/ SRT - GELFLEX, 36" x 80" x 7", 350 LB WEIGHT CAPACITY	\$1,095.00
AL3980	ALLAY w/ SRT - GELFLEX, 39" x 80" x 7", 350 LB WEIGHT CAPACITY	\$1,195.00
AL4280	ALLAY w/ SRT - GELFLEX, 42" x 80" x 7", 350 LB WEIGHT CAPACITY	\$1,295.00
AL4880	ALLAY w/ SRT - GELFLEX, 48" x 80" x 7", 350 LB WEIGHT CAPACITY	\$1,495.00

The **Align Bariatric** reactive HR foam multi-zone surface with SRT provides differing immersion levels that match the contours of the body. The solid torso, with 4" square cuts under the sacrum and 2" square cuts at the foot, release the surface tension and ensure excellent envelopment. *Custom sizing is available.*



Part Number	Description	Retail Price
A3680B	ALIGN BARIATRIC HIGH DENSITY FOAM w/ SRT, 36" x 80" x 7", 600 LB WEIGHT CAPACITY	\$1,095.00
A3980B	ALIGN BARIATRIC HIGH DENSITY FOAM w/ SRT, 39" x 80" x 7", 600 LB WEIGHT CAPACITY	\$1,195.00
A4280B	ALIGN BARIATRIC HIGH DENSITY FOAM w/ SRT, 42" x 80" x 7", 600 LB WEIGHT CAPACITY	\$1,295.00
A4880B	ALIGN BARIATRIC HIGH DENSITY FOAM w/ SRT, 48" x 80" x 7", 600 LB WEIGHT CAPACITY	\$1,395.00



Note: Full 5MM Product Guide is available upon request.

* **Additional 5% discount on 5 Minds Mobility mattresses only applicable when purchased with a Hälsa bed.**

5 Minds Mobility Key Terms

Shear Reduction Technology (SRT) is an innovative solution designed to address shearing forces present during client movement in bed or when changing the position of the bed frame, at the same time ensuring proper pressure redistribution.

A third party clinical trial was completed by Durham College Department of Nursing and Department of Engineering to validate SRT. The trial provided non-biased support for the efficacy and use in populations where shear forces can damage vulnerable skin.



GelFlex™ foam provides enhanced support while reducing motion transfer across the mattress surface for an undisturbed sleep. Polymer gel beads are infused in the body of the visco foam. The open cell structure promotes airflow allowing for excellent cooling. GelFlex™ particles also:

- promote a cool, comfortable sleep
- are non-toxic and hypo-allergenic
- ensure better heat conductivity than traditional visco elastic foam
- provide superior support

Hälsa Replacement Parts

Part Number	Description	Retail Price
HALSA-HFBC-NATURE	HALSA BED 36" HEAD & FOOT BOARD - NATURE (SET OF 2)	\$295.00
HALSA-HFBC-WINE	HALSA BED 36" HEAD & FOOT BOARD - WINE (SET OF 2)	\$295.00
HALSA-MC	HALSA BED MOTOR CABLE	\$15.00
HALSA-PC	HALSA BED POWER CABLE	\$15.00
HALSA-FD-STD	HALSA BED STANDARD FOOT DECK	\$280.00
HALSA-HD-PANS	HALSA BED DECK PANS FOR HEAD DECK	\$190.00
HALSA-HD-STD	HALSA BED STANDARD HEAD DECK	\$190.00
HALSA-HS	HALSA BED HAND PENDANT 10 BUTTON	\$60.00
HALSA-FD-PLUS	HALSA BED PLUS FOOT DECK	\$280.00
HALSA-HD-PLUS	HALSA BED PLUS HEAD DECK	\$280.00
HALSA-MOTOR-BACK	HALSA BED HEAD ACTUATOR - STANDARD ONLY (C108-00)	\$115.00
HALSA-MOTOR-BACK-PLUS	HALSA BED HEAD ACTUATOR - PLUS ONLY (20040)	\$115.00
HALSA-MOTOR-HI	HALSA BED HI-LOW HEAD ACTUATOR FOR STANDARD AND PLUS (C440-00)	\$145.00
HALSA-MOTOR-KNEE	HALSA BED FOOT ACTUATOR FOR STANDARD AND PLUS (270019-00)	\$115.00
HALSA-MOTOR-LOW	HALSA BED HI-LOW ACTUATOR FOR STANDARD AND PLUS (C243-00)	\$145.00
HALSA-CB	HALSA BED CONTROL BOX (CB6050-00)	\$260.00
HALSA-CASTERS-4"	HALSA CASTERS ASSEMBLY 4" (SET OF 4) STANDARD	\$255.00
HALSA-CASTERS-5"	HALSA CASTERS ASSEMBLY 5" (SET OF 4)	\$255.00

Warranty

Bed Warranty: 10 years on bed frame, 2 years on electrical components

5 Minds Mobility Mattress Warranty: 2 year limited warranty

Optima Mattress Guarantee: 3 year total surface and foam cover guarantee

Optima Symmetry Guarantee: 5 year total surface and foam cover guarantee

Discounts and Freight

Additional 5% discount on Optima and 5MM mattresses only applicable when purchased with a Hälsa bed.

NEW fixed freight pricing for Hälsa beds & packages

# of beds	ON	QC	Vancouver, BC	Other cities, BC	AB	SK	MB	NS/NB/PEI	NL
3 or more beds (each)	\$59.00	\$59.00	\$99.00	\$129.00	\$89.00	\$89.00	\$79.00	\$79.00	\$179.00
1 or 2 beds (total cost)	\$118.00	\$118.00	\$198.00	\$258.00	\$178.00	\$178.00	\$158.00	\$158.00	\$358.00

To Place Your Order

Permobil Canada is the exclusive distributor for Hälsa beds, and Optima and 5 Minds Mobility's mattresses. To order, please contact us at sales.ca@permobil.com or 800-265-9830. To learn more, visit our website at www.permobil.ca.