

Pride® Elegance Collection Lift Chairs are designed to be the ultimate in lift chair comfort, style and performance. Available in several beautiful fabrics and styles to fit any décor while utilizing a state-of-the-art lift system.



Easy-to-use  
hand control with  
USB charging port

ENERGY  
EFFICIENT  
WITH BUILT-IN  
BATTERY BACKUP  
(BATTERIES REQUIRED)



Scan here for more  
information on the  
**Elegance Collection**



FOLLOW US ON SOCIAL MEDIA

**PRIDE**®  
MOBILITY  
LIVE YOUR BEST®

## FEATURES & BENEFITS

- Industry-first Total Comfort Seating features a unique spring design and one-piece soft foam seat for outstanding relaxation (not available on LC-108)
- Premium furniture grade laminate/hardwood frames for added strength and durability
- Easy-to-use hand control with USB charger port (USB charging port not available when selecting the optional heat and massage)
- Standard head and arm covers
- Large, dual pockets for extra storage
- Back has access for adjusting the fill for comfort and fit (Back cushion of the LC-450 and LC-108 is fully removable)
- Quiet and smooth lift system
- Footrest extension standard (LC-108 and LC-570L only)
- Integrated battery backup ensures chair is operational when power goes out (Batteries not included)

## OPTIONS\*

- Cleaning kit
- Head pillow
- Footrest extension
- Wireless remote (LC-450C only)

\* Additional charges apply for all options. Extra lead time may apply for some options.



MODEL	LC-108	LC-450C	LC-570S	LC-570M	LC-570L
POSITION	Infinite-position Trendelenburg Chaise Lounger 5'4" - 6'0"	3-Position Chaise Lounger 5'4" - 6'0"	3-Position Chaise Lounger 5'3" and below	3-Position Chaise Lounger 5'4" - 6'0"	3-Position Chaise Lounger 5'9" - 6'2"
SIZE	Medium	Medium	Small	Medium	Large
WEIGHT CAPACITY	375 lbs.	375 lbs.	375 lbs.	375 lbs.	375 lbs.
SEAT-TO-FLOOR	19"	20"	19"	20"	21"
SEAT DEPTH	19"	20"	19"	20"	22"
SEAT WIDTH	20"	20"	19"	20"	22"
TOP OF BACK TO SEAT	29"	27"	24.5"	26.5"	27"
FULLY PADDED CHAISE	Yes (Removable and adjustable)	Yes (Removable and adjustable)	Yes	Yes	Yes
SEATING TYPE	Sinuous Spring	Total Comfort	Total Comfort	Total Comfort	Total Comfort
DISTANCE FROM WALL	24"	18"	16.5"	19"	24"
OVERALL WIDTH	31.5"	37"	35"	36"	38"
CHAIR WEIGHT	136 lbs.	119 lbs.	110 lbs.	113 lbs.	114 lbs.
WARRANTY	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated



**CLOUD 9 FABRICS**  
Standard

**LEXIS STA-KLEEN**  
(LC-108 only)

Beamsville, ON / 888-570-1113 / pridemobility.ca

Note: Measurements may vary +/- 1". The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.



**PRIDE**<sup>®</sup>  
MOBILITY  
LIVE YOUR BEST<sup>®</sup>