

The Essential Collection from Pride® Power Lift Recliners provides exceptional value and outstanding comfort. Several fabric options will match any décor while delivering rest and worry-free operation.



Easy-to-use hand control with USB charging port

ENERGY EFFICIENT WITH BUILT-IN BATTERY BACKUP (BATTERIES REQUIRED)



Scan here for more information on the **Essential Collection®**



FOLLOW US ON SOCIAL MEDIA

## FEATURES & BENEFITS

- Easy-to-use hand control adjusts your chair to find your perfect position
- Large dual pockets offer ample storage space that is easily accessible
- Smooth and durable lift and scissor mechanism for long-lasting and quiet performance
- Back has access for adjusting the fill for comfort and fit
- Integrated battery backup ensures chair is operational when the power goes out (Batteries not included)
- Energy-efficient power box equipped with self-testing electronics and 9' cord
- Lifetime warranty on steel lift frame
- Footrest extension standard (LC-107 only)

## OPTIONS\*

- Cleaning kit
- Non-skid leg levelers

\* Additional charges apply for all options. Extra lead time may apply for some options.



MODEL	LC-102	LC-106	LC-107
POSITION	3-Position Chaise Lounger 4'8" - 5'4"	3-Position Chaise Lounger 5'2" - 5'10"	Infinite Position Trendelenburg Chaise Lounger 5'4" - 6'0"
SIZE	Super Petite	Medium	Medium
WEIGHT CAPACITY	250 lbs.	325 lbs.	325 lbs.
SEAT-TO-FLOOR	16"	18"	20"
SEAT DEPTH	18"	19"	20"
SEAT WIDTH	16"	20"	20"
TOP OF BACK TO SEAT	23.5"	27.5"	27.5"
FULLY PADDED CHAISE	Yes	Yes	Yes
SEATING TYPE	Sinuous Spring	Sinuous Spring	Sinuous Spring
DISTANCE FROM WALL	16"	18"	26"
OVERALL WIDTH	25.5"	32.75"	35"
CHAIR WEIGHT	96 lbs.	96 lbs.	116 lbs.
WARRANTY	1-Year Limited	1-Year Limited	1-Year Limited



## CLOUD 9 FABRICS

Beamsville, ON / 888-570-1113 / pridemobility.ca

Note: Measurements may vary +/- 1". The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.



**PRIDE**<sup>®</sup>  
MOBILITY  
LIVE YOUR BEST<sup>®</sup>