



Pride Mobility Products Co.

5096 South Service Rd, Beamsville, ON, L3J 1V4

Phone: 888-570-1113 | Fax: 866-514-1303 | Email: productorders@pridemobility.com

# Q6 Edge® HD Order Form with TRU-Balance® 3 Positioning

Account Number: \_\_\_\_\_ Date: \_\_\_\_\_

Provider Name: \_\_\_\_\_

Contact: \_\_\_\_\_

Phone: \_\_\_\_\_ Fax: \_\_\_\_\_

Email: \_\_\_\_\_

PO Number: \_\_\_\_\_

Marked for: \_\_\_\_\_

Ship to Address: \_\_\_\_\_

City: \_\_\_\_\_ Prov: \_\_\_\_\_ PC: \_\_\_\_\_

## Q6 edge<sup>®</sup> HD



### INTRODUCTION

This form is interactive when viewed with [Adobe Acrobat Reader](#). Complete the form by placing checks in the desired boxes and provide information in the interactive fields. Buttons shown at the bottom of the form may be utilized to print or submit the order form through a desktop email application. To email via a web-based application, please 'Save As' and attach the PDF to your email.

This order form contains a large variety of options to fulfill various patient needs. Descriptions and section notes, such as optional or required, are included to help you complete your order. Other part numbers can be found on the quote acknowledgment after this order form is submitted. Please contact Quantum® Sales at 888-570-1113 if assistance is needed. Send the completed order form by fax (866-514-1303) or email ([productorders@pridemobility.com](mailto:productorders@pridemobility.com)). Incomplete forms may delay the quote or order. Customer service will contact you if the order is incomplete or if there are compatibility issues. If special order requests are needed, be sure that the Patient Information section is completed

All prices are MSRP. Prices, specifications, part numbers, and availability are subject to change without notice. Prices and part numbers as shown are only available when configured as a complete power chair. Please contact technical service for accurate parts ordering.

Thank you for choosing Quantum!


### PATIENT INFORMATION

Weight: \_\_\_\_\_ Height: \_\_\_\_\_

A. Shoulder Width: \_\_\_\_\_ D. Max Sitting Height: \_\_\_\_\_ G. Thigh Depth: \_\_\_\_\_

B. Chest Width: \_\_\_\_\_ E. Shoulder Height: \_\_\_\_\_ H. Lower Leg Length: \_\_\_\_\_

C. Hip Width: \_\_\_\_\_ F. Axilla Height: \_\_\_\_\_ I. Elbow Height: \_\_\_\_\_













### 1. BASE MODEL

**REQUIRED.** Select a group 3 power base. 450 lb. weight capacity with standard 4.7 mph programmable motors, 8 amp off board charger, and unoccupied transit loops. Dimensions of base without legrests: 26"Wx36.75"L.

Quantum Q6 Edge HD ..... \$8,815  
Part: Q6 EDGEHD

**2. BASE COLOR**

**REQUIRED.** Select one shroud color. Colors shown here may differ from actual product depending on monitor/printer calibration used.

- |  |   |   |   |
|--|---|---|---|
| <input type="checkbox"/>  Red Corvette<br>Part: DGN169178-PRD1  | <input type="checkbox"/>  Night Train<br>Part: DGN169178-PBK3      | <input type="checkbox"/>  Deep Purple<br>Part: DGN169178-PPR5    | <input type="checkbox"/>  Green Machine<br>Part: DGN169178-PGN11 |
| <input type="checkbox"/>  Electric Blue<br>Part: DGN169178-PBU9 | <input type="checkbox"/>  Pretty in Pink<br>Part: DGN169178-PPK3   | <input type="checkbox"/>  Grey Street<br>Part: DGN169178-PSL3    |   |
| <input type="checkbox"/>  Back in Black<br>Part: DGN169178-PBK1 | <input type="checkbox"/>  Evolution Orange<br>Part: DGN169178-POR7 | <input type="checkbox"/>  Midnight Blue<br>Part: DGN169178-PBU10 |   |

**3. FENDER LIGHTS**

**Optional.** Fender lights operate by pressing included single switch button. Fender lighting is not available with NE 4-Key Joystick or Single Actuator Control through Toggle Box.

- Switch Fender Lights without Lift System..... \$240  
Part: XRF000177
- Switch Fender Lights with Lift System ..... No Charge  
Part: DGN154117

**4. QUANTUM BACKUP CAMERA**

**Optional.** Requires Q-Logic electronics. The backup camera may interfere with armrest mounted accessories such as media device holder.

- Quantum Backup Camera ..... \$475  
Part: MEC1910365



**1. Optional Camera Mounting Bracket**

- Backup Camera Headrest Bracket..... \$30  
Part: MEC1910373  
Works with .5" x .5" vertical supports only.
- Backup Camera Track Mount Bracket ..... \$35  
Part: MEC2104227

**2. Optional Display Mounting Bracket**

- Backup Camera Q-Logic Joystick Mount Bracket ....\$145  
Part: MEC2106123
- Backup Camera NE/NE+ Joystick Mount Bracket ....\$145  
Part: MEC2106132

**5. DRIVE WHEELS**

**REQUIRED.** Select a type of drive wheels and accent color. Knobby tires increase overall base width by 1.5".

**5.1. Drive Wheels**

- |  |  |
|--|--|
| <input type="checkbox"/> 3" x 14" Run Flat Tires..... Standard<br>Part: PTOINDV3510    | <input type="checkbox"/> 4" x 14" Flat Free, Knobby, Extended Traction Tires..... \$120<br>Part: WHL121698 x2  |
| <input type="checkbox"/> 3" x 14" Flat Free Tires..... No Charge<br>Part: WHL124116 x2 | <input type="checkbox"/> 4" x 14" Pneumatic, Knobby, Extended Traction Tires ..... \$120<br>Part: WHL121703 x2 |
| <input type="checkbox"/> 3" x 14" Pneumatic Tires..... No Charge<br>Part: WHL124117 x2 |  |

**5.2. Drive Wheel Accent Color**

**Optional. Not available with knobby tires.**

- |  |  |  |
|--|--|--|
| <input type="checkbox"/> Back in Black Accent Color..... \$42.50 ea.<br>Part: XRF000179  | <input type="checkbox"/> Pretty in Pink Accent Color..... \$42.50 ea.<br>Part: XRF000185 | <input type="checkbox"/> Midnight Blue Accent Color ..... \$42.50 ea.<br>Part: XRF000187 |
| <input type="checkbox"/> Red Corvette Accent Color ..... \$42.50 ea.<br>Part: XRF000180  | <input type="checkbox"/> Evolution Orange Accent Color .. \$42.50 ea.<br>Part: XRF000186 | <input type="checkbox"/> Green Machine Accent Color ..... \$42.50 ea.<br>Part: XRF000190 |
| <input type="checkbox"/> Electric Blue Accent Color ..... \$42.50 ea.<br>Part: XRF000188 | <input type="checkbox"/> Deep Purple Accent Color..... \$42.50 ea.<br>Part: XRF000184    |  |
| <input type="checkbox"/> Night Train Accent Color ..... \$42.50 ea.<br>Part: XRF000181   | <input type="checkbox"/> Grey Street Accent Color..... \$42.50 ea.<br>Part: XRF000182    |  |

**6. BATTERIES**

**Optional.** The powerchair requires 2 batteries and will be installed if ordered.

- Group 24, 70 Amp Gel Cell, Interceptor 2 Battery ..... \$400 ea.  
Part: BAT2100315. HCPCS: E2342.50 ea.
- Group 24, 70 Amp Gel Cell, MK Battery .....\$510 ea.  
Part: BATGEL1009. HCPCS: E2342.50 ea.

**7. ELECTRONICS**

**REQUIRED.** Select a type of drive control from the following sections. The No Joystick option is used to omit the joystick and utilize another control, such as a specialty control at the end of this form or found on the [Alternative Controls and Electronics Order Forms](#).

**7.1a. Non-Expandable Controls**

Select a joystick below or omit selection if another controller type is desired. Non-expandable does not allow supplementary drive controllers.

- NE 4-Key Joystick.....Standard  
Part: CTLDC1554  
Available with Static Seat or Tilt Only through Toggle.
- NE+ 6-Key Joystick ..... No Charge  
Part: CTLDC1560  
Control of one actuator through Joystick on single power base model (3SPHD-SS) only.
- No Joystick, Setup for NE+ Joystick ..... No Charge  
Part: XRF000009



Electronics section continued on following page.

**7. ELECTRONICS CONTINUED**

**7.1b. Expandable Controls - Quantum Intelligent Systems™**

Select expandable controls or omit selection if another controller type is desired. Expandable controls allow multiple or specialized controllers.

**2. Q-Logic 3 Joystick**

- Q-Logic 3 EX Joystick .....\$900  
Part: CTL167724
- Q-Logic 3e Joystick .....\$620  
Part: CTL1910600  
Can operate full light package. Not available with alternative or specialty controls.
- Q-Logic 3 EX Joystick for Full Lights .....\$900  
Part: CTL167725
- No Joystick, Setup for Q-Logic  
Part: XRF000010  
Requires selection of Enhanced Display.



**7.2. Q-Logic 3 Joystick Armor**

**REQUIRED with Q-Logic 3 Joystick selections above.** Reference colors shown on page 2.

- Black, Soft Touch  
Part: DGN124307
- Red Corvette  
Part: DGN124336-PRD1
- Electric Blue  
Part: DGN124336-PBU9
- Back in Black  
Part: DGN124336-PBK1
- Night Train  
Part: DGN124336-PBK3
- Pretty in Pink  
Part: DGN124336-PPK3
- Evolution Orange  
Part: DGN124336-POR7
- Deep Purple  
Part: DGN124336-PPR5
- Grey Street  
Part: DGN124336-PSL3
- Midnight Blue  
Part: DGN124336-PBU10
- Green Machine  
Part: DGN124336-PGN11

**8. JOYSTICK HANDLES**

**Optional.** The NE/NE Plus joysticks are standard with a carrot handle (See A in image). The standard Q-Logic 3 joystick knob can be seen in the image in section 7.1b. The optional handle will be mounted to the joystick and the standard handle will be shipped with the unit.

- Mushroom Handle (B) ..... \$120  
Part: ACCASMB1860.
- Chin Cup (C) ..... \$85  
Part: ACCASMB186190
- T-Bar (D) ..... \$160  
Part: ACCASMB1858.
- Medium Ball Handle (E) ..... \$70  
Part: ACCASMB1859.
- Soft Ball Handle (F) ..... \$120  
Part: ACCASMB1856.
- Extended I-Handle (G) ..... \$140  
Part: ACCASMB1857.
- Bodypoint 3" Goal Post Handle ..... \$180  
Part: PLSKNCP1057.
- Bodypoint 4" Goal Post Handle ..... \$180  
Part: PLSKNCP1058.




**9. JOYSTICK MOUNTING BRACKETS**

**Optional.** Select a mount and side to be used with a joystick selected in previous section.

**9.1a. Fixed, In-line Joystick Bracket**


- Fixed, In-line, Bracket for Q-Logic 3 .....Standard  
Part: FRMASMB10617
- Fixed, In-line, Bracket for NE/NE+ .....Standard  
Part: FRMASMB12744



**9.1b. Swing-Away (Retractable) Joystick Bracket**

Allows the joystick to retract to the side of the armrest while maintaining the same joystick orientation.

- Swing-away, Q-Logic 3 .....\$265  
Part: ACC1908434. HCPCS: E1028
- Swing-away, NE/NE+ .....\$265  
Part: ACC1908437. HCPCS: E1028



**9.1c. Basic Swing Away Joystick Mount**

- Basic Swing Away, NE/NE+, Left Side .....\$265  
Part: ACC2201047. HCPCS: E1028
- Basic Swing Away, NE/NE+, Right Side .....\$265  
Part: ACC2201048. HCPCS: E1028
- Basic Swing Away, Q-Logic, Left Side .....\$265  
Part: ACC2201041. HCPCS: E1028
- Basic Swing Away, Q-Logic, Right Side .....\$265  
Part: ACC2201042. HCPCS: E1028

**9.1d. Offset Swing Away Joystick Mount**

- Offset Swing Away, NE/NE+, Left Side .....\$265  
Part: ACC142747. HCPCS: E1028
- Offset Swing Away, NE/NE+, Right Side .....\$265  
Part: ACC142748. HCPCS: E1028
- Offset Swing Away, Q-Logic, Left Side .....\$265  
Part: ACC142742. HCPCS: E1028
- Offset Swing Away, Q-Logic, Right Side .....\$265  
Part: ACC142743. HCPCS: E1028

**9.1e. Specialty Joystick Mounts**

- Retract4 Mount .....\$368  
Part: XRF000003. HCPCS: E1028  
Full 180° of movement that translates into multi-position use, straight, right offset, left offset, or fold away. Joystick orientation varies as bracket is adjusted.
- Stealth Swing-Away/Flip Down Bracket .....\$242.50 ea.  
Part: XRF000005. HCPCS: E1028  
Allows the joystick to pivot to side of armrest, or flip-down if bracket is adjusted.
- Multi-Axis Mount ..... \$126  
Part: XRF000004  
Allows height adjustment of joystick.



**9.2. Joystick Mount Location**

- Left Side Mount  
Part: PTOINDV1225
- Right Side Mount  
Part: PTOINDV1226

**10. TRU-BALANCE® 3 SEATING - QUANTUM INTELLIGENT SYSTEMS™**

**REQUIRED.** Maximum seat size is 32"W x 25"D.

**10.1. TRU-Balance 3 Configuration**

- TRU-Balance 3 Static Seat .....\$1,045  
Part: PTOINDV3442
- TRU-Balance 3 Heavy Duty, Power Recline ..... \$7,875  
Part: PTOINDV3645. HCPCS: E1004  
175° of power recline.
- TRU-Balance 3, Heavy Duty, Power Tilt .....\$6,245  
Part: PTOINDV3643. HCPCS: E1002  
50° Center of Gravity (CG) power tilt seating system.
- TRU-Balance 3 Heavy Duty, Power Tilt & Recline .....\$13,125  
Part: PTOINDV3644. HCPCS: E1007  
50° CG power tilt with 175° of power recline. Full drive lockout at 20°.

**10.2. 10" Power Lift**

**Optional.** ^Please verify the appropriate billing code to be used for this item with all third-party payors.

- TRU-Balance 3, Heavy Duty, 10" Power Lift with iLevel® and Fender Lights..... \$3,495  
Part: PTOINDV3797. HCPCS: E2298/E2300  
400 lb weight capacity. 24"W x 24"D maximum seat size. Includes fender lights found in section 5.



**10.3. Seat Width**

**REQUIRED.** Maximum non-lift seat width is 32"W. \*Widths over 24" are not available with 10" power lift. Part number is SETWDTH plus number following # below.

- 20" W (#1010)
- 22" W (#1012)
- 24" W (#1014)
- 26" W\* (#1054)
- 28" W\* (#1055)
- 30" W\* (#1059)
- 32" W\* (#1061)
- 21" W (#1011)
- 23" W (#1013)
- 25" W\* (#1056)
- 27" W\* (#1057)
- 29" W\* (#1058)
- 31" W\* (#1060)

**10.4. Seat Width Package**

**REQUIRED.**

- 20"-24" Width Package  
Part: XRF000158
- 24"-28" Width Package  
Part: XRF000192
- 28"-32" Width Package  
Part: XRF000193

**10.5. Seat Depth**

**REQUIRED.** Maximum seat depth is 25". Part number is SETDPTH plus number following # below.

- 20" D (#1010)
- 21" D (#1011)
- 22" D (#1012)
- 23" D (#1013)
- 24" D (#1014)
- 25" D (#1060)

**10.6. Seat Depth Package**

**REQUIRED.**

- 16"-20" Depth Package  
Part: XRF000013  
Minimum Seat Depth is 20"
- 20"-24" Depth Package  
Part: XRF000159  
Use with 25" Depth

Low	Med	High
Static Seat, Recline		
16.5"	17.5"	18.5"
Tilt, Tilt & Recline, Lift		
18.5"	19.5"	20.5"
TRU-Balance 3 Seat to Floor Heights		

**10.7. Seat to Floor Height**

**REQUIRED.**

- Lowest Seat to Floor  
Part: PTOINDV3558
- Middle Seat to Floor  
Part: PTOINDV3559
- Highest Seat to Floor  
Part: PTOINDV3560

**10.8. Matching Seat Cushion**

**Optional.** Specific cushion/component sizes available, part numbers, weight capacities, MSRP, HCPCS code, and special order cushions can be found on the respective complete [Stealth Cushion Order Forms](#). This section allows a quick selection of a single cushion that matches the seat pan width and depth. See specific order form for non-matching sizes, alternate covers and liners, or additional covers and liners.

**Stealth TRU-Comfort 2 SPP Cushion**

Skin protection and positioning with Coolcore, 4-way stretch, breathable cover.  
Part series: ST-TC2SPPXXYY. HCPCS: E2607/E2608


- TRU-Comfort 2 SPP, Matching Seat Size  
Reference part: XRF000116.



**Stealth Solution SPP Cushion**

Skin protection and positioning with Coolcore, 4-way stretch, breathable cover.  
Includes urethane liner.  
Part series: ST-SOLSPPXXYY. HCPCS: E2607/E2608


- Solution SPP, Matching Seat Size  
Reference part: XRF000113.



**Stealth Spectrum Foam SPP Cushion**

Skin protection and positioning with Coolcore, 4-way stretch, breathable cover.  
Includes urethane liner.  
Part series: ST-SPCFSPXXYY. HCPCS: E2607/E2608

- Spectrum Foam SPP, Matching Seat Size  
Reference part: XRF000112.



**Stealth Premiere P Cushion**

Positioning with Coolcore, 4-way stretch, breathable cover.  
Part series: ST-PREMPXXYY. HCPCS: E2605/E2606

- Premiere P, Matching Seat Size  
Reference part: XRF000122.



**Stealth Spectrum Gel SPP Cushion**

Skin protection and positioning with Coolcore, 4-way stretch, breathable cover.  
Includes urethane liner.  
Part series: ST-SPCGSPPXXYY. HCPCS: E2607/E2608

- Spectrum Gel SPP, Matching Seat Size  
Reference part: XRF000117.



**Stealth Essence SPP Cushion**

Skin protection and positioning with Coolcore, 4-way stretch, breathable cover.  
Part series: ST-ESSSPPXXYY. HCPCS: E2607/E2608


- Essence SPP, No Wedge, Matching Seat Size  
Reference part: XRF000118.
- Essence SPP, Small Wedge, Matching Seat Size  
Reference part: XRF000119.
- Essence SPP, Medium Wedge, Matching Seat Size  
Reference part: XRF000120.
- Essence SPP, Large Wedge, Matching Seat Size  
Reference part: XRF000121.



**Stealth Glacial SP Cushion**

Skin protection with Coolcore, 4-way stretch, breathable cover.  
Part series: ST-GLSPXXYY. HCPCS: E2603/E2604

- Glacial SP, Matching Seat Size  
Reference part: XRF000114.




TRU-Balance 3 Seating section continued on following pages.

**10. JOYSTICK MOUNTING BRACKETS CONTINUED**

**Stealth Zen SP Cushion**

Skin protection with Coolcore, 4-way stretch, breathable cover.  
Part series: ST-ZENSPXXYY. HCPCS: E2603/E2604

- Zen SP, Matching Seat Size  
Reference part: XRF000115.



**Stealth Simplicity G Cushion**

General use cushion with Stealth Tek 2, 2-way stretch, breathable cover.  
Part series: ST-SIMGXXYY. HCPCS: E2601/E2602

- Simplicity G, Matching Seat Size  
Reference part: XRF000111.



**10.9. Back Type Selection**

**REQUIRED.** Match choice in this section with the options in the following back cushion configuration section.

- TRU-Balance 3 Back Assembly (Static Back) ..... \$735  
Part: XRF000194  
Select a back cushion configuration in the subsequent section.
- Flip Forward Back..... \$735  
Part: XRF000220  
Available with Static Seat and Power Tilt (1 actuator). Only available with 2-post, flip back or pediatric armrests. Not available with 12"-14" width. Select a back cushion type in the subsequent section.
- Back Cane Option..... No Charge  
Part: XRF000015  
Available only with TRU-Balance 3 Static Seat and Power Tilt selections. Select back cane type, back cane height, and back cushion configuration.

**10.10. Back Cane Selection (Back Cane Option Only)**

**REQUIRED** when Back Cane Option is selected above. Select type and height of back canes.

**1. Back Cane Type**

- Angle Adjustable  
Part: MEC140266
- 1" Endo  
Part: XRF000269
- 2" Endo  
Part: XRF000271

**2. Back Cane Height**

- 16" High Cane  
Part: FRMASMB15834
- 17" High Cane  
Part: FRMASMB15835
- 20" High Cane  
Part: FRMASMB15838
- 21" High Cane  
Part: FRMASMB15839
- 24" High Cane  
Part: FRMASMB15842
- 25" High Cane  
Part: FRMASMB15843

**10.11. Back Cushion Configuration**

**REQUIRED.** Select a configuration type and applicable dimensions. See the Sport, Q-Back, and TRU-Comfort back availability chart for stock order back size ranges. Other size backs may be available with special order and extended lead time.

**1a. Sport Back/Omit/After-Market Interface**

- Sport Back..... No Charge  
Part: XRF000017  
Complete following width package and size.
- After-market Back Interface ..... \$210  
Part: XRF000018  
Select a back width within a range of 16"-24"W. Overall seat depth may need to be adjusted to accommodate certain backs. May interfere with gel/waterfall armpads and arm through. Does not include brackets to mount after-market backs.
- Omit Back Pan/Cushion..... No Charge  
Part: PTOINDV1287  
Select a following width package.

**1b. TRU-Comfort Back**

- TRU-Comfort Back, Back Cane Option ..... \$835  
Part: TC4V10W21WBC. HCPCS: E2620  
Complete following width package and size.
- TRU-Comfort Back (Static/Flip Forward), 22" Width or Greater.. \$1,235  
Part: TC4V22W32W. HCPCS: E2621  
Complete following width package and size.
- TRU-Comfort Back, Back Cane Option, 16" to 20" Width..... \$1,035  
Part: TC4V10W21WBC. HCPCS: E2620  
Complete following width package and size.

**1c. Q-Back 2.0, AL Series, Active Contour**

**Not available with Back Cane Option.** Select a back type and the specific overall back height in subsequent section.

- Active Contour, 18"Wx19-20"H Tilt/22"-23" H Recline ..... \$1,145  
Part: QB2AL18AS. HCPCS: E2613
- Active Contour, 20"Wx19-20"H Tilt/22"-23" H Recline..... \$1,145  
Part: QB2AL20AS. HCPCS: E2613
- Active Contour, 18"Wx21-22"H Tilt/24"-25" H Recline ..... \$1,145  
Part: QB2AL18AT. HCPCS: E2613
- Active Contour, 20"Wx21-22"H Tilt/24"-25" H Recline..... \$1,145  
Part: QB2AL20AT. HCPCS: E2613
- Active Contour, 18"Wx23-24"H Tilt/26"-27" H Recline ..... \$1,145  
Part: QB2AL18AX. HCPCS: E2613
- Active Contour, 20"Wx23-24"H Tilt/26"-27" H Recline..... \$1,145  
Part: QB2AL20AX. HCPCS: E2613
- Active Contour, 18"Wx25-26"H Tilt/28"-29" H Recline..... \$1,145  
Part: QB2AL18AXX. HCPCS: E2613
- Active Contour, 20"Wx25-26"H Tilt/28"-29" H Recline..... \$1,145  
Part: QB2AL20AXX. HCPCS: E2613

**1d. Q-Back 2.0, AL Series, Deep Contour**

**Not available with Back Cane Option.** Select a back type and the specific overall back height in subsequent section.

- Deep Contour, 18"Wx19-20"H Tilt/22"-23" H Recline ..... \$1,325  
Part: QB2AL18DS. HCPCS: E2620
- Deep Contour, 20"Wx19-20"H Tilt/22"-23" H Recline ..... \$1,325  
Part: QB2AL20DS. HCPCS: E2620
- Deep Contour, 18"Wx21-22"H Tilt/24"-25" H Recline ..... \$1,325  
Part: QB2AL18DT. HCPCS: E2620
- Deep Contour, 20"Wx21-22"H Tilt/24"-25" H Recline ..... \$1,325  
Part: QB2AL20DT. HCPCS: E2620
- Deep Contour, 18"Wx23-24"H Tilt/26"-27" H Recline ..... \$1,325  
Part: QB2AL18DX. HCPCS: E2620
- Deep Contour, 20"Wx23-24"H Tilt/26"-27" H Recline ..... \$1,325  
Part: QB2AL20DX. HCPCS: E2620
- Deep Contour, 18"Wx25-26"H Tilt/28"-29" H Recline ..... \$1,325  
Part: QB2AL18DXX. HCPCS: E2620
- Deep Contour, 20"Wx25-26"H Tilt/28"-29" H Recline ..... \$1,325  
Part: QB2AL20DXX. HCPCS: E2620

TRU-Balance 3 Seating section continued on following page.

**10. TRU-BALANCE 3 SEATING CONTINUED**

**1e. Q-Back 2.0, AL Series, Extra Deep Contour**

**Not available with Back Cane Option.** Select a back type and the specific overall back height in subsequent section.

- |   |   |
|---|---|
| <input type="checkbox"/> Extra Deep Contour, 18"Wx19-20"H Tilt/22"-23" H Recline .... \$1,410<br>Part: QB2AL18ES. HCPCS: E2620  | <input type="checkbox"/> Extra Deep Contour, 20"Wx19-20"H Tilt/22"-23" H Recline .... \$1,410<br>Part: QB2AL20ES. HCPCS: E2620  |
| <input type="checkbox"/> Extra Deep Contour, 18"Wx21-22"H Tilt/24"-25" H Recline.... \$1,410<br>Part: QB2AL18ET. HCPCS: E2620   | <input type="checkbox"/> Extra Deep Contour, 20"Wx21-22"H Tilt/24"-25" H Recline .... \$1,410<br>Part: QB2AL20ET. HCPCS: E2620  |
| <input type="checkbox"/> Extra Deep Contour, 18"Wx23-24"H Tilt/26"-27" H Recline .... \$1,410<br>Part: QB2AL18EX. HCPCS: E2620  | <input type="checkbox"/> Extra Deep Contour, 20"Wx23-24"H Tilt/26"-27" H Recline .... \$1,410<br>Part: QB2AL20EX. HCPCS: E2620  |
| <input type="checkbox"/> Extra Deep Contour, 18"Wx25-26"H Tilt/28"-29" H Recline .... \$1,410<br>Part: QB2AL18EXX. HCPCS: E2620 | <input type="checkbox"/> Extra Deep Contour, 20"Wx25-26"H Tilt/28"-29" H Recline .... \$1,410<br>Part: QB2AL20EXX. HCPCS: E2620 |

**1f. Q-Back 2.0, AFT Series, Active Contour**

**Not available with Back Cane Option.** Select a back type and the specific overall back height in subsequent section.

- |   |   |
|---|---|
| <input type="checkbox"/> AFT, Active, 16"Wx19-20"H Tilt/22"-23" H Recline ..... \$1,145<br>Part: QB2AFT16AS. HCPCS: E2613 | <input type="checkbox"/> AFT, Active, 18"Wx19-20"H Tilt/22"-23" H Recline ..... \$1,145<br>Part: QB2AFT18AS. HCPCS: E2613 |
|---|---|

**1g. Q-Back 2.0, AFT Series, Deep Contour**

**Not available with Back Cane Option.** Select a back type and the specific overall back height in subsequent section.

- |  |  |
|--|--|
| <input type="checkbox"/> AFT, Deep Contour, 16"Wx19-20"H Tilt/22"-23" H Recline..... \$1,325<br>Part: QB2AFT16DS. HCPCS: E2620 | <input type="checkbox"/> AFT, Deep Contour, 18"Wx21-22"H Tilt/24"-25" H Recline..... \$1,325<br>Part: QB2AFT18DT. HCPCS: E2620 |
| <input type="checkbox"/> AFT, Deep Contour, 16"Wx21-22"H Tilt/24"-25" H Recline..... \$1,325<br>Part: QB2AFT16DT. HCPCS: E2620 | <input type="checkbox"/> AFT, Deep Contour, 20"Wx21-22"H Tilt/24"-25" H Recline .... \$1,325<br>Part: QB2AFT20DT. HCPCS: E2620 |
| <input type="checkbox"/> AFT, Deep Contour, 18"Wx19-20"H Tilt/22"-23" H Recline .... \$1,325<br>Part: QB2AFT18DS. HCPCS: E2620 |  |

**1h. Q-Back 2.0, AFT Series, Extra Deep Contour**

**Not available with Back Cane Option.** Select a back type and the specific overall back height in subsequent section.

- |   |
|---|
| <input type="checkbox"/> AFT, Extra Deep Contour, 16"Wx25-26"H Tilt/28"-29" H Recline..... \$1,410<br>Part: QB2AFT16EXX. HCPCS: E2620 |
|---|

**10.12. TRU-Balance 3 Back Shroud**

**REQUIRED.** The back shroud is a feature of the TRU-Balance 3 seating system and is available in matching shroud colors. The back shroud is not available when Omit Back Pan/Cushion is selected. Reference colors shown in section 2.

- |   |   |  |  |   |
|---|---|--|--|---|
| <input type="checkbox"/> Red Corvette<br>Part: XRF000019  | <input type="checkbox"/> Night Train<br>Part: XRF000022 | <input type="checkbox"/> Deep Purple<br>Part: XRF000025    | <input type="checkbox"/> Evolution Orange<br>Part: XRF000028 | <input type="checkbox"/> Electric Blue<br>Part: XRF000030 |
| <input type="checkbox"/> Back in Black<br>Part: XRF000021 | <input type="checkbox"/> Grey Street<br>Part: XRF000023 | <input type="checkbox"/> Pretty in Pink<br>Part: XRF000027 | <input type="checkbox"/> Midnight Blue<br>Part: XRF000029    | <input type="checkbox"/> Green Machine<br>Part: XRF000033 |

**10.13. Lateral Positioning Wedge**

**Optional.**

- |  |
|--|
| <input type="checkbox"/> Lateral Positioning Wedge, pair ..... \$15 ea.<br>Part: SET123465<br>Available with Sport and TRU-Comfort 2 backs. Not available on power tilt backs 16"-17"H with limited adjustability on 18" H. Not available on power recline back 19"-20"H with limited adjustability on 21"H. |
|--|

**10.14. Back Width Package**

**REQUIRED** except when Omit Back Pan/Cushion is selected. Available with +/-4" of seat width.

- |   |   |   |
|---|---|---|
| <input type="checkbox"/> Match Seat Width Selections<br>Part: PTOINDV3468 | <input type="checkbox"/> 20"W-24"W Package<br>Part: XRF000040 | <input type="checkbox"/> 28"W-32"W Package<br>Part: XRF000196 |
| <input type="checkbox"/> 16"W-20"W Package<br>Part: XRF000038             | <input type="checkbox"/> 24"W-28"W Package<br>Part: XRF000195 |   |

**10.15. Back Width**

**REQUIRED** except when Omit Back Pan/Cushion is selected. Part number is PTOINDV plus number following # below. Available with +/-4" of seat width.

- |   |   |   |   |   |   |   |   |   |   |
|---|---|---|---|---|---|---|---|---|---|
| <input type="checkbox"/> 18" W<br>(#3566) | <input type="checkbox"/> 19" W<br>(#3567) | <input type="checkbox"/> 20" W<br>(#3568) | <input type="checkbox"/> 21" W<br>(#3569) | <input type="checkbox"/> 22" W<br>(#3570) | <input type="checkbox"/> 24" W<br>(#3571) | <input type="checkbox"/> 26" W<br>(#3572) | <input type="checkbox"/> 28" W<br>(#3573) | <input type="checkbox"/> 30" W<br>(#3574) | <input type="checkbox"/> 32" W<br>(#3575) |
|---|---|---|---|---|---|---|---|---|---|

TRU-Balance 3 Seating section continued on following page.

**10. TRU-BALANCE 3 SEATING CONTINUED**

**10.16. Overall Back Height**

Static Seat and Power Tilt available height range: 16"-25"H. Odd heights include 4" gap and even heights include 3" gap between seat pan and bottom of back cushion. Power Recline available height range: 19"-28"H. Odd heights include 6" gap and even heights include 7" gap. See the Sport and TRU-Comfort back availability chart for stock order back size ranges.

- 16" Overall Height  
Part: SETHGHT1029
- 17" Overall Height  
Part: SETHGHT1030
- 18" Overall Height  
Part: SETHGHT1031
- 19" Overall Height  
Part: SETHGHT1032
- 20" Overall Height  
Part: SETHGHT1033
- 21" Overall Height  
Part: SETHGHT1034
- 22" Overall Height  
Part: SETHGHT1035
- 23" Overall Height  
Part: SETHGHT1036
- 24" Overall Height  
Part: SETHGHT1037
- 25" Overall Height  
Part: SETHGHT1038
- 26" Overall Height  
Part: SETHGHT1039
- 27" Overall Height  
Part: SETHGHT1040
- 28" Overall Height  
Part: SETHGHT1041

Overall Back Height	Tilt or Static Back Width										Recline Back Width									
	18"	19"	20"	21"	22"	24"	26"	28"	30"	32"	18"	19"	20"	21"	22"	24"	26"	28"	30"	32"
16"	ST	ST	ST	S	ST	ST	S	S	S	S	-	-	-	-	-	-	-	-	-	-
17"	ST	ST	ST	S	ST	ST	S	S	S	S	-	-	-	-	-	-	-	-	-	-
18"	ST	ST	ST	S	ST	ST	ST	ST	S	S	-	-	-	-	-	-	-	-	-	-
19"	ST	ST	ST	S	ST	ST	ST	ST	S	S	ST	ST	ST	S	ST	S	S	S	S	S
20"	ST	ST	ST	S	ST	ST	ST	ST	S	S	ST	ST	ST	S	ST	ST	ST	ST	S	S
21"	ST	ST	ST	S	ST	ST	ST	ST	S	S	ST	ST	ST	S	ST	ST	ST	ST	S	S
22"	ST	ST	ST	S	ST	ST	S	S	S	S	ST	ST	ST	S	ST	ST	ST	ST	S	S
23"	ST	ST	ST	S	ST	ST	S	S	S	S	ST	ST	ST	S	ST	ST	ST	ST	S	S
24"	ST	S	ST	S	S	S	S	S	S	S	ST	ST	ST	S	ST	ST	ST	ST	S	S
25"	ST	S	ST	S	S	S	S	S	S	S	ST	ST	ST	S	ST	ST	S	S	S	S
26"	-	-	-	-	-	-	-	-	-	-	ST	ST	ST	S	ST	ST	S	S	S	S
27"	-	-	-	-	-	-	-	-	-	-	ST	S	ST	S	S	S	S	S	S	S
28"	-	-	-	-	-	-	-	-	-	-	ST	S	ST	S	S	S	S	S	S	S

Sport Back       Sport or TRU-Comfort 2 Back  
**Sport and TRU-Comfort 2 Back Availability Chart**

**10.17. Power Tilt Privacy Flap**

Optional. A privacy flap is standard with power recline systems. This selection is to add a flap for power tilt only systems.

- Privacy Flap .....\$48  
Part: XRF000351

**11. TRU-BALANCE 3 ARMRESTS**

REQUIRED. Select the desired type of armrest below. Please read notes to ensure seating compatibility. Be sure to also select an armrest size and armpad type.

**11.1a. Flip-Back, 8"-16" H Armrests**

Choose a left and right armrest and set height in the following section. Flip-back armrests include an armrest downpost on the Q6 Edge HD. Not available with waterfall armpads. Not available with Back Cane Option.

- Flip-Back, 22"-26"Widths, 8"-16"H Armrests .. \$145 ea. + \$115ea  
L R Left part: ARM2004694. Right part: ARM2004693  
\$115 charge is for downposts
- Flip-Back, 26"-32"Widths, 8"-16"H Armrests .. \$145 ea. + \$115ea  
L R Left part: ARM2004696. Right part: ARM2004695  
\$115 charge is for downposts

**11.1b. Height Adjustable, Removable, Heavy Duty Armrests**

Select a type of HD removable armrests and the set height in the following section. Only one height package is selectable as a normal configuration. Payer coding, coverage, and bundling guidelines may apply.

- 8"-12"H, Single Post, Height Adjustable, Removable, .....\$130 ea.  
Heavy Duty (Sold in pairs)  
Part: XRF000055. HCPCS w/ Static Seat: E0973
- 10"-14"H, Single Post, Height Adjustable, Removable, .....\$130 ea.  
Heavy Duty (Sold in pairs)  
Part: XRF000056. HCPCS w/ Static Seat: E0973

**11.2. Armrest Height**

Set left and right armrest height. Note the available height range in the option description.

**1. Left Armrest Height**

- 8" H PTOINDV3535       10" H PTOINDV3537       12" H PTOINDV3539       14" H PTOINDV3541       16" H PTOINDV3543
- 9" H PTOINDV3536       11" H PTOINDV3538       13" H PTOINDV3540       15" H PTOINDV3542

**2. Right Armrest Height**

- 8" H PTOINDV3450       10" H PTOINDV3452       12" H PTOINDV3454       14" H PTOINDV3456       16" H PTOINDV3458
- 9" H PTOINDV3451       11" H PTOINDV3453       13" H PTOINDV3455       15" H PTOINDV3457

**11.3. Armrest Size**

REQUIRED with armrest selection above. Select a left (L) and right (R) armrest size. This also affects armpad selection.

- 10" Length L part: XRF000057. R part: XRF000123        14" Length L part: XRF000058. R part: XRF000124        17" Length L part: XRF000059. R part: XRF000125  
L R      L R      L R  
Only available with Flip-Back Armrests and straight arm pads.

**12. TRU-BALANCE 3 ARMPADS**

**REQUIRED with armrest selection above except for Drop-In, T-Arm.** Please note that waterfall and gel armpads may cause interference with the seat back.

**12.1a. Omit Armpad**

Select if aftermarket armpads or handpads are required. See [Positioning Components Order Form](#).

- No Armpad Selection..... No Charge  
 Part: XRF000060  
 Not available with Flip-back, 8"-16" H armrests.

**12.1b. Foam Armpads**

- |   |   |
|---|---|
| <input type="checkbox"/> <input type="checkbox"/> Straight, Armpad ..... Standard                 | <input type="checkbox"/> <input type="checkbox"/> Waterfall, Armpad ..... No Charge               |
| <input type="checkbox"/> <input type="checkbox"/> L R Left part: XRF000061. Right part: XRF000126 | <input type="checkbox"/> <input type="checkbox"/> L R Left part: XRF000062. Right part: XRF000127 |
|   | Not available with Flip-back, 8"-16" armrests.  |

**12.1c. Gel Armpads**

Sold only in pairs. Armpads are 3" wide. Manufacturer part number for full length is 314G and desk length is 310G.

**a. Falcon Gel Armpads for Flip-Back Armrests**

- |   |   |
|---|---|
| <input type="checkbox"/> Falcon Gel, Full Length Left and Right ..... \$300<br>Part: ARM1806807         | <input type="checkbox"/> Falcon Gel, Desk Length Left, Full Length Right..... \$250<br>Part: ARM1806810 |
| <input type="checkbox"/> Falcon Gel, Full Length Left, Desk Length Right..... \$250<br>Part: ARM1806813 | <input type="checkbox"/> Falcon Gel, Desk Length Left and Right ..... \$250<br>Part: ARM1806806         |

**b. Falcon Gel Armpads for Tube Style Armrests**

- |  |  |
|--|--|
| <input type="checkbox"/> Falcon Gel, Full Length Left and Right ..... \$300<br>Part: ARM136438         | <input type="checkbox"/> Falcon Gel, Desk Length Left, Full Length Right..... \$250<br>Part: ARM136439 |
| <input type="checkbox"/> Falcon Gel, Full Length Left, Desk Length Right..... \$250<br>Part: ARM136442 | <input type="checkbox"/> Falcon Gel, Desk Length Left and Right ..... \$200<br>Part: ARM136437         |

**13. TRU-BALANCE 3 JOYSTICK BRACKET RECEIVER**

**REQUIRED.**

- |   |   |
|---|---|
| <input type="checkbox"/> Standard Fixed Receiver ..... Standard<br>Part: XRF000042.50 ea. | <input type="checkbox"/> Stealth Height Adjustable Receiver..... \$194<br>Part: XRF000064. HCPCS: K0108 |
|---|---|

**14. TRU-BALANCE 3 HEADRESTS**

**Optional.** Select a headrest pad and bracket.

**14.1a. Headrests, Comfort Plus with Cool Core**

**1. Headrest Pads**

- |   |   |
|---|---|
| <input type="checkbox"/> Stealth 8" Comfort Plus, Cool Core (CP08CCFPM) ..... \$279<br>Part: ST-CP08CCFPM-Q. HCPCS: E0955 | <input type="checkbox"/> Stealth 14" Comfort Plus, Cool Core (CP14CCFPM)..... \$279<br>Part: ST-CP14CCFPM-Q. HCPCS: E0955 |
| <input type="checkbox"/> Stealth 10" Comfort Plus, Cool Core (CP10CCFPM)..... \$279<br>Part: ST-CP10CCFPM-Q. HCPCS: E0955 | <input type="checkbox"/> Stealth 6", Comfort Plus, Cool Core (CP06CCFPM)..... \$279<br>Part: ST-CP06CCFPM-Q. HCPCS: E0955 |

**2. Headrest Brackets**

- |  |   |
|--|---|
| <input type="checkbox"/> Comfort Plus Headrest Mounting Bracket, UniLink..... \$485<br>Part: ST-UHB2D33200-Q. HCPCS: E1028 | <input type="checkbox"/> Comfort Plus Headrest Mounting Bracket, HMO475.....\$271<br>Part: ST-HMO475-17M-Q. |
|--|---|

**15. LAP BELTS**

**REQUIRED.** Select a lap belt below.

**15.1a. Standard Lap Belts**

- |  |  |  |  |  |  |
|--|--|--|--|--|--|
| <input type="checkbox"/> 38" Lap Belt<br>Part: ACCBELT1006 | <input type="checkbox"/> 50" Lap Belt<br>Part: ACCBELT1000 | <input type="checkbox"/> 60" Lap Belt<br>Part: ACCBELT1002 | <input type="checkbox"/> 70" Lap Belt<br>Part: ACCBELT1003 | <input type="checkbox"/> 80" Lap Belt<br>Part: ACCBELT1004 | <input type="checkbox"/> 90" Lap Belt<br>Part: ACCBELT1001 |
|--|--|--|--|--|--|

**15.1b. Optional Lap Belts**

- |   |   |
|---|---|
| <input type="checkbox"/> Stealth Retractable Lap Belt (IFRB5350).....\$175<br>Part: ST-IFRB5350   | <input type="checkbox"/> Stealth Positioning, 60" Single Pull, Push Button, Large ..... \$134<br>Padded (9.5"), 2" Belt (IFPB5314)<br>Part: ST-IFPB5314-Q   |
| <input type="checkbox"/> Stealth Positioning, 60" Single Pull, Push Button, Medium..... \$134<br>(7.5") Padded, 1.5" Belt (IFPB5213)<br>Part: ST-IFPB5213-Q | <input type="checkbox"/> Stealth Positioning, 60" Single Pull, Push Button, Extra Large. \$134<br>Padded (11.5"), 2" Belt (IFPB5315)<br>Part: ST-IFPB5315-Q |



**16. BASE MOUNTED FOOT PLATFORM**

**REQUIRED.** Select a foot platform option. A base mounted foot platform is not available with center mount foot platforms. May interfere with seat mounted legrests.

- Base Mounted Foot Platform ..... Standard  
Part: FRMASMB16077
- Omit Base Mounted Foot Platform ..... No Charge  
Part: PTOINDV3674
- Base Mounted, Bell Shaped Foot Platform ..... No Charge  
Part: RIGASMB7120041

**17. LEGRESTS**

**Optional.** Please note the set lower leg length may vary based on your seat type and seat to floor height selections. All efforts will be made to accommodate legrest length selection, but minor changes may need to be made to meet shipping guidelines. Some adjustments may be needed upon receipt of the unit. Reference the section's compatibility matrix for legrest interference.

**17.1a. Basic Center Mount Foot Platform**

Short	Medium	Long
8"-10"	10"-14"	14"-19"
Basic Center Mount Extension Range		

**1. Basic Center Mount Foot Platform Assembly**

- Basic Center Mount ..... \$200  
Part: XRF000323

**2. Basic Center Mount Extension and Angle Adjustable Footplates**

- Short Extension, 5"x9" Footplates  
Part: XRF000312. HCPCS: K0040
- Medium Lower Extension, 6.25"W Tapered 5"W x 12"D  
Part: XRF000349. HCPCS: K0040
- Medium Extension, 5"x9" Footplates  
Part: XRF000313. HCPCS: K0040
- Long Lower Extension, 6.25"W Tapered 5"W x 12"D  
Part: XRF000350. HCPCS: K0040
- Long Extension, 5"x9" Footplates  
Part: XRF000314. HCPCS: K0040

**3. Basic Center Mount Calf Pads**

- Small Calf Pads, Basic Center Mount ..... \$210  
Part: RIG2104369

**4. Basic Center Mount Additional Accessory**

- Basic Center Mount Footplate Connector ..... \$55  
Part: RIG2104240

**17.1b. Center Mount Foot Platform**

If power articulating foot platform is selected, be sure power legrest control has been selected in earlier section. Select legrest length in following section.

**1. Foot Platform Assembly**

- Power Articulating Foot Platform ..... \$3,145  
Part: PTOINDV3754. HCPCS: E1012+K0040
- Manually Adjustable Center Mount Foot Platform ..... \$265  
Part: XRF000155. HCPCS: K0040

**2. Lower Extension and Angle Adjustable Footplates**

- Medium Extension, 5"x6" Footplates  
Part: RIG150206
- Long Extension, 5"x6" Footplates  
Part: RIG151505
- Medium Extension, 5"x8" Footplates  
Part: RIG150203
- Long Extension, 5"x8" Footplates  
Part: RIG151506
- Medium Extension, 5"x11" Footplates  
Part: RIG150212
- Long Extension, 5"x11" Footplates  
Part: RIG151508

Medium	Long
11"-17"	14"-20"
Power AFP Extension Range	

Medium	Long
10"-16"	13.5"-19.5"
Manually Adjustable Extension Range	

**3a. Center Mount Calf Pads**

Calf pads are standard with a power articulating foot platform.

- Small Calf Pads, AFP Only, 6"Wx6"H (pair) ..... \$210  
Part: RIGASMB7120018
- Large Calf Pads, AFP Only, 8" taper to 6"Wx8"H (pair) ..... \$210  
Part: RIGASMB7120019
- Small Calf Pads, Adjustable Center Mount Only, ..... \$210  
6"Wx6"H (pair)  
Part: XRF000156
- Large Calf Pads, Adjustable Center Mount Only, ..... \$210  
8" taper to 6"Wx8"H (pair)  
Part: XRF000157

**3b. Center Mount Calf Panel**

Calf pads are standard with a power articulating foot platform.

- Small Calf Panel, 18"Wx8"H ..... \$210  
Part: FRMASMB7110007
- Large Calf Panel, 24"Wx8"H ..... \$265  
Part: RIG123824  
Only available on widths of 20" or greater.

**4. Modular Footplate Extension Components**

Optional. Assembly required. May require seating configuration adjustments to avoid front caster interference.

- Width Extension for 8"D Footplates ..... \$55  
Part: RIG152517
- 4" Outside Edge (pair) ..... \$55  
Part: RIG152510  
Only available with 6"D footplates and 8"D footplate extension.
- Width Extension for 11"D Footplates ..... \$55  
Part: RIG152518
- 6" Outside Edge (pair) ..... \$55  
Part: RIG152516  
Available with 8"D and 11"D footplates and 11"D footplate extension.
- 2" Bell-Shaped Length Extension ..... \$55  
Part: RIG152519  
Not compatible with width extensions.

Legrests section continued on following page.

**17. LEGRESTS CONTINUED**

**5. Center Mount Additional Accessories**

- |   |  |
|---|--|
| <input type="checkbox"/> Small Reversed Dartex Calf Pad Cover (pair) ..... \$150<br>Part: RIG133717 | <input type="checkbox"/> Large Foam Calf Pad Cover (pair) ..... \$150<br>Part: ST-CSCS0985-Q |
| <input type="checkbox"/> Large Reversed Dartex Calf Pad Cover (pair) ..... \$150<br>Part: RIG133719 | <input type="checkbox"/> Small Foam Panel Cover (pair) ..... \$150<br>Part: ST-CSCP0920-Q    |
| <input type="checkbox"/> Small Reversed Dartex Panel Cover (pair) ..... \$150<br>Part: RIG134673    | <input type="checkbox"/> Large Foam Panel Cover (pair) ..... \$150<br>Part: ST-CSCP0927-Q    |
| <input type="checkbox"/> Large Reversed Dartex Panel Cover (pair) ..... \$150<br>Part: RIG134744    | <input type="checkbox"/> Footplate Connector ..... \$55<br>Part: RIGASMB7120023              |
| <input type="checkbox"/> Small Foam Calf Pad Cover (pair) ..... \$150<br>Part: ST-CSCS0775-Q        |  |

**6. Stealth Gel Footplate Covers**

- |   |  |
|---|--|
| <input type="checkbox"/> 5"x8" Tapered Gel Footplate Covers ..... \$275<br>Part: ST-FPG5008T-Q. HCPCS: K0108  | <input type="checkbox"/> 5"x8" Tapered Gel Footplate Covers for Width Extension ..... \$275<br>Part: ST-FPG6508T-Q. HCPCS: K0108 |
| <input type="checkbox"/> 5"x11" Tapered Gel Footplate Covers ..... \$275<br>Part: ST-FPG5011T-Q. HCPCS: K0108 | <input type="checkbox"/> 5"x11" Tapered Gel Footplate Cover for Width Extension ..... \$275<br>Part: ST-FPG7011T-Q. HCPCS: K0108 |

**17.1c. Heavy Duty, Drop-In, Swing-Away Legrests**

Select a configuration, a legrest, extension, and type of footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes.

**1. Drop-In Legrest Selection**

- |  |   |
|--|---|
| <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> 70° HD Drop-in for 8.5"-12" Lower Extension ... \$222.50 ea.<br>Pair part: FRMASMB4762. Left: RIG137706. Right: RIG137707 | <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> 70° HD Drop-in 13"-18" or 16"-21" Lower ..... \$222.50 ea.<br>Extension<br>Pair part: FRMASMB4853. Left: RIG137708. Right: RIG137709 |
|--|---|

**2. Lower Extension**

- |   |  |
|---|--|
| <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> 8.5"-12" HD Drop-Ins/12.5"-16.5" ELR, Ext Length<br>Pair part: XRF000101. Left: XRF000102. Right: XRF000103. | <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> 16"-21" HD Drop-Ins/18.5"-22.5" ELR, Ext Length<br>Pair part: FRMASMB4754. L: FRMASMB4725. R: FRMASMB5393 |
| <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> 13"-18" HD Drop-Ins/15.5"-19.5" ELR, Ext Length<br>Pair part: XRF000104. Left: XRF000105. Right: XRF000106   |  |

**3. Angle Adjustable Footplate(s) with Heel Loops**

- |  |  |
|--|--|
| <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> 6"Wx8"D <sup>1</sup> ..... \$87.50 ea.<br>Pair part: FRMASMB4738. Left: RIG137793. Right: RIG137795.<br>HCPCS: K0040+E0951  | <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> 8"Wx11"D <sup>2</sup> ..... \$87.50 ea.<br>Pair part: FRMASMB4744. Left: RIG137798. Right: RIG137799.<br>HCPCS: K0040+E0951<br>May cause interference on seat widths less than 20". |
| <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> 6"Wx11"D <sup>1</sup> ..... \$87.50 ea.<br>Pair part: FRMASMB4741. Left: RIG137796. Right: RIG137797.<br>HCPCS: K0040+E0951<br>May cause interference on seat widths less than 18". |  |

**17.1d. Heavy Duty, Elevating Legrests**

Select a configuration, a legrest, extension, and type of footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes.

**1. Heavy Duty, Elevating Legrest Selection**

- |   |
|---|
| <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> Heavy Duty, Elevating Legrest(s) with Calf Pad(s) ..... \$275 ea.<br>Pair part: FRMASMB6702. Left: FRMASMB8444. Right: FRMASMB8445. HCPCS: E0990 |
|---|

**2. Lower Extension**

- |  |  |
|--|--|
| <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> 11"-14" Length<br>Part: FRMASMB4755. Left: RIG138026. Right: RIG138030<br>For reinforced footplates only.   | <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> 14"-17" Length<br>Part: FRMASMB4756. Left: RIG138033. Right: RIG138035<br>For reinforced footplates only.   |
| <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> 8.5"-12" HD Drop-Ins/12.5"-16.5" ELR, Ext Length<br>Pair part: XRF000101. Left: XRF000102. Right: XRF000103.<br>For angle adjustable footplates only. | <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> 13"-18" HD Drop-Ins/15.5"-19.5" ELR, Ext Length<br>Pair part: XRF000104. Left: XRF000105. Right: XRF000106<br>For angle adjustable footplates only. |

**3. Reinforced Footplates(s)**

- |   |  |
|---|--|
| <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> 7.25"Wx6"D Extruded Footplate ..... Standard<br>Part pair: FRMASMB4748. Left: RIG138055.<br>Right: RIG138056 | <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> 8.5"Wx6"D Extruded Footplate ..... Standard<br>Part pair: FRMASMB4747. Left: RIG138060.<br>Right: RIG138062 |
|---|--|

**4. Angle Adjustable Footplates(s) with Heel Loops**

- |   |  |
|---|--|
| <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> 6"Wx8"D ..... \$87.50 ea.<br>Pair part: FRMASMB4738. Left: RIG137793. Right: RIG137795.<br>HCPCS: K0040+E0951  | <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> 8"Wx11"D ..... \$87.50 ea.<br>Pair part: FRMASMB4744. Left: RIG137798. Right: RIG137799.<br>HCPCS: K0040+E0951<br>May cause interference with seat widths less than 20" |
| <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> 6"Wx11"D ..... \$87.50 ea.<br>Pair part: FRMASMB4741. Left: RIG137796. Right: RIG137797.<br>HCPCS: K0040+E0951 |  |

Legrests section continued on following page.

**17. LEGRESTS CONTINUED**

**17.2. Set Legrest Extension Length**

**REQUIRED** with center mounts and matching left and right rigging configurations. Part number is PTOINDV plus number following # below.

**1a. Centermount Extension Length**

- 8.5" L (#3681)    9.5" L (#3683)    10.5" L (#3685)    11.5" L (#3687)    12.5" L (#3689)    13.5" L (#3691)    14.5" L (#3693)    15.5" L (#3695)    16.5" L (#3697)    17.5" L (#3699)    18.5" L (#3701)
- 9" L (#3682)    10" L (#3684)    11" L (#3686)    12" L (#3688)    13" L (#3690)    14" L (#3692)    15" L (#3694)    16" L (#3696)    17" L (#3698)    18" L (#3700)    19" L (#3702)

**1b. Front Rigging Extension Length**

Non-matching or individual legrest configurations lengths are not set by the manufacturer. See previous section for center mount extension lengths.

- 8.5" L (#3473)    9.5" L (#3475)    10.5" L (#3477)    11.5" L (#3479)    12.5" L (#3481)    13.5" L (#3483)    14.5" L (#3485)    15.5" L (#3487)    16.5" L (#3489)    17.5" L (#3491)
- 9" L (#3474)    10" L (#3476)    11" L (#3478)    12" L (#3480)    13" L (#3482)    14" L (#3484)    15" L (#3486)    16" L (#3488)    17" L (#3490)    18" L (#3492)

**18. TRU-BALANCE 3 POWER POSITIONING ELECTRONICS**

**REQUIRED** with power seating or legrest options.

**18.1a. Tilt Only through Toggle**

- Single Actuator Control, Power Tilt Only, Through Toggle (NE 4-Key Joystick Only)..... No Charge  
Part: ELE1905907  
Select toggle side. May require extended lead time when selected with Fender Lights.

**1. Toggle Mount Side Selection**

- Left Side Toggle  
Part: PTOINDV1839
- Right Side Toggle  
Part: PTOINDV1838

**18.1b. NE+ 6 Key Controls through Joystick**

- Single or Dual Actuator Control Through NE+ Joystick..... No Charge  
Part: XRF000065

**18.1c. Q-Logic 3 Controls through Joystick**

- Single Actuator Through Q-Logic 3 Controller and Interface (AM3)..... No Charge
- 2 Actuators Through Q-Logic 3 Controller and Interface (AM3).....\$1,250  
Not funded by ADP
- 3 Actuators Through Q-Logic 3 Controller and Interface (AM3)..... \$1,500  
Not funded by ADP
- Up to 5 Actuators Through Q-Logic 3 Controller and Interface (AAM).....\$1,675  
ADP: SECX15005

**Optional Switch Control**

Available with Q-Logic 3 only. 4-Way switches will be mounted to same side of joystick.

- 4-Function Push Button.....\$225  
Part: CTL123218  
See next section for switch mounting.  
Only available with AAM.
- 4-Way Toggle.....\$575  
Part: ELEASMB5925  
Only available with AAM.
- 4-Way Button Toggle.....\$575  
Part: ELEASMB5926  
Only available with AAM.
- iAccess Flat Button.....\$575  
Part: CTL144009  
See next section for switch mounting.
- iAccess Toggle Button.....\$575  
Part: CTL144010  
See next section for switch mounting.



**18.2. 4-Function Push Button/iAccess Mount**

**REQUIRED** with 4-function push button and iAccess options. Under arm mount is not available with TRU-Balance 3 HD armrests.

- Joystick Mount ..... No Charge  
Part: MEC1907500
- Attendant Mount ..... \$35  
Part: MEC1906564
- Under Arm Mount, Flip-Back Armrest..... No Charge  
Part: MEC1906531
- Under Arm Mount, Flip-Back Armrest, plus iAccess Guard ..... \$35  
Part: MEC2001400
- Under Arm Mount, Tube Style Armrest ..... No Charge  
Part: MEC1906530
- Under Arm Mount, Tube Style Armrest, plus iAccess Guard ..... \$35  
Part: MEC2001401

**1. Under Arm Mount Side Selection**

- Left Side Under Arm Mount  
Part: PTOINDV3678
- Right Side Under Arm Mount  
Part: PTOINDV3679

**19. TRU-BALANCE 3 ATTENDANT CONTROL**

**Optional.** \*Attendant control will be at no charge when ordered with Tilt and Recline. Mounting bracket is included.

- Attendant Control ..... \$950  
Part: XRF000078. HCPCS: E2331

**20 FULL LIGHT PACKAGE**

**Optional.** Requires NE+ 8 Key or Q-Logic 3 Joystick for Full Light Package.


- Seat Mounted, LED, Full Light Package ..... \$770  
Part: XRF000130



**21. ACCESSORIES**

**Optional.** Multiple accessories can be mounted along the various tracks but some interference may occur. Accessories may need to be mounted by the provider.

- TRU-Balance 3 Cane/Crutch Holder ..... \$100  
Part: ACC152515. HCPCS: E2207
- Quantum Backpack ..... \$200  
Part: ACC1905853
- TRU-Balance 3 Transfer Handle (D) ..... \$152.50 ea.  
Part: MEC2205680. HCPCS: E0705  
Not available with Heavy Duty armrests.
- TRU-Balance 3 Cup Holder ..... \$20  
Part: ACC125003  
For use with TRU-Balance 3 (track style) armrests.
- Cup Holder (Tube-style Armrest) ..... \$20  
Part: ACCASMB242.50 ea.4  
For use with tube-style armrests.
- TRU-Balance 3 Personal Item Hook ..... \$20  
Part: ACC125008  
Not available with back cane option.
- TRU-Balance 3 Reflector Kit ..... \$35  
Part: ACC124670
- TRU-Balance 3 Backpack Holder (Pair) (A) ..... \$42  
Part: ACC147774  
Not available with back cane option.
- TRU-Balance 3 Front Mounted Cell Phone/ Media Holder, X-Grip (E) ..... \$95  
Part: ACC155835  
Not available with Heavy Duty armrests.
- TRU-Balance 3 Rear View Mirror ..... \$35  
Part: ACC155833  
Not available with Heavy Duty armrests.
- TRU-Balance 3 Oxygen Holder ..... \$110  
Part: ACC125006. HCPCS: E2208  
Includes accessory bar, a credit will be applied PTOIN-DV3769.
- TRU-Balance 3 Clothing Guard (Pair) ..... \$275  
Part: ACC125007  
Not available with Heavy Duty armrests.
- TRU-Balance 3 Glove Box (C) ..... \$85  
Part: ACC125011
- Stealth Large Tablet Holder ..... \$565  
Part: ACC141818  
Not available with Heavy Duty armrests.
- Trilogy Vent Tray ..... \$420  
Part: ACC136791. HCPCS: E1029  
Includes accessory bar, a credit will be applied PTOIN-
- Compact Vent Tray ..... \$1,395  
Part: ACC133889. HCPCS: E1029
- Push Handles (pair) ..... \$140  
Part: ACC2005027  
Not available with oxygen holder, back cane option, or vent tray. Required with attendant manual recline.
- Accessory Bar ..... \$90  
Part: XRF000128  
Horizontal rail attached to rear of seat. Standard with the Trilogy vent tray and oxygen holder. Not available with back cane option.
- Large, Black, Weather Cover ..... \$240  
Part: ACCCOVR1021
- Q-Logic EX/NE Programmer ..... \$625  
Part: ELEASMB6976
- Privacy Bag ..... \$80  
Part: ACC2005723
- VOCSN Vent Tray for Power Tilt, Seat Mount .... \$1,000  
Part: ACC2000162. HCPCS: E1029
- VOCSN Vent Tray for Power Lift, Seat Mount .... \$1,000  
Part: ACC2000159. HCPCS: E1029
- VOCSN Vent Tray for Accessory Bar Mount ..... \$1,000  
Part: ACC2000166. HCPCS: E1029  
Includes accessory bar, a credit will be applied PTOIN-DV3769.
- Flip-Back Armrest, Side Track Accessory Bracket .... \$25  
Part: MEC1809178  
Only available with flip-back, 8"-16" H armrests. This adds a 47mm channel to the side of the armrest to allow mounting of accessories with t-nuts.
- Personal Item Pouch ..... \$175  
Part: ACC2101017
- Depth Adjustable Vent Tray ..... \$1,000  
Part: ACC2302500. HCPCS: E1029  
Includes accessory bar, a credit will be applied PTOIN-DV3769.



**23. POSITIONING COMPONENTS**

Optional. More options are available on the [Positioning Components Order Form](#).

**23.1. Stealth Thoracic Laterals**

Swing-away laterals may interfere with TRU-Balance 3 armrests and 2-post, flip up armrests. These laterals are not available with the back cane option.

**1. Thoracic Lateral Bracket**

- Swing-away, Rail Mounts, pair (TWBTLTC) ..... \$714  
Part: INDPART2677. HCPCS: E1028  
Not available with Redesign back assemblies.
- Lateral Hardware, Q-Back Right ..... \$357  
Part: ST-TWBLT-QL-R. HCPCS: E1028
- Lateral Hardware, Q-Back Left ..... \$357  
Part: ST-TWBLT-QL-L. HCPCS: E1028

**2. Thoracic Lateral Pads**

- Curved Pads, 3"x5", pair (TWBL3X5)..... \$164  
Part: INDPART2663. HCPCS: E0956
- Flat Pads, 3"x5", pair (TWBLP3X5) ..... \$164  
Part: INDPART2676. HCPCS: E0956
- Curved Pads, 4"x6", pair (TWBL4X6)..... \$164  
Part: INDPART2664. HCPCS: E0956
- Flat Pads, 4"x6", pair (TWBLP4X6) ..... \$164  
Part: INDPART2675. HCPCS: E0956

**3. Thoracic Lateral Installation**

- Install Laterals  
Part: PTOINDV3499
- Do Not Install  
Part: PTOINDV3500

**23.2. Stealth UniLink Pelvic/Thigh Guide**

**1. UniLink Pelvic/Thigh Hardware**

- UniLink, 5" Height Adjustable, .....\$405 ea.  
Removable Mounting Hardware  
Part: ST-UTTP5-Q. HCPCS: E1028
- UniLink, 7" Height Adjustable, .....\$423 ea.  
Removable Mounting Hardware  
Part: ST-UTTP7-Q. HCPCS: E1028

**2. UniLink Pelvic/Thigh Pads**

- UniLink, 4" x 4" Pad ..... \$164 ea.  
Part: ST-UTP544056100-Q. HCPCS: E0953/E0956
- UniLink, 4" x 12" Pad ..... \$164 ea.  
Part: ST-UTP54D056100-Q. HCPCS: E0953/E0956
- UniLink, 4" x 8" Pad..... \$164 ea.  
Part: ST-UTP548056100-Q. HCPCS: E0953/E0956



**3. UniLink Pelvic/Thigh Additional Links**

**4. UniLink Pelvic/Thigh Guide Installation**

- Install on Left  
Part: PTOINDV3585
- Install on Right  
Part: PTOINDV3586
- Do Not Install  
Part: PTOINDV3719

**23.3. Stealth Pelvic/Thigh Guides**

Pelvic/thigh guides may interfere with the joystick receiver block and down posts. Please provide quantity next to check box.

**1. Pelvic/Thigh Guide Bracket**

- Flip Down Bracket (TWBTLADD-HW) ..... \$343 ea.  
Part: INDPART2804. HCPCS: E1028
- Flip Down with LEFT Full Surface Contact ..... \$459  
Feature (TWBADD-L & LPHW407)  
Part: POS140428. HCPCS: E1028
- Removable Bracket (LP430)..... \$383 ea.  
Part: INDPART2805. HCPCS: E1028
- Removable Bracket with 1" Offset ..... \$383 ea.  
(LP430 & LHW-121)  
Part: INDPART2854. HCPCS: E1028
- Flip Down with RIGHT Full Surface Contact..... \$459  
Feature (TWBADD-R & LPHW407)  
Part: POS140429. HCPCS: E1028

**2. Pelvic/Thigh Guide Adapter**

- 2" Height Adapter Plate (LHW-122)..... \$35  
Part: INDPART2821. HCPCS: K0108  
Height adapter is needed with TRU-Comfort seat cushions. Not available with the flip down bracket.
- 2 1/4" Adapter Plate (LHW-129) ..... \$82  
Part: POS140430. HCPCS: K0108  
Only available with the flip down full surface bracket.

**3. Pelvic/Thigh Guide Pad**

- Short Pad, 4"x4" (TWBADD-4X4)..... \$85 ea.  
Part: INDPART2801. HCPCS: E0953/E0956
- Long Pad, 4"x12" (TWBADD-4X12) ..... \$85 ea.  
Part: INDPART2803. HCPCS: E0953/E0956
- Medium Pad, 4"x8" (TWBADD-4X8)..... \$85 ea.  
Part: INDPART2802. HCPCS: E0953/E0956

**4. Pelvic/Thigh Guide Installation**

- Install on Left  
Part: PTOINDV3585
- Install on Right  
Part: PTOINDV3586
- Do Not Install  
Part: PTOINDV3719



**23.4. Stealth Positioning Chest Harness**

- Upper Classic, Medium (IFPS1103) ..... \$221  
Part: ST-IFPS1103-Q. HCPCS: E0960
- Upper Contour, Medium (IFPS2103)..... \$221  
Part: ST-IFPS2103-Q. HCPCS: E0960
- Upper Classic, Large (IFPS1104) ..... \$221  
Part: ST-IFPS1104-Q. HCPCS: E0960
- Upper Contour, Large (IFPS2104) ..... \$221  
Part: ST-IFPS2104-Q. HCPCS: E0960

Positioning components section continued on following page.

**23. POSITIONING COMPONENTS CONTINUED**

**23.5. Stealth Positioning Chest Straps**

- |   |  |
|---|--|
| <input type="checkbox"/> Chest Strap, Medium w/ Belt Anchor, ..... \$140<br>(IFCS4113)<br>Part: ST-IFCS4113-Q. HCPCS: E0960 | <input type="checkbox"/> Chest Strap, Large w/ Belt Anchor, ..... \$140<br>(IFCS4114)<br>Part: ST-IFCS4114-Q. HCPCS: E0960 |
|---|--|



**23.6. Stealth Swing-Away Headrest Mounts**

These swing-away headrest mounts are only available with the Stealth Comfort Series headrests.


**1. Facial Lateral Hardware**

- |  |   |
|--|---|
| <input type="checkbox"/> UniLink Facial Lateral, Swing Away, Right..... \$496<br>Part: ST-SAR2A-R-Q. HCPCS: E1028                                | <input type="checkbox"/> UniLink Facial Lateral, Dual Swing Away ..... \$765<br>2 Upper & 1 Lower Left<br>Part: ST-SAR4A-DBU-LL-Q. HCPCS: E1028 |
| <input type="checkbox"/> UniLink Facial Lateral, Swing Away, Left ..... \$496<br>Part: ST-SAR2A-L-Q. HCPCS: E1028                                | <input type="checkbox"/> UniLink Facial Lateral, Dual Swing Away ..... \$674<br>2 Upper<br>Part: ST-SAR4A-DBU-Q. HCPCS: E1028                   |
| <input type="checkbox"/> UniLink Facial Lateral, Dual Swing Away ..... \$765<br>2 Upper & 1 Lower Right<br>Part: ST-SAR4A-DBU-RL-Q. HCPCS: E1028 |   |



**2. Facial Lateral Pads and Accessories**


- |  |   |
|--|---|
| <input type="checkbox"/> Egg Switch Mount (ESM) (A) ..... \$70 ea.<br>Part: INDPART2815                                  | <input type="checkbox"/> Facial Pad, 6"x2", Cool Core (920M-Q)..... \$96 ea.<br>Part: ST-920M-Q. HCPCS: K0108   |
| <input type="checkbox"/> Figure 8 Bracket with 13" Rod..... \$366 ea.<br>and Egg Switch Mount (B)<br>Part: ST-F8675-13-Q | <input type="checkbox"/> Swing-Away Joystick Mount, 18" Rod (C)..... \$100 ea.<br>Part: ST-SUS18-Q<br>Allows mounting of facial pads and attaches to the headrest |
| <input type="checkbox"/> Facial Pad, 4".5x2", Cool Core (900M-Q) .... \$96 ea.<br>Part: ST-900M-Q. HCPCS: K0108          |   |



**23.7. Stealth Elbow Support**

**1. Elbow Support Bracket**

- |  |   |
|--|---|
| <input type="checkbox"/> Left, Swing-Away, Elbow Support ..... \$279<br>Bracket Tube Style Arms (ARMS-ES-L)<br>Part: POS124064. HCPCS: E1028               | <input type="checkbox"/> Right, Swing-Away, Elbow Support..... \$279<br>Bracket Tube Style Arms (ARMS-ES-R)<br>Part: POS124066. HCPCS: E1028                |
| <input type="checkbox"/> Left, UniLink, Fixed, Elbow Support Bracket..... \$309<br>for Flip-Back, 8"-16"H Armrests (UETB4021FL-Q)<br>Part: ST-UETB4021FL-Q | <input type="checkbox"/> Right, UniLink, Fixed, Elbow Support Bracket..... \$309<br>for Flip-Back, 8"-16"H Armrests (UETB4021FR-Q)<br>Part: ST-UETB4021FR-Q |



**2. Elbow Support Pad**

- |  |   |
|--|---|
| <input type="checkbox"/> 3"Wx3"H Elbow Pad (ARMS-ES-3X3) ..... \$85 ea.<br>Part: POS124065. HCPCS: K0108 | <input type="checkbox"/> 4"Wx4"H Elbow Pad (ARMS-ES-4X4)..... \$85 ea.<br>Part: POS124067. HCPCS: K0108 |
|--|---|

**23.8. Stealth Foot Boxes**

- |   |  |
|---|--|
| <input type="checkbox"/> Stealth Small Foot Boxes, 5"Wx8"Dx7"H..... \$725<br>(SPLE107) (pair)<br>Part: INDPART2786. HCPCS: E0954    | <input type="checkbox"/> Stealth Large Foot Boxes, 7"Wx12"Dx11"H..... \$725<br>(SPLE109) (pair)<br>Part: INDPART2788. HCPCS: E0954 |
| <input type="checkbox"/> Stealth Medium Foot Boxes, 6"Wx10"Dx9"H ..... \$725<br>(SPLE108) (pair)<br>Part: INDPART2787. HCPCS: E0954 |  |

**23.9. Stealth Positioning Shoe Holder**

- |  |   |
|--|---|
| <input type="checkbox"/> Small Shoe Holder, 3/4" Webbing (IFSH4211)..... \$93 ea.<br>Part: ST-IFSH4211-Q. HCPCS: E0951+E0952   | <input type="checkbox"/> Large Shoe Holder, 3/4" Webbing (IFSH4213) ..... \$93 ea.<br>Part: ST-IFSH4213-Q. HCPCS: E0951+E0952 |
| <input type="checkbox"/> Medium Shoe Holder, 3/4" Webbing (IFSH4212) ..... \$93 ea.<br>Part: ST-IFSH4212-Q. HCPCS: E0951+E0952 |   |

**24. ALTERNATIVE CONTROLS AND ELECTRONICS**

**Optional.** Select a type of specialty controls below. More options are available on the [Alternative Controls and Electronics order form](#).

**24.1. Enhanced Display**

- Q-Logic 3 Enhanced Display Kit..... \$2,045  
Part: CTL167726

Select a color option below. Includes enhanced display, flexible gooseneck and mounting bracket. Features: 3.5" LCD display; up to 5 profiles to control drive, seat, and auxiliary functions; Bluetooth connectivity for mobile devices and mouse emulation; integrated single switch scanner; IR controls.



Enhanced Display

**Q-Logic 3 Enhanced Display Shroud**

Required with Q-Logic 3 Enhanced Display selection above. Reference colors shown on page 2.

- |  |  |   |   |
|--|--|---|---|
| <input type="checkbox"/> Black, Soft Touch<br>Part: DGN137728  | <input type="checkbox"/> Night Train<br>Part: DGN137729-PBK3 | <input type="checkbox"/> Pretty in Pink<br>Part: DGN137729-PPK3 | <input type="checkbox"/> Green Machine<br>Part: DGN137729-PGN11 |
| <input type="checkbox"/> Red Corvette<br>Part: DGN137729-PRD1  | <input type="checkbox"/> Grey Street<br>Part: DGN137729-PSL3 | <input type="checkbox"/> Midnight Blue<br>Part: DGN137729-PBU10 |   |
| <input type="checkbox"/> Back in Black<br>Part: DGN137729-PBK1 | <input type="checkbox"/> Deep Purple<br>Part: DGN137729-PPR5 | <input type="checkbox"/> Electric Blue<br>Part: DGN137729-PBU9  |   |

**24.2. SCIM Module**

- SCIM Module Kit for Q-Logic 3..... \$850  
Part: CTL168055

Allows use of joystick display for specialty devices. Available with Q-Logic 3 EX Joystick selection above. Up to 5 profiles to control drive, seat, and auxiliary functions, excluding sip-n-puff and single switch functions. A switch is required to operate the input device attached to the SCIM module.



SCIM Module

**24.3. Q-Logic 3 SAJ**

- Q-Logic 3 EX SAJ..... \$850  
Part: CTL143980
- Must be used in conjunction with Enhanced Display or Q-Logic 3 Joystick.



EX SAJ

**24.4. Head Array**

- |  |   |
|--|---|
| <input type="checkbox"/> i-Drive 4.0 Bluetooth Tri-array Series,..... \$7,406<br>3 Switch with Egg Switch<br>Part: ST-IDHBT200-U-Q. HCPCS: E2330 | <input type="checkbox"/> Adult, Curved, i-Drive Tri-Array with ..... \$7,406<br>Mounting Hardware<br>Part: ST-IDHBT400-U-Q. HCPCS: E2330  |
| <input type="checkbox"/> i-Drive 4.0 Bluetooth Pro Series,..... \$7,486<br>3 Switch with Egg Switch<br>Part: ST-IDHBT300-U-Q. HCPCS: E2330       | <input type="checkbox"/> Junior, Curved, i-Drive Tri-Array with ..... \$7,406<br>Mounting Hardware<br>Part: ST-IDHBT450-U-Q. HCPCS: E2330 |



Adult Pro-Series

**24.5. Sip and Puff Drive**

**1. Sip and Puff Module**

- Sip and Puff Module for Q-Logic 3..... \$1,695  
Part: CTL168056. HCPCS: E2325  
Includes module and mounting hardware.



Sip & Puff Module

**2. Sip and Puff Straw and Mount**

The straw and mount do not include electrical components to operate the Sip&Puff with a head array. This component requires the selection of a Q-Logic Sip&Puff module or additional module from Stealth Products for operation.

- Stealth Sip and Puff for Stealth i-Drive, ..... \$496  
Head Array (SP625-U)  
Part: ST-SP625-U-Q. HCPCS: E2326  
Includes straw and mounting hardware. Use with Stealth i-Drive Head Arrays.



Stealth Sip & Puff

**3. Sip & Puff Accessories**

- Sip & Puff Straws, 12 pack (SP625-S12)..... \$47  
Part: ST-SP625-S12-Q

Alternative controls and electronics section continued on following page.

**24. ALTERNATIVE CONTROLS AND ELECTRONICS CONTINUED**

**24.6. Egg Switches**

- |  |  |
|--|--|
| <input type="checkbox"/> <u>QTY</u> Red Egg Switch (ESRED)..... \$116<br>Part: SWTMCRO1077   | <input type="checkbox"/> <u>QTY</u> Green Egg Switch (ESGRN) ..... \$116<br>Part: SWTMCRO1078  |
| <input type="checkbox"/> <u>QTY</u> Blue Egg Switch (ESBLU) ..... \$116<br>Part: SWTMCRO1079 | <input type="checkbox"/> <u>QTY</u> Yellow Egg Switch (ESYEL) ..... \$116<br>Part: SWTMCRO1080 |
| <input type="checkbox"/> <u>QTY</u> Black Egg Switch (ESBLK)..... \$116<br>Part: SWTMCRO1081 |  |



**24.7. Switch Adapter**

- 2 to 1 Switch Adapter ..... \$135  
Part: DWR1361H037  
Connects 2 single switches to a stereo jack to act as a dual switch.

**24.8. Environmental Controls**

- Q-Logic 3 ECU Module Kit ..... \$1,780  
Part: CTL168057



**24.9. USB Mobile Device Charger**

- USB Mobile Device Charger ..... \$185  
Part: ELE168069  
Can charge 1 device. A multiplier harness is required and must be added when a Static Seat or Power Tilt Only seating system is chosen. The selection of more than 3 modules requires an additional multiplier.