



Pride Mobility Products Co.

5096 South Service Rd, Beamsville, ON, L0R 1B3

Phone: 888-570-1113 | Fax: 866-514-1303 | Email: productorders@pridemobility.com

# Q6 Edge® 3 Order Form with TRU-Balance® 4 Positioning

Account Number: \_\_\_\_\_ Date: \_\_\_\_\_

Provider Name: \_\_\_\_\_

Contact: \_\_\_\_\_

Phone: \_\_\_\_\_ Fax: \_\_\_\_\_

Email: \_\_\_\_\_

PO Number: \_\_\_\_\_

Marked for: \_\_\_\_\_

Ship to Address: \_\_\_\_\_

City: \_\_\_\_\_ State: \_\_\_\_\_ Zip: \_\_\_\_\_

## EDGE 3



### INTRODUCTION

This form is interactive when viewed with [Adobe Acrobat Reader](#). Complete the form by placing checks in the desired boxes and provide information in the interactive fields. Buttons shown at the bottom of the form may be utilized to print or submit the order form through a desktop email application. To email via a web-based application, please 'Save As' and attach the PDF to your email.

This order form contains a large variety of options to fulfill various patient needs. Descriptions and section notes, such as optional or required, are included to help you complete your order. Other part numbers can be found on the quote acknowledgment after this order form is submitted. Please contact Quantum® Sales at 888-570-1113 if assistance is needed. Send the completed order form by fax (866-514-1303) or email ([productorders@pridemobility.com](mailto:productorders@pridemobility.com)). Incomplete forms may delay the quote or order. Customer service will contact you if the order is incomplete or if there are compatibility issues. If special order requests are needed, be sure that the Patient Information section is completed

All prices are MSRP. Prices, specifications, part numbers, and availability are subject to change without notice. Prices and part numbers as shown are only available when configured as a complete power chair. Please contact technical service for accurate parts ordering.

**Thank you for choosing Quantum!**


### PATIENT INFORMATION

Weight: \_\_\_\_\_ Height: \_\_\_\_\_

A. Shoulder Width: \_\_\_\_\_ D. Max Sitting Height: \_\_\_\_\_ G. Thigh Depth: \_\_\_\_\_

B. Chest Width: \_\_\_\_\_ E. Shoulder Height: \_\_\_\_\_ H. Lower Leg Length: \_\_\_\_\_

C. Hip Width: \_\_\_\_\_ F. Axilla Height: \_\_\_\_\_ I. Elbow Height: \_\_\_\_\_



\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

### 1. BASE MODEL

**REQUIRED.** 300 lb. weight capacity with standard 6mph programmable motors and 8 amp off board charger. XLR USB charger is included. Dimensions of base without legrests: 24"Wx35.5"L

Quantum Q6 Edge 3 3MP-SS ..... \$8,450  
Part: Q6 EDGE 3. ADP: WP0001575 (ADULT), WP0001576 (PEDIATRIC)

**2. DRIVE MOTORS**


















**REQUIRED.**

- |   |   |
|---|---|
| 6 mph 2 Pole Drive Motors..... Standard<br>Part: MOT147433  | 7 mph, 4 Pole Drive Motors..... \$700<br>Part: MOT1902455<br>Requires Q-Logic 3 electronics.  |
| 6 mph, 2 Pole, Accu-Trac Motor Technology..... \$995<br>Part: MOT147881.<br>Requires Q-Logic 3 electronics. Delivers enhanced traction and tracking, especially with switched drive controls. | 7 mph, 4 Pole, Accu-Trac Motor Technology..... \$1,695<br>Part: MOT2008501<br>Requires Q-Logic 3 electronics. Delivers enhanced traction and tracking, especially with switched drive controls. |

**3. BASE COLOR**

**REQUIRED.** Select one shroud color. Colors shown here may differ from actual product depending on monitor/printer calibration used.

**3.1a. Standard Shrouds**

- |   |  |  |  |
|---|--|--|--|
|  <b>Red Corvette</b><br>Part: DGN1801101-PRD1  |  <b>White Satin</b><br>Part: DGN1801101-PWT2      |  <b>Midnight Blue</b><br>Part: DGN1801101-PBU10 |  <b>Military Green</b><br>Part: DGN1801101-PGN8 |
|  <b>Cherry Bomb</b><br>Part: DGN1801101-PRD4   |  <b>Deep Purple</b><br>Part: DGN1801101-PPR5      |  <b>Electric Blue</b><br>Part: DGN1801101-PBU9  |  <b>Sugar Plum</b><br>Part: DGN1801101-PPR7     |
|  <b>Back in Black</b><br>Part: DGN1801101-PBK1 |  <b>Raspberry Beret</b><br>Part: DGN1801101-PPK4  |  <b>Blue Jean</b><br>Part: DGN1801101-PBU11     |  |
|  <b>Night Train</b><br>Part: DGN1801101-PBK3   |  <b>Pretty in Pink</b><br>Part: DGN1801101-PPK3   |  <b>Lemon Crush</b><br>Part: DGN1801101-PYL4    |  |
|  <b>Grey Street</b><br>Part: DGN1801101-PSL3   |  <b>Evolution Orange</b><br>Part: DGN1801101-POR7 |  <b>Green Machine</b><br>Part: DGN1801101-PGN11 |  |

**3.1b. Hydrographic Shrouds**

- |  |  |   |
|--|--|---|
| <input type="checkbox"/>  <b>Military Digital Camo</b> ..... \$325<br>Part: DGN1801101-WCM7 | <input type="checkbox"/>  <b>Faded Glory</b> ..... \$325<br>Part: DGN1801101-WFL2 | <input type="checkbox"/>  <b>Carbon Fiber</b> ..... \$325<br>Part: DGN1801101-WCB1 |
|--|--|---|

**4. BASE MOUNTED LIGHTS**

**Optional.** Select one option. Seat Mounted, LED, Full Light Package option can be found later on this form.

- |  |  |
|--|--|
| Switch Fender Lights..... No Charge<br>Part: DGN154117<br>Fender lights operate by pressing included single switch button. Fender lighting is not available with Single Actuator Control through Toggle Box. | Base Mounted, LED, Full Light Package ..... \$770<br>Part: XRF000129<br>Requires Q-Logic 3 EX Joystick for Full Lights or Q-Logic 3e Joystick. |
|--|--|

**5. QUANTUM BACKUP CAMERA**

**Optional.** The backup camera may interfere with armrest mounted accessories such as media device holder.

- Quantum Backup Camera ..... \$475  
Part: MEC1910365

**1. Optional Mounting Bracket**

- |  |  |
|--|--|
| Backup Camera Headrest Bracket..... \$30<br>Part: MEC1910373<br>Works with .5" x .5" vertical supports only. | Backup Camera Track Mount Bracket ..... \$35<br>Part: MEC2104227 |
|--|--|

**2. Optional Display Mounting Bracket**

- Backup Camera Q-Logic Joystick Mount Bracket ..... \$145  
Part: MEC2106123



**6. DRIVE WHEELS**

**REQUIRED.** Select a type of drive wheels and accent color. Knobby tires increase overall base width by 1.5".

**6.1. Drive Wheels**

- |  |  |
|--|--|
| 3"x14" Flat Free Black Tires ..... Standard<br>Part: PTOINDV3525 (pair)    | 3"x14" Pneumatic Black Tires ..... No Charge<br>Part: WHL1713925 x2        |
| 4"x14" Knobby, Flat Free Black Tires ..... \$60 ea.<br>Part: WHL2006647 x2 | 4"x14" Knobby, Pneumatic Black Tires ..... \$60 ea.<br>Part: WHL2006822 x2 |

**6.2. Drive Wheel Accent Color**

**Optional.** Not available with knobby tires.

- |   |   |  |
|---|---|--|
| <b>Back in Black</b> .....Standard<br>Part: DGN164245-PBK1 x2 | <b>White Satin</b> .....No Charge<br>Part: DGN164245-PWH2 x2      | <b>Midnight Blue</b> ..... No Charge<br>Part: DGN164245-PBU10 x2 |
| <b>Red Corvette</b> .....No Charge<br>Part: DGN164245-PRD1 x2 | <b>Deep Purple</b> .....No Charge<br>Part: DGN164245-PPR5 x2      | <b>Electric Blue</b> ..... No Charge<br>Part: DGN164245-PBU9 x2  |
| <b>Cherry Bomb</b> .....No Charge<br>Part: DGN164245-PRD4 x2  | <b>Raspberry Beret</b> .....No Charge<br>Part: DGN164245-PPK4 x2  | <b>Blue Jean</b> ..... No Charge<br>Part: DGN164245-PBU11 x2     |
| <b>Night Train</b> .....No Charge<br>Part: DGN164245-PBK3 x2  | <b>Pretty in Pink</b> ..... No Charge<br>Part: DGN164245-PPK3 x2  | <b>Lemon Crush</b> .....No Charge<br>Part: DGN164245-PYL4 x2     |
| <b>Grey Street</b> .....No Charge<br>Part: DGN164245-PSL3 x2  | <b>Evolution Orange</b> .....No Charge<br>Part: DGN164245-POR7 x2 | <b>Green Machine</b> ..... No Charge<br>Part: DGN164245-PGN11 x2 |

Drive wheels section continued on following page.

**6. DRIVE WHEELS CONTINUED**

Military Green ..... No Charge Part: DGN164245-PGN8 x2	Military Digital Camo ..... \$50 ea. Part: DGN164245-WCM7 x2	Carbon Fiber ..... \$50 ea. Part: DGN164245-WCB1 x2
Sugar Plum ..... No Charge Part: DGN164245-PPR7 x2	Faded Glory ..... \$50 ea. Part: DGN164245-WFL2 x2	

**7. BATTERIES**

**Optional.** The powerchair requires 2 batteries and will be installed if ordered.

NF-22, Gel Cell, Interceptor 2 Battery ..... \$300 ea. Part: BAT2004376 x2. HCPCS: E2361	NF-22, Gel Cell, MK Battery ..... \$400 ea. Part: BATGEL1008 x2. HCPCS: E2361
---	--

**8. PATIENT'S WEIGHT**

**REQUIRED.** The patient's weight is required to quote and order the power chair.

Up to 200 lb User Weight Part: PTOINDV3837	201 lb-300 lb User Weight Part: PTOINDV3838
---	--

**9. ELECTRONICS**

**REQUIRED.** Select a type of drive control from the following sections. The No Joystick option is used to omit the joystick and utilize another control, such as a specialty control at the end of this form or found on the [Alternative Controls and Electronics Order Form](#).

**9.1. Expandable Controls - Quantum Intelligent Systems™**

Select expandable controls or omit selection if another controller type is desired.

**2. Q-Logic 3 Joystick**

Q-Logic 3 EX Joystick ..... \$900 Part: CTL167724	Q-Logic 3 EX Joystick for Full Lights ..... \$900 Part: CTL167725
Q-Logic 3e Joystick ..... \$620 Part: CTL1910600 Can operate full light package. Not available with alternative or specialty controls.	No Joystick, Setup for Q-Logic Part: XRF000010



**9.2. Q-Logic 3 EX Joystick Armor**

**REQUIRED with Q-Logic 3 EX Joystick selections above. Not available with Q-Logic 3e Joystick.** Reference colors shown on page 2.

Red Corvette Part: DGN124336-PRD1	Grey Street Part: DGN124336-PSL3	Pretty in Pink Part: DGN124336-PPK3	Blue Jean Part: DGN124336-PBU11	Sugar Plum Part: DGN124336-PPR7
Cherry Bomb Part: DGN124336-PRD4	White Satin Part: DGN124336-PWH2	Evolution Orange Part: DGN124336-POR7	Lemon Crush Part: DGN124336-PYL4	Military Digital Camo Part: DGN124336-WCM7
Back in Black Part: DGN124336-PBK1	Deep Purple Part: DGN124336-PPR5	Midnight Blue Part: DGN124336-PBU10	Green Machine Part: DGN124336-PGN11	Faded Glory Part: DGN124336-WFL2
Night Train* See Black, Soft Touch	Raspberry Beret Part: DGN124336-PPK4	Electric Blue Part: DGN124336-PBU9	Military Green Part: DGN124336-PGN8	Carbon Fiber Part: DGN124336-WCB1

**10. JOYSTICK HANDLES**

**Optional.** The optional handle will be mounted to the joystick and the standard handle will be shipped with the unit.

Mushroom Handle (B) ..... \$120 Part: ACCASMB1860.	Extended I-Handle (G) ..... \$140 Part: ACCASMB1857.
Chin Cup (C) ..... \$85 Part: ACCASMB186190	Bodypoint 3" Goal Post Handle ..... \$180 Part: PLSKNCP1057.
T-Bar (D) ..... \$160 Part: ACCASMB1858.	Bodypoint 4" Goal Post Handle ..... \$180 Part: PLSKNCP1058.
Medium Ball Handle (E) ..... \$70 Part: ACCASMB1859.	
Soft Ball Handle (F) ..... \$120 Part: ACCASMB1856.	



## 11. JOYSTICK MOUNTING BRACKETS

**Optional.** Select a mount to be used with a joystick selected in previous section.

### 11.1a. Fixed, In-line Joystick Bracket


Fixed, In-line, Bracket for Q-Logic 3 ..... Standard  
 Part: FRMASMB10617



### 11.1b. Swing-Away (Retractable) Joystick Bracket

Allows the joystick to retract to the side of the armrest while maintaining the same joystick orientation.

Swing-away, Q-Logic 3 ..... \$265  
 Part: ACC1908434. HCPCS: E1028



### 11.1c. Basic Swing Away Joystick Mount

Basic Swing Away, Q-Logic, Left Side ..... \$265  
 Part: ACC2201041. HCPCS: E1028

Basic Swing Away, Q-Logic, Right Side ..... \$265  
 Part: ACC2201042. HCPCS: E1028

### 11.1d. Specialty Joystick Mounts

Retract4 Mount ..... \$450  
 Part: XRF000003. HCPCS: E1028

Full 180° of movement that translates into multi-position use, straight, right offset, left offset, or fold away. Joystick orientation varies as bracket is adjusted.

Multi-Axis Mount ..... \$100  
 Part: XRF000004  
 Allows height adjustment of joystick.

Stealth Swing-Away/Flip Down Bracket ..... \$300  
 Part: XRF000005. HCPCS: E1028

Allows the joystick to pivot to side of armrest, or flip-down if bracket is adjusted.



### 11.2. Joystick Mount Location

**Required except with basic swing away joystick mount.**

Left Side Mount  
 Part: PTOINDV1225

Right Side Mount  
 Part: PTOINDV1226

## 12. TRU-BALANCE® 4 SEATING

**REQUIRED.** Maximum seat size is 22"W x 22"D.

### 12.1. TRU-Balance 4 Configuration

**REQUIRED.**

✓ TRU-Balance 4 Power Tilt and Recline ..... \$14,995  
 Part: PTOINDV3954, ADP: WP0001925  
 50° of posterior tilt and 165° of recline. Contains positional feedback.

### 12.2. 12" Power Lift

**REQUIRED.**

✓ Power Lift with iLevel ..... \$2,495  
 Part: PTOINDV3956  
 12" scissor lift with positional feedback.

### 12.3. Anterior Tilt

Chest support is required with anterior tilt 20 and 30 degree functions. See chest support section. Omitting a chest support selection requires the provider to facilitate aftermarket chest supports that best fits consumer needs.

#### 1. Anterior Tilt Angle

**REQUIRED.**

10° Anterior Tilt Angle ..... Standard  
 Part: MEC2009356

20° Anterior Tilt Angle ..... \$2,695  
 Part: MEC2009399  
 Requires chest support and knee block.

30° Anterior Tilt Angle ..... \$3,195  
 Part: MEC2009400

Requires chest support and knee block.

#### 2. Disable Anterior Tilt

**Optional.**

Disable Anterior Tilt ..... No Charge  
 Part: PTOINDV3980  
 Only available with 10° Anterior Tilt Angle. Anterior tilt is restricted with programming.

### 12.4. Knee Block Assembly

**REQUIRED** with 20° or 30° anterior tilt. **Optional** with 10° anterior tilt.

Anterior Knee Block Assembly (ST-KBATB4ATC-1-Q) ..... \$700  
 Part: ST-KBATB4ATC-1-Q

### 12.5. Occupied Transit Kit

**Optional.** The power base is standard with Unoccupied Transit Loops. 60" Sure-Lok Lap Belt is required with Occupied Transit Option.

WC-19 Occupied Transit Option ..... \$400  
 Part: KIT152047

To determine if your chair is WC-19 compliant (occupied), please refer to this website: <http://wc-transportation-safety.umtri.umich.edu/crash-tested-product-lists/wheelchairs>

TRU-Balance 4 Seating section continued on following pages.

**12. TRU-BALANCE 4 SEATING CONTINUED**

**12.6. Seat Width**

**REQUIRED with seat selection.** Maximum seat width is 22". Part number is SETWDTH plus number following # below.

16" W (#1006)	17" W (#1007)	18" W (#1008)	19" W (#1009)	20" W (#1010)	21" W (#1011)	22" W (#1012)
------------------	------------------	------------------	------------------	------------------	------------------	------------------

**12.7. Seat Width Package**

**REQUIRED with seat selection.**

16"-20" Width Package  
Part: XRF000008

20"-24" Width Package  
Part: XRF000158

**12.8. Seat Depth**

**REQUIRED with seat selection.** Maximum seat depth is 22". Part number is SETDPHT plus number following # below.

14" D (#1004)	15" D (#1005)	16" D (#1006)	17" D (#1007)	18" D (#1008)	19" D (#1009)	20" D (#1010)	21" D (#1011)	22" D (#1012)
------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------

**12.9. Seat Depth Package**

**REQUIRED with seat selection.**

16"-22" Depth Package  
Part: XRF000309

**12.10. Seat to Floor Height**

**REQUIRED.**

Lowest Seat to Floor  
Part: PTOINDV3558

Low Seat to Floor  
Part: PTOINDV3976

Middle Seat to Floor  
Part: PTOINDV3559

Highest Seat to Floor  
Part: PTOINDV3560

Lowest	Low	Med	High
Q6 Edge 3			
17.50"	18.25"	18.75"	19.25"
TRU-Balance 4 Seat to Floor Heights			


**12.11. Matching Seat Cushion**

**Optional. Specific cushion/component sizes available, part numbers, weight capacities, MSRP, and special order cushions can be found on the respective order forms.** This section allows a quick selection of a single cushion that matches the seat pan width and depth. See specific order form for non-matching sizes, alternate covers and liners, or additional covers and liners.

**Stealth Simplicity G Cushion**

General use cushion with Stealth Tek 2, 2-way stretch, breathable cover.  
Part series: *ST-SIMGXXYY*

Simplicity G, Matching Seat Size  
Part: XRF000111



**Stealth Spectrum Foam SPP Cushion**

Skin protection and positioning with Coolcore, 4-way stretch, breathable cover. Includes urethane liner.  
Part series: *ST-SPCFSPXXYY*


Spectrum Foam SPP, Matching Seat Size  
Part: XRF000112



**Stealth Solution SPP Cushion**

Skin protection and positioning with Coolcore, 4-way stretch, breathable cover. Includes urethane liner.  
Part series: *ST-SOLSPPXXYY*


Solution SPP, Matching Seat Size  
Part: XRF000113



**Stealth Glacial SP Cushion**

Skin protection with Coolcore, 4-way stretch, breathable cover.  
Part series: *ST-GLSPXXYY*

Glacial SP, Matching Seat Size  
Part: XRF000114



**Stealth Zen SP Cushion**

Skin protection with Coolcore, 4-way stretch, breathable cover.  
Part series: *ST-ZENSPXXYY*


Zen SP, Matching Seat Size  
Part: XRF000115



**Stealth TRU-Comfort 2 SPP Cushion**

Skin protection and positioning with Coolcore, 4-way stretch, breathable cover.  
Part series: *ST-TC2SPPXXYY*

TRU-Comfort 2 SPP, Matching Seat Size  
Part: XRF000116



**Stealth Spectrum Gel SPP Cushion**

Skin protection and positioning with Coolcore, 4-way stretch, breathable cover. Includes urethane liner.  
Part series: *ST-SPCGSPXXYY*

Spectrum Gel SPP, Matching Seat Size  
Part: XRF000117



**Stealth Essence SPP Cushion**

Skin protection and positioning with Coolcore, 4-way stretch, breathable cover.  
Part series: *ST-ESSSPXXYY*

Essence SPP, No Wedge, Matching Seat Size  
Part: XRF000118

Essence SPP, Small Wedge, Matching Seat Size  
Part: XRF000119

Essence SPP, Medium Wedge, Matching Seat Size  
Part: XRF000120

Essence SPP, Large Wedge, Matching Seat Size  
Part: XRF000121

Essence Plus SPP, Large Wedge, Matching Seat Size  
Part: XRF000121



**Stealth Premiere P Cushion**

Positioning with Coolcore, 4-way stretch, breathable cover.  
Part series: *ST-PREMPXXYY*

Premiere P, Matching Seat Size  
Part: XRF000122



**12.12. Back Cushion Configuration**

**REQUIRED.** Select a configuration type and applicable dimensions. See the Sport, Q-Back, and TRU-Comfort back availability chart for stock order back size ranges. Other size backs may be available with special order and extended lead time.

**1a. Sport Back/Omit/After-Market Interface**

Sport Back..... No Charge  
Part: XRF000017  
Complete following width package and size.

Omit Back Pan/Cushion..... No Charge  
Part: PTOINDV1287

Select a following width package.

After-market Back Interface ..... No Charge  
Part: *ST-SSAPT4-Q*. Select a set width in following section, range of 12"-21". Overall seat depth may need to be adjusted to accommodate certain backs. May interfere with gel/waterfall armpads and arm

TRU-Balance 4 Seating section continued on following pages.

**12. TRU-BALANCE 4 SEATING CONTINUED**

**1b. TRU-Comfort Back**

TRU-Comfort Back (Static/Flip Forward), Up to 15" Width..... \$835  
 Part: SET123622. HCPCS: E2620  
 Complete following width package and size.

TRU-Comfort Back (Static/Flip Forward), 16" to 20" Width... \$1,035  
 Part: TC4V10W21W. HCPCS: E2620  
 Complete following width package and size.

TRU-Comfort Back (Static/Flip Forward), 22" Width .....\$1,235  
 Part: TC4V22W32W. HCPCS: E2621  
 Complete following width package and size.

**1c. Q-Back 2.0, AL Series, Active Contour**

Select a back type and the specific overall back height in subsequent section.

Active Contour, 14"Wx19-20"H Tilt/22"-23" H Recline ..... \$985  
 Part: QB2AL14AS. HCPCS: E2613

Active Contour, 14"Wx21-22"H Tilt/24"-25" H Recline ..... \$985  
 Part: QB2AL14AT. HCPCS: E2613

Active Contour, 16"Wx19-20"H Tilt/22"-23" H Recline ..... \$1,145  
 Part: QB2AL16AS. HCPCS: E2613

Active Contour, 16"Wx21-22"H Tilt/24"-25" H Recline ..... \$1,145  
 Part: QB2AL16AT. HCPCS: E2613

Active Contour, 16"Wx23-24"H Tilt/26"-27" H Recline ..... \$1,145  
 Part: QB2AL16AX. HCPCS: E2613

Active Contour, 16"Wx25-26"H Tilt/28"-29" H Recline ..... \$1,145  
 Part: QB2AL16AXX. HCPCS: E2613

Active Contour, 18"Wx19-20"H Tilt/22"-23" H Recline ..... \$1,145  
 Part: QB2AL18AS. HCPCS: E2613

Active Contour, 18"Wx21-22"H Tilt/24"-25" H Recline ..... \$1,145  
 Part: QB2AL18AT. HCPCS: E2613

Active Contour, 18"Wx23-24"H Tilt/26"-27" H Recline ..... \$1,145  
 Part: QB2AL18AX. HCPCS: E2613

Active Contour, 18"Wx25-26"H Tilt/28"-29" H Recline ..... \$1,145  
 Part: QB2AL18AXX. HCPCS: E2613

Active Contour, 20"Wx19-20"H Tilt/22"-23" H Recline ..... \$1,145  
 Part: QB2AL20AS. HCPCS: E2613

Active Contour, 20"Wx21-22"H Tilt/24"-25" H Recline ..... \$1,145  
 Part: QB2AL20AT. HCPCS: E2613

Active Contour, 20"Wx23-24"H Tilt/26"-27" H Recline ..... \$1,145  
 Part: QB2AL20AX. HCPCS: E2613

Active Contour, 20"Wx25-26"H Tilt/28"-29" H Recline ..... \$1,145  
 Part: QB2AL20AXX. HCPCS: E2613

**1d. Q-Back 2.0, AL Series, Deep Contour**

Select a back type and the specific overall back height in subsequent section.

Deep Contour, 14"Wx19-20"H Tilt/22"-23" H Recline ..... \$1,145  
 Part: QB2AL14DS. HCPCS: E2620

Deep Contour, 14"Wx21-22"H Tilt/24"-25" H Recline ..... \$1,145  
 Part: QB2AL14DT. HCPCS: E2620

Deep Contour, 16"Wx19-20"H Tilt/22"-23" H Recline ..... \$1,325  
 Part: QB2AL16DS. HCPCS: E2620

Deep Contour, 16"Wx21-22"H Tilt/24"-25" H Recline ..... \$1,325  
 Part: QB2AL16DT. HCPCS: E2620

Deep Contour, 16"Wx23-24"H Tilt/26"-27" H Recline ..... \$1,325  
 Part: QB2AL16DX. HCPCS: E2620

Deep Contour, 16"Wx25-26"H Tilt/28"-29" H Recline ..... \$1,325  
 Part: QB2AL16DXX. HCPCS: E2620

Deep Contour, 18"Wx19-20"H Tilt/22"-23" H Recline ..... \$1,325  
 Part: QB2AL18DS. HCPCS: E2620

Deep Contour, 18"Wx21-22"H Tilt/24"-25" H Recline ..... \$1,325  
 Part: QB2AL18DT. HCPCS: E2620

Deep Contour, 18"Wx23-24"H Tilt/26"-27" H Recline ..... \$1,325  
 Part: QB2AL18DX. HCPCS: E2620

Deep Contour, 18"Wx25-26"H Tilt/28"-29" H Recline ..... \$1,325  
 Part: QB2AL18DXX. HCPCS: E2620

Deep Contour, 20"Wx19-20"H Tilt/22"-23" H Recline ..... \$1,325  
 Part: QB2AL20DS. HCPCS: E2620

Deep Contour, 20"Wx21-22"H Tilt/24"-25" H Recline ..... \$1,325  
 Part: QB2AL20DT. HCPCS: E2620

Deep Contour, 20"Wx23-24"H Tilt/26"-27" H Recline ..... \$1,325  
 Part: QB2AL20DX. HCPCS: E2620

Deep Contour, 20"Wx25-26"H Tilt/28"-29" H Recline ..... \$1,325  
 Part: QB2AL20DXX. HCPCS: E2620

**1e. Q-Back 2.0, AL Series, Extra Deep Contour**

Select a back type and the specific overall back height in subsequent section.

Extra Deep Contour, 14"Wx19-20"H Tilt/22"-23" H Recline ..... \$1,225  
 Part: QB2AL14ES. HCPCS: E2620

Extra Deep Contour, 14"Wx21-22"H Tilt/24"-25" H Recline ..... \$1,225  
 Part: QB2AL14ET. HCPCS: E2620

Extra Deep Contour, 16"Wx19-20"H Tilt/22"-23" H Recline ..... \$1,410  
 Part: QB2AL16ES. HCPCS: E2620

Extra Deep Contour, 16"Wx21-22"H Tilt/24"-25" H Recline ..... \$1,410  
 Part: QB2AL16ET. HCPCS: E2620

Extra Deep Contour, 16"Wx23-24"H Tilt/26"-27" H Recline ..... \$1,410  
 Part: QB2AL16EX. HCPCS: E2620

Extra Deep Contour, 16"Wx25-26"H Tilt/28"-29" H Recline ..... \$1,410  
 Part: QB2AL16EXX. HCPCS: E2620

Extra Deep Contour, 18"Wx19-20"H Tilt/22"-23" H Recline ..... \$1,410  
 Part: QB2AL18ES. HCPCS: E2620

Extra Deep Contour, 18"Wx21-22"H Tilt/24"-25" H Recline ..... \$1,410  
 Part: QB2AL18ET. HCPCS: E2620

Extra Deep Contour, 18"Wx23-24"H Tilt/26"-27" H Recline ..... \$1,410  
 Part: QB2AL18EX. HCPCS: E2620

Extra Deep Contour, 18"Wx25-26"H Tilt/28"-29" H Recline ..... \$1,410  
 Part: QB2AL18EXX. HCPCS: E2620

Extra Deep Contour, 20"Wx19-20"H Tilt/22"-23" H Recline ..... \$1,410  
 Part: QB2AL20ES. HCPCS: E2620

Extra Deep Contour, 20"Wx21-22"H Tilt/24"-25" H Recline ..... \$1,410  
 Part: QB2AL20ET. HCPCS: E2620

Extra Deep Contour, 20"Wx23-24"H Tilt/26"-27" H Recline ..... \$1,410  
 Part: QB2AL20EX. HCPCS: E2620

Extra Deep Contour, 20"Wx25-26"H Tilt/28"-29" H Recline ..... \$1,410  
 Part: QB2AL20EXX. HCPCS: E2620

TRU-Balance 4 Seating section continued on following page.

**12. TRU-BALANCE 4 SEATING CONTINUED**

**12.13. TRU-Balance 4 Back Shroud**

**REQUIRED.** The back shroud is a feature of the TRU-Balance 4 seating system and is available in matching shroud colors. The back shroud is not available when Omit Back Pan/Cushion is selected. Reference colors shown on page 2.

<b>Red Corvette</b> Part: XRF000019	<b>Grey Street</b> Part: XRF000023	<b>Pretty in Pink</b> Part: XRF000027	<b>Blue Jean</b> Part: XRF000031	<b>Sugar Plum</b> Part: XRF000035
<b>Cherry Bomb</b> Part: XRF000020	<b>White Satin</b> Part: XRF000024	<b>Evolution Orange</b> Part: XRF000028	<b>Lemon Crush</b> Part: XRF000032	<b>Faded Glory</b> Part: DGN1714085-WFL2
<b>Back in Black</b> Part: XRF000021	<b>Deep Purple</b> Part: XRF000025	<b>Midnight Blue</b> Part: XRF000029	<b>Green Machine</b> Part: XRF000033	<b>Digital Camo</b> Part: DGN1714085-WCM7
<b>Night Train</b> Part: XRF000022	<b>Raspberry Beret</b> Part: XRF000026	<b>Electric Blue</b> Part: XRF000030	<b>Military Green</b> Part: XRF000034	

**12.14. Lateral Positioning Wedge**

**Optional.** Available with Sport and TRU-Comfort backs. Not available on power recline back 19"-20"H with limited adjustability on 21"H.

Lateral Positioning Wedge ..... \$30  
Part: SET123465 x2

**12.15. Back Width Package**

**REQUIRED** except when Omit Back Pan/Cushion is selected. Available with +/-4" of seat width.

Match Seat Width Selections  
Part: PTOINDV3468

14"W-16"W Package  
Part: XRF000037

16"W-20"W Package  
Part: XRF000038

18"W-22"W Package  
Part: XRF000039

- S Sport Back
- T TRU-Comfort Back
- ST Sport or TRU-Comfort Back
- SQT Sport, Q-Back, or TRU-Comfort Back

		Back Width												
		12"	14"	15"	16"	17"	18"	19"	20"	21"	22"			
Overall Height	17"	T	-	-	-	-	-	-	-	-	-	-	-	-
	18"	T	-	-	-	-	-	-	-	-	-	-	-	-
	19"	S	ST	ST	ST	ST	ST	ST	ST	ST	ST	ST	ST	ST
	20"	S	ST	ST	ST	ST	ST	ST	ST	ST	ST	ST	ST	ST
	21"	-	ST	ST	SQT	ST	ST	ST	ST	ST	ST	ST	ST	ST
	22"	-	ST	ST	SQT	ST	SQT	ST	SQT	ST	SQT	ST	ST	ST
	23"	-	ST	ST	SQT	ST	SQT	ST	SQT	ST	SQT	ST	ST	ST
	24"	-	ST	ST	SQT	ST	SQT	ST	SQT	ST	SQT	ST	ST	ST
	25"	-	ST	ST	SQT	ST	SQT	ST	SQT	ST	SQT	ST	ST	ST
	26"	-	ST	ST	ST	ST	ST	ST	ST	ST	ST	ST	ST	ST
	27"	-	ST	ST	ST	ST	ST	ST	ST	ST	ST	ST	ST	ST
	28"	-	ST	ST	ST	ST	ST	ST	ST	ST	ST	ST	ST	ST
	29"	-	ST	ST	ST	ST	ST	ST	ST	ST	ST	ST	ST	ST
	30"	-	ST	ST	ST	ST	ST	ST	ST	ST	ST	ST	ST	ST
31"	-	ST	ST	ST	ST	ST	ST	ST	ST	ST	ST	ST	ST	
32"	-	ST	ST	ST	ST	ST	ST	ST	ST	ST	ST	ST	ST	

Sport, Q-Back, and TRU-Comfort Back Availability Chart

**12.16. Back Width**

**REQUIRED** except when Omit Back Pan/Cushion is selected. Available with +/-4" of seat width.

<b>12" Back Width</b> Part: PTOINDV3562	<b>16" Back Width</b> Part: PTOINDV3564	<b>18" Back Width</b> Part: PTOINDV3566	<b>20" Back Width</b> Part: PTOINDV3568	<b>22" Back Width</b> Part: PTOINDV3570
<b>14" Back Width</b> Part: PTOINDV3563	<b>17" Back Width</b> Part: PTOINDV3565	<b>19" Back Width</b> Part: PTOINDV3567	<b>21" Back Width</b> Part: PTOINDV3569	

**12.17. Overall Back Height**

**Power Recline available height range: 19"-32"H. Odd heights include 6" gap and even heights include 7" gap between seat pan and bottom of back cushion. See the availability chart for stock order back size ranges.**

<b>19" Height</b> Part: SETHGHT1032	<b>22" Height</b> Part: SETHGHT1035	<b>25" Height</b> Part: SETHGHT1038	<b>28" Height</b> Part: SETHGHT1041	<b>31" Height</b> Part: SETHGHT1102
<b>20" Height</b> Part: SETHGHT1033	<b>23" Height</b> Part: SETHGHT1036	<b>26" Height</b> Part: SETHGHT1039	<b>29" Height</b> Part: SETHGHT1042	<b>32" Height</b> Part: SETHGHT1103
<b>21" Height</b> Part: SETHGHT1034	<b>24" Height</b> Part: SETHGHT1037	<b>27" Height</b> Part: SETHGHT1040	<b>30" Height</b> Part: SETHGHT1043	

**12.18. Optional Back Height Extension**

**Optional.** Adds 2" of seat back pan height. 25"H selection includes a back height extension.

Optional Back Height Extension ..... \$48  
Part: XRF000041

**13. ARMRESTS**

**REQUIRED.** Select the desired type of armrest below. Please read notes to ensure seating compatibility. Be sure to also select an armrest size and armpad type.

**13.1. Flip-Back, 8"-16" H Armrests**

Choose a left and right armrest and set height in the following section. Not available with waterfall armpads. Not available with Back Cane Option.

**L R** Flip-Back, 8"-16"H Armrests ..... \$145 ea.  
Left part: ARM1806662. Right part: ARM1806663

**Armrest Support**

**Optional.** Available up to 2" difference between seat and back width via endomorph or mesomorph configuration. Not available less than 16" width seat or back.

**L R** Flip-Back Armrest Downpost ..... \$115 ea.  
Left part: XRF000265. Right part: XRF000266

Armrests section continued on following page.

**13. ARMRESTS CONTINUED**

**13.2. Armrest Height**

**REQUIRED.** Set left and right armrest height. Note the available height range in the option description.

**1. Left Armrest Height**

8" H PTOINDV3535	10" H PTOINDV3537	12" H PTOINDV3539	14" H PTOINDV3541	16" H PTOINDV3543
9" H PTOINDV3536	11" H PTOINDV3538	13" H PTOINDV3540	15" H PTOINDV3542	

**2. Right Armrest Height**

8" H PTOINDV3450	10" H PTOINDV3452	12" H PTOINDV3454	14" H PTOINDV3456	16" H PTOINDV3458
9" H PTOINDV3451	11" H PTOINDV3453	13" H PTOINDV3455	15" H PTOINDV3457	

**13.3. Armrest Size**

**REQUIRED.** Select a left (L) and right (R) armrest size. This also affects armpad selection.

<b>L R</b> 10" Length L part: XRF000057. R part: XRF000123	<b>L R</b> 14" Length L part: XRF000058. R part: XRF000124	<b>L R</b> 17" Length L part: XRF000059. R part: XRF000125
---	---	---

**14. ARMPADS**

**REQUIRED.**

**14.1a. Aftermarket Armpad**

Not available with Chest Bar. Select one below if aftermarket armpads or handpads are required.

Interface for Aftermarket Pad ..... No Charge  
Part: XRF000308

**14.1b. Foam Armpads**

**L R** Straight, Armpad ..... Standard  
Left part: XRF000061. Right part: XRF000126

**14.1c. Gel Armpads**

Sold only in pairs. Armpads are 3" wide. Manufacturer part number for full length is 314G and desk length is 310G.

Falcon Gel, Full Length Left and Right ..... \$300 Part: ARM1806807	Falcon Gel, Desk Length Left, Full Length Right..... \$250 Part: ARM1806810
Falcon Gel, Full Length Left, Desk Length Right..... \$250 Part: ARM1806813	Falcon Gel, Desk Length Left and Right ..... \$200 Part: ARM1806806

**15. JOYSTICK BRACKETS**

**REQUIRED.**

Standard Fixed Receiver ..... Standard Part: XRF000063	Stealth Height Adjustable Receiver.....\$165 Part: XRF000064
---	---

**16. HEADRESTS**

**Optional.** Select a headrest pad and bracket.

**16.1a. Headrests, Comfort Plus with Cool Core**

**1. Headrest Pads**

Stealth 8" Comfort Plus, Cool Core (CP08CCFPM) ..... \$190 Part: ST-CP08CCFPM-Q. E0955	Stealth 14" Comfort Plus, Cool Core (CP14CCFPM) ..... \$190 Part: ST-CP14CCFPM-Q. E0955
Stealth 10" Comfort Plus, Cool Core (CP10CCFPM)..... \$190 Part: ST-CP10CCFPM-Q. E0955	Stealth 6", Comfort Plus, Cool Core (CP06CCFPM) ..... \$190 Part: ST-CP06CCFPM-Q. E0955

**2. Headrest Brackets**

Comfort Plus Headrest Mounting Bracket, UniLink.....\$255  
Part: ST-UHB2D33200-Q. E1028

**16.1b. Headrests, Comfort Plus**

**1. Headrest Pads**

Stealth 8" Comfort Plus (CP450)..... \$190 Part: POS148484. HCPCS: E0955	Stealth 14" Comfort Plus (CP150) ..... \$190 Part: INDPART2327. HCPCS: E0955
Stealth 10" Comfort Plus (CP250)..... \$190 Part: SETHEAD1057. HCPCS: E0955	Pediatric Comfort Plus (6"Wx4.5"H) (CPN350) ..... \$190 Part: SETHEAD1069. HCPCS: E0955

**2. Headrest Brackets**

Multi-Axis Removable Headrest Bracket (HMO475P) ..... \$200  
Part: HDWASMB1084. HCPCS: E1028



**17. LAP BELTS**

**REQUIRED.** Select a lap belt below.

**17.1a. Standard Lap Belts**

- 38" Lap Belt Part: ACCBELT1006     50" Lap Belt Part: ACCBELT1000     60" Lap Belt Part: ACCBELT1002     70" Lap Belt Part: ACCBELT1003     80" Lap Belt Part: ACCBELT1004     90" Lap Belt Part: ACCBELT1001

**17.1b. Optional Lap Belts**

- Stealth Positioning, 60" Single Pull, Push Button, Medium..... \$65  
(7.5") Padded, 1.5" Belt (IFPB5213)  
Part: ST-IFPB5213-Q
- Stealth Positioning, 60" Single Pull, Push Button, Extra Large, Padded (11.5"), 2" Belt (IFPB5315)  
Part: ST-IFPB5315-Q
- Stealth Positioning, 60" Single Pull, Push Button, Large ..... \$65  
Padded (9.5"), 2" Belt (IFPB5314)  
Part: ST-IFPB5314-Q
- 60" Sure-Lok Lap Belt (Non-retractable)..... \$200  
Part: ACCBELT1016  
Required with ISO 7176-19 Compliant Transit Option.

**18. LEGRESTS**

All efforts will be made to accommodate legrest length selection, but minor changes may need to be made to meet shipping guidelines. Some adjustments may be needed upon receipt of the unit.

**18.1. Center Mount Foot Platform**

**REQUIRED.**

**1. Foot Platform Assembly**

- ✓ Dual Actuator, Power Articulating Foot Platform..... \$4,295  
Part: RIG2009246. HCPCS: E1012+K0040, ADP: WPO001926  
Positional feedback with independent control of articulation.

**2. Lower Extension and Angle Adjustable Footplates**

Medium length extension range is 11"-17" and long extension range is 14"-20". See compatibility matrix for possible interference.

- Medium Extension, 5"x6" Footplates Part: RIG150206    Long Extension, 5"x6" Footplates Part: RIG151505
- Medium Extension, 5"x8" Footplates Part: RIG150203    Long Extension, 5"x8" Footplates Part: RIG151506
- Medium Extension, 5"x11" Footplates Part: RIG150212    Long Extension, 5"x11" Footplates Part: RIG151508

Legrest Compatibility Matrix Articulating Foot Platform TRU-Balance 4	Seat to Floor Option			
	Lowest	Low	Middle	High
Medium Length Extension	11"-15"	11"-16"	11"-16.5"	11"-17"
Long Length Extension	14"-15"	14"-16"	14"-16.5"	14"-17"

**3a. Center Mount Calf Pads**

Calf pads are standard with a power articulating foot platform.

- Small Calf Pads, AFP Only, 6"Wx6"H (pair) ..... Standard Part: RIGASMB7120018
- Large Calf Pads, AFP Only, 8" taper to 6"Wx8"H (pair) ..... Standard Part: RIGASMB7120019

**4. Center Mount Modular Footplate Extension Components**

Optional. Assembly required. May require seating configuration adjustments to avoid front caster interference.

- Width Extension for 8"D Footplates ..... \$55 Part: RIG152517
- Width Extension for 11"D Footplates ..... \$55 Part: RIG152518
- 2" Bell-Shaped Length Extension ..... \$55 Part: RIG152519  
Not compatible with width extensions. Can be used with 6" outside edge with 11" footplates.
- 4" Outside Edge (pair) ..... \$55 Part: RIG152510  
Only available with 6"D footplates and 8"D footplate extension.
- 6" Outside Edge (pair) ..... \$55 Part: RIG152516  
Available with 8"D and 11"D footplates and 11"D footplate extension.

**5. Center Mount Additional Accessories**

- Small Reversed Dartex Calf Pad Cover (pair) ..... \$150 Part: RIG133717
- Large Reversed Dartex Calf Pad Cover (pair) ..... \$150 Part: RIG133719
- Footplate Connector ..... \$55 Part: RIGASMB7120023

**18.2. Legrest Accessories**

- Wheel Bumper ..... \$15 ea. Part: XRF000100

**18.3. Legrest Length**

**REQUIRED with matching left and right legrest configuration. Part number is PTOINDV plus number following # below.** All efforts will be made to accommodate legrest length selection, but minor changes may need to be made to meet shipping guidelines. Some adjustments may be needed upon receipt of the unit. Reference the section's compatibility matrix for legrest interference. Non-matching or individual legrest configurations lengths are not set by the manufacturer.

11" L (#3478)	12" L (#3480)	13" L (#3482)	14" L (#3484)	15" L (#3486)	16" L (#3488)	17" L (#3490)
11.5" L (#3479)	12.5" L (#3481)	13.5" L (#3483)	14.5" L (#3485)	15.5" L (#3487)	16.5" L (#3489)	

**19. POWER POSITIONING ELECTRONICS**

**REQUIRED.**

**19.1a. Q-Logic 3 Controls through Joystick**

✓ Multiple Actuator Through Controller and Interface..... \$2,175  
 Part: XRF000336. HCPCS: E2311, ADP: SECX15005

**19.1b. Q-Logic 3 Controls with iAccess Switch Control**

iAccess Flat Button.....\$575 Part: CTL144009 See next section for switch mounting.	iAccess Toggle Button.....\$575 Part: CTL144010 See next section for switch mounting.
---	---

**19.2. iAccess Mount**

**REQUIRED with iAccess options.**

**1. iAccess Mount Selection**

Joystick Mount ..... No Charge Part: MEC1907500	Attendant Mount ..... \$35 Part: MEC1906564
Under Arm Mount, Flip-Back Armrest..... No Charge Part: MEC1906531	Under Arm Mount, Flip-Back Armrest, plus iAccess Guard ..... \$35 Part: MEC2001392

**2. Under Arm Mount Side**

Left Side Under Arm Mount Part: PTOINDV3678	Right Side Under Arm Mount Part: PTOINDV3679
--	---

**20. ATTENDANT CONTROL**

**Optional.** Mounting bracket is included.

Attendant Control ..... No Charge  
Part: XRF000078. HCPCS: E2331

**21. SEAT MOUNTED LED LIGHTS**


**Optional.** Requires Q-Logic 3 Joystick for Full Light Package or Q-Logic 3e Joystick.

Seat Mounted, LED, Full Light Package .....\$770  
 Part: XRF000130  
 Not available with Base Mounted, LED, Full Light Package.

**22. ACCESSORIES**

**Optional.** Multiple accessories can be mounted along the various tracks but some interference may occur.

Accessory Bar ..... \$95 Part: ACC1909634 Horizontal rail attached to rear of seat. Standard with the Trilogy vent tray and oxygen holder.	Personal Item Hook..... \$20 Part: ACC125008
Backpack Holder (Pair) (A)..... \$42 Part: ACC147774	Privacy Bag..... \$80 Part: ACC1907570
Cane/Crutch Holder ..... \$100 Part: ACC152515. HCPCS: E2207	Push Handles (pair)..... \$140 Part: ACC2005027 Not available with oxygen holder or vent tray.
Clothing Guard (Pair)..... \$275 Part: ACC125007	Q-Logic EX/NE Programmer ..... \$625 Part: ELEASMB6976
Compact Vent Tray .....\$1,395 Part: ACC133889. HCPCS: E1029	Quantum Backpack ..... \$200 Part: ACC1905853
Cup Holder ..... \$20 Part: ACC125003	Rear View Mirror ..... \$35 Part: ACC155833
Flag Kit ..... \$125 Part: ACC1700443	Reflector Kit ..... \$35 Part: ACC124670
Flip-Back Armrest, Side Track Accessory Bracket .. \$25 Part: MEC1809178	Stealth Large Tablet Holder ..... \$565 Part: ACC141818
Front Mounted, Cell Phone/ Media Holder, ..... \$95 X-Grip (E) Part: ACC155835	Transfer Handle (D)..... \$152.50 ea. Part: MEC1809930 x2. HCPCS: E0705
Glove Box (C)..... \$85 Part: ACC125011	Trilogy Vent Tray ..... \$420 Part: ACC136791. HCPCS: E1029 Includes accessory bar, a credit will be applied PTOIN-DV3769. Contact custom for models other than 100 and 200.
Large, Black, Weather Cover ..... \$240 Part: ACCCOVR1021	VOCSN Vent Tray for Accessory Bar Mount..... \$1,000 Part: ACC2000166. HCPCS: E1029 Includes accessory bar, a credit will be applied PTOIN-DV3769.
Oxygen Holder..... \$110 Part: ACC125006. HCPCS: E2208 Includes accessory bar, a credit will be applied PTOIN-DV3769.	



**24. POSITIONING COMPONENTS**

Optional. More options are available on the [Positioning Components Order Form](#).

**24.1. Stealth Unilink Pelvic/Thigh Guide**

**1. UniLink Pelvic/Thigh Hardware**

UniLink, 5" Height Adjustable, ..... \$320 ea.  
 Removable Mounting Hardware  
 Part: ST-UTTP5-Q  
 QTY \_\_\_\_\_

UniLink, 7" Height Adjustable, ..... \$320 ea.  
 Removable Mounting Hardware  
 Part: ST-UTTP7-Q  
 QTY \_\_\_\_\_



**2. UniLink Pelvic/Thigh Pads**

UniLink, 4" x 4" Pad ..... \$125 ea.  
 Part: ST-UTP544056100-Q  
 QTY \_\_\_\_\_

UniLink, 4" x 12" Pad ..... \$165 ea.  
 Part: ST-UTP54D056100-Q  
 QTY \_\_\_\_\_

UniLink, 4" x 8" Pad..... \$135 ea.  
 Part: ST-UTP548056100-Q  
 QTY \_\_\_\_\_

**3. UniLink Pelvic/Thigh Guide Installation**

Install on Left  
 Part: PTOINDV3585

Install on Right  
 Part: PTOINDV3586

Do Not Install  
 Part: PTOINDV3719

**24.2. Stealth Thoracic Laterals**

Swing-away laterals may interfere with flip back armrests. These laterals are not available with the back cane option.

**1. Bracket**

Swing-away, UniLink, Redesign Back, pair (ULS4042000-Q).... \$590  
 Part: ST-ULS4042000-Q. E1028  
 Only available with Redesigned back assemblies.

Swing-away, Rail Mounts, pair (TWBTLTC) ..... \$590  
 Part: INDPART2677. E1028  
 Not available with Redesign back assemblies.

**2. Pads**

Curved Pads, 3"x5", pair (TWBL3X5)..... \$140  
 Part: INDPART2663. E0956

Flat Pads, 3"x5", pair (TWBLP3X5) ..... \$140  
 Part: INDPART2676. E0956

Curved Pads, 4"x6", pair (TWBL4X6)..... \$140  
 Part: INDPART2664. E0956

Flat Pads, 4"x6", pair (TWBLP4X6) ..... \$140  
 Part: INDPART2675. E0956

**3. Installation**

Install Laterals  
 Part: PTOINDV3499

Do Not Install  
 Part: PTOINDV3500

**24.3. Stealth Pelvic/Thigh Guides**

Pelvic/thigh guides may interfere with the joystick receiver block and down posts. Please provide quantity next to check box.

**1. Pelvic/Thigh Guide Bracket**

Flip Down Bracket (TWBTLADD-HW) ..... \$300 ea.  
 Part: INDPART2804. : E1028  
 QTY \_\_\_\_\_

Flip Down with LEFT Full Surface Contact ..... \$400 ea.  
 Feature (TWBADD-L & LPHW407)  
 Part: POS140428.

Removable Bracket (LP430)..... \$330 ea.  
 Part: INDPART2805. : E1028  
 QTY \_\_\_\_\_

Flip Down with RIGHT Full Surface Contact..... \$400 ea.  
 Feature (TWBADD-R & LPHW407)  
 Part: POS140429.

Removable Bracket with 1" Offset (LP430 & LHW-121)... \$430 ea.  
 Part: INDPART2854. : E1028  
 QTY \_\_\_\_\_

**2. Pelvic/Thigh Guide Adapter**

2" Height Adapter Plate (LHW-122).....\$40  
 Part: INDPART2821.  
 Height adapter is needed with TRU-Comfort seat cushions. Not available with the flip down bracket.

2 1/4" Adapter Plate (LHW-129) ..... \$40  
 Part: POS140430.  
 Only available with the flip down full surface bracket.

**3. Pelvic/Thigh Guide Pad**

Short Pad, 4"x4" (TWBADD-4X4)..... \$70 ea.  
 Part: INDPART2801.  
 QTY \_\_\_\_\_

Long Pad, 4"x12" (TWBADD-4X12) ..... \$70 ea.  
 Part: INDPART2803.  
 QTY \_\_\_\_\_

Medium Pad, 4"x8" (TWBADD-4X8)..... \$70 ea.  
 Part: INDPART2802.  
 QTY \_\_\_\_\_

**4. Pelvic/Thigh Guide Installation**

Install on Left  
 Part: PTOINDV3585

Install on Right  
 Part: PTOINDV3586

Do Not Install  
 Part: PTOINDV3719

Positioning components section continued on following page.

**24. POSITIONING COMPONENTS CONTINUED**

**24.4. Chest Bar Assembly**

Small Chest Bar, 14" Foam Pad (ST-CBATB4ATT-S-Q)..... \$300      Medium Chest Bar, 16" Foam Pad (ST-CBATB4ATT-M-Q)..... \$300  
 Part: ST-IFPS1103-Q      Part: ST-CBATB4ATT-M-Q

**24.5. Stealth Positioning Chest Harness**

Upper Classic, Medium (IFPS1103) ..... \$190 ea.      Upper Contour, Medium (IFPS2103)..... \$190 ea.  
 Part: ST-IFPS1103-Q      Part: ST-IFPS2103-Q  
 Upper Classic, Large (IFPS1104) ..... \$190 ea.      Upper Contour, Large (IFPS2104) ..... \$190 ea.  
 Part: ST-IFPS1104-Q      Part: ST-IFPS2104-Q



**24.6. Stealth Positioning Chest Straps**

Chest Strap, Medium w/ Belt Anchor ..... \$120 ea.      Chest Strap, Large w/ Belt Anchor ..... \$120 ea.  
 (IFCS4113)      (IFCS4114)  
 Part: ST-IFCS4113-Q      Part: ST-IFCS4114-Q



**24.7. Stealth Swing-Away Headrest Mounts**

These swing-away headrest mounts are only available with the Stealth Comfort Series headrests.

**1. Facial Lateral Hardware**


Facial Lateral for Comfort Series..... \$430 ea.      Facial Laterals for Comfort Series ..... \$430 ea.  
 Right Side (CPS9-R-HW)      Left and Right Side (CPS9-B-HW)  
 Part: INDPART2647.      Part: INDPART2649.

Facial Lateral for Comfort Series..... \$430 ea.  
 Left Side (CPS9-L-HW)  
 Part: INDPART2648.



**2. Facial Lateral Pads and Accessories**


1" Ball Kit for Facial Laterals (TWB-1BK) (D) ..... \$100  
 Part: POS150532  
 Required with facial laterals to mount to comfort headrest brackets.  
 Egg Switch Mount (ESM) (A) ..... \$70 ea.  
 Part: INDPART2815  
 Figure 8 Bracket with 13" Rod..... \$315 ea.  
 and Egg Switch Mount (B)  
 Part: ST-F8675-13-Q  
 Facial Pad, 4"x2" (900-Q)..... \$70 ea.  
 Part: ST-900-Q : K0108  
 Facial Pad, 6"x2" (920-Q)..... \$70 ea.  
 Part: ST-920-Q : K0108  
 Swing-Away Joystick Mount, 18" Rod (C) .... \$295 ea.  
 Part: ST-SUS18-Q  
 Allows mounting of facial pads and attaches to the headrest



**24.8. Stealth Elbow Support**

**1. Elbow Support Bracket**

Left, Swing-Away, Elbow Support ..... \$230 ea.  
 Bracket Tube Style Arms (ARMS-ES-L)  
 Part: POS124064.  
 Left, UniLink, Fixed, Elbow Support Bracket ..\$360 ea.  
 for Flip-Back, 8"-16"H Armrests (UETB4021FL-Q)  
 Part: ST-UETB4021FL-Q  
 Right, Swing-Away, Elbow Support..... \$230 ea.  
 Bracket Tube Style Arms (ARMS-ES-R)  
 Part: POS124066.  
 Right, UniLink, Fixed, Elbow Support Bracket . \$360 ea.  
 for Flip-Back, 8"-16"H Armrest (UETB4021FR-Q)  
 Part: ST-UETB4021FR-Q



**2. Elbow Support Pad**

3"Wx3"H Elbow Pad (ARMS-ES-3X3) ..... \$70 ea.  
 Part: POS124065. HCPCS: K0108  
 4"Wx4"H Elbow Pad (ARMS-ES-4X4)..... \$70 ea.  
 Part: POS124067. HCPCS: K0108

**24.9. Stealth Foot Boxes**

Stealth Small Foot Boxes, 5"Wx8"Dx7"H..... \$625  
 (SPLE107) (pair)  
 Part: INDPART2786.  
 Stealth Medium Foot Boxes, 6"Wx10"Dx9"H ..... \$625  
 (SPLE108) (pair)  
 Part: INDPART2787.  
 Stealth Large Foot Boxes, 7"Wx12"Dx11"H..... \$625  
 (SPLE109) (pair)  
 Part: INDPART2788.

**24.10. Stealth Padded Shoe Holder**

Stealth Small Padded Shoe Holder 3.75"Wx6.25"D ..... \$220  
 (SPLE111) (pair)  
 Part: INDPART2783.  
 Stealth Medium Padded Shoe Holder, 4"Wx8.25"D..... \$220  
 (SPLE112) (pair)  
 Part: INDPART2784.  
 Stealth X-Large Padded Shoe Holder, 5"Wx10.75"D..... \$220  
 (SPLE114) (pair)  
 Part: INDPART2785.

**25. ALTERNATIVE CONTROLS AND ELECTRONICS**

**Optional.** Select a type of specialty controls below. More options are available on the [Alternative Controls and Electronics order form](#).

**25.1. Enhanced Display**

**Q-Logic 3 Enhanced Display Kit**..... \$2,045  
 Part: CTL167726  
 Select a color option below. Includes enhanced display, flexible gooseneck and mounting bracket. Features: 3.5" LCD display; up to 5 profiles to control drive, seat, and auxiliary functions; Bluetooth connectivity for mobile devices and mouse emulation; integrated single switch scanner; IR controls.



Enhanced Display

**Q-Logic 3 Enhanced Display Shroud**

Required with Q-Logic 3 Enhanced Display selection above. Reference colors shown on page 2.

<b>Black, Soft Touch</b> Part: DGN137728	<b>White Satin</b> Part: DGN137729-PWH2	<b>Electric Blue</b> Part: DGN137729-PBU9	<b>Military Digital Camo</b> Part: DGN137729-WCM7
<b>Red Corvette</b> Part: DGN137729-PRD1	<b>Deep Purple</b> Part: DGN137729-PPR5	<b>Blue Jean</b> Part: DGN137729-PBU11	<b>Faded Glory</b> Part: DGN137729-WFL2
<b>Cherry Bomb</b> Part: DGN137729-PRD4	<b>Raspberry Beret</b> Part: DGN137729-PPK4	<b>Lemon Crush</b> Part: DGN137729-PYL4	<b>Carbon Fiber</b> Part: DGN137729-WCB1
<b>Back in Black</b> Part: DGN137729-PBK1	<b>Pretty in Pink</b> Part: DGN137729-PPK3	<b>Green Machine</b> Part: DGN137729-PGN11	
<b>Night Train</b> Part: DGN137729-PBK3	<b>Evolution Orange</b> Part: DGN137729-POR7	<b>Military Green</b> Part: DGN137729-PGN8	
<b>Grey Street</b> Part: DGN137729-PSL3	<b>Midnight Blue</b> Part: DGN137729-PBU10	<b>Sugar Plum</b> Part: DGN137729-PPR7	

**25.2. SCIM Module**

**Optional.**  
**SCIM Module for Q-Logic 3** ..... \$850  
 Part: CTL2001963  
 Allows use of joystick display for specialty devices. Available with Q-Logic 3 EX Joystick selection above. Up to 5 profiles to control drive, seat, and auxiliary functions, excluding sip-n-puff and single switch functions. A switch is required to operate the input device attached to the SCIM module.



SCIM Module

**25.3. Q-Logic 3 SAJ**

**Q-Logic 3 EX SAJ** ..... \$850  
 Part: CTL143980  
 Must be used in conjunction with Enhanced Display or Q-Logic 3 Joystick.



EX SAJ

**25.4. Head Array**

<b>Adult, i-Drive Tri-Array Series, 3-Switch</b> ..... \$5,855 Head Array (IDHBT200-1) Part: ST-IDHBT200-1-Q. HCPCS: E2330 Features bluetooth, head support, 3 non-proportional proximity sensors for drive control, an egg switch for reset mode change, and i-Drive smart box.	<b>Adult, i-Drive Pro Series Head Array</b> ..... \$5,855 (IDHBT300-1) Part: ST-IDHBT300-1-Q. HCPCS: E2330 Features bluetooth, head support with occipital pad and 2 swing-away spot pads for left and right side, 3 non-proportional proximity sensors, an egg switch for reset mode change, and i-Drive smart box.
---	---



Adult Pro-Series

**25.5. Sip and Puff Drive**

**1. Sip and Puff Module**

**Sip and Puff Module for Q-Logic 3** ..... \$1,695  
 Part: CTL168056. HCPCS: E2325  
 Includes module and mounting hardware.



Sip & Puff Module

**2. Sip and Puff Straw and Mount**

**The straw and mount do not include electrical components to operate the Sip&Puff with a head array. This component requires the selection of a Q-Logic Sip&Puff module or additional module from Stealth Products for operation.**

<b>Stealth Sip and Puff for Comfort Series, Left Mount (SP625-CL)</b> ..... \$410 Part: ST-SP625-CL-Q. HCPCS: E2326 Includes straw and mounting hardware. Use with Stealth Comfort Plus headrests.	<b>Stealth Sip and Puff for Stealth i-Drive, Head Array (SP625-U)</b> ..... \$410 Part: ST-SP625-U-Q. HCPCS: E2326 Includes straw and mounting hardware. Use with Stealth i-Drive Head Arrays.
<b>Stealth Sip and Puff for Comfort Series, Right Mount (SP622-CR)</b> ..... \$410 Part: ST-SP625-CR-Q. HCPCS: E2326 Includes straw and mounting hardware. Use with Stealth Comfort Plus headrests.	



Stealth Sip & Puff

Alternative controls and electronics section continued on following page.

**25. ALTERNATIVE CONTROLS AND ELECTRONICS CONTINUED**

**25.6. Egg Switches**

Red Egg Switch (ESRED)..... \$100 Part: SWTMCRO1077	Green Egg Switch (ESGRN) ..... \$100 Part: SWTMCRO1078
Blue Egg Switch (ESBLU) ..... \$100 Part: SWTMCRO1079	Yellow Egg Switch (ESYEL) ..... \$100 Part: SWTMCRO1080
Black Egg Switch (ESBLK)..... \$100 Part: SWTMCRO1081	



**25.7. Switch Adapter**

2 to 1 Switch Adapter..... \$135  
Part: DWR1361H037  
Connects 2 single switches to a stereo jack to act as a dual switch.

**25.8. Environmental Controls**

Q-Logic 3 ECU Module Kit..... \$1,780  
Part: CTL168057

**25.9. USB Mobile Device Charger**

USB Mobile Device Charger ..... \$185  
Part: ELE168069  
Can charge 1 device. A multiplier harness is required and must be added when a Static Seat or Power Tilt Only seating system is chosen. The selection of more than 3 modules requires an additional multiplier.



**26. TRU-BALANCE 3 MULTIPLIER HARNESS**

**REQUIRED.** Multiple controllers or power accessories may require an additional multiplier harness or special order. Reference the chart below to determine the appropriate selection. The rows refer to the total amount of arm mounted devices and the columns refer to the total amount of back mounted devices. The online order configurator will automatically make this selection under specialty controls.

Secondary Multiplier (SM)..... \$250  
Part: XRF000107

Arm Mounts	Joystick (Hand Control)	1 Arm	2 Arm	3 Arm	4+ Arm
	Enhanced Display	NA	SM	SM	SOF
	iAccess	NA	SM	SM	SOF
	USB Charger	SM	SM	SM	SOF
	SAJ	SOF	SOF	SOF	SOF

Power Positioning Chart

Back Mounts	Attendant Control	1 Arm	2 Arm	3 Arm	4+ Arm
	SCIM	NA	SM	SOF	SOF
	Sip n Puff	NA	SM	SOF	SOF
	ECU	SM	SM	SOF	SOF
		SOF	SOF	SOF	SOF

Static Seat Chart

NA No additional Multiplier  
SM Secondary Multiplier  
SOF Special Order Required

**Q-Logic 3 Multiplier Charts**