



RAZ-AT

Attendant Tilt Mobile Shower Commode Chairs



Bathing and Toileting Products for Life

RAZ-AT FEATURES

The RAZ-AT (Attendant Tilt) tilts up to 40° to provide pressure reduction and maximum comfort. The unique frame design utilizes lockable gas springs that help lift the seat to an upright position, making the RAZ-AT one of the easiest tilting chairs to operate. The easy-to-clean, stainless steel frame fits over most toilets and is designed with unobstructed side and front access for the user or attendant.

EXTENSIVE RANGE OF OPTIONS AVAILABLE



| SPECIFICATIONS | RAZ-AT (Attendant Tilt) |
|--|---|
| WEIGHT CAPACITY | 400 lbs |
| FRAME WIDTH (seat and back frame) | Standard: 18" |
| OVERALL WIDTH (with Standard Armpad) | 22.5" |
| OVERALL LENGTH (with footrest) | 39" |
| SEAT WIDTH (Commode Seat) | 18" |
| SEAT DEPTH (front of back post to front of seat) | 16" or 19" |
| CLEARANCE HEIGHT UNDER SEAT | 17" – 21" with pan brackets (in 1" increments) 17.75" – 21.75" without pan brackets (in 1" increments) |
| ARMREST HEIGHT ABOVE SEAT (20" Adjusta-Back) | 7.5" - 10.5" in ½" increments |
| PRODUCT WEIGHT (with foot supports) | 61 lbs |
| SEAT TO FLOOR HEIGHT | 21.25" - 25.25" in 1" increments |
| FRAME AND HARDWARE MATERIAL | Stainless Steel – Life Time Warranty Against Corrosion |

All measurements taken with standard molded seat and standard arm pad.

A 2-year limited warranty covers the frame against manufacturing defects. A 1-year limited warranty covers all other components, subject to normal wear and tear exclusions. Specifications and options are subject to change without notice.