

# VivaLIFT!®

POWER RECLINER COLLECTION

THE MOST COMFORTABLE RECLINER IN THE WORLD!

THE 4-MOTOR RECLINER



New! Ergonomically designed toggle remote with two programmable memory positions



Shown in Canyon Walnut Fabric

NEW FEATURES:

**HEAT**  
TOGGLE REMOTE  
TRUE-INFINITE  
POSITIONING

## RADIANCE

Whether you cozy up with a hot beverage and your favorite book or kick back for a relaxing evening of television, this power recliner is right for you! Experience the warmth of the **Radiance Collection** from VivaLift!® Power Recliners. With three different heating pads in the shoulders, lumbar and seat, you can adjust your heat level accordingly and soothe your tired muscles. Keep your smart phone close at hand while you charge it with the wireless charging station. And if the power goes out, know you have peace of mind with the lithium battery backup.



Scan here for more information on **VivaLift!®**



**PRIDE**®  
MOBILITY  
LIVE YOUR BEST®

# RADIANCE

**VivaLIFT!**  
POWER RECLINER COLLECTION

COLLECTION	RADIANCE	RADIANCE	RADIANCE	RADIANCE
MODEL	PLR-3955LT	PLR-3955M	PLR-3955PW	PLR-3955S
POSITION	Infinite-Position Trendelenburg Zero Gravity Chaise Lounger Relaxed Reclined Upward Lift	Infinite-Position Trendelenburg Zero Gravity Chaise Lounger Relaxed Reclined Upward Lift	Infinite-Position Trendelenburg Zero Gravity Chaise Lounger Relaxed Reclined Upward Lift	Infinite-Position Trendelenburg Zero Gravity Chaise Lounger Relaxed Reclined Upward Lift
SIZE	Large/Tall 5'10" and above	Medium 5'4"-5'9"	Petite Wide 5'3" and below	Small 5'3" and below
WEIGHT CAPACITY	400 lbs.	400 lbs.	400 lbs.	400 lbs.
SEAT-TO-FLOOR	21.5"	19.5"	19.5"	18"
SEAT DEPTH	22"	21"	18"	18"
SEAT WIDTH	23.5"	22.5"	24"	20"
TOP OF BACK TO SEAT	31"	28"	23"	25"
SITTING DEPTH	39"	38"	34"	34"
FULLY PADDED CHAISE	Yes	Yes	Yes	Yes
FOOTREST EXTENSION	Standard	Standard	Standard	Standard
SEATING TYPE	Total Comfort	Total Comfort	Total Comfort	Total Comfort
BACK STYLE	Deluxe Buscuit Back	Deluxe Buscuit Back	Deluxe Buscuit Back	Deluxe Buscuit Back
HEAT	Standard	Standard	Standard	Standard
CHAIR HEIGHT	49"	43"	40"	41"
DISTANCE FROM WALL	21.5"	15"	14"	15"
OVERALL WIDTH	40"	39"	41"	37"
CHAIR WEIGHT	178 lbs.	158 lbs.	161 lbs.	156 lbs.
WARRANTY	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated







**NOTE:** The wireless charger is only compatible with smart phones that have wireless charging capabilities. Additionally, no other smart devices, (i.e. smart watches) should be used with the wireless charger.

Unit has a stowaway cup holder, and cell phone charging pad.

Unit has a Toggle, Memory hand control.

Measurements may vary +/- 1". Specifications subject to change without prior notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.


## INDEPENDENT FUNCTIONS

-  1. The power headrest lifts your head, neck and shoulders for maximum support.
-  2. The upward lift offers convenient stand assist for getting in and out of the recliner.
-  3. The power lumbar enables you to personalize your comfort.
-  4. The footrest extension allows you to relax and stretch out completely.
-  5. The backrest moves up and down independently to help you find the perfect position.
-  6. Infinite-Positioning recline enables the backrest and footrest to move independently from one another so that a variety of positions can be reached.

## STANDARD FEATURES

7. NEW! Ergonomically designed toggle remote has two programmable memory position buttons. Performs the up/down, recline functions, plus turns your heat on/off and allows you to adjust heat levels for your comfort. Convenient USB charger port allows you to charge your smart device.
8. Lithium battery backup ensures your recliner is operational when the power goes out.

## STANDARD FABRIC: CANYON



OCEAN

SILT

STEEL

WALNUT

Beamsville, Ontario | 888-570-1113  
pridemobility.ca

**PRIDE**<sup>®</sup>  
MOBILITY  
LIVE YOUR BEST<sup>®</sup>