

Dealer / Location: _____ T: _____ F: _____

PO# / Client: _____ Contact: _____ Email: _____

Raz-AT – Attendant Tilt Mobile Shower Commode Chair with 40° tilt range, head support fixture with small pad, height-adjustable stainless-steel frame, flip-up padded arm supports, foot supports, four dual-locking 5" (125mm) casters and commode pan. **User weight capacity 400 lbs.** See our [Evaluation and Medical Justification Form](#).

1. Chair Codes by size – Must Choose One

Z300S18	Raz-AT w/ 18"W seat and base frame – Standard size, fits over most toilets ¹	\$3,550
Z300S_ _*	Non-Standard Width Raz-AT – Must add upcharge for custom width: see below * add custom width to complete code (16" / 20" / 22" / 24") – Adjusta-Back must be selected for all custom widths [Sec.6]	\$3,550
Z3SW16	Upcharge for custom 16"W seat and base frame ² – will not fit over most toilets	\$645
Z3SW20	Upcharge for custom 20"W seat and base frame ³	\$645
Z3SW22	Upcharge for custom 22"W seat and base frame ⁴	\$750
Z3SW24	Upcharge for custom 24"W seat and base frame ⁴	\$850

¹ Fits all 18"W seats

² Visco-Foam Interface seats can be modified to fit.

³ 18"W seats can be selected but must also select Seat Support Adapter Kit ZSSKT [Sec.2]

⁴ ZVFI23FRON Visco Foam Interface HD seat can be modified to fit. Seat Support Kit (ZSSKT) required for 24"W seat frame.

Custom Seat / Base Frame Length and Seat Height

Z3LNGSF3	Upcharge for custom 3" Longer Seat Frame: required for seats ≥ 21"D	\$645
Z3LNGBF3	Upcharge for custom 3" Longer Base Frame: required for seats ≥ 22"D (typically used in conjunction with Z3SFL3)	\$645
Z3LOWF5	Upcharge for custom Super-Low Frame (5" lower) – height adjustment range (molded seat): 16"-20" ¹	\$645

¹ To lower height adjustment range by 1" (15"-19"), select 4" caster upgrade option ZCUT4 (will not fit over most toilets)

2. Seats – Standard

Molded Seats (thickness ~2.5"-3")			Visco-Foam Interface Seats (thickness ~3.5")		
ZMS17FRO	18"W x 17"D – Front Access	Standard	ZVFI17BRI	18"W x 17"D – Front Bridged	\$460
ZMS17BRI	18"W x 17"D – Front Bridged	N/C	ZVFI17FRO	18"W x 17"D – Front Access	\$460
ZMS20FRO	18"W x 20"D – Front Access	\$145	ZVFI17LTO	18"W x 17"D – Left Access	\$460
ZMS20BRI	18"W x 20"D – Front Bridged	\$145	ZVFI17RTO	18"W x 17"D – Right Access	\$460
Specialty Seats			ZVFI17RRO	18"W x 17"D – Rear Access	\$460
ZSOM	Shower-Only Seat – Molded 18"W x 16"-21"D	\$425	ZVFI20BRI	18"W x 20"D – Front Bridged	\$530
ZEJRM	E&J Replica Seat – Molded – 18"W x 17"D	\$210	ZVFI20FRO	18"W x 20"D – Front Access	\$530
ZVFI23FRON	23"W x 18"D HD – Visco Foam – Front Access Notched at rear to maintain fore/aft adjustment range (IPAS) – 22"W frame	\$905	ZVFI20LTO	18"W x 20"D – Left Access	\$530
			ZVFI20RTO	18"W x 20"D – Right Access	\$530
			ZVFI20RRO	18"W x 20"D – Rear Access	\$530

3. Seats – Custom

For custom seats, please specify seat width and depth. Indicate if the access opening is to be in the Front, Left, Right or Rear. Select Bridged if an access opening is not required.

ZCSEAT	Custom Seat (enter specs below) Select if both width and depth are ≤ 21"	\$1110
ZCEJR	Custom E&J Replica Seat (enter specs below) Select if both width and depth are ≤ 21"	\$1265
ZCSEATHD	Heavy Duty Custom Seats (enter specs below) Select if width and/or depth is ≥ 22"	\$1440

Custom Specifications (for above custom seats)

- Seat Frame Width: __ 16" / 18" / 20" / 22" / 24"
- Custom Seat Width: __ 16" / 18" / 20" / 22" / 24" (cannot be more than 2" wider than seat frame width)
- Custom Seat Depth: __ 15" / 16" / 17" / 18" / 19" / 20" / 21" / 22" / 23" / 24" [≥ 21" requires Longer Frame (Z3SFL3)]
- Access Opening: __ FR – Front / LT – Left / RT – Right / RR – Rear / BR – Bridged
- Access Opening Width: __ 2.5" / 3.0" / 3.5" / 4.0" / 4.5" / 5.0" (3.5"W is standard)
- Aperture Width: __ 3.0"-8.0" (6.5" recommended) – enter "NA" if aperture is not required

Unless otherwise requested, the width of the aperture (inside width dimension of the commode opening) will be 6.5" for custom seats.

Custom Seat width should not be narrower than frame width.

Custom Seat can be up to 2" wider than seat frame width to provide additional thigh support without ordering wider chair frame.

The Ischial Pelvic Alignment System (IPAS) allows 2" fore/aft adjustment of the seat (seat depth and aperture position) on all seats.

Special modifications on Custom Seats include: leg-length discrepancy, soft spots, pressure-relief areas, contouring, etc.

4. Seat Options

ZSSKT	Seat Support Kit for Raz-AT/CAT (to support seats that are narrower than seat frames)	\$155	ZGSB __ *	Custom Seat Bracket Spacers – for adding dump angle (pr) Specify additional height at front of seat to complete code (0.5"/1.0"/1.5")	\$315
-------	---	-------	-----------	--	-------

5. Fixed Back Frames (18"W only)

This is called a "Fixed" Back because the arm supports are not adjustable in height. Includes molded, vinyl back support, fixed-height flip-up padded arm supports, back angle adjustment (85°/90°/95°) with cane-style posts. Back height is measured from top of molded seat to top of back upholstery. Hip Adductors [ZHP1/2] are not compatible with Fixed Back.

ZFB1818	Fixed Back – 18"W x 18"H ^{1,2}	Standard
---------	---	----------




6. Adjusta-Back Frames

Includes height- and angle-adjustable arms with Armrest Lock Kit (ZALK), tension-adjustable back upholstery, flip-up padded arm supports, back angle adjustment (85°/90°/95°). Back heights of 15"/17"/20" include cane-style posts, on which back upholstery can be positioned 2.5" higher than nominal height indicated below. e.g. 20"H can be moved up to 22.5"H. 24"H Adjusta-Back includes straight posts, on which the upholstery cannot be raised.

Standard Width Adjusta-Backs (18"W)			Non-Standard Width Adjusta-Backs		
ZAB1815	Adjusta-Back – 15"H	\$265	* Complete Adjusta-Back order code by adding desired back frame width (16"/20"/22"/24") Note: 18" is the standard width.		
ZAB1817	Adjusta-Back – 17"H	\$265	ZAB __ 15*	Adjusta-Back – 15"H non-std width	\$265
ZAB1820	Adjusta-Back – 20"H	\$345	ZAB __ 17 *	Adjusta-Back – 17"H non-std width	\$265
ZAB1824	Adjusta-Back – 24"H	\$395	ZAB __ 20 *	Adjusta-Back – 20"H non-std width	\$345
ZDTL	Dual Triggers – both on left side ¹	N/C	ZAB __ 24 *	Adjusta-Back – 24"H non-std width	\$395
ZDTR	Dual Triggers – both on right side ¹	N/C	ZRBT __ 20*	Reclining Adjusta-Back ² 30° recline, 20" back height with 2-Point Reclining Arms (tilt is limited to 25°)	\$1,235
ZRBT1820	Reclining Adjusta-Back ² 30° recline, 20" back height with 2-Point Reclining Arms (tilt is limited to 25°)	\$1,235	ZRBT __ 24*	Reclining Adjusta-Back ² 30° recline, 24" back height with 2-Point Reclining Arms (tilt is limited to 25°)	\$1,290
ZRBT1824	Reclining Adjusta-Back ² 30° recline, 24" back height with 2-Point Reclining Arms (tilt is limited to 25°)	\$1,290			

¹ Dual Tilt Triggers option positions both triggers on one side, allowing for single-handed activation. Available for all back widths.
² On reclining backs, dual triggers for recline are mounted on the left side and for tilt on the right side.

7. Custom Back Frame Widths Different From Seat Frames – Must choose back frame [Sec.5] (configurations below are not available for recliner backs)

ZNB ___* Custom Back – 2" Narrower than Seat Frame \$370 * Enter seat and back frame widths to complete code. e.g. ZNB2018 specifies a 20"W seat frame with an 18"W back frame	Seat Frame Width	16" ¹	18" ²	20" ³	22" ⁴	24" ⁵	
	Back Frame Width	14"	16"	18"	20"	22"	
ZNB ___* Custom Back – 4" Narrower than Seat Frame – requires custom seat \$480 * Enter seat and back frame widths to complete code. e.g. ZNB2016 specifies a 20"W seat frame with a 16"W back frame	Seat Frame Width	N/A	18" ¹	20" ³	22" ⁴	24" ⁵	
	Back Frame Width	N/A	14"	16"	18"	20"	
ZWB ___* Custom Back – 2" Wider than Seat Frame \$425 * Enter seat and back frame widths to complete code. eg. ZWB1820 specifies an 18"W seat frame with a 20"W back frame	Seat Frame Width	16"	18"	20" ³	22" ⁶	N/A	
	Back Frame Width	18"	20"	22"	24"	N/A	

Consider also using Arm Spacer Kit (ZASK05/ZASK10) as another solution for providing more width between the arms.

¹ Requires Custom Seat
² Compatible with all 18"W seats but with reduced seat depth range. Also compatible with Custom or modified Visco-Foam Interface seats.
³ Custom Seat recommended. Also compatible with all 18"W seats, requires Seat Support Kit ZSSKT [Sec.2]
⁴ Custom Seat recommended. Also compatible with modified ZVFI23FRO (18" seat depth with limited adjustability).
⁵ Custom Seat recommended. Compatible with modified ZVFI23FRON (18" seat depth with limited adjustability, requires Seat Support Kit (ZSSKT) [Sec.2] and modified ZVFI26FRO (18" seat depth with limited adjustability).
⁶ Custom Seat recommended. Also compatible with ZVFI23FRO (18" seat depth).

8. Back Support Options – These back supports replace the standard, tension-adjustable upholstery and are recommended for infection-control environments.

ZSB14	14"W Symphony Back ¹ – recessed solid back	\$315	ZHB Harmony Back ² (contoured solid back – 18"W x 13"H) Fits 17" or 20" back heights, 18"W back frames only	\$315
ZSB16	16"W Symphony Back ¹ – recessed solid back	\$315		
ZSB18	18"W Symphony Back ¹ – recessed solid back	\$315		
ZSB20	20"W Symphony Back ¹ – recessed solid back	\$425		
ZSB22	22"W Symphony Back ¹ – recessed solid back	\$495		
ZSB24	24"W Symphony Back ¹ – recessed solid back	\$520		
ZSB26	26"W Symphony Back ¹ – recessed solid back	\$540		
			¹ Not compatible with 15"H back – all Symphony Back panels are 13"H	
			² Not recommended for reclining backs.	

9. Arm Supports and Accessories

ZARML	Flip-Up Arm Supports – 13"L pad (pr)	Std	ZPAML	Pivoting Arm Mount – Left (For ZFAPLL/ZTRL)	\$180
ZALK	Arm Support Lock Kit (for ZFB1818) (pr)	\$83	ZPAMR	Pivoting Arm Mount – Right (for ZFAPLR/ ZTRR)	\$180
Z2PLT	2-Point Flip-up Arm Supports ¹ (pr) (standard with reclining backs)	\$490	ZASK05	Arm Support Spacer Kit – 1/2" per side (pr)(For Adjusta-backs)	\$67
ZFAPLL	Long Molded Flat Armpad ² – Left (14"L x 3.25"W)	\$88	ZASK10	Arm Support Spacer Kit – 1" per side (pr)(For Adjusta-backs)	\$125
ZFAPLR	Long Molded Flat Armpad ² – Right (14"L x 3.25"W)	\$88	ZPHGL	Pivoting Hand Grips (pr)	\$425
ZTRL	Molded Arm Trough ² – Left (18.5"L x 4.5"W)	\$120	ZAPSB __ *	Anterior Posture Support Bar * Enter back frame width to complete code (14"/16"/18") ^{3,4}	\$210
ZTRR	Molded Arm Trough ² – Right (18.5"L x 4.5"W)	\$120	ZAPSB __ *	Anterior Posture Support Bar * Enter back frame width to complete code (20"/22"/24") ^{3,4}	\$245
¹ For chairs with same back width and frame width only ² May require ZPAM for arm support to flip up and/or for lateral supports to swing away			³ Not available with ZFL / ZTR ⁴ Special length required with ZASK05 (+1") and ZASK10 (+2").		

10. Foot and Leg Supports – Standard Foot Supports (ZFS) are included with the chair. The selection of any Foot / Leg Support Option below automatically replaces the relevant Standard Foot Support.

ZFS	Foot Supports 16"-20.5" (pr)	Std	ZELRL	Elevating Leg Support ³ – Left (includes ZADML)	\$530
ZFS16	Foot Supports for 16"W frames 16"-20.5" (pr)	N/C	ZELRR	Elevating Leg Support ³ – Right (includes ZADMR)	\$530
ZFET2	Footrest Extension Tubes – 2" longer ¹ (pr)	\$57	ZFOR	Forward Offset Foot Support Receiver (ea) Moves Foot support forward 2.5" to accommodate deeper seats (19" or deeper) Specify 1 or 2 ___	\$130
ZMFX8	MFX8 Foot Supports 8"-12" (pr)	\$110	ZRLSL	Universal Residual Limb Support – Left	\$210
ZMFX12	MFX12 Foot Supports 12"-16.5" (pr)	\$110	ZRLSR	Universal Residual Limb Support – Right	\$210
ZADM	Angle/Depth Adj Molded Footplates ² (pr)	\$210			
ZFET2A	Footrest Extension Tubes for ZADM (pr) (2" longer)	\$57			
ZADABS16	Angle/Depth Adj ABS Footplates (pr) For 16"W frames	\$210			
ZVS	V-Style Foot Support (1-piece) – for 18"W base frames	\$315			
ZVS __ *	V-Style Foot Support (1-piece) * Enter frame width to complete code (16" / 20" / 22" / 24")	\$430			
ZFBT	Flip-Back Foot Support (1-piece) with calf pads	\$490			

¹ Not compatible with ZMFX12 / ZADM / ZVS / ZFBAS
² Not compatible with Z300S16 or ZMFX8
³ ZELRL/R compatible with frame widths ≥ 18" and adjusted seat depths ≤ 18" or ≤ 21" w/ 3" longer seat frames (Z3LNGSF3)

11. Foot and Leg Support Accessories

ZFSLK	Foot Support Lock (pr)	\$102	ZAHL	Adjustable Heel Loops (pr) – not available for ZVS / ZFBT	\$67
ZCS __ *	Snap-On Calf Strap – Bodypoint Aeromesh * Enter seat frame width to complete code (16"/18"/20"/22"/24")	\$145	ZNC	Neoprene Covers for Standard Footplates (pr)	\$98
ZCSIC __ *	Calf Strap, Snap-On – Infection Control * Enter seat frame width to complete code (16" / 18" / 20")	\$155	ZNC16	Neoprene Covers for narrow molded footplates (pr)	\$98
ZCSIC22	Snap-On Infection Control Calf Strap – 22"W frames	\$165	ZNCA	Neoprene Covers for Angle/Depth Adj Molded Footplates (pr)	\$98
ZCPL	Calf Panel – Bodypoint ¹ Fits frame widths 16" / 18" / 20"	\$195	ZNCABS16	Neoprene Covers for ABS Adj Angle/Depth Footplates Z7ADABS (pr)	\$98
ZHS1820	H-Strap – for 18" / 20"W frames – hook and loop attachment	\$98	ZLMOR	Lateral / Medial Offset Foot Support Receiver (ea) – moves foot support 2" laterally or medially. Specify 1 or 2 ___	\$88
ZHS2224	H-Strap – for 22" / 24"W frames – hook and loop attachment	\$120			

¹ Not compatible with ZMFX12 / ZFBT / ZVS

12. Positioning Pads

ZHSP1	Head Support Pad – Small (10" x 4")	Std	ZLGT1L	Leg Adductor w/ Small Pad ¹ – Left	\$355
ZHSP2	Head Support Pad – Large (12" x 6")	\$57	ZLGT1R	Leg Adductor w/ Small Pad ¹ – Right	\$355
ZOHIP	Lateral Offset Head Support Interface Plate – up to 22° of adjustment	\$83	ZLGT2L	Leg Adductor w/ Large Pad ¹ – Left	\$355
ZLTS1L	Lateral Thoracic Support w/ Small Pad ^{1,2} – Left	\$355	ZLGT2R	Leg Adductor w/ Large Pad ¹ – Right	\$355
ZLTS1R	Lateral Thoracic Support w/ Small Pad ^{1,2} – Right	\$355	ZHP1L	Hip Adductor w/ Small Pad ^{1,3} – Left	\$355
ZLTS2L	Lateral Thoracic Support w/ Large Pad ^{1,2} – Left	\$355	ZHP1R	Hip Adductor w/ Small Pad ^{1,3} – Right	\$355
ZLTS2R	Lateral Thoracic Support w/ Large Pad ^{1,2} – Right	\$355	ZHP2L	Hip Adductor w/ Large Pad ^{1,3} – Left	\$355
ZLAEP	Lateral/Adductor Extension Plate	\$26	ZHP2R	Hip Adductor w/ Large Pad ^{1,3} – Right	\$355
			ZMTST1	Medial Thigh Support – Small	\$480
			ZMTST2	Medial Thigh Support – Large	\$480

¹ Adjustable Pad Sizes: Small – 5.5" x 3.5" / Large – 6.75" x 4.3".
² Pads on Lateral Thoracic Supports can be adjusted medially up to 3" per side. For additional medial offset, specify extension plate.
³ Fits Adjusta-Back only

13. Pelvic and Chest Belts – A chest belt should not be used without a pelvic belt

Pelvic Belts			Chest Belts		
ZPBS	Pelvic Belt – Small (frame widths 14"-16", 2-piece w/ side-release buckle)	\$46	ZCBS	Chest Belt ² – Small (back widths 14"-16", 1-piece w/ hook and loop closure)	\$36
ZPBM	Pelvic Belt – Medium (frame widths 16"-18", 2-piece w/ side-release buckle)	\$57	ZCBM	Chest Belt ² – Medium (back widths 16"-18", 1-piece w/ hook and loop closure)	\$41
ZPBL	Pelvic Belt – Large (frame widths 18"-22", 2-piece w/ side-release buckle)	\$52	ZCBL	Chest Belt ² – Large (back widths 18"-22", 1-piece w/ hook and loop closure)	\$52
ZPBXL	Pelvic Belt – XLarge (seat frame widths 22"-30", 2-piece w/ side-release buckle)	\$72	ZCBXL	Chest Belt – XLarge ² (back frame widths 22"-30", 1-piece w/ hook and loop closure)	\$67
ZPBICS	Infection Control Pelvic Belt ¹ – Small (frame widths 14"-16", 2-piece w/ side-release buckle)	\$102	ZCBICS	Infection Control Chest Belt ^{1,2} – Small (back widths 14"-16", 1-piece w/ side-release buckle)	\$83
ZPBICM	Infection Control Pelvic Belt ¹ – Medium (frame widths 16"-18", 2-piece w/ side-release buckle)	\$102	ZCBICM	Infection Control Chest Belt ^{1,2} – Medium (back widths 16"-18", 1-piece w/ side-release buckle)	\$83
ZPBICL	Infection Control Pelvic Belt ¹ – Large (frame widths 18"-22", 2-piece w/ side-release buckle)	\$102	ZCBICL	Infection Control Chest Belt ^{1,2} – Large (back widths 18"-22", 1-piece w/ side-release buckle)	\$93
ZPBICXL	Infection Control Pelvic Belt ¹ – XLarge (seat frame widths 22"-30", 2-piece w/ side-release buckle)	\$115	ZCBICXL	Infection Control Chest Belt ^{1,2} – XLarge (back frame widths 22"-30", 1-piece w/ side-release buckle)	\$115
ZPBBPS	Pelvic Belt – Bodypoint Aeromesh – Small	\$305	ZCBBPS	Chest Belt – Bodypoint Aeromesh – Small	\$340
ZPBBP__*	Pelvic Belt – Bodypoint Aeromesh * enter size to complete code (M / L)	\$275	ZCBBP__*	Chest Belt – Bodypoint Aeromesh * enter size to complete code (M / L / XL)	\$285
ZPBBPXL	Pelvic Belt – Bodypoint Aeromesh – XLarge	\$210			

Recommended sizes for Bodypoint belts	Pelvic belt	ZPBBPS	14"/16"W seat frame	ZPBBPM	16"/18"W seat frame	ZPBBPL	18"-22"W seat frame	ZPBBPXL	22"-26"W seat frame
	Chest belt	ZCBBPS	14"/16"W back frame	ZCBBPM	16"/18"W back frame	ZCBBPL	18"/20"W back frame	ZCBBPXL	22"-26"W back frame

¹ Non-absorbent belts – recommended for infection-control environments
² Tension-Adjustable Upholstery includes loops to retain a chest belt (ZCBM/ZCBL/ZCBXL). When the chair is supplied with a Harmony or Symphony Back, Chest Belt loops are provided.

14. Other

ZDCLT	Directional Caster Lock – fits on one rear caster (facilitates steering)	\$180	ZIHW	Whizard – winged urine deflector ²	\$72
ZATTRR	Rear anti-Tippers for 5" casters (pr)	\$370	ZBUDBH __ *	Bard Urine Drainage Bag Holder * Enter size to complete code (16"/18"/20"/22"/24")	\$175
ZCUT4	4" Caster Upgrade (4x)	\$130	ZSPKT22	AT Self-Propel Kit with 22" wheels – incl anti-tippers (tilt limited to 25°)	\$1650
ZCUT6	6" Caster Upgrade (4x)	\$375	ZSPKT24	AT Self-Propel Kit with 24" wheels – incl anti-tippers (tilt limited to 25°)	\$1650
ZPSK1	Commode Pan Spacer Kit ¹ – 1" lower	\$215			
ZPSK2	Commode Pan Spacer Kit ¹ – 2" lower	\$230			
ZPSK3	Commode Pan Spacer Kit ¹ – 3" lower	\$250			
ZPSK4	Commode Pan Spacer Kit ¹ – 4" lower	\$265			
ZTLS	Tilt Limit Spacers – 1/2" spacers snap onto gas struts to limit tilt (each reduces tilt range approx. 5°. Can be stacked for additional tilt limiting. Quantity ____	N/C			

¹ Commode Pan Spacer Kits provide extra clearance for self-care with pan in place. Reduces toilet bowl clearance height.
² Has patent-pending lateral wings that are positioned under the medial aspect of the thighs to keep the Whizard in place.

15. Credits and Omissions

Z905	Omit Commode Pan and Mounting Brackets	-\$20	Z911	Omit Molded Head Support Pad	-\$20
Z907	Omit Foot Supports - if no other Foot / Leg Support Option is selected	-\$45	Z912	Omit Complete Head Support Assembly	-\$50

SOME ASSEMBLY OF RAZ PRODUCTS IS REQUIRED.

Specifications	Raz-AT (Attendant Tilt)
Weight Capacity	400 lbs
Overall Width	Seat width plus 4.5" with Standard Arm Pads
Overall Length (with foot supports and headrest adjusted to middle positions)	
- at 0° tilt and standard foot supports	43"
- at 40°	54"
Seat Depth	16" or 21" depending on seat/cushion length
Toilet Clearance Height (under seat)	
- with pan brackets	17"-21" (1" increments)
- without pan brackets	17.75"-21.75" (1" increments)
Toilet Clearance Width	15.75" for standard 18"W seat frame. For other widths, seat frame width minus 2.25"
Arm Support Height - above molded seat (1/2" increments)	
- 18"H Fixed Back - standard	8.5"
- 15"H Adjusta-Back	5.5"-8.5" (1/2" increments)
- 17"H Adjusta-Back	6.5"-9.5" (1/2" increments)
- 20"H Adjusta-Back	7.5"-10.5" (1/2" increments)
- 24"H Adjusta-Back	7.5"-10.5" (1/2" increments)
- 20"H Reclining Adjusta-Back	7.5"-10.5" (1/2" increments)
- 24"H Reclining Adjusta-Back	7.5"-10.5" (1/2" increments)
Seat-to-Floor Height (with molded seat)	21"-25" in 1" increments
Product Weight (with 18"W seat)	61 lbs
Frame and Hardware Material	Stainless Steel (lifetime warranty against corrosion)