

Brochure

Rhythm Multi Mattress



Blake
Medical Distribution

1-866-662-3459
www.blakemedical.ca



The Rhythm Multi Mattress system provides progressive alternating pressure therapy with true low air loss to aid patients who have existing stage 1, 2, 3 or 4 ulcers. True low air loss provides microclimate, temperature management, and it is caregiver friendly. The Rhythm Multi is a comprehensive therapeutic solution for any clinical care setting.



Happy Heel provides further therapy and comfort by decreasing pressure in patient's heel area



Features

- 8" air cell mattress replacement
- 3-1 alternation
- Static low air loss therapy
- Pulsation therapy
- Happy Heel – dial softness control at foot section
- Auto Firm provides immediate firm surface
- High vapor permeable nylon cover
- 2" convoluted foam base
- Integrated glide sheet to base cover for easy transferring
- Visual and audio alarm for low pressure and power failure with alarm mute function
- Keypad lock out function
- 10 Digital comfort setting
- Modularized design for air mattress
- Integrated power cable retainer
- SWL 600lb
- Quick CPR deflation within 20 seconds

Control Unit

Size (LxWxH)	17.7" x 6.8" x 10.8"
Weight (lbs)	13.2
Cycle Time	3-95 min
Min Operating Pressure	10+/- 3mmHg
Max Operating Pressure	61 +/- 4mmHg
Flow-rate	N/A
Rated Voltage	AC 110-120V
Fuse Rating	5A 250V
Rated Frequency	60Hz
Classification	Class 1, Type BF, Not AP or AGP type
Max Current	5A
Operation Temperature	15°C (59°F)~35°C (95°F)
Operation Humidity	30%~75%
Mode of Operation	Continuous
Standard	IEC 60601-1, CAN/CSA C22.2 No. 601.1, IEC 60601-1-2

Mattress

Size (LxWxH)	36" (W) x 80" (L) x 10" (H)
Weight (lbs)	28.6
Cell Material	Nylon Coated with PU
Cell Number	18 cells
Cover Type	Zipper Cover
Cover Material	Nylon Woven Fabric w/PU coating finish
Base Material	Woven Polyester fabric w/PVC backing

