



FUTURE MOBILITY HEALTHCARE INC.

DIVISION OF ABF MFG GROUP INC.
WWW.FUTUREMOBILITY.CA



TEL: 905-671-1661
FAX: 905-671-3377
orderdesk@futuremobility.ca

Account Information

Account #: _____
Business Name: _____
Ship to Address: _____
Contact: _____
Tag For: _____

Order Date: _____
P.O.# _____
Phone: _____
City: _____
Province: _____ Postal Code: _____

1 FRAME TYPE

WM0001188

Stellato II (Growable) 116 \$1,395
Weight Capacity: 265 lbs.

2 FRAME FINISH

Wet Black 116C01 No Charge
Custom Colours
 Silver Lining 116C02 Net \$150
 Cobalt Blue 116C03 Net \$150
 Green Apple 116C04 Net \$150
 Burnt Orange 116C05 Net \$150
 Grape Purple 116C06 Net \$150
 Fire Red 116C07 Net \$150
 Yellow Ribbon 116C08 Net \$150
 Custom Colours 116C17 Net \$150

3 SEAT WIDTH (Chair Width = Seat Width + 8")

14" 116SW14 No Charge
 15" 116SW15 \$125
 16" 116SW16 No Charge
 17" 116SW17 \$125
 18" 116SW18 No Charge
 19" 116SW19 \$295
 20" 116SW20 No Charge
 21" 116SW21 \$295
 22" 116SW22 \$295

4 SEAT DEPTH

16" 116SD16 No Charge
 17" 116SD17 \$125
 18" 116SD18 No Charge
 19" 116SD19 \$225
 20" 116SD20 No Charge
 21" 116SD21 \$225
 22" 116SD21 \$125

5 REAR AXLE POSITION

Stable → Active Positions (P1 to P3)

P1 Position 116RA01 No Charge
 P2 Position (Standard) 116RA02 No Charge
 P3 Position 116RA03 No Charge
 Amputee 116RA07 \$110
Restricted to RP2-RP6 positions (see STF chart)

6 SEAT-TO-FLOOR HEIGHT (Before Cushion) (Flush Seat Slope)

13.5" 116SF03 No Charge
 14.5" 116SF05 No Charge
 15.5" 116SF07 No Charge
 16.5" 116SF09 No Charge
 17.5" 116SF11 No Charge
 18.5" 116SF13 No Charge
 19.5" 116SF15 No Charge
 20.5" 116SF17 No Charge

7 SEAT SLOPE

Flush 116SS01 No Charge
 1" Drop 116SS02 No Charge
 2" Drop 116SS03 No Charge

8 CASTER SIZE/OPTIONS

4" x 1" Urethane 116CW01 No Charge
 5" x 1.25" Urethane 116CW02 No Charge
 6" x 1.25" Urethane 116CW04 No Charge
 7" x 1.25" Urethane 116CW06 No Charge
 8" x 1.25" Urethane 116CW07 No Charge

9 REAR WHEEL OPTIONS

20" - Composite 116RW03 No Charge
 22" - Composite 116RW04 No Charge
 24" - Composite 116RW05 No Charge
 20" - Spoke (Standard) 116RW06 No Charge
 22" - Spoke (Standard) 116RW07 No Charge
 24" - Spoke (Standard) 116RW08 No Charge

10 REAR WHEEL TIRE


Airless (Solid Rubber) 116RT1 Standard
 Pneumatic 116RT2 No Charge

11 HANDRIM OPTIONS

Chrome 116HR1 Standard
 Plastic Coated 116HR2 \$160
 Omit Handrims 116HR3 No Charge


12 REAR AXLES

Fixed Bolt 116AL1 No Charge
 Quick Release 116AL3 \$110




STELLATO II (SEAT TO FLOOR HEIGHT) STF HEIGHT AVAILABILITY									
Flush Seat Slope		Wheel Size and Axle Position			Fork Position & Caster Wheel Size				
Front STF	Rear STF	20"	22"	24"	4"	5"	6"	7"	8"
13.5" FSTF	13.5" RSTF	RP1			FP0 SF1				
14.5" FSTF	14.5" RSTF	RP2	RP1			FP0 SF2			
15.5" FSTF	15.5" RSTF	RP3	RP2	RP1		FP0 MF2	FP0 MF1		
16.5" FSTF	16.5" RSTF	RP4	RP3	RP2		FP0 MF4	FP0 MF3	FP0 MF2	
17.5" FSTF	17.5" RSTF	RP5	RP4	RP3		FP0 LF4	FP0 LF3	FP0 MF4	FP0 LF1
18.5" FSTF	18.5" RSTF	RP6	RP5	RP4		FP1 LF4	FP1 LF3	FP1 MF4	FP0 LF3
19.5" FSTF	19.5" RSTF	RP7	RP6	RP5		FP2 LF4	FP2 LF3	FP1 LF4	FP1 LF3
20.5" FSTF	20.5" RSTF		RP7	RP6				FP2 LF4	FP2 LF3


POSITION FP0



POSITION FP1



POSITION FP2



STELLATO II (SEAT TO FLOOR HEIGHT) STF HEIGHT AVAILABILITY									
1" Drop Seat Slope		Wheel Size and Axle Position			Fork Position & Caster Wheel Size				
Front STF	Rear STF	20"	22"	24"	4"	5"	6"	7"	8"
14.5" FSTF	13.5" RSTF	RP1				FP0 SF2			
15.5" FSTF	14.5" RSTF	RP2	RP1			FP0 MF2	FP0 MF1		
16.5" FSTF	15.5" RSTF	RP3	RP2	RP1		FP0 MF4	FP0 MF3	FP0 MF2	
17.5" FSTF	16.5" RSTF	RP4	RP3	RP2		FP0 LF4	FP0 LF3	FP0 MF4	FP0 LF1
18.5" FSTF	17.5" RSTF	RP5	RP4	RP3		FP1 LF4	FP1 LF3	FP1 MF4	FP0 LF3
19.5" FSTF	18.5" RSTF	RP6	RP5	RP4		FP2 LF4	FP2 LF3	FP1 LF4	FP1 LF3
20.5" FSTF	19.5" RSTF	RP7	RP6	RP5				FP2 LF4	FP2 LF3

STELLATO II (SEAT TO FLOOR HEIGHT) STF HEIGHT AVAILABILITY									
2" Drop Seat Slope		Wheel Size and Axle Position			Fork Position & Caster Wheel Size				
Front STF	Rear STF	20"	22"	24"	4"	5"	6"	7"	8"
15.5" FSTF	13.5" RSTF	RP1				FP0 MF2	FP0 MF1		
16.5" FSTF	14.5" RSTF	RP2	RP1			FP0 MF4	FP0 MF3	FP0 MF2	
17.5" FSTF	15.5" RSTF	RP3	RP2	RP1		FP0 LF4	FP0 LF3	FP0 MF4	FP0 LF1
18.5" FSTF	16.5" RSTF	RP4	RP3	RP2		FP1 LF4	FP1 LF3	FP1 MF4	FP0 LF3
19.5" FSTF	17.5" RSTF	RP5	RP4	RP3		FP2 LF4	FP2 LF3	FP1 LF4	FP1 LF3
20.5" FSTF	18.5" RSTF	RP6	RP5	RP4				FP2 LF4	FP2 LF3

13 ARM TYPE

- NEW! Al. Flip Back Vinyl Arm Pad 10" 116AR30 \$260
- NEW! Al. Flip Back Vinyl Arm Pad 14" 116AR35 \$260
- NEW! Al. Flip Back Gel Arm Pad 16" 116AR36 \$300
- T-Style Height Adj. Vinyl Arm Pad 14" 116AR04 Standard
- T-Style Height Adj. Vinyl Arm Pad 10" 116AR03 No Charge
- T-Style Height Adj. Gel Arm Pad 16" 116AR40 \$91

14 BACK CANE

- Fixed Angle Adj 10° Bend w/Push Handle (low) 116PC3 \$90
- Fixed Angle Adj 10° Bend w/Push Handle (medium) 116PC13 Standard
- Fixed Angle Adj 10° Bend w/Push Handle (high) 116PC4 \$90
- NEW! Paddle Flip-Down Backcane (Height Adjustable) 116PC20 \$155

Set Seat to Back Angle at:

-21° ___ -17.5° ___ -14° ___ -10.5° ___ -7° ___ -3.5° ___ 0° ___ +3.5° ___ +7° ___

15 BACK OPTIONS (CHOOSE ONE)

- Slip-on Padded Back Upholstery 116BK48 \$100
- Adjustable Tension Back 116BK45 \$150
- Prism Basic Back 116BK01 \$150
- Prism Ultra Back 116BK06 \$250
- Prism Air Back 116BK11 \$400
- Prism TrueFitt Back 116BK16 \$650
- Prism Max Ultra Back 116BK21 \$300
- Prism Max Air Back 116BK26 \$450
- Prism Max TrueFitt Back 116BK31 \$700
- Lunar Back Mid Basic Lateral 116BK39 \$550
- Omit Back 116BK49 No Charge

W () x H ()

For other Backs, please fill in the required Backs Order Form and attach it to this order

16 BACK PINS

- Longer Pins for Back Hardware 1" Growth 116BP1 \$50
- Longer Pins for Back Hardware 2" Growth 116BP2 \$50

17 SEAT OPTIONS (CHOOSE ONE)

- Black Nylon Seat Upholstery 116SC19 Standard
- Prism Ideal Cushion 116SC01 \$100
- Prism Elite Cushion 116SC04 \$250
- Prism Supreme Cushion 116SC07 \$375
- Prism Supreme II Plus Cushion 116SC10 \$375
- Prism Supreme Air Cushion 116SC13 \$500
- Omit Seat Cushion 116SC20 No Charge

Add On Item

- 1" Drop Base on Cushion (Over Rail Only) 116SC37 \$170

Drop Base can only fit Future Mobility Cushions

W () x D ()

For other Cushions, please fill in the required Cushions Order Form and attach it to this order

18 FRONT RIGGING (PAIRS)

<input type="checkbox"/> 60° Swingaway Upgraded	116FR03	\$230
<input type="checkbox"/> 70° Swingaway Upgraded	116FR04	No Charge
<input type="checkbox"/> 80° Swingaway Upgraded	116FR45	\$230
<input type="checkbox"/> 90° Swingaway Upgraded	116FR05	\$230
<input type="checkbox"/> Elevating Legrest	116FR10	\$550
<small>Not available with composite/adjustable footplate with 6", 7" or 8" Casters</small>		
<input type="checkbox"/> Stump Support Swingaway Upgraded Pair	116FR32	\$240
Single Pick Option <small>Same limitations as above pairs</small>		
<input type="checkbox"/> 60° Swingaway Upgraded	<input type="checkbox"/> RT <input type="checkbox"/> LT 116FR35	\$230
<input type="checkbox"/> 70° Swingaway Upgraded	<input type="checkbox"/> RT <input type="checkbox"/> LT 116FR36	No Charge
<input type="checkbox"/> 80° Swingaway Upgraded	<input type="checkbox"/> RT <input type="checkbox"/> LT 116FR46	\$115
<input type="checkbox"/> 90° Swingaway Upgraded	<input type="checkbox"/> RT <input type="checkbox"/> LT 116FR37	\$115
<input type="checkbox"/> Elevating Legrest	<input type="checkbox"/> RT <input type="checkbox"/> LT 116FR42	\$275
<input type="checkbox"/> Stump Support Swingaway	<input type="checkbox"/> RT <input type="checkbox"/> LT 116FR44	\$120
<input type="checkbox"/> Omit Front Rigging	<input type="checkbox"/> RT <input type="checkbox"/> LT 116FR11	No Charge

19 FOOTPLATE TYPE (PAIRS)

<input type="checkbox"/> Composite Depth & Angle Adj. Flip-Up <small>Small</small>	116FP27	\$160
<input type="checkbox"/> Composite Depth & Angle Adj. Flip-Up <small>Medium</small>	116FP28	\$160
<input type="checkbox"/> Composite Depth & Angle Adj. Flip-Up <small>Large</small>	116FP29	\$160
<input type="checkbox"/> Aluminum Depth & Angle Adj. Flip-Up <small>Small</small>	116FP10	\$160
<input type="checkbox"/> Aluminum Depth & Angle Adj. Flip-Up <small>Medium</small>	116FP11	\$160
<input type="checkbox"/> Aluminum Depth & Angle Adj. Flip-Up <small>Large</small>	116FP12	\$160
<input type="checkbox"/> Composite Angle Adj. Flip-Up Cone Lock <small>Large</small>	116FP16	\$150
<input type="checkbox"/> Composite Angle Adj. Flip-Up Cone Lock <small>X-Large</small>	116FP17	\$170
<input type="checkbox"/> Composite Flip Up <small>Small</small>	116FP01	No Charge
<input type="checkbox"/> Composite Flip Up <small>Medium</small>	116FP02	No Charge
<input type="checkbox"/> Composite Flip Up <small>Large</small>	116FP03	No Charge
<input type="checkbox"/> Aluminum Flip Up <small>Small</small>	116FP05	Standard
<input type="checkbox"/> Aluminum Flip Up <small>Medium</small>	116FP06	Standard
<input type="checkbox"/> Aluminum Flip Up <small>Large</small>	116FP07	Standard
<input type="checkbox"/> Aluminum Flip Up <small>X-Large</small>	116FP08	Standard
<input type="checkbox"/> One Piece Flip-Up Footboard	116FP19	\$160
<input type="checkbox"/> Foam Pad One Piece Footboard	116FA06	\$210
<small>Not Available with 6", 7" or 8" Casters when footboard is facing rear on 90° Swingaway</small>		

20 FOOTREST/LEGREST ACCESSORIES

<input type="checkbox"/> Padded Calf Strap	116FA03	\$70
<input type="checkbox"/> Foam Guard on Elevating Legrest (pair)	116FA05	\$70
<input type="checkbox"/> Foam Pad on Footplate (pair)	116FA07	\$70
<input type="checkbox"/> Heel Loops	116FA01	Standard
<input type="checkbox"/> Omit Heel Loops	116FA02	No Charge

NOTES:**21 WHEEL LOCKS**

<input type="checkbox"/> Push-to-Lock	116WL5	Standard
<input type="checkbox"/> Pull-to-Lock	116WL6	No Charge
<small>Pull to Lock recommended for Seat Depth 18" and smaller</small>		
<input type="checkbox"/> Wheel Lock Extension	116WL7	No Charge
<input type="checkbox"/> Attendant Wheel Lock (Foot Activated) RT	116WL8	\$295
<input type="checkbox"/> Attendant Wheel Lock (Foot Activated) LT	116WL9	\$295
<input type="checkbox"/> Unilateral Push Wheel Lock RT Activated	116WL10	\$250
<input type="checkbox"/> Unilateral Push Wheel Lock LT Activated	116WL11	\$250
<input type="checkbox"/> Unilateral Pull Wheel Lock RT Activated	116WL12	\$250
<input type="checkbox"/> Unilateral Pull Wheel Lock LT Activated	116WL13	\$250

22 ANTI TIP

<input type="checkbox"/> Fixed	116MA22	No Charge
<input type="checkbox"/> Quick Flip	116MA23	\$100

23 MISCELLANEOUS ACCESSORIES

<input type="checkbox"/> Auto Seat Positioning Strap	116MA03	\$52
<input type="checkbox"/> Velcro Release Positioning Strap	116MA05	\$48
<input type="checkbox"/> IV Pole, telescopic	116MA20	\$168
<input type="checkbox"/> Cane Crutch Holder	116MA37	\$85
<input type="checkbox"/> Transit Tie Down (Unoccupied)	116MA74	\$50
<input type="checkbox"/> Spoke Guards	116MA49	\$130
<input type="checkbox"/> Seat Pouch	116MA75	\$50

WARRANTY:

Future Mobility Healthcare warrants the STELLATO II wheelchair to be free from defects in materials and workmanship for a period of two (2) year from the date of purchase. The side frame is warranted for the lifetime of the original purchaser and/or user.