

Poly-Fly High Strength, Lightweight Wheelchair/Flyweight Transport Chair Combo

HCPCS:



FIGURE A



FIGURE B



FIGURE C



FIGURE D

- Can be used as standard self-propelled wheelchair or a transport chair all in one
- Quick-release 24" wheels can be removed to transition self-propelled chair to a transport chair (Figure A - B)
- All-aluminum frame
- Comes with two sets of wheel locks: One for use with wheelchair, one for use with transport chair
- Back folds down with deluxe back release (Figure C)
- Nylon upholstery
- Comes standard with swing-away footrests with composite foot plates and heel loops
- 8" casters
- Comes with anti-tippers
- STDELR-TF - Swing-away Elevated Leg rest, Tool-free, (Optional front riggings)
- Weight Capacity: 250 lbs

PRODUCTS

Description

TR18	PolyFlyW/C/TrsPrt,w/SwgFt,18"AntiTp,1/cs
TR20	PolyFlyWC/TPrt,20w/SgwFt.AntiTp,1/cs

SPECIFICATIONS

Description

Warranty	Limited Lifetime
Weight Capacity	250 lbs.
Seat To Floor Height	19.5"
Carton Shipping Weight	38 lbs 39 lbs 41 lbs
Back Height	16"

PARTS

Description

STDS2A1621	RearAxle,QuickRelease,Polyfly,1/ea, EA
STDS1012	Rear Wheel Bearing,11.5 x 32 x 9.5mm,1/e, EA
STDS2A161R	Brake,Right,Polyfly,24"wheel,1/ea, EA
STDS2A161L	Brake,Left,Polyfly,24"wheel,1/ea, EA
STDS2A163R	Brake,Right, 8"Wheel,Polyfly,1/ea, EA
STDS2A163L	Brake,Left,8"Wheel,Polyfly,1/ea, EA
STDS2A1603	8"Rear Wheel,Polyfly,1/ea, EA
STDS1005	Rear Axle Bolts,Sentra,1/pkg, EA
STDS2A1611	RearAxle Receptacle,Polyfly,1/ea, EA
STDS2A1610	SlideTube w/Bushing,Polyfly,1/ea, EA
STDS2A1687	Left Side Frame,TR16 and TR18,1 ea, EA
STDS2A1612	Round Brake Pad,8",Polyfly,1/ea, EA

PARTS	Description
STDS2A161P	WheelBrakePad24"Polyfly,1 ea, EA
WCP 1287	Frame Plug,Polyfly,1/ea, EA
STDSSB	Spring Btn for Arm Rest,1/ea., EA
STDS004012	Handrim,Black,24"Wheel/1 ea.", EA
STDS2A1600	24" Wheel,BLK,HandRim TR18,1/ea, EA
STDS2A1658	Anti-tipper,Polyfly,1/pr, EA
STDS2A4204	Armrest SidePanel,Lt.Plastic,1/ea, EA
STDS2A162R	Armrest,Rt.Desk,Polyfly,1/ea, EA