

Bariatric Sentra Extra-Heavy-Duty Wheelchair

HCPCS: K0007,K0007



FIGURE A

- Weight Capacity: 500 lbs
- Carbon steel frame with triple-coated chrome for an attractive, chip-proof, maintainable finish
- Reinforced steel gussets at all weight bearing points provide additional strength
- Durable, heavy-gauge reinforced naugahyde upholstery
- Comes standard with carry pocket on backrest
- Reinforced side frame, caster journals and dual cross bars
- Dual-position front fork
- Dual-position front frame
- Composite Mag-style wheels with chrome hand rim are lightweight and maintenance free
- Hemi-adjustable caster is optional (Part #STDS1074)
- Precision sealed wheel bearings in front and rear ensure long-lasting performance and reliability
- Padded armrests provide added patient comfort
- Comes with swing-away footrests or elevating leg rests with tool-free adjustable length riggings
- Comes with push-to-lock wheel locks

PRODUCTS	Description
STD20DDA-ELR	SentraHD,20"DetDeskArmElvLeg,1/cs
STD20DDA-SF	SentraHD,20"DetDeskArm,Swg.Ft,1/cs
STD22DDA-ELR	SentraHD,22"DetDeskArmElvLeg,1/cs
STD22DDA-SF	SentraHD,22"DetDeskArm,Swg.Ft,1/cs
STD24DDA-ELR	SentraHD,24"DetDeskArmElvLeg,1/cs
STD24DDA-SF	SentraHD,24"DetDeskArm,Swg.Ft,1/cs

SPECIFICATIONS	Description
Back Height	18"
Carton Shipping Weight	60 lbs 62 lbs 72 lbs
Height	37"
Length with Rigging	42"
Product Weight with 20 Seat	54 lbs
Product Weight with 22 Seat	56 lbs
Product Weight with 24 Seat	66.25 lbs
Seat Depth	18"
Seat To Floor Height	17.5" - 19.5"
Warranty	Limited Lifetime
Weight Capacity	500 lbs
Width Closed	13"

SPECIFICATIONS	Description
Width Open	29" 31" 32.5"
ACCESSORIES	Description
WASL	Limb Support,Left,1/cs, CS
WASR	Limb Support, Right,1/cs, CS
STDS831	Heel Loops, 1pr, BX
STDS1034	Wheelchair Cane/Crutch Holder,1/ea, EA
STDS820	Telescoping IV Pole,W/C, 1/bx, BX
STDS803	Oxygn Tank Holdr for W/C,1/cs, CS
STDS834	OverheadAnti-Theft Bar,Single,1/ea, EA
STDS806	Anti-Fold/Theft Device for w/c,1/ea, BX
STDS855	Seat Belt,60"" Auto-type,1/ea, EA
STDELR-AL	ArticulatingElev.Legrest, Std,1pr/cs, CS
STDS856	Seat Belt,60"Velcro-type,1/ea, EA
STDS801	Brake Handle Ext.,8", 1/ea, EA
STDS860	Forwrd Stablzrs/Univ.Fit ,1pr/bx, BX
STDS802	Anti-Tipprs,w/Wheels,Adj,1pr/bx, BX
PARTS	Description
STDS1051	Dust Cap, fork stem, 1/ea, EA
STDS871	Hand Grip,Sentra,1/ea, EA
STDS1026	Front Bumper for Ft.Rst/LegRest,1/ea, EA
STDS1002-1	Axle for 8" caster, 1/ea, EA
STDS1060	Wheelchair Fork Bearing,1/ea, EA
STDS1011	Wheel Bearing, 8mm x 22mm 1/ea, EA
STDS8M4601	ForkStemBearing,SSP2,Cruiser,1/ea, EA
STDS1070	Seat Rail Guide, Removable Arm, 1/ea, EA
STDS4020	Tipping Sleeve 1/ea, EA
STDS005200	SlideTube,Sentra,1/ea, EA
STDS002001	Arm Pad Screw,1/ea, EA
STDS4011	BrakeTip, 1/ea, EA
STDS2329BLK	BackUphol black,SentraHD 20inx 18in,1/bx, BX
STDS2331BLK	Back Uphol 22"x18",Black,1/bx, BX
STDS1052	Rear AxleBolt,SentraHD,1/ea, EA
STDS2333BLK	Back Uphol 24"x18", Black,1/bx, BX
STDS2A5222	Crossframe,22"SentraHD/1ea, EA
STDS1074	6"Hemi Caster,SentraHD,1ea, EA
STDS4010R	Brake Asmb.Right Push to Lock,1/ea, EA
STDS1000HD	Wheel,24""Hvy.Dty.Hand Rl.Sentra,1/ea, EA