

## Half Length Bed Rail Bariatric Lightweight Beds

HCPCS:



- Durable steel construction in rich brown finish offers attractive and easy-to-maintain homecare style
- Easy to install half Rail designed specifically for Lightweight Bariatric Homecare Bed series
- Product Safety & Performance Standards: Incompliance with FDA Hospital Bed System Dimensional Assessment Guidance to Reduce Entrapment as well as CPSC ASTM standards
- Ideal for maximum protection and support when bed is elevated
- 
- 

PRODUCTS	Description
15021BV-LW	Half Side Rail,42in/48/54 LW Baribed,1pr

  

SPECIFICATIONS	Description
Color	Brown
Warranty	Steel Frame: Limited Lifetime   Knobs: 1 Year
Material Composition	Rails: Steel   Knobs: Molded plastic