

Med-Aire Assure 5" Air + 3" Foam Base Alternating Pressure and Low Air Loss Mattress System 14530



FIGURE A



FIGURE B



FIGURE C

- Active mattress replacement system provides both alternating pressure and low air loss to optimize pressure redistribution, shear/friction reduction, and microclimate control
- Designed to prevent, treat and heal pressure ulcers in the home or long term care setting
- Control unit offers 10 minute cycles when alternating and can also be used in static mode
- Adjustable patient weight settings allow for optimal immersion and patient comfort & compliance
- Audio and Visual alarm alerts when air pressure is low
- System includes 17 individual air cells over a 3" foam base for power outage protection
- Pillow feature provides head stability & comfort
- 5" deep air cells are easily removed for convenient servicing
- Removable quilted Nylon cover is fluid resistant, low shear and vapor permeable helping protect the skin from friction and moisture
- CPR valve allows for rapid deflation
- Tubing quickly and easily disconnects with EZ Lock/Quick Release design

PRODUCTS	Description
14530	8in Low air loss/APM w/foam base,1/ea

SPECIFICATIONS	Description
Control Unit Airflow	8 LPM
Control Unit Alarms	Visual / Audible
Control Unit Cycle Time	10 Minutes
Control Unit Type	Compressor
Control Unit Power	110 VAC \ 60 Hz
Mattress Cover Material	Fluid-Resistant Nylon
Warranty	1 Year Limited
Weight Capacity	350 lbs

ACCESSORIES	Description
14530P	Pump Only, 14530, 1/ea, EA
14530M	Mattress Only W/Foam Base, 14530, 1/ea, EA
14530C	Cover for 14530, 1/ea, EA
14530-cell	5 Air Cell, Non-Ventilated, 14530,1/ea, EA