



Pride Mobility Products Co.

5096 South Service Rd, Beamsville, ON, L0R 1B3

Phone: 888-570-1113 | Fax: 866-514-1303 | Email: productorders@pridemobility.com

Q6 Edge® 3 Order Form with Contoured Seating

Account Number: _____ Date: _____

Provider Name: _____

Contact: _____

Phone: _____ Fax: _____

Email: _____

PO Number: _____

Marked for: _____

Ship to Address: _____

City: _____ Prov: _____ PC: _____



INTRODUCTION

This form is interactive when viewed with [Adobe Acrobat Reader](#). Complete the form by placing checks in the desired boxes and provide information in the interactive fields. Buttons shown at the bottom of the form may be utilized to print or submit the order form through a desktop email application. To email via a web-based application, please 'Save As' and attach the PDF to your email.

This order form contains a large variety of options to fulfill various patient needs. Descriptions and section notes, such as optional or required, are included to help you complete your order. Other part numbers can be found on the quote acknowledgment after this order form is submitted. Please contact Quantum® Sales at 888-570-1113 if assistance is needed. Send the completed order form by fax (866-514-1303) or email (productorders@pridemobility.com). Incomplete forms may delay the quote or order. Customer service will contact you if the order is incomplete or if there are compatibility issues. If special order requests are needed, be sure that the Patient Information section is completed

All prices are MSRP. Prices, specifications, part numbers, and availability are subject to change without notice. Prices and part numbers as shown are only available when configured as a complete power chair. Please contact technical service for accurate parts ordering.

Thank you for choosing Quantum!


PATIENT INFORMATION

Weight: _____ Height: _____

A. Shoulder Width: _____ D. Max Sitting Height: _____ G. Thigh Depth: _____

B. Chest Width: _____ E. Shoulder Height: _____ H. Lower Leg Length: _____

C. Hip Width: _____ F. Axilla Height: _____ I. Elbow Height: _____



1. BASE MODEL

REQUIRED. Select a group 3 power base. 300 lb. weight capacity with standard 6 mph programmable motors and 8 amp off board charger. XLR USB charger is included. Dimensions of base without legrests: 24"Wx35.5"L

Quantum Q6 Edge 3 3S-SS \$8,450
Part: Q6 EDGE 3. ADP: WP0001575 (ADULT), WP0001576 (PEDIATRIC)

2. DRIVE MOTORS


















REQUIRED.

- | | |
|---|---|
| 6 mph 2 Pole Drive Motors..... Standard
Part: MOT147433 | 7 mph, 4 Pole Drive Motors..... \$700
Part: MOT1902455
Requires Q-Logic 3 electronics. |
| 6 mph, 2 Pole, Accu-Trac Motor Technology..... \$995
Part: MOT147881.
Requires Q-Logic 3 electronics. Delivers enhanced traction and tracking, especially with switched drive controls. | 7 mph, 4 Pole, Accu-Trac Motor Technology..... \$1,695
Part: MOT2008501
Requires Q-Logic 3 electronics. Delivers enhanced traction and tracking, especially with switched drive controls. |

3. BASE COLOR

REQUIRED. Select one shroud color. Colors shown here may differ from actual product depending on monitor/printer calibration used.

3.1a. Standard Shrouds

- | | | | |
|---|--|--|--|
|  Red Corvette
Part: DGN1801101-PRD1 |  White Satin
Part: DGN1801101-PWT2 |  Midnight Blue
Part: DGN1801101-PBU10 |  Military Green
Part: DGN1801101-PGN8 |
|  Cherry Bomb
Part: DGN1801101-PRD4 |  Deep Purple
Part: DGN1801101-PPR5 |  Electric Blue
Part: DGN1801101-PBU9 |  Sugar Plum
Part: DGN1801101-PPR7 |
|  Back in Black
Part: DGN1801101-PBK1 |  Raspberry Beret
Part: DGN1801101-PPK4 |  Blue Jean
Part: DGN1801101-PBU11 | |
|  Night Train
Part: DGN1801101-PBK3 |  Pretty in Pink
Part: DGN1801101-PPK3 |  Lemon Crush
Part: DGN1801101-PYL4 | |
|  Grey Street
Part: DGN1801101-PSL3 |  Evolution Orange
Part: DGN1801101-POR7 |  Green Machine
Part: DGN1801101-PGN11 | |

3.1b. Hydrographic Shrouds

- | | | |
|--|--|---|
| <input type="checkbox"/>  Military Digital Camo \$325
Part: DGN1801101-WCM7 | <input type="checkbox"/>  Faded Glory \$325
Part: DGN1801101-WFL2 | <input type="checkbox"/>  Carbon Fiber \$325
Part: DGN1801101-WCB1 |
|--|--|---|

4. POWER ADJUSTABLE SEAT HEIGHT

Optional. The seat to floor height will be increased. Lift/elevate will be operated through the joystick. Power elevating seat includes a 1/4 speed inhibit.

- 12" Power Lift for Captain Seat with iLevel..... \$2,495
Part: PTOINDV3857. HCPCS: E2300
For use with K0857 base code. Must select [Power Lift Captain Seat](#) in contoured seat section. Includes armrests and lap belt. Available with swing away legrests or manually adjustable center mount foot platform.

5. BASE MOUNTED LIGHTS

Optional. Select one option.

- | | |
|--|--|
| Switch Fender Lights..... No Charge
Part: DGN154117
Fender lights operate by pressing included single switch button. Fender lighting is not available with Single Actuator Control through Toggle Box. | Base Mounted, LED, Full Light Package \$770
Part: XRF000129
Requires Q-Logic 3 Joystick for Full Light Package or Q-Logic 3e Joystick. |
|--|--|

6. DRIVE WHEELS

REQUIRED. Select a type of drive wheels and accent color. Knobby tires increase overall base width by 1.5".

6.1. Drive Wheels

- | | |
|---|--|
| 3"x14" Flat Free Black Tires Standard
Part: PTOINDV3525 (pair) | 3"x14" Pneumatic Black Tires No Charge
Part: WHL1713925 x2 |
| 4"x14" Knobby, Flat Free Black Tires \$120
Part: WHL2006647 x2 | 4"x14" Knobby, Pneumatic Black Tires..... \$120
Part: WHL2006822 x2 |

6.2. Drive Wheel Accent Color

Not available with knobby tires.

- | | | |
|---|---|--|
| Back in BlackStandard
Part: DGN164245-PBK1 x2 | Raspberry BeretNo Charge
Part: DGN164245-PPK4 x2 | Green MachineNo Charge
Part: DGN164245-PGN11 x2 |
| Red CorvetteNo Charge
Part: DGN164245-PRD1 x2 | Pretty in PinkNo Charge
Part: DGN164245-PPK3 x2 | Military GreenNo Charge
Part: DGN164245-PGN8 x2 |
| Cherry BombNo Charge
Part: DGN164245-PRD4 x2 | Evolution OrangeNo Charge
Part: DGN164245-POR7 x2 | Sugar PlumNo Charge
Part: DGN164245-PPR7 x2 |
| Night TrainNo Charge
Part: DGN164245-PBK3 x2 | Midnight BlueNo Charge
Part: DGN164245-PBU10 x2 | Military Digital Camo \$50 ea.
Part: DGN164245-WCM7 x2 |
| Grey StreetNo Charge
Part: DGN164245-PSL3 x2 | Electric BlueNo Charge
Part: DGN164245-PBU9 x2 | Faded Glory \$50 ea.
Part: DGN164245-WFL2 x2 |
| White SatinNo Charge
Part: DGN164245-PWH2 x2 | Blue JeanNo Charge
Part: DGN164245-PBU11 x2 | Carbon Fiber \$50 ea.
Part: DGN164245-WCB1 x2 |
| Deep PurpleNo Charge
Part: DGN164245-PPR5 x2 | Lemon CrushNo Charge
Part: DGN164245-PYL4 x2 | |

7. BATTERIES

Optional. The powerchair requires 2 batteries and will be installed if ordered.

NF-22, Gel Cell, Interceptor Battery..... \$300 ea.
Part: BAT1704613.

NF-22, Gel Cell, MK Battery \$400 ea.
Part: BATGEL1008.

8. PATIENT'S WEIGHT

REQUIRED.

Up to 200 lb User Weight
Part: PTOINDV3837

201 lb-300 lb User Weight
Part: PTOINDV3838

10. ELECTRONICS

REQUIRED. Select a type of drive control from the following sections. The No Joystick option is used to omit the joystick and utilize another control, such as a specialty control at the end of this form or found on the [Alternative Controls and Electronics Order Form](#).

10.1a. Non-Expandable Controls

Select a joystick below or omit selection if another controller type is desired. Non-expandable does not allow supplementary drive controllers.

NE 4-Key Joystick.....Standard
Part: CTLDC1554

NE+ 8-Key Joystick (2 Actuators) No Charge
Part: CTLDC1561

NE+ 6-Key Joystick..... No Charge
Part: CTLDC1560

Control up to 2 actuators through joystick and only Seat Mounted, LED, full light package.

No Joystick, Setup for NE+ Joystick No Charge
Part: XRF000009



10.1b. Expandable Controls - Quantum Intelligent Systems™

Select expandable controls or omit selection if another controller type is desired. Expandable controls allow multiple or specialized controllers.

Q-Logic 3 Joystick

Q-Logic 3 EX Joystick \$900
Part: CTL167724

Q-Logic 3 EX Joystick for Full Lights \$900
Part: CTL167725

Q-Logic 3e Joystick \$620
Part: CTL1910600
Can operate full light package. Not available with alternative or specialty controls.

No Joystick, Setup for Q-Logic
Part: XRF000010



10.2. Q-Logic 3 EX Joystick Armor

REQUIRED with Q-Logic 3 EX Joystick selections above. Not available with Q-Logic 3e Joystick. Reference colors shown on page 2. *Night train joystick armor will be the black, soft touch.

Black, Soft Touch No Charge
Part: DGN124307
Also matches Night Train shrouds.

Non-matching Armor No Charge
Part: PTOINDV3422
Select color below.

Match Armor to Base Color Standard
Part: PTOINDV3421

Red Corvette Part: DGN124336-PRD1	Grey Street Part: DGN124336-PSL3	Pretty in Pink Part: DGN124336-PPK3	Blue Jean Part: DGN124336-PBU11	Sugar Plum Part: DGN124336-PPR7
Cherry Bomb Part: DGN124336-PRD4	White Satin Part: DGN124336-PWH2	Evolution Orange Part: DGN124336-POR7	Lemon Crush Part: DGN124336-PYL4	Military Digital Camo Part: DGN124336-WCM7
Back in Black Part: DGN124336-PBK1	Deep Purple Part: DGN124336-PPR5	Midnight Blue Part: DGN124336-PBU10	Green Machine Part: DGN124336-PGN11	Faded Glory Part: DGN124336-WFL2
Night Train* See Black, Soft Touch	Raspberry Beret Part: DGN124336-PPK4	Electric Blue Part: DGN124336-PBU9	Military Green Part: DGN124336-PGN8	Carbon Fiber Part: DGN124336-WCB1

11. JOYSTICK HANDLES

Optional. The NE/NE Plus Joysticks are standard with a carrot handle (See A in image). The optional handle will be mounted to the joystick and the standard handle will be shipped with the unit.

Mushroom Handle (B) \$120 Part: ACCASMB1860.	Extended I-Handle (G) \$140 Part: ACCASMB1857.
Chin Cup (C) \$85 Part: ACCASMB186190	Bodypoint 3" Goal Post Handle \$180 Part: PLSKNCP1057.
T-Bar (D) \$160 Part: ACCASMB1858.	Bodypoint 4" Goal Post Handle \$180 Part: PLSKNCP1058.
Medium Ball Handle (E) \$70 Part: ACCASMB1859.	
Soft Ball Handle (F) \$120 Part: ACCASMB1856.	



12. JOYSTICK MOUNTING BRACKETS

Optional. Select a mount to be used with a joystick selected in previous section.

12.1a. Fixed, In-line Joystick Bracket


Fixed, In-line, Bracket for Q-Logic 3Standard Part: FRMASMB10617	Fixed, In-line, Bracket for NE/NE+Standard Part: FRMASMB12744
---	--



12.1b. Swing-Away (Retractable) Joystick Bracket

Allows the joystick to retract to the side of the armrest while maintaining the same joystick orientation.

Swing-away, Q-Logic 3 \$265 Part: ACC1908434.	Swing-away, NE/NE+ \$265 Part: ACC1908437.
--	---



12.1c. Specialty Joystick Mounts

Select a bracket and side to mount it.

Retract4 Mount \$450 Part: XRF000003. Full 180° of movement that translates into multi-position use, straight, right offset, left offset, or fold away. Joystick orientation varies as bracket is adjusted.	Stealth Swing-Away/Flip Down Bracket \$300 Part: XRF000005. Allows the joystick to pivot to side of armrest, or flip-down if bracket is adjusted.
Multi-Axis Mount \$100 Part: XRF000004 Allows height adjustment of joystick.	



12.2. Joystick Mount Location

Left Side Mount Part: PTOINDV1225	Right Side Mount Part: PTOINDV1226
--------------------------------------	---------------------------------------

13. CONTOURED SEATING

Optional. Select only one seating system. Special order seat sizes can be requested by checking the seat type and writing in the seat size in the special order seat size.

13.2a. Quantum Captain Seat

16"W x 16"D No Charge Part: SET2005145	18"W x 20"D No Charge Part: SET2005151	20"W x 22"D No Charge Part: SET2005909
16"W x 18"D No Charge Part: SET2005146	18"W x 22"D No Charge Part: SET2005908	22"W x 18"D No Charge Part: SET2005155
16"W x 20"D No Charge Part: SET2005147	20"W x 16"D No Charge Part: SET2005152	22"W x 20"D No Charge Part: SET2005156
18"W x 16"D No Charge Part: SET2005149	20"W x 18"D No Charge Part: SET2005153	22"W x 22"D No Charge Part: SET2005154
18"W x 18"D No Charge Part: SET2005150	20"W x 20"D No Charge Part: SET2005154	

13.2b. Quantum Captain Seat with Solid Seat Pan

16"W x 16"D No Charge Part: SET2005394	18"W x 16"D No Charge Part: SET2005397	18"W x 22"D No Charge Part: SET2005910
16"W x 18"D No Charge Part: SET2005395	18"W x 18"D No Charge Part: SET2005398	20"W x 16"D No Charge Part: SET2005400
16"W x 20"D No Charge Part: SET2005396	18"W x 20"D No Charge Part: SET2005399	20"W x 18"D No Charge Part: SET2005401

Contoured seating section continued on following page.

13. CONTOURED SEATING CONTINUED

20"W x 20"D..... No Charge Part: SET2005402	22"W x 18"D..... No Charge Part: SET2005403	22"W x 22"D..... No Charge Part: SET2005405
20"W x 22"D..... No Charge Part: SET2005911	22"W x 20"D..... No Charge Part: SET2005404	


13.3. Matching Seat Cushion

Optional with Quantum Captain Seat with Solid Seat Pan. Specific cushion/component sizes available, part numbers, weight capacities, MSRP, HCPCS code, and special order cushions can be found on the respective complete [Stealth Cushion Order Forms](#). This section allows a quick selection of a single cushion that matches the seat pan width and depth. See specific order form for non-matching sizes, alternate covers and liners, or additional covers and liners.

Stealth Simplicity G Cushion

General use cushion with Stealth Tek 2, 2-way stretch, breathable cover.
Part series: ST-SIMGXXYY. HCPCS: E2601/E2602


Simplicity G, Matching Seat Size
Reference part: XRF000111. MSRP: \$110/\$165



Stealth Glacial SP Cushion

Skin protection with Coolcore, 4-way stretch, breathable cover.
Part series: ST-GLSPXXYY. HCPCS: E2603/E2604


Glacial SP, Matching Seat Size
Reference part: XRF000114. MSRP: \$330/\$440



Stealth Zen SP Cushion

Skin protection with Coolcore, 4-way stretch, breathable cover.
Part series: ST-ZENSPXXYY. HCPCS: E2603/E2604


Zen SP, Matching Seat Size
Reference part: XRF000115. MSRP: \$315/\$425



Stealth Premiere P Cushion

Positioning with Coolcore, 4-way stretch, breathable cover.
Part series: ST-PREMPXXYY. HCPCS: E2605/E2606


Premiere P, Matching Seat Size
Reference part: XRF000122. MSRP: \$325/\$475



Stealth Spectrum Foam SPP Cushion

Skin protection and positioning with Coolcore, 4-way stretch, breathable cover. Includes urethane liner.
Part series: ST-SPCFSPXXYY. HCPCS: E2607/E2608

Spectrum Foam SPP, Matching Seat Size
Reference part: XRF000112. MSRP: \$460/\$530



Stealth Essence SPP Cushion


Skin protection and positioning with Coolcore, 4-way stretch, breathable cover.
Part series: ST-ESSSPXXYY. HCPCS: E2607/E2608

Essence SPP, No Wedge, Matching Seat Size
Reference part: XRF000118. MSRP: \$400/\$500

Essence SPP, Small Wedge, Matching Seat Size
Reference part: XRF000119. MSRP: \$435/\$535

Essence SPP, Medium Wedge, Matching Seat Size
Reference part: XRF000120. MSRP: \$465/\$565


Essence SPP, Large Wedge, Matching Seat Size
Reference part: XRF000121. MSRP: \$495/\$595



Stealth Solution SPP Cushion

Skin protection and positioning with Coolcore, 4-way stretch, breathable cover. Includes urethane liner.
Part series: ST-SOLSPPXXYY. HCPCS: E2607/E2608

Solution SPP, Matching Seat Size
Reference part: XRF000113. MSRP: \$375/\$430



13.4. Armrests

Optional. Configured in pairs. Left (L) and right (R) boxes may be checked to configure different size arms within one section arm type.

1a. Drop-In Armrests

Includes flat armpads.

Full Length, Drop-In Armrests..... No Charge
L R Left part: ARM1803253. Right part: ARM156742

Desk Length, Drop-In Armrests No Charge
L R Left part: ARM1803254. Right part: ARM156743

1b. Flip-Up, Height Adjustable Armrests

Includes full length, flat armpads.

Full Length, Left Armrest No Charge
Part: MEC1906767

Full Length, Right Armrest..... No Charge
Part: MEC1906768

1c. Single Post, Height Adjustable, Removable Armrests

Height adjustment from 11"-17" in half-inch increments.

Full Length, Single Post Armrests\$145 ea.
L R Left part: ARM2000906. Right part: ARM2000907.

Desk Length, Single Post Armrests\$145 ea.
L R Left part: ARM2000908. Right part: ARM2000909.

1d. Bariatric, Height Adjustable, Removable Armrests

Select armrests and armpads. Armrests are sold in same-size pairs.

Full Length, Bariatric Armrests, 8"-12" H.....\$130 ea.
L R Left part: FRMASMB6728. Right part: FRMASMB6729

Full Length, Bariatric Armrests, 10"-14"H.....\$130 ea.
L R Left part: FRMASMB6730. Right part: FRMASMB6731

1e. Heavy Duty, Height Adjustable, Removable Armrests

Select armrests and armpads. Armrests are sold in same-size pairs.

Full Length, Heavy Duty Armrests, 8"-12" H\$130 ea.
L R Left part: XRF000253. Right part: XRF000254.

Full Length, Heavy Duty Armrests, 10"-14"H.....\$130 ea.
L R Left Part: XRF000257. Right part: XRF000258.

1f. Pediatric, Removable, Height Adjustable Armrests

Height adjustment from 7.75" to 8.25". Includes desk length, straight armpads.

Pediatric, Removable Armrests \$145 ea.
L R Left Part: ARM142876. Right part: ARM142877. HCPCS: E0973

Contoured seating section continued on following pages.

13. CONTOURED SEATING CONTINUED

13.5. Armrest Height

REQUIRED with Bariatric, Heavy Duty, and Pediatric Height Adjustable armrests.

2. Left Armrest Height

7.75" H PTOINDV3722	8.25" H PTOINDV3724	10" H PTOINDV3537	12" H PTOINDV3539	14" H PTOINDV3541	16" H PTOINDV3543
8" H PTOINDV3535	9" H PTOINDV3536	11" H PTOINDV3538	13" H PTOINDV3540	15" H PTOINDV3542	

3. Right Armrest Height

7.75" H PTOINDV3721	8.25" H PTOINDV3723	10" H PTOINDV3452	12" H PTOINDV3454	14" H PTOINDV3456	16" H PTOINDV3458
8" H PTOINDV3450	9" H PTOINDV3451	11" H PTOINDV3453	13" H PTOINDV3455	15" H PTOINDV3457	

13.6. Armpads

REQUIRED except with Flip-Up and Pediatric armrests.

1a. Straight Armpads for Drop-In, or Single Post Armrests

Straight, Full Length Left and Right No Charge Part: ARM134275	Straight, Desk Length Left, Full Length Right No Charge Part: ARM134280
Straight, Full Length Left, Desk Length Right No Charge Part: ARM134281	Straight, Desk Length Left and Right No Charge Part: ARM134254

1b. Waterfall Armpads for Drop-In, or Single Post Armrests

Waterfall, Full Length Left and Right No Charge Part: ARM134285	Waterfall, Desk Length Left, Full Length Right No Charge Part: ARM134288
Waterfall, Full Length Left, Desk Length Right No Charge Part: ARM134289	Waterfall, Desk Length Left and Right No Charge Part: ARM134282

1c. Falcon Gel Armpads for Drop-In, or Single Post Armrests

Falcon Gel, Full Length Left and Right \$300 Part: ARM136438	Falcon Gel, Desk Length Left, Full Length Right \$300 Part: ARM136439
Falcon Gel, Full Length Left, Desk Length Right \$250 Part: ARM136442	Falcon Gel, Desk Length Left and Right \$200 Part: ARM136437

1d. Armpads for Bariatric or Heavy Duty Armrests

Bariatric Armpads \$40 ea. Part: UPSASMB1027 Only for use with Bariatric Armrests.	Waterfall Armpads No Charge Part: UPSASMB1024 Only for use with Heavy Duty Armrests
Straight Armpads No Charge Part: UPSASMB1021 Only for use with Heavy Duty Armrests	

14. LAP BELTS

REQUIRED. Select one lap belt below.

14.1a. Standard Lap Belts

<input type="checkbox"/> 38" Lap Belt Part: ACCBELT1006	<input type="checkbox"/> 50" Lap Belt Part: ACCBELT1000	<input type="checkbox"/> 60" Lap Belt Part: ACCBELT1002	<input type="checkbox"/> 70" Lap Belt Part: ACCBELT1003	<input type="checkbox"/> 80" Lap Belt Part: ACCBELT1004	<input type="checkbox"/> 90" Lap Belt Part: ACCBELT1001
--	--	--	--	--	--

14.1b. Optional Lap Belts

Stealth Retractable Lap Belt (IFRB5350).....\$175 Part: ST-IFRB5350	Stealth Positioning, 60" Single Pull, Push Button, Extra Large... \$65 Padded (11.5"), 2" Belt (IFPB5315) Part: ST-IFPB5315-Q
Stealth Positioning, 60" Single Pull, Push Button, Medium..... \$65 (7.5") Padded, 1.5" Belt (IFPB5213) Part: ST-IFPB5213-Q	60" Sure-Lok Lap Belt (Non-retractable)..... \$200 Required with ISO 7176-19 Compliant Transit Option. Part: ACCBELT1016
Stealth Positioning, 60" Single Pull, Push Button, Large \$65 Padded (9.5"), 2" Belt (IFPB5314) Part: ST-IFPB5314-Q	

15. LEGRESTS

Optional. All efforts will be made to accommodate legrest length selection, but minor changes may need to be made to meet shipping guidelines. Some adjustments may be needed upon receipt of the unit. Reference the section's compatibility matrix for legrest interference.

15.1a. Base Mounted Foot Platform

May interfere with seat mounted legrests. Not available with center mount foot platforms.

Base Mounted Foot Platform Standard Part: FRMASMB16077	Omit base mounted foot platform No Charge Part: PTOINDV3674
Bell Shaped, Base Mounted Foot Platform \$235 Part: RIGASMB7120041	

15.1b. Center Mount Foot Platform

If power articulating foot platform is selected, be sure power legrest control has been selected in earlier section. Select legrest length in following section.

1. Foot Platform Assembly

Manually Adjustable Center Mount Foot Platform.....\$265 Part: XRF000080

2. Lower Extension and Angle Adjustable Footplates

10"-16" Extension, 5"x6" Footplates Part: RIG150206	10"-16" Extension, 5"x11" Footplates Part: RIG150212	13"-19" Extension, 5"x8" Footplates Part: RIG151506
10"-16" Extension, 5"x8" Footplates Part: RIG150203	13"-19" Extension, 5"x6" Footplates Part: RIG151505	13"-19" Extension, 5"x11" Footplates Part: RIG151508

3. Center Mount Calf Support

Small Calf Pads (pair) 6"Wx6"H.....\$210 Part: RIGASMB7120018
--

4. Center Mount Modular Footplate Extension Components

Optional. Assembly required. May require seating configuration adjustments to avoid front caster interference.

Width Extension for 8"D Footplates \$55 Part: RIG152517	4" Outside Edge (pair) \$55 Part: RIG152510 Only available with 6"D footplates and 8"D footplate extension.
Width Extension for 11"D Footplates \$55 Part: RIG152518	6" Outside Edge (pair) \$55 Part: RIG152516 Available with 8"D and 11"D footplates and 11"D footplate extension.
2" Bell-Shaped Length Extension \$55 Part: RIG152519 Not compatible with width extensions. Can be used with 6" outside edge with 11" footplates.	

6. Stealth Gel Footplate Covers

5"x8" Tapered Gel Footplate Covers..... \$275 Part: ST-FPG5008T-Q, K0108	6.5"x8" Tapered Gel Footplate Covers for Width Extension..... \$275 Part: ST-FPG6508T-Q, K0108
5"x11" Tapered Gel Footplate Covers..... \$275 Part: ST-FPG5011T-Q, K0108	7"x11" Tapered Gel Footplate Cover for Width Extension \$275 Part: ST-FPG7011T-Q, K0108

15.1c. Swing-Away Legrests

Select a legrest upper, a lower extension, and type of footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes.

1a. 60°-90° Swing-Away Legrest Uppers

Adjustable in 5° increments.

<input type="checkbox"/> P <input type="checkbox"/> L <input type="checkbox"/> R 60°-90° Swing-Away Legrest..... \$285 Pair part: RIG120062. Left: FRMASMB17558. Right: FRMASMB17559 Not available with 11"-13" extensions.	<input type="checkbox"/> P <input type="checkbox"/> L <input type="checkbox"/> R 60°-90° Swing-Away Legrest, Short \$285 Pair part: RIG136481. Left: RIG132607. Right: RIG132609. For use only with 11"-13" lower extensions.
--	--

2. Lower Extension

<input type="checkbox"/> P <input type="checkbox"/> L <input type="checkbox"/> R 11"-13" Length Pair part: XRF000088. Left: XRF000089. Right: XRF000090. N/A with Multi-axis.	<input type="checkbox"/> P <input type="checkbox"/> L <input type="checkbox"/> R 11"-16.5" Length Pair part: XRF000091. Left: XRF000092. Right: XRF000093. N/A with Angle Adjustable.	<input type="checkbox"/> P <input type="checkbox"/> L <input type="checkbox"/> R 13"-16.5" Length Pair part: XRF000094. Left: XRF000095. Right: XRF000096. N/A with Multi-Axis.	<input type="checkbox"/> P <input type="checkbox"/> L <input type="checkbox"/> R 15.5"-19.5" Length Pair part: XRF000097. Left: XRF000098. Right: XRF000099.
---	---	---	---

3a. Extruded Footplate(s)

<input type="checkbox"/> P <input type="checkbox"/> L <input type="checkbox"/> R 7.25"Wx6"D Extruded Footplate..... Standard Part pair: FRMASMB7847. Left: FRMASMB7861 Right: FRMASMB7862	<input type="checkbox"/> P <input type="checkbox"/> L <input type="checkbox"/> R 8.5"Wx6"D Extruded Footplate Standard Part pair: FRMASMB7848. Left: FRMASMB7863 Right: FRMASMB7864 May cause interference on seat widths less than 18."
---	---

Legrests section continued on following pages.

15. LEGRESTS CONTINUED

60°-90° Swing-Away Legrest Compatibility Matrix	Contoured/Synergy Seat		
	Low	Medium	High
Extruded, 13" - 16.5"	!	✓	✓
Extruded, 15.5" - 19.5"	!	!	!
Angle Adjustable, 13" - 16.5"	!	✓	✓
Angle Adjustable, 15.5" - 19.5"	!	!	!
Multi-Axis, 15.5" - 19.5"	!	!	!

This rigging configuration causes no interference issues at this seat-to-floor height.

This rigging configuration may have less than 1.75" clearance when fully extended and may interfere with the ground or reduce obstacle climbing ability.

This rigging configuration does not work with this seat-to-floor height.

3b. Angle Adjustable Footplate(s) with Heel Loops

¹The 6"x11" footplates may cause interference with seat widths less than 18". ²The 8"x11" footplates may cause interference with seat widths less than 20".

- | | |
|--|---|
| <p><input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> 4"Wx6"D \$87.50 ea.
Pair part: FRMASMB10564. Left: RIG133537. Right: RIG133560</p> <p><input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> 5"Wx6"D..... \$87.50 ea.
Pair part: FRMASMB10565. Left: RIG133538. Right: RIG133561</p> <p><input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> 5"Wx8"D \$87.50 ea.
Pair part: FRMASMB10566. Left: RIG133539. Right: RIG133562</p> | <p><input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> 6"Wx8"D¹ \$87.50 ea.
Pair part: FRMASMB10567. Left: RIG133540. Right: RIG133563</p> <p><input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> 6"Wx11"D¹ \$87.50 ea.
Pair part: FRMASMB10568. Left: RIG133541. Right: RIG133564</p> <p><input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> 8"Wx11"D² \$87.50 ea.
Pair part: FRMASMB10569. Left: RIG133542. Right: RIG133565</p> |
|--|---|

3c. Multi-Axis Footplate(s) with Heel Loops

¹The 6"x11" footplates may cause interference with seat widths less than 18". ²The 8"x11" footplates may cause interference with seat widths less than 20".

- | | |
|---|--|
| <p><input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> 4"Wx6"D \$260 ea.
Pair part: FRMASMB9456. Left: RIG133531. Right: RIG133552</p> <p><input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> 5"Wx6"D..... \$260 ea.
Pair part: FRMASMB9457. Left: RIG133532. Right: RIG133553</p> <p><input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> 5"Wx8"D \$260 ea.
Pair part: FRMASMB9458. Left: RIG133533. Right: RIG133554</p> | <p><input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> 6"Wx8"D¹ \$260 ea.
Pair part: FRMASMB9459. Left: RIG133534. Right: RIG133555</p> <p><input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> 6"Wx11"D¹ \$260 ea.
Pair part: FRMASMB9460. Left: RIG133535. Right: RIG133558</p> <p><input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> 8"Wx11"D² \$260 ea.
Pair part: FRMASMB9461. Left: RIG133536. Right: RIG133559</p> |
|---|--|

15.1d. Pediatric, Swing-Away Legrests

Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes. Select legrest length in following section.

1. Pediatric Legrests Upper

- 90°, Pediatric, Swing-Away Legrest (7"-10" Length Extension)** \$285
Pair part: FRMASMB8808. Left: FRMASMB8809
Right: FRMASMB8810

2. Angle Adjustable Footplate(s)

- | | |
|--|--|
| <p><input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> 4"Wx6"D Angle Adjustable \$87.50 ea.
Pair part: FRMASMB8812. Left: RIG133537. Right: RIG133560
HCPCS: K0040</p> | <p><input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> 5"Wx6"D Angle Adjustable \$87.50 ea.
Pair part: FRMASMB8811. Left: RIG133538. Right: RIG133561
HCPCS: K0040</p> |
|--|--|

Legrests section continued on following page.

15. LEGRESTS CONTINUED

15.1e. Elevating Legrests

Select a legrest upper, a legrest extension, and type of footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes. Select legrest length in following section.

1. Elevating Legrests Upper

P L R Elevating Legrest Assembly \$275 ea.
 Pair part: FRMASMB6483. Left: RIG138012
 Right: RIG138013. HCPCS: E0990
 Select lower extension and footplates

P L R **11"-14" Length**
 Pair: FRMASMB4755
 Left: RIG138026
 Right: RIG138030
 For reinforced only

P L R **14"-17" Length**
 Pair: FRMASMB4756
 Left: RIG138033
 Right: RIG138035
 For reinforced only

P L R **12"-14.5" Length**
 Pair part: XRF000101.
 Left: XRF000102.
 Right: XRF000103.
 For angle adjustable only

P L R **14.5"-18" Length**
 Pair part: XRF000104.
 Left: XRF000105.
 Right: XRF000106
 For angle adjustable only

Elevating Legrest Compatibility Matrix	Contoured/Synergy Seat		
	Low	Medium	High
Standard, 11"-14"	✓	✓	✓
Standard, 14"-17"	!	!	!
Angle Adjustable, 12"-14.5"	✓	✓	✓
Angle Adjustable, 14.5"-18"	!	!	!

This rigging configuration causes no interference issues at this seat-to-floor height.
 This rigging configuration may have less than 1.75" clearance when fully extended and may interfere with the ground or reduce obstacle climbing ability.

3a. Reinforced Cast Footplate(s)

P L R **7.25"Wx6"D Reinforced Footplate**..... Standard
 Pair part: FRMASMB4748. Left: RIG138055. Right: RIG138056

P L R **8.5"Wx6"D Reinforced Footplate** Standard
 Pair part: FRMASMB4747. Left: RIG138060. Right: RIG138062
 May cause interference on seat widths less than 18."

3b. Angle Adjustable Footplate(s)

¹The 6"x11" footplates may cause interference with seat widths less than 18". ²The 8"x11" footplates may cause interference with seat widths less than 20". 4"x6" and 5"x6" footplates are only available in pairs.

P L R **4"Wx6"D** \$87.50 ea.
 Pair parts: FRMASMB8812 x2.

P L R **6"Wx8"D¹ Footplates with Heel Loops** \$87.50 ea.
 Pair part: FRMASMB4738. Left: RIG137793. Right: RIG137795

P L R **5"Wx6"D**..... \$87.50 ea.
 Pair parts: FRMASMB8811 x2.

P L R **6"Wx11"D¹ Footplates with Heel Loops** \$87.50 ea.
 Pair part: FRMASMB4741. Left: RIG137796. Right: RIG137797

P L R **5"Wx8"D Footplates with Heel Loops** \$87.50 ea.
 Pair part: FRMASMB4735. Left: RIG137791. Right: RIG137792

P L R **8"Wx11"D² Footplates with Heel Loops** \$87.50 ea.
 Pair part: FRMASMB4744. Left: RIG137798. Right: RIG137799

15.2. Legrest Accessories

Heel Loops (pair)..... \$35
 Part: FRMASMB7873. HCPCS: E0951
 Standard on angle adjustable and multi axis footplates when noted.

Wheel Bumper \$15 ea.
 Part: XRF000100

15.3. Set Legrest Extension Length

REQUIRED with matching left and right legrest configuration. Part number is PTOINDV plus number following # below. All efforts will be made to accommodate legrest length selection, but minor changes may need to be made to meet shipping guidelines. Some adjustments may be needed upon receipt of the unit. Reference the section's compatibility matrix for legrest interference. Non-matching or individual legrest configurations lengths are not set by the manufacturer.

7" L (#3470)	8.5" L (#3473)	10" L (#3476)	11.5" L (#3479)	13" L (#3482)	14.5" L (#3485)	16" L (#3488)	17.5" L (#3491)	19" L (#3494)	20.5" L (#3497)	22" L (#3669)
7.5" L (#3471)	9" L (#3474)	10.5" L (#3477)	12" L (#3480)	13.5" L (#3483)	15" L (#3486)	16.5" L (#3489)	18" L (#3492)	19.5" L (#3495)	21" L (#3498)	22.5" L (#3670)
8" L (#3472)	9.5" L (#3475)	11" L (#3478)	12.5" L (#3481)	14" L (#3484)	15.5" L (#3487)	17" L (#3490)	18.5" L (#3493)	20" L (#3496)	21.5" L (#3600)	

16. POWER LIFT ELECTRONICS

REQUIRED with 12" Power Lift for Captain Seat with iLevel.

Single Actuator Control through NE+ Joystick..... No Charge
 Part: ELE2006331

Single Actuator Control through Q-Logic Joystick..... No Charge
 Part: ELE2006330

17. ACCESSORIES

Optional. Multiple accessories can be mounted along the various tracks but some interference may occur. Accessories may need to be installed by the provider.

Cup Holder, C Style \$20 <i>Part: ACCASMB2288</i> For use with Flip-Up armrests.	Quantum Backpack \$200 <i>Part: ACC1905853</i>
Cup Holder, Bariatric Arms or Gel Armpads \$20 <i>Part: ACCASMB2632</i>	Reflector Kit \$35 <i>Part: ACC124670</i>
Cup Holder, Tube Style Armrests \$20 <i>Part: XRF000221</i> For HD Removable, Single Post, or Pediatric Arms.	Clothing Guard (Pair) \$275 <i>Part: ACC125007</i> Available with Flip-Up armrests only. Not available with Elevating Legrests.
Oxygen Holder \$110 <i>Part: ACC2006353. HCPCS: E2208</i>	Transfer Handles \$305 <i>Part: MEC1809930. HCPCS: E0705</i> Available with Flip-Up armrests only. Not available with Elevating Legrests.
Cane/Crutch Holder \$195 <i>Part: ACC2006352. HCPCS: E2207</i>	Privacy Bag \$80 <i>Part: ACC1907570</i>
Rear Basket \$200 <i>Part: ACC2006351</i>	
Saddle Bag \$15 <i>Part: DWR1020D010</i>	
Large, Black, Weather Cover \$240 <i>Part: ACCCOVR1021</i>	

18. PROGRAMMER

Optional.

<input type="checkbox"/> Q-Logic EX/NE Programmer \$625 <i>Part: ELEASMB6976</i>
--