



Pride Mobility Products Co.

5096 South Service Rd, Beamsville, ON, L0R 1B3

Phone: 888-570-1113 | Fax: 866-514-1303 | Email: productorders@pridemobility.com

Q6 Edge® Z Order Form with Contoured Seating

Account Number: _____ Date: _____

Provider Name: _____

Contact: _____

Phone: _____ Fax: _____

Email: _____

PO Number: _____

Marked for: _____

Ship to Address: _____

City: _____ Prov: _____ PC: _____

Q6 edge[®] Z



INTRODUCTION

This form is interactive when viewed with [Adobe Acrobat Reader](#). Complete the form by placing checks in the desired boxes and provide information in the interactive fields. Buttons shown at the bottom of the form may be utilized to print or submit the order form through a desktop email application. To email via a web-based application, please 'Save As' and attach the PDF to your email.

This order form contains a large variety of options to fulfill various patient needs. Descriptions and section notes, such as optional or required, are included to help you complete your order. Other part numbers can be found on the quote acknowledgment after this order form is submitted. Please contact Quantum® Sales at 888-570-1113 if assistance is needed. Send the completed order form by fax (866-514-1303) or email (productorders@pridemobility.com). Incomplete forms may delay the quote or order. Customer service will contact you if the order is incomplete or if there are compatibility issues. If special order requests are needed, be sure that the Patient Information section is completed

All prices are MSRP. Prices, specifications, part numbers, and availability are subject to change without notice. Prices and part numbers as shown are only available when configured as a complete power chair. Please contact technical service for accurate parts ordering.


PATIENT INFORMATION

Weight: _____ Height: _____

A. Shoulder Width: _____ D. Max Sitting Height: _____ G. Thigh Depth: _____

B. Chest Width: _____ E. Shoulder Height: _____ H. Lower Leg Length: _____

C. Hip Width: _____ F. Axilla Height: _____ I. Elbow Height: _____



1. BASE MODEL

REQUIRED. Select a group 3 power base. 300 lb. weight capacity with standard 6mph programmable motors and 8 amp off board charger. Dimensions of base without legrests: 26.25"W x 35.5"L.

Quantum Q6 Edge Z 3S-SS.....\$8,080

Part: Q6 EDGE Z. ADP Adult: WP0001269, ADP Pediatric: WP0001270.

1A. FENDER LIGHTS

Optional. Fender lights operate by pressing included single switch button. Fender lighting is not available with Single Actuator Control through Toggle Box.

Switch Fender Lights..... \$240*
Part: DGN154117. *No charge with iLevel.

2. BASE COLOR








REQUIRED. Select one shroud color. Colors shown here may differ from actual product depending on monitor/printer calibration used.

2.1. Match Colors

Match All Color Options to Selection Below

Part: N/A. Any subsequent color selections, if available, will be matched to color chosen below when the order is placed. Please note that some color selections, such as joystick armor, may have a MSRP.

2.2. Power Base Shroud

- | | | | |
|---|--|--|--|
|  Red Corvette
Part: DGN1610108-PRD1 |  Night Train
Part: DGN1610108-PBK3 |  Evolution Orange
Part: DGN1610108-POR7 |  Green Machine
Part: DGN1610108-PGN11 |
|  Electric Blue
Part: DGN1610108-PBU9 |  Pretty in Pink
Part: DGN1610108-PPK3 |  Grey Street
Part: DGN1610108-PSL3 | To Be Determined
Part: N/A |
|  Back in Black
Part: DGN1610108-PBK1 |  Deep Purple
Part: DGN1610108-PPR5 |  Midnight Blue
Part: DGN1610108-PBU10 | |

3. DRIVE MOTORS

REQUIRED.

- | | |
|--|---|
| 6 mph Drive Motors Standard
Part: MOT150186 | Accu-Trac Motor Technology \$995
Part: MOT150181. |
| High Speed Motor Package..... \$700
Part: MOT150141
Programmable speed up to 8 mph. High speed motors decrease per-charge range. | Requires Q-Logic electronics. Delivers enhanced traction and tracking, especially with switched drive controls. Programmable speed up to 6 mph. |

4. BASE MOUNTED FOOT PLATFORM

REQUIRED. Select a foot platform option. May interfere with seat mounted legrests.

- | | |
|--|---|
| Base Mounted Foot Platform Standard
Part: FRMASMB16077
Not available with center mount foot platforms. | Bell Shaped, Base Mounted Foot Platform \$235
Part: RIGASMB7120041 |
| Omit Base Mounted Foot Platform No Charge
Part: Configuration selection | |

5. BATTERIES

REQUIRED. The powerchair requires 2 batteries and will be installed if ordered.

- | | |
|--|---|
| Group 24, 70 Amp Gel Cell, Interceptor Battery \$400 ea.
Part: BAT1704614. | Group 24, 70 Amp Gel Cell, MK Battery\$510 ea.
Part: BATGEL1009. |
| No batteries. Provider will assume responsibility of batteries.
Part: Configuration selection | |

6. DRIVE WHEELS

REQUIRED. Select a type of drive wheels. Knobby tires increase overall base width by 1.5".

6.1. Drive Wheels

- | | |
|--|--|
| 14" Run Flat Tires Standard
Part: PTOINDV3510 | 4"W Flat Free, Knobby, Extended Traction Tires \$120
Part: WHL121698 x2 |
| 14" Run Flat Tires for Drive Wheel Accents Standard
Part: WHL1703421 x2 | 4"W Pneumatic, Knobby, Extended Traction Tires \$120
Part: WHL153737 x2 |

6.2. Drive Wheel Accent Color


Optional. Not available with knobby tires. Select with WHL1703421 above. Sold in Pairs (parts below will be quantity of 2).

- | | | | | |
|--|---|--|--|---|
| Match to Base Color\$42.50 ea.
Part: Configuration selection, see specific part below | Non-matching Color\$42.50 ea.
Part: Configuration selection, see specific part below | | | |
| Red Corvette
Part: DGN164245-PRD1 | Night Train
Part: DGN164245-PBK3 | Deep Purple
Part: DGN164245-PPR5 | Grey Street
Part: DGN164245-PSL3 | Green Machine
Part: DGN164245-PGN11 |
| Electric Blue
Part: DGN164245-PBU9 | Pretty in Pink
Part: DGN164245-PPK3 | Evolution Orange
Part: DGN164245-POR7 | Midnight Blue
Part: DGN164245-PBU10 | To Be Determined
Follow-up required. |
| Back in Black
Part: DGN164245-PBK1 | | | | |

7. FULL LIGHT PACKAGE

Optional. Requires NE+ 8 Key or Q-Logic 3 Joystick for Full Light Package.

- | |
|--|
| Seat Mounted, LED, Full Light Package \$770
Part: ELE131477 |
|--|



8. ELECTRONICS

REQUIRED. Select a type of drive control from the following sections. The No Joystick option is used to omit the joystick and utilize another control, such as a specialty control at the end of this form or found on the [Alternative Controls and Electronics Order Forms](#).

8.1a. Non-Expandable Controls

Select a joystick below or omit selection if another controller type is desired. Non-expandable does not allow supplementary drive controllers.

NE+ 6-Key Joystick Standard

Part: CTLDC1560

Control up to 2 actuators through joystick.

No Joystick, Setup for NE+ Joystick No Charge

Part: Configuration selection

NE+ 8-Key Joystick (2 Actuators)..... No Charge

Part: CTLDC1561

Control up to 2 actuators and Seat Mounted LED Light Package through joystick.



NE+ 6 Key Joystick

8.1b. Expandable Controls - Quantum Intelligent Systems™

Select expandable controls or omit selection if another controller type is desired. Expandable controls allow multiple or specialized controllers.

2. Q-Logic 3 Joystick

Q-Logic 3 EX Joystick \$900

Part: CTL167724

No Joystick, Setup for Q-Logic

Part: Configuration selection.

Q-Logic 3 EX Joystick for Full Lights \$900

Part: CTL167725



Q-Logic 3 EX Joystick

8.2. Q-Logic 3 Joystick Armor

REQUIRED with Q-Logic 3 Joystick selections above. Reference colors shown on page 2.

Black, Soft Touch Standard

Part: DGN124307

Non-matching Armor (Select below) \$50

Part: PTOINDV3422

Match Armor to Base Color \$50

Part: PTOINDV3421

Red Corvette
Part: DGN124336-PRD1

Night Train
See Black, Soft Touch

Deep Purple
Part: DGN124336-PPR5

Midnight Blue
Part: DGN124336-PBU10

Electric Blue
Part: DGN124336-PBU9

Pretty in Pink
Part: DGN124336-PPK3

Grey Street
Part: DGN124336-PSL3

Green Machine
Part: DGN124336-PGN11

Back in Black
Part: DGN124336-PBK1

Evolution Orange
Part: DGN124336-POR7

9. JOYSTICK HANDLES

REQUIRED. The NE Plus are standard with a carrot handle (See A in image). The optional handle will be mounted to the joystick and the standard handle will be shipped with the unit.

Mushroom Handle (B) \$120

Part: ACCASMB1860.

Extended I-Handle (G) \$140

Part: ACCASMB1857.

Chin Cup (C) \$85

Part: ACCASMB186190

Bodypoint 3" Goal Post Handle \$180

Part: PLSKNCP1057.

T-Bar (D) \$160

Part: ACCASMB1858.

Bodypoint 4" Goal Post Handle \$180


Part: PLSKNCP1058.

Medium Ball Handle (E) \$70

Part: ACCASMB1859.

Soft Ball Handle (F) \$120

Part: ACCASMB1856.



Joystick Handles

10. JOYSTICK MOUNTING BRACKETS

Optional. Select a mount to be used with a joystick selected in previous section.

Omit Joystick Mounting Bracket


10.1a. Fixed, In-line Joystick Bracket

Fixed, In-line, Bracket Standard

Q-Logic Part: FRMASMB10617. NE/NE+ Part: FRMASMB12744

Left Side Mount
Part: PTOINDV1225

Right Side Mount
Part: PTOINDV1226



Swing-Away (Retractable) (right side)

10.1b. Swing-Away (Retractable) Joystick Bracket

Allows the joystick to retract to the side of the armrest while maintaining the same joystick orientation. "Inline bracket" will be offset 1" when used with tube style arms. Select "Offset Bracket" to have the joystick inline with armpad on tube style armrests.

Swing-away, Inline, Left Side Mount \$265

Part: See XRef. ADP: WEPM

Swing-away, 1" Offset, Left Side Mount \$265


Part: See XRef. ADP: WEPM

Swing-away, Inline, Right Side Mount \$265

Part: See XRef. ADP: WEPM

Swing-away, 1" Offset, Right Side Mount \$265

Part: See XRef. ADP: WEPM



Retract4

Joystick mounting brackets continued on following page.

10. JOYSTICK MOUNTING BRACKETS CONTINUED

10.1c. Specialty Joystick Mounts

Select a bracket and side to mount it.

Retract4 Mount \$450

NE/NE+ Part: MEC138876. ADP: WEPM
 Q-Logic Part: MEC138860. ADP: WEPM

Full 180° of movement that translates into multi-position use, straight, right offset, left offset, or fold away. Joystick orientation varies as bracket is adjusted.

Multi-axis Mount \$100

NE/NE+ Part: ACC142528. ADP: WEPM
 Q-Logic Part: ACC131216. ADP: WEPM
 SAJ Part: ACC142740. ADP: WEPM

Allows height adjustment of joystick.

Stealth Swing-Away/Flip Down Bracket \$300

NE/NE+ Part: ST-TWBM480QL-Q. ADP: WEPM
 Q-Logic Part: ST-TWBM480QL-Q. ADP: WEPM
 SAJ Part: ST-TWBM480CJ-Q. ADP: WEPM

Allows the joystick to pivot to side of armrest, or flip-down if bracket is adjusted.



Specialty Mount Side

Left Side Mount
 Part: PTOINDV1225

Right Side Mount
 Part: PTOINDV1226

11. CONTOURED SEATING

Optional. Select only one seating system. Special order seat sizes can be requested by checking the seat type and writing in the seat size in the special order seat size.

11.1a. Captain Seating

1a. Depth Adjustable High Back with Limited Recline

Full length armrests must be selected with a depth adjustable high back.

18"Wx18"-20"D with 60" Lap Belt No Charge
 Part: SETASMB6895

22"Wx20"-22"D with 70" Lap Belt No Charge
 Part: SETASMB6897

20"Wx18"-20"D with 60" Lap Belt No Charge
 Part: SETASMB6896

11.1c. Special Order Contoured Seat

Complete information below. Contact inside sales specialist for pricing.

Type of Seat: _____ Width: _____ Depth: _____

Contoured seating section continued on following page.

13. HEADRESTS

Optional. Select a headrest pad and bracket. **Not available with contoured seating.**

13.1. Headrest Pads

Stealth 8" Comfort Plus (CP450)..... \$190
 Part: POS148484.


Stealth 10" Comfort Plus (CP250)(A) \$190
 Part: SETHEAD1057.

Stealth 14" Comfort Plus (CP150) \$190
 Part: INDPART2327.

Pediatric Comfort Plus (6"Wx4.5"H) (CPN350) \$190
 Part: SETHEAD1069.

Stealth Small Contoured (9.75"Wx3"H) (STL-OTSM250) \$80
 Part: SETHEAD1052.

Stealth Large Contoured (11"Wx6"H) (B)..... \$100
 (STL-OTLG150)
 Part: SETHEAD1053.



13.2. Headrest Brackets

TWB Removable Headrest Bracket (TWB485-1) (C) \$255
 Part: FRMASMB15999.

Multi-Axis Removable Headrest Bracket (HMO475P) (D)..... \$200
 Part: HDWASMB1084.

14. LAP BELTS

REQUIRED. Select a lap belt below.

14.1a. Standard Lap Belts

38" Lap Belt Part: ACCBELT1006	50" Lap Belt Part: ACCBELT1000	60" Lap Belt Part: ACCBELT1002	70" Lap Belt Part: ACCBELT1003	80" Lap Belt Part: ACCBELT1004	90" Lap Belt Part: ACCBELT1001
--	--	--	--	--	--

14.1b. Optional Lap Belts

Stealth Retractable Lap Belt (IFRB5350).....\$175
 Part: ST-IFRB5350

i-Fit 60" Single Pull, Push Button, Medium (7.5") Padded, 1.5" Belt (IFPB5213) \$65
 Part: ST-IFPB5213-Q

i-Fit 60" Single Pull, Push Button, Large Padded (9.5"), 2" Belt (IFPB5314) \$65
 Part: ST-IFPB5314-Q

i-Fit 60" Single Pull, Push Button, Extra Large Padded (11.5"), 2" Belt (IFPB5315) \$65
 Part: ST-IFPB5315-Q

60" Sure-Lok Lap Belt (Non-retractable)..... \$200
 Required with ISO 7176-19 Compliant Transit Option.
 Part: ACCBELT1016

17. LEGRESTS CONTINUED

17.3. 60°-90° Swing-Away Legrests

Adjustable in 5° increments. Select a legrest upper, a legrest extension, and type of footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes. Select legrest length in following section.

1. Swing-Away Legrest Upper

- P L R 60°-90° Swing-Away Legrest..... \$142.50 ea.
Pair part: RIG120062. Left: FRMASMB17558.
Right: FRMASMB17559
- P L R 60°-90° Swing-Away Legrest, Short..... \$142.50 ea.
Pair part: RIG136481. Left: RIG132607.
Right: RIG132609.

2. Lower Extension

- P L R 11"-13" Length
Part: See XRef
N/A with Multi-axis.
- P L R 13"-16.5" Length
Part: See XRef
N/A with Multi-Axis.
- P L R 11"-16.5" Length
Part: See XRef
N/A with Angle Adjustable.
- P L R 15.5"-19.5" Length
Part: See XRef

60°-90° Swing-Away Legrest Compatibility Matrix	Contoured Seat			Synergy Seat		
	Low	Medium	High	Low	Medium	High
Extruded, 13" - 16.5"	✓	✓	✓	!	!	✓
Extruded, 15.5" - 19.5"	!	✓	✓	X	!	!
Angle Adjustable, 13" - 16.5"	✓	✓	✓	!	!	✓
Angle Adjustable, 15.5" - 19.5"	!	✓	✓	X	!	!
Multi-Axis, 11" - 16.5"	✓	✓	✓	!	!	✓
Multi-Axis, 15.5" - 19.5"	!	✓	✓	X	!	!

11"-13" Extensions cause no interference at all seat-to-floor heights.

- This rigging configuration causes no interference issues at this seat-to-floor height.
- This rigging configuration may have less than 1.75" clearance when fully extended and may interfere with the ground or reduce obstacle climbing ability.
- This rigging configuration does not work with this seat-to-floor height.

3a. Extruded Footplate(s)

- P L R 7.25"Wx6"D Extruded Footplate..... Standard
Part pair: FRMASMB7847. Left: FRMASMB7861
Right: FRMASMB7862
- P L R 8.5"Wx6"D Extruded Footplate Standard
Part pair: FRMASMB7848. Left: FRMASMB7863
Right: FRMASMB7864
May cause interference on seat widths less than 18."

3b. Angle Adjustable Footplate(s) with Heel Loop(s)

¹The 6"x11" footplates may cause interference with seat widths less than 18". ²The 8"x11" footplates may cause interference with seat widths less than 20".

- P L R 4"Wx6"D \$87.50 ea.
Pair part: FRMASMB10564. Left: RIG133537. Right: RIG133560
- P L R 5"Wx6"D..... \$87.50 ea.
Pair part: FRMASMB10565. Left: RIG133538. Right: RIG133561
- P L R 5"Wx8"D \$87.50 ea.
Pair part: FRMASMB10566. Left: RIG133539. Right: RIG133562
- P L R 6"Wx8"D¹ \$87.50 ea.
Pair part: FRMASMB10567. Left: RIG133540. Right: RIG133563
- P L R 6"Wx11"D¹ \$87.50 ea.
Pair part: FRMASMB10568. Left: RIG133541. Right: RIG133564
- P L R 8"Wx11"D² \$87.50 ea.
Pair part: FRMASMB10569. Left: RIG133542. Right: RIG133565

3c. Multi-Axis Footplate(s) with Heel Loops(s)

¹The 6"x11" footplates may cause interference with seat widths less than 18". ²The 8"x11" footplates may cause interference with seat widths less than 20".

- P L R 4"Wx6"D \$260 ea.
Pair part: FRMASMB9456. Left: RIG133531. Right: RIG133552
- P L R 5"Wx6"D..... \$260 ea.
Pair part: FRMASMB9457. Left: RIG133532. Right: RIG133553
- P L R 5"Wx8"D \$260 ea.
Pair part: FRMASMB9458. Left: RIG133533. Right: RIG133554
- P L R 6"Wx8"D¹ \$260 ea.
Pair part: FRMASMB9459. Left: RIG133534. Right: RIG133555
- P L R 6"Wx11"D¹ \$260 ea.
Pair part: FRMASMB9460. Left: RIG133535. Right: RIG133558
- P L R 8"Wx11"D² \$260 ea.
Pair part: FRMASMB9461. Left: RIG133536. Right: RIG133559

4. 60°-90° Swing-Away Legrest Accessories

- Heel Loops (pair)..... \$35
Part: FRMASMB7873.
Standard on angle adjustable and multi-axis footplates.
- Wheel Bumper..... \$15 ea.
Part: See XRef
- Gel Neoprene Footrest Sleeve (pair) \$130
Part: INDPART2243. : K0108

Legrests section continued on following page.

17. LEGRESTS CONTINUED

17.4. Elevating Legrests

Select a legrest upper, a legrest extension, and type of footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes. Select legrest length in following section.

1. Elevating Legrests Upper

P L R Elevating Legrest Assembly \$275 ea.
 Pair part: FRMASMB6483. Left: RIG138012
 Right: RIG138013.
 Select lower extension and footplates

2. Lower Extension

P L R **11"-14" Length**
 Pair: FRMASMB4755
 Left: RIG138026
 Right: RIG138030
 For reinforced only

P L R **14"-17" Length**
 Pair: FRMASMB4756
 Left: RIG138033
 Right: RIG138035
 For reinforced only

P L R **12"-14.5" Length**
 Part: See XRef
 For angle adjustable only

P L R **14.5"-18" Length**
 Part: See XRef
 For angle adjustable only

Elevating Legrest Compatibility Matrix	Contoured Seat			Synergy Seat		
	Low	Medium	High	Low	Medium	High
Standard, 11"-14"	✓	✓	✓	✓	✓	✓
Standard, 14"-17"	✓	✓	✓	!	!	!
Angle Adjustable, 12"-14.5"	✓	✓	✓	✓	✓	✓
Angle Adjustable, 14.5"-18"	✓	✓	✓	!	!	!

This rigging configuration causes no interference issues at this seat-to-floor height.
 This rigging configuration may have less than 1.75" clearance when fully extended and may interfere with the ground or reduce obstacle climbing ability.

3a. Reinforced Cast Footplate(s)

P L R **7.25"Wx6"D Reinforced Footplate**..... Standard
 Pair part: FRMASMB4748. Left: RIG138055. Right: RIG138056

P L R **8.5"Wx6"D Reinforced Footplate** Standard
 Pair part: FRMASMB4747. Left: RIG138060. Right: RIG138062
 May cause interference on seat widths less than 18."

3b. Angle Adjustable Footplate(s)

*The 6"x11" footplates may cause interference with seat widths less than 18". *The 8"x11" footplates may cause interference with seat widths less than 20". 4"x6" and 5"x6" footplates are only available in pairs.

P L R **4"Wx6"D** \$87.50 ea.
 Pair parts: FRMASMB8812 x2.

P L R **5"Wx6"D**..... \$87.50 ea.
 Pair parts: FRMASMB8811 x2.

P L R **5"Wx8"D Footplates with Heel Loops** \$87.50 ea.
 Pair part: FRMASMB4735. Left: RIG137791. Right: RIG137792

P L R **6"Wx8"D¹ Footplates with Heel Loops** \$87.50 ea.
 Pair part: FRMASMB4738. Left: RIG137793. Right: RIG137795

P L R **6"Wx11"D¹ Footplates with Heel Loops** \$87.50 ea.
 Pair part: FRMASMB4741. Left: RIG137796. Right: RIG137797

P L R **8"Wx11"D² Footplates with Heel Loops** \$87.50 ea.
 Pair part: FRMASMB4744. Left: RIG137798. Right: RIG137799

4. Elevating Legrest Accessories

Heel Loops (pair)..... \$35 Available on angle adjustable footplates.
 Part: FRMASMB7873

Wheel Bumper..... \$15 ea.
 Part: See XRef


18. SET LEGREST LENGTH


REQUIRED with matching left and right legrest configuration. Non-matching or individual legrest configurations lengths are not set by the manufacturer.

7" L (#3470)	8.5" L (#3473)	10" L (#3476)	11.5" L (#3479)	13" L (#3482)	14.5" L (#3485)	16" L (#3488)	17.5" L (#3491)	19" L (#3494)	20.5" L (#3497)	22" L (#3669)
7.5" L (#3471)	9" L (#3474)	10.5" L (#3477)	12" L (#3480)	13.5" L (#3483)	15" L (#3486)	16.5" L (#3489)	18" L (#3492)	19.5" L (#3495)	21" L (#3498)	22.5" L (#3670)
8" L (#3472)	9.5" L (#3475)	11" L (#3478)	12.5" L (#3481)	14" L (#3484)	15.5" L (#3487)	17" L (#3490)	18.5" L (#3493)	20" L (#3496)	21.5" L (#3600)	

Seat cushion section continued on following page.

19. SEAT CUSHION

Simplicity G\$310-\$370 


TRU-Comfort 2 SPP\$515-\$575 


Stealth Spectrum Foam SPP Cushion

Skin protection and positioning with Coolcore, 4-way stretch, breathable cover. Includes urethane liner.

Stealth Spectrum Gel SPP Cushion

Skin protection and positioning with Coolcore, 4-way stretch, breathable cover. Includes urethane liner.

Spectrum Foam SPP\$580-\$630 


Spectrum Gel SPP\$580-\$650 


Stealth Solution SPP Cushion

Skin protection and positioning with Coolcore, 4-way stretch, breathable cover. Includes urethane liner.

Stealth Essence SPP Cushion

Skin protection and positioning with Coolcore, 4-way stretch, breathable cover.

Solution SPP\$580-\$645 


Essence SPP\$545-\$705 
Essence Plus SPP\$800


Stealth Glacial SP Cushion

Skin protection with Coolcore, 4-way stretch, breathable cover.

Stealth Premiere P Cushion


Positioning with Coolcore, 4-way stretch, breathable cover.

Glacial SP\$265-\$530 

Premiere P\$585 

Stealth Zen SP Cushion

Skin protection and positioning with Coolcore, 4-way stretch, breathable cover.

Zen SP\$455-\$765 

19.2. Seat Cushion Dimensions

REQUIRED with cushion model selection. Select to match seating or pick non-matching dimensions. See respective cushion order form for pricing and available sizes.

Match Seat Width and Depth

Cushion Width

14" W 16" W 17" W 18" W 19" W 20" W 22" W

Cushion Depth

14" D 16" D 17" D 18" D 19" D 20" D 22" D

19.3. Seat Cushion Substitute Fluid Proof Liners

Fluid Proof Urethane Liner\$45-\$50
Part: URL1709070-XXYY
Not available with Simplicity, Essence, or TRU-Comfort 2.

Fluid Proof Polyurethane Liner\$70
Part: PUL1709030-XXYY
Not available with Simplicity.

19.4. Seat Cushion Substitute Cover

An additional cover may be ordered or the standard cover may be substituted with a cover by checking the option in subsection 2. The "*prefix" of the part number will be the part prefix of the cushion model selected above. For example, SPCGSPP-NBLK1816 will be an additional coolcore cover for a Spectrum Gel cushion.

Extra Stealth Tek 2 Cover\$90
Part: *prefix-TBLKXXYY
Not available with Simplicity, TRU-Comfort 2, Essence, or Premiere.

Extra Mesh Cover\$90
Part: *prefix-MBLKXXYY
Not available with Simplicity or TRU-Comfort 2.

20. ACCESSORIES

Optional. Multiple accessories can be mounted along the various tracks but some interference may occur. Accessories may need to be installed by the provider.

Contoured Seat Cup Holder\$20
Part: ACCASMB2287. Part with 12" lift: ACCASMB2288

Flag Holder\$125
Contoured seat part: ACCASMB1930
Synergy/Flex part: ACCASMB2415.

Synergy/Flex Seat Cup Holder\$20
Part: ACCASMB2634
Not available on TRU-Balance Flex seating without two-post armrests.

Saddle Bag\$15
Part: DWR1020D010

Oxygen Holder for Contoured Seat\$110
Part: ACCASMB2292.

Quantum Backpack\$65
Part: ACC190585

Oxygen Holder for Flex Seat\$110
Part: ACC125006

Accessory Bar\$90
12"W-16"W part: MEC145955. 16"W-20"W part: MEC143975.
20"W-22"W part: MEC146018

Walker Holder\$150
Contoured seat part: ACCASMB1912. Synergy/Flex part: FRMASMB7031

Large, Black, Weather Cover\$240
Part: ACCCOVR1021

Cane/Crutch Holder\$195
Contoured seat part: ACCASMB2294. HCPCS: E2207
Synergy part: FRMASMB4822. Flex part: ACC152515.

Rear Basket\$99
Contoured seat part: ACCASMB1914
Synergy part: FRMASMB4821