



**Pride Mobility Products Co.**

5096 South Service Rd, Beamsville, ON, L0R 1B3

Phone: 888-570-1113 | Fax: 866-514-1303 | Email: [productorders@pridemobility.com](mailto:productorders@pridemobility.com)

# Rival® Order Form with Contoured and TRU-Balance Flex Seating

Account Number: \_\_\_\_\_ Date: \_\_\_\_\_

Provider Name: \_\_\_\_\_

Contact: \_\_\_\_\_

Phone: \_\_\_\_\_ Fax: \_\_\_\_\_

Email: \_\_\_\_\_

PO Number: \_\_\_\_\_

Marked for: \_\_\_\_\_

Ship to Address: \_\_\_\_\_

City: \_\_\_\_\_ Prov: \_\_\_\_\_ PC: \_\_\_\_\_



## INTRODUCTION

This form is interactive when viewed with [Adobe Acrobat Reader](#). Complete the form by placing checks in the desired boxes and provide information in the interactive fields. Buttons shown at the bottom of the form may be utilized to print or submit the order form through a desktop email application. To email via a web-based application, please 'Save As' and attach the PDF to your email.

This order form contains a large variety of options to fulfill various patient needs. Descriptions and section notes, such as optional or required, are included to help you complete your order. Other part numbers can be found on the quote acknowledgment after this order form is submitted. Please contact Quantum® Sales at 888-570-1113 if assistance is needed. Send the completed order form by fax (866-514-1303) or email ([productorders@pridemobility.com](mailto:productorders@pridemobility.com)). Incomplete forms may delay the quote or order. Customer service will contact you if the order is incomplete or if there are compatibility issues. If special order requests are needed, be sure that the Patient Information section is completed

All prices are MSRP. Prices, specifications, part numbers, and availability are subject to change without notice. Prices and part numbers as shown are only available when configured as a complete power chair. Please contact technical service for accurate parts ordering.


## PATIENT INFORMATION

Weight: \_\_\_\_\_ Height: \_\_\_\_\_

A. Shoulder Width: \_\_\_\_\_ D. Max Sitting Height: \_\_\_\_\_ G. Thigh Depth: \_\_\_\_\_

B. Chest Width: \_\_\_\_\_ E. Shoulder Height: \_\_\_\_\_ H. Lower Leg Length: \_\_\_\_\_

C. Hip Width: \_\_\_\_\_ F. Axilla Height: \_\_\_\_\_ I. Elbow Height: \_\_\_\_\_



## 1. BASE MODEL

**REQUIRED.** Select a group 3 power base. 300 lb. weight capacity with standard 6 mph programmable motors and 8 amp off board charger. Dimensions of base without legrests: 25"Wx36"L.

Quantum Rival 3S-SS ..... \$7,240  
Part: QRIVAL, ADP Adult: WP0000968, ADP Pediatric: WP0000969

**2. BASE COLOR**











**REQUIRED.** Select one shroud color. Colors shown here may differ from actual product depending on monitor/printer calibration used.

**2.1. Match Colors**

Match All Color Options to Selection Below

Part: N/A. Any subsequent color selections, if available, will be matched to color chosen below when the order is placed. Please note that some color selections, such as joystick armor, may have a MSRP.

**2.2. Base Shroud Color**

- |  |  |   |  |
|--|--|---|--|
|  <b>Red Corvette</b><br>Part: DGN153201-PRD1  |  <b>Night Train</b><br>Part: DGN1801101-PBK3      |  <b>Deep Purple</b><br>Part: DGN153201-PPR5    |  <b>Green Machine</b><br>Part: DGN1801101-PGN11 |
|  <b>Electric Blue</b><br>Part: DGN153201-PBU9 |  <b>Pretty in Pink</b><br>Part: DGN153201-PPK3    |  <b>Grey Street</b><br>Part: DGN153201-PSL3    | <b>To Be Determined</b><br>Part: N/A   |
|  <b>Back in Black</b><br>Part: DGN153201-PBK1 |  <b>Evolution Orange</b><br>Part: DGN1801101-POR7 |  <b>Midnight Blue</b><br>Part: DGN153201-PBU10 |  |

**3. DRIVE MOTORS**

**REQUIRED.**

- |   |   |
|---|---|
| 6 mph Drive Motors ..... Standard<br>Part: MOT132462  | 8 mph Song Motors ..... \$160<br>Part: MOT132468<br>Not available with NE 4 Key joystick. |
| 6 mph, Accu-Trac Motor ..... \$630<br>Part: MOT134093.<br>Requires Q-Logic 2 electronics. Delivers enhanced traction and tracking, especially with switched drive controls. |   |

**4. DRIVE WHEELS**

**REQUIRED.** Select a type of drive wheels. Knobby tires increase overall base width by 1.5".

**4.1. Drive Wheels**

- |   |   |
|---|---|
| 3"x14" Flat Free Tires ..... Standard<br>Part: WHLASMB2246 x2                               | 3"x14" Pneumatic Tires ..... No Charge<br>Part: WHLASMB2247 x2                              |
| 4"x14" Knobby, Flat Free Tires ..... \$120<br>Extended Traction Tires<br>Part: WHL122174 x2 | 4"x14" Knobby, Pneumatic Tires ..... \$120<br>Extended Traction Tires<br>Part: WHL122175 x2 |

**4.2. Drive Wheel Accent Color**

Optional. Not available with knobby tires. Select with WHL1703421 above. Sold in Pairs (parts below will be quantity of 2).

- |   |  |
|---|--|
| Match to Base Color ..... \$42.50 ea.<br>Part: Configuration selection, see specific part below | Non-matching Color ..... \$42.50 ea.<br>Part: Configuration selection, see specific part below |
| <b>Red Corvette</b><br>Part: DGN164245-PRD1   | <b>Deep Purple</b><br>Part: DGN164245-PPR5   |
| <b>Electric Blue</b><br>Part: DGN164245-PBU9  | <b>Grey Street</b><br>Part: DGN164245-PSL3   |
| <b>Back in Black</b><br>Part: DGN164245-PBK1  | <b>Midnight Blue</b><br>Part: DGN164245-PBU10  |
| <b>Night Train</b><br>Part: DGN164245-PBK3  | <b>Green Machine</b><br>Part: DGN164245-PGN11  |
| <b>Pretty in Pink</b><br>Part: DGN164245-PPK3   | <b>To Be Determined</b><br>Follow-up required.   |
| <b>Evolution Orange</b><br>Part: DGN164245-POR7   |  |

**5. FRONT CASTER WHEELS**

**REQUIRED.** For optimal ride quality, pneumatic casters are recommended.

- |   |   |  |
|---|---|--|
| 9" Pneumatic ..... Standard<br>Part: FRMASMB17634 | 9" Flat Free ..... \$50<br>Part: FRMASMB17635 | 8" Flat Free ..... \$50<br>Part: MEC141453 |
|---|---|--|

**6. SUSPENSION**

**REQUIRED.**

- |   |   |
|---|---|
| Standard Suspension ..... Standard<br>Part: PTOINDV3545 | Heavy Suspension ..... No Charge<br>Part: MEC120220 |
|---|---|

**7. FULL LIGHT PACKAGE**

**Optional.** Requires NE+ 8 Key or Q-Logic 2 Joystick for Full Light Package.

- |   |
|---|
| Base Mounted, LED, Light Package ..... \$770<br>Part: ELEASMB7431 |
|---|

**8. CONTOURED SEAT, SEAT TO FLOOR HEIGHT**

**REQUIRED with contoured seating.**

- |   |   |  |
|---|---|--|
| <b>Lowest Height</b><br>Part: PTOINDV3558 | <b>Medium Height</b><br>Part: PTOINDV3559 | <b>Highest Height</b><br>Part: PTOINDV3560 |
|---|---|--|

**9. ELECTRONICS**

**REQUIRED.** Select a type of drive control from the following sections. The No Joystick option is used to omit the joystick and utilize another control, such as a specialty control at the end of this form or found on the [Alternative Controls and Electronics Order Forms](#).

**9.1a. Non-Expandable Controls**

Select a joystick below or omit selection if another controller type is desired. Non-expandable does not allow supplementary drive controllers.

NE 4-Key Joystick.....Standard  
 Part: CTLDC1554  
 Available with Tilt Only through Toggle.

No Joystick, Setup for NE+ Joystick ..... No Charge  
 Part: Configuration selection

NE+ 6-Key Joystick ..... No Charge  
 Part: CTLDC1560  
 Control up to 2 actuators through Joystick.



**9.1b. Expandable Controls - Quantum Intelligent Systems™**

Select expandable controls or omit selection if another controller type is desired. Expandable controls allow multiple or specialized controllers. This sections allows a selection of common configurations.

**2. Q-Logic 2 Joystick**

Q-Logic 2 EX Joystick ..... \$620  
 Part: CTL123217 ADP: WEX1

Stand Alone Joystick with ..... \$630  
 External Charger Port  
 Part: ELEASMB5205

Q-Logic 2 EX Joystick for Full Lights ..... \$620  
 Part: CTL124653 ADP: WEX1

Up to 3 drive profiles and 1 seating profile when used alone. Use in conjunction with Enhanced Display for up to 5 profiles and auxiliary functions such as Bluetooth Mouse and IR. Includes charger port mounting bracket and multiplier box harness.

No Joystick, Setup for Q-Logic  
 Part: Configuration selection.



**3. Q-Logic 2 Joystick Shroud**

**REQUIRED with Q-Logic 2 EX Joystick selection above.** Reference colors shown on page 2.

Black, Soft Touch ..... No Charge  
 Part: DGN124307

Match Joystick Shroud to Base Color ..... Standard  
 Part: PTOINDV3421

Non-matching Joystick Shroud..... No Charge  
 Part: PTOINDV3422  
 Select color below

Red Corvette  
 Part: DGN124336-PRD1

Night Train  
 Part: DGN124336-PBK3

Deep Purple  
 Part: DGN124336-PPR5

Green Machine  
 Part: DGN124336-PGN11

Electric Blue  
 Part: DGN124336-PBU9

Pretty in Pink  
 Part: DGN124336-PPK3

Grey Street  
 Part: DGN124336-PSL3

Back in Black  
 Part: DGN124336-PBK1

Evolution Orange  
 Part: DGN124336-POR7

Midnight Blue  
 Part: DGN124336-PBU10



**9.2. SCIM Module**

Optional.

SCIM Module Kit for Q-Logic 2..... \$730  
 Part: CTL168055

Available with Q-Logic 2 EX Joystick selection above. Allows use of joystick display for specialty devices and offers up to 5 profiles to control drive, seat, and auxiliary functions, excluding sip-n-puff and single switch functions. A switch is required to operate the input device attached to the SCIM module.



**10. JOYSTICK HANDLES**

**REQUIRED.** The NE Plus and Q-Logic 2 EX Joysticks are standard with a carrot handle (See A in image). The optional handle will be mounted to the joystick and the standard handle will be shipped with the unit.

Mushroom Handle (B) ..... \$120  
 Part: ACCASMB1860.

Extended I-Handle (G) ..... \$140  
 Part: ACCASMB1857.

Chin Cup (C)..... \$85  
 Part: ACCASMB186190

Bodypoint 3" Goal Post Handle ..... \$180  
 Part: PLSKNCP1057.

T-Bar (D) ..... \$160  
 Part: ACCASMB1858.

Bodypoint 4" Goal Post Handle ..... \$180  
 Part: PLSKNCP1058.

Medium Ball Handle (E)..... \$70  
 Part: ACCASMB1859.

Soft Ball Handle (F)..... \$120  
 Part: ACCASMB1856.



**11. JOYSTICK MOUNTING BRACKETS**

**Optional.** Select a mount to be used with a joystick selected in previous section.


Omit Joystick Mounting Bracket

**11.1a. Fixed, In-line Joystick Bracket**

Fixed, In-line, Bracket ..... Standard  
 Q-Logic Part: FRMASMB10617. NE/NE+ Part: FRMASMB12744

Left Side Mount  
 Part: PTOINDV1225

Right Side Mount  
 Part: PTOINDV1226



**11.1b. Swing-Away (Retractable) Joystick Bracket**


Allows the joystick to retract to the side of the armrest while maintaining the same joystick orientation. "Inline bracket" will be offset 1" when used with tube style arms. Select "Offset Bracket" to have the joystick inline with armpad on tube style armrests.

Swing-away, Inline, Left Side Mount ..... \$265  
 Part: See XRef. ADP: WEPM

Swing-away, 1" Offset, Left Side Mount ..... \$265  
 Part: See XRef. ADP: WEPM

Swing-away, Inline, Right Side Mount ..... \$265  
 Part: See XRef. ADP: WEPM

Swing-away, 1" Offset, Right Side Mount ..... \$265  
 Part: See XRef. ADP: WEPM



Retract4

**11.1c. Specialty Joystick Mounts**

Select a bracket and side to mount it.

Retract4 Mount ..... \$450  
 NE/NE+ Part: MEC138876. ADP: WEPM  
 Q-Logic Part: MEC138860. ADP: WEPM  
 Full 180° of movement that translates into multi-position use, straight, right offset, left offset, or fold away. Joystick orientation varies as bracket is adjusted.

Stealth Swing-Away/Flip Down Bracket ..... \$300  
 NE/NE+ Part: ST-TWBM480QL-Q. ADP: WEPM  
 Q-Logic Part: ST-TWBM480QL-Q. ADP: WEPM  
 SAJ Part: ST-TWBM480CJ-Q. ADP: WEPM  
 Allows the joystick to pivot to side of armrest, or flip-down if bracket is adjusted.

Multi-axis Mount ..... \$100  
 NE/NE+ Part: ACC142528. ADP: WEPM  
 Q-Logic Part: ACC131216. ADP: WEPM  
 SAJ Part: ACC142740. ADP: WEPM  
 Allows height adjustment of joystick.



Multi-Axis Joystick Mount



Swing-Away/Flip Down Bracket (Stand Alone Hardware Shown)

**Specialty Mount Side**

Left Side Mount  
 Part: PTOINDV1225

Right Side Mount  
 Part: PTOINDV1226

**13. CONTOURED SEATING**

**Optional.** Select only one seating system. Special order seat sizes can be requested by checking the seat type and writing in the seat size in the special order seat size.

**13.1a. Captain Style Seating**

**1a. Depth Adjustable High Back with Limited Recline**

Full length armrests must be selected with a depth adjustable high back.

18"Wx18"-20"D with 60" Lap Belt ..... No Charge  
 Part: SETASMB6895

22"Wx20"-22"D with 70" Lap Belt ..... No Charge  
 Part: SETASMB6897

20"Wx18"-20"D with 60" Lap Belt ..... No Charge  
 Part: SETASMB6896

**1b. High Back with Limited Recline**

16"Wx16"D ..... No Charge  
 Part: SETASMB2770

16"Wx20"D ..... No Charge  
 Part: SETASMB2772

18"Wx16"D ..... No Charge  
 Part: SETASMB2775

16"Wx18"D ..... No Charge  
 Part: SETASMB2771

16"Wx22"D ..... No Charge  
 Part: SETASMB4887

20"Wx16"D ..... No Charge  
 Part: SETASMB2780

**1c. Medium Back**

16"Wx16"D ..... No Charge  
 Part: SETASMB2488

16"Wx18"D ..... No Charge  
 Part: SETASMB2560

18"Wx16"D ..... No Charge  
 Part: SETASMB2159

18"Wx18"D ..... No Charge  
 Part: SETASMB2551

**13.1b. Contoured Seat with Solid Seat Pan**

**1a. High Back with Limited Recline and Solid Seat Pan**

16"Wx16"D ..... No Charge  
 Part: SETASMB3178

18"Wx16"D ..... No Charge  
 Part: SETASMB3181

20"Wx16"D ..... No Charge  
 Part: SETASMB3183

16"Wx18"D ..... No Charge  
 Part: SETASMB3179

18"Wx18"D ..... No Charge  
 Part: SETASMB2848

20"Wx18"D ..... No Charge  
 Part: SETASMB3184

16"Wx20"D ..... No Charge  
 Part: SETASMB3180

18"Wx20"D ..... No Charge  
 Part: SETASMB3182

20"Wx20"D ..... No Charge  
 Part: SETASMB3185

Contoured seating section continued on following page.

**13. CONTOURED SEATING CONTINUED**

**1b. Medium Back with Solid Seat Pan**

16"Wx16"D ..... No Charge  
Part: SETASMB3112

16"Wx18"D ..... No Charge  
Part: SETASMB3113

18"Wx16"D ..... No Charge  
Part: SETASMB2492

18"Wx18"D ..... No Charge  
Part: SETASMB3114

**13.1c. Special Order Contoured Seat**

Complete information below. Contact inside sales specialist for pricing.

Type of Seat: \_\_\_\_\_ Width: \_\_\_\_\_ Depth: \_\_\_\_\_

**15. TRU-BALANCE FLEX SEATING SYSTEM**

**Optional.** Seat width and depth can be adjusted in 1" increments from 15"W-22"W and 16"D-20"D. Maximum seat size 22"W x 20"D.

**15.1. TRU-Balance Flex Frame Selection**

TRU-Balance Flex Seating System ..... %650  
Part: PTOINDV3853

**15.2. TRU-Balance Flex Seat Dimension - Width**

**REQUIRED with TRU-Balance Flex seating.** Maximum seat width is 22". Part number is SETWDTH plus number following # below.

15" W (#1005)	16" W (#1006)	17" W (#1007)	18" W (#1008)	19" W (#1009)	20" W (#1010)	21" W (#1011)	22" W (#1012)
------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------

**15.3. TRU-Balance Flex Seat Dimension - Depth**

**REQUIRED with TRU-Balance Flex seating.** Maximum seat depth is 20". Part number is SETDPH plus number following # below.

15" D (#1005)	16" D (#1006)	17" D (#1007)	18" D (#1008)	19" D (#1009)	20" D (#1010)
------------------	------------------	------------------	------------------	------------------	------------------

**15.4. TRU-Balance Flex Back Options**

**REQUIRED with TRU-Balance Flex seating.** See the Synergy and TRU-Comfort back availability chart for stock order back size ranges. Other size backs may be available with special order and extended lead time.

**1. Back Selection**

Synergy Back Cushion ..... \$940  
Part: See XRef

Sport Back Cushion ..... No Charge  
Part: See XRef

TRU-Comfort 2, 20" Width or less ..... \$835  
Part: See XRef.

Omit Back ..... No Charge  
Part: PTOINDV1287

		Back Width										
		15"	16"	17"	18"	19"	20"	21"	22"			
Overall Height	16"	-	P	T	P	T	P	T	P	T	P	T
	17"	S	S	P	T	S	P	T	S	P	T	S
	18"	S	S	P	T	S	P	T	S	P	T	S
	19"	S	S	P	T	S	P	T	S	P	T	S
	20"	S	S	P	T	S	P	T	S	P	T	S
	21"	S	S	P	T	S	P	T	S	P	T	S
	22"	S	S	P	T	S	P	T	S	P	T	S
	23"	S	S	P	T	S	P	T	S	P	T	S
	24"	S	S	P	T	S	P	T	S	P	T	S
	25"	S	S	P	T	S	P	T	S	P	T	S

S Synergy Back   
 P Sport Back   
 T TRU-Comfort Back

**TRU-Balance Flex Seat Back Availability Chart**

**2. Back Cane Selection**

Angle Adjustable Back Cane ..... \$195 Part: MEC140266	4" Mesomorph (2" per side) Back Cane ..... \$205 Part: MEC1804739
Flip Forward Back Cane ..... \$195 Part: MEC140272 Not available with oxygen holder, or cane mounted and cantilever armrests.	2" Endomorph (1" per side) Back Cane ..... \$205 Part: MEC141455
2" Mesomorph (1" per side) Back Cane ..... \$205 Part: MEC1804736	4" Endomorph (2" per side) Back Cane ..... \$205 Part: MEC141456
	24"H Flip Forward Flex Cane ..... Standard

**3. TRU-Balance Flex Push Handles**

Push Handles ..... Standard Part: Configuration selection with back cane height	Omit Push Handles ..... \$125 Part: Configuration selection with back cane height Not available with Sling Back. Minimum cane height is 16".
--	--

**4. TRU-Balance Flex Accessory Bar**

Accessory Bar ..... \$90  
12"W-16"W Part: MEC145955. 16"W-20"W Part: MEC143975.

TRU-Balance Flex seating section continued on following page.



**15. TRU-BALANCE FLEX SEATING SYSTEM CONTINUED**

**5. TRU-Balance Flex Overall Back Height**

Refer to the Synergy and TRU-Comfort back availability chart for stock order back size ranges. *Parts: See XRef, configuration selection with back selection.*

12" High Back	15" High Back	18" High Back	21" High Back	24" High Back
13" High Back	16" High Back	19" High Back	22" High Back	25" High Back
14" High Back	17" High Back	20" High Back	23" High Back	

**6. TRU-Balance Flex Back Cane Height**

Omit this selection if a TRU-Comfort back is selected. 16" and 17" back canes are not available with a Synergy Back and cane mounted, flip up armrests. 14" and 15" back canes are not available without push handles. *Parts: See XRef, configuration selection with push handles selection.*

14" High Canes	17" High Canes	20" High Canes	23" High Canes
15" High Canes	18" High Canes	21" High Canes	24" High Canes
16" High Canes	19" High Canes	22" High Canes	25" High Canes

**15.5. TRU-Balance Flex Armrests**

**REQUIRED with TRU-Balance Flex seating.** Left (L) and right (R) boxes may be checked to configure different size arms within one section arm type.

**1a. TRU-Balance Flex Seat, Drop-In Armrests**

Full Length, Drop-In Armrests..... No Charge  
*Left part: ARM1803253. Right part: ARM156742*

Desk Length, Drop-In Armrests ..... No Charge  
*Left part: ARM1803254. Right part: ARM156743*

**1b. Flip Back, Cane Mounted Armrests**

May interfere with swing-away laterals.

Full Length, Flip Back Armrests (pair)..... No Charge  
*Left part: FRMASMB7038. Right part: FRMASMB7039*

Desk Length, Flip Back Armrests (pair) ..... No Charge  
*Left part: FRMASMB7040. Right part: FRMASMB7041*

**1c. Two-Post, Flip Back, Height Adjustable, Removable Armrests**

Not available with a depth less than 16". Required with flip forward back and laptrays.

Full Length, Two-Post Armrests.....\$145 ea.  
*Part: See XRef.*

Desk Length, Two-Post Armrests .....\$145 ea.  
*Part: See XRef.*

**1d. Single Post, Height Adjustable, Removable Armrests**

Height adjustment from 11"-17" in half-inch increments.

Full Length, Single Post Armrests .....\$145 ea.  
*Left part: FRMASMB8594. Right part: FRMASMB8595.*

Desk Length, Single Post Armrests .....\$145 ea.  
*Left part: FRMASMB8592. Right part: FRMASMB8593.*

**15.6. TRU-Balance Flex Armpads**

**REQUIRED with armrest selection above.** Waterfall armpads may cause interference with the seat back. Armpads are configured as pairs.

**1a. Straight Armpads**

Straight, Full Length Left and Right ..... No Charge  
*Part: ARM134275*

Straight, Desk Length Left, Full Length Right ..... No Charge  
*Part: ARM134275*

Straight, Full Length Left, Desk Length Right ..... No Charge  
*Part: ARM134281*

Straight, Desk Length Left and Right..... No Charge  
*Part: ARM134254*

**1b. Waterfall Armpads**

Waterfall, Full Length Left and Right..... No Charge  
*Part: ARM134275*

Waterfall, Desk Length Left, Full Length Right ..... No Charge  
*Part: ARM134275*

Waterfall, Full Length Left, Desk Length Right ..... No Charge  
*Part: ARM134281*

Waterfall, Desk Length Left and Right ..... No Charge  
*Part: ARM134254*

**1c. Falcon Gel Armpads**

Falcon Gel, Full Length Left and Right ..... \$300  
*Part: ARM136438*

Falcon Gel, Desk Length Left, Full Length Right..... \$250  
*Part: ARM136439*

Falcon Gel, Full Length Left, Desk Length Right..... \$250  
*Part: ARM136442*

Falcon Gel, Desk Length Left and Right ..... \$200  
*Part: ARM136437*

**16. TRU-BALANCE 3 HEADRESTS**

**Optional.** Select a headrest pad and bracket.

**16.1. Headrest Pads**

Stealth 8" Comfort Plus (CP450)..... \$190  
*Part: POS148484.*


Stealth Small Contoured (9.75"Wx3"H)..... \$80  
*(STL-OTS250)  
 Part: SETHEAD1052.*

Stealth 10" Comfort Plus (CP250)(A) ..... \$190  
*Part: SETHEAD1057.*

Stealth Large Contoured (11"Wx6"H) (B)..... \$100  
*(STL-OTLG150)  
 Part: SETHEAD1053.*

Stealth 14" Comfort Plus (CP150) ..... \$190  
*Part: INDPART2327.*

Pediatric Comfort Plus (6"Wx4.5"H) ..... \$190  
*(CPN350)  
 Part: SETHEAD1069.*



**16.2. Stealth Headrest Brackets**

TWB Removable Headrest Bracket (TWB485-1) (C) ..... \$255  
*Part: FRMASMB15999.*

Multi-Axis Removable Headrest Bracket (HMO475P) (D)..... \$200  
*Part: HDWASMB1084.*

**17. LAP BELTS**

**REQUIRED.** Select a lap belt below.

**17.1a. Standard Lap Belts**

<b>38" Lap Belt</b> Part: ACCBELT1006	<b>50" Lap Belt</b> Part: ACCBELT1000	<b>60" Lap Belt</b> Part: ACCBELT1002	<b>70" Lap Belt</b> Part: ACCBELT1003	<b>80" Lap Belt</b> Part: ACCBELT1004	<b>90" Lap Belt</b> Part: ACCBELT1001
--	--	--	--	--	--

**17.1b. Optional Lap Belts**

Not available with contoured or TRU-Balance Flex seating.

<b>i-Fit 60" Single Pull, Push Button, Medium (7.5") Padded, 1.5" Belt (IFPB5213)</b> Part: ST-IFPB5213-Q	\$65	<b>i-Fit 60" Single Pull, Push Button, Extra Large Padded (11.5"), 2" Belt (IFPB5315)</b> Part: ST-IFPB5315-Q	\$65
<b>i-Fit 60" Single Pull, Push Button, Large Padded (9.5"), 2" Belt (IFPB5314)</b> Part: ST-IFPB5314-Q	\$65		

**20. LEGRESTS**

**Optional.** All efforts will be made to accommodate legrest length selection, but minor changes may need to be made to meet shipping guidelines. Some adjustments may be needed upon receipt of the unit. Reference the section's compatibility matrix for legrest interference.

**20.1. Center Mount Foot Platform**

If power articulating foot platform is selected, be sure power legrest control has been selected in earlier section. Select legrest length in following section. **Not available with contoured or TRU-Balance Flex seating.**



**1. Foot Platform Assembly**

**Manually Adjustable Center Mount Foot Platform... \$265**  
12"-16" width part: FRMASMB17944. 14"-18" width part: FRMASMB17945. 16"-20" width part: FRMASMB17946.  
HCPCS: K0040.

**2. Lower Extension and Angle Adjustable Footplates**

<b>8"-12" Extension, 5"x6" Footplates</b> Part: RIG152577 Not available with power articulating platform.	<b>10"-16" Extension, 5"x6" Footplates</b> Part: RIG150206	<b>13"-19" Extension, 5"x6" Footplates</b> Part: RIG151505
<b>8"-12" Extension, 5"x8" Footplates</b> Part: RIG152578 Not available with power articulating platform.	<b>10"-16" Extension, 5"x8" Footplates</b> Part: RIG150203	<b>13"-19" Extension, 5"x8" Footplates</b> Part: RIG151506
	<b>10"-16" Extension, 5"x11" Footplates</b> Part: RIG150212	<b>13"-19" Extension, 5"x11" Footplates</b> Part: RIG151508

Center Mount Foot Platform Compatibility Matrix	Synergy Seat			<input checked="" type="checkbox"/> No interference issues at this seat-to-floor height. <input type="checkbox"/> May have less than 1.75" clearance when fully extended and may interfere with the ground or reduce obstacle climbing ability.
	Low	Medium	High	
10"-16" Length Extension	!	✓	✓	
13"-19" Length Extension	!	!	!	

**3. Center Mount Calf Support**

\*Calf pads are standard with a power articulating foot platform. Select calf pads covers for additional charge. Calf pad and panel covers include 1/2" of soft visco foam.

<b>Small Calf Pads (pair) 6"Wx6"H</b> ..... \$210* Part: RIGASMB7120018	<b>Small Calf Panel 18"Wx8"H</b> ..... \$210 Part: FRMASMB7110012
<b>Large Calf Pads (pair) 8" taper to 6"Wx8"H</b> ..... \$210* Part: RIGASMB7120019	<b>Large Calf Panel 24"Wx8"H</b> ..... \$265 Part: RIG123824 Only available on widths of 20" or greater.

**4. Modular Footplate Extension Components**

Optional. Assembly required. May require seating configuration adjustments to avoid front caster interference.

<b>Width Extension for 8"D Footplates</b> ..... \$55 Part: RIG152517	<b>4" Outside Edge (pair)</b> ..... \$55 Part: RIG152510 Only available with 6"D footplates and 8"D footplate extension.
<b>Width Extension for 11"D Footplates</b> ..... \$55 Part: RIG152518	<b>6" Outside Edge (pair)</b> ..... \$55 Part: RIG152516 Available with 8"D and 11"D footplates and 11"D footplate extension.
<b>2" Bell-Shaped Length Extension</b> ..... \$55 Part: RIG152519 Not compatible with width extensions. Can be used with 6" outside edge with 11" footplates.	

Legrests section continued on following page.

**20. LEGRESTS**

**5. Center Mount Foot Platform Additional Accessories**

- Reversed Dartex Calf Pad Cover (pair) ..... \$150  
Small part: RIG133717. Large part: RIG133719
- Reversed Dartex Panel Cover (pair) ..... \$150  
Small part: RIG134673. Large part: RIG134744
- Wheel Bumpers (pair)..... \$15 ea.  
Part: WHLASMB1578 x2
- Footplate Connector ..... \$55  
Part: RIGASMB7120023

**20.3. Pediatric, Swing-Away Legrests**

Not available with TRU-Balance Flex seating. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes. Select legrest length in following section.

**1. Pediatric Legrests Upper**

- 90°, Pediatric, Swing-Away Legrest (7"-10" Length Extension) ..... \$142.50 ea.  
Pair part: FRMASMB8808. Left: FRMASMB8809. Right: FRMASMB8810.

**2. Angle Adjustable Footplate(s) with Heel Loops**

- 4"Wx6"D Angle Adjustable ..... \$87.50 ea.  
Pair part: FRMASMB8812. Left: RIG133537. Right: RIG133560
- 5"Wx6"D Angle Adjustable ..... \$87.50 ea.  
Pair part: FRMASMB8811. Left: RIG133538. Right: RIG133561

**3. Accessories**

- Wheel Bumper(s) ..... \$15 ea.  
Pair part: WHL169508. Left: WHL169510. Right: WHL169511

**20.4. 60°-90° Swing-Away Legrests**

Adjustable in 5° increments. Select a legrest upper, a legrest extension, and type of footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes. Select legrest length in following section.

**1. Swing-Away Legrest Upper**

- 60°-90° Swing-Away Legrest ..... \$142.50 ea.  
Pair part: RIG120062. Left: FRMASMB17558.  
Right: FRMASMB17559  
Not available with 11"-13" extensions.
- 60°-90° Swing-Away Legrest, Short ..... \$142.50 ea.  
Pair part: RIG136481. Left: RIG132607.  
Right: RIG132609.  
For use only with 11"-13" lower extensions.

60°-90° Swing-Away Legrest Compatibility Matrix	Contoured Seat			Synergy Seat		
	Low	Medium	High	Low	Medium	High
Extruded, 13" - 16.5"	!	!	!	!	✓	✓
Extruded, 15.5" - 19.5"	X	X	!	!	!	!
Angle Adjustable, 13" - 16.5"	!	!	!	!	✓	✓
Angle Adjustable, 15.5" - 19.5"	X	X	!	!	!	!
Multi-Axis, 15.5" - 19.5"	X	X	!	!	!	!

This rigging configuration causes no interference issues at this seat-to-floor height.

This rigging configuration may have less than 1.75" clearance when fully extended and may interfere with the ground or reduce obstacle climbing ability.

This rigging configuration does not work with this seat-to-floor height.

Legrests section continued on following page.



**20. LEGRESTS**

**20.4. 60°-90° Swing-Away Legrests continued**

**2. Lower Extension**

<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>	<b>11"-13" Length</b> Part: See XRef N/A with Multi-axis.	<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>	<b>11"-16.5" Length</b> Part: See XRef N/A with Angle Adjustable.	<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>	<b>13"-16.5" Length</b> Part: See XRef N/A with Multi-Axis.	<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>	<b>15.5"-19.5" Length</b> Part: See XRef
--	---	--	---	--	---	--	---

**3a. Extruded Footplate(s)**

<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>	<b>7.25"Wx6"D Extruded Footplate</b> ..... Standard Part pair: FRMASMB7847. Left: FRMASMB7861 Right: FRMASMB7862	<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>	<b>8.5"Wx6"D Extruded Footplate</b> ..... Standard Part pair: FRMASMB7848. Left: FRMASMB7863 Right: FRMASMB7864 May cause interference on seat widths less than 18."
--	--	--	---

**3b. Angle Adjustable Footplate(s) with Heel Loops**

<sup>1</sup>The 6"x11" footplates may cause interference with seat widths less than 18". <sup>2</sup>The 8"x11" footplates may cause interference with seat widths less than 20".

<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>	<b>4"Wx6"D</b> ..... \$87.50 ea. Pair part: FRMASMB10564. Left: RIG133537. Right: RIG133560	<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>	<b>6"Wx8"D<sup>1</sup></b> ..... \$87.50 ea. Pair part: FRMASMB10567. Left: RIG133540. Right: RIG133563
<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>	<b>5"Wx6"D</b> ..... \$87.50 ea. Pair part: FRMASMB10565. Left: RIG133538. Right: RIG133561	<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>	<b>6"Wx11"D<sup>1</sup></b> ..... \$87.50 ea. Pair part: FRMASMB10568. Left: RIG133541. Right: RIG133564
<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>	<b>5"Wx8"D</b> ..... \$87.50 ea. Pair part: FRMASMB10566. Left: RIG133539. Right: RIG133562	<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>	<b>8"Wx11"D<sup>2</sup></b> ..... \$87.50 ea. Pair part: FRMASMB10569. Left: RIG133542. Right: RIG133565

**3c. Multi-Axis Footplate(s) with Heel Loops**

<sup>1</sup>The 6"x11" footplates may cause interference with seat widths less than 18". <sup>2</sup>The 8"x11" footplates may cause interference with seat widths less than 20".

<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>	<b>4"Wx6"D</b> ..... \$260 ea. Pair part: FRMASMB9456. Left: RIG133531. Right: RIG133552	<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>	<b>6"Wx8"D<sup>1</sup></b> ..... \$260 ea. Pair part: FRMASMB9459. Left: RIG133534. Right: RIG133555
<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>	<b>5"Wx6"D</b> ..... \$260 ea. Pair part: FRMASMB9457. Left: RIG133532. Right: RIG133553	<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>	<b>6"Wx11"D<sup>1</sup></b> ..... \$260 ea. Pair part: FRMASMB9460. Left: RIG133535. Right: RIG133558
<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>	<b>5"Wx8"D</b> ..... \$260 ea. Pair part: FRMASMB9458. Left: RIG133533. Right: RIG133554	<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>	<b>8"Wx11"D<sup>2</sup></b> ..... \$260 ea. Pair part: FRMASMB9461. Left: RIG133536. Right: RIG133559

**4. 60°-90° Swing-Away Legrest Accessories**

<b>Heel Loops (pair)</b> ..... \$35 Part: FRMASMB7873. Standard on angle adjustable and multi-axis footplates.	<b>Wheel Bumper</b> ..... \$15 ea. Part: See XRef
	<b>Gel Neoprene Footrest Sleeve (pair)</b> ..... \$130 Part: INDPART2243.

**20.5. Elevating Legrests**

Select a legrest upper, a legrest extension, and type of footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes. Select legrest length in following section.

**1. Elevating Legrests Upper**

<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>	<b>Elevating Legrest Assembly</b> ..... \$550 ea. Pair part: FRMASMB6483. Left: RIG138012 Right: RIG138013 Select lower extension and footplates
--	---

**2. Lower Extension**

<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>	<b>11"-14" Length</b> Pair: FRMASMB4755 Left: RIG138026 Right: RIG138030 For reinforced only	<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>	<b>14"-17" Length</b> Pair: FRMASMB4756 Left: RIG138033 Right: RIG138035 For reinforced only
<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>	<b>12"-14.5" Length</b> Part: See XRef For angle adjustable only	<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>	<b>14.5"-18" Length</b> Part: See XRef For angle adjustable only

**3a. Reinforced Cast Footplate(s)**

<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>	<b>7.25"Wx6"D Reinforced Footplate</b> ..... Standard Pair part: FRMASMB4748. Left: RIG138055. Right: RIG138056	<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>	<b>8.5"Wx6"D Reinforced Footplate</b> ..... Standard Pair part: FRMASMB4747. Left: RIG138060. Right: RIG138062 May cause interference on seat widths less than 18."
--	--	--	---

Elevating Legrest Compatibility Matrix	Contoured Seat			Synergy Seat		
	Low	Medium	High	Low	Medium	High
Standard, 11"-14"	!	✓	✓	✓	✓	✓
Standard, 14"-17"	X	!	!	!	!	!
Angle Adjustable, 12"-14.5"	!	✓	✓	✓	✓	✓
Angle Adjustable, 14.5"-18"	X	!	!	!	!	!

This rigging configuration causes no interference issues at this seat-to-floor height.  
 This rigging configuration may have less than 1.75" clearance when fully extended and may interfere with the ground or reduce obstacle climbing ability.  
 This rigging configuration does not work with this seat-to-floor height.

Legrests section continued on following page.

**20. LEGRESTS CONTINUED**

**20.5. Elevating Legrests continued**

**3b. Angle Adjustable Footplate(s)**

<sup>1</sup>The 6"x11" footplates may cause interference with seat widths less than 18". <sup>2</sup>The 8"x11" footplates may cause interference with seat widths less than 20". 4"x6" and 5"x6" footplates are only available in pairs.

<table border="0"> <tr> <td><input type="checkbox"/> P</td> <td><input type="checkbox"/> L</td> <td><input type="checkbox"/> R</td> <td>4"Wx6"D ..... \$87.50 ea. <i>Pair parts: FRMASMB8812 x2.</i></td> </tr> <tr> <td><input type="checkbox"/> P</td> <td><input type="checkbox"/> L</td> <td><input type="checkbox"/> R</td> <td>5"Wx6"D..... \$87.50 ea. <i>Pair parts: FRMASMB8811 x2.</i></td> </tr> <tr> <td><input type="checkbox"/> P</td> <td><input type="checkbox"/> L</td> <td><input type="checkbox"/> R</td> <td>5"Wx8"D Footplates with Heel Loops ..... \$87.50 ea. <i>Pair part: FRMASMB4735. Left: RIG137791. Right: RIG137792</i></td> </tr> </table>	<input type="checkbox"/> P	<input type="checkbox"/> L	<input type="checkbox"/> R	4"Wx6"D ..... \$87.50 ea. <i>Pair parts: FRMASMB8812 x2.</i>	<input type="checkbox"/> P	<input type="checkbox"/> L	<input type="checkbox"/> R	5"Wx6"D..... \$87.50 ea. <i>Pair parts: FRMASMB8811 x2.</i>	<input type="checkbox"/> P	<input type="checkbox"/> L	<input type="checkbox"/> R	5"Wx8"D Footplates with Heel Loops ..... \$87.50 ea. <i>Pair part: FRMASMB4735. Left: RIG137791. Right: RIG137792</i>	<table border="0"> <tr> <td><input type="checkbox"/> P</td> <td><input type="checkbox"/> L</td> <td><input type="checkbox"/> R</td> <td>6"Wx8"D<sup>1</sup> Footplates with Heel Loops ..... \$87.50 ea. <i>Pair part: FRMASMB4738. Left: RIG137793. Right: RIG137795</i></td> </tr> <tr> <td><input type="checkbox"/> P</td> <td><input type="checkbox"/> L</td> <td><input type="checkbox"/> R</td> <td>6"Wx11"D<sup>1</sup> Footplates with Heel Loops ..... \$87.50 ea. <i>Pair part: FRMASMB4741. Left: RIG137796. Right: RIG137797</i></td> </tr> <tr> <td><input type="checkbox"/> P</td> <td><input type="checkbox"/> L</td> <td><input type="checkbox"/> R</td> <td>8"Wx11"D<sup>2</sup> Footplates with Heel Loops ..... \$87.50 ea. <i>Pair part: FRMASMB4744. Left: RIG137798. Right: RIG137799</i></td> </tr> </table>	<input type="checkbox"/> P	<input type="checkbox"/> L	<input type="checkbox"/> R	6"Wx8"D <sup>1</sup> Footplates with Heel Loops ..... \$87.50 ea. <i>Pair part: FRMASMB4738. Left: RIG137793. Right: RIG137795</i>	<input type="checkbox"/> P	<input type="checkbox"/> L	<input type="checkbox"/> R	6"Wx11"D <sup>1</sup> Footplates with Heel Loops ..... \$87.50 ea. <i>Pair part: FRMASMB4741. Left: RIG137796. Right: RIG137797</i>	<input type="checkbox"/> P	<input type="checkbox"/> L	<input type="checkbox"/> R	8"Wx11"D <sup>2</sup> Footplates with Heel Loops ..... \$87.50 ea. <i>Pair part: FRMASMB4744. Left: RIG137798. Right: RIG137799</i>
<input type="checkbox"/> P	<input type="checkbox"/> L	<input type="checkbox"/> R	4"Wx6"D ..... \$87.50 ea. <i>Pair parts: FRMASMB8812 x2.</i>																						
<input type="checkbox"/> P	<input type="checkbox"/> L	<input type="checkbox"/> R	5"Wx6"D..... \$87.50 ea. <i>Pair parts: FRMASMB8811 x2.</i>																						
<input type="checkbox"/> P	<input type="checkbox"/> L	<input type="checkbox"/> R	5"Wx8"D Footplates with Heel Loops ..... \$87.50 ea. <i>Pair part: FRMASMB4735. Left: RIG137791. Right: RIG137792</i>																						
<input type="checkbox"/> P	<input type="checkbox"/> L	<input type="checkbox"/> R	6"Wx8"D <sup>1</sup> Footplates with Heel Loops ..... \$87.50 ea. <i>Pair part: FRMASMB4738. Left: RIG137793. Right: RIG137795</i>																						
<input type="checkbox"/> P	<input type="checkbox"/> L	<input type="checkbox"/> R	6"Wx11"D <sup>1</sup> Footplates with Heel Loops ..... \$87.50 ea. <i>Pair part: FRMASMB4741. Left: RIG137796. Right: RIG137797</i>																						
<input type="checkbox"/> P	<input type="checkbox"/> L	<input type="checkbox"/> R	8"Wx11"D <sup>2</sup> Footplates with Heel Loops ..... \$87.50 ea. <i>Pair part: FRMASMB4744. Left: RIG137798. Right: RIG137799</i>																						

**4. Elevating Legrest Accessories**

Heel Loops (pair)..... \$35 <i>Part: FRMASMB7873</i>	Falcon Gel Padded Knee Adductor Buttons (KBE) (pair) .....\$125 <i>Part: INDPART2242</i>
Wheel Bumper..... \$15 ea. <i>Part: See XRef</i> Available on angle adjustable footplates.	<sup>^</sup> Please verify the appropriate billing code to be used for this item with all third-party payors.

**21. SET LEGREST LENGTH**

**REQUIRED with matching left and right legrest configuration.** All efforts will be made to accommodate legrest length selection, but minor changes may need to be made to meet shipping guidelines. Some adjustments may be needed upon receipt of the unit. Reference the section's compatibility matrix for legrest interference. Non-matching or individual legrest configurations lengths are not set by the manufacturer.

7" L (#3470)	8.5" L (#3473)	10" L (#3476)	11.5" L (#3479)	13" L (#3482)	14.5" L (#3485)	16" L (#3488)	17.5" L (#3491)	19" L (#3494)	20.5" L (#3497)	22" L (#3669)
7.5" L (#3471)	9" L (#3474)	10.5" L (#3477)	12" L (#3480)	13.5" L (#3483)	15" L (#3486)	16.5" L (#3489)	18" L (#3492)	19.5" L (#3495)	21" L (#3498)	22.5" L (#3670)
8" L (#3472)	9.5" L (#3475)	11" L (#3478)	12.5" L (#3481)	14" L (#3484)	15.5" L (#3487)	17" L (#3490)	18.5" L (#3493)	20" L (#3496)	21.5" L (#3600)	


**22. SEAT CUSHION**

**Optional. Specific cushion and component sizes available, part numbers, and special order cushions can be found on the respective complete [Stealth Cushion Order Forms](#).** This section allows a quick selection of a single cushion. The XXYY portion of the part number directly relates to the dimensions of the cushion, where XX represents the width and YY represents the depth.

**22.1. Seat Cushion Model Selection**


**Stealth Simplicity G Cushion**

General use cushion with Stealth Tek 2, 2-way stretch, breathable cover.

Simplicity G ..... \$310-\$370 


**Stealth Spectrum Foam SPP Cushion**

Skin protection and positioning with Coolcore, 4-way stretch, breathable cover. Includes urethane liner.

Spectrum Foam SPP ..... \$580-\$630 


**Stealth Solution SPP Cushion**

Skin protection and positioning with Coolcore, 4-way stretch, breathable cover. Includes urethane liner.

Solution SPP ..... \$580-\$645 


**Stealth Glacial SP Cushion**

Skin protection with Coolcore, 4-way stretch, breathable cover.

Glacial SP ..... \$265-\$530 


**Stealth Zen SP Cushion**

Skin protection and positioning with Coolcore, 4-way stretch, breathable cover.

Zen SP ..... \$455-\$765 


**Stealth TRU-Comfort 2 SPP Cushion**

Skin protection and positioning with Coolcore, 4-way stretch, breathable cover.

TRU-Comfort 2 SPP ..... \$515-\$575 


**Stealth Spectrum Gel SPP Cushion**

Skin protection and positioning with Coolcore, 4-way stretch, breathable cover. Includes urethane liner.

Spectrum Gel SPP ..... \$580-\$650 


**Stealth Essence SPP Cushion**

Skin protection and positioning with Coolcore, 4-way stretch, breathable cover.

Essence SPP ..... \$545-\$705  
Essence Plus SPP ..... \$800 

**Stealth Premiere P Cushion**

Positioning with Coolcore, 4-way stretch, breathable cover.

Premiere P ..... \$585 

**22. SEAT CUSHION CONTINUED**

**22.2. Seat Cushion Dimensions**

**REQUIRED with cushion model selection.** Select to match or pick non-matching dimensions. See respective cushion order form for available sizes.

Match Seat Width and Depth

Cushion Width						Cushion Depth							
14" W	16" W	17" W	18" W	19" W	20" W	22" W	14" D	16" D	17" D	18" D	19" D	20" D	22" D

**22.3. Seat Cushion Substitute Fluid Proof Liners**

Fluid Proof Urethane Liner ..... \$45-\$50 Part: URL1709070-XXYY Not available with Simplicity, Essence, or TRU-Comfort 2.	Fluid Proof Polyurethane Liner ..... \$70 Part: PUL1709030-XXYY Not available with Simplicity.
--	--

**22.4. Seat Cushion Substitute Cover**

The standard cover may be substituted with a cover by checking the option below. The "\*"prefix" of the part number will be the part prefix of the cushion model selected above. For example, SPCGSP-NBLK1816 will be an additional coolcore cover for a Spectrum Gel cushion.

Extra Stealth Tek 2 Cover ..... \$90 Part: *prefix-TBLKXXYY Not available with Simplicity, TRU-Comfort 2, Essence, or Premiere.	Extra Mesh Cover ..... \$90 Part: *prefix-MBLKXXYY Not available with Simplicity or TRU-Comfort 2.
---	--

**23. PROVIDER INSTALLED ACCESSORIES**

**Optional.** Multiple accessories can be mounted along the various tracks but some interference may occur.

Contoured Seat Cup Holder ..... \$20 Part: ACCASMB2287. Part with 12" lift: ACCASMB2288	Flag Holder ..... \$125 Contoured seat part: ACCASMB1930 Synergy/Flex part: ACCASMB2415.
Synergy/Flex Seat Cup Holder ..... \$20 Part: ACCASMB2634 Not available on TRU-Balance Flex seating without two-post armrests.	Saddle Bag ..... \$15 Part: DWR1020D010
Oxygen Holder for Contoured Seat ..... \$110 Part: ACCASMB2292.	Quantum Backpack ..... \$65 Part: ACC190585
Oxygen Holder for Flex Seat ..... \$110 Part: ACC125006	Accessory Bar ..... \$90 12"W-16"W part: MEC145955. 16"W-20"W part: MEC143975. 20"W-22"W part: MEC146018
Walker Holder ..... \$150 Contoured seat part: ACCASMB1912. Synergy/Flex part: FRMASMB7031	Large, Black, Weather Cover ..... \$240 Part: ACCCOVR1021
Cane/Crutch Holder ..... \$195 Contoured seat part: ACCASMB2294. HCPCS: E2207 Synergy part: FRMASMB4822. Flex part: ACC152515.	
Rear Basket ..... \$99 Contoured seat part: ACCASMB1914 Synergy part: FRMASMB4821	

**24. BATTERIES**

**REQUIRED.** The powerchair requires 2 batteries and will be installed if ordered.

Group 24, Gel Cell (MK).....\$510 ea. Part: BATGEL1009.	No batteries. Provider will assume responsibility of batteries. Part: Configuration selection
--	--