



Pride Mobility Products Co.

5096 South Service Rd, Beamsville, ON, L0R 1B3

Phone: 888-570-1113 | Fax: 866-514-1303 | Email: productorders@pridemobility.com

Q6 Edge® 3 Stretto Order Form with TRU-Balance® Flex Seating

Account Number: _____ Date: _____

Provider Name: _____

Contact: _____

Phone: _____ Fax: _____

Email: _____

PO Number: _____

Marked for: _____

Ship to Address: _____

City: _____ Prov: _____ PC: _____



INTRODUCTION

This form is interactive when viewed with [Adobe Acrobat Reader](#). Complete the form by placing checks in the desired boxes and provide information in the interactive fields. Buttons shown at the bottom of the form may be utilized to print or submit the order form through a desktop email application. To email via a web-based application, please 'Save As' and attach the PDF to your email.

This order form contains a large variety of options to fulfill various patient needs. Descriptions and section notes, such as optional or required, are included to help you complete your order. Other part numbers can be found on the quote acknowledgment after this order form is submitted. Please contact Quantum® Sales at 888-570-1113 if assistance is needed. Send the completed order form by fax (866-514-1303) or email (productorders@pridemobility.com). Incomplete forms may delay the quote or order. Customer service will contact you if the order is incomplete or if there are compatibility issues. If special order requests are needed, be sure that the Patient Information section is completed

All prices are MSRP. Prices, specifications, part numbers, and availability are subject to change without notice. Prices and part numbers as shown are only available when configured as a complete power chair. Please contact technical service for accurate parts ordering.


PATIENT INFORMATION

Weight: _____ Height: _____

A. Shoulder Width: _____ D. Max Sitting Height: _____ G. Thigh Depth: _____

B. Chest Width: _____ E. Shoulder Height: _____ H. Lower Leg Length: _____

C. Hip Width: _____ F. Axilla Height: _____ I. Elbow Height: _____



1. BASE MODEL

REQUIRED. Select a group 3 power base. 300 lb. weight capacity with standard 6mph programmable motors and off board charger. XLR USB charger is included. Dimensions of base without legrests: 20.5"Wx35.5"L

Quantum Q6 Edge 3 Stretto\$8,450
Part: Q6 EDGE 3 STRETTO

DRIVE MOTORS

Optional. Requires Q-Logic 3 electronics. The standard 6 mph motor can be upgraded to Accu-Trac motors that delivers enhanced traction and tracking, especially with switched drive controls.

Accu-Trac Motor Technology\$995
Part: MOT147881.

Standard Motor Package..... Standard

2. BASE COLOR


















REQUIRED. Select one shroud color. Colors shown here may differ from actual product depending on monitor/printer calibration used.

2.1. Match Colors

Match All Color Options to Selection Below

Part: N/A. Any subsequent color selections, if available, will be matched to color chosen below when the order is placed. Please note that some color selections, such as joystick armor, may have a MSRP.

2.2a. Standard Shrouds

- | | | | |
|---|--|--|--|
|  Red Corvette
Part: DGN1900199-PRD1 |  White Satin
Part: DGN1900199-PWH2 |  Midnight Blue
Part: DGN1900199-PBU10 |  Military Green
Part: DGN1900199-PGN8 |
|  Cherry Bomb
Part: DGN1900199-PRD4 |  Deep Purple
Part: DGN1900199-PPR5 |  Electric Blue
Part: DGN1900199-PBU9 |  Sugar Plum
Part: DGN1900199-PPR7 |
|  Back in Black
Part: DGN1900199-PBK1 |  Raspberry Beret
Part: DGN1900199-PPK4 |  Blue Jean
Part: DGN1900199-PBU11 | To Be Determined
Follow-up required. |
|  Night Train
Part: DGN1900199-PBK3 |  Pretty in Pink
Part: DGN1801101-PPK3 |  Lemon Crush
Part: DGN1900199-PYL4 | |
|  Grey Street
Part: DGN1900199-PSL3 |  Evolution Orange
Part: DGN1900199-POR7 |  Green Machine
Part: DGN1900199-PGN11 | |

3. FENDER LIGHTS

Optional. Fender lights operate by pressing included single switch button. Fender lighting is not available with manual recline.

- Switch Fender Lights..... No Charge
Part: DGN1810041

4. FULL LIGHT PACKAGE

Optional. Requires NE+ 8 Key or Q-Logic 3 Joystick for Full Light Package.

- Full Light Package..... \$770
Part: ELE1905403

5. QUANTUM BACKUP CAMERA

Optional. The backup camera may interfere with armrest mounted accessories such as media device holder.

- Quantum Backup Camera \$475
Part: MEC1910365

1. Optional Mounting Bracket

- Backup Camera Headrest Bracket..... \$30
Part: MEC1910373
Works with .5" x .5" vertical supports only.



5. BATTERIES

REQUIRED. The powerchair requires 2 batteries and will be installed if ordered.

- | | |
|---|--|
| NF-22, Gel, 55AH, Interceptor Battery..... \$300 ea.
Part: BAT1704613.
Uses 8 AMP charger | Gel, 40AH Interceptor Battery..... \$200 ea.
Part: BAT1704612.
Uses 5 AMP charger. |
|---|--|
- No batteries. Provider will assume responsibility of batteries.
Part: Configuration selection

6. DRIVE WHEELS

REQUIRED. Select a type of drive wheel. Accent colors are available for 14" tires.

6.1. Drive Wheels

- | | |
|--|--|
| 2"x12.5" Pneumatic Tires No Charge
Part: WHL1810666 | 3"x14" Pneumatic Tires..... No Charge
Part: WHL1713925 |
| 2"x12.5" Flat Free Tires No Charge
Part: WHL1810665 | 3"x14" Flat Free Tires No Charge
Part: WHL1704058 |

6.2. Drive Wheel Accent Color for 14" Tires

Not available with 12.5" tires.

- | | | |
|---|--|--|
| Back in Black No Charge
Part: DGN164245-PBK1 x2 | Deep Purple No Charge
Part: DGN164245-PPR5 x2 | Blue Jean No Charge
Part: DGN164245-PBU11 x2 |
| Red Corvette No Charge
Part: DGN164245-PRD1 x2 | Raspberry Beret No Charge
Part: DGN164245-PPK4 x2 | Lemon Crush No Charge
Part: DGN164245-PYL4 x2 |
| Cherry Bomb No Charge
Part: DGN164245-PRD4 x2 | Pretty in Pink No Charge
Part: DGN164245-PPK3 x2 | Green Machine No Charge
Part: DGN164245-PGN11 x2 |
| Night Train No Charge
Part: DGN164245-PBK3 x2 | Evolution Orange No Charge
Part: DGN164245-POR7 x2 | Military Green No Charge
Part: DGN164245-PGN8 x2 |
| Grey Street No Charge
Part: DGN164245-PSL3 x2 | Midnight Blue No Charge
Part: DGN164245-PBU10 x2 | Sugar Plum No Charge
Part: DGN164245-PPR7 x2 |
| White Satin No Charge
Part: DGN164245-PWH2 x2 | Electric Blue No Charge
Part: DGN164245-PBU9 x2 | To Be Determined
Follow-up required. |

7. ELECTRONICS

REQUIRED. Select one type of drive control from the following sections. The No Joystick option is used to omit the joystick and utilize another control, such as a specialty control at the end of this form or found on the [Alternative Controls and Electronics Order Form](#).

7.1a. Non-Expandable Controls

Select a joystick below or omit selection if another controller type is desired. Non-expandable does not allow supplementary drive controllers.

NE 4-Key Joystick Standard
 Part: CTLDC1554
 Available with Tilt Only through Toggle.

NE+ 8-Key Joystick (2 Actuators) No Charge
 Part: CTLDC1561
 Control full light package.

NE+ 6-Key Joystick No Charge
 Part: CTLDC1560
 Control of one actuator through Joystick on single power base model (3SP-SS) only.

No Joystick, Setup for NE+ Joystick No Charge
 Part: Configuration selection



7.1b. Expandable Controls - Quantum Intelligent Systems™

Select expandable controls or omit selection if another controller type is desired. Expandable controls allow multiple or specialized controllers.

2. Q-Logic 3 Joystick



Q-Logic 3 EX Joystick \$900
 Part: CTL167724

Q-Logic 3 EX Joystick for Full Lights \$900
 Part: CTL167725

Q-Logic 3e Joystick \$620
 Part: CTL1910600
 Can operate full light package. Not available with alternative or specialty controls.

No Joystick, Setup for Q-Logic 3
 Part: XRF000010

7.2. Q-Logic 3 Joystick Armor

REQUIRED with Q-Logic 3 Joystick selections above. Reference colors shown on page 2. *Night train joystick armor will be the black, soft touch.

Black, Soft Touch No Charge
 Part: DGN124307
 Also matches Night Train shrouds.

Non-matching Armor No Charge
 Part: PTOINDV3422
 Select color below.

Match Armor to Base Color Standard
 Part: PTOINDV3421

Red Corvette
 Part: DGN124336-PRD1

Grey Street
 Part: DGN124336-PSL3

Pretty in Pink
 Part: DGN124336-PPK3

Blue Jean
 Part: DGN124336-PBU11

Sugar Plum
 Part: DGN124336-PPR7

Cherry Bomb
 Part: DGN124336-PRD4

White Satin
 Part: DGN124336-PWH2

Evolution Orange
 Part: DGN124336-POR7

Lemon Crush
 Part: DGN124336-PYL4

To Be Determined
 Follow-up required.

Back in Black
 Part: DGN124336-PBK1

Deep Purple
 Part: DGN124336-PPR5

Midnight Blue
 Part: DGN124336-PBU10

Green Machine
 Part: DGN124336-PGN11

Night Train*
 See Black, Soft Touch

Raspberry Beret
 Part: DGN124336-PPK4

Electric Blue
 Part: DGN124336-PBU9

Military Green
 Part: DGN124336-PGN8

8. JOYSTICK HANDLES

Optional. The NE/NE Plus Joysticks are standard with a carrot handle (See A in image). The optional handle will be mounted to the joystick and the standard handle will be shipped with the unit.

Mushroom Handle (B) \$120
 Part: ACCASMB1860.

Extended I-Handle (G) \$140
 Part: ACCASMB1857.

Chin Cup (C) \$85
 Part: ACCASMB186190

Bodypoint 3" Goal Post Handle \$180
 Part: PLSKNCP1057.

T-Bar (D) \$160
 Part: ACCASMB1858.

Bodypoint 4" Goal Post Handle \$180
 Part: PLSKNCP1058.

Medium Ball Handle (E) \$70
 Part: ACCASMB1859.

Soft Ball Handle (F) \$120
 Part: ACCASMB1856.




9. JOYSTICK MOUNTING BRACKETS

Optional. Select a mount to be used with a joystick selected in previous section.

Omit Joystick Mounting Bracket

9.1a. Fixed, In-line Joystick Bracket


- Fixed, In-line, Bracket Standard
Q-Logic Part: FRMASMB10617. NE/NE+ Part: FRMASMB12744
- Left Side Mount
Part: PTOINDV1225
- Right Side Mount
Part: PTOINDV1226



9.1b. Swing-Away (Retractable) Joystick Bracket

Allows the joystick to retract to the side of the armrest while maintaining the same joystick orientation. "Inline bracket" is used with TRU-Balance 3 arms when the joystick is inline with armpad. "Inline bracket" will be offset 1" when used with tube style arms. Select "Offset Bracket" to have the joystick inline with armpad on tube style armrests.


- Swing-away, Q-Logic 3 \$265
Part: ACC1908434.
- Swing-away, NE/NE+ \$265
Part: ACC1908437.
- Left Side Mount
Part: PTOINDV1225
- Right Side Mount
Part: PTOINDV1226



9.1c. Specialty Joystick Mounts

Select a bracket and side to mount it.

- Retract4 Mount \$450
NE/NE+ Part: MEC138876. ADP: WEPM
Q-Logic Part: MEC138860. ADP: WEPM
 Full 180° of movement that translates into multi-position use, straight, right offset, left offset, or fold away. Joystick orientation varies as bracket is adjusted.
- Stealth Swing-Away/Flip Down Bracket \$300
NE/NE+ Part: ST-TWBM480QL-Q. ADP: WEPM
Q-Logic Part: ST-TWBM480QL-Q. ADP: WEPM
SAJ Part: ST-TWBM480CJ-Q. ADP: WEPM
 Allows the joystick to pivot to side of armrest, or flip-down if bracket is adjusted.
- Multi-axis Mount \$100
NE/NE+ Part: ACC142528. ADP: WEPM
Q-Logic Part: ACC131216. ADP: WEPM
SAJ Part: ACC142740. ADP: WEPM
 Allows height adjustment of joystick.



Specialty Mount Side

- Left Side Mount
Part: PTOINDV1225
- Right Side Mount
Part: PTOINDV1226

10. TRU-BALANCE FLEX SEATING SYSTEM

Seat width and depth can be adjusted in 1" increments from 15"W-20"W and 16"D-20"D. Maximum seat size 20"W x 20"D.

10.1. TRU-Balance Flex Frame Selection

REQUIRED.

- TRU-Balance Flex Seating System \$650
Part: PTOINDV3853

10.2. TRU-Balance Flex Seat Dimension - Width

REQUIRED. Maximum seat width is 20". Part number is SETWDTH plus number following # below.

10" W (#1000)	11" W (#1001)	12" W (#1002)	13" W (#1003)	14" W (#1004)	
15" W (#1005)	16" W (#1006)	17" W (#1007)	18" W (#1008)	19" W (#1009)	20" W (#1010)

10.3. TRU-Balance Flex Seat Dimension - Depth

REQUIRED. Maximum seat depth is 20". Part number is SETDPTH plus number following # below.

10" D (#1000)	11" D (#1001)	12" D (#1002)	13" D (#1003)	14" D (#1004)	
15" D (#1005)	16" D (#1006)	17" D (#1007)	18" D (#1008)	19" D (#1009)	20" D (#1010)

10.4. TRU-Balance Flex Seat to Floor

REQUIRED.

- 16.5" at front and rear
Part: PTOINDV3459
- 17.5" at front with 2° posterior angle
Part: PTOINDV3460
- 18.25" at front with 4° posterior angle
Part: PTOINDV3805

10.5. TRU-Balance Flex Back Cane Selection

REQUIRED.

- Angle Adjustable Back Cane \$195
Part: MEC140266
- Flip Forward Back Cane \$195
Part: MEC140272
 Not available with oxygen holder, or cane mounted and cantilever armrests.
- 2" Mesomorph (1" per side) Back Cane \$205
Part: MEC1804736
- 4" Mesomorph (2" per side) Back Cane \$205
Part: MEC1804739
- 2" Endomorph (1" per side) Back Cane \$205
Part: MEC141455
- 4" Endomorph (2" per side) Back Cane \$205
Part: MEC141456
- 24"H Flip Forward Flex Cane Standard

TRU-Balance Flex seating section continued on following pages.

10. TRU-BALANCE FLEX SEATING SYSTEM CONTINUED

10.6. TRU-Balance Flex Back Cane Height Selection

REQUIRED. 16" and 17" back canes are not available with a Synergy Back and cane mounted, flip up armrests.

14" High Canes Part: FRMASMB15832	17" High Canes Part: FRMASMB15835	20" High Canes Part: FRMASMB15838	23" High Canes Part: FRMASMB15841
15" High Canes Part: FRMASMB15833	18" High Canes Part: FRMASMB15836	21" High Canes Part: FRMASMB15839	24" High Canes Part: FRMASMB15842
16" High Canes Part: FRMASMB15834	19" High Canes Part: FRMASMB15837	22" High Canes Part: FRMASMB15840	25" High Canes Part: FRMASMB15843

10.7. TRU-Balance Flex Back Options

REQUIRED. See the back availability chart for stock order back size ranges. Other size backs may be available with special order and extended lead time. Select a back type and overall back height for configured part number.

1. Back Selection

Synergy Back Cushion.....	\$940
Part: See XRef	
Sport Back Cushion.....	No Charge
Part: See XRef	
TRU-Comfort 2, 20" Width or less	\$835
Part: See XRef.	
Omit Back.....	No Charge
Part: PTOINDV1287	

		Back Width											
		15"	16"	17"	18"	19"	20"	21"	22"				
Overall Height	16"	-	P	T	P	T	P	T	P	T	P	T	
	17"	S	SPT	SPT	SPT	SPT	SPT	SPT	SPT	SPT	SPT	SPT	
	18"	S	SPT	SPT	SPT	SPT	SPT	SPT	SPT	SPT	SPT	SPT	
	19"	S	SPT	P	T	SPT	SPT	SPT	SPT	SPT	SPT	SPT	
	20"	S	SPT	P	T	SPT	SPT	SPT	SPT	SPT	SPT	SPT	
	21"	S	SPT	SPT	SPT	SPT	SPT	SPT	SPT	SPT	SPT	SPT	
	22"	S	SPT	SPT	SPT	SPT	SPT	SPT	SPT	SPT	SPT	SPT	
	23"	S	SPT	SPT	SPT	SPT	SPT	SPT	SPT	SPT	SPT	SPT	
	24"	S	S	P	S	P	SPT	S	P	SPT	S	P	S
	25"	S	S	P	S	P	SPT	S	P	SPT	S	P	S

S Synergy Back P Sport Back T TRU-Comfort Back

TRU-Balance Flex Seat Back Availability Chart

2. TRU-Balance Flex Overall Back Height

Refer to the Synergy and TRU-Comfort back availability chart for stock order back size ranges. Parts: See XRef, configuration selection with back selection.

14" Height Part: SETHGHT1027	17" Height Part: SETHGHT1030	20" Height Part: SETHGHT1033	23" Height Part: SETHGHT1036	26" Height Part: SETHGHT1039
15" Height Part: SETHGHT1028	18" Height Part: SETHGHT1031	21" Height Part: SETHGHT1034	24" Height Part: SETHGHT1037	27" Height Part: SETHGHT1040
16" Height Part: SETHGHT1029	19" Height Part: SETHGHT1032	22" Height Part: SETHGHT1035	25" Height Part: SETHGHT1038	28" Height Part: SETHGHT1041

10.8. TRU-Balance Flex Armrests

REQUIRED. Left (L) and right (R) boxes may be checked to configure different size arms within one section arm type.

1a. TRU-Balance Flex Seat, Drop-In Armrests

Includes full length, flat armpads.

Full Length, Drop-In Armrests.....	No Charge	Desk Length, Drop-In Armrests	No Charge
Left part: ARM1803253. Right part: ARM156742		Left part: ARM1803254. Right part: ARM156743	

1b. Flip Back, Cane Mounted Armrests

May interfere with swing-away laterals.

Full Length, Flip Back Armrests (pair).....	No Charge	Desk Length, Flip Back Armrests (pair)	No Charge
Left part: FRMASMB7038. Right part: FRMASMB7039		Left part: FRMASMB7040. Right part: FRMASMB7041	

1c. Two-Post, Flip Back, Height Adjustable, Removable Armrests

Not available with a depth less than 16". Required with flip forward back and laptrays.

Full Length, Two-Post Armrests.....	\$145 ea.	Desk Length, Two-Post Armrests	\$145 ea.
Part: See XRef.		Part: See XRef.	

1d. Single Post, Height Adjustable, Removable Armrests

Height adjustment from 11"-17" in half-inch increments.

Full Length, Single Post Armrests	\$145 ea.	Desk Length, Single Post Armrests	\$145 ea.
Left part: FRMASMB8594. Right part: FRMASMB8595.		Left part: FRMASMB8592. Right part: FRMASMB8593.	

10.9. TRU-Balance Flex Armpads

REQUIRED with armrest selection except Drop-In Armrests. Waterfall armpads may cause interference with the seat back. Armpads are configured as pairs.

1a. Straight Armpads

Straight, Full Length Left and Right	No Charge	Straight, Desk Length Left, Full Length Right	No Charge
Part: ARM134275		Part: ARM134275	
Straight, Full Length Left, Desk Length Right	No Charge	Straight, Desk Length Left and Right.....	No Charge
Part: ARM134281		Part: ARM134254	

TRU-Balance Flex seating section continued on following pages.

10. TRU-BALANCE FLEX SEATING SYSTEM CONTINUED

1b. Waterfall Armpads

Waterfall, Full Length Left and Right..... No Charge Part: ARM134275	Waterfall, Desk Length Left, Full Length Right No Charge Part: ARM134275
Waterfall, Full Length Left, Desk Length Right No Charge Part: ARM134281	Waterfall, Desk Length Left and Right No Charge Part: ARM134254

1c. Falcon Gel Armpads

Falcon Gel, Full Length Left and Right \$300 Part: ARM136438	Falcon Gel, Desk Length Left, Full Length Right..... \$250 Part: ARM136439
Falcon Gel, Full Length Left, Desk Length Right..... \$250 Part: ARM136442	Falcon Gel, Desk Length Left and Right \$200 Part: ARM136437

11. HEADRESTS

Optional. Select a headrest pad and bracket.

11.1. Headrest Pads

Stealth 8" Comfort Plus (CP450)..... \$190 Part: POS148484.	Pediatric Comfort Plus (6"Wx4.5"H) (CPN350) \$190 Part: SETHEAD1069.
Stealth 10" Comfort Plus (CP250) \$190 Part: SETHEAD1057.	Stealth Small Contoured (9.75"Wx3"H) (STL-OTSM250) \$80 Part: SETHEAD1052.
Stealth 14" Comfort Plus (CP150) \$190 Part: INDPART2327.	Stealth Large Contoured (11"Wx6"H) (STL-OTLG150) \$100 Part: SETHEAD1053.

11.2. Headrest Brackets

Stealth's Worlds Best Removable..... \$220.50 Headrest Bracket (TWB485-1) Part: FRMASMB15999. HCPCS: E1028	Pediatric Offset (HMO475P-17) \$160.50 Part: HDWASMB1166. HCPCS: E1028
--	---

12. LAP BELTS

REQUIRED. Select one lap belt below.

- 38" Lap Belt Part: ACCBELT1006
 50" Lap Belt Part: ACCBELT1000
 60" Lap Belt Part: ACCBELT1002
 70" Lap Belt Part: ACCBELT1003
 80" Lap Belt Part: ACCBELT1004
 90" Lap Belt Part: ACCBELT1001

13. LEGRESTS

Optional. All efforts will be made to accommodate legrest length selection, but minor changes may need to be made to meet shipping guidelines. Some adjustments may be needed upon receipt of the unit. Reference the section's compatibility matrix for legrest interference.

13.1. Center Mount Foot Platform

If power articulating foot platform is selected, be sure power legrest control has been selected in earlier section. Select legrest length in following section.

1. Foot Platform Assembly

Manually Adjustable Center Mount Foot Platform.....\$265 8"-12" extension part: RIG1802138. Other extensions: RIG1710345

2. Lower Extension and Angle Adjustable Footplates

10"-16" Extension, 5"x6" Footplates Part: RIG150206	10"-16" Extension, 5"x11" Footplates Part: RIG150212
10"-16" Extension, 5"x8" Footplates Part: RIG150203	

3. Center Mount Calf Support

Small Calf Pads (pair) 6"Wx6"H.....\$210 Part: RIGASMB7120018
--

4. Center Mount Modular Footplate Extension Components

Optional. Assembly required. May require seating configuration adjustments to avoid front caster interference.

Width Extension for 8"D Footplates \$55 Part: RIG152517	4" Outside Edge (pair) \$55 Part: RIG152510 Only available with 6"D footplates and 8"D footplate extension.
Width Extension for 11"D Footplates \$55 Part: RIG152518	6" Outside Edge (pair) \$55 Part: RIG152516 Available with 8"D and 11"D footplates and 11"D footplate extension.
2" Bell-Shaped Length Extension \$55 Part: RIG152519 Not compatible with width extensions. Can be used with 6" outside edge with 11" footplates.	

Legrests section continued on following pages.

13. LEGRESTS CONTINUED

5. Center Mount Additional Accessories

- Reversed Dartex Calf Pad Cover (pair) \$150
Small part: RIG133717. Large part: RIG133719
- Reversed Dartex Panel Cover (pair) \$150
Small part: RIG134673. Large part: RIG134744
- Wheel Bumpers (pair)..... \$15 ea.
Part: WHLASMB1578 x2
- Footplate Connector \$55
Part: RIGASMB7120023

6. Stealth Gel Footplate Covers

- 5"x8" Tapered Gel Footplate Covers..... \$275
Part: ST-FPG5008T-Q.
- 5"x11" Tapered Gel Footplate Covers..... \$275
Part: ST-FPG5011T-Q.
- 6.5"x8" Tapered Gel Footplate Covers for Width Extension..... \$275
Part: ST-FPG6508T-Q.
- 7"x11" Tapered Gel Footplate Cover for Width Extension \$275
Part: ST-FPG7011T-Q.

13.2. Swing-Away Legrests

Select a legrest upper, a legrest extension, and type of footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes.

1a. 60°-90° Swing-Away Legrest Upper

Adjustable in 5° increments.

- P L R 60°-90° Swing-Away Legrest..... \$142.50 ea.
Pair part: RIG120062. Left: FRMASMB17558.
Right: FRMASMB17559

- P L R 60°-90° Swing-Away Legrest, Short \$142.50 ea.
Pair part: RIG136481. Left: RIG132607.
Right: RIG132609.

2. Lower Extension

- P L R 11"-13" Length
Part: See XRef
N/A with Multi-axis.
- P L R 13"-16.5" Length
Part: See XRef
N/A with Multi-Axis.
- P L R 11"-16.5" Length
Part: See XRef
N/A with Angle Adjustable.
- P L R 15.5"-19.5" Length
Part: See XRef

3a. Extruded Footplate(s)

- P L R 7.25"Wx6"D Extruded Footplate..... Standard
Part pair: FRMASMB7847. Left: FRMASMB7861
Right: FRMASMB7862
- P L R 8.5"Wx6"D Extruded Footplate Standard
Part pair: FRMASMB7848. Left: FRMASMB7863
Right: FRMASMB7864
May cause interference on seat widths less than 18."

3b. Swing-Away Legrest, Angle Adjustable Footplate(s) with Heel Loops

¹The 6"x11" footplates may cause interference with seat widths less than 18". ²The 8"x11" footplates may cause interference with seat widths less than 20".

- P L R 4"Wx6"D \$87.50 ea.
Pair part: FRMASMB10564. Left: RIG133537. Right: RIG133560
- P L R 5"Wx6"D..... \$87.50 ea.
Pair part: FRMASMB10565. Left: RIG133538. Right: RIG133561
- P L R 5"Wx8"D \$87.50 ea.
Pair part: FRMASMB10566. Left: RIG133539. Right: RIG133562
- P L R 6"Wx8"D¹ \$87.50 ea.
Pair part: FRMASMB10567. Left: RIG133540. Right: RIG133563
- P L R 6"Wx11"D¹ \$87.50 ea.
Pair part: FRMASMB10568. Left: RIG133541. Right: RIG133564
- P L R 8"Wx11"D² \$87.50 ea.
Pair part: FRMASMB10569. Left: RIG133542. Right: RIG133565

3c. Swing-Away Legrest, Multi-Axis Footplate(s) with Heel Loops

¹The 6"x11" footplates may cause interference with seat widths less than 18". ²The 8"x11" footplates may cause interference with seat widths less than 20".

- P L R 4"Wx6"D \$260 ea.
Pair part: FRMASMB9456. Left: RIG133531. Right: RIG133552
- P L R 5"Wx6"D..... \$260 ea.
Pair part: FRMASMB9457. Left: RIG133532. Right: RIG133553
- P L R 5"Wx8"D \$260 ea.
Pair part: FRMASMB9458. Left: RIG133533. Right: RIG133554
- P L R 6"Wx8"D¹ \$260 ea.
Pair part: FRMASMB9459. Left: RIG133534. Right: RIG133555
- P L R 6"Wx11"D¹ \$260 ea.
Pair part: FRMASMB9460. Left: RIG133535. Right: RIG133558
- P L R 8"Wx11"D² \$260 ea.
Pair part: FRMASMB9461. Left: RIG133536. Right: RIG133559

4. Swing-Away Legrest Accessories

- Heel Loops (pair)..... \$35
Part: FRMASMB7873.
Standard on angle adjustable and multi-axis footplates.
- Wheel Bumper..... \$15 ea.
Part: See XRef
- Gel Neoprene Footrest Sleeve (pair) \$130
Part: INDPART2243.

60°-90° Swing-Away Legrest Compatibility Matrix	Power Recline or Static Seat			Power Tilt			Power Lift		
	Low	Med.	High	Low	Med.	High	Low	Med.	High
Extruded, 11"-16.5"	!	✓	✓	✓	✓	✓	✓	✓	✓
Extruded, 13"-16.5"	!	✓	✓	✓	✓	✓	✓	✓	✓
Extruded, 15.5"-19.5"	!	!	!	!	!	!	!	!	✓
Angle Adjustable, 13"-16.5"	!	✓	✓	✓	✓	✓	✓	✓	✓
Angle Adjustable, 15.5"-19.5"	!	!	!	!	!	!	!	!	✓
Multi-Axis, 11"-16.5"	!	✓	✓	✓	✓	✓	✓	✓	✓
Multi-Axis, 15.5"-19.5"	!	!	!	!	!	!	!	!	✓
11"-13" extension causes no interference at all seat-to-floor heights.									
<input checked="" type="checkbox"/> This rigging configuration causes no interference issues at this seat-to-floor height.									
<input type="checkbox"/> This rigging configuration may have less than 1.75" clearance when fully extended and may interfere with the ground or reduce obstacle climbing ability.									

Legrests section continued on following page.

13. LEGRESTS CONTINUED

13.3. Elevating Legrests

Select a legrest upper, a legrest extension, and type of footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes. Select legrest length in following section.

1. Elevating Legrests Upper

P L R Elevating Legrest Assembly \$275 ea.
 Pair part: FRMASMB6483. Left: RIG138012
 Right: RIG138013. HCPCS: E0990
 Select lower extension and footplates

2. Lower Extension

P L R **11"-14" Length** P L R **14"-17" Length**
 Pair: FRMASMB4755 Pair: FRMASMB4756
 Left: RIG138026 Left: RIG138033
 Right: RIG138030 Right: RIG138035
 For reinforced only For reinforced only

P L R **12"-14.5" Length** P L R **14.5"-18" Length**
 Part: See XRef Part: See XRef
 For angle adjustable only For angle adjustable only

Elevating Legrest Compatibility Matrix	Power Recline or Static Seat			Power Tilt			Power Lift		
	Low	Med.	High	Low	Med.	High	Low	Med.	High
Reinforced, 11"-14"	✓	✓	✓	✓	✓	✓	✓	✓	✓
Reinforced, 14"-17"	!	!	✓	!	✓	✓	✓	✓	✓
Angle Adjustable, 12"-14.5"	✓	✓	✓	✓	✓	✓	✓	✓	✓
Angle Adjustable, 14.5"-18"	!	!	!	!	!	✓	!	✓	✓

This rigging configuration causes no interference issues at this seat-to-floor height.
 This rigging configuration may have less than 1.75" clearance when fully extended and may interfere with the ground or reduce obstacle climbing ability.

3a. Reinforced Cast Footplate(s)

P L R **7.25"Wx6"D Reinforced Footplate**..... Standard
 Pair part: FRMASMB4748. Left: RIG138055. Right: RIG138056

P L R **8.5"Wx6"D Reinforced Footplate** Standard
 Pair part: FRMASMB4747. Left: RIG138060. Right: RIG138062
 May cause interference on seat widths less than 18."

3b. Angle Adjustable Footplate(s)

¹The 6"x11" footplates may cause interference with seat widths less than 18". ²The 8"x11" footplates may cause interference with seat widths less than 20". 4"x6" and 5"x6" footplates are only available in pairs.

P L R **4"Wx6"D** \$87.50 ea.
 Pair parts: FRMASMB8812 x2.

P L R **6"Wx8"D¹ Footplates with Heel Loops** \$87.50 ea.
 Pair part: FRMASMB4738. Left: RIG137793. Right: RIG137795

P L R **5"Wx6"D**..... \$87.50 ea.
 Pair parts: FRMASMB8811 x2.

P L R **6"Wx11"D¹ Footplates with Heel Loops** \$87.50 ea.
 Pair part: FRMASMB4741. Left: RIG137796. Right: RIG137797

P L R **5"Wx8"D Footplates with Heel Loops** \$87.50 ea.
 Pair part: FRMASMB4735. Left: RIG137791. Right: RIG137792

P L R **8"Wx11"D² Footplates with Heel Loops** \$87.50 ea.
 Pair part: FRMASMB4744. Left: RIG137798. Right: RIG137799

4. Elevating Legrest Accessories

Heel Loops (pair)..... \$35
 Part: FRMASMB7873

Falcon Gel Padded Knee Adductor Buttons (KBE) (pair)\$125
 Part: INDPART2242

Wheel Bumper..... \$15 ea.
 Part: See XRef
 Available on angle adjustable footplates.

[^]Please verify the appropriate billing code to be used for this item with all third-party payors.

14. SET LEGREST EXTENSION LENGTH

REQUIRED with matching left and right legrest configuration. All efforts will be made to accommodate legrest length selection, but minor changes may need to be made to meet shipping guidelines. Some adjustments may be needed upon receipt of the unit. Reference the section's compatibility matrix for legrest interference. Non-matching or individual legrest configurations lengths are not set by the manufacturer.

7" L (#3470)	8" L (#3472)	9" L (#3474)	10" L (#3476)	11" L (#3478)	12" L (#3480)	13" L (#3482)	14" L (#3484)	15" L (#3486)
7.5" L (#3471)	8.5" L (#3473)	9.5" L (#3475)	10.5" L (#3477)	11.5" L (#3479)	12.5" L (#3481)	13.5" L (#3483)	14.5" L (#3485)	

15. SEAT CUSHION


Optional. Specific cushion/component sizes available, part numbers, and special order cushions can be found on the respective complete [Stealth Cushion Order Forms](#). This section allows a quick selection of a single cushion. The XXYY portion of the part number directly relates to the dimensions of the cushion, where XX represents the width and YY represents the depth.

15.1. Seat Cushion Model Selection

Stealth Simplicity G Cushion

General use cushion with Stealth Tek 2, 2-way stretch, breathable cover.
Part series: *ST-SIMGXXYY*


Simplicity G\$310-\$370



Stealth TRU-Comfort 2 SPP Cushion

Skin protection and positioning with Coolcore, 4-way stretch, breathable cover.
Part series: *ST-TC2SPPXXYY*


TRU-Comfort 2 SPP\$515-\$575



Stealth Spectrum Foam SPP Cushion

Skin protection and positioning with Coolcore, 4-way stretch, breathable cover. Includes urethane liner.
Part series: *ST-SPCFSPXXYY*

Spectrum Foam SPP\$580-\$630



Stealth Spectrum Gel SPP Cushion

Skin protection and positioning with Coolcore, 4-way stretch, breathable cover. Includes urethane liner.
Part series: *ST-SPCGSPPXXYY*


Spectrum Gel SPP\$580-\$650



Stealth Solution SPP Cushion

Skin protection and positioning with Coolcore, 4-way stretch, breathable cover. Includes urethane liner.
Part series: *ST-SOLSPPXXYY*

Solution SPP\$580-\$645



Stealth Essence SPP Cushion

Skin protection and positioning with Coolcore, 4-way stretch, breathable cover.
Part series: *ST-ESSSPPXXYY*

Essence SPP\$545-\$705


Essence Plus SPP\$800



Stealth Glacial SP Cushion

Skin protection with Coolcore, 4-way stretch, breathable cover.
Part series: *ST-GLSPXXYY*


Glacial SP\$265-\$530



Stealth Zen SP Cushion

Skin protection and positioning with Coolcore, 4-way stretch, breathable cover.
Part series: *ST-ZENSPPXXYY*

Zen SP\$455-\$765



Stealth Premiere P Cushion

Positioning with Coolcore, 4-way stretch, breathable cover.
Part series: *ST-PREMPXXYY*

Premiere P\$585



15.3. Seat Cushion Substitute/Add Fluid Proof Liners

Fluid Proof Urethane Liner\$45-\$50
Part: *URL1709070-XXYY*
Not available with Simplicity, Essence, or TRU-Comfort 2.

Fluid Proof Polyurethane Liner\$70
Essence/TC2 part: *PUL1709080-XXYY*. Others part: *PUL1709030-XXYY*
Not available with Simplicity.

15.4. Seat Cushion Substitute Cover

The standard cover may be substituted by checking an option below. The "*"prefix" of the part number will be the part prefix of the cushion model selected above. For example, *SPCGSPP-NBLK1816* will be an additional coolcore cover for a Spectrum Gel cushion.

Stealth Tek 2 Cover\$90
Part: "*"prefix-TBLKXXYY
Not available with Simplicity, TRU-Comfort 2, Essence or Premiere.

Mesh Cover\$90
Part: "*"prefix-MBLKXXYY
Not available with Simplicity or TRU-Comfort 2.

16. ACCESSORIES

Optional. Multiple accessories can be mounted along the various tracks but some interference may occur. Accessories may need to be installed by the provider.

Flex Seat Cup Holder\$20
Part: *ACCASMB2634*
Not available on TRU-Balance Flex seating without two-post armrests.

Saddle Bag\$15
Part: *DWR1020D010*

Oxygen Holder for Flex Seat\$110
Part: *ACC125006*. HCPCS: *E2208*

Quantum Backpack\$65
Part: *ACC1905853*

Walker Holder\$150
Flex part: *FRMASMB7031*

Accessory Bar\$90
12"W-16"W part: *MEC145955*. 16"W-20"W part: *MEC143975*.
20"W-22"W part: *MEC146018*

Cane/Crutch Holder\$195
Flex part: *ACC152515*. HCPCS: *E2207*

Large, Black, Weather Cover\$240
Part: *ACCCOVR1021*

Flag Holder\$75
Part: *ACCASMB2415*.

Accessories section continued on following page.