Aluminum Transport Chair



HCPCS:


FIGURE A


FIGURE B


FIGURE C


FIGURE D


FIGURE E


FIGURE F

- Weighs only 19 lb
- Aluminum frame is lightweight and durable
- Composite 8" wheels are lightweight and maintenance free
- Both 17" and 19" come with swing-away footrests
- Seat belt included for added safety
- Padded armrests provide additional comfort
- Nylon upholstery is comfortable, lightweight and easy-to-clean
- Aluminum casted folding hinge allows back to fold down for easy storage and transport

PRODUCTS	Description
ATC17-BK	Transport Chair17" Black,1/cs
ATC17-BL	Transport Chair17" Blue 1/cs
ATC17-RD	Transport Chair,Alum.17"Red, 1/ea
ATC19-BK	Transport Chair,Blk,19" Black,1/cs
ATC19-BL	Transport Chair,Alum.19in Blue,1/cs
ATC19-RD	Transport Chair,Alum.19"Red,1/cv

SPECIFICATIONS	Description
Warranty	Limited Lifetime
Seat Dimensions	17" (W) x 16" (D)|19" (W) x 16" (D)
Weight Capacity	300 lbs
Seat To Floor Height	19"
Carton Shipping Weight	24 lbs
Weight	19 lbs
Width Between Arms	15"|17"

ACCESSORIES	Description
STDS1040S	Cup Holder, Swivel,1BX,6/CS, BX