

Rebel Wheelchair

HCPCS:


FIGURE A


FIGURE B


FIGURE C


FIGURE D



FIGURE E



FIGURE F

- Quick-release rear wheels allow for easy storage or transport
- Fold-down back and folding frame allows transport in almost any vehicle
- Easily removable swing-away footrests
- Tool-free removable front casters
- Desk-style armrests make for easy transfers and allows for comfortable seating at a table or desk
- Lightweight and extremely maneuverable
- Large accommodating seating surface
- Carton size makes the Rebel perfect for retail

PRODUCTS	Description
RTLREB18DDA-SF	Rebel W/C,18",FlameRd,DetDsk,SwgFt,1/ea

SPECIFICATIONS	Description
Warranty	Limited Lifetime
Product Weight	37 lb
Material	Steel frame
Seat Dimensions	18" (W) x 16" (D)
Weight Capacity	300 lb
Seat To Floor Height	19.5"
Carton Shipping Weight	42 lb
Overall Width	25.5" (open)|10.5" (closed)
Retail Box Dimensions	36.4" (L) x 11.3" (W) x 28" (H)
Back Height	17"

ACCESSORIES	Description
STDS1040S	Cup Holder, Swivel,1BX,6/CS, BX
STDS850	Seat Belt, Clasp Type,1/ea, EA
STDS855	Seat Belt,60"" Auto-type,1/ea, EA
STDS851	Seat Belt, Velcro Close,1/ea, EA
STDS856	Seat Belt,60"Velcro-type,1/ea, EA
STDS806	Anti-Fold/Theft Device for w/c,1/ea, BX
STDS801GT	WheelLockExtension,ViperGT,SSP2,1/ea, EA
STDS831	Heel Loops, 1pr, BX
STDS1034	Wheelchair Cane/Crutch Holder,1/ea, EA
STDS803	Oxygn Tank Holdr for W/C,1/cs, CS