Deluxe Bariatric Drop-Arm Commode

HCPCS: E0168




FIGURE A


FIGURE B


FIGURE C


FIGURE D


FIGURE E


FIGURE F


FIGURE G



FIGURE H

- Durable, heavy-duty, gray powder-coated steel tubing is sturdy and easy-to-maintain
- Width Between Legs: 26"
- Large, durable snap-on seat
- Easy-to-release arm mechanism allows for safe lateral patient transfers to and from commode
- Padded armrests provide comfort and stability
- Comes complete with 12 qt commode bucket and cover that removes from front
- Pail easily slides in from the front of the commode
- Weight Capacity 1000 lbs
- Made using recycled materials

PRODUCTS	Description
11135-1	Commode,Oversized DropArm,1/cs

SPECIFICATIONS	Description
Warranty	Limited Lifetime
Seat Dimensions	23.25" (W) x 18.5" (D)
Weight Capacity	1000 lbs
Carton Shipping Weight	26 lb
Height Adjustment	17.5" - 22"
Weight	21 lb