Quad Cane Tip



HCPCS:



FIGURE A



FIGURE B



FIGURE C

- Allows cane to stand upright on its own
- Provides better balance and stability with a 4-pronged base
- For use with ¾" diameter cane
- For use with Drive cane model #'s: 10302, 10303, 10304, 10306, 10307, 10309, 10317, 10331, 10335, 10339, 10342, 10345, 10347, 10350, 10351, 10370, 10372, and 10381 series

PRODUCTS	Description
RTL10353	Quad Tip, 1 blister, 10/CV