Adjustable Furniture Risers



HCPCS:




FIGURE A

- The adjustable furniture risers elevate the users bed to a height that allows them to get in and out of bed with greater ease
- The risers have a 4" (10cm) wide cut-out for your furniture legs
- 6298 (Medium) - 2 wood inserts adjustable height 2"-3" (5.08cm to 7.62cm)
- 6299 (Large) - 3 wood inserts adjustable height 3"-6" (7.6cm to 15.2cm)
- Made of a multi-laminated plywood
- Various heights are available by reversing the direction of the block or by adding and deleting inserts
- One block required per leg
- Block Height: 7" (17 8 cm); Inside Diameter: 4" (10 cm)

PRODUCTS	Description
6298	Adjustable Furniture Risers 2"-3"
6299	Adjustable Furniture Risers

SPECIFICATIONS	Description
Weight Capacity	500 lb (227 kg)