# Prime Care Bed Model P703



HCPCS:



FIGURE A




FIGURE B



FIGURE C



FIGURE D




FIGURE E

- 550 lbs. weight capacity, Resident Weight 500 lbs
- 7" to 30" height range
- Optional extension width kit adjusts to 39" - 42"
- Integrated 76", 80" and 84" sleep surface length
- Two way Red/Green color coded floor lock system
- Eight (8) Function Handheld Pendant Control
- Auto Contour
- 4 DC Motors Synchronized Hi/Lo control system keeps bed level
- Four sealed bearings secured in "G" track
- IV pole holders
- Extended Heel Lift
- Eight 3" Heavy Duty casters roll at any height
- Fixed adjustable wall bumper
- Bed does not include wood bed ends

| PRODUCTS | Description |
|---|---|
| P703 | Prime Care Model P703, 1ea |

| SPECIFICATIONS | Description |
|---|---|
| Warranty | Lifetime on Welds |
| Wheels | Accommodate up to 550 lbs Maximum Resident Weight 500 lbs |
| Bed Height Range | 7" - 30" |
| Warranty Limited | 15 Year on Frame |
| Warranty on Batteries | 5 years on Electrical |
| Weight without battery | 205 lbs. |

| ACCESSORIES | Description |
|---|---|
| PLTCDR | LTC Deluxe Rail,1PR, PR |
| PLTCAR | LTC Assist Rail,1PR, PR |
| PLTCAB | LTC Assist Bar, 1PR, PR |
| P703LE | P703 Length Extension 88",1ea, EA |
| P703-BBK | P703 Battery Backup,1ea, EA |
| P703SC04 | Staff Control Panel f/P703 ,1ea, EA |
| P703WE | P703 Width Extension 39-42",1pr, PR |
| 35BWBE-B-SC04 | 35" BRNWOOD BEW/SC04,OKOAK,EA, EA |