
Standard Telescoping Full-Length Side Rail, Brown-Vein Finish

HCPCS: E0310





FIGURE A


FIGURE B


FIGURE C


FIGURE D

- Attractive brown-vein finish (Figure A - B)
- Constructed of 1" steel
- Telescopes with movement of bed
- Easily adjusts up or down with spring-loaded release
- Spring-loaded crossbars are easy-to-install and remove (Figure C)

PRODUCTS	Description
15001ABV	Bed Rails,Full,BrownVein,4 Bar,1/cs

SPECIFICATIONS	Description
Warranty	Limited Lifetime

ACCESSORIES	Description
15050BV	Cross Bars, Universal,Brown Vein, 1/pr, CS