Prime Care Bed Model P903



HCPCS:


FIGURE A


FIGURE B


FIGURE C


FIGURE D


FIGURE E


FIGURE F



FIGURE G


FIGURE H

• 600 lbs. weight capacity*
• 7" super low height can help minimize falls (Figure A) while 30" high working height can facilitate staff &
• Integrated width with push button expansion from 36"-39"-42" provides extra measure of comfort and dignity (Figure I) 80", 84" and 88" sleep surface length helps comfortably accommodate the tallest residents
• 10 function Backlit Hand Pendant includes Auto Contour, Cardiac Comfort Chair (Figure D) and under bed night light provides an added measure of security at night
• Synchronized Hi/Lo control system keeps bed level for resident comfort and positioning
• Welded Steel construction, slat support deck is durable and easy to clean
• Roll at any height convenience with 3" heavy-duty dual footprint casters plus steer control facilitate easy maneuverability (Figure H)
• Color coded, dual pedal floor lock system keeps bed securely in place when engaged (Figure G)
• Adjustable wall bumper helps protect walls and bed from damage (Figure E)
• Standard features include: adj. mattress retainers, heel lift, IV pole holders
• Optional: bed ends (two styles and three wood grain laminate colors), Assist bars, rails, deluxe assist rails, trapeze, embedded staff control, battery backup, transporter
• Trendelenburg/Reverse Trendelenburg capable with staff control option (Figure C)
• Warranty: 15 yr. on frame, lifetime on welds, 5 yr. on electronics

PRODUCTS	Description
P903	Prime Care Bed Model 903

SPECIFICATIONS	Description
Warranty	Lifetime on Welds
Wheels	550 lbs
Optional Width Expansion	36"- 39"-42"
Sleep Surface	80", 84", 88"
Warranty Limited	15 Year on Frame
Warranty on Batteries	5 years on Electrical

ACCESSORIES	Description
PLTCAR	LTC Assist Rail,1PR, PR
PLTCAB	LTC Assist Bar, 1PR, PR
PLTCDR	LTC Deluxe Rail,1PR, PR
P903SC03	P903 Staff control for Linak Motors, EA
42BWBE-A	42" BRNWOOD BEDEND,CHERRY,EA, EA
42MCBE-FM	42" MILL CRK BED END,FIG MAHOG,EA, EA

ACCESSORIES	Description
35MCBE-A	35" MILL CRK BEDEND,CHERRY,EA, EA
42BWBE-B	42" BRNWOOD BEDEND,OKOAK,EA, EA
35BWBE-A	35" BRNWOOD BEDEND,CHERRY,EA, EA
35MCBE-B	35" MILL CRK BEDEND,OKOAK,EA, EA
35BWBE-B	35" BRNWOOD BEDEND,OKOAK,EA, EA