

AquaSense Adjustable Bath Seat with Backrest

HCPCS:


FIGURE A



FIGURE B


FIGURE C


FIGURE D



FIGURE E

• Provide comfort and safety while showering or bathing. The AquaSense® Adjustable Bath Seats with Backrest are made from lightweight aluminum, with blow-molded plastic for the seat and the backrest. They feature a textured surface, large drainage holes, handles in both seat and back, and slip resistant rubber safety tips. Adjustable legs raise or lower the seat height.
• Both lightweight and durable the rust-resistant aluminum frame supports up to 300 lbs.
• Each leg includes a non-skid, non-slip tip offering security on wet surfaces.
• Simply and quickly adjust the height with the help of a push pin to achieve the ideal seated position. Seat height adjustments on this model range from 14" to 21".
• Drainage holes help to prevent water from pooling, making the bench more comfortable, less slippery and ensuring transfer from bathing to exiting the tub is more secure.
• Textured seat surface helps to prevent slipping.
• Push button, snap together design requires no tools and assembles in just minutes.
• Feel good knowing that plastics parts of this product are made using recycled materials.

PRODUCTS	Description
770-510	AQS B.SEAT W/ BK WHITE, 1EA, 2CV, Ret.

SPECIFICATIONS	Description
Weight Capacity	300 lb