Folding Shower Chair

HCPCS: E0245




FIGURE A


FIGURE B


FIGURE C



FIGURE D

- Product available with back (main image) and without (Figure A)
- Conveniently folds flat (Figure B)
- Blow-molded bench provides comfort and strength
- Drainage holes in seat reduce slipping
- Aluminum frame is lightweight, durable and corrosion proof
- Angled legs provide additional stability
- Retail packaged
- Made using recycled materials

PRODUCTS	Description
12486	Folding Shower Bench

SPECIFICATIONS	Description
Warranty	Limited Lifetime
Seat Dimensions	19.75" (W) x 12" (D)|18" (W) x 11" (D)
Weight Capacity	300 lbs
Carton Shipping Weight	6.6 lbs.|7 lbs.
Weight	4.5 lbs.|5 lbs.

ACCESSORIES	Description
12486T	Replacement Tips for 12486 4/Bx, EA