GRAVIS Floor Lift, Power Base 600 lb



HCPCS: E0635


FIGURE A


FIGURE B


FIGURE C


FIGURE D


FIGURE E



FIGURE F


FIGURE G


FIGURE H

- Weight capacity of 600 lb
- Base option: Power
- Elegant, contemporary design with soft lines but still heavy duty
- Oversized sealed rear casters roll easily and lock securely
- Easily lifts resident from floor
- Redundant controls allow for use even if hand control is inoperable
- Smart & Safety Technology LCD display
- Visual battery fuel indicator on Hand Pendant and Lift Control Module
- Audible low battery alert fuel indicator, if maintenance is due, and if lift detects weight overload
- Lift safety stop if resident or caregiver accidentally gets in the way of the boom arm
- Maintenance monitoring: by months or lift cycles
- NEW Six-point cradle with spring clips offer a more secure and wider foot print reducing sling squeeze, shallow height offers greater height range
- Can be charged on, or off the lift
- Includes 1 battery and external charger stand

PRODUCTS	Description
FLP600	Gravis Floor lift, power base, 600 lbs

SPECIFICATIONS	Description
Battery Pack Weight	7.9 lb (5.0A), 4.9 lb (2.9A)
Boom Height	26" - 76"
Ground Clearance	2.9''
Warranty on Battery	1 yr

ACCESSORIES	Description
FLSCALE	Patient Lift Scale 660lbs capacity, EA